



Great Options Great Outcomes Great Value

The CR45 is a Complex Rotation chair with outstanding features at an excellent value. All the options you need and the quality you expect from Ki Mobility!

- Easier to tilt and manage weight with our patented Dual Path 45-degree tilt design.
- Low maintenance design incorporates inset sealed bearings.
- Easy and efficient setup and adjustments throughout.
 Industry leading durability (including lifetime footrest warranty.)

 Easy and efficient setup and adjustments throughout.
 Industry leading durability (including lifetime footrest warranty.)

Specifications

•		
Transport Weight:	38 lb.*	17 kg*
Seat Width:	14-22 in.	36-56 cm
Seat Depth:	14-22 in.	36-56 cm
Seat Height:	14-20.5 in.	36-52 cm
Weight Capacity:	300 lb.	136 kg
WC-19 Approved:	Yes	Yes

* 16"x16" with Stroller Handle Back Posts, 6"x1" Poly Casters. Armrests, Footrests and Rear Wheels Removed. Make it even easier to do business by ordering complete with Axiom and STAR on-chair seating.

