



LIBERTY IS INDEPENDENCE

Finally, one wheelchair that offers true portability, enables independent propulsion and provides the benefits of tilt-in-space!



Liberty is designed for independent mobility, especially foot propulsion. The front seat height remains the same regardless of the tilt angle. Simply choose the tilt angle that provides the best pelvic stability without compromising mobility.

REPOSITIONING

Everyone needs relief from sitting in a static position. Whether your primary concern is pressure relief, facilitating proper posture or just providing a break from gravity, Liberty makes it painless. Featuring an easy to adjust gas spring design, adjusting the system for effortless repositioning is simple.

PORTABILITY

Liberty Folding is the only adult tilt-in-space that can fold into a portable package ideal for transportation. Want an even smaller package? Choose the optional folding backrest. With a transport weight under 26 pounds, it's the lightest adult tilt-in-space available. Liberty gives you the independence to get out and go!



OPTION TOUR



Liberty is one of the most versatile solutions on the market providing the options needed to meet the end users changing needs. In most cases, the modular design allows for easy change of dimensional characteristics making it the ideal choice for short discharge notification in stroke rehabilitation.

LIBERTY BASES

Liberty Folding
Liberty Rigid

ACCESSORIES

O₂ Holder

IV Pole

Transit Wheels

ARMRESTS

Height Adjustable

Angle Adjustable Locking Extendable Flip-Up

SEAT & BACK

Growing Seat Pan

Standard Upholstery

Tension Adjustable

Axiom Back

BACK CANES

Stroller Handle

Height Adjustable

Height Adjustable with 8° Bend

Angle Adjustable (Folding)

Fixed (Non-Folding)

CASTERS

Widths: 1" | 1.5" | 2"

Heights: 4" | 5" | 6" | 7" | 8"

FOOTPLATES

Composite

Composite Angle Adjustable

Aluminum Angle Adjustable

Aluminum Locking Multi-

Angle Adjustable

One Piece Flip-Up Angle Adjustable

HANGERS

Front Mount

Extension Tube

Pro ELR

Pediatric Pro ELR

REAR WHEELS

12" | 16" | 20" | 22" | 24"

WHEEL LOCKS

Push to Lock

Pull to Lock

Hemi Wheel Lock

SPECIFICATIONS

Tilt Range:	0° - 20°	0° - 20°
Seat Width:	14 – 22"	36 - 56 cm
Seat Depth:	14 – 20"	36 - 51 cm
Seat Heights:	13 - 20"	33 - 51 cm
Liberty Folding Transport Weight *:	25.8 lbs.	11.7 kg
Liberty Rigid Transport Weight *:	24.5 lbs.	11.1 kg
Weight Capacity:	250 lbs.	113 kg
WC-19 Approved **:	Yes	Yes
US Pat#:	10,828,213	

^{* 16}x16 fixed back with armrests, footrests and rear wheels removed.
** Liberty is WC-19 approved with the folding backrest only.



Stevens Point, WI | 800-981-1540 | www.kimobility.com DCN0314.5