



Phone (800) 981-1540 | Fax (715) 254-0996  
sales@kimobility.com | www.kimobility.com



Company:

Ship To:

Account Number:

Date:

Address:

Quote: PO Number:

City:

State:

ATP:

Zip:

Phone:

Marked For:

Email/Fax:

Associated Quote/Order#:

Effective: 2/5/2024

CR45- Transit Included (E1161) \$3377

CR45 - Transit Not Included (E1161) \$3049

300 lb weight capacity Tilt range: 0° anterior to 45° posterior

## Frame

### Frame Width

Select box corresponding to desired Frame Width

Frame Width								
14"	15"	16"	17"	18"	19"	20" (E2201)	21" (E2201)	22" (E2201)
N/C	N/C	N/C	N/C	N/C	N/C	\$533	\$533	\$533

### Seating Setup

Make one selection from either Option 1 at right or Option 2 below based on desired seating setup

Setup By Ki Mobility (Option 1) N/C

Select Seat Depth corresponding to desired Frame Depth Range row and Backrest Offset column. Note: All Seat Depths in a column and within the same frame depth range are achievable with future adjustment.

Setup by Frame Depth (Option 2) N/C

Select Frame Depth within desired Frame Depth Range

### Option 2: Setup By Frame Depth

Frame Depth Range	Frame Depth				
Short	14" N/C	15" N/C	16" N/C	17" N/C	18" N/C
Medium	16" N/C	17" N/C	18" N/C	19" N/C	20" (E2203) \$646
Long	18" N/C	19" N/C	20" (E2203) \$646	21" (E2203) \$646	22" (E2204) \$1056

### Seat Pan

Growing Adjustable Aluminum Seat Pan (E2231) \$230  
Solid Seat Pan (E2231) \$165  
Omit Seat Pan N/C

### Transit


Transit System \$343  
Omit Transit N/C

### Option 1: Setup By Ki Mobility

Frame Depth Range	Backrest Offset A						
	-2"	-1"	Flush	+1"	+2"	+3"	+4"
Seat Depth							
Short	16"	15"	14"				
	17"	16"	15"	14"			
	18"	17"	16"	15"	14"		
	19"	18"	17"	16"	15"	14"	
	20" <sup>1</sup>	19"	18"	17"	16"	15"	14"
Medium	18"	17"	16"	15"	14"		
	19"	18"	17"	16"	15"	14"	
	20" <sup>1</sup>	19"	18"	17"	16"	15"	14"
	21" <sup>1</sup>	20" <sup>1</sup>	19"	18"	17"	16"	15"
	22" <sup>2</sup>	21" <sup>1</sup>	20" <sup>1</sup>	19"	18"	17"	16"
Long	20" <sup>1</sup>	19"	18"	17"	16"	15"	14"
	21" <sup>1</sup>	20" <sup>1</sup>	19"	18"	17"	16"	15"
	22" <sup>2</sup>	21" <sup>1</sup>	20" <sup>1</sup>	19"	18"	17"	16"
		22" <sup>2</sup>	21" <sup>1</sup>	20" <sup>1</sup>	19"	18"	17"
			22" <sup>2</sup>	21" <sup>1</sup>	20" <sup>1</sup>	19"	18"

1 Price is \$646 E2203

2 Price is \$1056 E2204





Phone (800) 981-1540 | Fax (715) 254-0996  
sales@kimobility.com | www.kimobility.com



## Tilt Mechanism

Hand Tilt - Dual	N/C
Foot Tilt - Right	N/C
Foot Tilt - Left	N/C

## Backrest

### Backrest

Stroller Handle Back Post (Short 20")	N/C
Stroller Handle Back Post (Tall 24")	N/C
Fixed Height Back with Adjustable Handle (Short 22")	N/C
Fixed Height Back with Adjustable Handle (Tall 26")	N/C
Adjustable Height Back with Adjustable Handle (Short 18" - 25")	N/C
<b>Not available with O2 Holder</b>	
Adjustable Height Back with Adjustable Handle (Tall 20" - 27")	N/C
Reclining Back with Adjustable Handle (Short 22") (Non-Folding)	\$1025
<b>Not available with O2 Holder. If no Gas Spring selection is made, Standard Spring will be used.</b>	
Reclining Back with Adjustable Handle (Tall 26") (Non-Folding)	\$1025
<b>Not available with O2 Holder. If no Gas Spring selection is made, Standard Spring will be used.</b>	

### Backrest Height

18" 19" 20" 21" 22" 23" 24" 25" 26" 27" N/C

### Stroller Handle Back Post Option

**Available only with Stroller Handle Back Post**

Removable Stroller Handle Extension	\$97
-------------------------------------	------

### Other Backrest Options

Non Folding Back	N/C
<b>Standard with Reclining Back</b>	
Adjustable Height Rigidizer Bar	\$82

### Reclining Back Options

**Select option below if Reclining Back is selected on chart above. If no selection is made, Standard Spring will be used.**

Standard Spring	N/C
Heavy Duty Spring	N/C

## Footrest

### Hanger / Front Frame Bend

Hanger	60°	70°	80°	90°		Retail
Extension Tube Mount			1			N/C
Front Mount			1			N/C
Pro Elevating Leg Rests E0990						\$318
Omit Hangers						(\$20)

**1 7" and 8" Casters may interfere with less than 5° tilt**

### PRO ELR Cover

**Available only on PRO ELR**

Gel Pad Covers (pair) (K0108)	\$425
-------------------------------	-------

### 3" Extension

**Available only with Front Mount Hangers**

3" Extension	N/C
3" Extension (Provided as Parts)	N/C



Phone (800) 981-1540 | Fax (715) 254-0996  
sales@kimobility.com | www.kimobility.com



## Seat to Footrest Length

Select Extension and Seat to Footrest Length in row corresponding to desired Footrest. Footrest Length of at least 2.5" shorter than Front Seat Height recommended. Extension Tube Hanger measurements based on 70° Hanger Bend.

Measurements based on 70° Hanger Bend.

Footrest	Extension	Seat to Footrest Length											
Extension Tube Hangers (60°,70°,80°)	Short (10" - 13.5")	10"	10.5"	11"	11.5"	12"	12.5"	13"	13.5"				
	Medium (11.5" - 15.5")	11.5"	12"	12.5"	13"	13.5"	14"	14.5"	15"	15.5"			
	Long (14" - 18")	14"	14.5"	15"	15.5"	16"	16.5"	17"	17.5"	18"			
Front Mount <sup>1</sup> Hangers		Standard Selections									Requires 3" Extension for \$35		
	Short Front Mount Hanger	4" <sup>2</sup>	4.5" <sup>2</sup>	5"	5.5"	6"					8.5"	9"	9.5"
		6.5"	7"	7.5"	8"					10"	10.5"	11"	
	Medium Front Mount Hanger	4" <sup>2</sup>	4.5" <sup>2</sup>	5"	5.5"	6"	6.5"	7"			10.5"	11"	11.5"
		7.5"	8"	8.5"	9"	9.5"	10"			12"	12.5"	13"	
	Long Front Mount Hanger	5"	5.5"	6"	6.5"	7"	7.5"	8"	8.5"		12.5"	13"	13.5"
9"		9.5"	10"	10.5"	11"	11.5"	12"			14"	14.5"	15"	
PRO ELR Low Pivot Setting	PRO ELR Low Pivot Setting Short (10.5" - 15.5")	10.5" - 15.5"											
	PRO ELR Low Pivot Setting Long (13" - 19")	13" - 19"											
PRO ELR High Pivot Setting	PRO ELR High Pivot Setting Short (8.5" - 17")	8.5" - 17"											
	PRO ELR High Pivot Setting Long (11.5" - 17.5")	11.5" - 17.5"											

<sup>1</sup> See Minimum Front Seat Height chart below

<sup>2</sup> Available only with 90° Front Mount Hangers

Minimum Front Seat Height with Front Mount Hangers						
Hanger Length	Short		Medium		Long	
Hanger Bend	70°	80/90°	70°	80/90°	70°	80°
Minimum Front Seat Height with Standard Selection	13"	13"	13"	13.5"	14.5"	15.5"
Minimum Front Seat Height with 3" Extension	13.5"	14"	15.5"	16"	17"	18"

## Hanger Release

If no selection is made 4-Way - In line Hanger Release will be used

Classic	N/C
4-Way - In-Line	N/C
Not available with Pro ELR	
4-Way - Out	N/C
Not available with Pro ELR	

## Footplates

Composite	N/C
Not available with Front Mount Hanger	
Composite Angle Adjustable (K0040)	\$195
Not available with 14" Frame Width	
Aluminum Angle Adjustable (K0040)	\$210
One Piece Flip Up Angle Adjustable (K0040)	\$205
Not available with Front Mount Hanger or PRO ELR	

## Heel Loops

Not available with One Piece Flip Up Angle Adjustable Footplate

Heel Loops (E0951)	\$46
Adjustable Heel Loops (E0951)	\$56
Omit Heel Loops	N/C

## Calf Straps

Velcro Adjustable (K0038)	\$41
Padded Velcro Adjustable (K0038)	\$72
Bodypoint Padded AeroMesh (K0038)	\$128

## Residual Limb Support

	Limb Support	Price	Omit	Cover		Pad Size					
Left	Left Limb Support (E1020)	\$431	Omit Left Hanger	Smooth Cover Left	Air Mesh Cover Left	5X8 Pad Left Side	7X10 Pad Left Side	9x10 Pad Left Side	9x14 Pad Left Side	11x10 Pad Left Side	11x14 Pad Left Side
Right	Right Limb Support (E1020)	\$431	Omit Right Hanger	Smooth Cover Right	Air Mesh Cover Right	5X8 Pad Right Side	7X10 Pad Right Side	9x10 Pad Right Side	9x14 Pad Right Side	11x10 Pad Right Side	11x14 Pad Right Side

## Seat Height, Casters and Rear Wheels

Select Seat Height corresponding to desired Rear Wheel and Caster Size

Rear Wheel Size	Caster Size					
	4"	5"	6"	6" x 2"	7"	8"
12" with Pneumatic or Pneumatic with Airless Insert Tire	Seat Height					
		15.5"				
	16"	16"	16"	16"		
	16.5"	16.5"	16.5"	16.5"		
		17"	17"	17"	17"	
12" with Full Poly Tire		17.5"	17.5"	17.5"	17.5"	
	15"	15"				
	15.5"	15.5"				
	16"	16"	16"	16"		
16" with Pneumatic, Pneumatic with Airless Insert or Full Poly Tire	16.5"	16.5"	16.5"	16.5"		
	14"					
	14.5"					
	15"	15"				
	15.5"	15.5" <sup>1</sup>				
		17"		17" <sup>2</sup>	17" <sup>2</sup>	
		17.5"	17.5"	17.5"	17.5"	
		18"	18"	18"	18"	18"
		18.5"	18.5"	18.5"	18.5"	18.5"
16" with Low Poly Tire		19"	19"	19"	19"	19"
	14"					
	14.5"					
		16"				
	16.5"	16.5"	16.5"	16.5"		
		17"	17"	17"	17"	
		17.5"	17.5"	17.5"	17.5"	
		18"	18"	18"	18"	18"

Rear Wheel Size	Caster Size					
	4"	5"	6"	6" x 2"	7"	8"
20"	Seat Height					
		16"	16"	16"	16"	
		16.5"	16.5"	16.5"	15.5"	
			17"	17"	17"	17"
			17.5"	17.5"	17.5"	17.5"
				19" <sup>1</sup>		
			19.5"	19.5"	19.5"	19.5"
			20"	20"	20"	20"
			20.5"	20.5"	20.5"	20.5"
22"			16.5"		16.5"	
			17"	17"	17"	17"
			17.5"	17.5"	17.5"	17.5"
			18"	18"	18"	18"
			18.5"	18.5"	18.5"	18.5"
			20.5"	20.5"	20.5"	20.5"
24"					17" <sup>1</sup>	17" <sup>1</sup>
			17.5"	17.5"	17.5"	17.5"
			18"	18"	18"	18"
			18.5"	18.5"	18.5"	18.5"
			19"	19"	19"	19"

<sup>1</sup> Available only with Full Poly tire

<sup>2</sup> Not available with Full Poly tire



Phone (800) 981-1540 | Fax (715) 254-0996  
sales@kimobility.com | www.kimobility.com



## Casters and Rear Wheels

### Casters

Caster Type		Caster Size					
		Price	4"	5"	6"	7"	8"
1" Poly		N/C					
1.5" Poly		\$95					
2" Poly		\$97					

### Axle

Fixed \_\_\_\_\_ N/C  
Quick Release \_\_\_\_\_ N/C

### Wheels

		Rear Wheel Size					
		Price	12"	16"	20"	22"	24"
Mag		N/C					

### Rear Tire Type

Full Poly (E2218) \_\_\_\_\_ N/C  
Pneumatic (E2211,E2212) \_\_\_\_\_ N/C  
Pneumatic with Airless Insert (E2211, E2213) \_\_\_\_\_ \$97  
Low Poly (E2218) \_\_\_\_\_ N/C  
*Available only with 16" wheel size*

### Handrims

*Handrims not available on 12" or 16" wheel sizes*

### Standard Handrims

Aluminum Anodized \_\_\_\_\_ N/C  
Plastic Coated \_\_\_\_\_ \$123  
Projection (8 Vertical Projections) \_\_\_\_\_ \$349  
*Available only with full profile tires in combination with 22" and 24" Mag wheels*  
Omit Handrim \_\_\_\_\_ N/C

## Wheel Locks

### Wheel Lock

Push to Lock \_\_\_\_\_ N/C  
Pull to Lock \_\_\_\_\_ N/C  
*Not available with 12" or 16" Full Poly tires*  
Attendant Foot Lock \_\_\_\_\_ \$87  
Omit Wheel Locks \_\_\_\_\_ N/C

### Wheel Lock Options

Extension Handles (E0961) \_\_\_\_\_ \$82  
*Available only with Push to Lock or Pull to Lock Wheel Locks*

## Armrests

Armrest Type	Length	Price
Height Adjustable T-Arm Tall (11.5"-15") E0973	Desk	\$323
	Full	
Height Adjustable T-Arm (8.5"-13.5") E0973	Desk	\$323
	Full	
Height Adjustable T-Arm Low (6.5"-10.5") E0973 <sup>1</sup>	Desk	\$323
	Full	

<sup>1</sup> If requested configuration is unavailable, Standard Height T-Arm will be used  
<sup>2</sup> Minimum armrest height with Reclining Backrest is 12"

Armrest Type	Length	Price
Angle Adjustable Locking Extendable Flip Up E0973 <sup>2</sup>	Desk	\$297
	Full	
Omit Armrests		N/C

### Arm Pads

Standard Arm Pad \_\_\_\_\_ N/C  
Waterfall Arm Pad \_\_\_\_\_ \$26  
*Not available with Angle Adjustable Locking Extendable Flip Up Armrests*  
Foam Grip \_\_\_\_\_ N/C  
*Available only with Angle Adjustable Locking Extendable Flip Up Armrest*

## Accessories

### Removable Side Guard

Composite Side Guards (K0108)	\$144
<i>Not available with T-Arm Armrests</i>	

### Anti-tips

Rear Anti-Tips (E0971)	\$168
<i>Required with Reclining Backrest</i>	
Omit	N/C

### Pelvic Positioning Belt

1 1/2" Auto Buckle Pelvic Positioning Belt (E0978)	\$87
1 1/2" Padded Auto Buckle Pelvic Positioning Belt (E0978)	\$179
2" Airline Buckle Pelvic Positioning Belt (E0978)	\$118
Small Bodypoint Evoflex (16-24") (E0978)	\$200
Medium Bodypoint Evoflex (19-29") (E0978)	\$200
Large Bodypoint Evoflex (25-37") (E0978)	\$200
Belt Mounting Band Clamp	\$36

### U.S. Military Service Patches

Army	Navy	Marines	Air Force	
Coast Guard	US Purple Heart			\$92
Combo Army/Purple Heart		Combo Navy/Purple Heart		
Combo Marines/Purple Heart		Combo Air Force/Purple Heart		
Combo Coast Guard/Purple Heart				\$149


### O2 Holder

O2 Holder (E2208)	\$179
-------------------	-------

### Other Accessories

Tool Kit	\$62
Removable Underseat Pouch	\$56
Ki Mobility Backpack	\$72
Black Neoprene Impact Guards	\$97
Tilt Stop (Pair) (K0108)	\$31
IV Pole	\$292
Canopy with Hardware (K0108)	\$297

*Not available on frame widths greater than 18". See image below.*



## Frame Color

### Standard Colors

Black	N/C
Military Green	N/C
Candy Blue	N/C
Charcoal Gray	N/C
Candy Purple	N/C
Iridescent Burgundy	N/C
Candy Red	N/C
Iridescent Navy	N/C
Copper	N/C
Shadow Black	N/C

### Special Order Colors

*Extended lead time*

Blue Skies	\$72
Emerald Green	\$72
Bubblegum Pink	\$72
Retro Red	\$72
Royal Blue	\$72
Granny Smith	\$72
Smooth White	\$72
Metallic Orange	\$72
Teal	\$72
Yellow	\$72