



Phone (800) 981-1540 | Fax (715) 254-0996
sales@kimobility.com | www.kimobility.com

CATALYST 4

Company:

Ship To:

Account Number:

Date:

Address:

Quote: PO Number:

City:

State:

ATP:

Zip:

Phone:

Marked For:

Email/Fax:

Associated Quote/Order#:

Effective: 1/10/2024

Catalyst 4 Swing Away Frame (K0004) \$1327

Catalyst 4 Super Low Frame (K0004) \$1327

15" Lowest Front Seat Height 250 lb weight capacity

Frame

Seat Width

14"	\$103
16"	N/C
18"	N/C
20" (E2201)	\$564
22" (E2201)	\$564

Weight Capacity

If no selection is made, 250 lb weight capacity will be used

250 lbs	N/C
350 lbs	\$564

Requires minimum 6" caster height

Seat Depth

16"	N/C
18"	N/C
20" (E2203)	\$564

Backrest

Back Post / Height Range / Back Height

Select Back Height in row corresponding to desired Back Post

Back Post	Price	Back Height			
Standard (16.5" - 19.5")	N/C	16.5"	17.5"	18.5"	19.5"
Angle Adjustable w/ Push Handle	\$200	16"	17"	18"	

Back Upholstery

If selecting both upholstery and Axiom Back, the upholstery will be installed.

Standard	N/C
Tension Adjustable with Standard Cover (E2611)	\$323
Tension Adjustable with Air Mesh Breathable Cover (E2611)	\$415
Omit	(\$10)

Footrest

Hanger

Frame	Hanger	60°	70°		Retail
Swing Away	Extension Tube Mount				N/C
	Pro Elevating Leg Rests (E0990)				\$359
	Omit Hangers				N/C

Footplates

Composite	N/C
Composite Angle Adjustable (K0040)	\$200
Not available with ELR or 14" Seat Width	
Aluminum Angle Adjustable (K0040)	\$205

PRO ELR Cover

Available only on PRO ELR

Gel Pad Covers (PAIR) (K0108) \$436

HCPCS codes provided are based on code verification by the PDAC or our interpretation of the code definitions.

Proper use of HCPCS codes does not ensure coverage or payment. For coverage information, verify the appropriate payer's coverage policy.

DCN0012.36 Pricing and specifications subject to change without notice.

Footrest / Extension / Seat to Footrest Length

Select Extension and Seat to Footrest Length in row corresponding to desired Footrest. Footrest Length of at least 2.5" shorter than Front Seat Height recommended. Extension Tube Hanger measurements are based on 70° Hanger

Measurements are based on 70° Hanger

Footrest	Extension	Seat to Footrest Length								
Extension Tube Hangers (60°,70°)	Short (11" - 15")	11"	11.5"	12"	12.5"	13"	13.5"	14"	14.5"	15"
	Medium (13" - 17")	13"	13.5"	14"	14.5"	15"	15.5"	16"	16.5"	17"
	Long (15.5" - 19.5")	15.5"	16"	16.5"	17"	17.5"	18"	18.5"	19"	19.5"
PRO ELR Low Pivot Setting	PRO ELR Low Pivot Setting Short (12" - 17.5")	12" - 17.5"								
	PRO ELR Low Pivot Setting Long (15" - 21")	15" - 21"								
PRO ELR High Pivot Setting	PRO ELR High Pivot Setting Short (10.5" - 15.5")	10.5" - 15.5"								
	PRO ELR High Pivot Setting Long (13" - 19")	13" - 19"								

Heel Loops

Heel Loops (E0951)	\$42
Adjustable Heel Loops (E0951)	\$51
Omit Heel Loops	N/C

Calf Straps

Velcro Adjustable (K0038)	\$41
Padded Velcro Adjustable (K0038)	\$72
Bodypoint Padded AeroMesh (K0038)	\$133

Residual Limb Support

	Limb Support	Price	Omit	Cover		Pad Size					
Left	Left Limb Support (E1020)	\$431	Omit Left Hanger	Smooth Cover Left	Air Mesh Cover Left	5X8 Pad Left Side	7X10 Pad Left Side	9x10 Pad Left Side	9x14 Pad Left Side	11x10 Pad Left Side	11x14 Pad Left Side
Right	Right Limb Support (E1020)	\$431	Omit Right Hanger	Smooth Cover Right	Air Mesh Cover Right	5X8 Pad Right Side	7X10 Pad Right Side	9x10 Pad Right Side	9x14 Pad Right Side	11x10 Pad Right Side	11x14 Pad Right Side

Front Seat Height and Casters

Front Seat Height

Select Front Seat Height in row corresponding to desired caster size. Note: Selections in first three Front Seat Height columns require Catalyst 4 Super Low.

Caster Size	Front Seat Height													
	Uses Super Low Front Frame													
4"	13.5"	14"	14.5"	15"	15.5"	16"	16.5"	17"	17.5"	18"	18.5"	19"	19.5"	20"
5"	14.5"	15"	15.5"	16"	16.5"	17"	17.5"	18"	18.5"	19"	19.5"	20"	20.5"	
6"	15.5"	16"	16.5"	17"	17.5"	18"	18.5"	19"	19.5"	20"	20.5"	21"		
7"	17"	17.5"		18"	18.5"	19"	19.5"	20"	20.5"	21"				
8"	17.5"	18"	18.5"	19"	19.5"	20"	20.5"	21"						

Selection for Catalyst 4 Super Low is part of base frame selection on page 1

Caster Type

Caster Type		Caster Size					
		Price	4"	5"	6"	7"	8"
Caster Type	1" Poly	N/C					
	1.5" Poly	\$67					
	2" Poly	\$97					

Maximum Seat Dump Reference

(Front Seat Height minus Rear Seat Height)

Seat Depth	16"	18"	20"
Maximum Dump	2"	2.5"	2.5"

Fork Stem

Standard Fork Stem	N/C
3/4" Longer Fork Stem	\$36
1 1/2" Longer Fork Stem	\$36

Rear Seat Height and Wheels

Wheels Size / Tire Profile / Rear Seat Height

Select Rear Seat Height in row corresponding to desired Wheel Size and Tire Profile

Wheel Size	Tire Profile	Rear Seat Height				
20"	Low	12.5" ¹	13.5"	14.5"	15.5"	16.5"
	Full	13" ¹	14"	15"	16"	17"
22"	Low	13.5" ¹	14.5"	15.5"	16.5"	17.5"
	Full	14" ¹	15"	16"	17"	18"
24"	Low	14.5" ¹	15.5"	16.5"	17.5"	18.5"
	Full	15" ¹	16"	17"	18"	19"

¹ Not available with Flip Back Armrests

Wheels

		Rear Wheel Size			
		Price	20"	22"	24"
Rear Wheel Type	Spoke	\$103			
	Mag	N/C			

Rear Tire Type

Low Profile

Low Poly 1" (E2218) _____ N/C
 Not available with 350 lb Weight Capacity

Full Profile

Full Poly 1 3/8" (E2218) _____ \$41
 Street Pneumatic 1 3/8" (E2211, E2212) _____ \$82
 Pneumatic with Airless Insert 1 3/8" (E2211, E2213) _____ \$210

Handrims

Standard Handrims

Aluminum Anodized _____ N/C
 Plastic Coated _____ \$110
 Projection (8 Vertical Projections) _____ \$359
Available only with full profile tires in combination with 22", 24", and 26" MAXX Spoke wheels, or 22" and 24" Mag wheels
 Omit Handrim _____ N/C

Ergonomic Handrims

Not available with 20" or 26" wheels. Tab mount only available with 24" wheels.

Lite Regular Grip _____ \$477
 Standard Regular Grip _____ \$456
 Lite Extra Grip _____ \$528
 Standard Extra Grip _____ \$507
 Lite Omit Grip _____ \$405
 Standard Omit Grip _____ \$354

Handrim Options

3M Micro Grip _____ \$72
 Not available on Plastic Coated or Projection Handrims

Axles

Axle Bracket

Offset Axle Bracket _____ N/C
 Reverse Mount Offset Axle Bracket (Amputee) _____ N/C
 Not available with Swing Away Armrests

Axle

Fixed Axle _____ N/C
 Quick Release _____ \$97

Wheel Locks

Wheel Lock

Push to Lock _____ N/C
 Push to Lock Flush Mount _____ N/C
 Pull to Lock _____ N/C
 Omit Wheel Locks _____ N/C

Wheel Lock Options

Extension Handles (E0961) _____ \$82
Available only with Push to Lock or Pull to Lock Wheel Locks
 Grade Aids (E0974) _____ \$159
Available only with Push to Lock Wheel Locks and Pneumatic or Pneumatic w/ Airless Insert tires



Phone (800) 981-1540 | Fax (715) 254-0996
sales@kimobility.com | www.kimobility.com

CATALYST 4

Armrests

Armrest Type	Length	Price
Height Adjustable T-Arm Tall (11.5"-15") E0973	Desk Full	\$323
Height Adjustable T-Arm (8.5"-13.5") E0973	Desk Full	\$323
Height Adjustable T-Arm Low (6.5"-10.5") E0973 ¹	Desk Full	\$323
Height Adjustable Flip Back T-Arm Tall (11.5"-15") E0973	Desk Full	\$350
Height Adjustable Flip Back T-Arm (9.5"-13.5") E0973	Desk Full	\$350

¹ If requested configuration is unavailable, Standard Height T-Arm will be used

² Not available with Angle Adjustable Back Posts

³ Not available with Arm Pads

⁴ Not Available with 15" Rear Seat Height and 24" Wheels

Armrest Type	Length	Price
Flip back (10.5") ^{2,4}	Desk Full	\$60
Swing Away (10.5"-12.5") E0973 ^{2,3}		\$251
Omit Armrests		N/C

Arm Pads

Not available with Swing Away Armrests

Waterfall Arm Pad \$26
Available only with Height Adjustable T-Arm Tall, Height Adjustable T-Arm, and Height Adjustable T-Arm Low

Accessories

Removable Side Guard

Composite Side Guards (K0108) \$159
Not available with T-Arm or Flip Back Armrests

Pelvic Positioning Belts

1 1/2" Auto Buckle Pelvic Positioning Belt (E0978) \$87
1 1/2" Padded Auto Buckle Pelvic Positioning Belt (E0978) \$179
2" Airline Buckle Pelvic Positioning Belt (E0978) \$118
Small Bodypoint Evoflex (16-24") (E0978) \$200
Medium Bodypoint Evoflex (19-29") (E0978) \$200
Large Bodypoint Evoflex (25-37") (E0978) \$200
Belt Mounting Band Clamp \$36

O2 Holder

O2 Holder - Left Side (E2208) \$313
O2 Holder - Right Side (E2208) \$313

Anti-tips

Rear Anti-Tips (E0971) \$159
Omit Rear Anti-Tips N/C

Other Accessories

Tool Kit \$46
Removable Underseat Pouch \$51
Ki Mobility Backpack \$80
Black Neoprene Impact Guards \$97
Clear Spoke Guards (K0065) \$123
Not available with Mag wheel
Black Spoke Guards (K0065) \$123
Not available with 20" wheel size or Mag wheel
Cane and Crutch Holder \$210
IV Pole (Straight Back Canes Will Be Used) \$292

Frame Color

Black Gloss N/C
Military Green N/C
Charcoal Gray N/C
Emerald Green N/C
Iridescent Burgundy N/C
Iridescent Navy N/C
Retro Red N/C
Royal Blue N/C
Shadow Black N/C
Smooth White N/C

Catalyst Front Seat Height Matrix																													
4" Caster				5" Caster				6" Caster				7" Caster				8" Caster													
FSH	Super Low	Stem	Fork	FSH	Super Low	Stem	Fork	FSH	Super Low	Stem	Fork	FSH	Super Low	Stem	Fork	FSH	Super Low	Stem	Fork										
15	13.5	Std	4"	16	14.5	Std	4"	17	15.5	Std	5"	18	16.5	Std	6"	19	17.5	Std	6"										
15.5	14			16.5	15.5	3/4		17.5	16.5	3/4		18.5	17			19.5	18.5	3/4											
15.5 ¹	14 ¹	17.5		16	1 1/2	18.5		17	1 1/2	18.5		17.5	3/4			20.5	n/a	1 1/2											
16	14.5	3/4		16	14.5	Std	5"	17	15.5	Std	6"	19	18	1 1/2		19	17.5	Std	7"										
16.5 ¹	15 ¹			16.5	15			17.5	16			19.5	18																
17	15.5	16.5		15.5	3/4			18	16.5			Std		20		18.5	Std			20	18.5								
15.5	14	Std	5"	17	16	1 1/2		17.5	16.5	3/4		6"	18.5	17		Std	20.5	19		3/4									
16	14.5			17.5	16			18	17				19	17.5			19.5	18.5											
16	15	3/4		18	16.5			Std	18.5				17.5	19.5		18	20	18.5			20.5	n/a							
16.5	15.5			17	15.5	19	17.5		19	18	20.5		n/a	21		n/a													
17	15.5	1 1/2		17.5	16	Std	19.5	18	19.5	18.5	3/4		21	n/a		1 1/2													
17.5	16			18	16.5		19	17.5	20	19			20.5	19		21	n/a												
16.5	15	Std	6"	17.5	16	3/4	6"	18	16.5	Std	7"	20	19	1 1/2															
17	15.5			18	16.5			18.5	17			20.5	n/a																
17	16	3/4		18.5	17			19	17.5			20	18.5																
17.5	16.5			18	16.5	19.5		18	20.5	19																			
18	16.5	1 1/2		18.5	17	Std		18.5	17.5	21		n/a																	
18.5	17			19	17.5																								
18	16.5	Std	7"	18	16.5	3/4	7"					19.5	18.5																
18.5	17			18	16.5							20	19																
18.5	17.5	3/4		19	17.5							19.5	18											20	18.5				
19	18			1 1/2	18.5	17.5						20.5	19																
19.5	18.5	19			18	20.5						19	21											n/a					
20	19	1 1/2		19.5	18.5																								
				19.5	18																								
				20	18.5																								
				20.5	19																								

¹ Not available with 4 x 1.5" Soft Roll Aluminum

Seat heights with "FSH" heading available with Standard Front Frame. Seat heights with "Super Low" heading available with Super Low Front Frame. Actual seat height may vary depending upon configuration.