

This form must be opened in Adobe Acrobat. The submit button will not work from most browser windows.



Phone (800) 981-1540 | Fax (715) 254-0996
 sales@kimobility.com | www.kimobility.com

Effective 03/13/2017

Company: _____

Ship To: _____

Account Number: _____ Date: _____

Address: _____

Quote PO Number: _____

City: _____ State: _____

RTS: _____

Zip: _____ Phone: _____

Marked For: _____

Email/Fax: _____

- Catalyst 4** K0004 _____ **\$1140**
 15" Lowest Front Seat Height
 250 lb weight capacity
- Catalyst 4 Super Low** K0004 _____ **\$1190**
 13.5" Lowest Front Seat Height
 250 lb weight capacity

*HCPCS codes provided are based on code verification by the PDAC or our interpretation of the code definitions.
 Proper use of HCPCS codes does not ensure coverage or payment. For coverage information, verify the appropriate payer's coverage policy.*

Frame

Seat Width / Seat Depth

Select box corresponding to desired Seat Width and Depth

		Seat Width					
		14"	16"	18"	20"	22"	
		Price HCPCS	\$100	N/C	N/C	\$420 E2201	\$495 E2201
Seat Depth	16"	N/C	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>
	18"	N/C	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>
	20"	\$450 E2203	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>

Heavy Duty Option

- 350 lb Weight Capacity _____ \$500
 Requires minimum 6" caster height

Backrest

Back Post / Height Range / Back Height

Select Back Height in row corresponding to desired Back Post

Back Post	Price	Back Height			
Standard	N/C	<input type="checkbox"/> 16.5	<input type="checkbox"/> 17.5	<input type="checkbox"/> 18.5	<input type="checkbox"/> 19.5
Angle Adjustable w/ Push Handle	\$185	<input type="checkbox"/> 16	<input type="checkbox"/> 17	<input type="checkbox"/> 18	

Back Upholstery

- Standard _____ Std
- Tension Adjustable E2611 _____ \$245
- Omit Back Upholstery _____ - \$10

Military Patch

- Attached only if ordered with Tension Adjustable Back Upholstery*
- US Military Service Patch _____ \$90
 - Army Navy Marines Air Force Coast Guard
 - US Purple Heart Patch _____ \$90
 - Purple Heart with Military Service Patch Combo _____ \$145
 - Army Navy Marines Air Force Coast Guard

Footrest

Hanger Bend	
<input type="checkbox"/> 60°	N/C
<input type="checkbox"/> 70°	N/C
<input type="checkbox"/> Elevating Leg Rests (ELR) E0990	\$295
<input type="checkbox"/> Pro Elevating Leg Rests (PRO ELR) E0990 New!	\$350
<input type="checkbox"/> Omit Hangers	N/C

Footplate	
<input type="checkbox"/> Composite	N/C <i>Not available with PRO ELR</i>
<input type="checkbox"/> Composite Angle Adjustable K0040	\$185 <i>Not available with ELR or 14" Seat Width</i>
<input type="checkbox"/> Aluminum Angle Adjustable K0040	\$185

Footrest / Extension / Seat to Footrest Length

Select Seat to Footrest Length in row corresponding to desired Footrest and Extension. Footrest Length of at least 2.5" shorter than Front Seat Height recommended. Extension Tube Hanger measurements are based on 70° Hanger / Front Frame Bend.

Footrest	Extension	Seat to Footrest Length									
Extension Tube Hangers (60°, 70°)	Short	<input type="checkbox"/> 11	<input type="checkbox"/> 11.5	<input type="checkbox"/> 12	<input type="checkbox"/> 12.5	<input type="checkbox"/> 13	<input type="checkbox"/> 13.5	<input type="checkbox"/> 14	<input type="checkbox"/> 14.5	<input type="checkbox"/> 15	
	Medium	<input type="checkbox"/> 13	<input type="checkbox"/> 13.5	<input type="checkbox"/> 14	<input type="checkbox"/> 14.5	<input type="checkbox"/> 15	<input type="checkbox"/> 15.5	<input type="checkbox"/> 16	<input type="checkbox"/> 16.5	<input type="checkbox"/> 17	
	Long	<input type="checkbox"/> 15.5	<input type="checkbox"/> 16	<input type="checkbox"/> 16.5	<input type="checkbox"/> 17	<input type="checkbox"/> 17.5	<input type="checkbox"/> 18	<input type="checkbox"/> 18.5	<input type="checkbox"/> 19	<input type="checkbox"/> 19.5	
Elevating Leg Rests (ELR)	Short	<input type="checkbox"/> 12	<input type="checkbox"/> 13	<input type="checkbox"/> 14	<input type="checkbox"/> 15	<input type="checkbox"/> 16					
	Long	<input type="checkbox"/> 13	<input type="checkbox"/> 14	<input type="checkbox"/> 15	<input type="checkbox"/> 16	<input type="checkbox"/> 17	<input type="checkbox"/> 18				
Extended Length ELR	Short	<input type="checkbox"/> 18	<input type="checkbox"/> 19	<input type="checkbox"/> 20	<input type="checkbox"/> 21	<input type="checkbox"/> 22					
	Long	<input type="checkbox"/> 18	<input type="checkbox"/> 19	<input type="checkbox"/> 20	<input type="checkbox"/> 21	<input type="checkbox"/> 22	<input type="checkbox"/> 23	<input type="checkbox"/> 24			
PRO ELR Low Pivot	Short	<input type="checkbox"/> 12 - 17.5									
	Long	<input type="checkbox"/> 15 - 21									
PRO ELR High Pivot	Short	<input type="checkbox"/> 10.5 - 15.5									
	Long	<input type="checkbox"/> 13 - 19									

Heel Loop	
<input type="checkbox"/> Heel Loops E0951	\$30
<input type="checkbox"/> Omit Heel Loops	N/C

Calf Strap	
<input type="checkbox"/> Velcro Adjustable K0038	\$35
<input type="checkbox"/> Padded Velcro Adjustable K0038	\$68
<input type="checkbox"/> Bodypoint Padded AeroMesh K0038	\$102

Residual Limb Support			
<input type="checkbox"/> Residual Limb Support E1020			
<input type="checkbox"/> Left	\$395	<input type="checkbox"/> Right	\$395
<input type="checkbox"/> Omit Left Hanger	N/C	<input type="checkbox"/> Omit Right Hanger	N/C
Cover			
<input type="checkbox"/> Smooth	<input type="checkbox"/> Air Mesh	<input type="checkbox"/> Smooth	<input type="checkbox"/> Air Mesh
Pad Size			
<input type="checkbox"/> 5 x 8	<input type="checkbox"/> 9 x 14	<input type="checkbox"/> 5 x 8	<input type="checkbox"/> 9 x 14
<input type="checkbox"/> 7 x 10	<input type="checkbox"/> 11 x 10	<input type="checkbox"/> 7 x 10	<input type="checkbox"/> 11 x 10
<input type="checkbox"/> 9 x 10	<input type="checkbox"/> 11 x 14	<input type="checkbox"/> 9 x 10	<input type="checkbox"/> 11 x 14

Front Seat Height and Casters

Caster Size / Front Frame Type / Front Seat Height

Select Front Seat Height in row corresponding to desired caster size.

Note: Selections in first three Front Seat Height columns require Catalyst 4 Super Low.

Caster Size	Front Seat Height													
	Uses Super Low Front Frame													
4"	<input type="checkbox"/> 13.5 ¹	<input type="checkbox"/> 14	<input type="checkbox"/> 14.5	<input type="checkbox"/> 15 ²	<input type="checkbox"/> 15.5	<input type="checkbox"/> 16	<input type="checkbox"/> 16.5	<input type="checkbox"/> 17	<input type="checkbox"/> 17.5	<input type="checkbox"/> 18	<input type="checkbox"/> 18.5	<input type="checkbox"/> 19	<input type="checkbox"/> 19.5	<input type="checkbox"/> 20
5"	<input type="checkbox"/> 14.5	<input type="checkbox"/> 15	<input type="checkbox"/> 15.5	<input type="checkbox"/> 16	<input type="checkbox"/> 16.5	<input type="checkbox"/> 17	<input type="checkbox"/> 17.5	<input type="checkbox"/> 18	<input type="checkbox"/> 18.5	<input type="checkbox"/> 19	<input type="checkbox"/> 19.5	<input type="checkbox"/> 20	<input type="checkbox"/> 20.5	
6"	<input type="checkbox"/> 15.5	<input type="checkbox"/> 16	<input type="checkbox"/> 16.5	<input type="checkbox"/> 17	<input type="checkbox"/> 17.5	<input type="checkbox"/> 18	<input type="checkbox"/> 18.5	<input type="checkbox"/> 19	<input type="checkbox"/> 19.5	<input type="checkbox"/> 20	<input type="checkbox"/> 20.5	<input type="checkbox"/> 21		
7"	<input type="checkbox"/> 17	<input type="checkbox"/> 17.5		<input type="checkbox"/> 18	<input type="checkbox"/> 18.5	<input type="checkbox"/> 19	<input type="checkbox"/> 19.5	<input type="checkbox"/> 20	<input type="checkbox"/> 20.5	<input type="checkbox"/> 21				
8"	<input type="checkbox"/> 17.5	<input type="checkbox"/> 18	<input type="checkbox"/> 18.5	<input type="checkbox"/> 19	<input type="checkbox"/> 19.5	<input type="checkbox"/> 20	<input type="checkbox"/> 20.5	<input type="checkbox"/> 21						

Selection for Catalyst 4 Super Low is part of base frame selection on page 1

¹ Front Seat Height not available with 4 x 1" Poly Caster

² Super Low Front Frame required when using 4 x 1" Poly Caster

Caster Type

- 1" Poly _____ N/C
- 1" Pneumatic _____ \$80
Not available with 4", 5" and 7" Caster Sizes
- 1.5" Poly _____ \$50
Not available with 7" and 8" Caster Sizes
- 2" Poly _____ \$70
Not available with 4", 5", 7" and 8" Caster Sizes
- 2" Pneumatic w/ Foam Insert _____ \$225
Not available with 4", 5", 6" and 7" Caster Sizes

Maximum Seat Dump Reference

(Front Seat Height minus Rear Seat Height)

Seat Depth	16"	18"	20"
Maximum Dump	2"	2.5"	2.5"

Fork Stem

- Standard Fork Stem _____ Std
- 3/4" Longer Fork Stem _____ \$35
- 1 1/2" Longer Fork Stem _____ \$35

Rear Seat Height and Wheels

Wheel Size / Tire Profile / Rear Seat Height

Select Rear Seat Height in row corresponding to desired Wheel Size and Tire Profile

Wheel Size	Tire Profile	Rear Seat Height				
20"	Low	<input type="checkbox"/> 12.5 ¹	<input type="checkbox"/> 13.5	<input type="checkbox"/> 14.5	<input type="checkbox"/> 15.5	<input type="checkbox"/> 16.5
	Full	<input type="checkbox"/> 13 ¹	<input type="checkbox"/> 14	<input type="checkbox"/> 15	<input type="checkbox"/> 16	<input type="checkbox"/> 17
22"	Low	<input type="checkbox"/> 13.5 ¹	<input type="checkbox"/> 14.5	<input type="checkbox"/> 15.5	<input type="checkbox"/> 16.5	<input type="checkbox"/> 17.5
	Full	<input type="checkbox"/> 14 ¹	<input type="checkbox"/> 15	<input type="checkbox"/> 16	<input type="checkbox"/> 17	<input type="checkbox"/> 18
24"	Low	<input type="checkbox"/> 14.5 ¹	<input type="checkbox"/> 15.5	<input type="checkbox"/> 16.5	<input type="checkbox"/> 17.5	<input type="checkbox"/> 18.5
	Full	<input type="checkbox"/> 15 ¹	<input type="checkbox"/> 16	<input type="checkbox"/> 17	<input type="checkbox"/> 18	<input type="checkbox"/> 19

¹ Not available with Flip Back Armrests

Rear Wheel Type

- Assembled 24" Six-Spoke Mag with Treaded Tire and Composite Handrim _____ - \$25
No substitutions available. Do not select Rear Tire Type or Handrim. Measures as Full Tire.
- Spoke _____ \$50
- Mag _____ N/C

Rear Tire Type

Low Profile

- Low Poly 1" E2218 _____ Std

Full Profile

- Full Poly 1 3/8" E2218 _____ \$30
- Street Pneumatic 1 3/8" E2211, E2212 _____ \$80
- Pneumatic with Airless Insert 1 3/8" E2211, E2213 _____ \$200

Handrim

Standard Handrims

- Aluminum Anodized _____ Std
- Plastic Coated _____ \$95
- Projection (8 Vertical Projections) _____ \$295
Available on 22" and 24" Spoke wheels and 24" Mag wheels with full tires
- Omit Handrim _____ N/C

Handrim Option

- 3M Micro Grip _____ \$65
Not available on Plastic Coated or Projection Handrims

Natural Fit Handrims

Not available on 20" rear wheel or 22" Mag wheels. Not recommended with Mag wheels.

- Natural Fit LT Regular Grip _____ \$445
- Natural Fit Standard Regular Grip _____ \$445
- Natural Fit LT Super Grip _____ \$495
- Natural Fit Standard Super Grip _____ \$495
- Natural Fit LT Omit Grip _____ \$345
- Natural Fit Standard Omit Grip _____ \$345

Axles

Axle Bracket

- Offset Axle Bracket _____ Std
- Reverse Mount Offset Axle Bracket _____ N/C
Not available with Swing Away Armrests

Axle

- Fixed Axle _____ Std
- Quick Release Axle _____ \$95

Wheel Locks

Wheel Lock

- Push to Lock _____ Std
- Push to Lock Flush Mount _____ N/C
- Pull to Lock _____ N/C
- Omit Wheel Locks _____ N/C

Wheel Lock Options

- Extension Handles E0961 _____ \$67
Available only with Push to Lock or Pull to Lock Wheel Locks
- Grade Aids E0974 _____ \$150
Available only with Push to Lock Wheel Locks and Pneumatic or Pneumatic w/ Airless Insert Rear Tire Types

Armrests

Armrest Please select Arm Pad option

- Height Adjustable T-Arm (8.5" - 13.5") E0973 _____ \$265
- Height Adjustable T-Arm Low (7.5" - 11.5") E0973 _____ \$265
- Height Adjustable Flip Back (9.5" - 13.5") E0973 _____ \$275
Not available with Angle Adjustable Back Posts
- Flip Back (10.5") _____ N/C
Not available with Angle Adjustable Back Posts

Arm Pad

Not available with Swing Away Armrests

- Standard Arm Pad _____ N/C Desk Full
- Waterfall Arm Pad _____ \$25 Desk Full
Available only with Height Adjustable T-Arms

Armrest continued

- Swing Away (10.5" - 12.5") E0973 _____ \$195
Not available with Angle Adjustable Back Posts or Arm Pads
- Omit Armrests _____ N/C

Accessories

Removable Side Guard

- Composite Side Guards K0108 _____ \$140
Not available with T-Arm or Flip Back Armrests

Anti-Tip

- Rear Anti-Tips E0971 _____ \$140

Positioning Belt

- 1 1/2" Auto Buckle E0978 _____ \$70
- 1 1/2" Padded Auto Buckle E0978 _____ \$155
- 2" Airline Buckle E0978 _____ \$110
- Bodypoint Evoflex - Pelvic Stabilizer E0978 _____ \$183
Size: Small (16-24") Medium (19-29") Large (25-37")
- Bodypoint Belt Mounting Kit _____ \$105
Style: Bracket Band Clamp 1"

Other Accessories

- Tool Kit _____ \$45
- Removable Underseat Pouch _____ \$50
- Ki Mobility Backpack _____ \$65
- Black Neoprene Impact Guards _____ \$85
- Clear Spoke Guards K0065 _____ \$115
Not available with Mag wheels
- Black Spoke Guards K0065 _____ \$115
Not available with 20" rear wheel or Mag wheels
- Cane and Crutch Holder _____ \$195
- IV Pole _____ \$265
Straight Back Canes will be used when this option is selected

O₂ Holder

- O₂ Holder E2208 _____ \$295
 Right Side Left Side

Frame Color

Matte Finish

- Army Green _____ N/C
- Charcoal Gray _____ N/C
- Emerald Green _____ N/C
- Iridescent Burgundy _____ N/C
- Iridescent Navy _____ N/C
- Retro Red _____ N/C
- Royal Blue _____ N/C
- Shadow Black _____ N/C
- Smooth White _____ N/C

Gloss Finish

- Black _____ N/C

Catalyst Front Seat Height Matrix																						
4" Caster				5" Caster				6" Caster				7" Caster				8" Caster						
FSH	Super Low	Stem	Fork	FSH	Super Low	Stem	Fork	FSH	Super Low	Stem	Fork	FSH	Super Low	Stem	Fork	FSH	Super Low	Stem	Fork			
15 ¹	13.5 ¹	Std	4"	16	14.5	Std	4"	17	15.5	Std	5"	18	17	Std	6"	19	18	Std	6"			
15.5	14			16.5	15.5	3/4		18	16.5	3/4		18.5	17.5			20	18.5	3/4				
16 ²	14.5 ²	17.5		16	1 1/2	18.5		17	1 1/2	19		17.5	20.5	n/a		1 1/2						
16.5	15	3/4		16	15	Std	17	16	Std	19.5	18	19	17.5	Std		19.5	18	Std	7"			
16.5 ²	15 ²	1 1/2		16.5	15.5		17.5	16.5		18	17	19.5	18.5			20	18.5					
17	15.5	Std		5"	17	15.5	3/4	5"	18	17	Std	6"	18.5	17.5		1 1/2	18.5	17.5	20.5	19	Std	
15.5	14.5				17.5	16	18		16.5	19			18	19.5			18.5	20	19			
16	15	1 1/2			17.5	16.5	1 1/2		18.5	17			3/4	19		18	19.5	18.5	20.5	19	19.5	18.5
16.5	15	3/4			18	17	1 1/2	19	17.5	Std	19.5	18.5	Std	20		19	3/4	20.5	n/a	3/4	7"	
17	15.5	16.5			15.5	18.5		17.5	19		18	19.5		18.5		20		19				
17	16	1 1/2	6"		17	16	Std	6"	19	18	1 1/2	7"	19	18	3/4	21	n/a	1 1/2	7"			
17.5	16.5				17.5	16.5	18.5		18.5	19.5			18.5	20.5		19	21			n/a		
17	15.5	Std			17.5	16	3/4		18	16.5			Std	20	19	1 1/2	20.5	n/a	8" x 2" caster not available with 6" fork			
17.5	16	3/4			18.5	17	19	17.5	19.5	18.5	20	18.5										
17.5	16.5	Std			7"	18	17	1 1/2	7"	19.5	18.5	3/4	7"	20.5	19	1 1/2	21	n/a				
18	17			18.5		17.5	1 1/2			19	18			19.5	18		20	18.5				
18.5	17	1 1/2		18.5		17.5	1 1/2	19		18	19.5	18		20	18.5							
19	17.5	Std		7"		19	17.5	Std	7"	20	19	1 1/2	7"	20.5	19	1 1/2	21	n/a				
19.5	18					3/4	18.5			17.5	3/4			19	18		20.5	19				
19.5	18.5	1 1/2				19	18	3/4		20.5	19	21		n/a								
20	19	1 1/2	7"			19.5	18.5	1 1/2	7"	6" x 2" caster not available with 5" fork				7"	20.5	19						
						19.5	18			20	18.5	20.5	19									
						20.5	19															

¹ Front Seat Height not available with 4 x 1" Poly Caster

² Not available with 4 x 1" Poly, 4 x 1.5" Poly or 4 x 1.5" Soft Roll Aluminum

Seat heights with "FSH" heading available with Standard Front Frame. Seat heights with "Super Low" heading available with Super Low Front Frame. Actual seat height may vary depending upon configuration.

Cushions

General Use Cushions

STEP 1: Choose Cushion

- 2" General Use Foam Cushion (14"x16" to 20"x20") E2601 ____ \$99
- 2" General Use Foam Cushion (Over 20"x20") E2601 ____ \$250

STEP 2: Choose Cushion Size

- Match to Chair Size
- Other Size _____

Matrix Seat Cushions

Cushion size is determined by chair dimensions. Certain sizes will have higher price and longer lead time; visit www.motionconcepts.com/matrix_international.html for size specifications.

STEP 1: Choose Cushion

- PS Cushion _____ \$414
Cushions 14" wide have a weight capacity of 150 lbs. Cushions 16" through 20" have a weight capacity of 300 lbs.
- Vi Cushion _____ \$510
Weight capacity of 300 lbs
- Libra Cushion _____ \$582
Weight capacity of 300 lbs

STEP 2: Choose Cushion Size

- Match to Chair Size
- Other Size - Provide Part Number _____

Other Matrix Products

Must provide part number

- _____

ROHO Seat Cushions

Cushion size is determined by chair dimensions. Certain sizes will have higher price and longer lead time; visit therohogroup.com for size specifications.

STEP 1: Choose Cushion

Cushions

For chairs under 22" wide

- HIGH PROFILE® Single Compartment E2622 ____ \$603
- MID PROFILE™ Single Compartment E2622 ____ \$603
- LOW PROFILE® Single Compartment E2622 ____ \$603
- QUADTRO SELECT® High Profile® E2624 ____ \$603
- QUADTRO SELECT® Mid Profile® E2624 ____ \$603
- HYBRID ELITE® Single Comp Smart Check Ready E2622 _ \$623

X Large Cushions

For 22" wide chairs

- HIGH PROFILE® Single Compartment X Large E2623 ____ \$670
- LOW PROFILE® Single Compartment X Large E2623 ____ \$670
- QUADTRO SELECT® High Profile® X Large E2625 ____ \$670

Cushions with Heavy Duty Cover

- HIGH PROFILE® Single Compartment w/ HD Cover E2622 _ \$629
- QUADTRO SELECT® High Profile® w/ HD Cover E2624 ____ \$629

STEP 2: Choose Cushion Size

- Match to Chair Size
- Other Size - Provide Part Number _____

Backs

Matrix Backs		Therafin Headrest for MATRIX Back Systems	
<input type="checkbox"/> Elite Back _____ \$750 <i>Not available with 22" Seat Width. Back Height of 14" not available with 20" Seat Width.</i>	<input type="checkbox"/> Headrest, Multi-Axis Offset Hardware with Smooth/Soft Pad (E1028, E0955) 34187 _____ \$420		
Back Height: <input type="checkbox"/> 10" <input type="checkbox"/> 12" <input type="checkbox"/> 14" <input type="checkbox"/> 16" <input type="checkbox"/> 18" <input type="checkbox"/> 20"	PAD SIZE: <input type="checkbox"/> 5 x 8 <input type="checkbox"/> 5 x 10 <input type="checkbox"/> 5 x 14		
<input type="checkbox"/> Elite Deep Back _____ \$912 <i>Not available with 22" Seat Width. Back Height of 10" not available with 14" Seat Width.</i>	Other Matrix Mounting Hardware		
Back Height: <input type="checkbox"/> 10" <input type="checkbox"/> 12" <input type="checkbox"/> 16" <input type="checkbox"/> 18" <input type="checkbox"/> 20"	Must provide part number		
<input type="checkbox"/> Mini Back _____ \$750 <i>Available only with 14" Seat Width. Comes standard with Kidabra cover.</i>	<input type="checkbox"/> _____ <i>For mounting hardware other than quick release</i>		
Back Height: <input type="checkbox"/> 7" <input type="checkbox"/> 10" <input type="checkbox"/> 12" <input type="checkbox"/> 14"	Other Matrix Products		
	Must provide part number		
	<input type="checkbox"/> _____		

ROHO Back Systems (Available for chairs up to 20" wide)

STEP 1: Choose Back

- | | |
|---|--|
| <input type="checkbox"/> AGILITY™ Min. Contour w/ Foam Only (E2611) _____ \$549
BACKREST HEIGHT: <input type="checkbox"/> 10" <input type="checkbox"/> 13" <input type="checkbox"/> 16" <input type="checkbox"/> 18" | <input type="checkbox"/> AGILITY™ Mid Contour w/ Air Only (E2620) _____ \$781
BACKREST HEIGHT: <input type="checkbox"/> 10" <input type="checkbox"/> 13" <input type="checkbox"/> 16" <input type="checkbox"/> 18" <input type="checkbox"/> 20" |
| <input type="checkbox"/> AGILITY™ Min. Contour w/ Foam & Lumbar Pad (E2613) _____ \$659
BACKREST HEIGHT: <input type="checkbox"/> 10" <input type="checkbox"/> 13" <input type="checkbox"/> 16" <input type="checkbox"/> 18" | <input type="checkbox"/> AGILITY™ Mid Contour w/ Air & Lumbar Pad (E2615) _____ \$781
BACKREST HEIGHT: <input type="checkbox"/> 10" <input type="checkbox"/> 13" <input type="checkbox"/> 16" <input type="checkbox"/> 18" <input type="checkbox"/> 20" |
| <input type="checkbox"/> AGILITY™ Min. Contour w/ Air Only (E2620) _____ \$769
BACKREST HEIGHT: <input type="checkbox"/> 10" <input type="checkbox"/> 13" <input type="checkbox"/> 16" <input type="checkbox"/> 18" | <input type="checkbox"/> AGILITY™ Max Contour w/ Air Only (E2620) _____ \$920
BACKREST HEIGHT: <input type="checkbox"/> 16" <input type="checkbox"/> 18" <input type="checkbox"/> 20" |
| <input type="checkbox"/> AGILITY™ Min. Contour w/ Air & Lumbar Pad (E2615) _____ \$769
BACKREST HEIGHT: <input type="checkbox"/> 10" <input type="checkbox"/> 13" <input type="checkbox"/> 16" <input type="checkbox"/> 18" | <input type="checkbox"/> AGILITY™ Max Contour w/ Air & Lumbar Pad (E2615) _____ \$920
BACKREST HEIGHT: <input type="checkbox"/> 16" <input type="checkbox"/> 18" <input type="checkbox"/> 20" |

STEP 2: Choose Back Hardware

- Fixed _____ N/C
Not compatible with 18" and 20" tall backs
- Quick Release _____ N/C

ROHO Laterals for AGILITY™ Back Systems

Sold in pairs

- Swing Away with Curved Pad (E1028/E0956) _____ \$1014
 PAD SIZE: Medium Large X-Large
 4 T x 4 L 4 T x 5 L 5 T x 5 L
- OFFSET TYPE: Flush Offset 1" Offset 2" Offset
- Fixed with Flat Pad (E0956) _____ \$452
 PAD SIZE: Medium Large X-Large
 4 T x 4 L 4 T x 5 L 5 T x 5 L
- OFFSET TYPE: Flush Offset

Therafin Headrest for AGILITY Back Systems

- Headrest, Multi-Axis Offset Hardware with Smooth/Soft Pad (E1028, E0955) 34187 _____ \$420
 PAD SIZE: 5 x 8 5 x 10 5 x 14



Vest

Not compatible with Back Upholstery

- Full Shape, Comfort Strap (E0960) _____ \$167
 - Small 30589 Medium 30590 Large 30591
 - X-Large 30592 2X-Large 30593
- Full Shape, Extended Strap (E0960) _____ \$167
 - Small 30594 Medium 30595 Large 30596
 - X-Large 30597 2X-Large 30598
- Trim Shape, Comfort Strap (E0960) _____ \$169
 - Small 30618 Medium 30619 Large 30620
 - X-Large 30621 2X-Large 30622
- Trim Shape, Extended Strap (E0960) _____ \$169
 - Small 30604 Medium 30624 Large 30625
 - X-Large 30626 2X-Large 30627

Shoulder Harness

With metal clips and cam locks. Not compatible with Back Upholstery.

- Small (E0960) 61403 _____ \$157
- Medium (E0960) 61503 _____ \$164
- Large (E0960) 61703 _____ \$170
- X-Large (E0960) 61803 _____ \$177

Hip Belt

- 4 Point Y-Style, Metal Buckle (E0978) _____ \$173
 - Small 60024 Medium 60025 Large 60026
- Dual Rear Pull, Metal Attaching Clips (E0978)
 - 1" Straps, 6.25"W x 1.75"H Pads 60161 _____ \$97
 - 1.5" Straps, 7.25"W x 2.25"H Pads 60171 _____ \$106
 - 1.5" Straps, 9.25"W x 2.5"H Pads 60181 _____ \$106
- Single Pull, Metal Buckle (E0978)
 - 1" Straps, 4.25"W x 1.75"H Pads 61001 _____ \$86
 - 1" Straps, 5.25"W x 1.75"H Pads 61101 _____ \$86
 - 1" Straps, 6.25"W x 1.75"H Pads 61201 _____ \$86
 - 1.5" Straps, 5.25"W x 2.25"H Pads 62001 _____ \$95
 - 1.5" Straps, 9.25"W x 2.5"H Pads 62201 _____ \$95
 - 2" Straps, 9.25"W x 3"H Pads 64001 _____ \$102
 - 2" Straps, 11.25"W x 3"H Pads 64101 _____ \$102

Cane Mounted Swing Away Lateral Thoracic Support

Sold as pair. Not compatible with Back Upholstery.

- Swing Away Kits w/ 4"x6" Gel Pads (E0956, E1028) 50021 _____ \$478
- Swing Away Kits, No Pads (E0956, E1028) 50022 _____ \$309

Knee Adductor

- Flip Down Hardware (E1028) _____ \$634
Sold in pairs, one right side and one left side
 - Short Arm, Fits 3/4", 7/8" and 1" Tubing 30248
 - Long Arm, Fits 3/4", 7/8" and 1" Tubing 30284
- Knee Adductor Pad with Neoprene Cover (E0956) _____ \$117
Sold single
 - Small Pad, 4.25" x 3.25" 50490 Quantity: _____
 - Medium Pad, 6.25" x 4.25" 50491 Quantity: _____

Foot / Shoe Holder

Compatible only with Aluminum Angle Adjustable and Composite Angle Adjustable Footplates

- FootSure Ankle Support (K0108) _____ \$178
Sold in pairs
 - X-Small w/ Buckle Closure, 6"-8" Circum. 60250
 - Small w/ Buckle Closure, 7"-8.5" Circum. 60251
 - Medium w/ Buckle Closure, 8"-11" Circum. 60252
 - Large w/ Buckle Closure, 9"-12" Circum. 60253
 - X-Large w/ Buckle Closure, 12.5"-15" Circum. 60254
 - X-Small w/ Hook-N-Loop Closure, 5"-6" Circum. 60255
 - Small w/ Hook-N-Loop Closure, 6"-8" Circum. 60256
 - Med w/ Hook-N-Loop Closure, 8"-10.5" Circum. 60257
 - Large w/ Hook-N-Loop Closure, 9"-12" Circum. 60258
 - X-Lrg w/ Hook-N-Loop Closure, 12"-15" Circum. 60259
- Padded Footbox (K0108)
 - Small, Full Design, Fits 14" Chair 30111 _____ \$287
 - Medium, Full Design, Fits 16" Chair 30112 _____ \$293
 - Large, Full Design, Fits 18" Chair 30113 _____ \$298
 - Small, Split Design, Fits 14" Chair 30152 _____ \$367
 - Medium, Split Design, Fits 16" Chair 30153 _____ \$373
 - Large, Split Design, Fits 18" Chair 30154 _____ \$378
- Molded Shoe Holders (K0108) _____ \$162
Sold in pairs
 - Small with Plain Straps, 3.5" x 5.25" 30293
 - Medium with Plain Straps, 4" x 6.5" 30303
 - Large with Plain Straps, 4.75" x 8" 30309
 - X-Large with Plain Straps, 5" x 10" 30316



Half Tray

Not compatible with Swing Away Armrests

- Standard Shape, Clear, Flip-Away (E0950, E1028) _____ \$246
 - Right Side 41310 Left Side 41311
- Standard Shape, Padded, Flip-Away (E0950, E1028) _____ \$269
 - Right Side 41350 Left Side 41353
- Trimline Shape, Padded, Flip-Away (E0950, E1028) _____ \$269
 - Right Side 41352 Left Side 41351

Clear Tray with Border and Pre-Drilled Holes

STEP 1: Choose Tray

Chair Width	Tray Dimension (W x D)	Body Cutout (W x D)	HCPCS	Part Number	Price	<input type="checkbox"/>
14"	20 x 18	10.5 x 6	E0950	33210	\$220	<input type="checkbox"/>
	20 x 18	11.5 x 6.5	E0950	33211	\$220	<input type="checkbox"/>
16"	22 x 21	12.5 x 7.5	E0950	33214	\$255	<input type="checkbox"/>
	22 x 21	13.5 x 8.5	E0950	33215	\$255	<input type="checkbox"/>
18"	24.5 x 22	14.5 x 9.5	E0950	33218	\$285	<input type="checkbox"/>
	24.5 x 22	15.5 x 10.5	E0950	33219	\$285	<input type="checkbox"/>
20"	26.5 x 22	17 x 11	E0950	33222	\$309	<input type="checkbox"/>
	26.5 x 22	18 x 11	E0950	33223	\$309	<input type="checkbox"/>
22"	28.5 x 23.5	18.5 x 11	E0950	33226	\$309	<input type="checkbox"/>
	28.5 x 23.5	20 x 11	E0950	33227	\$309	<input type="checkbox"/>

STEP 2: Choose Tray Hardware

- E-Z Lock Clamps, Tube Arm (E1028) 31072 _____ \$150
Compatible only with Swing Away Armrests. Sold in Pairs.
- E-Z Lock Clamps, Standard (E1028) 31098 _____ \$142
Not compatible with Swing Away Armrests. Sold in Pairs.

Communication Mount Hardware

Choose Device Holder AND Mount

STEP 1: What DEVICE will be used?

- Smart Phone 32359 _____ \$150
- Large Smart Phone 32354 _____ \$150
- I-Pad 32355 _____ \$184
- I-Pad Mini 32356 _____ \$178
- 7" Tablet 32357 _____ \$173
- 10" Tablet 32358 _____ \$395

STEP 2: Where does it need to MOUNT?

- Frame 32350 _____ \$512
- Armrest 32351 _____ \$395
Not compatible with Swing Away Armrests
- Tray 32352 _____ \$212

Cup Holder

- Heavy Duty 30889 _____ \$44
Compatible only with Flip Back Armrests