



Phone (800) 981-1540 | Fax (715) 254-0996
sales@kimobility.com | www.kimobility.com

CATALYST 4C

Back Upholstery

If selecting both upholstery and Axiom Back, the upholstery will be installed.

Standard	N/C
<i>Not available with Short or Medium Back Height Ranges on 21" and 22" Seat Widths</i>	
Tension Adjustable with Standard Cover (E2611)	\$323
Tension Adjustable with Air Mesh Breathable Cover (E2611)	\$415
Omit Back Upholstery	(\$10)

U.S. Military Service Patches

Attached only if ordered with Tension Adjustable Back Upholstery

Army	Navy	Marines	Air Force	
Coast Guard	US Purple Heart			\$92
Combo Army/Purple Heart		Combo Navy/Purple Heart		
Combo Marines/Purple Heart		Combo Air Force/Purple Heart		
Combo Coast Guard/Purple Heart				\$149

Footrest

Frame Type

Swing Away	N/C
Super Low - Swing Away	N/C

Footrest Hanger Type

Hanger	60°	70°	80°	90°		Retail
Extension Tube Mount			1			N/C
Front Mount			1	1		N/C
Pro Elevating Leg Rests E0990						\$359
Omit Hangers						(\$20)

1 Not available with 7" or 8" casters

PRO ELR Cover

Available only on PRO ELR

Gel Pad Covers (PAIR)	\$436
-----------------------	-------

Seat to Footrest Length

Select Extension and Seat to Footrest Length in row corresponding to desired Footrest. Footrest Length of at least 2.5" shorter than Front Seat Height recommended. Extension Tube Hanger measurements are based on 70° Hanger.

Measurements are based on 70" Hanger.														
Footrest	Extension	Seat to Footrest Length												
Extension Tube Hangers (60°,70°,80°)	Short (11" - 15")	11"	11.5"	12"	12.5"	13"	13.5"	14"	14.5"	15"				
	Medium (13" - 17")	13"	13.5"	14"	14.5"	15"	15.5"	16"	16.5"	17"				
	Long (15.5" - 19.5")	15.5"	16"	16.5"	17"	17.5"	18"	18.5"	19"	19.5"				
Front Mount ¹ Hangers 1			Standard Selections									Requires 3" Extension		
	Short Front Mount Hanger	4" ²	4.5" ²	5"	5.5"	6"	6.5"				10"	10.5"	11"	
		7"	7.5"	8"	8.5"	9"	9.5"				11.5"	12"	12.5"	
	Medium Front Mount Hanger	4" ²	4.5" ²	5"	5.5"	6"	6.5"	7"	7.5"		12"	12.5"	13"	
		8"	8.5"	9"	9.5"	10"	10.5"	11"	11.5"		13.5"	14"	14.5"	
	Long Front Mount Hanger	5"	5.5"	6"	6.5"	7"	7.5"	8"	8.5"	9"	14"	14.5"	15"	
		9.5"	10"	10.5"	11"	11.5"	12"	12.5"	13"	13.5"	15.5"	16"	16.5"	
PRO ELR Low Pivot Setting	PRO ELR Low Pivot Setting Short (12" - 17.5")	12" - 17.5"												
	PRO ELR Low Pivot Setting Long (15" - 21")	15" - 21"												
PRO ELR High Pivot Setting	PRO ELR High Pivot Setting Short (10.5" - 15.5")	10.5" - 15.5"												
	PRO ELR High Pivot Setting Long (13" - 19")	13" - 19"												

¹ See Minimum Front Seat Height chart above.

² Available with 90° only

3" Extension

Available only with Front Mount Hangers

3" Extension	\$36
3" Extension (Provided as Parts)	\$36

Footplates

Composite	N/C
<i>Not available with Front Mount Hanger</i>	
Composite Angle Adjustable (K0040)	\$205
<i>Not available with Front Mount Hanger or 14" Seat Width</i>	
Aluminum Angle Adjustable (K0040)	\$205
One Piece Flip Up Angle Adjustable (K0040)	\$205
<i>Not available with Front Mount Hanger or PRO ELR</i>	



Phone (800) 981-1540 | Fax (715) 254-0996
sales@kimobility.com | www.kimobility.com

CATALYST 4C

Hanger Release

If no selection is made Classic Hanger Release will be used

Classic	N/C
4-Way - In-Line	N/C
4-Way - Out	N/C

Calf Straps

Velcro Adjustable (K0038)	\$41
Padded Velcro Adjustable (K0038)	\$72
Bodypoint Padded AeroMesh (K0038)	\$133

Heel Loops

Not available with One Piece Flip Up Angle Adjustable Footplate

Heel Loops (E0951)	\$42
Adjustable Heel Loops (E0951)	\$51
Omit Heel Loops	N/C

Residual Limb Support

	Limb Support	Price	Omit	Cover		Pad Size					
Left	Left Limb Support (E1020)	\$431	Omit Left Hanger	Smooth Cover Left	Air Mesh Cover Left	5X8 Pad Left Side	7X10 Pad Left Side	9x10 Pad Left Side	9x14 Pad Left Side	11x10 Pad Left Side	11x14 Pad Left Side
Right	Right Limb Support (E1020)	\$431	Omit Right Hanger	Smooth Cover Right	Air Mesh Cover Right	5X8 Pad Right Side	7X10 Pad Right Side	9x10 Pad Right Side	9x14 Pad Right Side	11x10 Pad Right Side	11x14 Pad Right Side

Front Seat Height and Casters

Front Seat Height

Select Front Seat Height in row corresponding to desired caster size.

Note: Selections in first three Front Seat Height columns use Super Low Front Frame.

Caster Size	Front Seat Height													
	Uses Super Low Front Frame													
4"	13.5"	14"	14.5"	15"	15.5"	16"	16.5"	17"	17.5"	18"	18.5"	19"	19.5"	20"
5"	14.5"	15"	15.5"	16"	16.5"	17"	17.5"	18"	18.5"	19"	19.5"	20"	20.5"	
6"	15.5"	16"	16.5"	17"	17.5"	18"	18.5"	19"	19.5"	20"	20.5"	21"		
7"	17"	17.5"		18"	18.5"	19"	19.5"	20"	20.5"	21"				
8"	17.5"	18"	18.5"	19"	19.5"	20"	20.5"	21"						

Castors

		Caster Size					
		Price	4"	5"	6"	7"	8"
Caster Type	1" Poly	N/C					
	1" Pneumatic	\$92					
	1.5" Poly	\$62					
	1.5" Soft Roll Aluminum	\$205					
	2" Poly	\$97					
	2" Pneumatic w/Foam Insert	\$256					

Maximum Seat Dump Reference

(Front Seat Height minus Rear Seat Height)

Seat Depth	14"	15"	16"	17"	18"	19"	20"
Maximum Dump	1.5"	2"	2"	2"	2.5"	2.5"	2.5"

Fork / Stem Options

User Weight Under 80 lbs	N/C
User Weight 80-120 lbs	N/C
User Weight 120-160 lbs	N/C
User Weight 160-220 lbs	N/C
User Weight 220-260 lbs	N/C
User Weight 260-300 lbs	N/C

Caster Pin Locks (K0073) \$220

Not available with 1 1/2" Longer Fork Stem or Frog Legs Fork

Fork Stem

Standard Fork Stem	N/C
3/4" Longer Fork Stem	\$36
1 1/2" Longer Fork Stem	\$36
Frog Legs Suspension Fork	\$487

Available only with 4", 5" or 6" Caster Sizes. Not available with Caster Pin Locks, 350 lb Weight Capacity, 2" Poly caster or 2" Pneumatic w/ Foam Insert caster.



Phone (800) 981-1540 | Fax (715) 254-0996
sales@kimobility.com | www.kimobility.com

CATALYST 4C

Wheel Locks

Push to Lock	N/C
Push to Lock Flush Mount	N/C
Pull to Lock	N/C
Omit	N/C

Wheel Lock Options

Extension Handles (E0961)	\$82
<i>Available only with Push to Lock or Pull to Lock Wheel Locks</i>	
Grade Aids (E0974)	\$159
<i>Available only with Push to Lock Wheel Locks and Pneumatic or Pneumatic with Airless Insert tires</i>	

Rear Seat Height and Wheels

Wheel Size / Tire Profile / Rear Seat Height

Select Rear Seat Height in row corresponding to desired Wheel Size and Tire Profile. 14" Seat Depth not available with 24", 25" or 26" wheel sizes.

Wheel Size	Tire Profile	Rear Seat Height				
20"	Low	12.5" ¹	13.5"	14.5"	15.5"	16.5"
	Full	13" ¹	14"	15"	16"	17"
22"	Low	13.5" ¹	14.5"	15.5"	16.5"	17.5"
	Full	14" ¹	15"	16"	17"	18"
24"	Low	14.5" ¹	15.5"	16.5"	17.5"	18.5"
	Full	15" ¹	16"	17"	18"	19"
25"	Low	15" ¹	16"	17"	18"	19"
26"	Low	15.5" ¹	16.5"	17.5"	18.5"	19.5"
	Full	16" ¹	17"	18"	19"	20"

¹ Not available with Flip Back Armrests

Rear Wheel

		Rear Wheel Size					
		Price	20"	22"	24"	25"	26"
Rear Wheel Type	Spoke	N/C					
	Mag	N/C					
	MAXX Spoke	\$154					
	Superlight Spoke ¹	\$313				²	
	Spinergy Spox	\$938					

¹ Not available with 350 lb Weight Capacity

² Available only with Superlight and Plastic Coated Handrims

Spinergy Options

Select spoke color and handrim tab for Spinergy wheels

Black Spokes Red Spokes Blue Spokes White Spokes _____ N/C

Rear Tire Type

Full Profile

Street Pneumatic 1 3/8" (E2211, E2212)	N/C
<i>Not available on 25" wheel size</i>	
Full Poly 1 3/8" (E2218)	N/C
<i>Not available on Superlight wheel or 25" or 26" wheel sizes</i>	
Pneumatic with Airless Insert 1 3/8" (E2211, E2213)	\$108
<i>Not available on Superlight wheel or 25" wheel size</i>	

Low Profile

Low Poly 1" (E2218)	N/C
<i>Not available on Superlight wheel, 25" or 26" wheel sizes or 350 lb Weight Cap.</i>	
Sentinel Solid Tire 1" (E2218)	\$103
<i>Not available on Superlight wheel or 20" and 22" Mag wheels</i>	
High Pressure 1" (E2211, E2212)	\$103
<i>Not available on Mag wheel</i>	

Handrims

Standard Handrims

Aluminum Anodized	N/C
Superlight Aluminum Anodized	\$56
<i>Available only on 24" and 25" wheel sizes. Available only on MAXX Spoke, Superlight or Spinergy wheels. Not available with 350 lb Weight Capacity.</i>	
Plastic Coated	\$113
Projection (8 Vertical Projections)	\$359
<i>Available only with full profile tires in combination with 22", 24", and 26" MAXX Spoke wheels, or 22" and 24" Mag wheels</i>	
Omit Handrim	N/C

Ergonomic Handrims

Not available with 20" or 26" wheels. Tab mount only available with 24" wheels.

Ergo Lite Regular Grip	\$477
Ergo Standard Regular Grip	\$456
Ergo Lite Extra Grip	\$528
Ergo Standard Extra Grip	\$507
Ergo Lite Omit Grip	\$405
Ergo Standard Omit Grip	\$354

Handrim Options

Micro-Grip

3M Micro Grip	\$72
<i>Not available on Plastic Coated or Projection Handrims</i>	

Handrim Mount

Close Mount	\$31
<i>Chair will come with Standard Mount unless Close Mount is selected</i>	

Axles

Axle Bracket		Axle	
Offset Axle Bracket	N/C	Quick Release	N/C
Reverse Mount Offset Axle Bracket (Amputee)	N/C	Quad Quick Release	\$200
<i>Not available with Swing Away Armrests</i>		<i>Not available with Mag or Spinergy wheels</i>	

Armrests

Armrest Type	Length		Price
Height Adjustable T-Arm Tall (11.5"-15") E0973	Desk		\$323
	Full		
Height Adjustable T-Arm (8.5"-13.5") E0973	Desk		\$323
	Full		
Height Adjustable T-Arm Low (6.5"-10.5") E0973 ¹	Desk		\$323
	Full		
Height Adjustable Flip Back T-Arm Tall (11.5"-15") E0973	Desk		\$338
	Full		
Height Adjustable Flip Back T-Arm (9.5"-13.5") E0973	Desk		\$338
	Full		

Armrest Type	Length		Price
Flip back (10.5") ^{2,3}	Desk		N/C
	Full		
Swing Away (10.5"-12.5") E0973 ^{3,4}			\$251
Omit Armrests			N/C

Arm Pads

Not available with Swing Away Armrests

Waterfall Arm Pad

¹ If requested configuration is unavailable, Standard Height T-Arm will be used
² Not Available with 15" Rear Seat Height and 24" Wheels or Seat Depths less than 16"
³ Not Available with Angle Adjustable Back Posts
⁴ Not available with Arm Pads

Arm Pads

Not available with Swing Away Armrests

Available only with Height Adjustable T-Arm Tall, Height Adjustable T-Arm, and Height Adjustable T-Arm Low

Accessories

Side Guards		Anti-tips	
Composite Side Guards (K0108)	\$159	Rear Anti-Tips (E0971)	\$160
<i>Not available with T-Arm or Flip Back Armrests</i>		Omit	N/C

Accessories

Pelvic Positioning Belts

1 1/2" Auto Buckle Pelvic Positioning Belt (E0978)	\$87
1 1/2" Padded Auto Buckle Pelvic Positioning Belt (E0978)	\$179
2" Airline Buckle Pelvic Positioning Belt (E0978)	\$118
Small Bodypoint Evoflex (16-24") (E0978)	\$200
Medium Bodypoint Evoflex (19-29") (E0978)	\$200
Large Bodypoint Evoflex (25-37") (E0978)	\$200

Belt Mounting Band Clamp

Belt Mounting Band Clamp	\$36
--------------------------	------

O2 Holder

O2 Holder - Left Side (E2208)	\$313
O2 Holder - Right Side (E2208)	\$313

Other Accessories

Tool Kit	\$46
Removable Underseat Pouch	\$51
Ki Mobility Backpack	\$80
Black Neoprene Impact Guards	\$97
Clear Spoke Guards (K0065)	\$123
<i>Not available with Mag wheel</i>	
Black Spoke Guards (K0065)	\$123
<i>Not available with 20" wheel size or Mag wheel</i>	
Cane and Crutch Holder	\$210
IV Pole	\$292
<i>Available only with Standard Straight Back Posts</i>	



Phone (800) 981-1540 | Fax (715) 254-0996
sales@kimobility.com | www.kimobility.com

CATALYST 4C

Frame Color

Gloss Finish

Black	N/C
Blue Skies	N/C
Bubblegum Pink	N/C
Candy Blue	N/C
Candy Purple	N/C
Candy Red	N/C
Copper	N/C
Granny Smith	N/C
Metallic Orange	N/C
Yellow	N/C

Matte Finish

Military Green	N/C
Charcoal Gray	N/C
Emerald Green	N/C
Iridescent Burgundy	N/C
Iridescent Navy	N/C
Retro Red	N/C
Royal Blue	N/C
Shadow Black	N/C
Smooth White	N/C

Hydrographics

1 Select Hydrographics Base Color. Zebra not available in Black.

AOR1 (Desert Digital Camo)	
NEXT Camo	
Pink Camo	
Carbon Fiber 1	
Muddy Girl	
Zebra 1	\$825

Hydrographics Base Color

Black	N/C
Blue	N/C
Granny Smith	N/C
Pink	N/C
Purple	N/C
Red	N/C

Catalyst Front Seat Height Matrix																																		
4" Caster				5" Caster				6" Caster				7" Caster				8" Caster																		
FSH	Super Low	Stem	Fork	FSH	Super Low	Stem	Fork	FSH	Super Low	Stem	Fork	FSH	Super Low	Stem	Fork	FSH	Super Low	Stem	Fork															
15	13.5	Std	4"	16	14.5	Std	4"	17	15.5	Std	5"	18	16.5	Std	6"	19	17.5	Std	6"															
15.5	14			16.5	15.5	3/4		17.5	16.5	3/4		18.5	17			19.5	18.5	3/4																
15.5 ¹	14 ¹			17.5	16	1 1/2		18.5	17	1 1/2		18.5	17.5	20.5		n/a	1 1/2																	
16	14.5	3/4		16	14.5	Std	5"	17	15.5	Std	6"	19	18	1 1/2		19	17.5	Std	7"															
16.5 ¹	15 ¹			16.5	15			17.5	16			19.5	18																					
17	15.5	16.5		15.5	18			16.5	20			18.5	20			18.5																		
15.5	14	Std	17	16	3/4	17.5		16.5	3/4	18.5		17	Std	20.5	19	3/4																		
16	14.5		17.5	16		18		17		19		17.5		19.5	18.5																			
16	15	3/4	18	16.5		Std		18.5		17.5		1 1/2		19.5	18		3/4	20		19														
16.5	15.5		16.5	15	18.5			17	20	18.5			20.5	n/a																				
17	15.5	1 1/2	17	15.5	Std			19	17.5	1 1/2			19	18	3/4	21		n/a		1 1/2														
17.5	16		17.5	16		19.5		18	20.5			19																						
16.5	15	Std	17.5	16		3/4	18	16.5	Std		20	19	1 1/2																					
17	15.5		18	16.5	18.5		17	20.5		n/a																								
17	16	3/4	18.5	17	1 1/2		19	17.5		3/4	20	18.5						1 1/2																
17.5	16.5		18	16.5		19.5	18	20.5	19																									
18	16.5	1 1/2	18.5	17		3/4	18.5	17.5	1 1/2		21	n/a																						
18.5	17		19	17.5	19		18																											
18	16.5	Std	18	16.5	Std		19.5	18.5		3/4																								
18.5	17		18.5	17		20	19																											
18.5	17.5	3/4	19	17.5		1 1/2	19.5	18	1 1/2																									
19	18		18.5	17.5	20		18.5																											
19.5	18.5	1 1/2	19	18	3/4		20.5	19		1 1/2																								
20	19		19.5	18.5		21	n/a																											
				19.5		18	6" x 2" caster not available with 5" fork																											
				20	18.5																													
				20.5	19																													

¹ Not available with 4 x 1.5" Soft Roll Aluminum

Seat heights with "FSH" heading available with Standard Front Frame. Seat heights with "Super Low" heading available with Super Low Front Frame. Actual seat height may vary depending upon configuration.