



Phone (800) 981-1540 | Fax (715) 254-0996
 sales@kimobility.com | www.kimobility.com

CATALYST 4C

Effective 05/23/2017

Company: _____

Account Number: _____ Date: _____

Quote PO Number: _____

RTS: _____

Marked For: _____

Ship To: _____

Address: _____

City: _____ State: _____

Zip: _____ Phone: _____

Email/Fax: _____

Catalyst 4C K0004 **\$1525**
 250 lb weight capacity

*HCPCS codes provided are based on code verification by the PDAC or our interpretation of the code definitions.
 Proper use of HCPCS codes does not ensure coverage or payment. For coverage information, verify the appropriate payer's coverage policy.*

Frame

Seat Width / Seat Depth

Select box corresponding to desired Seat Width and Depth

		Seat Width									
		14"	15"	16"	17"	18"	19"	20"	21"	22"	
Seat Depth	Price HCPCS	\$45	\$45	N/C	\$45	N/C	\$150	\$395 E2201	\$495 E2201	\$495 E2201	
	14"	N/C	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	
	15"	\$70	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	
	16"	N/C	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	
	17"	\$90	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	
	18"	N/C	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	
	19"	\$195	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	
	20"	\$450 E2203	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	

Heavy Duty Option

350 lb Weight Capacity _____ \$500
 Requires minimum 6" caster height

Cross Brace

Standard Cross Braces _____ Std
 Seating System Cross Braces _____ N/C
 Allows for easy attachment of seating system mounting hardware. Not available with Seat Sling.

Backrest

Back Post / Height Range / Back Height

Select Back Height in row corresponding to desired Back Post and Height Range

Back Post	Price	Height Range	Back Height			
Standard	N/C	Short	<input type="checkbox"/> 8.5	<input type="checkbox"/> 9.5	<input type="checkbox"/> 10.5	<input type="checkbox"/> 11.5
		Medium	<input type="checkbox"/> 12.5	<input type="checkbox"/> 13.5	<input type="checkbox"/> 14.5	<input type="checkbox"/> 15.5
		Tall	<input type="checkbox"/> 16.5	<input type="checkbox"/> 17.5	<input type="checkbox"/> 18.5	<input type="checkbox"/> 19.5
Angle Adjustable w/ Push Handle	\$185	Short	<input type="checkbox"/> 14	<input type="checkbox"/> 15	<input type="checkbox"/> 16	
		Medium	<input type="checkbox"/> 16	<input type="checkbox"/> 17	<input type="checkbox"/> 18	
		Tall	<input type="checkbox"/> 18	<input type="checkbox"/> 19	<input type="checkbox"/> 20	

Standard Back Post

Select option below only if a Standard Back Post option is selected on chart at left

Straight with Push Handle _____ Std
 Straight without Push Handle _____ N/C
 8° Bend with Push Handle _____ \$25
 Not available with Short Back Height Range

Backrest continued

Back Upholstery

- Standard _____ Std
Not available with Short or Medium Back Height Ranges on 20" - 22" Seat Widths
- Tension Adjustable E2611 _____ \$245
- Omit Back Upholstery _____ - \$10

Footrest

Hanger Bend

- 60° _____ N/C
Available only with Extension Tube Hangers
- 70° _____ N/C
- 80° _____ N/C
Not available with 7" or 8" Casters
- 90° _____ N/C
Available only with Front Mount Hangers. Not available with 7" or 8" Casters.
- Elevating Leg Rests (ELR) E0990 _____ \$295
- Pro Elevating Leg Rests (PRO ELR) E0990 **New!** _____ \$350
- Omit Hangers _____ - \$20

Military Patch

- Attached only if ordered with Tension Adjustable Back Upholstery**
- US Military Service Patch _____ \$90
 Army Navy Marines Air Force Coast Guard
 - US Purple Heart Patch _____ \$90
 - Purple Heart with Military Service Patch Combo _____ \$145
 Army Navy Marines Air Force Coast Guard

Footplate

- Composite _____ N/C
Not available with Front Mount Hangers or PRO ELR
- Composite Angle Adjustable K0040 _____ \$185
Not available with Front Mount Hangers, ELR or 14" Seat Width
- Aluminum Angle Adjustable K0040 _____ \$185
- One Piece Flip Up Angle Adjustable K0040 _____ \$190
Not available with Front Mount Hangers, ELR or PRO ELR

Footrest / Extension / Seat to Footrest Length

Select Seat to Footrest Length in row corresponding to desired Footrest and Extension. Footrest Length of at least 2.5" shorter than Front Seat Height recommended. Extension Tube Hanger measurements are based on 70° Hanger / Front Frame Bend.

Footrest	Extension	Seat to Footrest Length													
Extension Tube Hangers (60°, 70°, 80°)	Short	<input type="checkbox"/> 11	<input type="checkbox"/> 11.5	<input type="checkbox"/> 12	<input type="checkbox"/> 12.5	<input type="checkbox"/> 13	<input type="checkbox"/> 13.5	<input type="checkbox"/> 14	<input type="checkbox"/> 14.5	<input type="checkbox"/> 15					
	Medium	<input type="checkbox"/> 13	<input type="checkbox"/> 13.5	<input type="checkbox"/> 14	<input type="checkbox"/> 14.5	<input type="checkbox"/> 15	<input type="checkbox"/> 15.5	<input type="checkbox"/> 16	<input type="checkbox"/> 16.5	<input type="checkbox"/> 17					
	Long	<input type="checkbox"/> 15.5	<input type="checkbox"/> 16	<input type="checkbox"/> 16.5	<input type="checkbox"/> 17	<input type="checkbox"/> 17.5	<input type="checkbox"/> 18	<input type="checkbox"/> 18.5	<input type="checkbox"/> 19	<input type="checkbox"/> 19.5					
Front Mount Hangers ¹	Standard Selections											Requires 3" Extension for \$35			
	Short (70°, 80°, 90°)	<input type="checkbox"/> 5	<input type="checkbox"/> 5.5	<input type="checkbox"/> 6	<input type="checkbox"/> 6.5	<input type="checkbox"/> 7							<input type="checkbox"/> 10	<input type="checkbox"/> 10.5	<input type="checkbox"/> 11
		<input type="checkbox"/> 7.5	<input type="checkbox"/> 8	<input type="checkbox"/> 8.5	<input type="checkbox"/> 9	<input type="checkbox"/> 9.5							<input type="checkbox"/> 11.5	<input type="checkbox"/> 12	<input type="checkbox"/> 12.5
	Medium (70°, 80°, 90°)	<input type="checkbox"/> 5	<input type="checkbox"/> 5.5	<input type="checkbox"/> 6	<input type="checkbox"/> 6.5	<input type="checkbox"/> 7	<input type="checkbox"/> 7.5	<input type="checkbox"/> 8					<input type="checkbox"/> 12	<input type="checkbox"/> 12.5	<input type="checkbox"/> 13
<input type="checkbox"/> 8.5		<input type="checkbox"/> 9	<input type="checkbox"/> 9.5	<input type="checkbox"/> 10	<input type="checkbox"/> 10.5	<input type="checkbox"/> 11	<input type="checkbox"/> 11.5					<input type="checkbox"/> 13.5	<input type="checkbox"/> 14	<input type="checkbox"/> 14.5	
Long (70°, 80°)	<input type="checkbox"/> 5	<input type="checkbox"/> 5.5	<input type="checkbox"/> 6	<input type="checkbox"/> 6.5	<input type="checkbox"/> 7	<input type="checkbox"/> 7.5	<input type="checkbox"/> 8	<input type="checkbox"/> 8.5	<input type="checkbox"/> 9			<input type="checkbox"/> 14	<input type="checkbox"/> 14.5	<input type="checkbox"/> 15	
	<input type="checkbox"/> 9.5	<input type="checkbox"/> 10	<input type="checkbox"/> 10.5	<input type="checkbox"/> 11	<input type="checkbox"/> 11.5	<input type="checkbox"/> 12	<input type="checkbox"/> 12.5	<input type="checkbox"/> 13	<input type="checkbox"/> 13.5	<input type="checkbox"/> 15.5	<input type="checkbox"/> 16	<input type="checkbox"/> 16.5			
Elevating Leg Rests (ELR)	Short	<input type="checkbox"/> 12	<input type="checkbox"/> 13	<input type="checkbox"/> 14	<input type="checkbox"/> 15	<input type="checkbox"/> 16									
	Long	<input type="checkbox"/> 13	<input type="checkbox"/> 14	<input type="checkbox"/> 15	<input type="checkbox"/> 16	<input type="checkbox"/> 17	<input type="checkbox"/> 18								
Extended Length ELR	Short	<input type="checkbox"/> 18	<input type="checkbox"/> 19	<input type="checkbox"/> 20	<input type="checkbox"/> 21	<input type="checkbox"/> 22									
	Long	<input type="checkbox"/> 18	<input type="checkbox"/> 19	<input type="checkbox"/> 20	<input type="checkbox"/> 21	<input type="checkbox"/> 22	<input type="checkbox"/> 23	<input type="checkbox"/> 24							
PRO ELR Low Pivot	Short	<input type="checkbox"/> 12 - 17.5													
	Long	<input type="checkbox"/> 15 - 21													
PRO ELR High Pivot	Short	<input type="checkbox"/> 10.5 - 15.5													
	Long	<input type="checkbox"/> 13 - 19													

¹ Available only with Aluminum Angle Adjustable Footplates. See Minimum Front Seat Height chart below.

Minimum Front Seat Height with Front Mount Hangers						
Hanger Length	Short		Medium		Long	
	70°	80/90°	70°	80/90°	70°	80°
Minimum Front Seat Height with Standard Selection	13.5"	13.5"	14"	15"	16"	17"
Minimum Front Seat Height with 3" Extension	14.5"	16"	16.5"	18"	18.5"	20"

3" Extension


- Available only with Front Mount Hangers. Provided as parts.
- 3" Extension for Front Mount Hanger _____ \$35

Footrest continued

Hanger Release

If no selection is made Classic Hanger Release will be used

- Classic _____ Std
- 4-Way - In line **New!** _____ N/C
See reference A to the right. Not available with ELR or PRO ELR.
- 4-Way - Out **New!** _____ N/C
See reference B to the right



Heel Loop

- Heel Loops E0951 _____ \$30
Not available with One Piece Flip Up Angle Adjustable Footplate
- Omit Heel Loops _____ N/C

Calf Strap

- Velcro Adjustable K0038 _____ \$35
- Padded Velcro Adjustable K0038 _____ \$68
- Bodypoint Padded AeroMesh K0038 _____ \$102

Residual Limb Support

- Residual Limb Support E1020
 - Left _____ \$395
 - Right _____ \$395
 - Omit Left Hanger _ - \$10
 - Omit Right Hanger _ - \$10
- Cover**
- Smooth Air Mesh
 - Smooth Air Mesh
- Pad Size**
- 5 x 8 9 x 14
 - 7 x 10 11 x 10
 - 9 x 10 11 x 14
 - 5 x 8 9 x 14
 - 7 x 10 11 x 10
 - 9 x 10 11 x 14

Front Seat Height and Casters

Caster Size / Front Frame Type / Front Seat Height

Select Front Seat Height in row corresponding to desired caster size.
Note: Selections in first three Front Seat Height columns use Super Low Front Frame.

Caster Size	Front Seat Height													
	Uses Super Low Front Frame													
4"	<input type="checkbox"/> 13.5 ¹	<input type="checkbox"/> 14	<input type="checkbox"/> 14.5	<input type="checkbox"/> 15 ²	<input type="checkbox"/> 15.5	<input type="checkbox"/> 16	<input type="checkbox"/> 16.5	<input type="checkbox"/> 17	<input type="checkbox"/> 17.5	<input type="checkbox"/> 18	<input type="checkbox"/> 18.5	<input type="checkbox"/> 19	<input type="checkbox"/> 19.5	<input type="checkbox"/> 20
5"	<input type="checkbox"/> 14.5	<input type="checkbox"/> 15	<input type="checkbox"/> 15.5	<input type="checkbox"/> 16	<input type="checkbox"/> 16.5	<input type="checkbox"/> 17	<input type="checkbox"/> 17.5	<input type="checkbox"/> 18	<input type="checkbox"/> 18.5	<input type="checkbox"/> 19	<input type="checkbox"/> 19.5	<input type="checkbox"/> 20	<input type="checkbox"/> 20.5	
6"	<input type="checkbox"/> 15.5	<input type="checkbox"/> 16	<input type="checkbox"/> 16.5	<input type="checkbox"/> 17	<input type="checkbox"/> 17.5	<input type="checkbox"/> 18	<input type="checkbox"/> 18.5	<input type="checkbox"/> 19	<input type="checkbox"/> 19.5	<input type="checkbox"/> 20	<input type="checkbox"/> 20.5	<input type="checkbox"/> 21		
7"	<input type="checkbox"/> 17	<input type="checkbox"/> 17.5		<input type="checkbox"/> 18	<input type="checkbox"/> 18.5	<input type="checkbox"/> 19	<input type="checkbox"/> 19.5	<input type="checkbox"/> 20	<input type="checkbox"/> 20.5	<input type="checkbox"/> 21				
8"	<input type="checkbox"/> 17.5	<input type="checkbox"/> 18	<input type="checkbox"/> 18.5	<input type="checkbox"/> 19	<input type="checkbox"/> 19.5	<input type="checkbox"/> 20	<input type="checkbox"/> 20.5	<input type="checkbox"/> 21						

¹ Front Seat Height not available with 4 x 1" Poly Caster
² Super Low Front Frame required when using 4 x 1" Poly Caster

Caster Type

- 1" Poly _____ N/C
- 1" Poly Aluminum _____ \$160
Not available with 6", 7" and 8" Caster Sizes
- 1" Pneumatic _____ \$80
Not available with 4", 5" and 7" Caster Sizes
- 1.5" Poly _____ \$50
Not available with 7" and 8" Caster Sizes
- 1.5" Soft Roll Aluminum _____ \$175
Not available with 7" and 8" Caster Sizes
- 2" Poly _____ \$70
Not available with 4", 5", 7" and 8" Caster Sizes
- 2" Pneumatic w/ Foam Insert _____ \$225
Not available with 4", 5", 6" and 7" Caster Sizes

Maximum Seat Dump Reference

(Front Seat Height minus Rear Seat Height)

Seat Depth	14"	15"	16"	17"	18"	19"	20"
Maximum Dump	1.5"	2"	2"	2"	2.5"	2.5"	2.5"

Fork Stem

- Standard Fork Stem _____ Std
- 3/4" Longer Fork Stem _____ \$35
- 1 1/2" Longer Fork Stem _____ \$35

Fork Options

- Frog Legs Ultra-Sport Fork - Black _____ \$425
Available only with 4", 5" or 6" Casters. Not available with Caster Pin Locks, Heavy Duty Option, 2" Poly caster or 2" Pneumatic w/ Foam Insert caster.

End User Weight Required: _____ lbs

- Frog Legs PHASE2 Carbon Fork **New!** _____ \$425
Lightweight shock absorbing fork. Available only with 4", 5" or 6" Casters. Not available with Caster Pin Locks, Heavy Duty Option, 2" Poly caster or 2" Pneumatic w/ Foam Insert caster.

- Caster Pin Locks K0073 _____ \$195
Not available with 1 1/2" Longer Fork Stem or Frog Legs Fork

Wheel Locks

- Wheel Lock**
- Push to Lock _____ Std
 - Push to Lock Flush Mount _____ N/C
 - Pull to Lock _____ N/C
 - Omit Wheel Locks _____ N/C

- Wheel Lock Options**
- Extension Handles E0961 _____ \$67
Available only with Push to Lock or Pull to Lock Wheel Locks
 - Grade Aids E0974 _____ \$150
Available only with Push to Lock Wheel Locks and Pneumatic or Pneumatic with Airless Insert Rear Tire Types

Rear Seat Height and Wheels

Wheel Size / Tire Profile / Rear Seat Height

Select Rear Seat Height in row corresponding to desired Wheel Size and Tire Profile. 14" Seat Depth not available with 24", 25" or 26" Rear Wheels.

Wheel Size	Tire Profile	Rear Seat Height				
20"	Low	<input type="checkbox"/> 12.5 ¹	<input type="checkbox"/> 13.5	<input type="checkbox"/> 14.5	<input type="checkbox"/> 15.5	<input type="checkbox"/> 16.5
	Full	<input type="checkbox"/> 13 ¹	<input type="checkbox"/> 14	<input type="checkbox"/> 15	<input type="checkbox"/> 16	<input type="checkbox"/> 17
22"	Low	<input type="checkbox"/> 13.5 ¹	<input type="checkbox"/> 14.5	<input type="checkbox"/> 15.5	<input type="checkbox"/> 16.5	<input type="checkbox"/> 17.5
	Full	<input type="checkbox"/> 14 ¹	<input type="checkbox"/> 15	<input type="checkbox"/> 16	<input type="checkbox"/> 17	<input type="checkbox"/> 18
24"	Low	<input type="checkbox"/> 14.5 ¹	<input type="checkbox"/> 15.5	<input type="checkbox"/> 16.5	<input type="checkbox"/> 17.5	<input type="checkbox"/> 18.5
	Full	<input type="checkbox"/> 15 ¹	<input type="checkbox"/> 16	<input type="checkbox"/> 17	<input type="checkbox"/> 18	<input type="checkbox"/> 19
25"	Low	<input type="checkbox"/> 15 ¹	<input type="checkbox"/> 16	<input type="checkbox"/> 17	<input type="checkbox"/> 18	<input type="checkbox"/> 19
26"	Low	<input type="checkbox"/> 15.5 ¹	<input type="checkbox"/> 16.5	<input type="checkbox"/> 17.5	<input type="checkbox"/> 18.5	<input type="checkbox"/> 19.5
	Full	<input type="checkbox"/> 16 ¹	<input type="checkbox"/> 17	<input type="checkbox"/> 18	<input type="checkbox"/> 19	<input type="checkbox"/> 20

¹ Not available with Flip Back Armrests

Rear Wheel Type

- Assembled 24" Six-Spoke Mag with Treaded Tire and Composite Handrim _____ - \$25
No substitutions available. Do not select Rear Tire Type or Handrim. Measures as Full Tire.
 - 24" Spoke _____ N/C
 - Spoke (20", 22", 25" or 26") _____ \$100
 - Mag _____ N/C
Not available in 25" or 26" rear wheels
 - Superlight Spoke _____ \$295
Available only in 24" and 25" rear wheels. Not available with Heavy Duty Option.
 - Spinergy Spox _____ \$795
Not available in 20" rear wheel
- Color: Black Red White Blue Yellow
- Handrim: Short Tab Long Tab
Available only with Aluminum Anodized and Plastic Coated Handrims

Rear Tire Type

- Full Profile**
- Street Pneumatic 1 3/8" E2211, E2212 _____ Std
Not available on 25" rear wheel
 - Full Poly 1 3/8" E2218 _____ N/C
Not available on Superlight or 25" rear wheel
 - Pneumatic with Airless Insert 1 3/8" E2211, E2213 _____ \$95
Not available on Superlight or 25" rear wheel
 - Schwalbe HS-127 1 3/8" E2211, E2212 **New!** _____ \$125
Available only on 24" rear wheel. Not available on Mag wheels.
 - Euro Trek Knobby 1 3/8" E2211, E2212 _____ \$95
Available only on 24" rear wheel

- Low Profile**
- Low Poly 1" E2218 _____ N/C
Not available on Superlight, 25" or 26" rear wheels
 - Sentinel Solid Tire 1" E2218 _____ \$95
Not available on Superlight wheels
 - High Pressure 1" E2211, E2212 _____ \$95
Not available on Mag wheels

Handrim

- Standard Handrims**
- Aluminum Anodized _____ Std
 - Superlight Aluminum Anodized _____ \$50
Available only on 24" and 25" Superlight and Spinergy wheels. Not available with Heavy Duty Option.
 - Plastic Coated _____ \$95
 - Projection (8 Vertical Projections) _____ \$295
Available on 22", 24", 25" and 26" Spoke wheels & 24" Mag wheels w/full tires
 - Omit Handrim _____ N/C

- Natural Fit Handrims**
- Not available on 20" rear wheel or 22" Mag wheels. Not recommended with Mag wheels.*
- Natural Fit LT Regular Grip _____ \$445
 - Natural Fit Standard Regular Grip _____ \$445
 - Natural Fit LT Super Grip _____ \$495
 - Natural Fit Standard Super Grip _____ \$495
 - Natural Fit LT Omit Grip _____ \$345
 - Natural Fit Standard Omit Grip _____ \$345

Handrim Option

- 3M Micro Grip _____ \$65
Not available on Plastic Coated or Projection Handrims

Axles

Axle Bracket	
<input type="checkbox"/> Offset Axle Bracket _____ Std	
<input type="checkbox"/> Reverse Mount Offset Axle Bracket _____ N/C	<i>Not available with Swing Away Armrests</i>

Axle	
<input type="checkbox"/> Quick Release _____ Std	
<input type="checkbox"/> Quad Quick Release _____ \$195	<i>Not available with Mag wheels</i>

Armrests

Armrest Please select Arm Pad option	
<input type="checkbox"/> Height Adjustable T-Arm (8.5" - 13.5") E0973 _____ \$265	
<input type="checkbox"/> Height Adjustable T-Arm Low (7.5" - 11.5") E0973 _____ \$265	
<input type="checkbox"/> Height Adjustable Flip Back (9.5" - 13.5") E0973 _____ \$275	<i>Not available with Seat Depths less than 16" or Angle Adjustable Back Posts</i>
<input type="checkbox"/> Flip Back (10.5") _____ N/C	<i>Not available with Seat Depths less than 16" or Angle Adjustable Back Posts</i>

Armrest continued	
<input type="checkbox"/> Swing Away (10.5" - 12.5") E0973 _____ \$195	<i>Not available with Angle Adjustable Back Posts or Arm Pads</i>
<input type="checkbox"/> Omit Armrests _____ N/C	

Arm Pad	
<i>Not available with Swing Away Armrests</i>	
Standard Arm Pad _____ N/C	<input type="checkbox"/> Desk <input type="checkbox"/> Full
Waterfall Arm Pad _____ \$25	<input type="checkbox"/> Desk <input type="checkbox"/> Full
<i>Available only with Height Adjustable T-Arms</i>	

Accessories

Removable Side Guard	
<input type="checkbox"/> Composite Side Guards K0108 _____ \$140	<i>Not available with T-Arm or Flip Back Armrests</i>

Seat Pan	
<input type="checkbox"/> Removable Aluminum Seat Pan E2231 _____ \$245	<i>Pan depth 1/2" shorter than Seat Depth</i>

Anti-Tip	
<input type="checkbox"/> Rear Anti-Tips E0971 _____ \$140	

Other Accessories	
<input type="checkbox"/> Tool Kit _____ \$45	
<input type="checkbox"/> Removable Underseat Pouch _____ \$50	
<input type="checkbox"/> Ki Mobility Backpack _____ \$65	
<input type="checkbox"/> Black Neoprene Impact Guards _____ \$85	
<input type="checkbox"/> Clear Spoke Guards K0065 _____ \$115	<i>Not available with Mag wheels</i>
<input type="checkbox"/> Black Spoke Guards K0065 _____ \$115	<i>Not available with 20" rear wheel or Mag wheels</i>
<input type="checkbox"/> Cane and Crutch Holder _____ \$195	
<input type="checkbox"/> IV Pole _____ \$265	<i>Available only with Standard Straight Back Posts</i>

Positioning Belt	
<input type="checkbox"/> 1 1/2" Auto Buckle E0978 _____ \$70	
<input type="checkbox"/> 1 1/2" Padded Auto Buckle E0978 _____ \$155	
<input type="checkbox"/> 2" Airline Buckle E0978 _____ \$110	
<input type="checkbox"/> Bodypoint Evoflex - Pelvic Stabilizer E0978 _____ \$183	Size: <input type="checkbox"/> Small (16-24") <input type="checkbox"/> Medium (19-29") <input type="checkbox"/> Large (25-37")
<input type="checkbox"/> Bodypoint Belt Mounting Kit _____ \$105	Style: <input type="checkbox"/> Bracket <input type="checkbox"/> Band Clamp 1"

O ₂ Holder	
<input type="checkbox"/> O ₂ Holder E2208 _____ \$295	<input type="checkbox"/> Right Side <input type="checkbox"/> Left Side

Frame Color

Gloss Finish	
<input type="checkbox"/> Black _____ N/C	
<input type="checkbox"/> Blue Skies _____ N/C	
<input type="checkbox"/> Bubblegum Pink _____ N/C	
<input type="checkbox"/> Candy Blue _____ N/C	
<input type="checkbox"/> Candy Purple _____ N/C	
<input type="checkbox"/> Candy Red _____ N/C	
<input type="checkbox"/> Copper _____ N/C	
<input type="checkbox"/> Fuchsia _____ N/C	
<input type="checkbox"/> Granny Smith _____ N/C	
<input type="checkbox"/> Yellow _____ N/C	

Matte Finish	
<input type="checkbox"/> Army Green _____ N/C	
<input type="checkbox"/> Charcoal Gray _____ N/C	
<input type="checkbox"/> Emerald Green _____ N/C	
<input type="checkbox"/> Iridescent Burgundy _____ N/C	
<input type="checkbox"/> Iridescent Navy _____ N/C	
<input type="checkbox"/> Retro Red _____ N/C	
<input type="checkbox"/> Royal Blue _____ N/C	
<input type="checkbox"/> Shadow Black _____ N/C	
<input type="checkbox"/> Smooth White _____ N/C	

Graphic Image Transfer	
<i>Visit www.kimobility.com for current graphic image transfer options</i>	
<i>Note selection below</i>	
<input type="checkbox"/> _____	\$595

Catalyst Front Seat Height Matrix																						
4" Caster				5" Caster				6" Caster				7" Caster				8" Caster						
FSH	Super Low	Stem	Fork	FSH	Super Low	Stem	Fork	FSH	Super Low	Stem	Fork	FSH	Super Low	Stem	Fork	FSH	Super Low	Stem	Fork			
15 ¹	13.5 ¹	Std	4"	16	14.5	Std	4"	17	15.5	Std	5"	18	17	Std	6"	19	18	Std	6"			
15.5	14			16.5	15.5	3/4		18	16.5	3/4		18.5	17.5			20	18.5	3/4				
16 ²	14.5 ²	17.5		16	1 1/2	18.5		17	1 1/2	19		17.5	20.5	n/a		1 1/2						
16.5	15	3/4		16	15	Std	5"	17	16	Std	6"	19.5	18	3/4		19	17.5	Std	7"			
16.5 ²	15 ²	1 1/2		16.5	15.5			17.5	16.5			18.5	18.5			19.5	18.5			19.5	18	
17	15.5	17		15.5	3/4	18		17	18			17	20	19		20	18.5					
15.5	14.5	Std		5"	17.5	16		3/4	6"	18	16.5	3/4	7"	18.5		17.5	Std	20.5	19	Std	7"	
16	15				17.5	16.5				1 1/2	18.5			17		19		18	19.5			18.5
16.5	15	3/4			18	17				19	17.5			20		19	20.5	19	20.5			n/a
17	15.5	1 1/2			5"	16.5		15.5	Std	6"	18.5	17.5	1 1/2	7"		19	18	3/4	20	19	1 1/2	7"
17	16		17			16		19			18	20			19							
17.5	16.5	17.5	16.5			19.5		18.5			21	n/a										
17	15.5	Std	6"			17.5		16	3/4	6"	18	16.5	Std	7"	20	19	1 1/2	8" x 2" caster not available with 6" fork				
17.5	16					18	16.5	18.5			17	20.5			n/a							
17.5	16.5	3/4				18.5	17	19			17.5	20			18.5							
18	17	1 1/2				6"	18	17	1 1/2	6"	19.5	18.5	3/4	7"	20.5	19	1 1/2					
18.5	17			18.5			17.5	18.5			17.5	21			n/a							
19	17.5	19		18			19	18														
18	17	Std		7"			18	16.5	Std	7"	19.5	18.5	1 1/2	7"								
18.5	17.5				18.5		17	20			19											
19	17.5	19			17.5		19.5	18														
19.5	18	3/4			7"		18.5	17.5	3/4	7"	20	18.5	1 1/2	7"								
19.5	18.5	19	18				20.5	19														
20	19	1 1/2	19.5				18.5	21			n/a											
							19.5	18	1 1/2	7"	6" x 2" caster not available with 5" fork				7"	7"						
						20	18.5	20.5			19											
						20.5	19															

¹ Front Seat Height not available with 4 x 1" Poly Caster

² Not available with 4 x 1" Poly, 4 x 1.5" Poly or 4 x 1.5" Soft Roll Aluminum

Seat heights with "FSH" heading available with Standard Front Frame. Seat heights with "Super Low" heading available with Super Low Front Frame. Actual seat height may vary depending upon configuration.

Cushions

General Use Cushions

STEP 1: Choose Cushion

- 2" General Use Foam Cushion (14"x14" to 20"x20") E2601 ____ \$99
- 2" General Use Foam Cushion (Over 20"x20") E2601 ____ \$250

STEP 2: Choose Cushion Size

- Match to Chair Size
- Other Size _____

Matrix Seat Cushions

Cushion size is determined by chair dimensions. Certain sizes will have higher price and longer lead time; visit www.motionconcepts.com/matrix_international.html for size specifications.

STEP 1: Choose Cushion

- PS Cushion _____ \$414
Cushions 14" wide have a weight capacity of 150 lbs. Cushions 15" through 20" have a weight capacity of 300 lbs.
- Vi Cushion _____ \$510
Weight capacity of 300 lbs
- Libra Cushion _____ \$582
Weight capacity of 300 lbs

STEP 2: Choose Cushion Size

- Match to Chair Size
- Other Size - Provide Part Number _____

Other Matrix Products

Must provide part number

- _____

ROHO Seat Cushions

Cushion size is determined by chair dimensions. Certain sizes will have higher price and longer lead time; visit therohogroup.com for size specifications.

STEP 1: Choose Cushion

Standard Cushions

For chairs under 22" wide

- HIGH PROFILE® Single Compartment (E2622) ____ \$603
- MID PROFILE™ Single Compartment (E2622) ____ \$603
- LOW PROFILE® Single Compartment (E2622) ____ \$603
- QUADTRO SELECT® High Profile® (E2624) ____ \$603
- QUADTRO SELECT® Mid Profile® (E2624) ____ \$603
- HYBRID ELITE® Single Comp Smart Check Ready (E2622) ____ \$623
- CONTOUR SELECT® (E2624) ____ \$603

X Large Cushions

For 22" wide chairs

- HIGH PROFILE® Single Compartment X Large (E2623) ____ \$670
- LOW PROFILE® Single Compartment X Large (E2623) ____ \$670
- QUADTRO SELECT® High Profile® X Large (E2625) ____ \$670

Cushions with Heavy Duty Cover

- HIGH PROFILE® Single Compartment w/ HD Cover (E2622) ____ \$629
- QUADTRO SELECT® High Profile® w/ HD Cover (E2624) ____ \$629

STEP 2: Choose Cushion Size

- Match to Chair Size
- Other Size - Provide Part Number _____

Backs

Matrix Backs

Elite Back _____ \$750
Not available with 22" Seat Width. Back Height of 14" not available with 19", 20" and 21" Seat Widths.

Back Height:

10" 12" 14" 16" 18" 20"

Elite Deep Back _____ \$912
Not available with 22" Seat Width. Back Height of 10" not available with 14" Seat Width.

Back Height:

10" 12" 16" 18" 20"

Mini Back _____ \$750
Available only with 14" and 15" Seat Widths. Comes standard with Kidabra cover.

Back Height:

7" 10" 12" 14"

Therafin Headrest for MATRIX Back Systems

Headrest, Multi-Axis Offset Hardware with Smooth/Soft Pad (E1028, E0955) 34187 _____ \$420

PAD SIZE: 5 x 8 5 x 10 5 x 14

Other Matrix Mounting Hardware

Must provide part number

For mounting hardware other than quick release

Other Matrix Products

Must provide part number

ROHO Back Systems (Available for chairs up to 20" wide)

STEP 1: Choose Back

AGILITY™ Min. Contour w/ Foam Only (E2611) _____ \$549
 BACKREST HEIGHT: 10" 13" 16" 18"

AGILITY™ Min. Contour w/ Foam & Lumbar Pad (E2613) _____ \$659
 BACKREST HEIGHT: 10" 13" 16" 18"

AGILITY™ Min. Contour w/ Air Only (E2620) _____ \$769
 BACKREST HEIGHT: 10" 13" 16" 18"

AGILITY™ Min. Contour w/ Air & Lumbar Pad (E2615) _____ \$769
 BACKREST HEIGHT: 10" 13" 16" 18"

AGILITY™ Mid Contour w/ Air Only (E2620) _____ \$781
 BACKREST HEIGHT: 10" 13" 16" 18" 20"

AGILITY™ Mid Contour w/ Air & Lumbar Pad (E2615) _____ \$781
 BACKREST HEIGHT: 10" 13" 16" 18" 20"

AGILITY™ Max Contour w/ Air Only (E2620) _____ \$920
 BACKREST HEIGHT: 16" 18" 20"

AGILITY™ Max Contour w/ Air & Lumbar Pad (E2615) _____ \$920
 BACKREST HEIGHT: 16" 18" 20"

STEP 2: Choose Back Hardware

Fixed _____ N/C
Not compatible with 18" and 20" tall backs

Quick Release _____ N/C

ROHO Laterals for AGILITY™ Back Systems

Sold in pairs

Swing Away with Curved Pad (E1028/E0956) _____ \$1014
 PAD SIZE: Medium Large X-Large
 4 T x 4 L 4 T x 5 L 5 T x 5 L

OFFSET TYPE: Flush Offset 1" Offset 2" Offset

Fixed with Flat Pad (E0956) _____ \$452
 PAD SIZE: Medium Large X-Large
 4 T x 4 L 4 T x 5 L 5 T x 5 L

OFFSET TYPE: Flush Offset

Therafin Headrest for AGILITY Back Systems

Headrest, Multi-Axis Offset Hardware with Smooth/Soft Pad (E1028, E0955) 34187 _____ \$420

PAD SIZE: 5 x 8 5 x 10 5 x 14



Vest

Not compatible with Back Upholstery

- Full Shape, Comfort Strap (E0960) _____ \$167
 - Small 30589 Medium 30590 Large 30591
 - X-Large 30592 2X-Large 30593
- Full Shape, Extended Strap (E0960) _____ \$167
 - Small 30594 Medium 30595 Large 30596
 - X-Large 30597 2X-Large 30598
- Trim Shape, Comfort Strap (E0960) _____ \$169
 - Small 30618 Medium 30619 Large 30620
 - X-Large 30621 2X-Large 30622
- Trim Shape, Extended Strap (E0960) _____ \$169
 - Small 30604 Medium 30624 Large 30625
 - X-Large 30626 2X-Large 30627

Shoulder Harness

With metal clips and cam locks. Not compatible with Back Upholstery.

- Small (E0960) 61403 _____ \$157
- Medium (E0960) 61503 _____ \$164
- Large (E0960) 61703 _____ \$170
- X-Large (E0960) 61803 _____ \$177

Hip Belt

- 4 Point Y-Style, Metal Buckle (E0978) _____ \$173
 - Small 60024 Medium 60025 Large 60026
- Dual Rear Pull, Metal Attaching Clips (E0978)
 - 1" Straps, 6.25"W x 1.75"H Pads 60161 _____ \$97
 - 1.5" Straps, 7.25"W x 2.25"H Pads 60171 _____ \$106
 - 1.5" Straps, 9.25"W x 2.5"H Pads 60181 _____ \$106
- Single Pull, Metal Buckle (E0978)
 - 1" Straps, 4.25"W x 1.75"H Pads 61001 _____ \$86
 - 1" Straps, 5.25"W x 1.75"H Pads 61101 _____ \$86
 - 1" Straps, 6.25"W x 1.75"H Pads 61201 _____ \$86
 - 1.5" Straps, 5.25"W x 2.25"H Pads 62001 _____ \$95
 - 1.5" Straps, 9.25"W x 2.5"H Pads 62201 _____ \$95
 - 2" Straps, 9.25"W x 3"H Pads 64001 _____ \$102
 - 2" Straps, 11.25"W x 3"H Pads 64101 _____ \$102

Cane Mounted Swing Away Lateral Thoracic Support

Sold as pair. Not compatible with Back Upholstery.

- Swing Away Kits w/ 4"x6" Gel Pads (E0956, E1028) 50021 _____ \$478
- Swing Away Kits, No Pads (E0956, E1028) 50022 _____ \$309

Knee Adductor

- Flip Down Hardware (E1028) _____ \$634
Sold in pairs, one right side and one left side
 - Short Arm, Fits 3/4", 7/8" and 1" Tubing 30248
 - Long Arm, Fits 3/4", 7/8" and 1" Tubing 30284
- Knee Adductor Pad with Neoprene Cover (E0956) _____ \$117
Sold single
 - Small Pad, 4.25" x 3.25" 50490 Quantity: _____
 - Medium Pad, 6.25" x 4.25" 50491 Quantity: _____

Foot / Shoe Holder

- Compatible only with Aluminum Angle Adjustable and Composite Angle Adjustable Footplates**
- FootSure Ankle Support (K0108) _____ \$178
Sold in pairs
 - X-Small w/ Buckle Closure, 6"-8" Circum. 60250
 - Small w/ Buckle Closure, 7"-8.5" Circum. 60251
 - Medium w/ Buckle Closure, 8"-11" Circum. 60252
 - Large w/ Buckle Closure, 9"-12" Circum. 60253
 - X-Large w/ Buckle Closure, 12.5"-15" Circum. 60254
 - X-Small w/ Hook-N-Loop Closure, 5"-6" Circum. 60255
 - Small w/ Hook-N-Loop Closure, 6"-8" Circum. 60256
 - Med w/ Hook-N-Loop Closure, 8"-10.5" Circum. 60257
 - Large w/ Hook-N-Loop Closure, 9"-12" Circum. 60258
 - X-Lrg w/ Hook-N-Loop Closure, 12"-15" Circum. 60259
 - Padded Footbox (K0108)
 - Small, Full Design, Fits 14" Chair 30111 _____ \$287
 - Medium, Full Design, Fits 16" Chair 30112 _____ \$293
 - Large, Full Design, Fits 18" Chair 30113 _____ \$298
 - Small, Split Design, Fits 14" Chair 30152 _____ \$367
 - Medium, Split Design, Fits 16" Chair 30153 _____ \$373
 - Large, Split Design, Fits 18" Chair 30154 _____ \$378
 - Molded Shoe Holders (K0108) _____ \$162
Sold in pairs
 - Small with Plain Straps, 3.5" x 5.25" 30293
 - Medium with Plain Straps, 4" x 6.5" 30303
 - Large with Plain Straps, 4.75" x 8" 30309
 - X-Large with Plain Straps, 5" x 10" 30316



Half Tray

Not compatible with Swing Away Armrests

- Standard Shape, Clear, Flip-Away (E0950, E1028) _____ \$246
 - Right Side 41310 Left Side 41311
- Standard Shape, Padded, Flip-Away (E0950, E1028) _____ \$269
 - Right Side 41350 Left Side 41353
- Trimline Shape, Padded, Flip-Away (E0950, E1028) _____ \$269
 - Right Side 41352 Left Side 41351

Clear Tray with Border and Pre-Drilled Holes

STEP 1: Choose Tray

Chair Width	Tray Dimension (W x D)	Body Cutout (W x D)	HCPCS	Part Number	Price	
14"	20 x 18	10.5 x 6	E0950	33210	\$220	<input type="checkbox"/>
	20 x 18	11.5 x 6.5	E0950	33211	\$220	<input type="checkbox"/>
15"	21 x 18	11.5 x 6.5	E0950	33212	\$220	<input type="checkbox"/>
	21 x 18	12.5 x 7	E0950	33213	\$220	<input type="checkbox"/>
16"	22 x 21	12.5 x 7.5	E0950	33214	\$255	<input type="checkbox"/>
	22 x 21	13.5 x 8.5	E0950	33215	\$255	<input type="checkbox"/>
17"	23.5 x 22	13.5 x 8.5	E0950	33216	\$255	<input type="checkbox"/>
	23.5 x 22	14.5 x 9.5	E0950	33217	\$285	<input type="checkbox"/>
18"	24.5 x 22	14.5 x 9.5	E0950	33218	\$285	<input type="checkbox"/>
	24.5 x 22	15.5 x 10.5	E0950	33219	\$285	<input type="checkbox"/>
19"	26.5 x 22	15.5 x 10.5	E0950	33220	\$285	<input type="checkbox"/>
	26.5 x 22	16.5 x 11	E0950	33221	\$285	<input type="checkbox"/>
20"	26.5 x 22	17 x 11	E0950	33222	\$309	<input type="checkbox"/>
	26.5 x 22	18 x 11	E0950	33223	\$309	<input type="checkbox"/>
21"	27.5 x 23.5	17.5 x 11	E0950	33224	\$309	<input type="checkbox"/>
	27.5 x 23.5	19 x 11	E0950	33225	\$309	<input type="checkbox"/>
22"	28.5 x 23.5	18.5 x 11	E0950	33226	\$309	<input type="checkbox"/>
	28.5 x 23.5	20 x 11	E0950	33227	\$309	<input type="checkbox"/>

STEP 2: Choose Tray Hardware

- E-Z Lock Clamps, Tube Arm (E1028) 31072 _____ \$150
Compatible only with Swing Away Armrests. Sold in Pairs.
- E-Z Lock Clamps, Standard (E1028) 31098 _____ \$142
Not compatible with Swing Away Armrests. Sold in Pairs.

Communication Mount Hardware

Choose Device Holder AND Mount

STEP 1: What DEVICE will be used?

- Smart Phone 32359 _____ \$150
- Large Smart Phone 32354 _____ \$150
- I-Pad 32355 _____ \$184
- I-Pad Mini 32356 _____ \$178
- 7" Tablet 32357 _____ \$173
- 10" Tablet 32358 _____ \$395

STEP 2: Where does it need to MOUNT?

- Frame 32350 _____ \$512
- Armrest 32351 _____ \$395
Not compatible with Swing Away Armrests
- Tray 32352 _____ \$212

Cup Holder

- Heavy Duty 30889 _____ \$44
Compatible only with Flip Back Armrests