



Phone (800) 981-1540 | Fax (715) 254-0996  
 sales@kimobility.com | www.kimobility.com

# CATALYST 4TTL

Effective 05/23/2017

Company: \_\_\_\_\_

Ship To: \_\_\_\_\_

Account Number: \_\_\_\_\_ Date: \_\_\_\_\_

Address: \_\_\_\_\_

Quote  PO Number: \_\_\_\_\_

City: \_\_\_\_\_ State: \_\_\_\_\_

RTS: \_\_\_\_\_

Zip: \_\_\_\_\_ Phone: \_\_\_\_\_

Marked For: \_\_\_\_\_

Email/Fax: \_\_\_\_\_

**Catalyst 4TTL K0004** **\$1850**  
 250 lb weight capacity

*HCPCS codes provided are based on code verification by the PDAC or our interpretation of the code definitions.*

*Proper use of HCPCS codes does not ensure coverage or payment. For coverage information, verify the appropriate payer's coverage policy.*

## Frame

### Seat Width / Seat Depth

Select box corresponding to desired Seat Width and Depth

		Seat Width									
		14"	15"	16"	17"	18"	19"	20"	21"	22"	
Seat Depth	14"	N/C	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>
	15"	N/C	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>
	16"	N/C	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>
	17"	N/C	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>
	18"	N/C	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>
	19"	N/C	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>
	20"	\$450 E2203	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>

### Heavy Duty Option

350 lb Weight Capacity \_\_\_\_\_ \$500  
 Requires minimum 6" caster height

### Cross Brace

Standard Cross Braces \_\_\_\_\_ Std  
 Seating System Cross Braces \_\_\_\_\_ N/C  
 Allows for easy attachment of seating system mounting hardware. Not available with Seat Sling.

## Backrest

### Back Post / Height Range / Back Height

Select Back Height in row corresponding to desired Back Post and Height Range

Back Post	Price	Height Range	Back Height			
Standard	N/C	Short	<input type="checkbox"/> 8.5	<input type="checkbox"/> 9.5	<input type="checkbox"/> 10.5	<input type="checkbox"/> 11.5
		Medium	<input type="checkbox"/> 12.5	<input type="checkbox"/> 13.5	<input type="checkbox"/> 14.5	<input type="checkbox"/> 15.5
		Tall	<input type="checkbox"/> 16.5	<input type="checkbox"/> 17.5	<input type="checkbox"/> 18.5	<input type="checkbox"/> 19.5
Angle Adjustable w/ Push Handle	\$185	Short	<input type="checkbox"/> 14	<input type="checkbox"/> 15	<input type="checkbox"/> 16	
		Medium	<input type="checkbox"/> 16	<input type="checkbox"/> 17	<input type="checkbox"/> 18	
		Tall	<input type="checkbox"/> 18	<input type="checkbox"/> 19	<input type="checkbox"/> 20	

### Standard Back Post

Select option below only if a Standard Back Post option is selected on chart at left

Straight with Push Handle \_\_\_\_\_ Std  
 Straight without Push Handle \_\_\_\_\_ N/C  
 8° Bend with Push Handle \_\_\_\_\_ N/C  
 Not available with Short Back Height Range

**Backrest continued**

- Back Upholstery**
- Standard \_\_\_\_\_ Std  
Not available with Short or Medium Back Height Ranges on 20" - 22" Seat Widths
  - Tension Adjustable E2611 \_\_\_\_\_ \$245
  - Omit Back Upholstery \_\_\_\_\_ - \$10

- Military Patch**
- Attached only if ordered with Tension Adjustable Back Upholstery*
- US Military Service Patch \_\_\_\_\_ \$90  
 Army  Navy  Marines  Air Force  Coast Guard
  - US Purple Heart Patch \_\_\_\_\_ \$90
  - Purple Heart with Military Service Patch Combo \_\_\_\_\_ \$145  
 Army  Navy  Marines  Air Force  Coast Guard

**Footrest**

- Hanger Bend**
- 60° \_\_\_\_\_ N/C  
Available only with Extension Tube Hangers
  - 70° \_\_\_\_\_ N/C
  - 80° \_\_\_\_\_ N/C  
Not available with 7" or 8" Casters
  - 90° \_\_\_\_\_ N/C  
Available only with Front Mount Hangers. Not available with 7" or 8" Casters.
  - Elevating Leg Rests (ELR) E0990 \_\_\_\_\_ \$295
  - Pro Elevating Leg Rests (PRO ELR) E0990 **New!** \_\_\_\_\_ \$350
  - Omit Hangers \_\_\_\_\_ N/C

- Footplate**
- Composite \_\_\_\_\_ N/C  
Not available with Front Mount Hangers or PRO ELR
  - Composite Angle Adjustable K0040 \_\_\_\_\_ \$185  
Not available with Front Mount Hangers, ELR or 14" Seat Width
  - Aluminum Angle Adjustable K0040 \_\_\_\_\_ \$185
  - One Piece Flip Up Angle Adjustable K0040 \_\_\_\_\_ \$190  
Not available with Front Mount Hangers, ELR or PRO ELR

**Footrest / Extension / Seat to Footrest Length**

Select Seat to Footrest Length in row corresponding to desired Footrest and Extension. Footrest Length of at least 2.5" shorter than Front Seat Height recommended. Extension Tube Hangers measurements are based on 70° Hanger / Front Frame Bend.

Footrest	Extension	Seat to Footrest Length													
Extension Tube Hangers (60°, 70°, 80°)	Short	<input type="checkbox"/> 11	<input type="checkbox"/> 11.5	<input type="checkbox"/> 12	<input type="checkbox"/> 12.5	<input type="checkbox"/> 13	<input type="checkbox"/> 13.5	<input type="checkbox"/> 14	<input type="checkbox"/> 14.5	<input type="checkbox"/> 15					
	Medium	<input type="checkbox"/> 13	<input type="checkbox"/> 13.5	<input type="checkbox"/> 14	<input type="checkbox"/> 14.5	<input type="checkbox"/> 15	<input type="checkbox"/> 15.5	<input type="checkbox"/> 16	<input type="checkbox"/> 16.5	<input type="checkbox"/> 17					
	Long	<input type="checkbox"/> 15.5	<input type="checkbox"/> 16	<input type="checkbox"/> 16.5	<input type="checkbox"/> 17	<input type="checkbox"/> 17.5	<input type="checkbox"/> 18	<input type="checkbox"/> 18.5	<input type="checkbox"/> 19	<input type="checkbox"/> 19.5					
Front Mount Hangers <sup>1</sup>	<b>Standard Selections</b>											<b>Requires 3" Extension for N/C</b>			
	Short (70°, 80°, 90°)	<input type="checkbox"/> 5	<input type="checkbox"/> 5.5	<input type="checkbox"/> 6	<input type="checkbox"/> 6.5	<input type="checkbox"/> 7							<input type="checkbox"/> 10	<input type="checkbox"/> 10.5	<input type="checkbox"/> 11
		<input type="checkbox"/> 7.5	<input type="checkbox"/> 8	<input type="checkbox"/> 8.5	<input type="checkbox"/> 9	<input type="checkbox"/> 9.5							<input type="checkbox"/> 11.5	<input type="checkbox"/> 12	<input type="checkbox"/> 12.5
	Medium (70°, 80°, 90°)	<input type="checkbox"/> 5	<input type="checkbox"/> 5.5	<input type="checkbox"/> 6	<input type="checkbox"/> 6.5	<input type="checkbox"/> 7	<input type="checkbox"/> 7.5	<input type="checkbox"/> 8					<input type="checkbox"/> 12	<input type="checkbox"/> 12.5	<input type="checkbox"/> 13
<input type="checkbox"/> 8.5		<input type="checkbox"/> 9	<input type="checkbox"/> 9.5	<input type="checkbox"/> 10	<input type="checkbox"/> 10.5	<input type="checkbox"/> 11	<input type="checkbox"/> 11.5					<input type="checkbox"/> 13.5	<input type="checkbox"/> 14	<input type="checkbox"/> 14.5	
Long (70°, 80°)	<input type="checkbox"/> 5	<input type="checkbox"/> 5.5	<input type="checkbox"/> 6	<input type="checkbox"/> 6.5	<input type="checkbox"/> 7	<input type="checkbox"/> 7.5	<input type="checkbox"/> 8	<input type="checkbox"/> 8.5	<input type="checkbox"/> 9			<input type="checkbox"/> 14	<input type="checkbox"/> 14.5	<input type="checkbox"/> 15	
	<input type="checkbox"/> 9.5	<input type="checkbox"/> 10	<input type="checkbox"/> 10.5	<input type="checkbox"/> 11	<input type="checkbox"/> 11.5	<input type="checkbox"/> 12	<input type="checkbox"/> 12.5	<input type="checkbox"/> 13	<input type="checkbox"/> 13.5	<input type="checkbox"/> 15.5	<input type="checkbox"/> 16	<input type="checkbox"/> 16.5			
Elevating Leg Rests (ELR)	Short	<input type="checkbox"/> 12	<input type="checkbox"/> 13	<input type="checkbox"/> 14	<input type="checkbox"/> 15	<input type="checkbox"/> 16									
	Long	<input type="checkbox"/> 13	<input type="checkbox"/> 14	<input type="checkbox"/> 15	<input type="checkbox"/> 16	<input type="checkbox"/> 17	<input type="checkbox"/> 18								
Extended Length ELR	Short	<input type="checkbox"/> 18	<input type="checkbox"/> 19	<input type="checkbox"/> 20	<input type="checkbox"/> 21	<input type="checkbox"/> 22									
	Long	<input type="checkbox"/> 18	<input type="checkbox"/> 19	<input type="checkbox"/> 20	<input type="checkbox"/> 21	<input type="checkbox"/> 22	<input type="checkbox"/> 23	<input type="checkbox"/> 24							
PRO ELR Low Pivot	Short	<input type="checkbox"/> 12 - 17.5													
	Long	<input type="checkbox"/> 15 - 21													
PRO ELR High Pivot	Short	<input type="checkbox"/> 10.5 - 15.5													
	Long	<input type="checkbox"/> 13 - 19													

<sup>1</sup> Available only with Aluminum Angle Adjustable Footplates. See Minimum Front Seat Height chart below.

**Minimum Front Seat Height with Front Mount Hangers**

Hanger Length	Short		Medium		Long	
	70°	80/90°	70°	80/90°	70°	80°
Minimum Front Seat Height with Standard Selection	13.5"	13.5"	14"	15"	16"	17"
Minimum Front Seat Height with 3" Extension	14.5"	16"	16.5"	18"	18.5"	20"

- 3" Extension**
- Available only with Front Mount Hangers. Provided as parts.*
- 3" Extension for Front Mount Hanger \_\_\_\_\_ N/C

**Footrest continued**

- Heel Loop**
- Heel Loops E0951 \_\_\_\_\_ \$30  
*Not available with One Piece Flip Up Angle Adjustable Footplate*
  - Omit Heel Loops \_\_\_\_\_ N/C
- Calf Strap**
- Velcro Adjustable K0038 \_\_\_\_\_ \$35
  - Padded Velcro Adjustable K0038 \_\_\_\_\_ \$68
  - Bodypoint Padded AeroMesh K0038 \_\_\_\_\_ \$102

- Residual Limb Support**
- Residual Limb Support E1020
    - Left \_\_\_\_\_ \$395
    - Right \_\_\_\_\_ \$395
    - Omit Left Hanger \_\_\_\_\_ N/C
    - Omit Right Hanger \_\_\_\_\_ N/C
- Cover**
- Smooth  Air Mesh
  - Smooth  Air Mesh
- Pad Size**
- 5 x 8  9 x 14
  - 7 x 10  11 x 10
  - 9 x 10  11 x 14

**Front Seat Height and Casters**

**Caster Size / Front Frame Type / Front Seat Height**

Select Front Seat Height in row corresponding to desired caster size.  
Note: Selections in first three Front Seat Height columns use Super Low Front Frame.

Caster Size	Front Seat Height													
	Uses Super Low Front Frame													
4"	<input type="checkbox"/> 13.5 <sup>1</sup>	<input type="checkbox"/> 14	<input type="checkbox"/> 14.5	<input type="checkbox"/> 15 <sup>2</sup>	<input type="checkbox"/> 15.5	<input type="checkbox"/> 16	<input type="checkbox"/> 16.5	<input type="checkbox"/> 17	<input type="checkbox"/> 17.5	<input type="checkbox"/> 18	<input type="checkbox"/> 18.5	<input type="checkbox"/> 19	<input type="checkbox"/> 19.5	<input type="checkbox"/> 20
5"	<input type="checkbox"/> 14.5	<input type="checkbox"/> 15	<input type="checkbox"/> 15.5	<input type="checkbox"/> 16	<input type="checkbox"/> 16.5	<input type="checkbox"/> 17	<input type="checkbox"/> 17.5	<input type="checkbox"/> 18	<input type="checkbox"/> 18.5	<input type="checkbox"/> 19	<input type="checkbox"/> 19.5	<input type="checkbox"/> 20	<input type="checkbox"/> 20.5	
6"	<input type="checkbox"/> 15.5	<input type="checkbox"/> 16	<input type="checkbox"/> 16.5	<input type="checkbox"/> 17	<input type="checkbox"/> 17.5	<input type="checkbox"/> 18	<input type="checkbox"/> 18.5	<input type="checkbox"/> 19	<input type="checkbox"/> 19.5	<input type="checkbox"/> 20	<input type="checkbox"/> 20.5	<input type="checkbox"/> 21		
7"	<input type="checkbox"/> 17	<input type="checkbox"/> 17.5		<input type="checkbox"/> 18	<input type="checkbox"/> 18.5	<input type="checkbox"/> 19	<input type="checkbox"/> 19.5	<input type="checkbox"/> 20	<input type="checkbox"/> 20.5	<input type="checkbox"/> 21				
8"	<input type="checkbox"/> 17.5	<input type="checkbox"/> 18	<input type="checkbox"/> 18.5	<input type="checkbox"/> 19	<input type="checkbox"/> 19.5	<input type="checkbox"/> 20	<input type="checkbox"/> 20.5	<input type="checkbox"/> 21						

<sup>1</sup> Front Seat Height not available with 4 x 1" Poly Caster  
<sup>2</sup> Super Low Front Frame required when using 4 x 1" Poly Caster

- Caster Type**
- 1" Poly \_\_\_\_\_ N/C
  - 1" Poly Aluminum \_\_\_\_\_ \$110  
*Not available with 6", 7" and 8" Caster Sizes*
  - 1" Pneumatic \_\_\_\_\_ \$80  
*Not available with 4", 5" and 7" Caster Sizes*
  - 1.5" Poly \_\_\_\_\_ N/C  
*Not available with 7" and 8" Caster Sizes*
  - 1.5" Soft Roll Aluminum \_\_\_\_\_ \$50  
*Not available with 7" and 8" Caster Sizes*
  - 2" Poly \_\_\_\_\_ N/C  
*Not available with 4", 5", 7" and 8" Caster Sizes*
  - 2" Pneumatic w/ Foam Insert \_\_\_\_\_ \$125  
*Not available with 4", 5", 6" and 7" Caster Sizes*

Maximum Seat Dump Reference							
(Front Seat Height minus Rear Seat Height)							
Seat Depth	14"	15"	16"	17"	18"	19"	20"
Maximum Dump	1.5"	2"	2"	2"	2.5"	2.5"	2.5"

- Fork Stem**
- Best Fit \_\_\_\_\_ Std  
*Ki Mobility will pick optimal Fork Stem based on Caster Size and Front Seat Height*
  - Standard Fork Stem \_\_\_\_\_ N/C
  - 3/4" Longer Fork Stem \_\_\_\_\_ N/C
  - 1 1/2" Longer Fork Stem \_\_\_\_\_ N/C

- Fork Options**
- Frog Legs Ultra-Sport Fork - Black \_\_\_\_\_ \$350  
*Available only with 4", 5" or 6" Casters. Not available with Caster Pin Locks, Heavy Duty Option, 2" Poly caster or 2" Pneumatic w/ Foam Insert caster.*
- End User Weight Required:** \_\_\_\_\_ lbs
- Caster Pin Locks K0073 \_\_\_\_\_ \$195  
*Not available with 1 1/2" Longer Fork Stem or Frog Legs Fork*

**Wheel Locks**

- Wheel Lock**
- Push to Lock \_\_\_\_\_ Std
  - Push to Lock Flush Mount \_\_\_\_\_ N/C
  - Pull to Lock \_\_\_\_\_ N/C
  - Omit Wheel Locks \_\_\_\_\_ N/C

- Wheel Lock Options**
- Extension Handles E0961 \_\_\_\_\_ \$67  
*Available only with Push to Lock or Pull to Lock Wheel Locks*
  - Grade Aids E0974 \_\_\_\_\_ \$150  
*Available only with Push to Lock Wheel Locks and Pneumatic or Pneumatic with Airless Insert Rear Tire Types*

## Rear Seat Height and Wheels

### Wheel Size / Tire Profile / Rear Seat Height

Select Rear Seat Height in row corresponding to desired Wheel Size and Tire Profile. 14" Seat Depth not available with 24", 25" or 26" Rear Wheels.

Wheel Size	Tire Profile	Rear Seat Height				
20"	Low	<input type="checkbox"/> 12.5 <sup>1</sup>	<input type="checkbox"/> 13.5	<input type="checkbox"/> 14.5	<input type="checkbox"/> 15.5	<input type="checkbox"/> 16.5
	Full	<input type="checkbox"/> 13 <sup>1</sup>	<input type="checkbox"/> 14	<input type="checkbox"/> 15	<input type="checkbox"/> 16	<input type="checkbox"/> 17
22"	Low	<input type="checkbox"/> 13.5 <sup>1</sup>	<input type="checkbox"/> 14.5	<input type="checkbox"/> 15.5	<input type="checkbox"/> 16.5	<input type="checkbox"/> 17.5
	Full	<input type="checkbox"/> 14 <sup>1</sup>	<input type="checkbox"/> 15	<input type="checkbox"/> 16	<input type="checkbox"/> 17	<input type="checkbox"/> 18
24"	Low	<input type="checkbox"/> 14.5 <sup>1</sup>	<input type="checkbox"/> 15.5	<input type="checkbox"/> 16.5	<input type="checkbox"/> 17.5	<input type="checkbox"/> 18.5
	Full	<input type="checkbox"/> 15 <sup>1</sup>	<input type="checkbox"/> 16	<input type="checkbox"/> 17	<input type="checkbox"/> 18	<input type="checkbox"/> 19
25"	Low	<input type="checkbox"/> 15 <sup>1</sup>	<input type="checkbox"/> 16	<input type="checkbox"/> 17	<input type="checkbox"/> 18	<input type="checkbox"/> 19
26"	Low	<input type="checkbox"/> 15.5 <sup>1</sup>	<input type="checkbox"/> 16.5	<input type="checkbox"/> 17.5	<input type="checkbox"/> 18.5	<input type="checkbox"/> 19.5
	Full	<input type="checkbox"/> 16 <sup>1</sup>	<input type="checkbox"/> 17	<input type="checkbox"/> 18	<input type="checkbox"/> 19	<input type="checkbox"/> 20

<sup>1</sup> Not available with Flip Back Armrests

### Rear Wheel Type

- Assembled 24" Six-Spoke Mag with Treaded Tire and Composite Handrim \_\_\_\_\_ - \$25  
*No substitutions available. Do not select Rear Tire Type or Handrim. Measures as Full Tire.*
  - Spoke \_\_\_\_\_ N/C
  - Mag \_\_\_\_\_ N/C  
*Not available in 25" or 26" rear wheels*
  - Superlight Spoke \_\_\_\_\_ \$265  
*Available only in 24" and 25" rear wheels. Not available with Heavy Duty Option.*
  - Spinergy Spox \_\_\_\_\_ \$795  
*Not available in 20" rear wheel*
- Color:  Black  Red  White  Blue  Yellow
- Handrim:  Short Tab  Long Tab  
*Available only with Aluminum Anodized and Plastic Coated Handrims*

### Rear Tire Type

#### Full Profile

- Street Pneumatic 1 3/8" E2211, E2212 \_\_\_\_\_ Std  
*Not available on 25" rear wheel*
- Full Poly 1 3/8" E2218 \_\_\_\_\_ N/C  
*Not available on Superlight or 25" rear wheel*
- Pneumatic with Airless Insert 1 3/8" E2211, E2213 \_\_\_\_\_ \$95  
*Not available on Superlight or 25" rear wheel*
- Schwalbe HS-127 1 3/8" E2211, E2212 **New!** \_\_\_\_\_ \$125  
*Available only on 24" rear wheel. Not available on Mag wheels.*
- Euro Trek Knobby 1 3/8" E2211, E2212 \_\_\_\_\_ \$95  
*Available only on 24" rear wheel*

#### Low Profile

- Low Poly 1" E2218 \_\_\_\_\_ N/C  
*Not available on Superlight, 25" or 26" rear wheels*
- Sentinel Solid Tire 1" E2218 \_\_\_\_\_ \$95  
*Not available on Superlight wheels*
- High Pressure 1" E2211, E2212 \_\_\_\_\_ \$95  
*Not available on Mag wheels*

### Handrim

#### Standard Handrims

- Aluminum Anodized \_\_\_\_\_ Std
- Superlight Aluminum Anodized \_\_\_\_\_ \$50  
*Available only on 24" and 25" Superlight and Spinergy wheels. Not available with Heavy Duty Option.*
- Plastic Coated \_\_\_\_\_ \$50
- Projection (8 Vertical Projections) \_\_\_\_\_ \$295  
*Available on 22", 24", 25" and 26" Spoke wheels & 24" Mag wheels w/full tires*
- Omit Handrim \_\_\_\_\_ N/C

#### Natural Fit Handrims

- Not available on 20" rear wheel or 22" Mag wheels. Not recommended with Mag wheels.**
- Natural Fit LT Regular Grip \_\_\_\_\_ \$445
  - Natural Fit Standard Regular Grip \_\_\_\_\_ \$445
  - Natural Fit LT Super Grip \_\_\_\_\_ \$495
  - Natural Fit Standard Super Grip \_\_\_\_\_ \$495
  - Natural Fit LT Omit Grip \_\_\_\_\_ \$345
  - Natural Fit Standard Omit Grip \_\_\_\_\_ \$345

#### Handrim Option

- 3M Micro Grip \_\_\_\_\_ \$65  
*Not available on Plastic Coated or Projection Handrims*

## Axles

### Axle Bracket

- Offset Axle Bracket \_\_\_\_\_ Std
- Reverse Mount Offset Axle Bracket \_\_\_\_\_ N/C  
*Not available with Swing Away Armrests*

### Axle

- Quick Release \_\_\_\_\_ Std
- Quad Quick Release \_\_\_\_\_ \$195  
*Not available with Mag wheels*

## Armrests

### Armrest Please select Arm Pad option

- Height Adjustable T-Arm (8.5" - 13.5") E0973 \$265
- Height Adjustable T-Arm Low (7.5" - 11.5") E0973 \$265
- Height Adjustable Flip Back (9.5" - 13.5") E0973 \$275  
*Not available with Seat Depths less than 16" or Angle Adjustable Back Posts*
- Flip Back (10.5") N/C  
*Not available with Seat Depths less than 16" or Angle Adjustable Back Posts*

### Arm Pad

**Not available with Swing Away Armrests**

- Standard Arm Pad N/C  Desk  Full
- Waterfall Arm Pad \$25  Desk  Full  
*Available only with Height Adjustable T-Arms*

### Armrest continued

- Swing Away (10.5" - 12.5") E0973 \$195  
*Not available with Angle Adjustable Back Posts or Arm Pads*
- Omit Armrests N/C

## Accessories

### Removable Side Guard

- Composite Side Guards K0108 \$140  
*Not available with T-Arm or Flip Back Armrests*

### Anti-Tip

- Rear Anti-Tips E0971 \$140

### Positioning Belt

- 1 1/2" Auto Buckle E0978 \$70
- 1 1/2" Padded Auto Buckle E0978 \$155
- 2" Airline Buckle E0978 \$110
- Bodypoint Evoflex - Pelvic Stabilizer E0978 \$183  
Size:  Small (16-24")  Medium (19-29")  Large (25-37")
- Bodypoint Belt Mounting Kit \$105  
Style:  Bracket  Band Clamp 1"

### O<sub>2</sub> Holder

- O<sub>2</sub> Holder E2208 \$295  
 Right Side  Left Side

## Frame Color

### Gloss Finish

- Black N/C
- Blue Skies N/C
- Bubblegum Pink N/C
- Candy Blue N/C
- Candy Purple N/C
- Candy Red N/C
- Copper N/C
- Fuchsia N/C
- Granny Smith N/C
- Yellow N/C

### Matte Finish

- Army Green N/C
- Charcoal Gray N/C
- Emerald Green N/C
- Iridescent Burgundy N/C
- Iridescent Navy N/C
- Retro Red N/C
- Royal Blue N/C
- Shadow Black N/C
- Smooth White N/C

### Graphic Image Transfer

Visit [www.kimobility.com](http://www.kimobility.com) for current graphic image transfer options  
 Note selection below

- \_\_\_\_\_ \$595

Catalyst Front Seat Height Matrix																																																			
4" Caster				5" Caster				6" Caster				7" Caster				8" Caster																																			
FSH	Super Low	Stem	Fork	FSH	Super Low	Stem	Fork	FSH	Super Low	Stem	Fork	FSH	Super Low	Stem	Fork	FSH	Super Low	Stem	Fork																																
15 <sup>1</sup>	13.5 <sup>1</sup>	Std	4"	16	14.5	Std	4"	17	15.5	Std	5"	18	17	Std	6"	19	18	Std	6"																																
15.5	14			16.5	15.5	3/4		18	16.5	3/4		18.5	17.5			20	18.5	3/4																																	
16 <sup>2</sup>	14.5 <sup>2</sup>	17.5		16	1 1/2	18.5		17	1 1/2	19		17.5	20.5	n/a		1 1/2																																			
16.5	15	3/4		Std	16	15	Std	5"	17	16	Std	6"	19.5	18		3/4	7"	19	17.5	Std	7"																														
16.5 <sup>2</sup>	15 <sup>2</sup>	1 1/2			16.5	15.5			17.5	16.5			Std	19.5				18.5	1 1/2			19.5	18																												
17	15.5	1 1/2		Std	17	15.5	3/4		5"	18	17		Std	6"		20		19	1 1/2	7"		20	18.5	Std	7"																										
15.5	14.5				17.5	16				18	16.5					18.5		17.5				20.5	19																												
16	15	3/4		5"	17.5	16.5	1 1/2		6"	18.5	17		3/4	7"		19		18	Std	8"		19.5	18.5	3/4	7"																										
16.5	15				18	17				19	17.5					20		19																																	
17	15.5	1 1/2		Std	16.5	15.5	Std		5"	18.5	17.5		1 1/2	6"		20		19	3/4	7"		20.5	n/a	3/4	7"																										
17	16		17		16	19				18	21				n/a																																				
17.5	16.5	1 1/2	Std	17.5	16.5	3/4	6"		19.5	18.5	1 1/2		7"	20.5	19	1 1/2		8"	20.5	19		1 1/2	7"																												
17	15.5			18	16.5				20	19				21	n/a																																				
17.5	16.5	3/4	6"	18.5	17.5	1 1/2	7"	19.5	18.5	3/4	8"	8" x 2" caster not available with 6" fork																																							
18	17			18	17			19.5	18.5											20.5	19																														
18.5	17	1 1/2	Std	18.5	17.5	1 1/2	6"	18.5	17.5	3/4	7"									8" x 2" caster not available with 6" fork																															
19	17.5			19	18			19	18																			20.5	19																						
18	17	Std	7"	18	16.5	Std	6"	19.5	18.5	3/4	7"																	8" x 2" caster not available with 6" fork																							
18.5	17.5			18.5	17			20	19																											20.5	18.5														
19	17.5	3/4	7"	19	17.5	3/4	7"	19.5	18	1 1/2	8"																									8" x 2" caster not available with 6" fork															
19.5	18			18.5	17.5			20	18.5																																			20.5	19						
19.5	18.5	1 1/2	Std	19	18	3/4	7"	20.5	19	1 1/2	8"																																	8" x 2" caster not available with 6" fork							
20	19			19.5	18.5			21	n/a																																										
				19.5	18	1 1/2	7"	6" x 2" caster not available with 5" fork				8" x 2" caster not available with 6" fork																																							
				20	18.5																																														
				20.5	19																																														

<sup>1</sup> Front Seat Height not available with 4 x 1" Poly Caster

<sup>2</sup> Not available with 4 x 1" Poly, 4 x 1.5" Poly or 4 x 1.5" Soft Roll Aluminum

Seat heights with "FSH" heading available with Standard Front Frame. Seat heights with "Super Low" heading available with Super Low Front Frame. Actual seat height may vary depending upon configuration.

## Cushions

### General Use Cushions

#### STEP 1: Choose Cushion

- 2" General Use Foam Cushion (14"x14" to 20"x20") E2601 \_\_\_\_ \$99
- 2" General Use Foam Cushion (Over 20"x20") E2601 \_\_\_\_ \$250

#### STEP 2: Choose Cushion Size

- Match to Chair Size
- Other Size \_\_\_\_\_

### Matrix Seat Cushions

*Cushion size is determined by chair dimensions. Certain sizes will have higher price and longer lead time; visit [www.motionconcepts.com/matrix\\_international.html](http://www.motionconcepts.com/matrix_international.html) for size specifications.*

#### STEP 1: Choose Cushion

- PS Cushion \_\_\_\_\_ \$414  
*Cushions 14" wide have a weight capacity of 150 lbs. Cushions 15" through 20" have a weight capacity of 300 lbs.*
- Vi Cushion \_\_\_\_\_ \$510  
*Weight capacity of 300 lbs*
- Libra Cushion \_\_\_\_\_ \$582  
*Weight capacity of 300 lbs*

#### STEP 2: Choose Cushion Size

- Match to Chair Size
- Other Size - Provide Part Number \_\_\_\_\_

### Other Matrix Products

*Must provide part number*

- \_\_\_\_\_

### ROHO Seat Cushions

*Cushion size is determined by chair dimensions. Certain sizes will have higher price and longer lead time; visit [therohogroup.com](http://therohogroup.com) for size specifications.*

#### STEP 1: Choose Cushion

##### Standard Cushions

*For chairs under 22" wide*

- HIGH PROFILE® Single Compartment (E2622) \_\_\_\_\_ \$603
- MID PROFILE™ Single Compartment (E2622) \_\_\_\_\_ \$603
- LOW PROFILE® Single Compartment (E2622) \_\_\_\_\_ \$603
- QUADTRO SELECT® High Profile® (E2624) \_\_\_\_\_ \$603
- QUADTRO SELECT® Mid Profile® (E2624) \_\_\_\_\_ \$603
- HYBRID ELITE® Single Comp Smart Check Ready (E2622) \_\_\_\_\_ \$623
- CONTOUR SELECT® (E2624) \_\_\_\_\_ \$603

##### X Large Cushions

*For 22" wide chairs*

- HIGH PROFILE® Single Compartment X Large (E2623) \_\_\_\_\_ \$670
- LOW PROFILE® Single Compartment X Large (E2623) \_\_\_\_\_ \$670
- QUADTRO SELECT® High Profile® X Large (E2625) \_\_\_\_\_ \$670

##### Cushions with Heavy Duty Cover

- HIGH PROFILE® Single Compartment w/ HD Cover (E2622) \_\_\_\_\_ \$629
- QUADTRO SELECT® High Profile® w/ HD Cover (E2624) \_\_\_\_\_ \$629

#### STEP 2: Choose Cushion Size

- Match to Chair Size
- Other Size - Provide Part Number \_\_\_\_\_

**Backs**

**Matrix Backs**

- Elite Back \_\_\_\_\_ \$750  
*Not available with 22" Seat Width. Back Height of 14" not available with 19", 20" and 21" Seat Widths.*  
**Back Height:**  
 10"  12"  14"  16"  18"  20"
- Elite Deep Back \_\_\_\_\_ \$912  
*Not available with 22" Seat Width. Back Height of 10" not available with 14" Seat Width.*  
**Back Height:**  
 10"  12"  16"  18"  20"
- Mini Back \_\_\_\_\_ \$750  
*Available only with 14" and 15" Seat Widths. Comes standard with Kidabra cover.*  
**Back Height:**  
 7"  10"  12"  14"

**Therafin Headrest for MATRX Back Systems**

- Headrest, Multi-Axis Offset Hardware with Smooth/Soft Pad (E1028, E0955) 34187 \_\_\_\_\_ \$420  
 PAD SIZE:  5 x 8  5 x 10  5 x 14

**Other Matrix Mounting Hardware**

- Must provide part number**
- \_\_\_\_\_  
*For mounting hardware other than quick release*

**Other Matrix Products**

- Must provide part number**
- \_\_\_\_\_

**ROHO Back Systems (Available for chairs up to 20" wide)**

**STEP 1: Choose Back**

- AGILITY™ Min. Contour w/ Foam Only (E2611) \_\_\_\_\_ \$549  
 BACKREST HEIGHT:  10"  13"  16"  18"
- AGILITY™ Min. Contour w/ Foam & Lumbar Pad (E2613) \_\_\_\_\_ \$659  
 BACKREST HEIGHT:  10"  13"  16"  18"
- AGILITY™ Min. Contour w/ Air Only (E2620) \_\_\_\_\_ \$769  
 BACKREST HEIGHT:  10"  13"  16"  18"
- AGILITY™ Min. Contour w/ Air & Lumbar Pad (E2615) \_\_\_\_\_ \$769  
 BACKREST HEIGHT:  10"  13"  16"  18"
- AGILITY™ Mid Contour w/ Air Only (E2620) \_\_\_\_\_ \$781  
 BACKREST HEIGHT:  10"  13"  16"  18"  20"
- AGILITY™ Mid Contour w/ Air & Lumbar Pad (E2615) \_\_\_\_\_ \$781  
 BACKREST HEIGHT:  10"  13"  16"  18"  20"
- AGILITY™ Max Contour w/ Air Only (E2620) \_\_\_\_\_ \$920  
 BACKREST HEIGHT:  16"  18"  20"
- AGILITY™ Max Contour w/ Air & Lumbar Pad (E2615) \_\_\_\_\_ \$920  
 BACKREST HEIGHT:  16"  18"  20"

**STEP 2: Choose Back Hardware**

- Fixed \_\_\_\_\_ N/C  
*Not compatible with 18" and 20" tall backs*
- Quick Release \_\_\_\_\_ N/C

**ROHO Laterals for AGILITY™ Back Systems**

**Sold in pairs**

- Swing Away with Curved Pad (E1028/E0956) \_\_\_\_\_ \$1014  
 PAD SIZE:  Medium 4 T x 4 L  Large 4 T x 5 L  X-Large 5 T x 5 L  
 OFFSET TYPE:  Flush  Offset  1" Offset  2" Offset
- Fixed with Flat Pad (E0956) \_\_\_\_\_ \$452  
 PAD SIZE:  Medium 4 T x 4 L  Large 4 T x 5 L  X-Large 5 T x 5 L  
 OFFSET TYPE:  Flush  Offset

**Therafin Headrest for AGILITY Back Systems**

- Headrest, Multi-Axis Offset Hardware with Smooth/Soft Pad (E1028, E0955) 34187 \_\_\_\_\_ \$420  
 PAD SIZE:  5 x 8  5 x 10  5 x 14





**Vest**

*Not compatible with Back Upholstery*

- Full Shape, Comfort Strap (E0960) \_\_\_\_\_ \$167
  - Small 30589     Medium 30590     Large 30591
  - X-Large 30592     2X-Large 30593
- Full Shape, Extended Strap (E0960) \_\_\_\_\_ \$167
  - Small 30594     Medium 30595     Large 30596
  - X-Large 30597     2X-Large 30598
- Trim Shape, Comfort Strap (E0960) \_\_\_\_\_ \$169
  - Small 30618     Medium 30619     Large 30620
  - X-Large 30621     2X-Large 30622
- Trim Shape, Extended Strap (E0960) \_\_\_\_\_ \$169
  - Small 30604     Medium 30624     Large 30625
  - X-Large 30626     2X-Large 30627

**Shoulder Harness**

*With metal clips and cam locks. Not compatible with Back Upholstery.*

- Small (E0960) 61403 \_\_\_\_\_ \$157
- Medium (E0960) 61503 \_\_\_\_\_ \$164
- Large (E0960) 61703 \_\_\_\_\_ \$170
- X-Large (E0960) 61803 \_\_\_\_\_ \$177

**Hip Belt**

- 4 Point Y-Style, Metal Buckle (E0978) \_\_\_\_\_ \$173
  - Small 60024     Medium 60025     Large 60026
- Dual Rear Pull, Metal Attaching Clips (E0978)
  - 1" Straps, 6.25"W x 1.75"H Pads 60161 \_\_\_\_\_ \$97
  - 1.5" Straps, 7.25"W x 2.25"H Pads 60171 \_\_\_\_\_ \$106
  - 1.5" Straps, 9.25"W x 2.5"H Pads 60181 \_\_\_\_\_ \$106
- Single Pull, Metal Buckle (E0978)
  - 1" Straps, 4.25"W x 1.75"H Pads 61001 \_\_\_\_\_ \$86
  - 1" Straps, 5.25"W x 1.75"H Pads 61101 \_\_\_\_\_ \$86
  - 1" Straps, 6.25"W x 1.75"H Pads 61201 \_\_\_\_\_ \$86
  - 1.5" Straps, 5.25"W x 2.25"H Pads 62001 \_\_\_\_\_ \$95
  - 1.5" Straps, 9.25"W x 2.5"H Pads 62201 \_\_\_\_\_ \$95
  - 2" Straps, 9.25"W x 3"H Pads 64001 \_\_\_\_\_ \$102
  - 2" Straps, 11.25"W x 3"H Pads 64101 \_\_\_\_\_ \$102

**Cane Mounted Swing Away Lateral Thoracic Support**

*Sold as pair. Not compatible with Back Upholstery.*

- Swing Away Kits w/ 4"x6" Gel Pads (E0956, E1028) 50021 \_\_\_\_\_ \$478
- Swing Away Kits, No Pads (E0956, E1028) 50022 \_\_\_\_\_ \$309

**Knee Adductor**

- Flip Down Hardware (E1028) \_\_\_\_\_ \$634  
*Sold in pairs, one right side and one left side*
  - Short Arm, Fits 3/4", 7/8" and 1" Tubing 30248
  - Long Arm, Fits 3/4", 7/8" and 1" Tubing 30284
- Knee Adductor Pad with Neoprene Cover (E0956) \_\_\_\_\_ \$117  
*Sold single*
  - Small Pad, 4.25" x 3.25" 50490    Quantity: \_\_\_\_\_
  - Medium Pad, 6.25" x 4.25" 50491    Quantity: \_\_\_\_\_

**Foot / Shoe Holder**

*Compatible only with Aluminum Angle Adjustable and Composite Angle Adjustable Footplates*

- FootSure Ankle Support (K0108) \_\_\_\_\_ \$178  
*Sold in pairs*
  - X-Small w/ Buckle Closure, 6"-8" Circum. 60250
  - Small w/ Buckle Closure, 7"-8.5" Circum. 60251
  - Medium w/ Buckle Closure, 8"-11" Circum. 60252
  - Large w/ Buckle Closure, 9"-12" Circum. 60253
  - X-Large w/ Buckle Closure, 12.5"-15" Circum. 60254
  - X-Small w/ Hook-N-Loop Closure, 5"-6" Circum. 60255
  - Small w/ Hook-N-Loop Closure, 6"-8" Circum. 60256
  - Med w/ Hook-N-Loop Closure, 8"-10.5" Circum. 60257
  - Large w/ Hook-N-Loop Closure, 9"-12" Circum. 60258
  - X-Lrg w/ Hook-N-Loop Closure, 12"-15" Circum. 60259
- Padded Footbox (K0108)
  - Small, Full Design, Fits 14" Chair 30111 \_\_\_\_\_ \$287
  - Medium, Full Design, Fits 16" Chair 30112 \_\_\_\_\_ \$293
  - Large, Full Design, Fits 18" Chair 30113 \_\_\_\_\_ \$298
  - Small, Split Design, Fits 14" Chair 30152 \_\_\_\_\_ \$367
  - Medium, Split Design, Fits 16" Chair 30153 \_\_\_\_\_ \$373
  - Large, Split Design, Fits 18" Chair 30154 \_\_\_\_\_ \$378
- Molded Shoe Holders (K0108) \_\_\_\_\_ \$162  
*Sold in pairs*
  - Small with Plain Straps, 3.5" x 5.25" 30293
  - Medium with Plain Straps, 4" x 6.5" 30303
  - Large with Plain Straps, 4.75" x 8" 30309
  - X-Large with Plain Straps, 5" x 10" 30316



**Half Tray**

*Not compatible with Swing Away Armrests*

- Standard Shape, Clear, Flip-Away (E0950, E1028) \_\_\_\_\_ \$246
  - Right Side 41310       Left Side 41311
- Standard Shape, Padded, Flip-Away (E0950, E1028) \_\_\_\_\_ \$269
  - Right Side 41350       Left Side 41353
- Trimline Shape, Padded, Flip-Away (E0950, E1028) \_\_\_\_\_ \$269
  - Right Side 41352       Left Side 41351

**Clear Tray with Border and Pre-Drilled Holes**

**STEP 1: Choose Tray**

Chair Width	Tray Dimension (W x D)	Body Cutout (W x D)	HCPCS	Part Number	Price	
14"	20 x 18	10.5 x 6	E0950	33210	\$220	<input type="checkbox"/>
	20 x 18	11.5 x 6.5	E0950	33211	\$220	<input type="checkbox"/>
15"	21 x 18	11.5 x 6.5	E0950	33212	\$220	<input type="checkbox"/>
	21 x 18	12.5 x 7	E0950	33213	\$220	<input type="checkbox"/>
16"	22 x 21	12.5 x 7.5	E0950	33214	\$255	<input type="checkbox"/>
	22 x 21	13.5 x 8.5	E0950	33215	\$255	<input type="checkbox"/>
17"	23.5 x 22	13.5 x 8.5	E0950	33216	\$255	<input type="checkbox"/>
	23.5 x 22	14.5 x 9.5	E0950	33217	\$285	<input type="checkbox"/>
18"	24.5 x 22	14.5 x 9.5	E0950	33218	\$285	<input type="checkbox"/>
	24.5 x 22	15.5 x 10.5	E0950	33219	\$285	<input type="checkbox"/>
19"	26.5 x 22	15.5 x 10.5	E0950	33220	\$285	<input type="checkbox"/>
	26.5 x 22	16.5 x 11	E0950	33221	\$285	<input type="checkbox"/>
20"	26.5 x 22	17 x 11	E0950	33222	\$309	<input type="checkbox"/>
	26.5 x 22	18 x 11	E0950	33223	\$309	<input type="checkbox"/>
21"	27.5 x 23.5	17.5 x 11	E0950	33224	\$309	<input type="checkbox"/>
	27.5 x 23.5	19 x 11	E0950	33225	\$309	<input type="checkbox"/>
22"	28.5 x 23.5	18.5 x 11	E0950	33226	\$309	<input type="checkbox"/>
	28.5 x 23.5	20 x 11	E0950	33227	\$309	<input type="checkbox"/>

**STEP 2: Choose Tray Hardware**

- E-Z Lock Clamps, Tube Arm (E1028) 31072 \_\_\_\_\_ \$150  
*Compatible only with Swing Away Armrests. Sold in Pairs.*
- E-Z Lock Clamps, Standard (E1028) 31098 \_\_\_\_\_ \$142  
*Not compatible with Swing Away Armrests. Sold in Pairs.*

**Communication Mount Hardware**

*Choose Device Holder AND Mount*

**STEP 1: What DEVICE will be used?**

- Smart Phone 32359 \_\_\_\_\_ \$150
- Large Smart Phone 32354 \_\_\_\_\_ \$150
- I-Pad 32355 \_\_\_\_\_ \$184
- I-Pad Mini 32356 \_\_\_\_\_ \$178
- 7" Tablet 32357 \_\_\_\_\_ \$173
- 10" Tablet 32358 \_\_\_\_\_ \$395

**STEP 2: Where does it need to MOUNT?**

- Frame 32350 \_\_\_\_\_ \$512
- Armrest 32351 \_\_\_\_\_ \$395  
*Not compatible with Swing Away Armrests*
- Tray 32352 \_\_\_\_\_ \$212

**Cup Holder**

- Heavy Duty 30889 \_\_\_\_\_ \$44  
*Compatible only with Flip Back Armrests*