

This form must be opened in Adobe Acrobat. The submit button will not work from most browser windows.



Phone (800) 981-1540 | Fax (715) 254-0996
 sales@kimobility.com | www.kimobility.com

Effective 03/13/2017

Company: _____

Ship To: _____

Account Number: _____ Date: _____

Address: _____

Quote PO Number: _____

City: _____ State: _____

RTS: _____

Zip: _____ Phone: _____

Marked For: _____

Email/Fax: _____

Catalyst 5Ti K0005 **\$3595**
 300 lb weight capacity

HCPCS codes provided are based on code verification by the PDAC or our interpretation of the code definitions.
 Proper use of HCPCS codes does not ensure coverage or payment. For coverage information, verify the appropriate payer's coverage policy.

Frame

Seat Width / Seat Depth

Select box corresponding to desired Seat Width and Depth

		Seat Width									
		14"	15"	16"	17"	18"	19"	20"	21"	22"	
Seat Depth	Price HCPCS	N/C	N/C	N/C	N/C	N/C	\$195	\$395 E2201	\$495 E2201	\$495 E2201	
	14"	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	
	15"	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	
	16"	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	
	17"	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	
	18"	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	
	19"	\$195	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	
	20"	\$450 E2203	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	

Heavy Duty Option

350 lb Weight Capacity _____ \$500
 Requires minimum 6" caster height

Transit Option

Transit System K0108 _____ \$320
 Not available with Removable Seat Pan or Heavy Duty Option

Cross Brace

Aluminum _____ Std
 Titanium _____ \$495
 Not available with 21" and 22" Seat Widths

New! Wheelchair Drive Ready

Chair will ship from Ki Mobility with adaption kit installed for selected Alber system. Order system separately from Alber USA. Select Rear Seat Height (full profile selections only) in drive wheel size on Rear Seat Height chart (page 5). Requires Rear Anti-Tips. Not available with Spinergy wheels.

e-motion Ready _____ N/C
 twion Ready _____ N/C

Backrest

Back Post / Height Range / Back Height

Select Back Height in row corresponding to desired Back Post and Height Range

Back Post	Price	Height Range	Back Height			
Standard	N/C	Short	<input type="checkbox"/> 8.5	<input type="checkbox"/> 9.5	<input type="checkbox"/> 10.5	<input type="checkbox"/> 11.5
		Medium	<input type="checkbox"/> 12.5	<input type="checkbox"/> 13.5	<input type="checkbox"/> 14.5	<input type="checkbox"/> 15.5
		Tall	<input type="checkbox"/> 16.5	<input type="checkbox"/> 17.5	<input type="checkbox"/> 18.5	<input type="checkbox"/> 19.5
Angle Adjustable w/ Push Handle	\$185	Short	<input type="checkbox"/> 14	<input type="checkbox"/> 15	<input type="checkbox"/> 16	
		Medium	<input type="checkbox"/> 16	<input type="checkbox"/> 17	<input type="checkbox"/> 18	
		Tall	<input type="checkbox"/> 18	<input type="checkbox"/> 19	<input type="checkbox"/> 20	

Standard Back Post

Select option below only if a Standard Back Post option is selected on chart at left

Straight without Push Handle _____ N/C
 Straight with Push Handle _____ \$50
 8° Bend with Push Handle _____ \$80
 Not available with Short Back Height Range
 Straight with Fold Down Push Handle _____ \$125
 Not available with Heavy Duty option

Push Handle Options

Not available with Straight Back Post with Fold Down Push Handle or Heavy Duty Option
 Low Profile Push Handle _____ \$65
 2" shorter than standard push handles

Backrest continued

Back Upholstery

Select option below only if Standard or Angle Adjustable Back Post option is selected on chart on previous page

- Standard _____ Std
Not available with Short or Medium Back Height Ranges on 20" - 22" Seat Widths
- Tension Adjustable E2611 _____ \$245
- Omit Back Upholstery _____ - \$10

Military Patch

Attached only if ordered with Tension Adjustable Back Upholstery

- US Military Service Patch _____ \$90
 Army Navy Marines Air Force Coast Guard
- US Purple Heart Patch _____ \$90
- Purple Heart with Military Service Patch Combo _____ \$145
 Army Navy Marines Air Force Coast Guard

Footrest

Hanger / Front Frame Bend

- 70° _____ N/C
- 80° _____ N/C
Not available with 7" or 8" Casters
- Elevating Leg Rests (ELR) E0990 _____ \$295
- Pro Elevating Leg Rests (PRO ELR) E0990 **New!** _____ \$350
- Omit Hangers _____ - \$20

Footplate

- Composite _____ N/C
Not available with PRO ELR
- Composite Angle Adjustable K0040 _____ \$185
Not available with ELR or 14" Seat Width
- Aluminum Angle Adjustable K0040 _____ \$185
- One Piece Flip Up Angle Adjustable K0040 _____ \$190
Not available with ELR or PRO ELR

Footrest / Extension / Seat to Footrest Length


Select Seat to Footrest Length in row corresponding to desired Footrest and Extension. Footrest Length of at least 2.5" shorter than Front Seat Height recommended. Extension Tube Hangers and Fixed Front Frame measurements are based on 70° Hanger / Front Frame Bend.

Footrest	Extension	Seat to Footrest Length									
Extension Tube Hangers (60°, 70°, 80°)	Short	<input type="checkbox"/> 11	<input type="checkbox"/> 11.5	<input type="checkbox"/> 12	<input type="checkbox"/> 12.5	<input type="checkbox"/> 13	<input type="checkbox"/> 13.5	<input type="checkbox"/> 14	<input type="checkbox"/> 14.5	<input type="checkbox"/> 15	
	Medium	<input type="checkbox"/> 13	<input type="checkbox"/> 13.5	<input type="checkbox"/> 14	<input type="checkbox"/> 14.5	<input type="checkbox"/> 15	<input type="checkbox"/> 15.5	<input type="checkbox"/> 16	<input type="checkbox"/> 16.5	<input type="checkbox"/> 17	
	Long	<input type="checkbox"/> 15.5	<input type="checkbox"/> 16	<input type="checkbox"/> 16.5	<input type="checkbox"/> 17	<input type="checkbox"/> 17.5	<input type="checkbox"/> 18	<input type="checkbox"/> 18.5	<input type="checkbox"/> 19	<input type="checkbox"/> 19.5	
Elevating Leg Rests (ELR)	Short	<input type="checkbox"/> 12	<input type="checkbox"/> 13	<input type="checkbox"/> 14	<input type="checkbox"/> 15	<input type="checkbox"/> 16					
	Long	<input type="checkbox"/> 13	<input type="checkbox"/> 14	<input type="checkbox"/> 15	<input type="checkbox"/> 16	<input type="checkbox"/> 17	<input type="checkbox"/> 18				
Extended Length ELR	Short	<input type="checkbox"/> 18	<input type="checkbox"/> 19	<input type="checkbox"/> 20	<input type="checkbox"/> 21	<input type="checkbox"/> 22					
	Long	<input type="checkbox"/> 18	<input type="checkbox"/> 19	<input type="checkbox"/> 20	<input type="checkbox"/> 21	<input type="checkbox"/> 22	<input type="checkbox"/> 23	<input type="checkbox"/> 24			
PRO ELR Low Pivot	Short	<input type="checkbox"/> 12 - 17.5									
	Long	<input type="checkbox"/> 15 - 21									
PRO ELR High Pivot	Short	<input type="checkbox"/> 10.5 - 15.5									
	Long	<input type="checkbox"/> 13 - 19									
Fixed Front Frame (70°, 80°)	Short	<input type="checkbox"/> 11.5	<input type="checkbox"/> 12	<input type="checkbox"/> 12.5	<input type="checkbox"/> 13	<input type="checkbox"/> 13.5	<input type="checkbox"/> 14	<input type="checkbox"/> 14.5	<input type="checkbox"/> 15		
	Medium	<input type="checkbox"/> 13	<input type="checkbox"/> 13.5	<input type="checkbox"/> 14	<input type="checkbox"/> 14.5	<input type="checkbox"/> 15	<input type="checkbox"/> 15.5	<input type="checkbox"/> 16	<input type="checkbox"/> 16.5	<input type="checkbox"/> 17	
	Long	<input type="checkbox"/> 15.5	<input type="checkbox"/> 16	<input type="checkbox"/> 16.5	<input type="checkbox"/> 17	<input type="checkbox"/> 17.5	<input type="checkbox"/> 18	<input type="checkbox"/> 18.5	<input type="checkbox"/> 19	<input type="checkbox"/> 19.5	

Hanger Release

If no selection is made Classic Hanger Release will be used

- Classic _____ Std
- 4-Way - In line **New!** _____ N/C
See reference A to the right. Not available with ELR or PRO ELR.
- 4-Way - Out **New!** _____ N/C
See reference B to the right



Footrest continued

- Heel Loop**
- Heel Loops E0951 _____ \$30
Not available with One Piece Flip Up Angle Adjustable Footplate
 - Omit Heel Loops _____ N/C
- Calf Strap**
- Velcro Adjustable K0038 _____ \$35
 - Padded Velcro Adjustable K0038 _____ \$68
 - Bodypoint Padded AeroMesh K0038 _____ \$102

- Residual Limb Support**
- Residual Limb Support E1020
 - Left _____ \$395
 - Right _____ \$395
 - Omit Left Hanger _ - \$10
 - Omit Right Hanger _ - \$10
- Cover**
- Smooth Air Mesh
 - Smooth Air Mesh
- Pad Size**
- 5 x 8 9 x 14
 - 7 x 10 11 x 10
 - 9 x 10 11 x 14

Front Seat Height and Casters

Caster Size / Front Frame Type / Front Seat Height

Select Front Seat Height in row corresponding to desired caster size

Caster Size	Front Seat Height										
4"	<input type="checkbox"/> 15 ¹	<input type="checkbox"/> 15.5	<input type="checkbox"/> 16	<input type="checkbox"/> 16.5	<input type="checkbox"/> 17	<input type="checkbox"/> 17.5	<input type="checkbox"/> 18	<input type="checkbox"/> 18.5	<input type="checkbox"/> 19	<input type="checkbox"/> 19.5	<input type="checkbox"/> 20
5"	<input type="checkbox"/> 16	<input type="checkbox"/> 16.5	<input type="checkbox"/> 17	<input type="checkbox"/> 17.5	<input type="checkbox"/> 18	<input type="checkbox"/> 18.5	<input type="checkbox"/> 19	<input type="checkbox"/> 19.5	<input type="checkbox"/> 20	<input type="checkbox"/> 20.5	
6"	<input type="checkbox"/> 17	<input type="checkbox"/> 17.5	<input type="checkbox"/> 18	<input type="checkbox"/> 18.5	<input type="checkbox"/> 19	<input type="checkbox"/> 19.5	<input type="checkbox"/> 20	<input type="checkbox"/> 20.5	<input type="checkbox"/> 21		
7"	<input type="checkbox"/> 18	<input type="checkbox"/> 18.5	<input type="checkbox"/> 19	<input type="checkbox"/> 19.5	<input type="checkbox"/> 20	<input type="checkbox"/> 20.5	<input type="checkbox"/> 21				
8"	<input type="checkbox"/> 19	<input type="checkbox"/> 19.5	<input type="checkbox"/> 20	<input type="checkbox"/> 20.5	<input type="checkbox"/> 21						

¹ Not available with 4 x 1" Poly Caster

- Caster Type**
- .75" Rollerblade _____ N/C
Not available with 5", 6", 7" and 8" Caster Sizes
 - .75" Lighted Rollerblade _____ \$70
Not available with 5", 6", 7" and 8" Caster Sizes
 - 1" Poly _____ N/C
 - 1" Poly Aluminum _____ \$160
Not available with 6", 7" and 8" Caster Sizes
 - 1" Pneumatic _____ \$80
Not available with 4", 5" and 7" Caster Sizes
 - 1.5" Poly _____ \$50
Not available with 7" and 8" Caster Sizes
 - 1.5" Soft Roll Aluminum _____ \$175
Not available with 7" and 8" Caster Sizes
 - 2" Poly _____ \$70
Not available with 4", 5", 7" and 8" Caster Sizes
 - 2" Pneumatic w/ Foam Insert _____ \$225
Not available with 4", 5", 6" and 7" Caster Sizes

- Caster Housing**
- Standard _____ N/C
 - Performance _____ \$35

Maximum Seat Dump Reference
(Front Seat Height minus Rear Seat Height)

Seat Depth	14"	15"	16"	17"	18"	19"	20"
Maximum Dump Standard Caster Housing	1.5"	2"	2"	2"	2.5"	2.5"	2.5"
Maximum Dump Performance Caster Housing	3"	3"	3"	3.5"	3.5"	4"	4"

Depending on configuration, chairs with Fixed Front Frame Type may achieve up to an additional .75"

- Fork Stem**
- Best Fit _____ Std
Ki Mobility will pick optimal Fork Stem based on Caster Size and Front Seat Height
 - Standard Fork Stem _____ N/C
 - 3/4" Longer Fork Stem _____ N/C
 - 1 1/2" Longer Fork Stem _____ N/C

- Fork Options**
- Frog Legs Ultra-Sport Fork - Black _____ \$425
Available only with 4", 5" or 6" casters. Not available with Caster Pin Locks, Heavy Duty Option, 2" Poly caster or 2" Pneumatic w/ Foam Insert caster.
 - End User Weight Required: _____ lbs**
 - Frog Legs PHASE 2 Carbon Fork **New!** _____ \$425
Lightweight shock absorbing fork. Available only with 4", 5" or 6" casters. Not available with Caster Pin Locks, Heavy Duty Option, 2" Poly caster or 2" Pneumatic w/ Foam Insert caster.
 - Caster Pin Locks K0073 _____ \$195
Not available with Performance Caster Housing, 1 1/2" Longer Fork Stem or Frog Legs Fork

Wheel Locks

- Wheel Lock**
- Push to Lock _____ Std
 - Push to Lock Flush Mount _____ N/C
 - Pull to Lock _____ N/C
 - Low Profile Scissor Lock _____ \$45
Not available with Low Poly tire
 - Short Thro Scissor Lock _____ \$110
 - Omit Wheel Locks _____ N/C

- Wheel Lock Options**
- Extension Handles E0961 _____ \$67
Available only with Push to Lock or Pull to Lock Wheel Locks
 - Grade Aids E0974 _____ \$150
Available only with Push to Lock Wheel Locks and Pneumatic or Pneumatic with Airless Insert Rear Tire Types

Rear Seat Height and Wheels

Wheel Size / Tire Profile / Rear Seat Height
Select Rear Seat Height in row corresponding to Wheel Size and Tire Profile. 14" Seat Depth not available with 24", 25" or 26" Rear Wheels.
 Requires curved axle plate Any axle plate can be used (Curved, Standard, Amputee, One Arm Drive)

Wheel Size	Tire Profile	Rear Seat Height									
20"	Low	<input type="radio"/> 11.5	<input type="checkbox"/> 12 ¹	<input type="radio"/> 12.5	<input type="checkbox"/> 13 ¹	<input type="radio"/> 13.5	<input type="checkbox"/> 14	<input type="radio"/> 14.5	<input type="checkbox"/> 15	<input type="radio"/> 15.5	<input type="radio"/> 16
	Full	<input type="radio"/> 12	<input type="checkbox"/> 12.5 ¹	<input type="radio"/> 13	<input type="checkbox"/> 13.5 ¹	<input type="radio"/> 14	<input type="checkbox"/> 14.5	<input type="radio"/> 15	<input type="checkbox"/> 15.5	<input type="radio"/> 16	<input type="radio"/> 16.5
22"	Low	<input type="radio"/> 12.5	<input type="checkbox"/> 13 ¹	<input type="radio"/> 13.5	<input type="checkbox"/> 14 ¹	<input type="radio"/> 14.5	<input type="checkbox"/> 15	<input type="radio"/> 15.5	<input type="checkbox"/> 16	<input type="radio"/> 16.5	<input type="radio"/> 17
	Full	<input type="radio"/> 13	<input type="checkbox"/> 13.5 ¹	<input type="radio"/> 14	<input type="checkbox"/> 14.5 ¹	<input type="radio"/> 15	<input type="checkbox"/> 15.5	<input type="radio"/> 16	<input type="checkbox"/> 16.5	<input type="radio"/> 17	<input type="radio"/> 17.5
24"	Low	<input type="radio"/> 13.5	<input type="checkbox"/> 14 ¹	<input type="radio"/> 14.5	<input type="checkbox"/> 15 ¹	<input type="radio"/> 15.5	<input type="checkbox"/> 16	<input type="radio"/> 16.5	<input type="checkbox"/> 17	<input type="radio"/> 17.5	<input type="radio"/> 18
	Full	<input type="radio"/> 14	<input type="checkbox"/> 14.5 ¹	<input type="radio"/> 15	<input type="checkbox"/> 15.5 ¹	<input type="radio"/> 16	<input type="checkbox"/> 16.5	<input type="radio"/> 17	<input type="checkbox"/> 17.5	<input type="radio"/> 18	<input type="radio"/> 18.5
25"	Low	<input type="radio"/> 14	<input type="checkbox"/> 14.5 ¹	<input type="radio"/> 15	<input type="checkbox"/> 15.5 ¹	<input type="radio"/> 16	<input type="checkbox"/> 16.5	<input type="radio"/> 17	<input type="checkbox"/> 17.5	<input type="radio"/> 18	<input type="radio"/> 18.5
26"	Low	<input type="radio"/> 14.5	<input type="checkbox"/> 15 ¹	<input type="radio"/> 15.5	<input type="checkbox"/> 16 ¹	<input type="radio"/> 16.5	<input type="checkbox"/> 17	<input type="radio"/> 17.5	<input type="checkbox"/> 18	<input type="radio"/> 18.5	<input type="radio"/> 19
	Full	<input type="radio"/> 15	<input type="checkbox"/> 15.5 ¹	<input type="radio"/> 16	<input type="checkbox"/> 16.5 ¹	<input type="radio"/> 17	<input type="checkbox"/> 17.5	<input type="radio"/> 18	<input type="checkbox"/> 18.5	<input type="radio"/> 19	<input type="radio"/> 19.5

¹ Not available with Amputee Axle Plate

- One Arm Drive**
- One Arm Drive E0958 _____ \$1095
Available only with 24" Mag wheels. Not available with 14" Rear Seat Height in combination with 24" rear wheel or 4" camber. Standard Handrim will be Aluminum Anodized. Axle not quick release.
 - Drive Handrim Side: Left Right
 - Drive Handrim Type: Aluminum _____ N/C
 Poly Coated _____ \$200

- Axle Plate**
- No selection needed with One Arm Drive*
- Standard _____ Std
 - Curved _____ \$85
 - Amputee _____ \$175
Not available with Swing Away Armrests or Wheelchair Drive Ready Option

- Rear Wheel Type**
- 24" Spoke _____ N/C
 - Spoke (20", 22", 25" or 26") _____ \$200
 - Mag _____ \$50
Not available in 25" or 26" rear wheels
 - Superlight Spoke _____ \$295
Available only in 24" and 25" rear wheels. Not available with Heavy Duty Option.
 - Spinery Spox _____ \$795
Not available in 20" rear wheel
 - Color: Black Red White Blue Yellow
 - Handrim: Short Tab Long Tab
Available only with Aluminum Anodized and Plastic Coated Handrims
 - Spinery LX _____ \$995
Not available in 20" rear wheel
 - Color: Black Red White Blue Yellow
 Pink Green Orange
 - Handrim: Short Tab Long Tab
Available only with Aluminum Anodized and Plastic Coated Handrims

- Rear Tire Type**
- Full Profile**

 - Street Pneumatic 1 3/8" E2211, E2212 _____ Std
Not available on 25" rear wheel
 - Full Poly 1 3/8" E2218 _____ \$95
Not available on Superlight or 25" rear wheel
 - Pneumatic with Airless Insert 1 3/8" E2211, E2213 _____ \$125
Not available on Superlight or 25" rear wheel
 - Schwalbe HS-127 1 3/8" E2211, E2212 **New!** _____ \$125
Available only on 24" rear wheel. Not available on Mag wheels.
 - Euro Trek Knobby 1 3/8" E2211, E2212 _____ \$95
Available only on 24" rear wheel

Low Profile

 - Low Poly 1" E2218 _____ \$75
Not available on Superlight, 25" or 26" rear wheels
 - Sentinel Solid Tire 1" E2218 _____ \$95
Not available on Superlight wheels
 - High Pressure 1" E2211, E2212 _____ \$95
Not available on Mag wheels
 - Schwalbe One 1" E2211, E2212 **New!** _____ \$125
Not available on Mag, 20", 22" or 26" rear wheels
 - Schwalbe Marathon Plus Evolution 1" E2211, E2212 _____ \$125
Not available on Mag, 20" or 22" rear wheels

Rear Seat Height and Wheels continued

Handrim

No selection needed with One Arm Drive

Standard Handrims

- Aluminum Anodized _____ Std
- Superlight Aluminum Anodized _____ \$50
Available only on 24" and 25" Superlight and Spinergy wheels. Not available with Heavy Duty Option.
- Plastic Coated _____ \$95
- Projection (8 Vertical Projections) _____ \$295
Available on 22", 24", 25" and 26" Spoke wheels & 24" Mag wheels w/full tires
- Spinergy Flex Rim _____ \$600
Available only on 24" and 25" Spinergy LX wheels
- Omit Handrim _____ N/C

Handrim Option

- 3M Micro Grip _____ \$65
Not available on Plastic Coated or Projection Handrims

Axles

Axle

No selection needed with Wheelchair Drive Ready Option and One Arm Drive

- Quick Release _____ Std
- Quad Quick Release _____ \$175
Not available with Mag wheels

Armrests

Armrest Please select Arm Pad option

- Height Adjustable T-Arm (8.5" - 13.5") E0973 _____ \$375
- Height Adjustable T-Arm Low (7.5" - 11.5") E0973 _____ \$375
- Height Adjustable Flip Back T-Arm (9.5"-13.5") E0973 **New!** _____ \$275
- Flip Back (10.5") _____ \$425
Not available with Seat Depths less than 16" or Angle Adjustable Back Posts

Accessories

Removable Side Guard

- Composite Side Guards K0108 _____ \$140
Not available with T-Arm or Flip Back Armrests

Anti-Tip

- Rear Anti-Tips E0971 _____ \$140

Positioning Belt

- 1 1/2" Auto Buckle E0978 _____ \$70
- 1 1/2" Padded Auto Buckle E0978 _____ \$155
- 2" Airline Buckle E0978 _____ \$110
- Bodypoint Evoflex - Pelvic Stabilizer E0978 _____ \$183
Size: Small (16-24") Medium (19-29") Large (25-37")
- Bodypoint Belt Mounting Kit _____ \$105
Style: Bracket Band Clamp 1"

O₂ Holder

- O₂ Holder E2208 _____ \$295
 Right Side Left Side

Natural Fit Handrims

Not available on 20" rear wheel or 22" Mag wheels. Not recommended with Mag wheels.

- Natural Fit LT Regular Grip _____ \$445
- Natural Fit Standard Regular Grip _____ \$445
- Natural Fit LT Super Grip _____ \$495
- Natural Fit Standard Super Grip _____ \$495
- Natural Fit LT Omit Grip _____ \$345
- Natural Fit Standard Omit Grip _____ \$345

Axle Position

- Rear Wheel Camber: 0° 2° 4°
- Center of Gravity: Rear Center Forward

Arm Pad

Not available with Swing Away Armrests

- Standard Arm Pad _____ N/C Desk Full
- Waterfall Arm Pad _____ \$25 Desk Full
Available only with Height Adjustable T-Arms

Armrest continued

- Swing Away (10.5" - 12.5") E0973 _____ \$195
Not available with Angle Adjustable Back Posts or Arm Pads
- Omit Armrests _____ N/C

Seat Pan

- Removable Aluminum Seat Pan E2231 _____ \$245
Not available with Transit System. Pan depth 1/2" shorter than Seat Depth.

Other Accessories

- Tool Kit _____ \$45
- Removable Underseat Pouch _____ \$50
- Ki Mobility Backpack _____ \$65
- Luggage Carrier **New!** _____ \$85
Not available with ELR, PRO ELR or Residual Limb Support
- Black Neoprene Impact Guards _____ \$85
- Clear Spoke Guards K0065 _____ \$115
Not available with Mag wheels
- Black Spoke Guards K0065 _____ \$115
Not available with 20" rear wheel or Mag wheels
- Cane and Crutch Holder _____ \$195
- IV Pole _____ \$265
Available only with Standard Straight Back Posts

Catalyst Front Seat Height Matrix																										
4" Caster			5" Caster			6" Caster			7" Caster			8" Caster														
FSH	Stem	Fork	FSH	Stem	Fork	FSH	Stem	Fork	FSH	Stem	Fork	FSH	Stem	Fork												
15 ¹	Std	4"	16	Std	4"	17	Std	5"	18	Std	6"	19	Std	6"												
15.5			16.5	3/4		18	3/4		18.5			20	3/4													
16 ²	17.5		1 1/2	18.5		1 1/2	19		1 1/2	20.5		1 1/2														
16.5	3/4		16	Std	5"	17	Std	6"	19.5	3/4		7"	19	Std	7"											
16.5 ²	16.5		17.5			18			19.5				20			19.5	20									
17	1 1/2		17	3/4		18	1 1/2		18.5	1 1/2			19.5	1 1/2		20	Std	20.5	Std							
15.5	Std	5"	17.5	3/4		6"	18		3/4	7"	18.5		Std	7"		19.5	3/4	7"								
16			17.5				18				19					19.5			20	20.5						
16.5	3/4		5"	18			1 1/2		6"		18.5		3/4			7"	19		Std	7"	19.5	Std	7"			
17				16.5	18.5			19			19.5	20			20.5											
17	1 1/2			5"	17		Std	6"			19	1 1/2	7"		19		3/4		7"		21	1 1/2		7"		
17.5					17.5						19.5				19.5						20.5				21	
17	Std	6"			17.5	3/4	6"			18	Std	7"		20	1 1/2		7"	21			1 1/2	7"				
17.5					18					18.5				19				20.5							21	
17.5	3/4		6"		18.5	1 1/2			6"	19	Std			7"	20	1 1/2		7"		8" x 2" caster not available with 6" fork						
18					18					19.5					20								20.5			
18.5	1 1/2			6"	18.5	1 1/2		6"		18.5	3/4		7"		21	1 1/2			7"				8" x 2" caster not available with 6" fork			
19					19					19					20											20.5
18	Std	7"			18	Std	7"			19.5	3/4	7"			8" x 2" caster not available with 6" fork											
18.5					18.5					20																20.5
19	3/4		7"		19	Std			7"	19.5	1 1/2			7"				8" x 2" caster not available with 6" fork								
19.5					18.5					20											20.5					
19.5	1 1/2			7"	19	3/4		7"		20.5	1 1/2		7"								8" x 2" caster not available with 6" fork					
20					19.5					21														21		
					19.5	1 1/2	7"			6" x 2" caster not available with 5" fork					8" x 2" caster not available with 6" fork											
					20					20.5																
			20.5		21																					

¹ Front Seat Height not available with 4 x 1" Poly Caster
² Not available with 4 x 1" Poly, 4 x 1.5" Poly or 4 x 1.5" Soft Roll Aluminum

Actual seat height may vary depending upon configuration.

Cushions

General Use Cushions

STEP 1: Choose Cushion

- 2" General Use Foam Cushion (14"x14" to 20"x20") E2601 ____ \$99
- 2" General Use Foam Cushion (Over 20"x20") E2601 ____ \$250

STEP 2: Choose Cushion Size

- Match to Chair Size
- Other Size _____

Matrx Seat Cushions

Cushion size is determined by chair dimensions. Certain sizes will have higher price and longer lead time; visit www.motionconcepts.com/matrx_international.html for size specifications.

STEP 1: Choose Cushion

- PS Cushion _____ \$414
Cushions 14" wide have a weight capacity of 150 lbs. Cushions 15" through 20" have a weight capacity of 300 lbs.
- Vi Cushion _____ \$510
Weight capacity of 300 lbs
- Libra Cushion _____ \$582
Weight capacity of 300 lbs

STEP 2: Choose Cushion Size

- Match to Chair Size
- Other Size - Provide Part Number _____

Other Matrx Products

Must provide part number

- _____

ROHO Seat Cushions

Cushion size is determined by chair dimensions. Certain sizes will have higher price and longer lead time; visit therohogroup.com for size specifications.

STEP 1: Choose Cushion

Standard Cushions

For chairs under 22" wide

- HIGH PROFILE® Single Compartment (E2622) ____ \$603
- MID PROFILE™ Single Compartment (E2622) ____ \$603
- LOW PROFILE® Single Compartment (E2622) ____ \$603
- QUADTRO SELECT® High Profile® (E2624) ____ \$603
- QUADTRO SELECT® Mid Profile® (E2624) ____ \$603
- HYBRID ELITE® Single Comp Smart Check Ready (E2622) ____ \$623
- CONTOUR SELECT® (E2624) ____ \$603

X Large Cushions

For 22" wide chairs

- HIGH PROFILE® Single Compartment X Large (E2623) ____ \$670
- LOW PROFILE® Single Compartment X Large (E2623) ____ \$670
- QUADTRO SELECT® High Profile® X Large (E2625) ____ \$670

Cushions with Heavy Duty Cover

- HIGH PROFILE® Single Compartment w/ HD Cover (E2622) ____ \$629
- QUADTRO SELECT® High Profile® w/ HD Cover (E2624) ____ \$629

STEP 2: Choose Cushion Size

- Match to Chair Size
- Other Size - Provide Part Number _____

Backs

<p>Matrix Backs</p> <p><input type="checkbox"/> Elite Back _____ \$750 <i>Not available with 22" Seat Width. Back Height of 14" not available with 19", 20" and 21" Seat Widths.</i> Back Height: <input type="checkbox"/> 10" <input type="checkbox"/> 12" <input type="checkbox"/> 14" <input type="checkbox"/> 16" <input type="checkbox"/> 18" <input type="checkbox"/> 20"</p> <p><input type="checkbox"/> Elite Deep Back _____ \$912 <i>Not available with 22" Seat Width. Back Height of 10" not available with 14" Seat Width.</i> Back Height: <input type="checkbox"/> 10" <input type="checkbox"/> 12" <input type="checkbox"/> 16" <input type="checkbox"/> 18" <input type="checkbox"/> 20"</p> <p><input type="checkbox"/> Mini Back _____ \$750 <i>Available only with 14" and 15" Seat Widths. Comes standard with Kidabra cover.</i> Back Height: <input type="checkbox"/> 7" <input type="checkbox"/> 10" <input type="checkbox"/> 12" <input type="checkbox"/> 14"</p>	<p>Therafin Headrest for MATRX Back Systems</p> <p><input type="checkbox"/> Headrest, Multi-Axis Offset Hardware with Smooth/Soft Pad (E1028, E0955) 34187 _____ \$420 PAD SIZE: <input type="checkbox"/> 5 x 8 <input type="checkbox"/> 5 x 10 <input type="checkbox"/> 5 x 14</p>
	<p>Other Matrix Mounting Hardware</p> <p><i>Must provide part number</i></p> <p><input type="checkbox"/> _____ For mounting hardware other than quick release</p>
	<p>Other Matrix Products</p> <p><i>Must provide part number</i></p> <p><input type="checkbox"/> _____</p>

ROHO Back Systems (Available for chairs up to 20" wide)

STEP 1: Choose Back

- | | |
|---|---|
| <p><input type="checkbox"/> AGILITY™ Min. Contour w/ Foam Only (E2611) _____ \$549
 BACKREST HEIGHT: <input type="checkbox"/> 10" <input type="checkbox"/> 13" <input type="checkbox"/> 16" <input type="checkbox"/> 18"</p> <p><input type="checkbox"/> AGILITY™ Min. Contour w/ Foam & Lumbar Pad (E2613) _____ \$659
 BACKREST HEIGHT: <input type="checkbox"/> 10" <input type="checkbox"/> 13" <input type="checkbox"/> 16" <input type="checkbox"/> 18"</p> <p><input type="checkbox"/> AGILITY™ Min. Contour w/ Air Only (E2620) _____ \$769
 BACKREST HEIGHT: <input type="checkbox"/> 10" <input type="checkbox"/> 13" <input type="checkbox"/> 16" <input type="checkbox"/> 18"</p> <p><input type="checkbox"/> AGILITY™ Min. Contour w/ Air & Lumbar Pad (E2615) _____ \$769
 BACKREST HEIGHT: <input type="checkbox"/> 10" <input type="checkbox"/> 13" <input type="checkbox"/> 16" <input type="checkbox"/> 18"</p> | <p><input type="checkbox"/> AGILITY™ Mid Contour w/ Air Only (E2620) _____ \$781
 BACKREST HEIGHT: <input type="checkbox"/> 10" <input type="checkbox"/> 13" <input type="checkbox"/> 16" <input type="checkbox"/> 18" <input type="checkbox"/> 20"</p> <p><input type="checkbox"/> AGILITY™ Mid Contour w/ Air & Lumbar Pad (E2615) _____ \$781
 BACKREST HEIGHT: <input type="checkbox"/> 10" <input type="checkbox"/> 13" <input type="checkbox"/> 16" <input type="checkbox"/> 18" <input type="checkbox"/> 20"</p> <p><input type="checkbox"/> AGILITY™ Max Contour w/ Air Only (E2620) _____ \$920
 BACKREST HEIGHT: <input type="checkbox"/> 16" <input type="checkbox"/> 18" <input type="checkbox"/> 20"</p> <p><input type="checkbox"/> AGILITY™ Max Contour w/ Air & Lumbar Pad (E2615) _____ \$920
 BACKREST HEIGHT: <input type="checkbox"/> 16" <input type="checkbox"/> 18" <input type="checkbox"/> 20"</p> |
|---|---|

STEP 2: Choose Back Hardware

- Fixed _____ N/C
Not compatible with 18" and 20" tall backs
- Quick Release _____ N/C

ROHO Laterals for AGILITY™ Back Systems

Sold in pairs

- Swing Away with Curved Pad (E1028/E0956) _____ \$1014
 PAD SIZE: Medium Large X-Large
 4 T x 4 L 4 T x 5 L 5 T x 5 L
- OFFSET TYPE: Flush Offset 1" Offset 2" Offset
- Fixed with Flat Pad (E0956) _____ \$452
 PAD SIZE: Medium Large X-Large
 4 T x 4 L 4 T x 5 L 5 T x 5 L
- OFFSET TYPE: Flush Offset

Therafin Headrest for AGILITY Back Systems

- Headrest, Multi-Axis Offset Hardware with Smooth/Soft Pad (E1028, E0955) 34187 _____ \$420
 PAD SIZE: 5 x 8 5 x 10 5 x 14



Vest

Not compatible with Back Upholstery

- Full Shape, Comfort Strap (E0960) _____ \$167
 - Small 30589 Medium 30590 Large 30591
 - X-Large 30592 2X-Large 30593
- Full Shape, Extended Strap (E0960) _____ \$167
 - Small 30594 Medium 30595 Large 30596
 - X-Large 30597 2X-Large 30598
- Trim Shape, Comfort Strap (E0960) _____ \$169
 - Small 30618 Medium 30619 Large 30620
 - X-Large 30621 2X-Large 30622
- Trim Shape, Extended Strap (E0960) _____ \$169
 - Small 30604 Medium 30624 Large 30625
 - X-Large 30626 2X-Large 30627

Shoulder Harness

With metal clips and cam locks. Not compatible with Back Upholstery.

- Small (E0960) 61403 _____ \$157
- Medium (E0960) 61503 _____ \$164
- Large (E0960) 61703 _____ \$170
- X-Large (E0960) 61803 _____ \$177

Hip Belt

- 4 Point Y-Style, Metal Buckle (E0978) _____ \$173
 - Small 60024 Medium 60025 Large 60026
- Dual Rear Pull, Metal Attaching Clips (E0978)
 - 1" Straps, 6.25"W x 1.75"H Pads 60161 _____ \$97
 - 1.5" Straps, 7.25"W x 2.25"H Pads 60171 _____ \$106
 - 1.5" Straps, 9.25"W x 2.5"H Pads 60181 _____ \$106
- Single Pull, Metal Buckle (E0978)
 - 1" Straps, 4.25"W x 1.75"H Pads 61001 _____ \$86
 - 1" Straps, 5.25"W x 1.75"H Pads 61101 _____ \$86
 - 1" Straps, 6.25"W x 1.75"H Pads 61201 _____ \$86
 - 1.5" Straps, 5.25"W x 2.25"H Pads 62001 _____ \$95
 - 1.5" Straps, 9.25"W x 2.5"H Pads 62201 _____ \$95
 - 2" Straps, 9.25"W x 3"H Pads 64001 _____ \$102
 - 2" Straps, 11.25"W x 3"H Pads 64101 _____ \$102

Cane Mounted Swing Away Lateral Thoracic Support

Sold as pair. Not compatible with Back Upholstery.

- Swing Away Kits w/ 4"x6" Gel Pads (E0956, E1028) 50021 _____ \$478
- Swing Away Kits, No Pads (E0956, E1028) 50022 _____ \$309

Knee Adductor

- Flip Down Hardware (E1028) _____ \$634
Sold in pairs, one right side and one left side
 - Short Arm, Fits 3/4", 7/8" and 1" Tubing 30248
 - Long Arm, Fits 3/4", 7/8" and 1" Tubing 30284
- Knee Adductor Pad with Neoprene Cover (E0956) _____ \$117
Sold single
 - Small Pad, 4.25" x 3.25" 50490 Quantity: _____
 - Medium Pad, 6.25" x 4.25" 50491 Quantity: _____

Foot / Shoe Holder

Compatible only with Aluminum Angle Adjustable and Composite Angle Adjustable Footplates

- FootSure Ankle Support (K0108) _____ \$178
Sold in pairs
 - X-Small w/ Buckle Closure, 6"-8" Circum. 60250
 - Small w/ Buckle Closure, 7"-8.5" Circum. 60251
 - Medium w/ Buckle Closure, 8"-11" Circum. 60252
 - Large w/ Buckle Closure, 9"-12" Circum. 60253
 - X-Large w/ Buckle Closure, 12.5"-15" Circum. 60254
 - X-Small w/ Hook-N-Loop Closure, 5"-6" Circum. 60255
 - Small w/ Hook-N-Loop Closure, 6"-8" Circum. 60256
 - Med w/ Hook-N-Loop Closure, 8"-10.5" Circum. 60257
 - Large w/ Hook-N-Loop Closure, 9"-12" Circum. 60258
 - X-Lrg w/ Hook-N-Loop Closure, 12"-15" Circum. 60259
- Padded Footbox (K0108)
 - Small, Full Design, Fits 14" Chair 30111 _____ \$287
 - Medium, Full Design, Fits 16" Chair 30112 _____ \$293
 - Large, Full Design, Fits 18" Chair 30113 _____ \$298
 - Small, Split Design, Fits 14" Chair 30152 _____ \$367
 - Medium, Split Design, Fits 16" Chair 30153 _____ \$373
 - Large, Split Design, Fits 18" Chair 30154 _____ \$378
- Molded Shoe Holders (K0108) _____ \$162
Sold in pairs
 - Small with Plain Straps, 3.5" x 5.25" 30293
 - Medium with Plain Straps, 4" x 6.5" 30303
 - Large with Plain Straps, 4.75" x 8" 30309
 - X-Large with Plain Straps, 5" x 10" 30316



Half Tray

Not compatible with Swing Away Armrests

- Standard Shape, Clear, Flip-Away (E0950, E1028) _____ \$246
 - Right Side 41310 Left Side 41311
- Standard Shape, Padded, Flip-Away (E0950, E1028) _____ \$269
 - Right Side 41350 Left Side 41353
- Trimline Shape, Padded, Flip-Away (E0950, E1028) _____ \$269
 - Right Side 41352 Left Side 41351

Clear Tray with Border and Pre-Drilled Holes

STEP 1: Choose Tray

Chair Width	Tray Dimension (W x D)	Body Cutout (W x D)	HCPCS	Part Number	Price	
14"	20 x 18	10.5 x 6	E0950	33210	\$220	<input type="checkbox"/>
	20 x 18	11.5 x 6.5	E0950	33211	\$220	<input type="checkbox"/>
15"	21 x 18	11.5 x 6.5	E0950	33212	\$220	<input type="checkbox"/>
	21 x 18	12.5 x 7	E0950	33213	\$220	<input type="checkbox"/>
16"	22 x 21	12.5 x 7.5	E0950	33214	\$255	<input type="checkbox"/>
	22 x 21	13.5 x 8.5	E0950	33215	\$255	<input type="checkbox"/>
17"	23.5 x 22	13.5 x 8.5	E0950	33216	\$255	<input type="checkbox"/>
	23.5 x 22	14.5 x 9.5	E0950	33217	\$285	<input type="checkbox"/>
18"	24.5 x 22	14.5 x 9.5	E0950	33218	\$285	<input type="checkbox"/>
	24.5 x 22	15.5 x 10.5	E0950	33219	\$285	<input type="checkbox"/>
19"	26.5 x 22	15.5 x 10.5	E0950	33220	\$285	<input type="checkbox"/>
	26.5 x 22	16.5 x 11	E0950	33221	\$285	<input type="checkbox"/>
20"	26.5 x 22	17 x 11	E0950	33222	\$309	<input type="checkbox"/>
	26.5 x 22	18 x 11	E0950	33223	\$309	<input type="checkbox"/>
21"	27.5 x 23.5	17.5 x 11	E0950	33224	\$309	<input type="checkbox"/>
	27.5 x 23.5	19 x 11	E0950	33225	\$309	<input type="checkbox"/>
22"	28.5 x 23.5	18.5 x 11	E0950	33226	\$309	<input type="checkbox"/>
	28.5 x 23.5	20 x 11	E0950	33227	\$309	<input type="checkbox"/>

STEP 2: Choose Tray Hardware

- E-Z Lock Clamps, Tube Arm (E1028) 31072 _____ \$150
Compatible only with Swing Away Armrests. Sold in Pairs.
- E-Z Lock Clamps, Standard (E1028) 31098 _____ \$142
Not compatible with Swing Away Armrests. Sold in Pairs.

Communication Mount Hardware

Choose Device Holder AND Mount

STEP 1: What DEVICE will be used?

- Smart Phone 32359 _____ \$150
- Large Smart Phone 32354 _____ \$150
- I-Pad 32355 _____ \$184
- I-Pad Mini 32356 _____ \$178
- 7" Tablet 32357 _____ \$173
- 10" Tablet 32358 _____ \$395

STEP 2: Where does it need to MOUNT?

- Frame 32350 _____ \$512
- Armrest 32351 _____ \$395
Not compatible with Swing Away Armrests
- Tray 32352 _____ \$212

Cup Holder

- Heavy Duty 30889 _____ \$44
Compatible only with Flip Back Armrests