



Phone (800) 981-1540 | Fax (715) 254-0996
sales@kimobility.com | www.kimobility.com

CATALYST 5 TTL

Company:

Ship To:

Account Number:

Date:

Address:

Quote: PO Number:

City:

State:

ATP:

Zip:

Phone:

Marked For:

Email/Fax:

Associated Quote/Order#:

Effective: 4/1/2024

Catalyst 5TTL (K0005) _____ \$2696
300 lb weight capacity

Frame

Seat Width

14"	N/C
15"	N/C
16"	N/C
17"	N/C
18"	N/C
19"	N/C
20" (E2201)	\$518
21" (E2201)	\$518
22" (E2201)	\$518

Seat Depth

14"	N/C
15"	N/C
16"	N/C
17"	N/C
18"	N/C
19"	N/C
20" (E2203)	\$523

Transit Option

Transit System (Weight Limit 300lbs)	\$333
Not available with Removable Seat Pan or 350 lb Weight Capacity	
Omit	N/C

Wheelchair Drive Ready

Chair will ship from Ki Mobility with adaption kit installed for selected Alber system. Order system separately from Alber USA. Select Rear Seat Height (full profile selections only) in drive wheel size on Rear Seat Height chart (page 5). Requires Rear Anti-Tips. Not available with Spinergy wheels or Height Adjustable Flip Back T-Arm.

e-motion Ready	N/C
twion Ready	N/C

Titanium Option

Titanium Frame	\$1819
Includes hangers and base frame. Does not include Cross Braces.	

Weight Capacity

300 lbs	N/C
350 lb (Heavy Duty)	\$513
Requires minimum 6" caster height	

Cross Brace

Standard Cross Braces	N/C
Seating System Cross Braces	N/C
Allows for easy attachment of seating system mounting hardware. Not available with Seat Sling.	
Titanium Cross Braces	\$725
Not available with 21" and 22" Seat Widths. 250 lb maximum weight capacity.	

Seat Pan

Removable Aluminum Seat Pan (1/2" Shorter than Depth) (E2231)	\$300
Pan depth 1/2" shorter than Seat Depth	

Rear Frame Type

No selection needed with Titanium Frame

Standard (Curved)	N/C
Select from Standard or Angle Adjustable Back Posts on chart on the next page	
Depth Adjustable (K0108)	\$108
Select from depth adjustable back posts in chart on next page. Provides 2" forward and 2" rearward adjustment with minimum Seat Depth of 14" and maximum Seat Depth of 22". Available only with T-Arm or Flip Up Armrests.	
Reclining Frame	\$1122
Select from Reclining Back Post in chart on the next page. Not available with Dynamic Back, O2 Holder, Transit Option, 350 lb Weight Capacity, Swing Away Armrests, Flip Back Armrests, Height Adjustable Flip Back Armrests, Back Upholstery or IV Pole.	

HCPCS Codes provided are based on code verification by the PDAC or our interpretation of the code definitions.

Proper use of HCPCS Codes does not ensure coverage or payment. For coverage information, verify the appropriate payer's coverage policy.

DCN0021.41 Pricing and specifications subject to change without notice.



Phone (800) 981-1540 | Fax (715) 254-0996
sales@kimobility.com | www.kimobility.com

CATALYST 5 TTL

Backrest

Back Post / Height Range / Back Height

Depth Adjustable not available with Titanium Frame

Back Post	Price	Height Range	Back Height				
Standard	N/C	Short ^{3,4}	8.5"	9.5"	10.5"	11.5"	
		Medium ⁴	12.5" ⁶	13.5"	14.5"	15.5"	
		Tall	16.5"	17.5"	18.5"	19.5"	
Angle Adjustable w/ Push Handle	\$190	Short ⁴	14"	15"	16"		
		Medium ⁴	16"	17"	18"		
		Tall	18"	19"	20"		
Stroller Handle Back Post (Depth Adjustable)	N/C	Short	20"				
		Tall	24"				
Height Adjustable Straight with Push Handle (Depth Adjustable)	\$87	Short ¹	13"	14"	15"	16"	17"
		Medium	16"	17"	18"	19"	20"
		Tall	20" ⁸	21" ⁵	22" ⁵	23" ⁵	24" ⁵
Reclining w/ Adjustable Handle (E1225) ^{1,7}	N/C	Short	22"				
		Tall	26"				

¹ Not available with Dynamic Back

² Not available with Standard Back Upholstery

³ Not available with Tension Adjustable Back Upholstery

⁴ Not available with Standard Back Upholstery in combination with 21"-22" Seat Widths

⁵ Not available with Tension Adjustable Back Upholstery in combination with 14" and 15" Seat Widths

⁶ Not available with Tension Adjustable Back Upholstery in combination with 22" Seat Width

⁷ Not available with O2 Holder, Transit Option, 350 lb Weight Capacity, Swing Away Armrests, Flip Back Armrest, Height Adjustable Flip Back T-Arm or IV Pole

⁸ Not available with Angle Adjustable Locking Extendable Flip Up Armrest

⁹ Not available with Tension Adjustable Back Upholstery in combination with Angle Adjustable Locking Extendable Flip Up Armrests

Reclining Back Options

Select options below if Reclining Back is selected on chart above. Recline angle range is 0° - 65°.

Standard Spring	N/C
Heavy Duty Spring	N/C

Standard Back Post

Select option below only if a Standard Back Post option is selected on chart above

Straight with Push Handle	N/C
Straight without Push Handle	N/C
8 Degree Bend with Push Handle	N/C
<i>Not available with Short Back Height Range</i>	
Straight with Fold Down Push Handle	\$128
<i>Not available with 350 lb Weight Capacity</i>	

Push Handle Options

Not available with Stroller Handle Back Post, Straight Back Post with Fold Down Push Handle, Reclining Back Post or 350 Weight Capacity

Low Profile Push Handle	\$87
<i>2" shorter than standard push handles</i>	

Depth Adjustable Back Post Options

Options below may be selected only if a Depth Adjustable back post option is selected on the chart above

Stroller Handle Extensions	\$87
<i>Not available with Height Adjustable Straight with Push Handle Back Post</i>	
Dynamic Back (E2398)	\$1395
<i>Available only with T-Arms.</i>	

Back Upholstery Options

Not available with Reclining or Stroller Handle back post. If selecting both upholstery and Axiom Back, the upholstery will be installed.

Standard	N/C
Tension Adjustable with Standard Cover (E2611)	\$330
Tension Adjustable with Air Mesh Breathable Cover (E2611)	\$400
Omit	(\$10)

U.S.

Attached only if ordered with Tension Adjustable Back Upholstery

Army	Navy	Marines	Air Force	
Coast Guard	US Purple Heart			\$92
Combo Army/Purple Heart		Combo Navy/Purple Heart		
Combo Marines/Purple Heart		Combo Air Force/Purple Heart		
Combo Coast Guard/Purple Heart				\$149

Frame Type / Footrest

Frame Type

Swing Away	N/C
Super Low - Swing Away	N/C
70° Fixed Front	N/C
80° Fixed Front	N/C

Hanger / Front Frame Bend

Frame	Hanger	60°	70°	80°	90°		Retail
Swing Away	Extension Tube Mount	2		1			N/C
	Front Mount 2			1	1		N/C
	Pro Elevating Leg Rests E0990						\$302
	Omit Hangers						N/C
Fixed	Extension Tube Mount			1			N/C

1 Not available with 7" or 8" casters

2 Not available with Titanium Frame

PRO ELR Cover

Available only on PRO ELR

Gel Pad Covers (PAIR) (K0108) \$400

Footrest / Extension / Seat to Footrest Length

Select Extension and Seat to Footrest Length in row corresponding to desired Footrest. Footrest Length of at least 2.5" shorter than Front Seat Height recommended. Extension Tube Hangers and Fixed Front Frame measurements are based on 70° Hanger / Front Frame Bend.

Footrest	Extension	Seat to Footrest Length									
		11"	11.5"	12"	12.5"	13"	13.5"	14"	14.5"	15"	
Extension Tube Hangers (60°, 70°, 80°)	Short (11" - 15")	11"	11.5"	12"	12.5"	13"	13.5"	14"	14.5"	15"	
	Medium (13" - 17")	13"	13.5"	14"	14.5"	15"	15.5"	16"	16.5"	17"	
	Long (15.5" - 19.5")	15.5"	16"	16.5"	17"	17.5"	18"	18.5"	19"	19.5"	
Front Mount Hangers	Standard Selections										Requires 3" Extension
	Short Front Mount Hanger	4" 2	4.5" 2	5"	5.5"	6"	6.5"	7"			10" 10.5" 11"
		7.5"	8"	8.5"	9"	9.5"					11.5" 12" 12.5"
	Medium Front Mount Hanger	4" 2	4.5" 2	5"	5.5"	6"	6.5"	7"	7.5"	8"	12" 12.5" 13"
		8.5"	9"	9.5"	10"	10.5"	11"	11.5"			13.5" 14" 14.5"
	Long Front Mount Hanger	5"	5.5"	6"	6.5"	7"	7.5"	8"	8.5"	9"	14" 14.5" 15"
		9.5"	10"	10.5"	11"	11.5"	12"	12.5"	13"	13.5"	15.5" 16" 16.5"
PRO ELR Low Pivot Setting	PRO ELR Low Pivot Setting Short (12" - 17.5")	12" - 17.5"									
	PRO ELR Low Pivot Setting Long (15" - 21")	15" - 21"									
PRO ELR High Pivot Setting	PRO ELR High Pivot Setting Short (10.5" - 15.5")	10.5" - 15.5"									
	PRO ELR High Pivot Setting Long (13" - 19")	13" - 19"									
Fixed Front Frame (70°, 80°)	Short (11" - 15")	11.5"	12"	12.5"	13"	13.5"	14"	14.5"	15"		
	Medium (13" - 17")	13"	13.5"	14"	14.5"	15"	15.5"	16"	16.5"	17"	
	Long (15.5" - 19.5")	15.5"	16"	16.5"	17"	17.5"	18"	18.5"	19"	19.5"	

1 Not available with Titanium Frame. See Minimum Front Seat Height chart above.

2 Available with 90° only

3" Extension

Available only with Front Mount Hangers

3" Extension	N/C
3" Extension (Provided as Parts)	N/C

Hanger Release

If no selection is made Classic Hanger Release will be used

Classic	N/C
4-Way - In-Line	N/C
4-Way - Out	N/C

Footplates

Composite	N/C
Not available with Front Mount Hanger	
Composite Angle Adjustable (K0040)	\$195
Not available with Front Mount Hanger or 14" Seat Width	
Aluminum Angle Adjustable (K0040)	\$200
Aluminum Locking Multi-Angle Adjustable Footplate (K0108)	\$308
Available only with Front Mount Hanger	
One Piece Flip Up Angle Adjustable (K0040)	\$195
Not available with Front Mount Hanger or PRO ELR	



Phone (800) 981-1540 | Fax (715) 254-0996
sales@kimobility.com | www.kimobility.com

CATALYST 5 TTL

Heel Loops

Not available with One Piece Flip Up Angle Adjustable Footplate

Heel Loops (E0951)	\$42
Adjustable Heel Loops (E0951)	\$51
Omit Heel Loops	N/C

Calf Straps

Velcro Adjustable (K0038)	\$41
Padded Velcro Adjustable (K0038)	\$70
Bodypoint Padded AeroMesh (K0038)	\$121

Residual Limb Support

	Limb Support	Price	Omit	Cover		Pad Size					
Left	Left Limb Support (E1020)	\$415	Omit Left Hanger	Smooth Cover Left	Air Mesh Cover Left	5x8 Pad Left Side	7x10 Pad Left Side	9x10 Pad Left Side	9x14 Pad Left Side	11x10 Pad Left Side	11x14 Pad Left Side
Right	Right Limb Support (E1020)	\$415	Omit Right Hanger	Smooth Cover Right	Air Mesh Cover Right	5x8 Pad Right Side	7x10 Pad Right Side	9x10 Pad Right Side	9x14 Pad Right Side	11x10 Pad Right Side	11x14 Pad Right Side

Front Seat Height and Casters

Caster Size / Front Frame Type / Front Seat Height

Select Front Seat Height in row corresponding to desired caster size. Note: Selections in first three Front Seat Height columns use Super Low Front Frame. Super Low Front Frame not available with Titanium Frame or Fixed Front Frame.

Caster Size	Front Seat Height													
	Uses Super Low Front Frame													
4"	13.5"	14"	14.5"	15"	15.5"	16"	16.5"	17"	17.5"	18"	18.5"	19"	19.5"	20"
5"	14.5"	15"	15.5"	16"	16.5"	17"	17.5"	18"	18.5"	19"	19.5"	20"	20.5"	
6"	15.5"	16"	16.5"	17"	17.5"	18"	18.5"	19"	19.5"	20"	20.5"	21"		
7"	17"	17.5"		18"	18.5"	19"	19.5"	20"	20.5"	21"				
8"	17.5"	18"	18.5"	19"	19.5"	20"	20.5"	21"						

Maximum Seat Dump Reference

(Front Seat Height minus Rear Seat Height)

Seat Depth	14"	15"	16"	17"	18"	19"	20"
Maximum Dump Standard Caster Housing	1.5"	2"	2"	2"	2.5"	2.5"	2.5"
Maximum Dump Performance Caster Housing	3"	3"	3"	3.5"	3.5"	4"	4"

Depending on configuration, chairs with Fixed Front Frame Type may achieve up to an additional .75"

Casters

		Caster Size					
		Price	4"	5"	6"	7"	8"
Caster Type	.75" Rollerblade	N/C					
	.75" Lighted Rollerblade	\$26					
	1" Poly	N/C					
	1" Poly Aluminum	\$113					
	1" Pneumatic	\$87					
	1.5" Poly	N/C					
	1.5" Soft Roll Aluminum	\$67					
	2" Poly	N/C					
	2" Pneumatic w/Foam Insert	\$128					

Caster Pin Locks

Caster Pin Locks (K0073)	\$200
<i>Not available with Performance Caster Housing, 1 1/2" Longer Fork Stem or Frog Legs Fork</i>	

Caster Housing

Standard Housing	N/C
Performance Housing	\$36
<i>Not available with 350 lb Weight Capacity</i>	

Fork / Stem

Best Fit	N/C
<i>Ki Mobility will pick Fork and Stem combination to optimize performance based on Caster Size and Front Seat Height</i>	
Standard Fork Stem	N/C
3/4" Longer Fork Stem	N/C
1 1/2" Longer Fork Stem	N/C
Frog Legs II Suspension System	\$461
<i>Lightweight shock absorbing fork. Available only with 4", 5" or 6" Caster Sizes. Not available with Caster Pin Locks, 350 lb Weight Capacity, 2" Poly caster or 2" Pneumatic w/ Foam Insert caster.</i>	
Frog Legs Suspension Fork	\$461
<i>Available only with 4", 5" or 6" Caster Sizes. Not available with Caster Pin Locks, 350 lb Weight Capacity, 2" Poly caster or 2" Pneumatic w/ Foam Insert caster.</i>	

Frog Legs Option

User Weight Under 80 lbs	N/C
User Weight 80-120 lbs	N/C
User Weight 120-160 lbs	N/C
User Weight 160-220 lbs	N/C
User Weight 220-260 lbs	N/C
User Weight 260-300 lbs	N/C

Wheel Locks

Wheel Lock	
Push to Lock	N/C
Push to Lock Flush Mount	N/C
Pull to Lock	N/C
Low Profile Scissor Lock	N/C
Not available with Low Poly, Full Poly or Sentinel Solid tires	
Drum Brake (E2228)	\$574
Available only with 20", 22" and 24" MAXX Spoke wheels. Not available with Reclining Backrest, Straight without Push Handle or Straight with Fold Down Push Handle Back Posts, or Low Profile Push Handle. Axle is Quick Release.	
Dual - Drum Brake with Push to Lock (E2228)	\$677
Available only with 20", 22" and 24" MAXX Spoke wheels. Not available with Reclining Backrest, Straight without Push Handle or Straight with Fold Down Push Handle Back Posts, or Low Profile Push Handle. Axle is Quick Release.	
Dual - Drum Brake with Pull to Lock (E2228)	\$677
Available only with 20", 22" and 24" MAXX Spoke wheels. Not available with Reclining Backrest, Straight without Push Handle or Straight with Fold Down Push Handle Back Posts, or Low Profile Push Handle. Axle is Quick Release.	
Omit	N/C

Wheel Lock Options	
Extension Handles (E0961)	\$77
Available only with Push to Lock or Pull to Lock Wheel Locks	
Grade Aids (E0974)	\$154
Available only with Push to Lock Wheel Locks and Pneumatic or Pneumatic with Airless Insert tires	

Rear Seat Height and Wheels

Wheel Size / Tire Profile / Rear Seat Height											
<i>Select Rear Seat Height in row corresponding to Wheel Size and Tire Profile. 14" Seat Depth not available with 24", 25" or 26" wheel sizes.</i>											
Wheel	Tire	Rear Seat Height									
20"	Low	11.5" ³	12" ¹	12.5" ³	13" ¹	13.5" ³	14"	14.5" ³	15"	15.5" ³	16" ³
	Full	12" ³	12.5" ¹	13" ³	13.5" ¹	14" ³	14.5"	15" ³	15.5"	16" ³	16.5" ³
22"	Low	12.5" ³	13" ¹	13.5" ³	14" ¹	14.5" ³	15"	15.5" ³	16"	16.5" ³	17" ³
	Full	13" ³	13.5" ¹	14" ³	14.5" ¹	15" ³	15.5"	16" ³	16.5"	17" ³	17.5" ³
24"	Low	13.5" ³	14" ¹	14.5" ³	15" ¹	15.5" ³	16"	16.5" ³	17"	17.5" ³	18" ³
	Full	14" ³	14.5" ¹	15" ³	15.5" ¹	16" ³	16.5"	17" ³	17.5"	18" ³	18.5" ³
25"	Low	14" ³	14.5" ¹	15" ³	15.5" ¹	16" ³	16.5"	17" ³	17.5"	18" ³	18.5" ³
26"	Low		15" ^{1,2}	15.5" ³	16" ¹	16.5" ³	17"	17.5" ³	18"	18.5" ³	19" ³
	Full	15" ^{2,3}	15.5" ^{1,2}	16" ³	16.5" ¹	17" ³	17.5"	18" ³	18.5"	19" ³	19.5" ³

1 Not available with Amputee Axle Plate

2 Not available with Height Adjustable Flip Back T-Arm

3 Requires Universal Axle Plate

One Arm Drive		Wheels					
One Arm Drive - Quick Release (E0958)		<i>No selection needed with One Arm Drive</i>					
Available only with 20", 22" and 24" spoke wheels and Quick Release Axles. Spoke wheel style is exclusive to OAD. Not available with 4 degrees of camber. Standard handrim will be Aluminum Anodized.							
\$1122							
Drive Side							
Left Side Drive							
Right Side Drive							
Drive Handrim							
Aluminum Drive Handrim							
Poly Coated Drive Handrim							
\$205							
Spinergy Options							
<i>Select Spoke Color and Handrim Tab for Spinergy wheels</i>							
Black Spokes		Red Spokes		Blue Spokes		White Spokes	N/C

Wheels

No selection needed with One Arm Drive

		Rear Wheel Size					
		Price	20"	22"	24"	25"	26"
Rear Wheel Type	Maxx Spoke	N/C					
	Mag	N/C					
	Superlight Spoke ¹	\$302				²	
	Spinergy Spox	\$866					
	Spinergy LX	\$1076					

1 Not available with 350 lb Weight Capacity

2 Available only with Superlight and Plastic Coated Handrims

Axle Plates

No selection needed with One Arm Drive

Universal	N/C
Amputee (E0959)	\$179
<i>Not available with Swing Away Armrests or Wheelchair Drive Ready Option</i>	

Handrimms

Standard Handrim

Aluminum Anodized	N/C
Superlight Aluminum Anodized	\$51
<i>Available only on 24" and 25" wheel sizes. Available only on MAXX Spoke, Superlight or Spinergy wheels.</i>	
Plastic Coated	\$51
Projection (8 Vertical Projections)	\$333
<i>Available only with full profile tires in combination with 22", 24", and 26" MAXX Spoke wheels, or 22" and 24" Mag wheels</i>	

Ergonomic Handrim

Not available with 20" or 26" wheels. Tab mount only available with 24" wheels.

Lite Regular Grip	\$477
Standard Regular Grip	\$456
Lite Extra Grip	\$528
Standard Extra Grip	\$507
Lite Omit Grip	\$405
Standard Omit Grip	\$354

Handrim Options

3M Micro Grip	\$26
<i>Not available on Plastic Coated or Projection Handrimms</i>	

Handrim Mount

Close Mount	\$26
<i>If no selection is made, Standard Mount will be used</i>	

Rear Tire Type

Full Profile

Street Pneumatic 1 3/8" (E2211, E2212)	N/C
<i>Not available on 25" wheel size</i>	
Full Poly 1 3/8" (E2218)	N/C
<i>Not available on Superlight wheel or 25" or 26" wheel sizes</i>	
Pneumatic with Airless Insert 1 3/8" (E2211, E2213)	\$87
<i>Not available on Superlight wheel or 25" wheel size</i>	
Schwalbe HS-127 1 3/8" (E2211, E2212)	\$128
<i>Available only on 24" wheel size. Not available on Mag wheel.</i>	

Low Profile

Low Poly 1" (E2218)	N/C
<i>Not available on Superlight wheel, 25" or 26" wheel sizes or 350 lb Weight Cap.</i>	
Sentinel Solid Tire 1" (E2218)	\$26
<i>Not available on Superlight wheels or 20" and 22" Mag wheels</i>	
High Pressure 1" (E2211, E2212)	N/C
<i>Not available on Mag wheel</i>	
Schwalbe One 1" (E2211, E2212)	\$118
<i>Not available on Mag wheel, 20", 22" or 26" wheel sizes</i>	
Schwalbe Marathon Plus Evolution 1" (E2211, E2212)	\$128
<i>Not available on Mag wheel, 20" or 22" wheel sizes</i>	

Axles

No selection needed with One Arm Drive

Quick Release	N/C
Quad Quick Release	\$128
<i>Not available with Mag or Spinergy wheels</i>	

Axle Position

Rear Wheel Camber

0° 2° 4°	N/C
----------	-----

No selection needed with Amputee Axle Plate

Rear Center Forward	N/C
---------------------	-----

Armrests

Armrest Type	Length	Price
Height Adjustable T-Arm Tall (11.5"-15") E0973	Desk	\$318
	Full	
Height Adjustable T-Arm (8.5"-13.5") E0973	Desk	\$318
	Full	
Height Adjustable T-Arm Low (6.5"-10.5") E0973 ¹	Desk	\$318
	Full	
Height Adjustable Flip Back T-Arm Tall (11.5"-15") E0973 ²	Desk	\$333
	Full	
Height Adjustable Flip Back T-Arm (9.5"-13.5") E0973 ²	Desk	\$333
	Full	

¹ If requested configuration is unavailable, Standard Height T-Arm will be used

² Available only with Curved Rear Frame

³ Not available with Seat Depths less than 16". Not Available with 15" Rear Seat Height and 24" Wheels.

⁴ Available only with Standard Back Posts

⁵ Available only with Depth Adjustable and Reclining Back Posts

⁶ Not available with Arm Pads

Armrest Type	Length	Price
Flip back (10.5") ^{3,4}	Desk	N/C
	Full	
Angle Adjustable Locking Extendable Flip Up E0973 ⁵	Desk	\$272
	Full	
Swing Away (10.5"-12.5") E0973 ^{4,6}		\$231
Swing Away with Side Guard (10.5"-12.5") E0973 ^{4,6}		\$302
Omit Armrests		N/C

Arm Pads

Not available with Swing Away Armrests

Standard Arm Pad	N/C
Waterfall Arm Pad	\$26
<i>Available only with Height Adjustable T-Arm Tall, Height Adjustable T-Arm, and Height Adjustable T-Arm Low</i>	
Foam Grip	N/C
<i>Available only with Angle Adjustable Locking Extendable Flip Up Armrest</i>	



Phone (800) 981-1540 | Fax (715) 254-0996
sales@kimobility.com | www.kimobility.com

CATALYST 5 TTL

Accessories

Removable Side Guard

Not available with T-Arms, Flip Back or Swing Away with Side Guard Armrests

Adult Composite Side Guards (K0108)	\$144
Aluminum Straight - Black Matte (K0108)	\$200
Aluminum Straight - Match Frame Color (K0108)	\$251
<i>Not available with Hydrographics or Titanium Frame</i>	
Aluminum Adjustable Fender w/Foam Pad - Black Matte (K0108)	\$282
Aluminum Adjustable Fender w/Foam Pad - Match Frame Color (K0108)	\$333
<i>Not available with Hydrographics or Titanium Frame</i>	
Carbon Fiber Adjustable Fender (K0108)	\$497

Anti-Tips

Rear Anti-Tips (E0971)	\$149
Swing-In Anti-Tip Single - Left (E0971)	\$159
<i>Not available with 350 lb Weight Capacity, Reclining Back Post or Depth Adjustable Rear Frame. Not available with left side O2 Holder.</i>	
Swing-In Anti-Tip Single - Right (E0971)	\$159
<i>Not available with 350 lb Weight Capacity, Reclining Back Post or Depth Adjustable Rear Frame. Not available with right side O2 Holder.</i>	
Swing-In Anti-Tip Pair (E0971)	\$318
<i>Not available with 350 lb Weight Capacity, Reclining Back Post or Depth Adjustable Rear Frame. Not available with Cane and Crutch Holder or O2 Holder.</i>	
Omit	N/C

Pelvic Positioning Belts

1 1/2" Auto Buckle Pelvic Positioning Belt (E0978)	\$87
1 1/2" Padded Auto Buckle Pelvic Positioning Belt (E0978)	\$179
2" Airline Buckle Pelvic Positioning Belt (E0978)	\$118
Small Bodypoint Evoflex (16-24") (E0978)	\$200
Medium Bodypoint Evoflex (19-29") (E0978)	\$200
Large Bodypoint Evoflex (25-37") (E0978)	\$200
Bodypoint Belt Mounting Bracket Kit	\$113
Bodypoint 1" Band Clamp Belt Mounting Kit	\$113

O2 Holder

Not available with Depth Adjustable Rear Frame

O2 Holder - Left Side (E2208)	\$313
O2 Holder - Right Side (E2208)	\$313

Other Accessories

Tool Kit	\$46
Removable Underseat Pouch	\$51
Pediatric Ki Mobility Backpack	\$51
Ki Mobility Backpack	\$80
Luggage Carrier	\$87
<i>Not available with PRO ELR or Residual Limb Support</i>	
Black Neoprene Impact Guards	\$97
Clear Spoke Guards (K0065)	\$123
<i>Not available with Mag wheel</i>	
Black Spoke Guards (K0065)	\$123
<i>Not available with 20" wheel size or Mag wheel</i>	
Cane and Crutch Holder	\$210
IV Pole	\$292
<i>Available only with Curved Rear Frame Type and Straight Back Posts</i>	

Frame Color

Gloss Finish

Black	N/C
Blue Skies	N/C
Bubblegum Pink	N/C
Candy Blue	N/C
Candy Purple	N/C
Candy Red	N/C
Copper	N/C
Granny Smith	N/C
Metallic Orange	N/C
Yellow	N/C

Matte Finish

Military Green	N/C
Charcoal Gray	N/C
Emerald Green	N/C
Iridescent Burgundy	N/C
Iridescent Navy	N/C
Retro Red	N/C
Royal Blue	N/C
Shadow Black	N/C
Smooth White	N/C

Hydrographics

1 Select Hydrographics Base Color. Zebra not available in Black.

AOR1 (Desert Digital Camo)	NEXT Camo	
Pink Camo	Carbon Fiber ¹	
Muddy Girl	Zebra ¹	\$794

Hydrographics Base Color

Black Blue Granny Smith Pink Purple Red N/C

Colored Anodizing Package

	Option	Standard Caster Housing ²	Caster Fork ³	Axle Plate ⁴
	Retail	\$77	\$77	\$77
Color	Black ¹			
	Blue			
	Bright Green			
	Orange			
	Pewter			
	Pink			
	Purple			
	Red			

1 Black is no charge

2 Available only if corresponding option selection was made

3 Not available with 2" wide casters, Frog Legs Suspension Fork or Frog Legs II Suspension System

Catalyst Front Seat Height Matrix																			
4" Caster				5" Caster				6" Caster				7" Caster				8" Caster			
FSH	Super Low	Stem	Fork	FSH	Super Low	Stem	Fork	FSH	Super Low	Stem	Fork	FSH	Super Low	Stem	Fork	FSH	Super Low	Stem	Fork
15	13.5	Std	4"	16	14.5	Std	4"	17	15.5	Std	5"	18	16.5	Std	6"	19	17.5	Std	6"
15.5	14			16.5	15.5	3/4		17.5	16.5	3/4		18.5	17			19.5	18.5	3/4	
15.5 ¹	14 ¹	17.5		16	1 1/2	18.5		17	1 1/2	18.5		17.5	20.5	n/a		1 1/2			
16	14.5	3/4		16	14.5	Std	5"	17	15.5	Std	19	18	3/4	19		17.5	Std	7"	
16.5 ¹	15 ¹			16.5	15			17.5	16		19.5	18							
17	15.5	16.5		15.5	18	16.5		20	18.5		20	18.5							
15.5	14	Std	17	16	3/4	17.5		16.5	Std	18.5	17	1 1/2	20.5	19	3/4				
16	14.5		17.5	16		19		17.5		19.5	18.5								
16	15	3/4	18	16.5	1 1/2	18.5		17.5	Std	19.5	18	Std	20	19		3/4			
16.5	15.5		16.5	15		20	18.5	20.5		n/a									
17	15.5	1 1/2	17	15.5	Std	19	17.5	1 1/2	19	18	3/4	21	n/a	1 1/2					
17.5	16		17.5	16		19.5	18		20.5	19									
16.5	15	Std	17.5	16	3/4	18	16.5	Std	20	19	1 1/2								
17	15.5		18	16.5		18.5	17		20.5	n/a									
17	16	3/4	18.5	17	1 1/2	19	17.5	3/4	20	18.5	1 1/2								
17.5	16.5		18	16.5		19.5	18		20.5	19									
18	16.5	1 1/2	18.5	17	Std	18.5	17.5	3/4	21	n/a									
18.5	17		19	17.5		19	18												
18	16.5	Std	18	16.5	Std	19.5	18.5	1 1/2											
18.5	17		18.5	17		20	19												
18.5	17.5	3/4	19	17.5	3/4	19.5	18	1 1/2											
19	18		18.5	17.5		20	18.5												
19.5	18.5	1 1/2	19	18	3/4	20.5	19	1 1/2											
20	19		19.5	18.5		21	n/a												
				19.5	18	6" x 2" caster not available with 5" fork													
				20	18.5												1 1/2		
				20.5	19														

¹ Not available with 4 x 1.5" Soft Roll Aluminum

Seat heights with "FSH" heading available with Standard Front Frame. Seat heights with "Super Low" heading available with Super Low Front Frame. Actual seat height may vary depending upon configuration.