



Phone (800) 981-1540 | Fax (715) 254-0996
 sales@kimobility.com | www.kimobility.com

CATALYST 5 TTL

Effective 05/23/2017

Company: _____

Ship To: _____

Account Number: _____ Date: _____

Address: _____

Quote PO Number: _____

City: _____ State: _____

RTS: _____

Zip: _____ Phone: _____

Marked For: _____

Email/Fax: _____

Catalyst 5TTL K0005 **\$2475**
 300 lb weight capacity

HCPCS codes provided are based on code verification by the PDAC or our interpretation of the code definitions.

Proper use of HCPCS codes does not ensure coverage or payment. For coverage information, verify the appropriate payer's coverage policy.

Frame

Seat Width / Seat Depth

Select box corresponding to desired Seat Width and Depth

		Seat Width									
		14"	15"	16"	17"	18"	19"	20"	21"	22"	
Seat Depth	14"	N/C	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>
	15"	N/C	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>
	16"	N/C	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>
	17"	N/C	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>
	18"	N/C	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>
	19"	N/C	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>
	20"	\$450 E2203	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>

Transit Option

Transit System K0108 \$320
 Not available with Removable Seat Pan or Heavy Duty Option

Wheelchair Drive Ready

Chair will ship from Ki Mobility with adaption kit installed for selected Alber system. Order system separately from Alber USA. Select Rear Seat Height (full profile selections only) in drive wheel size on Rear Seat Height chart (page 5). Requires Rear Anti-Tips. Not available with Spinergy wheels.

e-motion Ready _____ N/C
 twion Ready _____ N/C

Titanium Option

Titanium Frame K0108 \$1650
 Includes hangers and base frame. Does not include Cross Braces.

Heavy Duty Option

350 lb Weight Capacity \$500
 Requires minimum 6" caster height

Cross Brace

Standard Cross Braces _____ Std
 Seating System Cross Braces _____ N/C
Allows for easy attachment of seating system mounting hardware. Not available with Seat Sling.
 Titanium Cross Braces \$600
 Not available with 21" and 22" Seat Widths. 250 lb maximum weight capacity.

Rear Frame Type

No selection needed with Titanium Frame

Curved _____ Std
 Select from Standard or Angle Adjustable Back Posts on chart on the next page
 Depth Adjustable K0108 \$105
 Select from Standard Depth Adjustable or Height Adjustable Depth Adjustable Back Posts in chart on the next page. Provides 2" forward and 2" rearward adjustment with minimum Seat Depth of 14" and maximum Seat Depth of 22". Available only with T-Arm or Flip Up Armrests. Not available with Back Upholstery.
 Reclining \$1095
 Select from Reclining Back Post in chart on the next page. Not available with Seating Dynamics Rocker Back, O₂ Holder, Transit Option, Heavy Duty Option, Swing Away Armrests, Flip Back Armrests, Height Adjustable Flip Back Armrests, Back Upholstery or IV Pole.

Backrest

Back Post / Height Range / Back Height

Select Back Height in row corresponding to desired Back Post and Height Range. Depth Adjustable not available with Titanium Frame.

Back Post	Price	Height Range	Back Height				
Standard	N/C	Short	<input type="checkbox"/> 8.5	<input type="checkbox"/> 9.5	<input type="checkbox"/> 10.5	<input type="checkbox"/> 11.5	
		Medium	<input type="checkbox"/> 12.5	<input type="checkbox"/> 13.5	<input type="checkbox"/> 14.5	<input type="checkbox"/> 15.5	
		Tall	<input type="checkbox"/> 16.5	<input type="checkbox"/> 17.5	<input type="checkbox"/> 18.5	<input type="checkbox"/> 19.5	
Angle Adjustable w/ Push Handle	\$185	Short	<input type="checkbox"/> 14	<input type="checkbox"/> 15	<input type="checkbox"/> 16		
		Medium	<input type="checkbox"/> 16	<input type="checkbox"/> 17	<input type="checkbox"/> 18		
		Tall	<input type="checkbox"/> 18	<input type="checkbox"/> 19	<input type="checkbox"/> 20		
Standard Depth Adjustable	N/C	Short	<input type="checkbox"/> 20				
		Tall	<input type="checkbox"/> 24				
Height Adjustable Depth Adjustable	\$75	Short ¹	<input type="checkbox"/> 13	<input type="checkbox"/> 14	<input type="checkbox"/> 15	<input type="checkbox"/> 16	<input type="checkbox"/> 17
		Medium	<input type="checkbox"/> 16	<input type="checkbox"/> 17	<input type="checkbox"/> 18	<input type="checkbox"/> 19	<input type="checkbox"/> 20
		Tall	<input type="checkbox"/> 20	<input type="checkbox"/> 21	<input type="checkbox"/> 22	<input type="checkbox"/> 23	<input type="checkbox"/> 24
Reclining w/ Adjustable Handle ^{1, 2} <small>(Removable)</small>	N/C	Short	<input type="checkbox"/> 22				
		Tall	<input type="checkbox"/> 26				

¹ Not available with Seating Dynamics Rocker Back
² Not available with O₂ Holder, Transit Option, Heavy Duty Option, Swing Away Armrests, Flip Back Armrests, Height Adjustable Flip Back Armrests, Back Upholstery or IV Pole

Reclining Back Post Option

Select options below if Reclining Back is selected on chart above. Recline angle range is 0° - 65°.

- Gas Spring: Standard Spring
 Heavy Duty Spring

Standard Back Post

Select option below only if a Standard Back Post option is selected on chart at left

- Straight with Push Handle _____ Std
 Straight without Push Handle _____ N/C
 8° Bend with Push Handle _____ N/C
Not available with Short Back Height Range
 Straight with Fold Down Push Handle _____ \$125
Not available with Heavy Duty option

Push Handle Options

Not available with Standard Depth Adjustable Back Post, Straight Back Post with Fold Down Push Handle, Reclining Back Post or Heavy Duty Option

- Low Profile Push Handle _____ \$65
2" shorter than standard push handles

Depth Adjustable Back Post Options

Options below may be selected only if Standard Depth Adjustable or Height Adjustable Depth Adjustable option is selected on chart at left

- Stroller Handle Extensions _____ N/C
Not available with Height Adjustable Depth Adjustable Back Posts
 Seating Dynamics Rocker Back K0108 _____ \$1995
Available only with T-Arms

Back Upholstery

Select option below only if Standard or Angle Adjustable Back Post option is selected on chart at left

- Standard _____ Std
Not available with Short or Medium Back Height Ranges on 20" - 22" Seat Widths
 Tension Adjustable E2611 _____ \$245
 Omit Back Upholstery _____ - \$10

Military Patch

Attached only if ordered with Tension Adjustable Back Upholstery

- US Military Service Patch _____ \$90
 Army Navy Marines Air Force Coast Guard
 US Purple Heart Patch _____ \$90
 Purple Heart with Military Service Patch Combo _____ \$145
 Army Navy Marines Air Force Coast Guard

Footrest

Hanger / Front Frame Bend

- 60° _____ N/C
Available only with Extension Tube Hangers. Not available with Titanium Frame.
 70° _____ N/C
 80° _____ N/C
Not available with 7" or 8" Casters
 90° _____ N/C
Available only with Front Mount Hangers. Not available with Titanium Frame or 7" and 8" Casters.
 Elevating Leg Rests (ELR) E0990 _____ \$295
 Pro Elevating Leg Rests (PRO ELR) E0990 **New!** _____ \$350
 Omit Hangers _____ N/C

Footplate

- Composite _____ N/C
Not available with Front Mount Hangers or PRO ELR
 Composite Angle Adjustable K0040 _____ \$185
Not available with Front Mount Hangers, ELR or 14" Seat Width
 Aluminum Angle Adjustable K0040 _____ \$185
 One Piece Flip Up Angle Adjustable K0040 _____ \$190
Not available with Front Mount Hangers, ELR or PRO ELR

Footrest continued

Footrest / Extension / Seat to Footrest Length

Select Seat to Footrest Length in row corresponding to desired Footrest and Extension. Footrest Length of at least 2.5" shorter than Front Seat Height recommended. Extension Tube Hangers and Fixed Front Frame measurements are based on 70° Hanger / Front Frame Bend.

Footrest	Extension	Seat to Footrest Length													
Extension Tube Hangers (60°, 70°, 80°)	Short	<input type="checkbox"/> 11	<input type="checkbox"/> 11.5	<input type="checkbox"/> 12	<input type="checkbox"/> 12.5	<input type="checkbox"/> 13	<input type="checkbox"/> 13.5	<input type="checkbox"/> 14	<input type="checkbox"/> 14.5	<input type="checkbox"/> 15					
	Medium	<input type="checkbox"/> 13	<input type="checkbox"/> 13.5	<input type="checkbox"/> 14	<input type="checkbox"/> 14.5	<input type="checkbox"/> 15	<input type="checkbox"/> 15.5	<input type="checkbox"/> 16	<input type="checkbox"/> 16.5	<input type="checkbox"/> 17					
	Long	<input type="checkbox"/> 15.5	<input type="checkbox"/> 16	<input type="checkbox"/> 16.5	<input type="checkbox"/> 17	<input type="checkbox"/> 17.5	<input type="checkbox"/> 18	<input type="checkbox"/> 18.5	<input type="checkbox"/> 19	<input type="checkbox"/> 19.5					
Front Mount Hangers ¹	Standard Selections											Requires 3" Extension for N/C			
	Short (70°, 80°, 90°)	<input type="checkbox"/> 5	<input type="checkbox"/> 5.5	<input type="checkbox"/> 6	<input type="checkbox"/> 6.5	<input type="checkbox"/> 7						<input type="checkbox"/> 10	<input type="checkbox"/> 10.5	<input type="checkbox"/> 11	
		<input type="checkbox"/> 7.5	<input type="checkbox"/> 8	<input type="checkbox"/> 8.5	<input type="checkbox"/> 9	<input type="checkbox"/> 9.5						<input type="checkbox"/> 11.5	<input type="checkbox"/> 12	<input type="checkbox"/> 12.5	
	Medium (70°, 80°, 90°)	<input type="checkbox"/> 5	<input type="checkbox"/> 5.5	<input type="checkbox"/> 6	<input type="checkbox"/> 6.5	<input type="checkbox"/> 7	<input type="checkbox"/> 7.5	<input type="checkbox"/> 8					<input type="checkbox"/> 12	<input type="checkbox"/> 12.5	<input type="checkbox"/> 13
		<input type="checkbox"/> 8.5	<input type="checkbox"/> 9	<input type="checkbox"/> 9.5	<input type="checkbox"/> 10	<input type="checkbox"/> 10.5	<input type="checkbox"/> 11	<input type="checkbox"/> 11.5					<input type="checkbox"/> 13.5	<input type="checkbox"/> 14	<input type="checkbox"/> 14.5
	Long (70°, 80°)	<input type="checkbox"/> 5	<input type="checkbox"/> 5.5	<input type="checkbox"/> 6	<input type="checkbox"/> 6.5	<input type="checkbox"/> 7	<input type="checkbox"/> 7.5	<input type="checkbox"/> 8	<input type="checkbox"/> 8.5	<input type="checkbox"/> 9				<input type="checkbox"/> 14	<input type="checkbox"/> 14.5
<input type="checkbox"/> 9.5		<input type="checkbox"/> 10	<input type="checkbox"/> 10.5	<input type="checkbox"/> 11	<input type="checkbox"/> 11.5	<input type="checkbox"/> 12	<input type="checkbox"/> 12.5	<input type="checkbox"/> 13	<input type="checkbox"/> 13.5	<input type="checkbox"/> 15.5	<input type="checkbox"/> 16	<input type="checkbox"/> 16.5			
Elevating Leg Rests (ELR)	Short	<input type="checkbox"/> 12	<input type="checkbox"/> 13	<input type="checkbox"/> 14	<input type="checkbox"/> 15	<input type="checkbox"/> 16									
	Long	<input type="checkbox"/> 13	<input type="checkbox"/> 14	<input type="checkbox"/> 15	<input type="checkbox"/> 16	<input type="checkbox"/> 17	<input type="checkbox"/> 18								
Extended Length ELR	Short	<input type="checkbox"/> 18	<input type="checkbox"/> 19	<input type="checkbox"/> 20	<input type="checkbox"/> 21	<input type="checkbox"/> 22									
	Long	<input type="checkbox"/> 18	<input type="checkbox"/> 19	<input type="checkbox"/> 20	<input type="checkbox"/> 21	<input type="checkbox"/> 22	<input type="checkbox"/> 23	<input type="checkbox"/> 24							
PRO ELR Low Pivot	Short	<input type="checkbox"/> 12 - 17.5													
	Long	<input type="checkbox"/> 15 - 21													
PRO ELR High Pivot	Short	<input type="checkbox"/> 10.5 - 15.5													
	Long	<input type="checkbox"/> 13 - 19													
Fixed Front Frame (70°, 80°)	Short	<input type="checkbox"/> 11.5	<input type="checkbox"/> 12	<input type="checkbox"/> 12.5	<input type="checkbox"/> 13	<input type="checkbox"/> 13.5	<input type="checkbox"/> 14	<input type="checkbox"/> 14.5	<input type="checkbox"/> 15						
	Medium	<input type="checkbox"/> 13	<input type="checkbox"/> 13.5	<input type="checkbox"/> 14	<input type="checkbox"/> 14.5	<input type="checkbox"/> 15	<input type="checkbox"/> 15.5	<input type="checkbox"/> 16	<input type="checkbox"/> 16.5	<input type="checkbox"/> 17					
	Long	<input type="checkbox"/> 15.5	<input type="checkbox"/> 16	<input type="checkbox"/> 16.5	<input type="checkbox"/> 17	<input type="checkbox"/> 17.5	<input type="checkbox"/> 18	<input type="checkbox"/> 18.5	<input type="checkbox"/> 19	<input type="checkbox"/> 19.5					


¹ Available only with Aluminum Angle Adjustable Footplates. Not available with Titanium Frame. See Minimum Front Seat Height chart below.

Minimum Front Seat Height with Front Mount Hangers						
Hanger Length	Short		Medium		Long	
	70°	80/90°	70°	80/90°	70°	80°
Minimum Front Seat Height with Standard Selection	13.5"	13.5"	14"	15"	16"	17"
Minimum Front Seat Height with 3" Extension	14.5"	16"	16.5"	18"	18.5"	20"

Hanger Release

If no selection is made Classic Hanger Release will be used

- Classic _____ Std
- 4-Way - In line **New!** _____ N/C
See reference A below. Not available with ELR or PRO ELR.
- 4-Way - Out **New!** _____ N/C
See reference B below



3" Extension

Available only with Front Mount Hangers. Provided as parts.

- 3" Extension for Front Mount Hanger _____ N/C

Residual Limb Support

- Residual Limb Support E1020
 - Left _____ \$395
 - Right _____ \$395
 - Omit Left Hanger ___ N/C
 - Omit Right Hanger ___ N/C

Cover

- Smooth
- Air Mesh

Pad Size

- | | | | |
|---------------------------------|----------------------------------|---------------------------------|----------------------------------|
| <input type="checkbox"/> 5 x 8 | <input type="checkbox"/> 9 x 14 | <input type="checkbox"/> 5 x 8 | <input type="checkbox"/> 9 x 14 |
| <input type="checkbox"/> 7 x 10 | <input type="checkbox"/> 11 x 10 | <input type="checkbox"/> 7 x 10 | <input type="checkbox"/> 11 x 10 |
| <input type="checkbox"/> 9 x 10 | <input type="checkbox"/> 11 x 14 | <input type="checkbox"/> 9 x 10 | <input type="checkbox"/> 11 x 14 |

Footrest continued

Heel Loop

- Heel Loops E0951 _____ \$30
Not available with One Piece Flip Up Angle Adjustable Footplate
- Omit Heel Loops _____ N/C

Calf Strap

- Velcro Adjustable K0038 _____ \$35
- Padded Velcro Adjustable K0038 _____ \$68
- Bodypoint Padded AeroMesh K0038 _____ \$102

Front Seat Height and Casters

Caster Size / Front Frame Type / Front Seat Height

Select Front Seat Height in row corresponding to desired caster size.

Note: Selections in first three Front Seat Height columns use Super Low Front Frame. Super Low Front Frame not available with Titanium Frame.

Caster Size	Front Seat Height													
	Uses Super Low Front Frame													
4"	<input type="checkbox"/> 13.5 ¹	<input type="checkbox"/> 14	<input type="checkbox"/> 14.5	<input type="checkbox"/> 15 ²	<input type="checkbox"/> 15.5	<input type="checkbox"/> 16	<input type="checkbox"/> 16.5	<input type="checkbox"/> 17	<input type="checkbox"/> 17.5	<input type="checkbox"/> 18	<input type="checkbox"/> 18.5	<input type="checkbox"/> 19	<input type="checkbox"/> 19.5	<input type="checkbox"/> 20
5"	<input type="checkbox"/> 14.5	<input type="checkbox"/> 15	<input type="checkbox"/> 15.5	<input type="checkbox"/> 16	<input type="checkbox"/> 16.5	<input type="checkbox"/> 17	<input type="checkbox"/> 17.5	<input type="checkbox"/> 18	<input type="checkbox"/> 18.5	<input type="checkbox"/> 19	<input type="checkbox"/> 19.5	<input type="checkbox"/> 20	<input type="checkbox"/> 20.5	
6"	<input type="checkbox"/> 15.5	<input type="checkbox"/> 16	<input type="checkbox"/> 16.5	<input type="checkbox"/> 17	<input type="checkbox"/> 17.5	<input type="checkbox"/> 18	<input type="checkbox"/> 18.5	<input type="checkbox"/> 19	<input type="checkbox"/> 19.5	<input type="checkbox"/> 20	<input type="checkbox"/> 20.5	<input type="checkbox"/> 21		
7"	<input type="checkbox"/> 17	<input type="checkbox"/> 17.5		<input type="checkbox"/> 18	<input type="checkbox"/> 18.5	<input type="checkbox"/> 19	<input type="checkbox"/> 19.5	<input type="checkbox"/> 20	<input type="checkbox"/> 20.5	<input type="checkbox"/> 21				
8"	<input type="checkbox"/> 17.5	<input type="checkbox"/> 18	<input type="checkbox"/> 18.5	<input type="checkbox"/> 19	<input type="checkbox"/> 19.5	<input type="checkbox"/> 20	<input type="checkbox"/> 20.5	<input type="checkbox"/> 21						

¹ Front Seat Height not available with 4 x 1" Poly Caster

² Super Low Front Frame required when using 4 x 1" Poly Caster

Caster Type

- .75" Rollerblade _____ N/C
Not available with 5", 6", 7" and 8" Caster Sizes
- .75" Lighted Rollerblade _____ \$15
Not available with 5", 6", 7" and 8" Caster Sizes
- 1" Poly _____ N/C
- 1" Poly Aluminum _____ \$110
Not available with 6", 7" and 8" Caster Sizes
- 1" Pneumatic _____ \$80
Not available with 4", 5" and 7" Caster Sizes
- 1.5" Poly _____ N/C
Not available with 7" and 8" Caster Sizes
- 1.5" Soft Roll Aluminum _____ \$50
Not available with 7" and 8" Caster Sizes
- 2" Poly _____ N/C
Not available with 4", 5", 7" and 8" Caster Sizes
- 2" Pneumatic w/ Foam Insert _____ \$125
Not available with 4", 5", 6" and 7" Caster Sizes

Caster Housing

- Standard _____ N/C
- Performance _____ \$35

Maximum Seat Dump Reference

(Front Seat Height minus Rear Seat Height)

Seat Depth	14"	15"	16"	17"	18"	19"	20"
Maximum Dump Standard Caster Housing	1.5"	2"	2"	2"	2.5"	2.5"	2.5"
Maximum Dump Performance Caster Housing	3"	3"	3"	3.5"	3.5"	4"	4"

Depending on configuration, chairs with Fixed Front Frame Type may achieve up to an additional .75"

Fork Stem

- Best Fit _____ Std
Ki Mobility will pick optimal Fork Stem based on Caster Size and Front Seat Height
- Standard Fork Stem _____ N/C
- 3/4" Longer Fork Stem _____ N/C
- 1 1/2" Longer Fork Stem _____ N/C

Fork Options

- Frog Legs Ultra-Sport Fork - Black _____ \$350
Available only with 4", 5" or 6" casters. Not available with Caster Pin Locks, Heavy Duty Option, 2" Poly caster or 2" Pneumatic w/ Foam Insert caster.
End User Weight Required: _____ lbs
- Frog Legs PHASE 2 Carbon Fork **New!** _____ \$425
Lightweight shock absorbing fork. Available only with 4", 5" or 6" casters. Not available with Caster Pin Locks, Heavy Duty Option, 2" Poly caster or 2" Pneumatic w/ Foam Insert caster.
- Caster Pin Locks K0073 _____ \$195
Not available with Performance Caster Housing, 1 1/2" Longer Fork Stem or Frog Legs Fork

Wheel Locks

Wheel Lock	
<input type="checkbox"/> Push to Lock	Std
<input type="checkbox"/> Push to Lock Flush Mount	N/C
<input type="checkbox"/> Pull to Lock	N/C
<input type="checkbox"/> Low Profile Scissor Lock	N/C
<i>Not available with Low Poly tire</i>	
<input type="checkbox"/> Short Thro Scissor Lock	\$110
<input type="checkbox"/> Omit Wheel Locks	N/C

Wheel Lock Options	
<input type="checkbox"/> Extension Handles E0961	\$67
<i>Available only with Push to Lock or Pull to Lock Wheel Locks</i>	
<input type="checkbox"/> Grade Aids E0974	\$150
<i>Available only with Push to Lock Wheel Locks and Pneumatic or Pneumatic with Airless Insert Tire Types</i>	

Rear Seat Height and Wheels

Wheel Size / Tire Profile / Rear Seat Height

Select Rear Seat Height in row corresponding to Wheel Size and Tire Profile. 14" Seat Depth not available with 24", 25" or 26" Rear Wheels.

Requires curved axle plate Any axle plate can be used (Curved, Standard, Amputee, One Arm Drive)

Wheel Size	Tire Profile	Rear Seat Height									
20"	Low	<input type="radio"/> 11.5	<input type="checkbox"/> 12 ¹	<input type="radio"/> 12.5	<input type="checkbox"/> 13 ¹	<input type="radio"/> 13.5	<input type="checkbox"/> 14	<input type="radio"/> 14.5	<input type="checkbox"/> 15	<input type="radio"/> 15.5	<input type="radio"/> 16
	Full	<input type="radio"/> 12	<input type="checkbox"/> 12.5 ¹	<input type="radio"/> 13	<input type="checkbox"/> 13.5 ¹	<input type="radio"/> 14	<input type="checkbox"/> 14.5	<input type="radio"/> 15	<input type="checkbox"/> 15.5	<input type="radio"/> 16	<input type="radio"/> 16.5
22"	Low	<input type="radio"/> 12.5	<input type="checkbox"/> 13 ¹	<input type="radio"/> 13.5	<input type="checkbox"/> 14 ¹	<input type="radio"/> 14.5	<input type="checkbox"/> 15	<input type="radio"/> 15.5	<input type="checkbox"/> 16	<input type="radio"/> 16.5	<input type="radio"/> 17
	Full	<input type="radio"/> 13	<input type="checkbox"/> 13.5 ¹	<input type="radio"/> 14	<input type="checkbox"/> 14.5 ¹	<input type="radio"/> 15	<input type="checkbox"/> 15.5	<input type="radio"/> 16	<input type="checkbox"/> 16.5	<input type="radio"/> 17	<input type="radio"/> 17.5
24"	Low	<input type="radio"/> 13.5	<input type="checkbox"/> 14 ¹	<input type="radio"/> 14.5	<input type="checkbox"/> 15 ¹	<input type="radio"/> 15.5	<input type="checkbox"/> 16	<input type="radio"/> 16.5	<input type="checkbox"/> 17	<input type="radio"/> 17.5	<input type="radio"/> 18
	Full	<input type="radio"/> 14	<input type="checkbox"/> 14.5 ¹	<input type="radio"/> 15	<input type="checkbox"/> 15.5 ¹	<input type="radio"/> 16	<input type="checkbox"/> 16.5	<input type="radio"/> 17	<input type="checkbox"/> 17.5	<input type="radio"/> 18	<input type="radio"/> 18.5
25"	Low	<input type="radio"/> 14	<input type="checkbox"/> 14.5 ¹	<input type="radio"/> 15	<input type="checkbox"/> 15.5 ¹	<input type="radio"/> 16	<input type="checkbox"/> 16.5	<input type="radio"/> 17	<input type="checkbox"/> 17.5	<input type="radio"/> 18	<input type="radio"/> 18.5
26"	Low	<input type="radio"/> 14.5	<input type="checkbox"/> 15 ¹	<input type="radio"/> 15.5	<input type="checkbox"/> 16 ¹	<input type="radio"/> 16.5	<input type="checkbox"/> 17	<input type="radio"/> 17.5	<input type="checkbox"/> 18	<input type="radio"/> 18.5	<input type="radio"/> 19
	Full	<input type="radio"/> 15	<input type="checkbox"/> 15.5 ¹	<input type="radio"/> 16	<input type="checkbox"/> 16.5 ¹	<input type="radio"/> 17	<input type="checkbox"/> 17.5	<input type="radio"/> 18	<input type="checkbox"/> 18.5	<input type="radio"/> 19	<input type="radio"/> 19.5

¹ Not available with Amputee Axle Plate

One Arm Drive

<input type="checkbox"/> One Arm Drive E0958	\$1095
<i>Available only with 24" Mag wheels. Not available with 14" Rear Seat Height in combination with 24" rear wheel or 4° camber. Standard Handrim will be Aluminum Anodized. Axle not quick release.</i>	
Drive Handrim Side:	<input type="checkbox"/> Left <input type="checkbox"/> Right
Drive Handrim Type:	<input type="checkbox"/> Aluminum N/C <input type="checkbox"/> Poly Coated \$200

Axle Plate

No selection needed with One Arm Drive

<input type="checkbox"/> Standard	Std
<input type="radio"/> Curved	\$85
<input type="checkbox"/> Amputee	\$175
<i>Not available with Swing Away Armrests or Wheelchair Drive Ready Option. Required with Reclining Backrest.</i>	

Rear Wheel Type

<input type="checkbox"/> Spoke	N/C
<input type="checkbox"/> Mag	N/C
<i>Not available in 25" or 26" rear wheels</i>	
<input type="checkbox"/> Superlight Spoke	\$295
<i>Available only in 24" and 25" rear wheels. Not available with Heavy Duty Option.</i>	
<input type="checkbox"/> Spinergy Spox	\$795
<i>Not available in 20" rear wheel</i>	
Color:	<input type="checkbox"/> Black <input type="checkbox"/> Red <input type="checkbox"/> White <input type="checkbox"/> Blue <input type="checkbox"/> Yellow
Handrim:	<input type="checkbox"/> Short Tab <input type="checkbox"/> Long Tab
<i>Available only with Aluminum Anodized and Plastic Coated Handrims</i>	
<input type="checkbox"/> Spinergy LX	\$995
<i>Not available in 20" rear wheel</i>	
Color:	<input type="checkbox"/> Black <input type="checkbox"/> Red <input type="checkbox"/> White <input type="checkbox"/> Blue <input type="checkbox"/> Yellow <input type="checkbox"/> Pink <input type="checkbox"/> Green <input type="checkbox"/> Orange
Handrim:	<input type="checkbox"/> Short Tab <input type="checkbox"/> Long Tab
<i>Available only with Aluminum Anodized and Plastic Coated Handrims</i>	

Rear Tire Type

Full Profile

<input type="checkbox"/> Street Pneumatic 1 3/8" E2211, E2212	Std
<i>Not available on 25" rear wheel</i>	
<input type="checkbox"/> Full Poly 1 3/8" E2218	N/C
<i>Not available on Superlight or 25" rear wheel</i>	
<input type="checkbox"/> Pneumatic with Airless Insert 1 3/8" E2211, E2213	\$80
<i>Not available on Superlight or 25" rear wheel</i>	
<input type="checkbox"/> Schwalbe HS-127 1 3/8" E2211, E2212 New!	\$125
<i>Available only on 24" rear wheel. Not available on Mag wheels.</i>	
<input type="checkbox"/> Euro Trek Knobby 1 3/8" E2211, E2212	\$95
<i>Available only on 24" rear wheel</i>	

Low Profile

<input type="checkbox"/> Low Poly 1" E2218	N/C
<i>Not available on Superlight, 25" or 26" rear wheels</i>	
<input type="checkbox"/> Sentinel Solid Tire 1" E2218	\$25
<i>Not available on Superlight wheels</i>	
<input type="checkbox"/> High Pressure 1" E2211, E2212	N/C
<i>Not available on Mag wheels</i>	
<input type="checkbox"/> Schwalbe One 1" E2211, E2212 New!	\$95
<i>Not available on Mag, 20", 22" or 26" rear wheels</i>	
<input type="checkbox"/> Schwalbe Marathon Plus Evolution 1" E2211, E2212	\$95
<i>Not available on Mag, 20" or 22" rear wheels</i>	

Rear Seat Height and Wheels continued

Handrim

No selection needed with One Arm Drive

Standard Handrims

- Aluminum Anodized _____ Std
- Superlight Aluminum Anodized _____ \$50
Available only on 24" and 25" Superlight and Spinergy wheels. Not available with Heavy Duty Option.
- Plastic Coated _____ \$50
- Projection (8 Vertical Projections) _____ \$295
Available on 22", 24", 25" and 26" Spoke wheels & 24" Mag wheels w/full tires
- Spinergy Flex Rim _____ \$600
Available only on 24" and 25" Spinergy LX wheels
- Omit Handrim _____ N/C

Natural Fit Handrims

Not available on 20" rear wheel or 22" Mag wheels. Not recommended with Mag wheels.

- Natural Fit LT Regular Grip _____ \$445
- Natural Fit Standard Regular Grip _____ \$445
- Natural Fit LT Super Grip _____ \$495
- Natural Fit Standard Super Grip _____ \$495
- Natural Fit LT Omit Grip _____ \$345
- Natural Fit Standard Omit Grip _____ \$345

Handrim Option

- 3M Micro Grip _____ \$25
Not available on Plastic Coated or Projection Handrims

Axles

Axle

No selection needed with Wheelchair Drive Ready Option and One Arm Drive

- Quick Release _____ Std
- Quad Quick Release _____ \$125
Not available with Mag wheels

Axle Position

- Rear Wheel Camber: 0° 2° 4°
 Center of Gravity: Rear Center Forward

Armrests

Armrest Please select Arm Pad Option

- Height Adjustable T-Arm (8.5" - 13.5") E0973 _____ \$265
- Height Adjustable T-Arm Low (7.5" - 11.5") E0973 _____ \$265
- Height Adjustable Flip Back T-Arm (9.5"-13.5") E0973 **New!** _____ \$275
Not available with Depth Adjustable Back Posts
- Flip Back (10.5") _____ N/C
Not available with Seat Depths less than 16" or Angle Adjustable or Depth Adjustable Back Posts
- Angle Adjustable Locking Extendable Flip Up E0973 _____ \$265
Available only with Depth Adjustable Back Posts

Arm Pad

Not available with Swing Away Armrests

- Standard Arm Pad _____ N/C Desk Full
 Waterfall Arm Pad _____ \$25 Desk Full
Available only with Height Adjustable T-Arms
 Foam Grip _____ N/C Desk Full
Available only with Angle Adjustable Locking Extendable Flip Up Armrest

Armrest continued

- Swing Away (10.5" - 12.5") E0973 _____ \$195
Not available with Angle Adjustable or Depth Adjustable Back Posts or Arm Pads
- Swing Away with Side Guard (10.5" - 12.5") E0973 _____ \$285
Not available with Angle Adjustable or Depth Adjustable Back Posts or Arm Pads
- Omit Armrests _____ N/C

Accessories

Removable Side Guard

- Composite Side Guards K0108 _____ \$140
Not available with T-Arm, Flip Back or Swing Away with Side Guard Armrests

Anti-Tip

- Rear Anti-Tips E0971 _____ \$140

Positioning Belt

- 1 1/2" Auto Buckle E0978 _____ \$70
- 1 1/2" Padded Auto Buckle E0978 _____ \$155
- 2" Airline Buckle E0978 _____ \$110
- Bodypoint Evoflex - Pelvic Stabilizer E0978 _____ \$183
 Size: Small (16-24") Medium (19-29") Large (25-37")
- Bodypoint Belt Mounting Kit _____ \$105
 Style: Bracket Band Clamp 1"

Seat Pan

- Removable Aluminum Seat Pan E2231 _____ \$245
Not available with Transit System. Pan depth 1/2" shorter than Seat Depth.

Frame Color

Frame Color not available with Titanium Frame Option

Gloss Finish

- Black _____ N/C
- Blue Skies _____ N/C
- Bubblegum Pink _____ N/C
- Candy Blue _____ N/C
- Candy Purple _____ N/C
- Candy Red _____ N/C
- Copper _____ N/C
- Fuchsia _____ N/C
- Granny Smith _____ N/C
- Yellow _____ N/C

Other Accessories

- Tool Kit _____ \$45
- Removable Underseat Pouch _____ \$50
- Ki Mobility Backpack _____ \$65
- Luggage Carrier **New!** _____ \$85
Not available with ELR, PRO ELR or Residual Limb Support
- Black Neoprene Impact Guards _____ \$85
- Clear Spoke Guards K0065 _____ \$115
Not available with Mag wheels
- Black Spoke Guards K0065 _____ \$115
Not available with 20" rear wheel or Mag wheels
- Cane and Crutch Holder _____ \$195
- IV Pole _____ \$265
Available only with Curved Rear Frame Type and Straight Back Posts

O₂ Holder

Not available with Depth Adjustable Rear Frame

- O₂ Holder E2208 _____ \$295
 Right Side Left Side

Matte Finish

- Army Green _____ N/C
- Charcoal Gray _____ N/C
- Emerald Green _____ N/C
- Iridescent Burgundy _____ N/C
- Iridescent Navy _____ N/C
- Retro Red _____ N/C
- Royal Blue _____ N/C
- Shadow Black _____ N/C
- Smooth White _____ N/C

Graphic Image Transfer

*Visit www.kimobility.com for current graphic image transfer options
 Note selection below*

- _____ \$595

Catalyst Front Seat Height Matrix																				
4" Caster				5" Caster				6" Caster				7" Caster				8" Caster				
FSH	Super Low	Stem	Fork	FSH	Super Low	Stem	Fork	FSH	Super Low	Stem	Fork	FSH	Super Low	Stem	Fork	FSH	Super Low	Stem	Fork	
15 ¹	13.5 ¹	Std	4"	16	14.5	Std	4"	17	15.5	Std	5"	18	17	Std	6"	19	18	Std	6"	
15.5	14			16.5	15.5	3/4		18	16.5	3/4		18.5	17.5			20	18.5	3/4		
16 ²	14.5 ²	17.5		16	1 1/2	18.5		17	1 1/2	19		17.5	20.5	n/a		1 1/2				
16.5	15	3/4		16	15	Std	17	16	Std	19.5	18	19	17.5	Std						
16.5 ²	15 ²	1 1/2		16.5	15.5		17.5	16.5		19.5	18.5	19.5	18			20	18.5			
17	15.5	1 1/2		17	15.5	3/4	5"	18	17	Std	20	19	1 1/2	20		18.5	Std			
15.5	14.5	Std		17.5	16	18		16.5	18.5		17.5	20.5		19						
16	15	3/4		5"	17.5	16.5	1 1/2	6"	18.5	17	3/4	6"	19	18		Std	19.5	18.5	3/4	7"
16.5	15				18	17	19		17.5	19.5	18.5		20	19						
17	15.5	16.5			15.5	18.5	17.5		20	19	20.5		n/a							
17	16	1 1/2	17		16	Std	19	18	1 1/2	19	18	3/4	21	n/a	1 1/2					
17.5	16.5	17.5	16.5		19.5	18.5	20.5	19	21	n/a										
17	15.5	Std	6"		17.5	16	3/4	6"	18	16.5	Std	7"	20	19	1 1/2	8" x 2" caster not available with 6" fork				
17.5	16				18	16.5			18.5	17			20.5	n/a						
17.5	16.5	18.5			17	19			17.5	20			18.5							
18	17	3/4			18	17	1 1/2	19.5	18.5	20.5	19	1 1/2								
18.5	17	1 1/2			18.5	17.5		18.5	17.5	21	n/a									
19	17.5	1 1/2		19	18	3/4	19	18	7"	3/4	7"									
18	17	Std		18	16.5		19.5	18.5												
18.5	17.5	3/4		7"	18.5	17	Std	7"	20	19	1 1/2	7"	19.5	18	1 1/2					
19	17.5				19	17.5			19.5	18			20	18.5						
19.5	18	18.5			17.5	20			18.5	20.5			19							
19.5	18.5	1 1/2	19		18	3/4	7"	20.5	19	1 1/2	7"									
20	19		1 1/2		19.5	18.5	21	n/a												
					19.5	18	1 1/2	6" x 2" caster not available with 5" fork												
					20	18.5														
					20.5	19														

¹ Front Seat Height not available with 4 x 1" Poly Caster

² Not available with 4 x 1" Poly, 4 x 1.5" Poly or 4 x 1.5" Soft Roll Aluminum

Seat heights with "FSH" heading available with Standard Front Frame. Seat heights with "Super Low" heading available with Super Low Front Frame. Actual seat height may vary depending upon configuration.

Cushions

General Use Cushions

STEP 1: Choose Cushion

- 2" General Use Foam Cushion (14"x14" to 20"x20") E2601 ____ \$99
- 2" General Use Foam Cushion (Over 20"x20") E2601 ____ \$250

STEP 2: Choose Cushion Size

- Match to Chair Size
- Other Size _____

Matrx Seat Cushions

Cushion size is determined by chair dimensions. Certain sizes will have higher price and longer lead time; visit www.motionconcepts.com/matrx_international.html for size specifications.

STEP 1: Choose Cushion

- PS Cushion _____ \$414
Cushions 14" wide have a weight capacity of 150 lbs. Cushions 15" through 20" have a weight capacity of 300 lbs.
- Vi Cushion _____ \$510
Weight capacity of 300 lbs
- Libra Cushion _____ \$582
Weight capacity of 300 lbs

STEP 2: Choose Cushion Size

- Match to Chair Size
- Other Size - Provide Part Number _____

Other Matrx Products

Must provide part number

- _____

ROHO Seat Cushions

Cushion size is determined by chair dimensions. Certain sizes will have higher price and longer lead time; visit therohogroup.com for size specifications.

STEP 1: Choose Cushion

Standard Cushions

For chairs under 22" wide

- HIGH PROFILE® Single Compartment (E2622) ____ \$603
- MID PROFILE™ Single Compartment (E2622) ____ \$603
- LOW PROFILE® Single Compartment (E2622) ____ \$603
- QUADTRO SELECT® High Profile® (E2624) ____ \$603
- QUADTRO SELECT® Mid Profile® (E2624) ____ \$603
- HYBRID ELITE® Single Comp Smart Check Ready (E2622) ____ \$623
- CONTOUR SELECT® (E2624) ____ \$603

X Large Cushions

For 22" wide chairs

- HIGH PROFILE® Single Compartment X Large (E2623) ____ \$670
- LOW PROFILE® Single Compartment X Large (E2623) ____ \$670
- QUADTRO SELECT® High Profile® X Large (E2625) ____ \$670

Cushions with Heavy Duty Cover

- HIGH PROFILE® Single Compartment w/ HD Cover (E2622) ____ \$629
- QUADTRO SELECT® High Profile® w/ HD Cover (E2624) ____ \$629

STEP 2: Choose Cushion Size

- Match to Chair Size
- Other Size - Provide Part Number _____

Backs

<p>Matrix Backs</p> <p><input type="checkbox"/> Elite Back _____ \$750 <i>Not available with 22" Seat Width. Back Height of 14" not available with 19", 20" and 21" Seat Widths.</i></p> <p>Back Height: <input type="checkbox"/> 10" <input type="checkbox"/> 12" <input type="checkbox"/> 14" <input type="checkbox"/> 16" <input type="checkbox"/> 18" <input type="checkbox"/> 20"</p> <p><input type="checkbox"/> Elite Deep Back _____ \$912 <i>Not available with 22" Seat Width. Back Height of 10" not available with 14" Seat Width.</i></p> <p>Back Height: <input type="checkbox"/> 10" <input type="checkbox"/> 12" <input type="checkbox"/> 16" <input type="checkbox"/> 18" <input type="checkbox"/> 20"</p> <p><input type="checkbox"/> Mini Back _____ \$750 <i>Available only with 14" and 15" Seat Widths. Comes standard with Kidabra cover.</i></p> <p>Back Height: <input type="checkbox"/> 7" <input type="checkbox"/> 10" <input type="checkbox"/> 12" <input type="checkbox"/> 14"</p>	<p>Therafin Headrest for MATRX Back Systems</p> <p><input type="checkbox"/> Headrest, Multi-Axis Offset Hardware with Smooth/Soft Pad (E1028, E0955) 34187 _____ \$420</p> <p>PAD SIZE: <input type="checkbox"/> 5 x 8 <input type="checkbox"/> 5 x 10 <input type="checkbox"/> 5 x 14</p>
	<p>Other Matrix Mounting Hardware</p> <p><i>Must provide part number</i></p> <p><input type="checkbox"/> _____ <i>For mounting hardware other than quick release</i></p>
	<p>Other Matrix Products</p> <p><i>Must provide part number</i></p> <p><input type="checkbox"/> _____</p>

ROHO Back Systems (Available for chairs up to 20" wide)

STEP 1: Choose Back

<p><input type="checkbox"/> AGILITY™ Min. Contour w/ Foam Only (E2611) _____ \$549 BACKREST HEIGHT: <input type="checkbox"/> 10" <input type="checkbox"/> 13" <input type="checkbox"/> 16" <input type="checkbox"/> 18"</p> <p><input type="checkbox"/> AGILITY™ Min. Contour w/ Foam & Lumbar Pad (E2613) _____ \$659 BACKREST HEIGHT: <input type="checkbox"/> 10" <input type="checkbox"/> 13" <input type="checkbox"/> 16" <input type="checkbox"/> 18"</p> <p><input type="checkbox"/> AGILITY™ Min. Contour w/ Air Only (E2620) _____ \$769 BACKREST HEIGHT: <input type="checkbox"/> 10" <input type="checkbox"/> 13" <input type="checkbox"/> 16" <input type="checkbox"/> 18"</p> <p><input type="checkbox"/> AGILITY™ Min. Contour w/ Air & Lumbar Pad (E2615) _____ \$769 BACKREST HEIGHT: <input type="checkbox"/> 10" <input type="checkbox"/> 13" <input type="checkbox"/> 16" <input type="checkbox"/> 18"</p>	<p><input type="checkbox"/> AGILITY™ Mid Contour w/ Air Only (E2620) _____ \$781 BACKREST HEIGHT: <input type="checkbox"/> 10" <input type="checkbox"/> 13" <input type="checkbox"/> 16" <input type="checkbox"/> 18" <input type="checkbox"/> 20"</p> <p><input type="checkbox"/> AGILITY™ Mid Contour w/ Air & Lumbar Pad (E2615) _____ \$781 BACKREST HEIGHT: <input type="checkbox"/> 10" <input type="checkbox"/> 13" <input type="checkbox"/> 16" <input type="checkbox"/> 18" <input type="checkbox"/> 20"</p> <p><input type="checkbox"/> AGILITY™ Max Contour w/ Air Only (E2620) _____ \$920 BACKREST HEIGHT: <input type="checkbox"/> 16" <input type="checkbox"/> 18" <input type="checkbox"/> 20"</p> <p><input type="checkbox"/> AGILITY™ Max Contour w/ Air & Lumbar Pad (E2615) _____ \$920 BACKREST HEIGHT: <input type="checkbox"/> 16" <input type="checkbox"/> 18" <input type="checkbox"/> 20"</p>
---	---

STEP 2: Choose Back Hardware

Fixed _____ N/C
Not compatible with 18" and 20" tall backs

Quick Release _____ N/C

ROHO Laterals for AGILITY™ Back Systems

Sold in pairs

Swing Away with Curved Pad (E1028/E0956) _____ \$1014

PAD SIZE: Medium Large X-Large
 4 T x 4 L 4 T x 5 L 5 T x 5 L

OFFSET TYPE: Flush Offset 1" Offset 2" Offset

Fixed with Flat Pad (E0956) _____ \$452

PAD SIZE: Medium Large X-Large
 4 T x 4 L 4 T x 5 L 5 T x 5 L

OFFSET TYPE: Flush Offset

Therafin Headrest for AGILITY Back Systems

Headrest, Multi-Axis Offset Hardware with Smooth/Soft Pad (E1028, E0955) 34187 _____ \$420

PAD SIZE: 5 x 8 5 x 10 5 x 14



Vest

Not compatible with Back Upholstery

- Full Shape, Comfort Strap (E0960) _____ \$167
 - Small 30589 Medium 30590 Large 30591
 - X-Large 30592 2X-Large 30593
- Full Shape, Extended Strap (E0960) _____ \$167
 - Small 30594 Medium 30595 Large 30596
 - X-Large 30597 2X-Large 30598
- Trim Shape, Comfort Strap (E0960) _____ \$169
 - Small 30618 Medium 30619 Large 30620
 - X-Large 30621 2X-Large 30622
- Trim Shape, Extended Strap (E0960) _____ \$169
 - Small 30604 Medium 30624 Large 30625
 - X-Large 30626 2X-Large 30627

Shoulder Harness

With metal clips and cam locks. Not compatible with Back Upholstery.

- Small (E0960) 61403 _____ \$157
- Medium (E0960) 61503 _____ \$164
- Large (E0960) 61703 _____ \$170
- X-Large (E0960) 61803 _____ \$177

Hip Belt

- 4 Point Y-Style, Metal Buckle (E0978) _____ \$173
 - Small 60024 Medium 60025 Large 60026
- Dual Rear Pull, Metal Attaching Clips (E0978)
 - 1" Straps, 6.25"W x 1.75"H Pads 60161 _____ \$97
 - 1.5" Straps, 7.25"W x 2.25"H Pads 60171 _____ \$106
 - 1.5" Straps, 9.25"W x 2.5"H Pads 60181 _____ \$106
- Single Pull, Metal Buckle (E0978)
 - 1" Straps, 4.25"W x 1.75"H Pads 61001 _____ \$86
 - 1" Straps, 5.25"W x 1.75"H Pads 61101 _____ \$86
 - 1" Straps, 6.25"W x 1.75"H Pads 61201 _____ \$86
 - 1.5" Straps, 5.25"W x 2.25"H Pads 62001 _____ \$95
 - 1.5" Straps, 9.25"W x 2.5"H Pads 62201 _____ \$95
 - 2" Straps, 9.25"W x 3"H Pads 64001 _____ \$102
 - 2" Straps, 11.25"W x 3"H Pads 64101 _____ \$102

Cane Mounted Swing Away Lateral Thoracic Support

Sold as pair. Not compatible with Back Upholstery.

- Swing Away Kits w/ 4"x6" Gel Pads (E0956, E1028) 50021 _____ \$478
- Swing Away Kits, No Pads (E0956, E1028) 50022 _____ \$309

Knee Adductor

- Flip Down Hardware (E1028) _____ \$634
Sold in pairs, one right side and one left side
 - Short Arm, Fits 3/4", 7/8" and 1" Tubing 30248
 - Long Arm, Fits 3/4", 7/8" and 1" Tubing 30284
- Knee Adductor Pad with Neoprene Cover (E0956) _____ \$117
Sold single
 - Small Pad, 4.25" x 3.25" 50490 Quantity: _____
 - Medium Pad, 6.25" x 4.25" 50491 Quantity: _____

Foot / Shoe Holder

Compatible only with Aluminum Angle Adjustable and Composite Angle Adjustable Footplates

- FootSure Ankle Support (K0108) _____ \$178
Sold in pairs
 - X-Small w/ Buckle Closure, 6"-8" Circum. 60250
 - Small w/ Buckle Closure, 7"-8.5" Circum. 60251
 - Medium w/ Buckle Closure, 8"-11" Circum. 60252
 - Large w/ Buckle Closure, 9"-12" Circum. 60253
 - X-Large w/ Buckle Closure, 12.5"-15" Circum. 60254
 - X-Small w/ Hook-N-Loop Closure, 5"-6" Circum. 60255
 - Small w/ Hook-N-Loop Closure, 6"-8" Circum. 60256
 - Med w/ Hook-N-Loop Closure, 8"-10.5" Circum. 60257
 - Large w/ Hook-N-Loop Closure, 9"-12" Circum. 60258
 - X-Lrg w/ Hook-N-Loop Closure, 12"-15" Circum. 60259
- Padded Footbox (K0108)
 - Small, Full Design, Fits 14" Chair 30111 _____ \$287
 - Medium, Full Design, Fits 16" Chair 30112 _____ \$293
 - Large, Full Design, Fits 18" Chair 30113 _____ \$298
 - Small, Split Design, Fits 14" Chair 30152 _____ \$367
 - Medium, Split Design, Fits 16" Chair 30153 _____ \$373
 - Large, Split Design, Fits 18" Chair 30154 _____ \$378
- Molded Shoe Holders (K0108) _____ \$162
Sold in pairs
 - Small with Plain Straps, 3.5" x 5.25" 30293
 - Medium with Plain Straps, 4" x 6.5" 30303
 - Large with Plain Straps, 4.75" x 8" 30309
 - X-Large with Plain Straps, 5" x 10" 30316



Half Tray

Not compatible with Swing Away Armrests

- Standard Shape, Clear, Flip-Away (E0950, E1028) _____ \$246
 - Right Side 41310 Left Side 41311
- Standard Shape, Padded, Flip-Away (E0950, E1028) _____ \$269
 - Right Side 41350 Left Side 41353
- Trimline Shape, Padded, Flip-Away (E0950, E1028) _____ \$269
 - Right Side 41352 Left Side 41351

Clear Tray with Border and Pre-Drilled Holes

STEP 1: Choose Tray

Chair Width	Tray Dimension (W x D)	Body Cutout (W x D)	HCPCS	Part Number	Price	<input type="checkbox"/>
14"	20 x 18	10.5 x 6	E0950	33210	\$220	<input type="checkbox"/>
	20 x 18	11.5 x 6.5	E0950	33211	\$220	<input type="checkbox"/>
15"	21 x 18	11.5 x 6.5	E0950	33212	\$220	<input type="checkbox"/>
	21 x 18	12.5 x 7	E0950	33213	\$220	<input type="checkbox"/>
16"	22 x 21	12.5 x 7.5	E0950	33214	\$255	<input type="checkbox"/>
	22 x 21	13.5 x 8.5	E0950	33215	\$255	<input type="checkbox"/>
17"	23.5 x 22	13.5 x 8.5	E0950	33216	\$255	<input type="checkbox"/>
	23.5 x 22	14.5 x 9.5	E0950	33217	\$285	<input type="checkbox"/>
18"	24.5 x 22	14.5 x 9.5	E0950	33218	\$285	<input type="checkbox"/>
	24.5 x 22	15.5 x 10.5	E0950	33219	\$285	<input type="checkbox"/>
19"	26.5 x 22	15.5 x 10.5	E0950	33220	\$285	<input type="checkbox"/>
	26.5 x 22	16.5 x 11	E0950	33221	\$285	<input type="checkbox"/>
20"	26.5 x 22	17 x 11	E0950	33222	\$309	<input type="checkbox"/>
	26.5 x 22	18 x 11	E0950	33223	\$309	<input type="checkbox"/>
21"	27.5 x 23.5	17.5 x 11	E0950	33224	\$309	<input type="checkbox"/>
	27.5 x 23.5	19 x 11	E0950	33225	\$309	<input type="checkbox"/>
22"	28.5 x 23.5	18.5 x 11	E0950	33226	\$309	<input type="checkbox"/>
	28.5 x 23.5	20 x 11	E0950	33227	\$309	<input type="checkbox"/>

STEP 2: Choose Tray Hardware

- E-Z Lock Clamps, Tube Arm (E1028) 31072 _____ \$150
Available only with Swing Away, Swing Away with Side Guard and Angle Adjustable Locking Extendable Flip Up with Foam Grip Armrests. Sold in Pairs.
- E-Z Lock Clamps, Standard (E1028) 31098 _____ \$142
Not available with Swing Away, Swing Away with Side Guard and Angle Adjustable Locking Extendable Flip Up with Foam Grip Armrests. Sold in Pairs.

Communication Mount Hardware

Choose Device Holder AND Mount

STEP 1: What DEVICE will be used?

- Smart Phone 32359 _____ \$150
- Large Smart Phone 32354 _____ \$150
- I-Pad 32355 _____ \$184
- I-Pad Mini 32356 _____ \$178
- 7" Tablet 32357 _____ \$173
- 10" Tablet 32358 _____ \$395

STEP 2: Where does it need to MOUNT?

- Frame 32350 _____ \$512
- Armrest 32351 _____ \$395
Not compatible with Swing Away Armrests
- Tray 32352 _____ \$212

Cup Holder

- Heavy Duty 30889 _____ \$44
Compatible only with Flip Back Armrests