



Phone (800) 981-1540 | Fax (715) 254-0996
sales@kimobility.com | www.kimobility.com

CATALYST 5Ti

Back Upholstery Options

Select option below only if Standard or Angle Adjustable Back Post option is selected on chart on previous page. If selecting both upholstery and Axiom Back, the upholstery will be installed.

Standard	N/C
<i>Not available with Short or Medium Back Height Ranges on 21" - 22" Seat Widths</i>	
Tension Adjustable with Standard Cover (E2611)	\$345
Tension Adjustable with Air Mesh Breathable Cover (E2611)	\$405
Omit	(\$15)

U.S. Military Patches

Attached only if ordered with Tension Adjustable Back Upholstery

Army	Navy	Marines	Air Force	
Coast Guard	US Purple Heart			\$92
Combo Army/Purple Heart		Combo Navy/Purple Heart		
Combo Marines/Purple Heart		Combo Air Force/Purple Heart		
Combo Coast Guard/Purple Heart				\$149

Frame Type / Footrest

Hanger / Front Frame Bend

Frame	Hanger	70°	80°		Retail
Swing Away	Extension Tube Mount		1		N/C
	Pro ELR				\$320
	Omit Hangers				(\$20)
Fixed	Extension Tube Mount		1		N/C

1 Not available with 7" or 8" Casters

PRO ELR Cover

Available only on PRO ELR

Gel Pad Covers (PAIR)	\$436
-----------------------	-------

Footplates

Composite	N/C
Composite Angle Adjustable (K0040)	\$200
<i>Not available with 14" Seat Width</i>	
Aluminum Angle Adjustable (K0040)	\$205
One Piece Flip Up Angle Adjustable (K0040)	\$200
<i>Not available with PRO ELR</i>	

Footrest / Extension / Seat to Footrest Length

Select Extension and Seat to Footrest Length in row corresponding to desired Footrest. Footrest Length of at least 2.5" shorter than Front Seat Height recommended. Extension Tube Hangers and Fixed Front Frame measurements are based on 70° Hanger / Front Frame Bend.

Footrest	Extension	Seat to Footrest Length								
Extension Tube Hangers (70°,80°)	Short (11" - 15")	11"	11.5"	12"	12.5"	13"	13.5"	14"	14.5"	15"
	Medium (13" - 17")	13"	13.5"	14"	14.5"	15"	15.5"	16"	16.5"	17"
	Long (15.5" - 19.5")	15.5"	16"	16.5"	17"	17.5"	18"	18.5"	19"	19.5"
PRO ELR Low Pivot Setting	PRO ELR Low Pivot Setting Short (12" - 17.5")	12" - 17.5"								
	PRO ELR Low Pivot Setting Long (15" - 21")	15" - 21"								
PRO ELR High Pivot Setting	PRO ELR High Pivot Setting Short (10.5" - 15.5")	10.5" - 15.5"								
	PRO ELR High Pivot Setting Long (13" - 19")	13" - 19"								
Fixed Front Frame (70°,80°)	Short (11" - 15")	11.5"	12"	12.5"	13"	13.5"	14"	14.5"	15"	
	Medium (13" - 17")	13"	13.5"	14"	14.5"	15"	15.5"	16"	16.5"	17"
	Long (15.5" - 19.5")	15.5"	16"	16.5"	17"	17.5"	18"	18.5"	19"	19.5"

Hanger Release

If no selection is made Classic Hanger Release will be used

Classic	N/C
4-Way - In-Line	N/C
<i>Not available with PRO ELR.</i>	
4-Way - Out	N/C
<i>Not available with PRO ELR.</i>	

Heel Loops

Not available with One Piece Flip Up Angle Adjustable Footplate

Heel Loops (E0951)	\$42
Adjustable Heel Loops (E0951)	\$51
Omit Heel Loops	N/C

Calf Straps

Velcro Adjustable (K0038)	\$41
Padded Velcro Adjustable (K0038)	\$72
Bodypoint Padded AeroMesh (K0038)	\$133

Residual Limb Support

	Limb Support	Price	Omit	Cover		Pad Size					
Left	Left Limb Support (E1020)	\$441	Omit Left Hanger	Smooth Cover Left	Air Mesh Cover Left	5X8 Pad Left Side	7X10 Pad Left Side	9x10 Pad Left Side	9x14 Pad Left Side	11x10 Pad Left Side	11x14 Pad Left Side
Right	Right Limb Support (E1020)	\$441	Omit Right Hanger	Smooth Cover Right	Air Mesh Cover Right	5X8 Pad Right Side	7X10 Pad Right Side	9x10 Pad Right Side	9x14 Pad Right Side	11x10 Pad Right Side	11x14 Pad Right Side



Phone (800) 981-1540 | Fax (715) 254-0996
sales@kimobility.com | www.kimobility.com

CATALYST 5Ti

Front Seat Height and Casters

Caster Size / Front Seat Height

Select Front Seat Height in row corresponding to desired caster size

Caster Size	Front Seat Height										
4"	15"	15.5"	16"	16.5"	17"	17.5"	18"	18.5"	19"	19.5"	20"
5"	16"	16.5"	17"	17.5"	18"	18.5"	19"	19.5"	20"	20.5"	
6"	17"	17.5"	18"	18.5"	19"	19.5"	20"	20.5"	21"		
7"	18"	18.5"	19"	19.5"	20"	20.5"	21"				
8"	19"	19.5"	20"	20.5"	21"						

Maximum Seat Dump Reference

(Front Seat Height minus Rear Seat Height)

Seat Depth	14"	15"	16"	17"	18"	19"	20"
Maximum Dump Standard Caster Housing	1.5"	2"	2"	2"	2.5"	2.5"	2.5"
Maximum Dump Performance Caster Housing	3"	3"	3"	3.5"	3.5"	4"	4"

Depending on configuration, chairs with Fixed Front Frame Type may achieve up to an additional .75"

Casters

		Caster Size					
		Price	4"	5"	6"	7"	8"
Caster Type	.75" Rollerblade	N/C					
	.75" Lighted Rollerblade	\$77					
	1" Poly	N/C					
	1" Poly Aluminum	\$169					
	1" Pneumatic	\$92					
	1.5" Poly	\$60					
	1.5" Soft Roll Aluminum	\$212					
	2" Poly	\$97					
	2" Pneumatic w/Foam Insert	\$246					

Caster Pin Locks

Caster Pin Locks (K0073) \$210

Not available with Performance Caster Housing, 1 1/2" Longer Fork Stem or Frog Legs Fork

Caster Housing

Standard Housing N/C
Performance Housing \$46
Not available with 350 lb Weight Capacity

Fork / Stem

Best Fit N/C
Ki Mobility will pick Fork and Stem combination to optimize performance based on Caster Size and Front Seat Height
Standard Fork Stem N/C
3/4" Longer Fork Stem N/C
1 1/2" Longer Fork Stem N/C
Frog Legs II Suspension System \$461
Lightweight shock absorbing fork. Available only with 4", 5" or 6" Caster Sizes. Not available with Caster Pin Locks, 350 lb Weight Capacity, 2" Poly caster or 2" Pneumatic w/ Foam Insert caster.
Frog Legs Suspension Fork \$461
Available only with 4", 5" or 6" Caster Sizes. Not available with Caster Pin Locks, 350 lb Weight Capacity, 2" Poly caster or 2" Pneumatic w/ Foam Insert caster.

Frog Legs Option

User Weight Under 80 lbs N/C
User Weight 80-120 lbs N/C
User Weight 120-160 lbs N/C
User Weight 160-220 lbs N/C
User Weight 220-260 lbs N/C
User Weight 260-300 lbs N/C



Phone (800) 981-1540 | Fax (715) 254-0996
sales@kimobility.com | www.kimobility.com

CATALYST 5Ti

Wheel Locks

Wheel locks

Pull to Lock	N/C
Short Thro Scissor Lock	\$118
Drum Brake (E2228)	\$595
<i>Available only with 20", 22" and 24" MAXX Spoke wheels. Not available with Straight without Push Handle or Straight with Fold Down Push Handle Back Posts, or Low Profile Push Handle. Axle is Quick Release.</i>	
Hemi Wheel Lock - Right (K0108)	\$359
<i>Includes Extension Handle</i>	

Wheel Locks continued

Dual - Drum Brake with Pull to Lock (E2228)	\$702
<i>Available only with 20", 22" and 24" MAXX Spoke wheels. Not available with Straight without Push Handle or Straight with Fold Down Push Handle Back Posts, or Low Profile Push Handle. Axle is Quick Release.</i>	
Omit	N/C

Wheel Lock Options

Extension Handles (E0961)	\$81
<i>Available only with Push to Lock or Pull to Lock Wheel Locks</i>	
Grade Aids (E0974)	\$159
<i>Available only with Push to Lock Wheel Locks and Pneumatic or Pneumatic with Airless Insert tires</i>	

Rear Seat Height and Wheels

Wheel Size / Tire Profile / Rear Seat Height

Select Rear Seat Height in row corresponding to Wheel Size and Tire Profile. 14" Seat Depth not available with 24", 25" or 26" wheel sizes.

Wheel	Tire	Rear Seat Height									
		11.5" ³	12" ¹	12.5" ³	13" ¹	13.5" ³	14"	14.5" ³	15"	15.5" ³	16" ³
20"	Low	11.5" ³	12" ¹	12.5" ³	13" ¹	13.5" ³	14"	14.5" ³	15"	15.5" ³	16" ³
	Full	12" ³	12.5" ¹	13" ³	13.5" ¹	14" ³	14.5"	15" ³	15.5"	16" ³	16.5" ³
22"	Low	12.5" ³	13" ¹	13.5" ³	14" ¹	14.5" ³	15"	15.5" ³	16"	16.5" ³	17" ³
	Full	13" ³	13.5" ¹	14" ³	14.5" ¹	15" ³	15.5"	16" ³	16.5"	17" ³	17.5" ³
24"	Low	13.5" ³	14" ¹	14.5" ³	15" ¹	15.5" ³	16"	16.5" ³	17"	17.5" ³	18" ³
	Full	14" ³	14.5" ¹	15" ³	15.5" ¹	16" ³	16.5"	17" ³	17.5"	18" ³	18.5" ³
25"	Low	14" ³	14.5" ¹	15" ³	15.5" ¹	16" ³	16.5"	17" ³	17.5"	18" ³	18.5" ³
26"	Low		15" ^{1,2}	15.5" ³	16" ¹	16.5" ³	17"	17.5" ³	18"	18.5" ³	19" ³
	Full	15" ^{2,3}	15.5" ^{1,2}	16" ³	16.5" ¹	17" ³	17.5"	18" ³	18.5"	19" ³	19.5" ³

¹ Not available with Amputee Axle Plate

² Not available with Height Adjustable Flip Back T-Arm

³ Requires Universal Axle Plate

One Arm Drive

One Arm Drive - Quick Release (E0958)	\$1122
<i>Available only with 20", 22" and 24" spoke wheels and Quick Release Axles. Spoke wheel style is exclusive to OAD. Not available with 4 degrees of camber. Standard handrim will be Aluminum Anodized.</i>	

Drive Side

Left Side Drive	N/C
Right Side Drive	N/C

Drive Handrim

Aluminum Drive Handrim	N/C
Poly Coated Drive Handrim	\$210

Axle Plates

No selection needed with One Arm Drive

Universal	N/C
Amputee (E0959)	\$187
<i>Not available with Swing Away Armrests or Wheelchair Drive Ready Option</i>	

Wheels

		Rear Wheel Size					
		Price	20"	22"	24"	25"	26"
Rear Wheel Type	Mag	N/C					
	Maxx Spoke	\$236					
	Maxx Spoke	N/C					
	Superlight Spoke	\$313				¹	
	Spinergy Spox	\$938					
	Spinergy LX	\$1177					

¹ Available only with Superlight and Plastic Coated Handrims

Wheel Color

Select spoke color and handrim tab for Spinergy wheels

Black Spokes Red Spokes Blue Spokes White Spokes _____ N/C

Rear Tire Type

Full Profile

Street Pneumatic 1 3/8" (E2211, E2212)	N/C
<i>Not available on 25" wheel size</i>	
Full Poly 1 3/8" (E2218)	\$103
<i>Not available on Superlight wheel or 25" or 26" wheel sizes</i>	
Pneumatic with Airless Insert 1 3/8" (E2211, E2213)	\$142
<i>Not available on Superlight wheel or 25" wheel size</i>	
Schwalbe HS-127 1 3/8" (E2211, E2212)	\$133
<i>Available only on 24" wheel size. Not available on Mag wheel</i>	
Schwalbe Downtown 1 3/8" (E2211, E2212)	\$103
<i>Available only on 24" wheel size</i>	

Low Profile

Low Poly 1" (E2218)	\$82
<i>Not available on Superlight wheel, 25" or 26" wheel sizes or 350 lb Weight Cap</i>	
Sentinel Solid Tire 1"	\$103
<i>Not available on Superlight wheel or 20" and 22" Mag wheels</i>	
High Pressure 1" (E2211, E2212)	\$103
<i>Not available on Mag wheel</i>	
Schwalbe Marathon Plus Evolution 1" (E2211, E2212)	\$149
<i>Not available on Mag wheel, 20" or 22" wheel sizes</i>	
Schwalbe One 1" (E2211, E2212)	\$133
<i>Not available on Mag wheel, 20", 22" or 26" wheel sizes</i>	

Handrims

No selection needed with One Arm Drive

Standard Handrim

Aluminum Anodized	N/C
Superlight Aluminum Anodized	\$56
<i>Available only on 24" and 25" wheel sizes. Available only on MAXX Spoke, Superlight or Spinergy wheels.</i>	
Plastic Coated	\$110
Projection (8 Vertical Projections)	\$359
<i>Available only with full profile tires in combination with 22", 24", and 26" MAXX Spoke wheels, or 22" and 24" Mag wheels</i>	
Omit Handrim	N/C

Ergonomic Handrim

Not available with 20" or 26" wheels. Tab mount only available with 24" wheels.

Ergo Lite Regular Grip	\$477
Ergo Standard Regular Grip	\$456
Ergo Lite Extra Grip	\$528
Ergo Standard Extra Grip	\$507
Ergo Lite Omit Grip	\$405
Ergo Standard Omit Grip	\$354

Handrim Options

Micro-Grip

3M Micro Grip	\$72
<i>Not available on Plastic Coated or Projection Handrims</i>	

Handrim Mount

Available only with Spinergy, Mag or Maxx Spoke and Aluminum Anodized or Plastic Coated Handrims

Close Mount	\$31
-------------	------

Axles

Axles

No selection needed with One Arm Drive

Quick Release	N/C
Quad Quick Release	\$185
<i>Not available with Mag or Spinergy wheels</i>	

Axle Position

Rear Wheel Camber

0°	2°	4°	N/C
----	----	----	-----

Center of Gravity

Rear	Center	Forward	N/C
------	--------	---------	-----

No selection needed with Amputee Axle Plate

Armrests

Armrest Type	Length	Price
Height Adjustable T-Arm Tall (11.5"-15") E0973	Desk	\$323
	Full	
Height Adjustable T-Arm (8.5"-13.5") E0973	Desk	\$323
	Full	
Height Adjustable T-Arm Low (6.5"-10.5") E0973 ¹	Desk	\$323
	Full	
Height Adjustable Flip Back T-Arm Tall (11.5"-15") E0973	Desk	\$338
	Full	
Height Adjustable Flip Back T-Arm (9.5"-13.5") E0973	Desk	\$338
	Full	

¹ If requested configuration is unavailable, Standard Height T-Arm will be used

² Not available with Seat Depths less than 16". Not Available with 15" Rear Seat Height and 24" Wheels.

³ Not available with Angle Adjustable Back Posts

⁴ Not available with Arm Pads

Armrest Type	Length	Price
Flip back (10.5") ^{2,3}	Desk	\$436
	Full	
Swing Away (10.5"-12.5") E0973 ^{3,4}		\$241
Omit Armrests		N/C

Arm Pads

Not available with Swing Away Armrests

Waterfall Arm Pad	\$26
<i>Available only with Height Adjustable T-Arm Tall, Height Adjustable T-Arm, and Height Adjustable T-Arm Low</i>	



Phone (800) 981-1540 | Fax (715) 254-0996
sales@kimobility.com | www.kimobility.com

CATALYST 5Ti

Accessories

Side Guards

Available only with Swing Away Armrests

Adult Composite Side Guards (K0108)	\$154
Pediatric Composite Side Guards (K0108)	\$154
Aluminum Straight - Black Matte (K0108)	\$210
Aluminum Adjustable Fender w/Foam Pad - Black Matte (K0108)	\$302
Carbon Fiber Adjustable Fender (K0108)	\$513

Anti-Tip

Rear Anti-Tips (E0971)	\$165
Swing-In Anti-Tip Single - Left (E0971)	\$164
<i>Not available with 350 lb Weight Capacity, Reclining Back Post or Depth Adjustable Rear Frame. Not available with left side O2 Holder.</i>	
Swing-In Anti-Tip Single - Right (E0971)	\$164
<i>Not available with 350 lb Weight Capacity, Reclining Back Post or Depth Adjustable Rear Frame. Not available with right side O2 Holder.</i>	
Swing-In Anti-Tip Pair (E0971)	\$328
<i>Not available with 350 lb Weight Capacity, Reclining Back Post or Depth Adjustable Rear Frame. Not available with Cane and Crutch Holder or O2 Holder.</i>	
Omit	N/C

Pelvic Positioning Belts

1 1/2" Auto Buckle Pelvic Positioning Belt (E0978)	\$87
1 1/2" Padded Auto Buckle Pelvic Positioning Belt (E0978)	\$179
2" Airline Buckle Pelvic Positioning Belt (E0978)	\$118
Small Bodypoint Evoflex (16-24") (E0978)	\$200
Medium Bodypoint Evoflex (19-29") (E0978)	\$200
Large Bodypoint Evoflex (25-37") (E0978)	\$200
Belt Mounting Band Clamp	\$36

O2 Holder

O2 Holder - Left Side (E2208)	\$313
O2 Holder - Right Side (E2208)	\$313

Other Accessories

Tool Kit	\$46
Removable Underseat Pouch	\$51
Pediatric Ki Mobility Backpack	\$51
Ki Mobility Backpack	\$80
Luggage Carrier	\$87
<i>Not available with PRO ELR or Residual Limb Support</i>	
Black Neoprene Impact Guards	\$97
Clear Spoke Guards (K0065)	\$123
<i>Not available with Mag wheel</i>	
Black Spoke Guards (K0065)	\$123
<i>Not available with 20" wheel size or Mag wheel</i>	
Cane and Crutch Holder	\$210
IV Pole	\$292
<i>Available only with Standard Straight Back Posts</i>	

Catalyst Front Seat Height Matrix																																		
4" Caster				5" Caster				6" Caster				7" Caster				8" Caster																		
FSH	Super Low	Stem	Fork	FSH	Super Low	Stem	Fork	FSH	Super Low	Stem	Fork	FSH	Super Low	Stem	Fork	FSH	Super Low	Stem	Fork															
15	13.5	Std	4"	16	14.5	Std	4"	17	15.5	Std	5"	18	16.5	Std	6"	19	17.5	Std	6"															
15.5	14			16.5	15.5	3/4		17.5	16.5	3/4		18.5	17			19.5	18.5	3/4																
15.5 ¹	14 ¹	17.5		16	1 1/2	18.5		17	1 1/2	18.5		17.5	20.5	n/a		1 1/2																		
16	14.5	3/4		16	14.5	Std	5"	17	15.5	Std	6"	19	18	3/4		19	17.5	Std	7"															
16.5 ¹	15 ¹			16.5	15			17.5	16			19.5	18			19.5	18																	
17	15.5	16.5		15.5	18	16.5		20	18.5			20	18.5																					
15.5	14	Std		17	16	17.5		16.5	18.5	17		19	17.5	20.5		19	3/4																	
16	14.5			17.5	16	18		17	19.5	18		20.5	18.5	21		n/a																		
16	15	3/4		18	16.5	18.5		17.5	19	18		20	19	21		n/a		1 1/2																
16.5	15.5			16.5	15	19		17.5	19.5	18.5		20	18.5	21	n/a																			
17	15.5	1 1/2		17	15.5	Std		6"	18	16.5		1 1/2	7"	19	18	3/4																		
17.5	16			17.5	16				19.5	18				20	19							21	n/a											
16.5	15	Std	17.5	16	18	16.5			20.5	n/a				20	18.5	20.5						19												
17	15.5		18	16.5	19.5	18			21	n/a																								
17.5	16	1 1/2	18	16.5	3/4	19			17.5	1 1/2		20		18.5	21	n/a																		
18	16.5		19	17.5		20			19			21		n/a																				
18.5	17	1 1/2	18	16.5	Std	7"		19.5	18.5	3/4																								
18	16.5		18	16.5				20	19																									
18.5	17	19	17.5	19.5	18			20	18.5																									
18.5	17.5	3/4	18.5	17.5	1 1/2			20	19	20.5												19												
19	18		19	18				21	n/a																									
19.5	18.5	1 1/2	19	18	3/4			20.5	19																									
20	19		19.5	18.5		21	n/a																											
				19.5	18	1 1/2	6" x 2" caster not available with 5" fork																											
				20	18.5																													
				20.5	19																													

¹ Not available with 4 x 1.5" Soft Roll Aluminum

Seat heights with "FSH" heading available with Standard Front Frame. Seat heights with "Super Low" heading available with Super Low Front Frame. Actual seat height may vary depending upon configuration.