



Phone (800) 981-1540 | Fax (715) 254-0996
sales@kimobility.com | www.kimobility.com

FOCUS CR

Company:

Ship To:

Account Number:

Date:

Address:

Quote: PO Number:

City:

State:

ATP:

Zip:

Phone:

Marked For:

Email/Fax:

Associated Quote/Order#:

Effective: 4/1/2024

Focus CR - Transit Included (E1161) \$3731

Focus CR - Transit Not Included (E1161) \$3385

300 lb weight capacity Tilt range: -5° anterior to 50° posterior

Frame

Frame Width

Select Frame Width within desired Frame Width Range

Frame Width Range	Frame Width				
Small (14" - 18")	14" N/C	15" N/C	16" N/C	17" N/C	18" N/C
Medium (16" - 20")	16" N/C	17" N/C	18" N/C	19" N/C	20" (E2201)
Large (18" - 22")	18" N/C	19" N/C	20" (E2201)	21" (E2201)	22" (E2201)
*FTFRW008 ¹	20" (E2201)	21" (E2201)	22" (E2201)	23" (E2201)	24" (E2202)

¹ 400 lb Weight Capacity is required

Seating Setup

Make one selection from either Option 1 at right or Option 2 below based on desired seating setup

Setup By Ki Mobility (Option 1) N/C

Select Seat Depth corresponding to desired Frame Depth Range row and Backrest Offset column. Note: All Seat Depths in a column and within the same frame depth range are achievable with future adjustment.

Setup by Frame Depth (Option 2) N/C

Select Frame Depth within desired Frame Depth Range

Option 2: Setup By Frame Depth

Frame Depth Range	Frame Depth				
Short	14" N/C	15" N/C	16" N/C	17" N/C	18" N/C
Medium	16" N/C	17" N/C	18" N/C	19" N/C	20" (E2203) \$630
Long	18" N/C	19" N/C	20" (E2203) \$630	21" (E2203) \$630	22" (E2204) \$1056


Seat Options

Growing Adjustable Aluminum Seat Pan \$231

Omit Seat Pan N/C

Option 1: Setup By Ki Mobility

Frame Depth Range	Backrest Offset A						
	-2"	-1"	Flush	+1"	+2"	+3"	+4"
	Seat Depth						
Short	16"	15"	14"				
	17"	16"	15"	14"			
	18"	17"	16"	15"	14"		
	19"	18"	17"	16"	15"	14"	
	20" ¹	19"	18"	17"	16"	15"	14"
Medium	18"	17"	16"	15"	14"		
	19"	18"	17"	16"	15"	14"	
	20" ¹	19"	18"	17"	16"	15"	14"
	21" ¹	20" ¹	19"	18"	17"	16"	15"
	22" ²	21" ¹	20" ¹	19"	18"	17"	16"
Long	20" ¹	19"	18"	17"	16"	15"	14"
	21" ¹	20" ¹	19"	18"	17"	16"	15"
	22" ²	21" ¹	20" ¹	19"	18"	17"	16"
		22" ²	21" ¹	20" ¹	19"	18"	17"
			22" ²	21" ¹	20" ¹	19"	18"

¹ Price is \$630 E2203² Price is \$1056 E2204

HCPCS Codes provided are based on code verification by the PDAC or our interpretation of the code definitions.

Proper use of HCPCS Codes does not ensure coverage or payment. For coverage information, verify the appropriate payer's coverage policy.

DCN0031.42 Pricing and specifications subject to change without notice.



Phone (800) 981-1540 | Fax (715) 254-0996
sales@kimobility.com | www.kimobility.com

FOCUS CR

Frame Options

Weight Capacity		Offset Option	
300 lbs	N/C	Offset Hanger with Tilt Stops	\$36
400 lbs (Heavy Duty)	\$451	<i>Tilt Stop included. No anterior tilt.</i>	
<i>Available only with Long Frame Depth Range and 8" x 2" or 6" x 2" Casters</i>			
Transit		Vent and Battery Tray	
Transit System	\$343	<i>Requires minimum 16" Seat Height. Not available on 14" or 15" Frame Widths, Attendant Foot Lock, 12" wheels, One Arm Drive, 5" casters</i>	
Omit Transit	N/C	Vent Tray (E1029)	\$282
		<i>Not available with Self Propel Axle Plate or Foot Tilt. Requires Anti-Tip.</i>	
		Battery Box Holder (E1029)	\$282
		<i>Available only with Dual Hand Tilt</i>	
		Vent Tray and Battery Box Holder (E1029)	\$564
		<i>Requires Anti-Tip. Available only with Dual Hand Tilt. Not available with Self Propel Axle Plate.</i>	

Tilt Mechanism

Tilt Mechanism	Power Tilt Options
<i>Power Tilt not available with Vent Tray, Battery Box Holder, or 400 lb Weight Capacity. Power Tilt includes Anterior and Posterior Tilt Stops.</i>	
Hand Tilt - Right	N/C
Hand Tilt - Left	N/C
Hand Tilt - Dual	N/C
Foot Tilt - Right	N/C
Foot Tilt - Left	N/C
Power Tilt - Right (E1002)	\$2250
Power Tilt - Left (E1002)	\$2250
	<i>Power Tilt selection required</i>
	Toggle Switch
	Pendant Control
	<i>Sent as accessory</i>
	Interface for Adaptive Switching
	<i>Sent as accessory</i>
	\$282

Backrest

Backrest	Other Backrest Options
Stroller Handle Back Post (Short 20")	N/C
Stroller Handle Back Post (Tall 24")	N/C
Fixed Height Back with Adjustable Handle (Short 22")	N/C
<i>Not available with Dynamic Back</i>	
Fixed Height Back with Adjustable Handle (Tall 26")	N/C
<i>Not available with Dynamic Back</i>	
Adjustable Height Back with Adjustable Handle (Short 18" - 25")	N/C
<i>Not available with O2 Holder or Dynamic Back</i>	
Adjustable Height Back with Adjustable Handle (Tall 20" - 27")	N/C
<i>20"-22" Heights not available with Dynamic Back</i>	
Reclining Back with Adjustable Handle (Short 22") (Non-Folding)	\$1025
<i>Not available with O2 Holder or Dynamic Back. If no Recline Angle and Gas Spring selections are made, 0° - 65° and Standard Spring will be used.</i>	
Reclining Back with Adjustable Handle (Tall 26") (Non-Folding)	\$1025
<i>Not available with O2 Holder or Dynamic Back. If no Recline Angle and Gas Spring selections are made, 0° - 65° and Standard Spring will be used.</i>	
Backrest Height	Reclining Back Options
<i>Required when Adjustable Height Back with Adjustable Handle Back Post is selected</i>	
18" ¹	0-65° Recline Angle (E1225)
19" ¹	10-75° Recline Angle (E1225)
20"	20-85° Recline Angle (E1226)
21"	N/C
22"	
23"	
24"	
25"	
26" ²	
27" ²	N/C
<i>1 Available only with Short</i>	
<i>2 Available only with Tall</i>	
	Gas Spring
	Standard Spring
	Heavy Duty Spring
	N/C
	N/C

Footrest

Hanger / Front Frame Bend						
Hanger	60°	70°	80°	90°		Retail
Extension Tube Mount			1			N/C
Front Mount			1	2		N/C
Pro Elevating Leg Rests E0990						\$318
Contracture Footrest Bilateral Mount ³						\$382
Contracture Footrest Center Mount ³						\$382
Omit Hangers						(\$20)

1 7" and 8" Casters may interfere with less than 5° tilt
2 Not available with 7" or 8" Casters
3 Aluminum Angle Adjustable footplates included. See image to right.

PRO ELR Cover

Available only on PRO ELR


Gel Pad Covers (PAIR) (K0108) \$425

3" Extension

Available only with Front Mount Hangers.

3" Extension N/C

3" Extension (Provided as Parts) N/C



Footrest / Extension / Seat to Footrest Length

Select Extension and Seat to Footrest Length in row corresponding to desired Footrest. Footrest Length of at least 2.5" shorter than Front Seat Height recommended. Extension Tube Hanger measurements based on 70° Hanger Bend.

Footrest	Extension	Seat to Footrest Length									
		10"	10.5"	11"	11.5"	12"	12.5"	13"	13.5"	14"	14.5"
Extension Tube Hangers (60°, 70°, 80°)	Short (10" - 13.5")	10"	10.5"	11"	11.5"	12"	12.5"	13"	13.5"	14"	14.5"
	Medium (11.5" - 15.5")	11.5"	12"	12.5"	13"	13.5"	14"	14.5"	15"	15.5"	16"
	Long (14" - 18")	14"	14.5"	15"	15.5"	16"	16.5"	17"	17.5"	18"	18.5"
Front Mount Hangers ¹	Short Front Mount Hanger	Standard Selections									Requires 3" Extension for \$35
		4" ²	4.5" ²	5"	5.5"	6"	6.5"	7"	7.5"	8"	8.5"
	Medium Front Mount Hanger	6.5"	7"	7.5"	8"	8.5"	9"	9.5"	10"	10.5"	11"
		4" ²	4.5" ²	5"	5.5"	6"	6.5"	7"	7.5"	8"	8.5"
	Long Front Mount Hanger	7.5"	8"	8.5"	9"	9.5"	10"	10.5"	11"	11.5"	12"
		5"	5.5"	6"	6.5"	7"	7.5"	8"	8.5"	9"	9.5"
Contracture Footrest	Short (10" - 13.5")	10.5"	11"	11.5"	12"	12.5"	13"	13.5"	14"	14.5"	15"
	Medium (11.5" - 15.5")	13.5"	14"	14.5"	15"	15.5"	16"	16.5"	17"	17.5"	18"
	Long (14" - 18")	16.5"	17"	17.5"	18"	18.5"	19"	19.5"	20"	20.5"	21"
PRO ELR Low Pivot Setting	PRO ELR Low Pivot Setting Short (10.5" - 15.5")	10.5" - 15.5"									
	PRO ELR Low Pivot Setting Long (13" - 19")	13" - 19"									
PRO ELR High Pivot Setting	PRO ELR High Pivot Setting Short (8.5" - 17")	8.5" - 17"									
	PRO ELR High Pivot Setting Long (11.5" - 17.5")	11.5" - 17.5"									

1 See Minimum Front Seat Height chart below

2 Available with 90° only

Minimum Front Seat Height with Front Mount Hangers						
Hanger Length	Short		Medium		Long	
Hanger Bend	70°	80/90°	70°	80/90°	70°	80°
Minimum Front Seat Height with Standard Selection	13"	13"	13"	13.5"	14.5"	15.5"
Minimum Front Seat Height with 3" Extension	13.5"	14"	15.5"	16"	17"	18"



Phone (800) 981-1540 | Fax (715) 254-0996
sales@kimobility.com | www.kimobility.com

FOCUS CR

Footplates

No selection needed with Contracture Footrest

Composite	N/C
<i>Not available with Front Mount Hanger</i>	
Composite Angle Adjustable (K0040)	\$195
<i>Not available with Front Mount Hanger or 14" Seat Width without Offset Hanger</i>	
Aluminum Angle Adjustable (K0040)	\$210
Aluminum Locking Multi-Angle Adjustable (K0108)	\$318
<i>Available only with Front Mount Hanger</i>	
One Piece Flip Up Angle Adjustable (K0040)	\$205
<i>Not available with Front Mount Hanger or PRO ELR</i>	

Hanger Release

No selection needed with Contracture Footrest. If no selection is made 4-Way - In line Hanger Release will be used.

Classic	N/C
4-Way - In-Line	N/C
4-Way - Out	N/C

Heel Loops

Not available with One Piece Flip Up Angle Adjustable Footplate

Heel Loops (E0951)	\$46
Adjustable Heel Loops (E0951)	\$56
Omit Heel Loops	N/C

Calf Straps

Velcro Adjustable (K0038)	\$41
Padded Velcro Adjustable (K0038)	\$72
Bodypoint Padded AeroMesh (K0038)	\$128

Residual Limb Support

	Limb Support	Price	Omit	Cover		Pad Size					
Left	Left Limb Support (E1020)	\$431	Omit Left Hanger	Smooth Cover Left	Air Mesh Cover Left	5X8 Pad Left Side	7X10 Pad Left Side	9x10 Pad Left Side	9x14 Pad Left Side	11x10 Pad Left Side	11x14 Pad Left Side
Right	Right Limb Support (E1020)	\$431	Omit Right Hanger	Smooth Cover Right	Air Mesh Cover Right	5X8 Pad Right Side	7X10 Pad Right Side	9x10 Pad Right Side	9x14 Pad Right Side	11x10 Pad Right Side	11x14 Pad Right Side

Seat Height, Casters and Rear Wheels

Select Seat Height corresponding to desired Rear Wheel and Caster Size

Rear Wheel Size	Caster Size					
	5"	6"	6" x 2"	7"	8"	8" x 2"
	Seat Height					
12" with Pneumatic or Pneumatic with Airless Insert Tire ³	13.5"					
	14"					
	14.5"	14.5"	14.5"		17"	
	17.5"	17.5"	17.5"	17.5"	17.5"	17.5"
	18"	18"	18"	18"	18"	18"
	18.5"	18.5"	18.5"	18.5"	18.5"	18.5"
	19"	19"	19"	19"	19"	19"
	19.5"	19.5"	19.5"	19.5"	19.5"	19.5"
12" with Full Poly ³	16.5"	16.5"	16.5"			
	17"	17"	17"	17"	17"	17"
	17.5"	17.5"	17.5"	17.5"	17.5"	17.5"
	18"	18"	18"	18"	18"	18"
	18.5"	18.5"	18.5"	18.5"	18.5"	18.5"
	19"	19"	19"	19"		
12" with Drum Brake	13"					
	13.5"					
	14"	14"	14"			
	14.5"	14.5"	14.5"			
	17.5"	17.5"	17.5"	17.5"	17.5"	17.5"
	18"	18"	18"	18"	18"	18"
	18.5"	18.5"	18.5"	18.5"	18.5"	18.5"
	19"	19"	19"	19"	19"	19"
	19.5"	19.5"	19.5"	19.5"	19.5"	19.5"
	20" ²		20" ²	20" ²	20" ²	
16" with Pneumatic, Pneumatic with Airless Insert or Full Poly Tire ⁴	13"					
	13.5"					
	14"	14"	14"			
	14.5"	14.5"	14.5"			
	15"	15"	15"	15"		
	15.5"	15.5"	15.5"	15.5"		
	16"	16"	16"	16"	16"	16"
	16.5" ¹	16.5" ¹	16.5" ¹	16.5" ¹		
	19"	19"	19"	19" ²	19"	19"
	19.5"	19.5"	19.5"	19.5"	19.5"	19.5"
	20"	20"	20"	20"	20"	20"
	20.5"	20.5"	20.5"	20.5"	20.5"	20.5"
16" with Low Poly ⁴	13"					
	13.5"					
	14"	14"	14"			
	14.5"	14.5"	14.5"			
	15"	15"	15"	15"		
	18"	18"	18"	18"	18"	18"
	18.5"	18.5"	18.5"	18.5"	18.5"	18.5"
	19"	19"	19"	19"	19"	19"
	19.5"	19.5"	19.5"	19.5"	19.5"	19.5"
	20"	20"	20"	20"	20"	20"
			20.5"	20.5"	20.5"	

Rear Wheel Size	Caster Size					
	5"	6"	6" x 2"	7"	8"	8" x 2"
	Seat Height					
16" with Drum ⁴	13.5"					
	14"	14"	14"			
	14.5"	14.5"	14.5"			
	15"	15"	15"	15"		
	15.5"	15.5"	15.5"	15.5"		
	16"	16"	16"			16"
		18.5"	18.5"	18.5"	18.5"	18.5"
	19"	19"	19"	19"	19"	19"
	19.5"	19.5"	19.5"	19.5"	19.5"	19.5"
	20"	20"	20"	20"	20"	20"
	20.5"	20.5"	20.5"	20.5"	20.5"	20.5"
20"	15"	15"	15"	15"		
	15.5"	15.5"	15.5"	15.5"		
	16"	16"	16"	16"	16"	16"
	16.5"	16.5"	16.5"	16.5"	16.5"	16.5"
	17"	17"	17"	17"	17"	17"
	17.5"	17.5"	17.5"	17.5"	17.5"	17.5"
	18"	18"	18"	18"	18"	18"
	18.5" ²	18.5" ²	18.5" ²	18.5" ²		
22"	16"	16"	16"	16"	16"	
	16.5"	16.5"	16.5"	16.5"	16.5"	16.5"
	17"	17"	17"	17"	17"	17"
	17.5"	17.5"	17.5"	17.5"	17.5"	17.5"
	18"	18"	18"	18"	18"	18"
	18.5"	18.5"	18.5"	18.5"	18.5"	18.5"
	19"	19"	19"	19"	19"	19"
	19.5"	19.5"	19.5"	19.5"	19.5"	
24"	16.5" ⁵	16.5" ^{1.5}	16.5" ^{1.5}			
	17" ⁵	17" ⁵	17" ⁵	17" ⁵	17" ⁵	17" ⁵
	17.5" ⁶	17.5" ⁶	17.5" ⁶	17.5"	17.5"	17.5" ⁶
	18"	18"	18"	18"	18"	18"
	18.5"	18.5"	18.5"	18.5"	18.5"	18.5"
	19"	19"	19"	19"	19"	19"
	19.5"	19.5"	19.5"	19.5"	19.5"	19.5"
	20"	20"	20"	20"	20"	20" ²

¹ Available only with Full Poly tire

² Not available with Full Poly tire

³ Not available with Vent and Battery Tray combination

⁴ Seat heights higher than 16.5" with 16" wheel size not available with Vent and Battery Tray combination

⁵ Not available with Drum Brake

⁶ Not available with Drum Brake in combination with Pneumatic tires



Phone (800) 981-1540 | Fax (715) 254-0996
sales@kimobility.com | www.kimobility.com

FOCUS CR

Casters and Rear Wheels

One Arm Drive

One Arm Drive - Quick Release (E0958) \$1379

*Available only with 24" Spoke wheels. Spoke wheel style is exclusive to OAD.
Not available with Vent Tray, Battery Box Holder, Self Propel Axle Plate or
Angle Adjustable Locking Extendable Flip Up. Standard Handrim will be
Aluminum Anodized. Axle is Quick Release.*

Left Side Drive N/C
Right Side Drive N/C
Aluminum Drive Handrim N/C
Poly Coated Drive Handrim \$205

Casters

		Caster Size				
		Price	5"	6"	7"	8"
Caster Type	1" Poly	N/C				
	1" Poly Aluminum	\$118				
	1" Pneumatic	\$118				
	1.5" Poly	\$95				
	1.5" Soft Roll Aluminum	\$190				
	2" Poly	\$97				
	2" Pneumatic w/Foam Insert	\$246				

Axle

No selection needed with One Arm Drive or Drum Brake

Fixed N/C
Quick Release N/C

Self Propel Axle Plate

Self Propel Axle Plate \$103
Not available with Vent Tray or 12" or 16" wheel sizes

Wheel Locks

Wheel locks

Push to Lock N/C
Pull to Lock N/C
Not available with 12" or 16" Full Poly tires
Attendant Foot Lock \$87
Combination Attendant Foot Lock \$159
Available only with 12" or 16" wheel sizes. Includes Extension Handles. Push to Lock only.
Drum Brake (E2228) \$595
Available only with 12" Mag, 16" Mag, or MAXX Spoke wheels. Not available with 16" Low Poly tire, Reclining Backrest, One Arm Drive or Self Propel Axle Plate. Axle is Quick Release.
Dual - Drum Brake with Push to Lock (E2228) \$702
Available only with 12" Mag, 16" Mag, or MAXX Spoke wheels. Not available with 16" Low Poly tire, Reclining Backrest, One Arm Drive or Self Propel Axle Plate. Axle is Quick Release.
Dual - Drum Brake with Pull to Lock (E2228) \$702
Available only with 12" Mag, 16" Mag, or MAXX Spoke wheels. Not available with 16" Low Poly tire, Reclining Backrest, One Arm Drive or Self Propel Axle Plate. Axle is Quick Release.
Omit Wheel Locks N/C

Rear Wheel

No selection needed with One Arm Drive

		Rear Wheel Size					
		Price	12"	16"	20"	22"	24"
Rear Wheel Type	Mag	N/C					
	Spoke	\$205					

Rear Tire Type

Full Poly (E2218) N/C
Pneumatic (E2211,E2212) N/C
Pneumatic with Airless Insert (E2211, E2213) \$97
Low Poly (E2218) N/C
Available only with 16" wheel size

Handrims

Handrims not available on 12" or 16" wheel sizes. No selection needed with One Arm Drive

Standard Handrims

Aluminum Anodized N/C
Plastic Coated \$123
Projection (8 Vertical Projections) \$349
Available only with full profile tires in combination with 22", 24", and 26" MAXX Spoke wheels, or 22" and 24" Mag wheels
Omit Handrim N/C

Handrim Options

Not available on Plastic Coated or Projection Handrims

3M Micro Grip \$72

Armrests

Armrest Type	Length	Price
Height Adjustable T-Arm Tall (11.5"-15") E0973	Desk	\$323
	Full	
Height Adjustable T-Arm (8.5"-13.5") E0973	Desk	\$323
	Full	
Height Adjustable T-Arm Low (6.5"-10.5") E0973 ¹	Desk	\$323
	Full	

¹ If requested configuration is unavailable, Standard Height T-Arm will be used. Not available with 14" seat depth and Offset Hangers"

² Not available with IV Pole

³ Not available with Dynamic Back

⁴ Minimum armrest height with Reclining Backrest is 12"

Armrest Type	Length	Price
Height Adj. Dual Post Flip Back (9.5"-14.5") 2 ²	Desk	\$320
	Full	
Angle Adjustable Locking Extendable Flip Up 3,4 ^{3,4}	Desk	\$297
	Full	
Omit Armrests		N/C

Arm Pads

Standard Arm Pad	N/C
Waterfall Arm Pad	\$26
Not available with Angle Adjustable Locking Extendable Flip Up Armrests	
Foam Grip	N/C
Available only with Angle Adjustable Locking Extendable Flip Up Armrest	

Accessories

Removable Side Guard

Composite Side Guards (K0108)	\$149
Not available with T-Arm or Dual Post Flip Back Armrests	

Anti-Tip

Rear Anti-Tips (E0971)	\$168
Required with Reclining Backrest or Vent Tray	
Omit	N/C

Pelvic Positioning Belt

1 1/2" Auto Buckle Pelvic Positioning Belt (E0978)	\$87
1 1/2" Padded Auto Buckle Pelvic Positioning Belt (E0978)	\$179
2" Airline Buckle Pelvic Positioning Belt (E0978)	\$118
Small Bodypoint Evoflex (16-24") (E0978)	\$200
Medium Bodypoint Evoflex (19-29") (E0978)	\$200
Large Bodypoint Evoflex (25-37") (E0978)	\$200

Belt Mounting Kit

Ki Mobility Belt Mounting Brackets	\$51
Belt Mounting Band Clamp	\$36

U.S. Military Service Patches


Army	Navy	Marines	Air Force	
Coast Guard	US Purple Heart			\$92
Combo Army/Purple Heart		Combo Navy/Purple Heart		
Combo Marines/Purple Heart		Combo Air Force/Purple Heart		
Combo Coast Guard/Purple Heart				\$149

O2 Holder

O2 Holder (E2208)	\$179
Not available with Reclining Backrest or Dynamic Back	

Other Accessories

Tool Kit	\$62
Removable Underseat Pouch	\$56
Ki Mobility Backpack	\$72
Black Neoprene Impact Guards	\$97
Tilt Stop (Pair)	\$31
Included with Offset Hanger	
IV Pole	\$292
Canopy with Hardware (K0108)	\$297
Not available on frame widths greater than 18". See image below.	





Phone (800) 981-1540 | Fax (715) 254-0996
sales@kimobility.com | www.kimobility.com



Frame Color

Gloss Finish		Matte Finish	
Black	N/C	Military Green	N/C
Blue Skies	N/C	Charcoal Gray	N/C
Bubblegum Pink	N/C	Emerald Green	N/C
Candy Blue	N/C	Iridescent Burgundy	N/C
Candy Purple	N/C	Iridescent Navy	N/C
Candy Red	N/C	Retro Red	N/C
Copper	N/C	Royal Blue	N/C
Granny Smith	N/C	Shadow Black	N/C
Metallic Orange	N/C	Smooth White	N/C
Yellow	N/C		