



Phone (800) 981-1540 | Fax (715) 254-0996
 sales@kimobility.com | www.kimobility.com

Effective 05/23/2017

Company: _____

Ship To: _____

Account Number: _____ Date: _____

Address: _____

Quote PO Number: _____

City: _____ State: _____

RTS: _____

Zip: _____ Phone: _____

Marked For: _____

Email/Fax: _____

Focus CR E1161 _____ **\$2980**
 300 lb weight capacity. Tilt range: -5° anterior to 50° posterior.

HCPCS codes provided are based on code verification by the PDAC or our interpretation of the code definitions.

Proper use of HCPCS codes does not ensure coverage or payment. For coverage information, verify the appropriate payer's coverage policy.

Frame Width

Frame Width Range / Frame Width					
<i>Select Frame Width within desired Frame Width Range</i>					
Frame Width Range	Frame Width				
Small 14" - 18"	<input type="checkbox"/> 14 N/C	<input type="checkbox"/> 15 N/C	<input type="checkbox"/> 16 N/C	<input type="checkbox"/> 17 N/C	<input type="checkbox"/> 18 N/C
Medium 16" - 20"	<input type="checkbox"/> 16 N/C	<input type="checkbox"/> 17 N/C	<input type="checkbox"/> 18 N/C	<input type="checkbox"/> 19 N/C	<input type="checkbox"/> 20 \$470 E2201
Large 18" - 22"	<input type="checkbox"/> 18 N/C	<input type="checkbox"/> 19 N/C	<input type="checkbox"/> 20 \$470 E2201	<input type="checkbox"/> 21 \$470 E2201	<input type="checkbox"/> 22 \$470 E2201

OPTION 1: Setup By Ki Mobility

Select Seat Depth corresponding to desired Frame Depth Range row and Backrest Offset column.

Note: All Seat Depths in a column and within the same frame depth range are achievable with future adjustment.

Frame Depth Range	Backrest Offset A						
	-2	-1	Flush	+1	+2	+3	+4
Short	Seat Depth						
	<input type="checkbox"/> 16	<input type="checkbox"/> 15	<input type="checkbox"/> 14				
	<input type="checkbox"/> 17	<input type="checkbox"/> 16	<input type="checkbox"/> 15	<input type="checkbox"/> 14			
	<input type="checkbox"/> 18	<input type="checkbox"/> 17	<input type="checkbox"/> 16	<input type="checkbox"/> 15	<input type="checkbox"/> 14		
	<input type="checkbox"/> 19	<input type="checkbox"/> 18	<input type="checkbox"/> 17	<input type="checkbox"/> 16	<input type="checkbox"/> 15	<input type="checkbox"/> 14	
Medium	<input type="checkbox"/> 20 ¹	<input type="checkbox"/> 19	<input type="checkbox"/> 18	<input type="checkbox"/> 17	<input type="checkbox"/> 16	<input type="checkbox"/> 15	<input type="checkbox"/> 14
	<input type="checkbox"/> 21 ¹	<input type="checkbox"/> 20 ¹	<input type="checkbox"/> 19	<input type="checkbox"/> 18	<input type="checkbox"/> 17	<input type="checkbox"/> 16	<input type="checkbox"/> 15
	<input type="checkbox"/> 22 ²	<input type="checkbox"/> 21 ¹	<input type="checkbox"/> 20 ¹	<input type="checkbox"/> 19	<input type="checkbox"/> 18	<input type="checkbox"/> 17	<input type="checkbox"/> 16
	<input type="checkbox"/> 23 ¹	<input type="checkbox"/> 22 ²	<input type="checkbox"/> 21 ¹	<input type="checkbox"/> 20 ¹	<input type="checkbox"/> 19	<input type="checkbox"/> 18	<input type="checkbox"/> 17
	<input type="checkbox"/> 24 ¹	<input type="checkbox"/> 23 ¹	<input type="checkbox"/> 22 ²	<input type="checkbox"/> 21 ¹	<input type="checkbox"/> 20 ¹	<input type="checkbox"/> 19	<input type="checkbox"/> 18
Long	<input type="checkbox"/> 25 ¹	<input type="checkbox"/> 24 ¹	<input type="checkbox"/> 23 ¹	<input type="checkbox"/> 22 ²	<input type="checkbox"/> 21 ¹	<input type="checkbox"/> 20 ¹	<input type="checkbox"/> 19
	<input type="checkbox"/> 26 ¹	<input type="checkbox"/> 25 ¹	<input type="checkbox"/> 24 ¹	<input type="checkbox"/> 23 ¹	<input type="checkbox"/> 22 ²	<input type="checkbox"/> 21 ¹	<input type="checkbox"/> 20 ¹
	<input type="checkbox"/> 27 ¹	<input type="checkbox"/> 26 ¹	<input type="checkbox"/> 25 ¹	<input type="checkbox"/> 24 ¹	<input type="checkbox"/> 23 ¹	<input type="checkbox"/> 22 ²	<input type="checkbox"/> 21 ¹
	<input type="checkbox"/> 28 ¹	<input type="checkbox"/> 27 ¹	<input type="checkbox"/> 26 ¹	<input type="checkbox"/> 25 ¹	<input type="checkbox"/> 24 ¹	<input type="checkbox"/> 23 ¹	<input type="checkbox"/> 22 ²
	<input type="checkbox"/> 29 ¹	<input type="checkbox"/> 28 ¹	<input type="checkbox"/> 27 ¹	<input type="checkbox"/> 26 ¹	<input type="checkbox"/> 25 ¹	<input type="checkbox"/> 24 ¹	<input type="checkbox"/> 23 ¹

¹ Price is \$605 E2203 ² Price is \$1025 E2204


Seating Setup

Make one selection from either Option 1 at right or Option 2 below based on desired seating setup

OPTION 2: Setup By Frame Depth

Select Frame Depth within desired Frame Depth Range

Frame Depth Range	Frame Depth				
Short 14" - 18"	<input type="checkbox"/> 14 N/C	<input type="checkbox"/> 15 N/C	<input type="checkbox"/> 16 N/C	<input type="checkbox"/> 17 N/C	<input type="checkbox"/> 18 N/C
Medium 16" - 20"	<input type="checkbox"/> 16 N/C	<input type="checkbox"/> 17 N/C	<input type="checkbox"/> 18 N/C	<input type="checkbox"/> 19 N/C	<input type="checkbox"/> 20 \$605 E2203
Long 18" - 22"	<input type="checkbox"/> 18 N/C	<input type="checkbox"/> 19 N/C	<input type="checkbox"/> 20 \$605 E2203	<input type="checkbox"/> 21 \$605 E2203	<input type="checkbox"/> 22 \$1025 E2204



Frame Options

Heavy Duty Option

- 400 lb Weight Capacity _____ \$425
Available only with Long Frame Depth Range and 8" x 2" or 6" x 2" Casters

Transit Option

- Transit System K0108 _____ \$320

Offset Hanger Option

- Offset Hanger _____ \$30
Tilt Stop included. No anterior tilt.

Vent and Battery Tray

Not available on 14" or 15" Frame Widths, Foot Tilt Mechanism, Attendant Foot Lock or 5" casters. Requires minimum 16" Seat Height.

- Vent Tray E1029 _____ \$250
Requires Anti-Tip
- Battery Box Holder E1029 _____ \$250

Tilt Mechanism

Hand

- Hand Tilt - Right _____ N/C
- Hand Tilt - Left _____ N/C
- Hand Tilt - Dual _____ N/C

Foot

- Foot Tilt - Right _____ N/C
- Foot Tilt - Left _____ N/C

Power Tilt

Not available with Vent Tray, Battery Box Holder or Heavy Duty Option

- Power Tilt - Right _____ \$1995
- Power Tilt - Left _____ \$1995

Backrest

Back Post / Height Range / Back Height

Select Back Height in row corresponding to desired Back Post and Height Range

Back Post	Price	Height Range	Back Height							
Stroller Handle Back Post (Folding)	N/C	Short	<input type="checkbox"/> 20							
		Tall	<input type="checkbox"/> 24							
Fixed Height with Adjustable Handle ² (Folding) New!	N/C	Short	<input type="checkbox"/> 22							
		Tall	<input type="checkbox"/> 26							
Adjustable Height with Adjustable Handle (Folding)	\$25	Short ^{1,2}	<input type="checkbox"/> 18	<input type="checkbox"/> 19	<input type="checkbox"/> 20	<input type="checkbox"/> 21	<input type="checkbox"/> 22	<input type="checkbox"/> 23	<input type="checkbox"/> 24	<input type="checkbox"/> 25
		Tall	<input type="checkbox"/> 20 ²	<input type="checkbox"/> 21 ²	<input type="checkbox"/> 22 ²	<input type="checkbox"/> 23	<input type="checkbox"/> 24	<input type="checkbox"/> 25	<input type="checkbox"/> 26	<input type="checkbox"/> 27
Reclining Back with Adjustable Handle ^{1,2} (Non-folding)	\$885	Short	<input type="checkbox"/> 22							
		Tall	<input type="checkbox"/> 26							

¹ Not available with O₂ Holder

² Not available with Seating Dynamics Rocker Back

Stroller Handle Back Post Option

Available only with Stroller Handle Back Post

- Removable Stroller Handle Extension _____ \$75

Other Backrest Options

- Non Folding Back _____ N/C
Standard with Reclining Back and Seating Dynamics Rocker Back
- Seating Dynamics Rocker Back K0108 _____ \$1995
Not available with Angle Adjustable Locking Extendable Flip Up Armrests, Heavy Duty Option or Reclining Backrest. Available only with Non Folding Back.
- Adjustable Height Rigidizer Bar _____ \$75

Reclining Backrest Options

Select options below if Reclining Back is selected on chart above

- Recline Angle: 0° - 65° E1225
 10° - 75° E1225
 20° - 85° E1226

- Gas Spring: Standard Spring
 Heavy Duty Spring

Footrest

Hanger / Front Frame Bend

- 60° _____ N/C
Available only with Extension Tube Hangers
- 70° _____ N/C
- 80° _____ N/C
Not available with 7" or 8" Casters
- 90° _____ N/C
Available only with Front Mount Hangers. Not available with 7" or 8" Casters.

- Elevating Leg Rests (ELR) E0990 _____ \$295
- Pro Elevating Leg Rests (PRO ELR) E0990 **New!** _____ \$350
- Omit Hangers _____ - \$20

Footrest continued

- Footplate**
- Composite _____ N/C
Not available with Front Mount Hangers or PRO ELR
 - Composite Angle Adjustable K0040 _____ \$185
Not available with Front Mount Hangers, ELR or 14" Seat Width without Offset Hangers
 - Aluminum Angle Adjustable K0040 _____ \$190
 - One Piece Flip Up Angle Adjustable K0040 _____ \$190
Not available with Front Mount Hangers, ELR, PRO ELR or 21" or 22" Seat Widths in combination with Offset Hangers

Footrest / Extension / Seat to Footrest Length

Select Seat to Footrest Length in row corresponding to desired Footrest and Extension. Footrest Length of at least 2.5" shorter than Front Seat Height recommended. Extension Tube Hanger measurements based on 70° Hanger Bend.

Footrest	Extension	Seat to Footrest Length													
Extension Tube Hangers (60°, 70°, 80°)	Short	<input type="checkbox"/> 10	<input type="checkbox"/> 10.5	<input type="checkbox"/> 11	<input type="checkbox"/> 11.5	<input type="checkbox"/> 12	<input type="checkbox"/> 12.5	<input type="checkbox"/> 13	<input type="checkbox"/> 13.5						
	Medium	<input type="checkbox"/> 11.5	<input type="checkbox"/> 12	<input type="checkbox"/> 12.5	<input type="checkbox"/> 13	<input type="checkbox"/> 13.5	<input type="checkbox"/> 14	<input type="checkbox"/> 14.5	<input type="checkbox"/> 15	<input type="checkbox"/> 15.5					
	Long	<input type="checkbox"/> 14	<input type="checkbox"/> 14.5	<input type="checkbox"/> 15	<input type="checkbox"/> 15.5	<input type="checkbox"/> 16	<input type="checkbox"/> 16.5	<input type="checkbox"/> 17	<input type="checkbox"/> 17.5	<input type="checkbox"/> 18					
Front Mount Hangers ¹	Standard Selections											Requires 3" Extension			
	Short (70°, 80°, 90°)	<input type="checkbox"/> 5	<input type="checkbox"/> 5.5	<input type="checkbox"/> 6	<input type="checkbox"/> 6.5						<input type="checkbox"/> 8.5	<input type="checkbox"/> 9	<input type="checkbox"/> 9.5		
		<input type="checkbox"/> 7	<input type="checkbox"/> 7.5	<input type="checkbox"/> 8						<input type="checkbox"/> 10	<input type="checkbox"/> 10.5	<input type="checkbox"/> 11			
	Medium (70°, 80°, 90°)	<input type="checkbox"/> 5	<input type="checkbox"/> 5.5	<input type="checkbox"/> 6	<input type="checkbox"/> 6.5	<input type="checkbox"/> 7	<input type="checkbox"/> 7.5					<input type="checkbox"/> 10.5	<input type="checkbox"/> 11	<input type="checkbox"/> 11.5	
		<input type="checkbox"/> 8	<input type="checkbox"/> 8.5	<input type="checkbox"/> 9	<input type="checkbox"/> 9.5	<input type="checkbox"/> 10						<input type="checkbox"/> 12	<input type="checkbox"/> 12.5	<input type="checkbox"/> 13	
	Long (70°, 80°)	<input type="checkbox"/> 5	<input type="checkbox"/> 5.5	<input type="checkbox"/> 6	<input type="checkbox"/> 6.5	<input type="checkbox"/> 7	<input type="checkbox"/> 7.5	<input type="checkbox"/> 8	<input type="checkbox"/> 8.5				<input type="checkbox"/> 12.5	<input type="checkbox"/> 13	<input type="checkbox"/> 13.5
<input type="checkbox"/> 9		<input type="checkbox"/> 9.5	<input type="checkbox"/> 10	<input type="checkbox"/> 10.5	<input type="checkbox"/> 11	<input type="checkbox"/> 11.5	<input type="checkbox"/> 12				<input type="checkbox"/> 14	<input type="checkbox"/> 14.5	<input type="checkbox"/> 15		
Elevating Leg Rests (ELR)	Short	<input type="checkbox"/> 10.5	<input type="checkbox"/> 11.5	<input type="checkbox"/> 12.5	<input type="checkbox"/> 13.5	<input type="checkbox"/> 14.5									
	Long	<input type="checkbox"/> 11.5	<input type="checkbox"/> 12.5	<input type="checkbox"/> 13.5	<input type="checkbox"/> 14.5	<input type="checkbox"/> 15.5	<input type="checkbox"/> 16.5								
Extended Length ELR	Short	<input type="checkbox"/> 16.5	<input type="checkbox"/> 17.5	<input type="checkbox"/> 18.5	<input type="checkbox"/> 19.5	<input type="checkbox"/> 20.5									
	Long	<input type="checkbox"/> 16.5	<input type="checkbox"/> 17.5	<input type="checkbox"/> 18.5	<input type="checkbox"/> 19.5	<input type="checkbox"/> 20.5	<input type="checkbox"/> 21.5	<input type="checkbox"/> 22.5							
PRO ELR Low Pivot	Short	<input type="checkbox"/> 10.5 - 15.5													
	Long	<input type="checkbox"/> 13 - 19													
PRO ELR High Pivot	Short	<input type="checkbox"/> 8.5 - 17													
	Long	<input type="checkbox"/> 11.5 - 17.5													

¹ Available only with Aluminum Angle Adjustable Footplates. See Minimum Front Seat Height chart below.

Hanger Length	Short		Medium		Long	
	70°	80/90°	70°	80/90°	70°	80°
Hanger Bend	70°	80/90°	70°	80/90°	70°	80°
Minimum Front Seat Height with Standard Selection	13"	13"	13"	13.5"	14.5"	15.5"
Minimum Front Seat Height with 3" Extension	13.5"	14"	15.5"	16"	17"	18"

- Heel Loop**
- Heel Loops E0951 _____ \$30
Not available with One Piece Flip Up Angle Adjustable Footplate
 - Omit Heel Loops _____ N/C
- Calf Strap**
- Velcro Adjustable K0038 _____ \$35
 - Padded Velcro Adjustable K0038 _____ \$68
 - Bodypoint Padded AeroMesh K0038 _____ \$102

3" Extension

Available only with Front Mount Hangers. Provided as parts.

- 3" Extension for Front Mount Hanger _____ N/C

Residual Limb Support

- Residual Limb Support E1020
 - Left _____ \$395
 - Right _____ \$395
 - Omit Left Hanger ___ - \$10
 - Omit Right Hanger ___ - \$10

Cover

- Smooth
- Air Mesh

Pad Size


- | | | | |
|---------------------------------|----------------------------------|---------------------------------|----------------------------------|
| <input type="checkbox"/> 5 x 8 | <input type="checkbox"/> 9 x 14 | <input type="checkbox"/> 5 x 8 | <input type="checkbox"/> 9 x 14 |
| <input type="checkbox"/> 7 x 10 | <input type="checkbox"/> 11 x 10 | <input type="checkbox"/> 7 x 10 | <input type="checkbox"/> 11 x 10 |
| <input type="checkbox"/> 9 x 10 | <input type="checkbox"/> 11 x 14 | <input type="checkbox"/> 9 x 10 | <input type="checkbox"/> 11 x 14 |

Footrest continued

Hanger Release

If no selection is made Classic Hanger Release will be used

- Classic _____ Std
- 4-Way - In line **New!** _____ N/C
See reference A to the right. Not available with ELR or PRO ELR.
- 4-Way - Out **New!** _____ N/C
See reference B to the right



Seat Height, Casters and Rear Wheels

Seat Height / Caster Size / Rear Wheel Size

Select Seat Height corresponding to desired Rear Wheel and Caster Size

Rear Wheel Size	Caster Size					
	5"	6"	6" x 2"	7"	8"	8" x 2"
12" ⁵ with Pneumatic or Pneumatic with Airless Insert Tire	Seat Height					
	<input type="checkbox"/> 13.5					
	<input type="checkbox"/> 14					
	<input type="checkbox"/> 14.5	<input type="checkbox"/> 14.5	<input type="checkbox"/> 14.5		<input type="checkbox"/> 17	
	<input type="checkbox"/> 17.5	<input type="checkbox"/> 17.5	<input type="checkbox"/> 17.5	<input type="checkbox"/> 17.5	<input type="checkbox"/> 17.5	<input type="checkbox"/> 17.5
	<input type="checkbox"/> 18	<input type="checkbox"/> 18	<input type="checkbox"/> 18	<input type="checkbox"/> 18	<input type="checkbox"/> 18	<input type="checkbox"/> 18
	<input type="checkbox"/> 18.5	<input type="checkbox"/> 18.5	<input type="checkbox"/> 18.5	<input type="checkbox"/> 18.5	<input type="checkbox"/> 18.5	<input type="checkbox"/> 18.5 ⁷
	<input type="checkbox"/> 19	<input type="checkbox"/> 19	<input type="checkbox"/> 19	<input type="checkbox"/> 19	<input type="checkbox"/> 19	<input type="checkbox"/> 19
	<input type="checkbox"/> 19.5	<input type="checkbox"/> 19.5	<input type="checkbox"/> 19.5	<input type="checkbox"/> 19.5	<input type="checkbox"/> 19.5	<input type="checkbox"/> 19.5
	12" ⁵ with Full Poly Tire	<input type="checkbox"/> 16.5	<input type="checkbox"/> 16.5	<input type="checkbox"/> 16.5		
<input type="checkbox"/> 17		<input type="checkbox"/> 17	<input type="checkbox"/> 17	<input type="checkbox"/> 17	<input type="checkbox"/> 17	<input type="checkbox"/> 17
<input type="checkbox"/> 17.5		<input type="checkbox"/> 17.5	<input type="checkbox"/> 17.5	<input type="checkbox"/> 17.5	<input type="checkbox"/> 17.5	<input type="checkbox"/> 17.5
<input type="checkbox"/> 18		<input type="checkbox"/> 18	<input type="checkbox"/> 18	<input type="checkbox"/> 18	<input type="checkbox"/> 18	<input type="checkbox"/> 18
<input type="checkbox"/> 18.5		<input type="checkbox"/> 18.5	<input type="checkbox"/> 18.5	<input type="checkbox"/> 18.5	<input type="checkbox"/> 18.5	<input type="checkbox"/> 18.5 ⁷
<input type="checkbox"/> 19		<input type="checkbox"/> 19	<input type="checkbox"/> 19	<input type="checkbox"/> 19		
16" ⁶ with Pneumatic, Pneumatic with Airless Insert or Full Poly Tire	<input type="checkbox"/> 13					
	<input type="checkbox"/> 13.5					
	<input type="checkbox"/> 14	<input type="checkbox"/> 14	<input type="checkbox"/> 14			
	<input type="checkbox"/> 14.5	<input type="checkbox"/> 14.5	<input type="checkbox"/> 14.5			
	<input type="checkbox"/> 15	<input type="checkbox"/> 15	<input type="checkbox"/> 15	<input type="checkbox"/> 15		
	<input type="checkbox"/> 15.5	<input type="checkbox"/> 15.5	<input type="checkbox"/> 15.5	<input type="checkbox"/> 15.5		
	<input type="checkbox"/> 16	<input type="checkbox"/> 16	<input type="checkbox"/> 16	<input type="checkbox"/> 16	<input type="checkbox"/> 16	<input type="checkbox"/> 16
	<input type="checkbox"/> 16.5 ¹	<input type="checkbox"/> 16.5	<input type="checkbox"/> 16.5 ¹	<input type="checkbox"/> 16.5 ¹		<input type="checkbox"/> 16.5 ¹
	<input type="checkbox"/> 19	<input type="checkbox"/> 19	<input type="checkbox"/> 19	<input type="checkbox"/> 19 ²	<input type="checkbox"/> 19	<input type="checkbox"/> 19
	<input type="checkbox"/> 19.5	<input type="checkbox"/> 19.5	<input type="checkbox"/> 19.5	<input type="checkbox"/> 19.5	<input type="checkbox"/> 19.5	<input type="checkbox"/> 19.5
	<input type="checkbox"/> 20	<input type="checkbox"/> 20	<input type="checkbox"/> 20	<input type="checkbox"/> 20	<input type="checkbox"/> 20	<input type="checkbox"/> 20
<input type="checkbox"/> 20.5	<input type="checkbox"/> 20.5	<input type="checkbox"/> 20.5	<input type="checkbox"/> 20.5	<input type="checkbox"/> 20.5	<input type="checkbox"/> 20.5	
16" ⁶ with Low Poly Tire	<input type="checkbox"/> 13					
	<input type="checkbox"/> 13.5					
	<input type="checkbox"/> 14	<input type="checkbox"/> 14	<input type="checkbox"/> 14			
	<input type="checkbox"/> 14.5	<input type="checkbox"/> 14.5	<input type="checkbox"/> 14.5			
	<input type="checkbox"/> 15	<input type="checkbox"/> 15	<input type="checkbox"/> 15	<input type="checkbox"/> 15		
	<input type="checkbox"/> 18	<input type="checkbox"/> 18	<input type="checkbox"/> 18	<input type="checkbox"/> 18	<input type="checkbox"/> 18	<input type="checkbox"/> 18
	<input type="checkbox"/> 18.5	<input type="checkbox"/> 18.5	<input type="checkbox"/> 18.5	<input type="checkbox"/> 18.5	<input type="checkbox"/> 18.5	<input type="checkbox"/> 18.5
	<input type="checkbox"/> 19	<input type="checkbox"/> 19	<input type="checkbox"/> 19	<input type="checkbox"/> 19	<input type="checkbox"/> 19	<input type="checkbox"/> 19
	<input type="checkbox"/> 19.5	<input type="checkbox"/> 19.5	<input type="checkbox"/> 19.5	<input type="checkbox"/> 19.5	<input type="checkbox"/> 19.5	<input type="checkbox"/> 19.5
	<input type="checkbox"/> 20	<input type="checkbox"/> 20	<input type="checkbox"/> 20	<input type="checkbox"/> 20	<input type="checkbox"/> 20	<input type="checkbox"/> 20
			<input type="checkbox"/> 20.5	<input type="checkbox"/> 20.5	<input type="checkbox"/> 20.5	

Rear Wheel Size	Caster Size					
	5"	6"	6" x 2"	7"	8"	8" x 2"
16" ⁶ with Drum Brake	Seat Height					
	<input type="checkbox"/> 13.5					
	<input type="checkbox"/> 14	<input type="checkbox"/> 14	<input type="checkbox"/> 14			
	<input type="checkbox"/> 14.5	<input type="checkbox"/> 14.5	<input type="checkbox"/> 14.5			
	<input type="checkbox"/> 15	<input type="checkbox"/> 15	<input type="checkbox"/> 15	<input type="checkbox"/> 15		
	<input type="checkbox"/> 15.5	<input type="checkbox"/> 15.5	<input type="checkbox"/> 15.5	<input type="checkbox"/> 15.5		
	<input type="checkbox"/> 16	<input type="checkbox"/> 16	<input type="checkbox"/> 16			<input type="checkbox"/> 16
		<input type="checkbox"/> 18.5	<input type="checkbox"/> 18.5 ⁷	<input type="checkbox"/> 18.5	<input type="checkbox"/> 18.5	<input type="checkbox"/> 18.5 ⁷
	<input type="checkbox"/> 19	<input type="checkbox"/> 19	<input type="checkbox"/> 19	<input type="checkbox"/> 19	<input type="checkbox"/> 19	<input type="checkbox"/> 19
	<input type="checkbox"/> 19.5	<input type="checkbox"/> 19.5	<input type="checkbox"/> 19.5	<input type="checkbox"/> 19.5	<input type="checkbox"/> 19.5	<input type="checkbox"/> 19.5
20"	<input type="checkbox"/> 15	<input type="checkbox"/> 15	<input type="checkbox"/> 15	<input type="checkbox"/> 15		
	<input type="checkbox"/> 15.5	<input type="checkbox"/> 15.5	<input type="checkbox"/> 15.5	<input type="checkbox"/> 15.5		
	<input type="checkbox"/> 16	<input type="checkbox"/> 16	<input type="checkbox"/> 16	<input type="checkbox"/> 16	<input type="checkbox"/> 16	<input type="checkbox"/> 16
	<input type="checkbox"/> 16.5	<input type="checkbox"/> 16.5	<input type="checkbox"/> 16.5	<input type="checkbox"/> 16.5	<input type="checkbox"/> 16.5	<input type="checkbox"/> 16.5
	<input type="checkbox"/> 17	<input type="checkbox"/> 17	<input type="checkbox"/> 17	<input type="checkbox"/> 17	<input type="checkbox"/> 17	<input type="checkbox"/> 17
	<input type="checkbox"/> 17.5	<input type="checkbox"/> 17.5	<input type="checkbox"/> 17.5	<input type="checkbox"/> 17.5	<input type="checkbox"/> 17.5	<input type="checkbox"/> 17.5
	<input type="checkbox"/> 18	<input type="checkbox"/> 18	<input type="checkbox"/> 18	<input type="checkbox"/> 18	<input type="checkbox"/> 18	<input type="checkbox"/> 18
	<input type="checkbox"/> 18.5	<input type="checkbox"/> 18.5	<input type="checkbox"/> 18.5	<input type="checkbox"/> 18.5		
22"	<input type="checkbox"/> 16	<input type="checkbox"/> 16	<input type="checkbox"/> 16	<input type="checkbox"/> 16	<input type="checkbox"/> 16	
	<input type="checkbox"/> 16.5	<input type="checkbox"/> 16.5	<input type="checkbox"/> 16.5	<input type="checkbox"/> 16.5	<input type="checkbox"/> 16.5	<input type="checkbox"/> 16.5
	<input type="checkbox"/> 17	<input type="checkbox"/> 17	<input type="checkbox"/> 17	<input type="checkbox"/> 17	<input type="checkbox"/> 17	<input type="checkbox"/> 17
	<input type="checkbox"/> 17.5	<input type="checkbox"/> 17.5	<input type="checkbox"/> 17.5	<input type="checkbox"/> 17.5	<input type="checkbox"/> 17.5	<input type="checkbox"/> 17.5
	<input type="checkbox"/> 18	<input type="checkbox"/> 18	<input type="checkbox"/> 18	<input type="checkbox"/> 18	<input type="checkbox"/> 18	<input type="checkbox"/> 18
	<input type="checkbox"/> 18.5	<input type="checkbox"/> 18.5	<input type="checkbox"/> 18.5	<input type="checkbox"/> 18.5	<input type="checkbox"/> 18.5	<input type="checkbox"/> 18.5 ⁷
	<input type="checkbox"/> 19	<input type="checkbox"/> 19	<input type="checkbox"/> 19	<input type="checkbox"/> 19	<input type="checkbox"/> 19	<input type="checkbox"/> 19
	<input type="checkbox"/> 19.5	<input type="checkbox"/> 19.5	<input type="checkbox"/> 19.5	<input type="checkbox"/> 19.5	<input type="checkbox"/> 19.5	
24"	<input type="checkbox"/> 16.5 ¹	<input type="checkbox"/> 16.5 ¹	<input type="checkbox"/> 16.5 ¹	<input type="checkbox"/> 16.5 ³		
	<input type="checkbox"/> 17	<input type="checkbox"/> 17	<input type="checkbox"/> 17	<input type="checkbox"/> 17	<input type="checkbox"/> 17	<input type="checkbox"/> 17
	<input type="checkbox"/> 17.5	<input type="checkbox"/> 17.5	<input type="checkbox"/> 17.5	<input type="checkbox"/> 17.5	<input type="checkbox"/> 17.5	<input type="checkbox"/> 17.5
	<input type="checkbox"/> 18	<input type="checkbox"/> 18	<input type="checkbox"/> 18	<input type="checkbox"/> 18	<input type="checkbox"/> 18	<input type="checkbox"/> 18
	<input type="checkbox"/> 18.5	<input type="checkbox"/> 18.5	<input type="checkbox"/> 18.5	<input type="checkbox"/> 18.5	<input type="checkbox"/> 18.5	<input type="checkbox"/> 18.5 ⁷
	<input type="checkbox"/> 19	<input type="checkbox"/> 19	<input type="checkbox"/> 19	<input type="checkbox"/> 19	<input type="checkbox"/> 19	<input type="checkbox"/> 19
	<input type="checkbox"/> 19.5	<input type="checkbox"/> 19.5	<input type="checkbox"/> 19.5	<input type="checkbox"/> 19.5	<input type="checkbox"/> 19.5	<input type="checkbox"/> 19.5
	<input type="checkbox"/> 20	<input type="checkbox"/> 20 ⁴	<input type="checkbox"/> 20	<input type="checkbox"/> 20	<input type="checkbox"/> 20	<input type="checkbox"/> 20 ²

¹ Available only with Full Poly tire
² Not available with Full Poly tire
³ Available only with spoke wheel and Full Poly tire combination
⁴ Not available with spoke wheel and Full Poly tire combination
⁵ Not available with Vent and Battery Tray combination
⁶ Seat heights higher than 16.5" with 16" rear wheel are not available with Vent and Battery Tray combination
⁷ Not available with Heavy Duty Option

Seat Height, Casters and Rear Wheels continued

Caster Type	
<input type="checkbox"/> 1" Poly	N/C
<input type="checkbox"/> 1" Poly Aluminum	\$110
<i>Not available with 6", 7" or 8" Caster Sizes</i>	
<input type="checkbox"/> 1" Pneumatic	\$110
<i>Not available with 5" and 7" Caster Sizes</i>	
<input type="checkbox"/> 1.5" Poly	\$50
<i>Not available with 7" and 8" Caster Sizes</i>	
<input type="checkbox"/> 1.5" Soft Roll Aluminum	\$175
<i>Not available with 7" and 8" Caster Sizes</i>	
<input type="checkbox"/> 2" Poly	\$50
<i>Not available with 5", 7" or 8" Caster Sizes</i>	
<input type="checkbox"/> 2" Pneumatic w/ Foam Insert	\$225
<i>Not available with 5", 6" or 7" Caster Sizes</i>	

One Arm Drive	
<input type="checkbox"/> One Arm Drive E0958	\$1295
<i>Available only with 24" Mag wheels. Not available with Vent Tray, Battery Box Holder or Self Propel Axle Plate. Standard Handrim will be Aluminum Anodized. Axle not quick release.</i>	
Drive Handrim Side:	<input type="checkbox"/> Left <input type="checkbox"/> Right
Drive Handrim Type:	<input type="checkbox"/> Aluminum N/C <input type="checkbox"/> Poly Coated \$200

Axle	
No selection needed with One Arm Drive or Drum Brake	
<input type="checkbox"/> Fixed	N/C
<input type="checkbox"/> Quick Release	N/C

Self Propel Axle Plate	
<input type="checkbox"/> Self Propel Axle Plate	\$95
<i>Not available with Vent Tray or 12" or 16" rear wheels</i>	

Wheel Locks

Wheel Lock	
<input type="checkbox"/> Push to Lock	N/C
<input type="checkbox"/> Pull to Lock	N/C
<i>Not available with 12" or 16" Full Poly tires</i>	
<input type="checkbox"/> Attendant Foot Lock	\$75
<input type="checkbox"/> Drum Brake E2228	\$1200
<i>Not available with Hand Tilt Mechanism or Reclining Backrest. Axle is Quick Release.</i>	
<input type="checkbox"/> Omit Wheel Locks	N/C

Armrests

Armrest Please select Arm Pad option	
<input type="checkbox"/> Height Adjustable T-Arm (8"-13") E0973	\$265
<i>Not available with Heavy Duty option</i>	
<input type="checkbox"/> Height Adjustable T-Arm Low (7.5" - 11.5") E0973	\$265
<i>Not available with Heavy Duty option</i>	
<input type="checkbox"/> Height Adj. Dual Post Flip Back (9.5"-14.5") E0973	\$275
<input type="checkbox"/> Angle Adjustable Locking Extendable Flip Up E0973	\$265
<i>Not available with Seating Dynamics Rocker Back or Heavy Duty option. Minimum armrest height with Reclining Backrest is 12".</i>	
<input type="checkbox"/> Omit Armrests	N/C

Rear Wheel Type	
<input type="checkbox"/> Mag	N/C
<input type="checkbox"/> 24" Spoke	\$200

Rear Tire Type	
<input type="checkbox"/> Full Poly E2218	N/C
<input type="checkbox"/> Low Poly E2218 New!	N/C
<input type="checkbox"/> Pneumatic E2211, E2212	N/C
<input type="checkbox"/> Pneumatic with Airless Insert E2211, E2213	\$85

Handrim
Handrims not available on 12" or 16" rear wheels. No selection needed with One Arm Drive

Standard Handrims	
<input type="checkbox"/> Aluminum Anodized	Std
<input type="checkbox"/> Plastic Coated	\$95
<input type="checkbox"/> Projection (8 Vertical Projections)	\$295
<i>Available on 24" rear wheel</i>	
<input type="checkbox"/> Omit Handrims	N/C

Natural Fit Handrims	
<i>Not available on 20" rear wheel or 22" Mag wheels. Not recommended with Mag wheels.</i>	
<input type="checkbox"/> Natural Fit LT Regular Grip	\$445
<input type="checkbox"/> Natural Fit Standard Regular Grip	\$445
<input type="checkbox"/> Natural Fit LT Super Grip	\$495
<input type="checkbox"/> Natural Fit Standard Super Grip	\$495
<input type="checkbox"/> Natural Fit LT Omit Grip	\$345
<input type="checkbox"/> Natural Fit Standard Omit Grip	\$345

Handrim Option	
<i>Not available on Plastic Coated or Projection Handrims</i>	
<input type="checkbox"/> 3M Micro Grip	\$65

Wheel Lock continued	
<input type="checkbox"/> Combination Attendant Foot Lock	\$145
<i>Available only with 12" or 16" rear wheels. Includes Extension Handles. Push to Lock only.</i>	
<input type="checkbox"/> Dual - Drum Brake	\$1270
<i>Not available with Hand Tilt Mechanism or Reclining Backrest. Axle is Quick Release.</i>	
<input type="checkbox"/> Push to Lock <input type="checkbox"/> Pull to Lock	

Wheel Lock Options	
<input type="checkbox"/> Extension Handles E0961	\$70
<i>Available only with Push to Lock or Pull to Lock Wheel Locks.</i>	

Arm Pad			
Standard Arm Pad	N/C	<input type="checkbox"/> Desk	<input type="checkbox"/> Full
Waterfall Arm Pad	\$25	<input type="checkbox"/> Desk	<input type="checkbox"/> Full
<i>Not available with Angle Adjustable Locking Extendable Flip Up Armrests</i>			
Foam Grip	N/C	<input type="checkbox"/> Desk	<input type="checkbox"/> Full
<i>Available only with Angle Adjustable Locking Extendable Flip Up Armrest</i>			

Accessories

Removable Side Guard

- Composite Side Guards K0108 _____ \$140
Not available with T-Arm or Dual Post Flip Back Armrests

Anti-Tip

- Rear Anti-Tips E0971 _____ \$140
Required with Reclining Backrest and Vent Tray

Positioning Belt

- 1 1/2" Auto Buckle E0978 _____ \$70
 1 1/2" Padded Auto Buckle E0978 _____ \$155
 2" Airline Buckle E0978 _____ \$110
 Bodypoint Evoflex - Pelvic Stabilizer E0978 _____ \$183
 Size: Small (16-24") Medium (19-29") Large (25-37")

Belt Mounting Kit

- Ki Mobility Belt Mounting Brackets _____ \$50
 Bodypoint 1" Band Clamp Belt Mounting Kit _____ \$105

Seat Pan

- Growing Adjustable Aluminum Seat Pan E2231 _____ \$200

O₂ Holder

- O₂ Holder E2208 _____ \$170
Not available with Reclining Backrest or Rocker Back

Other Accessories

- Tool Kit _____ \$55
 Removable Underseat Pouch _____ \$50
 Ki Mobility Backpack _____ \$65
 Black Neoprene Impact Guards _____ \$85
 Clear Spoke Guards K0065 _____ \$115
Available only with Spoke wheels
 Black Spoke Guards K0065 _____ \$115
Available only with Spoke wheels
 Tilt Stop (Pair) K0108 _____ \$30
Included with Offset Hangers
 IV Pole _____ \$265

Frame Color

Gloss Finish

- Black _____ N/C
 Blue Skies _____ N/C
 Bubblegum Pink _____ N/C
 Candy Blue _____ N/C
 Candy Purple _____ N/C
 Candy Red _____ N/C
 Copper _____ N/C
 Fuchsia _____ N/C
 Granny Smith _____ N/C
 Yellow _____ N/C

Matte Finish

- Army Green _____ N/C
 Charcoal Gray _____ N/C
 Emerald Green _____ N/C
 Iridescent Burgundy _____ N/C
 Iridescent Navy _____ N/C
 Retro Red _____ N/C
 Royal Blue _____ N/C
 Shadow Black _____ N/C
 Smooth White _____ N/C

Cushions

Matrix Seat Cushions

Cushion size is determined by chair dimensions. Certain sizes will have higher price and longer lead time; visit www.motionconcepts.com/matrix_international.html for size specifications.

STEP 1: Choose Cushion

- PS Cushion _____ \$414
Cushions 14" wide have a weight capacity of 150 lbs
- Vi Cushion _____ \$510

STEP 2: Choose Cushion Size

- Match to Chair Size
- Other Size - Provide Part Number _____

Other Matrix Products

Must provide part number

- _____

ROHO Seat Cushions

Cushion size is determined by chair dimensions. Certain sizes will have higher price and longer lead time; visit therohogroup.com for size specifications.

Cushions

For chairs under 22" wide

- HIGH PROFILE® Single Compartment (E2622) _____ \$603
- MID PROFILE™ Single Compartment (E2622) _____ \$603
- LOW PROFILE® Single Compartment (E2622) _____ \$603
- QUADTRO SELECT® High Profile® (E2624) _____ \$603
- QUADTRO SELECT® Mid Profile® (E2624) _____ \$603
- QUADTRO SELECT® Low Profile® (E2624) _____ \$603
- HYBRID ELITE® Single Comp Smart Check Ready E2622 _ \$623
- CONTOUR SELECT® (E2624) _____ \$603

X Large Cushions

For 22" wide chairs

- QUADTRO SELECT® High Profile® X Large (E2625) _____ \$670

Cushions with Heavy Duty Cover

- HIGH PROFILE® Single Compartment w/ HD Cover (E2622) _ \$629
- QUADTRO SELECT® High Profile® w/ HD Cover (E2624) _ \$629

Backs

Matrix Backs

Matrix Backs come standard with quick release mounting hardware

Elite TR Back _____ \$780
Not available with 14", 15" and 22" Seat Widths

Back Height:
 16" 18" 20"

Elite Deep Back _____ \$912
Not available with 22" Seat Width. Back Height of 10" not available with 14" Seat Width.

Back Height:
 10" 12" 16" 18" 20"

Mini Back _____ \$750
Available only with 14" and 15" Seat Widths. Comes standard with Kidabra cover.

Back Height:
 7" 10" 12" 14"

Therafin Headrest for MATRIX Back Systems

Headrest, Multi-Axis Offset Hardware with Smooth/Soft Pad (E1028, E0955) 34187 _____ \$420
 PAD SIZE: 5 x 8 5 x 10 5 x 14

Other Matrix Mounting Hardware

Must provide part number

For mounting hardware other than quick release

Other Matrix Products

Must provide part number

ROHO Back Systems (Available for chairs up to 20" wide)

STEP 1: Choose Back

AGILITY™ Min. Contour w/ Foam Only (E2611) _____ \$549
 BACKREST HEIGHT: 10" 13" 16" 18"

AGILITY™ Min. Contour w/ Foam & Lumbar Pad (E2613) _____ \$659
 BACKREST HEIGHT: 10" 13" 16" 18"

AGILITY™ Min. Contour w/ Air Only (E2620) _____ \$769
 BACKREST HEIGHT: 10" 13" 16" 18"

AGILITY™ Min. Contour w/ Air & Lumbar Pad (E2615) _____ \$769
 BACKREST HEIGHT: 10" 13" 16" 18"

AGILITY™ Mid Contour w/ Air Only (E2620) _____ \$781
 BACKREST HEIGHT: 10" 13" 16" 18" 20"

AGILITY™ Mid Contour w/ Air & Lumbar Pad (E2615) _____ \$781
 BACKREST HEIGHT: 10" 13" 16" 18" 20"

AGILITY™ Max Contour w/ Air Only (E2620) _____ \$920
 BACKREST HEIGHT: 16" 18" 20"

AGILITY™ Max Contour w/ Air & Lumbar Pad (E2615) _____ \$920
 BACKREST HEIGHT: 16" 18" 20"

STEP 2: Choose Back Hardware

Fixed _____ N/C
Not compatible with 18" and 20" tall backs

Quick Release _____ N/C

ROHO Laterals for AGILITY™ Back Systems

Sold in pairs

Swing Away with Curved Pad (E1028/E0956) _____ \$1014
 PAD SIZE: Medium Large X-Large
 4 T x 4 L 4 T x 5 L 5 T x 5 L
 OFFSET TYPE: Flush Offset 1" Offset 2" Offset

Fixed with Flat Pad (E0956) _____ \$452
 PAD SIZE: Medium Large X-Large
 4 T x 4 L 4 T x 5 L 5 T x 5 L
 OFFSET TYPE: Flush Offset

Therafin Headrest for AGILITY Back Systems

Headrest, Multi-Axis Offset Hardware with Smooth/Soft Pad (E1028, E0955) 34187 _____ \$420
 PAD SIZE: 5 x 8 5 x 10 5 x 14



Vest

- Full Shape, Comfort Strap (E0960) _____ \$167
 - Small 30589 Medium 30590 Large 30591
 - X-Large 30592 2X-Large 30593
- Full Shape, Extended Strap (E0960) _____ \$167
 - Small 30594 Medium 30595 Large 30596
 - X-Large 30597 2X-Large 30598
- Trim Shape, Comfort Strap (E0960) _____ \$169
 - Small 30618 Medium 30619 Large 30620
 - X-Large 30621 2X-Large 30622
- Trim Shape, Extended Strap (E0960) _____ \$169
 - Small 30604 Medium 30624 Large 30625
 - X-Large 30626 2X-Large 30627

Shoulder Harness

With metal clips and cam locks

- Small (E0960) 61403 _____ \$157
- Medium (E0960) 61503 _____ \$164
- Large (E0960) 61703 _____ \$170
- X-Large (E0960) 61803 _____ \$177

Hip Belt

- 4 Point Y-Style, Metal Buckle (E0978) _____ \$173
 - Small 60024 Medium 60025 Large 60026
- Dual Rear Pull, Metal Attaching Clips (E0978)
 - 1" Straps, 6.25"W x 1.75"H Pads 60161 _____ \$97
 - 1.5" Straps, 7.25"W x 2.25"H Pads 60171 _____ \$106
 - 1.5" Straps, 9.25"W x 2.5"H Pads 60181 _____ \$106
- Single Pull, Metal Buckle (E0978)
 - 1" Straps, 4.25"W x 1.75"H Pads 61001 _____ \$86
 - 1" Straps, 5.25"W x 1.75"H Pads 61101 _____ \$86
 - 1" Straps, 6.25"W x 1.75"H Pads 61201 _____ \$86
 - 1.5" Straps, 5.25"W x 2.25"H Pads 62001 _____ \$95
 - 1.5" Straps, 9.25"W x 2.5"H Pads 62201 _____ \$95
 - 2" Straps, 9.25"W x 3"H Pads 64001 _____ \$102
 - 2" Straps, 11.25"W x 3"H Pads 64101 _____ \$102

Cane Mounted Swing Away Lateral Thoracic Support

Sold as pair

- Swing Away Kits w/ 4"x6" Gel Pads (E0956, E1028) 50021 _____ \$478
- Swing Away Kits, No Pads (E0956, E1028) 50022 _____ \$309

Knee Adductor

- Flip Down Hardware (E1028) _____ \$634
Sold in pairs, one right side and one left side
 - Short Arm, Fits 3/4", 7/8" and 1" Tubing 30248
 - Long Arm, Fits 3/4", 7/8" and 1" Tubing 30284
- Knee Adductor Pad with Neoprene Cover (E0956) _____ \$117
Sold single
 - Small Pad, 4.25" x 3.25" 50490 Quantity: _____
 - Medium Pad, 6.25" x 4.25" 50491 Quantity: _____

Foot / Shoe Holder

Compatible only with Aluminum Angle Adjustable and Composite Angle Adjustable Footplates

- FootSure Ankle Support (K0108) _____ \$178
Sold in pairs
 - X-Small w/ Buckle Closure, 6"-8" Circum. 60250
 - Small w/ Buckle Closure, 7"-8.5" Circum. 60251
 - Medium w/ Buckle Closure, 8"-11" Circum. 60252
 - Large w/ Buckle Closure, 9"-12" Circum. 60253
 - X-Large w/ Buckle Closure, 12.5"-15" Circum. 60254
 - X-Small w/ Hook-N-Loop Closure, 5"-6" Circum. 60255
 - Small w/ Hook-N-Loop Closure, 6"-8" Circum. 60256
 - Med w/ Hook-N-Loop Closure, 8"-10.5" Circum. 60257
 - Large w/ Hook-N-Loop Closure, 9"-12" Circum. 60258
 - X-Lrg w/ Hook-N-Loop Closure, 12"-15" Circum. 60259
- Padded Footbox (K0108)
 - Small, Full Design, Fits 14" Chair 30111 _____ \$287
 - Medium, Full Design, Fits 16" Chair 30112 _____ \$293
 - Large, Full Design, Fits 18" Chair 30113 _____ \$298
 - Small, Split Design, Fits 14" Chair 30152 _____ \$367
 - Medium, Split Design, Fits 16" Chair 30153 _____ \$373
 - Large, Split Design, Fits 18" Chair 30154 _____ \$378
- Molded Shoe Holders (K0108) _____ \$162
Sold in pairs
 - Small with Plain Straps, 3.5" x 5.25" 30293
 - Medium with Plain Straps, 4" x 6.5" 30303
 - Large with Plain Straps, 4.75" x 8" 30309
 - X-Large with Plain Straps, 5" x 10" 30316



Half Tray

Power tilt switch may interfere with standard trays

- Standard Shape, Clear, Flip-Away (E0950, E1028) _____ \$246
 - Right Side 41310 Left Side 41311
- Standard Shape, Padded, Flip-Away (E0950, E1028) _____ \$269
 - Right Side 41350 Left Side 41353
- Trimline Shape, Padded, Flip-Away (E0950, E1028) _____ \$269
 - Right Side 41352 Left Side 41351

Clear Tray with Border and Pre-Drilled Holes

Power tilt switch may interfere with standard trays

STEP 1: Choose Tray

Chair Width	Tray Dimension (W x D)	Body Cutout (W x D)	HCPCS	Part Number	Price	
14"	20 x 18	10.5 x 6	E0950	33210	\$220	<input type="checkbox"/>
	20 x 18	11.5 x 6.5	E0950	33211	\$220	<input type="checkbox"/>
15"	21 x 18	11.5 x 6.5	E0950	33212	\$220	<input type="checkbox"/>
	21 x 18	12.5 x 7	E0950	33213	\$220	<input type="checkbox"/>
16"	22 x 21	12.5 x 7.5	E0950	33214	\$255	<input type="checkbox"/>
	22 x 21	13.5 x 8.5	E0950	33215	\$255	<input type="checkbox"/>
17"	23.5 x 22	13.5 x 8.5	E0950	33216	\$255	<input type="checkbox"/>
	23.5 x 22	14.5 x 9.5	E0950	33217	\$285	<input type="checkbox"/>
18"	24.5 x 22	14.5 x 9.5	E0950	33218	\$285	<input type="checkbox"/>
	24.5 x 22	15.5 x 10.5	E0950	33219	\$285	<input type="checkbox"/>
19"	26.5 x 22	15.5 x 10.5	E0950	33220	\$285	<input type="checkbox"/>
	26.5 x 22	16.5 x 11	E0950	33221	\$285	<input type="checkbox"/>
20"	26.5 x 22	17 x 11	E0950	33222	\$309	<input type="checkbox"/>
	26.5 x 22	18 x 11	E0950	33223	\$309	<input type="checkbox"/>
21"	27.5 x 23.5	17.5 x 11	E0950	33224	\$309	<input type="checkbox"/>
	27.5 x 23.5	19 x 11	E0950	33225	\$309	<input type="checkbox"/>
22"	28.5 x 23.5	18.5 x 11	E0950	33226	\$309	<input type="checkbox"/>
	28.5 x 23.5	20 x 11	E0950	33227	\$309	<input type="checkbox"/>

STEP 2: Choose Tray Hardware

- E-Z Lock Clamps, Tube Arm (E1028) 31072 _____ \$150
Available only with Angle Adjustable Locking Extendable Flip Up with Foam Grip Armrest. Sold in Pairs.
- E-Z Lock Clamps, Standard (E1028) 31098 _____ \$142
Not available with Angle Adjustable Locking Extendable Flip Up with Foam Grip Armrest. Sold in Pairs.

Communication Mount Hardware

Choose Device Holder AND Mount

STEP 1: What DEVICE will be used?

- Smart Phone 32359 _____ \$150
- Large Smart Phone 32354 _____ \$150
- I-Pad 32355 _____ \$184
- I-Pad Mini 32356 _____ \$178
- 7" Tablet 32357 _____ \$173
- 10" Tablet 32358 _____ \$395

STEP 2: Where does it need to MOUNT?

- Frame 32350 _____ \$512
- Armrest 32351 _____ \$395
- Tray 32352 _____ \$212

Cup Holder

- Heavy Duty 30889 _____ \$44
Compatible only with Height Adjustable Dual Post Flip Back Armrests