



Phone (800) 981-1540 | Fax (715) 254-0996  
sales@kimobility.com | www.kimobility.com

# FOCUS CRe

Company:

Ship To:

Account Number:

Date:

Address:

Quote: PO Number:

City:

State:

ATP:

Zip:

Phone:

Marked For:

Email/Fax:

Associated Quote/Order#:

Effective: 1/10/2024

Focus CRe - Transit Included (E1161) \$3606

Focus CRe - Transit Not Included (E1161) \$3260

300 lb weight capacity Tilt range: -5° anterior to 50° posterior

## Frame

### Frame Width Range / Frame Width

Select Frame Width within desired Frame Width Range

Frame Width Range	Frame Width				
Small (14" - 18")	14" N/C	15" N/C	16" N/C	17" N/C	18" N/C
Medium (16" - 20")	16" N/C	17" N/C	18" N/C	19" N/C	20" (E2201) \$533
Large (18" - 22")	18" N/C	19" N/C	20" (E2201) \$533	21" (E2201) \$533	22" (E2201) \$533
X-Large (20" - 24") <sup>1</sup>	20" (E2201) \$533	21" (E2201) \$533	22" (E2201) \$533	23" (E2201) \$743	24" (E2202) \$743

<sup>1</sup> 400 lb Weight Capacity is required

### Setup

Make one selection from either Option 1 at right or Option 2 below based on desired seating setup

Setup By Ki Mobility (Option 1) N/C

Select Seat Depth corresponding to desired Frame Depth Range row and Backrest Offset column. Note: All Seat Depths in a column and within the same frame depth range are achievable with future adjustment.

Setup by Frame Depth (Option 2) N/C


Select Frame Depth within desired Frame Depth Range

### Option 2: Setup By Frame Depth

Frame Depth Range	Frame Depth				
Short	14" N/C	15" N/C	16" N/C	17" N/C	18" N/C
Medium	16" N/C	17" N/C	18" N/C	19" N/C	20" (E2203) \$630
Long	18" N/C	19" N/C	20" (E2203) \$630	21" (E2203) \$630	22" (E2204) \$1056

### Option 1: Setup By Ki Mobility

Frame Depth Range	Backrest Offset A						
	-2"	-1"	Flush	+1"	+2"	+3"	+4"
Seat Depth							
Short	16"	15"	14"				
	17"	16"	15"	14"			
	18"	17"	16"	15"	14"		
	19"	18"	17"	16"	15"	14"	
	20" <sup>1</sup>	19"	18"	17"	16"	15"	14"
Medium	18"	17"	16"	15"	14"		
	19"	18"	17"	16"	15"	14"	
	20" <sup>1</sup>	19"	18"	17"	16"	15"	14"
	21" <sup>1</sup>	20" <sup>1</sup>	19"	18"	17"	16"	15"
	22" <sup>2</sup>	21" <sup>1</sup>	20" <sup>1</sup>	19"	18"	17"	16"
Long	20" <sup>1</sup>	19"	18"	17"	16"	15"	14"
	21" <sup>1</sup>	20" <sup>1</sup>	19"	18"	17"	16"	15"
	22" <sup>2</sup>	21" <sup>1</sup>	20" <sup>1</sup>	19"	18"	17"	16"
		22" <sup>2</sup>	21" <sup>1</sup>	20" <sup>1</sup>	19"	18"	17"
			22" <sup>2</sup>	21" <sup>1</sup>	20" <sup>1</sup>	19"	18"

<sup>1</sup> Price is \$630 E2203<sup>2</sup> Price is \$1056 E2204

### Seat Pan

Growing Adjustable Aluminum Seat Pan (E2231) \$231  
Omit Seat Pan N/C

HCPCS codes provided are based on code verification by the PDAC or our interpretation of the code definitions.

Proper use of HCPCS codes does not ensure coverage or payment. For coverage information, verify the appropriate payer's coverage policy.

DCN0049.35 Pricing and specifications subject to change without notice.



Phone (800) 981-1540 | Fax (715) 254-0996  
sales@kimobility.com | www.kimobility.com

# FOCUS CRe

## Frame Options

Weight Capacity	Offset Hanger Option
300 lbs N/C	Offset Hanger with Tilt Stops
400 lbs (Heavy Duty) \$451	Tilt Stop included. No anterior tilt. \$36
Available only with Long Frame Depth Range and 8" x 2" or 6" x 2" Casters	
Transit Option	Vent and Battery Tray
Omit Transit N/C	Requires minimum 16" Seat Height. Not available on 14" or 15" Frame Widths, Attendant Foot Lock, 12" wheels, One Arm Drive, 5" casters.
Transit System \$343	Vent Tray (E1029) \$282
	Not available with Self Propel Axle Plate or Foot Tilt. Requires Anti-Tip.
	Battery Box Holder (E1029) \$282
	Available only with Dual Hand Tilt
	Vent Tray and Battery Box Holder (E1029) \$564
	Requires Anti-Tip. Available only with Dual Hand Tilt. Not available with Self Propel Axle Plate.

## Tilt Mechanism

Tilt Mechanism	Power Tilt Options
Power Tilt not available with Vent Tray, Battery Box Holder, or 400 lb Weight Capacity. Power Tilt includes Anterior and Posterior Tilt Stops.	Power Tilt selection required
Hand Tilt - Right N/C	Toggle Switch N/C
Hand Tilt - Left N/C	Pendant Control N/C
Hand Tilt - Dual N/C	Sent as accessory
Foot Tilt - Right N/C	Interface for Adaptive Switching \$282
Foot Tilt - Left N/C	Sent as accessory
Power Tilt - Right (E1002) \$2250	
Power Tilt - Left (E1002) \$2250	

## Backrest

Backrest	Reclining Back Options
Stroller Handle Back Post (Short 20") N/C	Required when Reclining Back is selected
Stroller Handle Back Post (Tall 24") N/C	Recline Angle
Fixed Height Back with Adjustable Handle (Short 22") N/C	0-65° Recline Angle (E1225) N/C
Not available with Dynamic Back	10-75° Recline Angle (E1225) N/C
Fixed Height Back with Adjustable Handle (Tall 26") N/C	20-85° Recline Angle (E1226) N/C
Not available with Dynamic Back	
Adjustable Height Back with Adjustable Handle (Short 18" - 25") N/C	Gas Spring
Not available with O2 Holder or Dynamic Back	Standard Spring N/C
Adjustable Height Back with Adjustable Handle (Tall 20" - 27") N/C	Heavy Duty Spring N/C
20"-22" Heights not available with Dynamic Back	
Reclining Back with Adjustable Handle (Short 22") (Non-Folding) \$1025	Stroller Handle Extensions
Not available with O2 Holder or Dynamic Back. If no Recline Angle and Gas Spring selections are made, 0° - 65° and Standard Spring will be used.	Available only with Stroller Handle Back Post
Reclining Back with Adjustable Handle (Tall 26") (Non-Folding) \$1025	Removable Stroller Handle Extension \$97
Not available with O2 Holder or Dynamic Back. If no Recline Angle and Gas Spring selections are made, 0° - 65° and Standard Spring will be used.	
Backrest Height	Backrest Options
Required when Adjustable Height Back with Adjustable Handle Back Post is selected	Non Folding Back N/C
18" <sup>1</sup> 19" <sup>1</sup> 20" 21" 22" 23" 24" 25" 26" <sup>2</sup>	Standard with Reclining Back and Dynamic Back
27" <sup>2</sup> N/C	Dynamic Back (E2398) \$2045
	Available only with Non Folding Back. Not available with Angle Adjustable Locking Extendable Flip Up Armrest or 400 lb Weight Capacity.
	Adjustable Height Rigidizer Bar \$82
	Not available with 23" or 24" frame widths

1 Available only with Short  
2 Available only with Tall

## Footrest

### Hanger / Front Frame Bend

Hanger	60°	70°	80°	90°		Retail
Extension Tube Mount			1			N/C
Front Mount			1	2		N/C
Pro Elevating Leg Rests E0990						\$318
Pediatric Pro Elevating Leg Rests E0990						\$318
Contracture Footrest Bilateral Mount <sup>3</sup>						\$382
Contracture Footrest Center Mount <sup>3</sup>						\$382
Omit Hangers						(\$20)

1 7" or 8" Casters may interfere with less than 5° tilt

2 Not available with 7" or 8" Casters

3 Aluminum Angle Adjustable footplates included. See image to right.

### PRO ELR Cover

Available only on PRO ELR


Gel Pad Covers (PAIR) (K0108) \$425

### 3" Extension

Available only with Front Mount Hangers

3" Extension N/C

3" Extension (Provided as Parts) N/C



### Footrest / Extension / Seat to Footrest Length

Select Seat to Footrest Length in row corresponding to desired Footrest and Extension. Footrest Length of at least 2.5" shorter than Front Seat Height recommended. Extension Tube Hanger measurements based on 70° Hanger Bend.

Measurements based on 70° Hanger Bent.														
Footrest	Extension	Seat to Footrest Length												
Extension Tube Hangers (60°,70°,80°)	Short (10" - 13.5")	10"	10.5"	11"	11.5"	12"	12.5"	13"	13.5"					
	Medium (11.5" - 15.5")	11.5"	12"	12.5"	13"	13.5"	14"	14.5"	15"	15.5"				
	Long (14" - 18")	14"	14.5"	15"	15.5"	16"	16.5"	17"	17.5"	18"				
Front Mount Hangers <sup>1</sup>			Standard Selections									Requires 3" Extension for \$35		
	Short Front Mount Hanger	4" <sup>2</sup>	4.5" <sup>2</sup>	5"	5.5"	6"						8.5"	9"	9.5"
		6.5"	7"	7.5"	8"						10"	10.5"	11"	
	Medium Front Mount Hanger	4" <sup>2</sup>	4.5" <sup>2</sup>	5"	5.5"	6"	6.5"	7"				10.5"	11"	11.5"
		7.5"	8"	8.5"	9"	9.5"	10"				12"	12.5"	13"	
	Long Front Mount Hanger	5"	5.5"	6"	6.5"	7"	7.5"	8"	8.5"			12.5"	13"	13.5"
		9"	9.5"	10"	10.5"	11"	11.5"	12"				14"	14.5"	15"
Contracture Footrest	Short (10" - 13.5")	10.5"	11"	11.5"	12"	12.5"	13"	13.5"	14"	14.5"				
	Medium (11.5" - 15.5")	13.5"	14"	14.5"	15"	15.5"	16"	16.5"	17"	17.5"				
	Long (14" - 18")	16.5"	17"	17.5"	18"	18.5"	19"	19.5"	20"	20.5"				
PRO ELR Low Pivot Setting	PRO ELR Low Pivot Setting Short (10.5" - 15.5")	10.5" - 15.5"												
	PRO ELR Low Pivot Setting Long (13" - 19")	13" - 19"												
PRO ELR High Pivot Setting	PRO ELR High Pivot Setting Short (8.5" - 17")	8.5" - 17"												
	PRO ELR High Pivot Setting Long (11.5" - 17.5")	11.5" - 17.5"												

1 See Minimum Front Seat Height chart below

2 Available with 90° only

Minimum Front Seat Height with Front Mount Hangers					
Hanger Length	Short		Medium		Long
Hanger Bend	70°	80/90°	70°	80/90°	70° 80°
Minimum Front Seat Height with Standard Selection	13"	13"	13"	13.5"	14.5"
Minimum Front Seat Height with 3" Extension	13.5"	14"	15.5"	16"	17"



Phone (800) 981-1540 | Fax (715) 254-0996  
sales@kimobility.com | www.kimobility.com

# FOCUS CRe

## Footplates

*No selection needed with Contracture Footrest*

Composite	N/C
<i>Not available with Front Mount Hanger</i>	
Composite Angle Adjustable (K0040)	\$195
<i>Not available with Front Mount Hanger or 14" Seat Width without Offset Hanger</i>	
Aluminum Angle Adjustable (K0040)	\$210
Aluminum Locking Multi-Angle Adjustable (K0108)	\$318
<i>Available only with Front Mount Hanger</i>	
One Piece Flip Up Angle Adjustable (K0040)	\$205
<i>Not available with Front Mount Hanger or PRO ELR</i>	

## Heel Loops

*Not available with One Piece Flip Up Angle Adjustable Footplate*

Heel Loops (E0951)	\$46
Adjustable Heel Loops (E0951)	\$56
Omit Heel Loops	N/C

## Calf Straps

Velcro Adjustable (K0038)	\$41
Padded Velcro Adjustable (K0038)	\$72
Bodypoint Padded AeroMesh (K0038)	\$128

## Hanger Release

*No selection needed with Contracture Footrest. If no selection is made 4-Way - In line Hanger Release will be used.*

Classic	N/C
4-Way - In-Line	N/C
<i>Not available with PRO ELR or Peds PRO ELR</i>	
4-Way - Out	N/C
<i>Not available with PRO ELR or Peds PRO ELR</i>	

## Residual Limb Support

	Limb Support	Price	Omit	Cover		Pad Size					
Left	Left Limb Support (E1020)	\$431	Omit Left Hanger	Smooth Cover Left	Air Mesh Cover Left	5X8 Pad Left Side	7X10 Pad Left Side	9x10 Pad Left Side	9x14 Pad Left Side	11x10 Pad Left Side	11x14 Pad Left Side
Right	Right Limb Support (E1020)	\$431	Omit Right Hanger	Smooth Cover Right	Air Mesh Cover Right	5X8 Pad Right Side	7X10 Pad Right Side	9x10 Pad Right Side	9x14 Pad Right Side	11x10 Pad Right Side	11x14 Pad Right Side

## Seat Height, Casters and Rear Wheels

Select Seat Height corresponding to desired Rear Wheel and Caster Size

Rear Wheel Size	Caster Size					
	5"	6"	6" x 2"	7"	8"	8" x 2"
12" with Pneumatic or Pneumatic with Airless Insert Tire <sup>3</sup>	Seat Height					
	13.5"					
	14"					
	14.5"	14.5"	14.5"		17"	
	17.5"	17.5"	17.5"	17.5"	17.5"	17.5"
	18"	18"	18"	18"	18"	18"
	18.5"	18.5"	18.5"	18.5"	18.5"	18.5"
	19"	19"	19"	19"	19"	19"
12" with Full Poly Tire <sup>3</sup>	16.5"	16.5"	16.5"			
	17"	17"	17"	17"	17"	17"
	17.5"	17.5"	17.5"	17.5"	17.5"	17.5"
	18"	18"	18"	18"	18"	18"
	18.5"	18.5"	18.5"	18.5"	18.5"	18.5"
	19"	19"	19"	19"		
12" with Drum Brake	13"					
	13.5"					
	14"	14"	14"			
	14.5"	14.5"	14.5"			
	17.5"	17.5"	17.5"	17.5"	17.5"	17.5"
	18"	18"	18"	18"	18"	18"
	18.5"	18.5"	18.5"	18.5"	18.5"	18.5"
	19"	19"	19"	19"	19"	19"
	19.5"	19.5"	19.5"	19.5"	19.5"	19.5"
	20" <sup>2</sup>		20" <sup>2</sup>	20" <sup>2</sup>	20" <sup>2</sup>	
16" with Pneumatic, Pneumatic with Airless Insert or Full Poly Tire <sup>4</sup>	13"					
	13.5"					
	14"	14"	14"			
	14.5"	14.5"	14.5"			
	15"	15"	15"	15"		
	15.5"	15.5"	15.5"	15.5"		
	16"	16"	16"	16"	16"	16"
	16.5" <sup>1</sup>	16.5" <sup>1</sup>	16.5" <sup>1</sup>	16.5" <sup>1</sup>		16.5" <sup>1</sup>
	19"	19"	19"	19" <sup>2</sup>	19"	19"
	19.5"	19.5"	19.5"	19.5"	19.5"	19.5"
	20"	20"	20"	20"	20"	20"
	20.5"	20.5"	20.5"	20.5"	20.5"	20.5"
16" with Low Poly Tire <sup>4</sup>	13"					
	13.5"					
	14"	14"	14"			
	14.5"	14.5"	14.5"			
	15"	15"	15"	15"		
	18"	18"	18"	18"	18"	18"
	18.5"	18.5"	18.5"	18.5"	18.5"	18.5"
	19"	19"	19"	19"	19"	19"
	19.5"	19.5"	19.5"	19.5"	19.5"	19.5"
	20"	20"	20"	20"	20"	20"
			20.5"	20.5"	20.5"	

Rear Wheel Size	Caster Size					
	5"	6"	6" x 2"	7"	8"	8" x 2"
16" with Drum Brake <sup>4</sup>	Seat Height					
	13.5"					
	14"	14"	14"			
	14.5"	14.5"	14.5"			
	15"	15"	15"	15"		
	15.5"	15.5"	15.5"	15.5"		
	16"	16"	16"			16"
		18.5"	18.5"	18.5"	18.5"	18.5"
	19"	19"	19"	19"	19"	19"
	19.5"	19.5"	19.5"	19.5"	19.5"	19.5"
	20"	20"	20"	20"	20"	20"
	20.5"	20.5"	20.5"	20.5"	20.5"	20.5"
20"	15"	15"	15"	15"		
	15.5"	15.5"	15.5"	15.5"		
	16"	16"	16"	16"	16"	16"
	16.5"	16.5"	16.5"	16.5"	16.5"	16.5"
	17"	17"	17"	17"	17"	17"
	17.5"	17.5"	17.5"	17.5"	17.5"	17.5"
	18"	18"	18"	18"	18"	18"
	18.5" <sup>2</sup>	18.5" <sup>2</sup>	18.5" <sup>2</sup>	18.5" <sup>2</sup>		
22"	16"	16"	16"	16"	16"	
	16.5"	16.5"	16.5"	16.5"	16.5"	16.5"
	17"	17"	17"	17"	17"	17"
	17.5"	17.5"	17.5"	17.5"	17.5"	17.5"
	18"	18"	18"	18"	18"	18"
	18.5"	18.5"	18.5"	18.5"	18.5"	18.5"
	19"	19"	19"	19"	19"	19"
24"	16.5" <sup>5</sup>	16.5" <sup>1,5</sup>	16.5" <sup>1,5</sup>			
	17" <sup>5</sup>	17" <sup>5</sup>	17" <sup>5</sup>	17" <sup>5</sup>	17" <sup>5</sup>	17" <sup>5</sup>
	17.5" <sup>6</sup>	17.5" <sup>6</sup>	17.5" <sup>6</sup>	17.5"	17.5"	17.5" <sup>6</sup>
	18"	18"	18"	18"	18"	18"
	18.5"	18.5"	18.5"	18.5"	18.5"	18.5"
	19"	19"	19"	19"	19"	19"
	19.5"	19.5"	19.5"	19.5"	19.5"	19.5"
	20"	20"	20"	20"	20"	20" <sup>2</sup>

<sup>1</sup> Available only with Full Poly tire

<sup>2</sup> Not available with Full Poly tire

<sup>3</sup> Not available with Vent and Battery Tray combination

<sup>4</sup> Seat heights higher than 16.5" with 16" wheel size are not available with Vent and Battery Tray combination

<sup>5</sup> Not available with Drum Brake

<sup>6</sup> Not available with Drum Brake in combination with Pneumatic tires



Phone (800) 981-1540 | Fax (715) 254-0996  
sales@kimobility.com | www.kimobility.com

# FOCUS CRe

## Casters and Rear Wheels

### One Arm Drive

One Arm Drive - Quick Release (E0958) \$1379

*Available only with 24" Spoke wheels. Spoke wheel style is exclusive to OAD.  
Not available with Vent Tray, Battery Box Holder, Self Propel Axle Plate or  
Angle Adjustable Locking Extendable Flip Up. Standard Handrim will be  
Aluminum Anodized. Axle is Quick Release*

Left Side Drive N/C  
Right Side Drive N/C  
Aluminum Drive Handrim N/C  
Poly Coated Drive Handrim \$205

### Casters

		Caster Size				
		Price	5"	6"	7"	8"
Caster Type	1" Poly	N/C				
	1" Poly Aluminum	\$118				
	1" Pneumatic	\$118				
	1.5" Poly	\$95				
	1.5" SoFE Roll Aluminum	\$190				
	2" Poly	\$97				
	2" Pneumatic w/Foam Insert	\$246				

### Axle

*No selection needed with One Arm Drive or Drum Brake*

Fixed N/C  
Quick Release N/C

### Self Propel Axle Plate

Self Propel Axle Plate \$103  
*Not available with Vent Tray or 12" or 16" wheel sizes*

## Wheel Locks

### Wheel Lock

Push to Lock N/C  
Pull to Lock N/C  
*Not available with 12" or 16" Full Poly tires*  
Attendant Foot Lock \$87  
Combination Attendant Foot Lock with Push to Lock with Extension Handles \$159  
*Available only with 12" or 16" wheel sizes. Includes Extension Handles. Push to Lock only.*  
Drum Brake (E2228) \$595  
*Available only with 12" Mag, 16" Mag, or MAXX Spoke wheels. Not available with 16" Low Poly tire, Reclining Backrest, One Arm Drive or Self Propel Axle Plate. Axle is Quick Release.*  
Dual - Drum Brake with Push to Lock (E2228) \$702  
*Available only with 12" Mag, 16" Mag, or MAXX Spoke wheels. Not available with 16" Low Poly tire, Reclining Backrest, One Arm Drive or Self Propel Axle Plate. Axle is Quick Release.*  
Dual - Drum Brake with Pull to Lock (E2228) \$702  
*Available only with 12" Mag, 16" Mag, or MAXX Spoke wheels. Not available with 16" Low Poly tire, Reclining Backrest, One Arm Drive or Self Propel Axle Plate. Axle is Quick Release.*  
Omit Wheel Locks N/C

### Rear Wheel

*No selection needed with One Arm Drive*

		Rear Wheel Size					
		Price	12"	16"	20"	22"	24"
Rear Wheel Type	Mag	N/C					
	MAXX Spoke	\$200					

### Rear Tire Type

Full Poly (E2218) N/C  
Pneumatic (E2211,E2212) N/C  
Pneumatic with Airless Insert (E2211, E2213) \$97  
Low Poly (E2218) N/C  
*Available only with 16" wheel size*

### Handrims

*Handrims not available on 12" or 16" wheel sizes. No selection needed with One Arm Drive*

#### Standard Handrims

Aluminum Anodized N/C  
Plastic Coated \$123  
Projection (8 Vertical Projections) \$349  
*Available only with full profile tires in combination with 22", 24", and 26" MAXX Spoke wheels, or 22" and 24" Mag wheels*  
Omit Handrim N/C

### Handrim Option

*Not available on Plastic Coated or Projection Handrims*

3M Micro Grip \$67

## Armrests

Armrest Type	Length	Price
Height Adjustable T-Arm Tall (11.5"-15") E0973	Desk Full	\$323
Height Adjustable T-Arm (8.5"-13.5") E0973	Desk Full	\$323
Height Adjustable T-Arm Low (6.5"-10.5") E0973 <sup>1</sup>	Desk Full	\$323

<sup>1</sup> If requested configuration is unavailable, Standard Height T-Arm will be used. Not available with 14" depth and offset hangers.

<sup>2</sup> Not available with IV Pole

<sup>3</sup> Not available with Dynamic Back

<sup>4</sup> Minimum armrest height with Reclining Backrest is 12"

Armrest Type	Length	Price
Height Adj. Dual Post Flip Back (9.5"-14.5") E0973 <sup>2</sup>	Desk Full	\$320
Angle Adjustable Locking Extendable Flip Up E0973 <sup>3,4</sup>	Desk Full	\$297
Omit Armrests		N/C

## Arm Pads

Standard Arm Pad	N/C
Waterfall Arm Pad	\$26
Not available with Angle Adjustable Locking Extendable Flip Up Armrests	
Foam Grip	N/C
Available only with Angle Adjustable Locking Extendable Flip Up Armrest	

## Accessories

### Removable Side Guard

Composite Side Guards (K0108)	\$149
Not available with T-Arm or Dual Post Flip Back Armrests	

### Anti-tips

Rear Anti-Tips (E0971)	\$168
Required with Reclining Backrest or Vent Tray	
Omit	N/C

### Pelvic Positioning Belt

1 1/2" Auto Buckle Pelvic Positioning Belt (E0978)	\$87
1 1/2" Padded Auto Buckle Pelvic Positioning Belt (E0978)	\$179
2" Airline Buckle Pelvic Positioning Belt (E0978)	\$118
Small Bodypoint Evoflex (16-24") (E0978)	\$200
Medium Bodypoint Evoflex (19-29") (E0978)	\$200
Large Bodypoint Evoflex (25-37") (E0978)	\$200

### Belt Mounting Kit

Ki Mobility Belt Mounting Brackets	\$51
Belt Mounting Band Clamp	\$36

### U.S. Military Service Patches


Army	Navy	Marines	Air Force	
Coast Guard	US Purple Heart			\$92
Combo Army/Purple Heart		Combo Navy/Purple Heart		
Combo Marines/Purple Heart		Combo Air Force/Purple Heart		
Combo Coast Guard/Purple Heart				\$149

### O2 Holder

O2 Holder (E2208)	\$179
Not available with Reclining Backrest or Dynamic Back	

### Other Accessories

Tool Kit	\$62
Removable Underseat Pouch	\$56
Ki Mobility Backpack	\$72
Black Neoprene Impact Guards	\$97
Tilt Stop (Pair) (K0108)	\$31
Included with Offset Hanger	
IV Pole	\$292
Canopy with Hardware (K0108)	\$297
Not available on frame widths greater than 18"	





Phone (800) 981-1540 | Fax (715) 254-0996  
sales@kimobility.com | www.kimobility.com



Frame Color

Gloss Finish		Matte Finish	
Black	N/C	Military Green	N/C
Blue Skies	N/C	Charcoal Gray	N/C
Bubblegum Pink	N/C	Emerald Green	N/C
Candy Blue	N/C	Iridescent Burgundy	N/C
Candy Purple	N/C	Iridescent Navy	N/C
Candy Red	N/C	Retro Red	N/C
Copper	N/C	Royal Blue	N/C
Granny Smith	N/C	Shadow Black	N/C
Metallic Orange	N/C	Smooth White	N/C
Yellow	N/C		