



Phone (800) 981-1540 | Fax (715) 254-0996
sales@kimobility.com | www.kimobility.com

LIBERTY

Company:

Ship To:

Account Number:

Date:

Address:

Quote: PO Number:

City:

State:

ATP:

Zip:

Phone:

Marked For:

Email/Fax:

Associated Quote/Order#:

Effective: 1/10/2024

Liberty - Transit Included (E1161)	\$3629
Liberty - Transit Not Included (E1161)	\$3383
Liberty - Rigid, Transit Included (E1161)	\$3475
Liberty - Rigid, Transit Not Included (E1161)	\$3229
250 lb weight capacity Tilt range: 0° - 20°	

Frame

Frame Width

Select Seat Width within desired Seat Width Range

Frame Width Range	Frame Width				
Small (14" - 18")	14" N/C	15" N/C	16" N/C	17" N/C	18" N/C
Medium (16" - 20")	16" N/C	17" N/C	18" N/C	19" N/C	20" (E2201) \$520
Large (18" - 22")	18" N/C	19" N/C	20" (E2201) \$520	21" (E2201) \$510	22" (E2201) \$502

Seat Depth

Select Seat Depth within desired Seat Depth Range

Seat Depth Range	Seat Depth		
Short (14" - 16")	14" N/C	15" N/C	16" N/C
Medium (16" - 18")	16" N/C	17" N/C	18" N/C
Long (18" - 20")	18" N/C	19" N/C	20" (E2203) \$630

Assistive Gas Strut Set-Up

Client Weight > 200 lbs (Setting 1, 450N Spring)	N/C
Client Weight 176-200 lbs (Setting 2, 450N Spring)	N/C
Client Weight 151-175 lbs (Setting 3, 450N Spring)	N/C
Standard on chairs less than 16" deep	
Client Weight 151-175 lbs (Setting 1, 300N Spring)	N/C
Standard on chairs 16" deep or larger	
Client Weight 125-150 lbs (Setting 2, 300N Spring)	N/C
Client Weight < 125 lbs (Setting 3, 300N Spring)	N/C

Seat Option

Standard Seat Sling	N/C
Growing Adjustable Seat Pan (E2231)	\$231
Omit Seat	N/C

Transit Option

Transit System	\$250
No charge when selected with Liberty - Transit Included	
Omit	N/C



Phone (800) 981-1540 | Fax (715) 254-0996
sales@kimobility.com | www.kimobility.com

LIBERTY

Backrest

Back Post / Height Range / Back Height

Back Post	Price	Height Range	Back Height				
Stroller Handle Back Post ¹	N/C	Short	20"				
		Tall	24"				
Height Adjustable Straight with Push Handle	\$26	Short	14"	15"	16"	17"	18"
		Medium	16"	17"	18"	19"	20"
		Tall	18"	19"	20"	21"	22"
Height Adjustable 8° Bend with Push Handle ²	\$26	Short	14"	15"	16"	17"	18"
		Medium	16"	17"	18"	19"	20"
		Tall	18"	19"	20"	21"	22"

¹ Not available with Back Upholstery

² Not available with IV Pole

Stroller Handle Back Post Option

Available only with Stroller Handle Back Post

Removable Stroller Handle Extension _____ \$97

Backrest Options

Fixed (Non-Folding) _____ N/C

Angle Adjustable (Folding) _____ \$103

Back Upholstery

If selecting both upholstery and Axiom Back, the upholstery will be installed.

Standard _____ N/C

Tension Adjustable (E2611) _____ \$297

Not available on Seat Widths less than 16" combined with 21" and 22" Backrest Heights

Omit _____ (\$15)

U.S. Military Service Patches

Attached only if ordered with Tension Adjustable Back Upholstery

Army Navy Marines Air Force
Coast Guard US Purple Heart _____ \$92

Combo Army/Purple Heart Combo Navy/Purple Heart
Combo Marines/Purple Heart Combo Air Force/Purple Heart
Combo Coast Guard/Purple Heart _____ \$149

Footrest

Hanger

Hanger	60°	70°	80°	90°		Retail
Extension Tube Mount			1			N/C
Front Mount			1	2		N/C
Pro Elevating Leg Rests E0990						\$318
Omit Hangers						(\$20)

¹ Not available with 8" casters

² Not available with 8" casters in combination with 3" extension

PRO ELR Cover

Available only on PRO ELR

Gel Pad Covers (PAIR) (K0108) _____ \$425

Footplates

Composite _____ N/C

Not available with Front Mount Hanger

Composite Angle Adjustable (K0040) _____ \$195

Not available with Front Mount Hanger or 14" Seat Width

Aluminum Angle Adjustable (K0040) _____ \$210

Aluminum Locking Multi-Angle Adjustable (K0108) _____ \$318

Available only with Front Mount Hanger

One Piece Flip Up Angle Adjustable (K0040) _____ \$205

Not available with Front Mount Hanger or PRO ELR

Footrest / Extension / Seat to Footrest Length

Select Extension and Seat to Footrest Length in row corresponding to desired Footrest. Footrest Length of at least 2.5" shorter than Seat Height recommended. Extension Tube Hanger measurements based on 70° Hanger Bend.

Measurements Based on 70° Hanger Bend.														
Footrest	Extension	Seat to Footrest Length												
Extension Tube Hangers (60°,70°,80°)	Short (10" - 13.5")	10"	10.5"	11"	11.5"	12"	12.5"	13"	13.5"					
	Medium (11.5" - 15.5")	11.5"	12"	12.5"	13"	13.5"	14"	14.5"	15"	15.5"				
	Long (14" - 18")	14"	14.5"	15"	15.5"	16"	16.5"	17"	17.5"	18"				
Front Mount Hangers ¹		Standard Selections									Requires 3" Extension for \$35			
	Short Front Mount Hanger	4" ²	4.5" ²	5"	5.5"	6"					8.5"	9"	9.5"	
		6.5"	7"	7.5"	8"					10"	10.5"	11"		
	Medium Front Mount Hanger	4" ²	4.5" ²	5"	5.5"	6"	6.5"	7"			10.5"	11"	11.5"	
		7.5"	8"	8.5"	9"	9.5"	10"			12"	12.5"	13"		
	Long Front Mount Hanger	5"	5.5"	6"	6.5"	7"	7.5"	8"	8.5"			12.5"	13"	13.5"
		9"	9.5"	10"	10.5"	11"	11.5"	12"			14"	14.5"	15"	
PRO ELR Low Pivot Setting	PRO ELR Low Pivot Setting Short (10.5" - 15.5")	10.5" - 15.5"												
	PRO ELR Low Pivot Setting Long (13" - 19")	13" - 19"												
PRO ELR High Pivot Setting	PRO ELR High Pivot Setting Short (8.5" - 17")	8.5" - 17"												
	PRO ELR High Pivot Setting Long (11.5" - 17.5")	11.5" - 17.5"												

¹ See Minimum Seat Height chart below

² Available with 90° only

Minimum Seat Height with Front Mount Hangers						
Hanger Length	Short		Medium		Long	
Hanger Bend	70°	80/90°	70°	80/90°	70°	80°
Minimum Seat Height with Standard Selection	13"	13"	13"	13.5"	14.5"	15.5"
Minimum Seat Height with 3" Extension	13.5"	14"	15.5"	16"	17"	18"

Hanger Extension Tubes

Available only with Front Mount Hangers

3" Extension	N/C
3" Extension (Provided as Parts)	N/C

Heel Loops

Not available with One Piece Flip Up Angle Adjustable Footplate

Heel Loops (E0951)	\$45
Omit Heel Loops	N/C

Calf Straps

Velcro Adjustable (K0038)	\$45
Padded Velcro Adjustable (K0038)	\$72
Bodypoint Padded AeroMesh (K0038)	\$128

Residual Limb Support

	Limb Support	Price	Omit	Cover		Pad Size					
Left	Left Limb Support (E1020)	\$431	Omit Left Hanger	Smooth Cover Left	Air Mesh Cover Left	5X8 Pad Left Side	7X10 Pad Left Side	9x10 Pad Left Side	9x14 Pad Left Side	11x10 Pad Left Side	11x14 Pad Left Side
Right	Right Limb Support (E1020)	\$431	Omit Right Hanger	Smooth Cover Right	Air Mesh Cover Right	5X8 Pad Right Side	7X10 Pad Right Side	9x10 Pad Right Side	9x14 Pad Right Side	11x10 Pad Right Side	11x14 Pad Right Side

Hanger Release

4-Way - In-Line	N/C
Not available with PRO ELR	
4-Way - Out	N/C



Phone (800) 981-1540 | Fax (715) 254-0996
sales@kimobility.com | www.kimobility.com

LIBERTY

Seat Height, Casters and Rear Wheels

Select Seat Height corresponding to desired Rear Wheel and Caster Size

Rear Wheel Size	Caster Size			
	4"	5"	6"	8"
12"	Seat Height			
	13"	14"	15"	
	13.5"	14.5"	15.5"	17"
	14"	15"	16"	17.5"
	14.5"	15.5"	16.5"	18"
	15"	16"	17"	18.5"
	15.5"	16.5"	17.5"	19" ²
	16"	17"	18"	
		17.5"	18.5"	
		18"	19" ²	
16"	13.5" ¹	14"	15"	
	14" ¹	14.5"	15.5"	17"
	14.5"	15"	16"	17.5"
	15"	15.5"	16.5"	18"
	15.5"	16"	17"	18.5"
	16"	16.5"	17.5"	19"
		17"	18"	19.5"
		17.5"	18.5"	20"
		18"	19"	
		18.5"		
20"	13"	14"	15"	
	13.5"	14.5"	15.5"	17"
	14"	15"	16"	17.5"
	14.5"	15.5"	16.5"	18"
	15"	16"	17"	18.5"
	15.5"	16.5"	17.5"	19"
	16"	17"	18"	19.5"
		17.5"	18.5"	20"
		18"	19"	
		18.5"		

¹ Available only with Low Profile tires

² Not available with Full Poly tire

³ Not available with Low Profile tires in combination with Mag wheel

Rear Wheel Size	Caster Size			
	4"	5"	6"	8"
22" ³	Seat Height			
	13.5" ¹	14"	15"	
	14"	14.5"	15.5"	17"
	14.5"	15"	16"	17.5"
	15"	15.5"	16.5"	18"
	15.5"	16"	17"	18.5"
	16"	16.5"	17.5"	19"
		17"	18"	19.5"
		17.5"	18.5"	20"
		18"	19" ²	
24" ³	14.5" ¹	14.5" ¹	15"	
	15"	15"	15.5"	17"
	15.5"	15.5"	16"	17.5"
	16"	16"	16.5"	18"
		16.5"	17"	18.5"
		17"	17.5"	19"
		17.5"	18"	19.5"
		18"	18.5"	20"
		18.5"	19"	
		19"		



Phone (800) 981-1540 | Fax (715) 254-0996
sales@kimobility.com | www.kimobility.com

LIBERTY

Casters

		Caster Size				
		Price	4"	5"	6"	8"
Caster Type	1" Poly	N/C				
	1" Poly Aluminum	\$118				
	1" Pneumatic	\$118				
	1.5" Poly	\$95				
	1.5" Soft Roll Aluminum	\$190				
	2" Poly	\$100				
	2" Pneumatic w/Foam Insert	\$256				

Rear Wheel

		Rear Wheel Size					
		Price	12"	16"	20"	22"	24"
Rear Wheel Type	Mag	N/C					
	Spoke	N/C					

Rear Tire Type

Full Profile

Street Pneumatic (E2211, E2212)	N/C
Full Poly 1 3/8" (E2218)	N/C
Pneumatic with Airless Insert (E2211, E2213)	\$97

Low Profile

Low Poly 1" (E2218)	N/C
<i>Available only with 16" wheel size</i>	
Sentinel Solid Tire 1" (E2218)	\$103
<i>Not available on 12" or 16" wheel sizes or with 20" and 22" Mag wheels</i>	
Schwalbe Marathon Plus Evolution 1" (E2211, E2212)	\$133
<i>Available only on 24" Spoke</i>	

Wheel Locks

Wheel Lock

Push to Lock	N/C
Pull to Lock	N/C
<i>Not available with 12" or 16" Full Poly tires</i>	
Hemi Wheel Lock - Left (K0108)	\$318
<i>Extension Handle included</i>	
Hemi Wheel Lock - Right (K0108)	\$318
<i>Extension Handle included</i>	
Drum Brake (E2228)	\$595
<i>Available only with Spoke wheels. Axle is Quick Release.</i>	
Dual - Drum Brake with Push to Lock (E2228)	\$702
<i>Available only with Spoke wheels. Axle is Quick Release.</i>	
Dual - Drum Brake with Pull to Lock (E2228)	\$702
<i>Available only with Spoke wheels. Axle is Quick Release.</i>	
Omit Wheel Locks	N/C

Handrims

Standard Handrims

Handrims not available on 12" or 16" wheel sizes

Aluminum Anodized	N/C
Plastic Coated	\$123
Projection (8 Vertical Projections)	\$349
<i>Available only with full profile tires in combination with 22", 24", and 26" Spoke wheels, or 22" and 24" Mag wheels</i>	
Omit Handrim	N/C

Standard Handrim Mount

If no selection is made, Standard will be used

Close Mount	N/C
<i>Reduces overall width by .75"</i>	

Handrim Options

Not available on Plastic Coated or Projection Handrims

3M Micro Grip	\$72
---------------	------

Axles

No selection needed with Drum Brake

Fixed	N/C
Quick Release	N/C

Wheel Lock Options

Extension Handles (E0961)	\$82
<i>Available only with Push to Lock or Pull to Lock Wheel Locks</i>	
Grade Aids (E0974)	\$159
<i>Available only with Push to Lock Wheel Locks and Pneumatic or Pneumatic with Airless Insert tires</i>	



Phone (800) 981-1540 | Fax (715) 254-0996
sales@kimobility.com | www.kimobility.com

LIBERTY

Armrests

Armrest Type	Length	Price
Height Adjustable T-Arm Tall (11.5"-15") E0973	Desk	\$323
	Full	
Height Adjustable T-Arm (8.5"-13.5") E0973	Desk	\$323
	Full	
Angle Adjustable Locking Extendable Flip Up E0973 ¹	Desk	\$297
	Full	
Omit Armrests		N/C

¹ Available only with 12" and 16" wheel sizes and Stroller Handle Back Post

Accessories

Removable Side Guard

Composite Side Guards (K0108)	\$149
<i>Not available with T-Arms</i>	

Pelvic Positioning Belt

1 1/2" Auto Buckle Pelvic Positioning Belt (E0978)	\$87
1 1/2" Padded Auto Buckle Pelvic Positioning Belt (E0978)	\$179
2" Airline Buckle Pelvic Positioning Belt (E0978)	\$118
Small Bodypoint Evoflex (16-24") (E0978)	\$200
Medium Bodypoint Evoflex (19-29") (E0978)	\$200
Large Bodypoint Evoflex (25-37") (E0978)	\$200

Belt Mounting Kit

Seat Belt Mounting Hardware	\$16
Belt Mounting Band Clamp	\$36

O2 Holder

O2 Holder (E2208)	\$313
<i>Mounted to the left side. Not available with Fold Down Back.</i>	

Arm Pads

Standard Arm Pad	N/C
<i>Not available with Angle Adjustable Locking Extendable Flip Up Armrests in combination with 20"-24" wheel sizes</i>	
Waterfall Arm Pad	\$26
<i>Not available with Angle Adjustable Locking Extendable Flip Up Armrests</i>	
Foam Grip	N/C
<i>Available only with Angle Adjustable Locking Extendable Flip Up Armrest</i>	

Anti-Tip

Rear Anti-Tips (E0971)	\$165
Omit	N/C

Transit Wheels

Transit Wheels (K0108)	\$318
<i>Not available with Long Depth Seat Frame, Seat Heights less than 14" or greater than 18", 12" or 16" rear wheels, Cane and Crutch Holder, and O2 Holder</i>	

Other Accessories

Tool Kit	\$62
Removable Underseat Pouch	\$56
Ki Mobility Backpack	\$72
Black Neoprene Impact Guards	\$97
Cane and Crutch Holder	\$210
IV Pole - Mounted Right Side	\$292
<i>Not available with Transit Option or 8° Bend with Push Handle Back Post. Left side mount not available with O2 Holder.</i>	
IV Pole - Mounted Left Side	\$292
<i>Not available with Transit Option or 8° Bend with Push Handle Back Post. Left side mount not available with O2 Holder.</i>	

Frame Color

Standard Colors

Black	N/C
Military Green	N/C
Candy Blue	N/C
Charcoal Gray	N/C
Candy Purple	N/C
Iridescent Burgundy	N/C
Candy Red	N/C
Iridescent Navy	N/C
Copper	N/C
Shadow Black	N/C

Optional Colors

Extended Lead Time

Blue Skies	\$72
Emerald Green	\$72
Bubblegum Pink	\$72
Retro Red	\$72
Royal Blue	\$72
Granny Smith	\$72
Smooth White	\$72
Metallic Orange	\$72
Teal	\$72
Yellow	\$72