



Phone (800) 981-1540 | Fax (715) 254-0996
sales@kimobility.com | www.kimobility.com

LIBERTY Basic

Company:

Ship To:

Account Number:

Date:

Address:

Quote: PO Number:

City:

State:

ATP:

Zip:

Phone:

Marked For:

Email/Fax:

Associated Quote/Order#:

Effective: 1/10/2024

Liberty Basic (E1161) \$3383

Liberty Basic - Rigid (E1161) \$3229

250 lb weight capacity Tilt range: 0° - 20° Standard seat upholstery

Frame

Frame Width

14"-16" Seat Widths ship with 300 N Assistive Gas Strut, 18"-22" Seat Widths ship with 450 N Assistive Gas Strut

14"	N/C
16"	N/C
18"	N/C
20" (E2201)	\$510
22" (E2201)	\$510

Seat Depth

14"	N/C
16"	N/C
18"	N/C
20" (E2203)	\$630

Seat Option

Standard Seat Sling	N/C
Growing Adjustable Seat Pan (E2231)	\$231
Omit Seat	N/C

Transit Option

Transit System	\$250
Omit	N/C

Backrest

Back Post / Height Range / Back Height

Back Post	Price	Height Range	Back Height				
Height Adjustable 8° Bend with Push Handle	\$26	Short	14"	15"	16"	17"	18"
		Medium	16"	17"	18"	19"	20"
		Tall	18"	19"	20"	21"	22"

Backrest Options

Fixed (Non-Folding)	N/C
Angle Adjustable (Folding)	\$103

Back Upholstery

If selecting both upholstery and Axiom Back, the upholstery will be installed

Standard	N/C
Tension Adjustable (E2611)	\$265
Not available on Seat Widths less than 16" combined with 21" and 22" Backrest Heights	
Omit	(\$15)

Footrest

Hanger / Front Frame Bend

Hanger	60°	70°		Retail
Extension Tube Mount				N/C
Pro Elevating Leg Rests E0990				\$318
Omit Hangers				N/C

Footplates

Composite	N/C
Aluminum Angle Adjustable (K0040)	\$210

HCPCS codes provided are based on code verification by the PDAC or our interpretation of the code definitions.

Proper use of HCPCS codes does not ensure coverage or payment. For coverage information, verify the appropriate payer's coverage policy.

DCN0524.10 Pricing and specifications subject to change without notice.



Phone (800) 981-1540 | Fax (715) 254-0996
sales@kimobility.com | www.kimobility.com

LIBERTY Basic

Footrest / Extension / Seat to Footrest Length

Select Extension and Seat to Footrest Length in row corresponding to desired Footrest. Footrest Length of at least 2.5" shorter than Seat Height recommended. Extension Tube Hanger measurements based on 70° Hanger Bend.

Footrest	Extension	Seat to Footrest Length								
Extension Tube Hangers (60°, 70°)	Short (10" - 13.5")	10"	10.5"	11"	11.5"	12"	12.5"	13"	13.5"	
	Medium (11.5" - 15.5")	11.5"	12"	12.5"	13"	13.5"	14"	14.5"	15"	15.5"
	Long (14" - 18")	14"	14.5"	15"	15.5"	16"	16.5"	17"	17.5"	18"
PRO ELR Low Pivot Setting	PRO ELR Low Pivot Setting Short (10.5" - 15.5")	10.5" - 15.5"								
	PRO ELR Low Pivot Setting Long (13" - 19")	13" - 19"								
PRO ELR High Pivot Setting	PRO ELR High Pivot Setting Short (8.5" - 17")	8.5" - 17"								
	PRO ELR High Pivot Setting Long (11.5" - 17.5")	11.5" - 17.5"								

PRO ELR Cover

Available only on PRO ELR

Gel Pad Cover (K0108) \$425

Calf Straps

Velcro Adjustable (K0038) \$45

Hanger Release

4-Way-In-Line N/C

Not available with PRO ELR

4-Way - Out N/C

Heel Loops

Heel Loops (E0951) \$45

Omit Heel Loops N/C

Seat Height, Casters and Rear Wheels

Select Seat Height corresponding to desired Rear Wheel and Caster Size

Rear Wheel Size	Caster Size		
	4"	5"	6"
16"	Seat Height		
	13.5" ¹	14"	15"
	14" ¹	14.5"	15.5"
	14.5"	15"	16"
	15"	15.5"	16.5"
	15.5"	16"	17"
	16"	16.5"	17.5"
		17"	18"
		17.5"	18.5"
		18"	19"
		18.5"	
		19"	
20"	13"	14"	15"
	13.5"	14.5"	15.5"
	14"	15"	16"
	14.5"	15.5"	16.5"
	15"	16"	17"
	15.5"	16.5"	17.5"
	16"	17"	18"
		17.5"	18.5"
		18"	19"
		18.5"	
		19"	

¹ Available only with Low Profile tires

Rear Wheel Size	Caster Size		
	4"	5"	6"
22"	Seat Height		
	14"	14"	15"
	14.5"	14.5"	15.5"
	15"	15"	16"
	15.5"	15.5"	16.5"
	16"	16"	17"
		16.5"	17.5"
		17"	18"
		17.5"	18.5"
		18"	19"
		18.5"	
		19"	
24"	15"	14.5"	15"
	15.5"	15"	15.5"
	16"	15.5"	16"
		16"	16.5"
		16.5"	17"
		17"	17.5"
		17.5"	18"
		18"	18.5"
		18.5"	19"
		19"	



Phone (800) 981-1540 | Fax (715) 254-0996
sales@kimobility.com | www.kimobility.com

LIBERTY Basic

Casters

		Caster Size			
		Price	4"	5"	6"
Caster Type	1" Poly	N/C			
	1.5" Poly	\$95			

Rear Tire Type

Full Profile

Street Pneumatic (E2211, E2212) _____ N/C
Full Poly 1 3/8" (E2218) _____ N/C
Pneumatic with Airless Insert (E2211, E2213) _____ \$97

Low Profile

Low Poly 1" (E2218) _____ N/C
Available only with 16" Mag Wheels

Wheels

		Rear Wheel Size				
		Price	16"	20"	22"	24"
Rear Wheel Type	Mag	N/C				
	MAXX Spoke	N/C				

Handrims

Not Available with 16" Rear Wheels

Aluminum Anodized _____ N/C
Plastic Coated _____ \$123

Handrim Mount

If no selection is made, Standard Mount will be used

Close Mount _____ N/C
Reduces overall width by .75"

Axles

Quick Release _____ N/C

Wheel Locks

Wheel Lock

Push to Lock _____ N/C
Hemi Wheel Lock - Left (K0108) _____ \$318
Extension Handle included
Hemi Wheel Lock - Right (K0108) _____ \$318
Extension Handle included
Omit Wheel Locks _____ N/C

Wheel Lock Options

Extension Handles (E0961) _____ \$82

Armrests

Armrest

Selection Required

Height Adjustable T-Arm Desk (8.5" - 13.5") (E0973) _____ \$323
Height Adjustable T-Arm Full (8.5" - 13.5") (E0973) _____ \$323

Arm Pad

Standard Arm Pad _____ N/C

Accessories

Anti-tips

Rear Anti-Tips (E0971) _____ \$165
Omit _____ N/C

Transit Wheels

Transit Wheels (K0108) _____ \$318
Not available with Long Depth Seat Frame, Seat Heights less than 14" or greater than 18", 12" or 16" rear wheels, Cane and Crutch Holder, and O2 Holder

Other Accessories

1 1/2" Auto Buckle Pelvic Positioning Belt (E0978) _____ \$87
Seat Belt Mounting Hardware _____ \$16
Belt Mounting Band Clamp _____ \$36
O2 Holder (E2208) _____ \$313
Mounted to the left side. Not available with Fold Down Back.
Cane and Crutch Holder _____ \$210

Frame Color

Black _____ Charcoal Gray _____ Candy Blue _____ Shadow Black _____ Candy Red _____ N/C