



Phone (800) 981-1540 | Fax (715) 254-0996  
 sales@kimobility.com | www.kimobility.com

Effective 05/23/2017

Company: \_\_\_\_\_

Ship To: \_\_\_\_\_

Account Number: \_\_\_\_\_ Date: \_\_\_\_\_

Address: \_\_\_\_\_

Quote  PO Number: \_\_\_\_\_

City: \_\_\_\_\_ State: \_\_\_\_\_

RTS: \_\_\_\_\_

Zip: \_\_\_\_\_ Phone: \_\_\_\_\_

Marked For: \_\_\_\_\_

Email/Fax: \_\_\_\_\_

- Little Wave Flip XP-Rigid Base** E1233, E1161 **\$2750**  
 165 lb weight capacity. Price includes one time free growth kit.  
 Tilt range: 0° - 45°.
- Little Wave Flip XP-Folding Base** E1234, E1161 **\$2850**  
 165 lb weight capacity. Price includes one time free growth kit.  
 Tilt range: 0° - 45°.

*HCPCS codes provided are based on code verification by the PDAC or our interpretation of the code definitions.  
 Proper use of HCPCS codes does not ensure coverage or payment. For coverage information, verify the appropriate payer's coverage policy.*

**Frame**

**Frame Width**

Select box corresponding to desired Frame Width

Frame Width								
<input type="checkbox"/> 10	<input type="checkbox"/> 11	<input type="checkbox"/> 12	<input type="checkbox"/> 13	<input type="checkbox"/> 14	<input type="checkbox"/> 15	<input type="checkbox"/> 16	<input type="checkbox"/> 17	<input type="checkbox"/> 18

**Frame Depth**

Select Frame Depth within desired Frame Depth Range

Frame Depth Range	Frame Depth				
Short 12" - 16"	<input type="checkbox"/> 12	<input type="checkbox"/> 13	<input type="checkbox"/> 14	<input type="checkbox"/> 15	<input type="checkbox"/> 16
Medium 14" - 18"	<input type="checkbox"/> 14	<input type="checkbox"/> 15	<input type="checkbox"/> 16	<input type="checkbox"/> 17	<input type="checkbox"/> 18
Long 16" - 20"	<input type="checkbox"/> 16	<input type="checkbox"/> 17	<input type="checkbox"/> 18	<input type="checkbox"/> 19	<input type="checkbox"/> 20 \$605 E2203

**Transit Option**

Transit System K0108 \_\_\_\_\_ N/C

**Offset Hanger Option**

Offset Hanger \_\_\_\_\_ \$30  
 Available only with Rigid Base. Not available with 12" Frame Depth.

**Tilt Mechanism**

Hand \_\_\_\_\_ N/C  
 Foot \_\_\_\_\_ N/C

**Backrest**

**Back Post / Height Range / Back Height**

Select Back Height in row corresponding to desired Back Post and Height Range

Back Post	Price	Height Range	Back Height								
Stroller Handle Back Post (Folding)	N/C	Short	<input type="checkbox"/> 20								
		Tall	<input type="checkbox"/> 24								
Fixed Height with Adjustable Handle <sup>1</sup> (Folding) <b>New!</b>	N/C	Short	<input type="checkbox"/> 22								
		Tall	<input type="checkbox"/> 26								
Adjustable Height with Adjustable Handle (Folding)	\$25	Short <sup>1</sup>	<input type="checkbox"/> 18	<input type="checkbox"/> 19	<input type="checkbox"/> 20	<input type="checkbox"/> 21	<input type="checkbox"/> 22	<input type="checkbox"/> 23	<input type="checkbox"/> 24	<input type="checkbox"/> 25	
		Tall	<input type="checkbox"/> 20 <sup>1</sup>	<input type="checkbox"/> 21 <sup>1</sup>	<input type="checkbox"/> 22 <sup>1</sup>	<input type="checkbox"/> 23	<input type="checkbox"/> 24	<input type="checkbox"/> 25	<input type="checkbox"/> 26	<input type="checkbox"/> 27	

<sup>1</sup> Not available with Seating Dynamics Rocker Back

**Stroller Handle Back Post Option**

Available only with Stroller Handle Back Post

Removable Stroller Handle Extension \_\_\_\_\_ \$85

**Other Backrest Options**

Non Folding Back \_\_\_\_\_ N/C  
 Standard with Seating Dynamics Rocker Back

Seating Dynamics Rocker Back K0108 \_\_\_\_\_ \$1995  
 Not available with Flip Up Armrests. Available only with Non Folding Back.

Adjustable Height Rigidizer Bar \_\_\_\_\_ \$75

## Footrest

### Hanger / Front Frame Bend

- 60° \_\_\_\_\_ N/C  
*Available only with Extension Tube Hangers*
- 70° \_\_\_\_\_ N/C
- 80° \_\_\_\_\_ N/C  
*Not available with 8" Casters*
- 90° \_\_\_\_\_ N/C  
*Available only with Front Mount Hangers. Not available with 8" Casters.*
- Elevating Leg Rests (ELR) E0990 \_\_\_\_\_ \$295  
*Not available with Frame Widths less than 14" without Offset Hangers or Frame Widths less than 12" with Offset Hangers*
- Pro Elevating Leg Rests (PRO ELR) E0990 **New!** \_\_\_\_\_ \$350  
*Not available with Frame Widths less than 13" without Offset Hangers or Frame Widths less than 12" with Offset Hangers*
- Omit Hangers \_\_\_\_\_ - \$20

### Footplate

- Composite \_\_\_\_\_ N/C  
*Not available with Front Mount Hangers, Frame Widths less than 14" without Offset Hangers, Frame Widths less than 12" with Offset Hangers or PRO ELR*
- Composite Angle Adjustable K0040 \_\_\_\_\_ \$185  
*Not available with Front Mount Hangers, ELR or 10"-12" Frame Widths. 13" and 14" Frame Widths require Offset Hangers.*
- Aluminum Angle Adjustable K0040 \_\_\_\_\_ \$190
- One Piece Flip Up Angle Adjustable K0040 \_\_\_\_\_ \$190  
*Not available with Front Mount Hangers, ELR, PRO ELR, Frame Widths less than 14" without Offset Hangers or Frame Widths less than 12" with Offset Hangers*

### Footrest / Extension / Seat to Footrest Length

Select Seat to Footrest Length in row corresponding to desired Footrest and Extension. Footrest Length of at least 2.5" shorter than Front Seat Height recommended. Extension Tube Hanger measurements based on 70° Hanger Bend.

Footrest	Extension	Seat to Footrest Length										
Extension Tube Hangers (60°, 70°, 80°)	Short	<input type="checkbox"/> 10	<input type="checkbox"/> 10.5	<input type="checkbox"/> 11	<input type="checkbox"/> 11.5	<input type="checkbox"/> 12	<input type="checkbox"/> 12.5	<input type="checkbox"/> 13	<input type="checkbox"/> 13.5			
	Medium	<input type="checkbox"/> 11.5	<input type="checkbox"/> 12	<input type="checkbox"/> 12.5	<input type="checkbox"/> 13	<input type="checkbox"/> 13.5	<input type="checkbox"/> 14	<input type="checkbox"/> 14.5	<input type="checkbox"/> 15	<input type="checkbox"/> 15.5		
	Long	<input type="checkbox"/> 14	<input type="checkbox"/> 14.5	<input type="checkbox"/> 15	<input type="checkbox"/> 15.5	<input type="checkbox"/> 16	<input type="checkbox"/> 16.5	<input type="checkbox"/> 17	<input type="checkbox"/> 17.5	<input type="checkbox"/> 18		
Front Mount Hangers <sup>1</sup>	<b>Standard Selections</b>											
	Short (70°, 80°, 90°)	<input type="checkbox"/> 5	<input type="checkbox"/> 5.5	<input type="checkbox"/> 6	<input type="checkbox"/> 6.5					<input type="checkbox"/> 8.5	<input type="checkbox"/> 9	<input type="checkbox"/> 9.5
		<input type="checkbox"/> 7	<input type="checkbox"/> 7.5	<input type="checkbox"/> 8						<input type="checkbox"/> 10	<input type="checkbox"/> 10.5	<input type="checkbox"/> 11
	Medium (70°, 80°, 90°)	<input type="checkbox"/> 5	<input type="checkbox"/> 5.5	<input type="checkbox"/> 6	<input type="checkbox"/> 6.5	<input type="checkbox"/> 7	<input type="checkbox"/> 7.5			<input type="checkbox"/> 10.5	<input type="checkbox"/> 11	<input type="checkbox"/> 11.5
		<input type="checkbox"/> 8	<input type="checkbox"/> 8.5	<input type="checkbox"/> 9	<input type="checkbox"/> 9.5	<input type="checkbox"/> 10				<input type="checkbox"/> 12	<input type="checkbox"/> 12.5	<input type="checkbox"/> 13
	Long (70°, 80°)	<input type="checkbox"/> 5	<input type="checkbox"/> 5.5	<input type="checkbox"/> 6	<input type="checkbox"/> 6.5	<input type="checkbox"/> 7	<input type="checkbox"/> 7.5	<input type="checkbox"/> 8	<input type="checkbox"/> 8.5	<input type="checkbox"/> 12.5	<input type="checkbox"/> 13	<input type="checkbox"/> 13.5
<input type="checkbox"/> 9		<input type="checkbox"/> 9.5	<input type="checkbox"/> 10	<input type="checkbox"/> 10.5	<input type="checkbox"/> 11	<input type="checkbox"/> 11.5	<input type="checkbox"/> 12		<input type="checkbox"/> 14	<input type="checkbox"/> 14.5	<input type="checkbox"/> 15	
Elevating Leg Rests (ELR)	Short	<input type="checkbox"/> 10.5	<input type="checkbox"/> 11.5	<input type="checkbox"/> 12.5	<input type="checkbox"/> 13.5	<input type="checkbox"/> 14.5						
	Long	<input type="checkbox"/> 11.5	<input type="checkbox"/> 12.5	<input type="checkbox"/> 13.5	<input type="checkbox"/> 14.5	<input type="checkbox"/> 15.5	<input type="checkbox"/> 16.5					
Extended Length ELR	Short	<input type="checkbox"/> 16.5	<input type="checkbox"/> 17.5	<input type="checkbox"/> 18.5	<input type="checkbox"/> 19.5	<input type="checkbox"/> 20.5						
	Long	<input type="checkbox"/> 16.5	<input type="checkbox"/> 17.5	<input type="checkbox"/> 18.5	<input type="checkbox"/> 19.5	<input type="checkbox"/> 20.5	<input type="checkbox"/> 21.5	<input type="checkbox"/> 22.5				
PRO ELR Low Pivot	Short	<input type="checkbox"/> 10.5 - 15.5										
	Long	<input type="checkbox"/> 13 - 19										
PRO ELR High Pivot	Short	<input type="checkbox"/> 8.5 - 17										
	Long	<input type="checkbox"/> 11.5 - 17.5										

<sup>1</sup> Available only with Aluminum Angle Adjustable Footplates. See Minimum Front Seat Height chart below.

Minimum Front Seat Height with Front Mount Hangers						
Hanger Length	Short		Medium		Long	
	70°	80/90°	70°	80/90°	70°	80°
Minimum Front Seat Height with Standard Selection	13"	13"	13"	13.5"	14.5"	15.5"
Minimum Front Seat Height with 3" Extension	13.5"	14"	15.5"	16"	17"	18"

### 3" Extension

Available only with Front Mount Hangers. Provided as parts.

- 3" Extension for Front Mount Hanger \_\_\_\_\_ \$55

**Footrest continued**

**Heel Loop**

- Heel Loops E0951 \_\_\_\_\_ \$30  
*Not available with One Piece Flip Up Angle Adjustable Footplates, Frame Widths less than 14" without Offset Hangers or Frame Widths less than 12" with Offset Hangers*
- Omit Heel Loops \_\_\_\_\_ N/C

**Calf Strap**

- Velcro Adjustable K0038 \_\_\_\_\_ \$35
- Padded Velcro Adjustable K0038 \_\_\_\_\_ \$68
- Bodypoint Padded AeroMesh K0038 \_\_\_\_\_ \$102



**Residual Limb Support**

- Residual Limb Support E1020
    - Left \_\_\_\_\_ \$395
    - Right \_\_\_\_\_ \$395
    - Omit Left Hanger \_ - \$10
    - Omit Right Hanger \_ - \$10
- |                                 |                                   |                                 |                                   |
|---------------------------------|-----------------------------------|---------------------------------|-----------------------------------|
| <b>Cover</b>                    |                                   | <b>Cover</b>                    |                                   |
| <input type="checkbox"/> Smooth | <input type="checkbox"/> Air Mesh | <input type="checkbox"/> Smooth | <input type="checkbox"/> Air Mesh |
| <b>Pad Size</b>                 |                                   | <b>Pad Size</b>                 |                                   |
| <input type="checkbox"/> 5 x 8  | <input type="checkbox"/> 9 x 14   | <input type="checkbox"/> 5 x 8  | <input type="checkbox"/> 9 x 14   |
| <input type="checkbox"/> 7 x 10 | <input type="checkbox"/> 11 x 10  | <input type="checkbox"/> 7 x 10 | <input type="checkbox"/> 11 x 10  |
| <input type="checkbox"/> 9 x 10 | <input type="checkbox"/> 11 x 14  | <input type="checkbox"/> 9 x 10 | <input type="checkbox"/> 11 x 14  |

**Hanger Release**

*If no selection is made 4-Way - In line will be used with hangers or Classic will be used with ELR and PRO ELR*

- Classic \_\_\_\_\_ Std
- 4-Way - In line **New!** \_\_\_\_\_ N/C  
*See reference A to the right. Allows chair to be folded without removing footrest hanger. Not available with ELR or PRO ELR.*
- 4-Way - Out **New!** \_\_\_\_\_ N/C  
*See reference B to the right*



## Seat Height, Casters and Rear Wheels

### Seat Height / Caster Size / Rear Wheel Size

Select Seat Height corresponding to desired Rear Wheel and Caster Size

Rear Wheel Size	Caster Size					
	5"	6"	6" x 2"	New! 7"	8"	8" x 2"
Seat Height						
12" with Pneumatic or Pneumatic with Airless Insert Tire	<input type="checkbox"/> 13.5					
	<input type="checkbox"/> 14	<input type="checkbox"/> 14 <sup>4</sup>	<input type="checkbox"/> 14 <sup>4</sup>			
	<input type="checkbox"/> 14.5	<input type="checkbox"/> 14.5	<input type="checkbox"/> 14.5			
	<input type="checkbox"/> 15	<input type="checkbox"/> 15	<input type="checkbox"/> 15	<input type="checkbox"/> 15 <sup>4</sup>		
	<input type="checkbox"/> 15.5	<input type="checkbox"/> 15.5	<input type="checkbox"/> 15.5	<input type="checkbox"/> 15.5		
	<input type="checkbox"/> 16	<input type="checkbox"/> 16	<input type="checkbox"/> 16	<input type="checkbox"/> 16		
	<input type="checkbox"/> 16.5	<input type="checkbox"/> 16.5	<input type="checkbox"/> 16.5	<input type="checkbox"/> 16.5	<input type="checkbox"/> 16.5	<input type="checkbox"/> 16.5
	<input type="checkbox"/> 17	<input type="checkbox"/> 17	<input type="checkbox"/> 17	<input type="checkbox"/> 17 <sup>3</sup>	<input type="checkbox"/> 17 <sup>3</sup>	<input type="checkbox"/> 17 <sup>3</sup>
	<input type="checkbox"/> 17.5	<input type="checkbox"/> 17.5	<input type="checkbox"/> 17.5	<input type="checkbox"/> 17.5	<input type="checkbox"/> 17.5	<input type="checkbox"/> 17.5
12" with Full Poly Tire	<input type="checkbox"/> 13.5					
	<input type="checkbox"/> 14					
	<input type="checkbox"/> 14.5	<input type="checkbox"/> 14.5	<input type="checkbox"/> 14.5			
	<input type="checkbox"/> 15	<input type="checkbox"/> 15	<input type="checkbox"/> 15			
	<input type="checkbox"/> 15.5	<input type="checkbox"/> 15.5	<input type="checkbox"/> 15.5	<input type="checkbox"/> 15.5		
	<input type="checkbox"/> 16	<input type="checkbox"/> 16	<input type="checkbox"/> 16	<input type="checkbox"/> 16	<input type="checkbox"/> 16 <sup>4</sup>	<input type="checkbox"/> 16 <sup>4</sup>
	<input type="checkbox"/> 16.5	<input type="checkbox"/> 16.5	<input type="checkbox"/> 16.5	<input type="checkbox"/> 16.5	<input type="checkbox"/> 16.5	<input type="checkbox"/> 16.5
	<input type="checkbox"/> 17	<input type="checkbox"/> 17	<input type="checkbox"/> 17	<input type="checkbox"/> 17 <sup>3</sup>	<input type="checkbox"/> 17 <sup>3</sup>	<input type="checkbox"/> 17 <sup>3</sup>
	<input type="checkbox"/> 17.5	<input type="checkbox"/> 17.5	<input type="checkbox"/> 17.5	<input type="checkbox"/> 17.5	<input type="checkbox"/> 17.5	<input type="checkbox"/> 17.5
16" with Pneumatic, Pneumatic with Airless Insert or Full Poly Tire	<input type="checkbox"/> 15	<input type="checkbox"/> 15	<input type="checkbox"/> 15	<input type="checkbox"/> 15 <sup>4</sup>		
	<input type="checkbox"/> 15.5	<input type="checkbox"/> 15.5	<input type="checkbox"/> 15.5	<input type="checkbox"/> 15.5		
	<input type="checkbox"/> 16	<input type="checkbox"/> 16	<input type="checkbox"/> 16	<input type="checkbox"/> 16		<input type="checkbox"/> 16 <sup>4</sup>
	<input type="checkbox"/> 16.5	<input type="checkbox"/> 16.5	<input type="checkbox"/> 16.5	<input type="checkbox"/> 16.5	<input type="checkbox"/> 16.5	<input type="checkbox"/> 16.5
	<input type="checkbox"/> 17	<input type="checkbox"/> 17	<input type="checkbox"/> 17	<input type="checkbox"/> 17	<input type="checkbox"/> 17	<input type="checkbox"/> 17
	<input type="checkbox"/> 17.5	<input type="checkbox"/> 17.5	<input type="checkbox"/> 17.5	<input type="checkbox"/> 17.5	<input type="checkbox"/> 17.5	<input type="checkbox"/> 17.5
	<input type="checkbox"/> 18	<input type="checkbox"/> 18	<input type="checkbox"/> 18	<input type="checkbox"/> 18	<input type="checkbox"/> 18	<input type="checkbox"/> 18
	<input type="checkbox"/> 18.5 <sup>4</sup>	<input type="checkbox"/> 18.5	<input type="checkbox"/> 18.5	<input type="checkbox"/> 18.5	<input type="checkbox"/> 18.5	<input type="checkbox"/> 18.5
	<input type="checkbox"/> 19	<input type="checkbox"/> 19	<input type="checkbox"/> 19	<input type="checkbox"/> 19	<input type="checkbox"/> 19	<input type="checkbox"/> 19
<input type="checkbox"/> 19.5 <sup>2</sup>	<input type="checkbox"/> 19.5 <sup>4</sup>	<input type="checkbox"/> 19.5 <sup>4</sup>			<input type="checkbox"/> 19.5 <sup>3</sup>	
16" with Low Poly Tire	<input type="checkbox"/> 14	<input type="checkbox"/> 14 <sup>4</sup>	<input type="checkbox"/> 14 <sup>4</sup>			
	<input type="checkbox"/> 14.5	<input type="checkbox"/> 14.5	<input type="checkbox"/> 14.5			
	<input type="checkbox"/> 15	<input type="checkbox"/> 15	<input type="checkbox"/> 15			
	<input type="checkbox"/> 15.5	<input type="checkbox"/> 15.5	<input type="checkbox"/> 15.5	<input type="checkbox"/> 15.5		
	<input type="checkbox"/> 16	<input type="checkbox"/> 16	<input type="checkbox"/> 16	<input type="checkbox"/> 16		<input type="checkbox"/> 16 <sup>4</sup>
	<input type="checkbox"/> 16.5	<input type="checkbox"/> 16.5	<input type="checkbox"/> 16.5	<input type="checkbox"/> 16.5	<input type="checkbox"/> 16.5	<input type="checkbox"/> 16.5
	<input type="checkbox"/> 17	<input type="checkbox"/> 17	<input type="checkbox"/> 17	<input type="checkbox"/> 17	<input type="checkbox"/> 17	<input type="checkbox"/> 17
	<input type="checkbox"/> 17.5	<input type="checkbox"/> 17.5	<input type="checkbox"/> 17.5	<input type="checkbox"/> 17.5	<input type="checkbox"/> 17.5	<input type="checkbox"/> 17.5
	<input type="checkbox"/> 18	<input type="checkbox"/> 18	<input type="checkbox"/> 18	<input type="checkbox"/> 18	<input type="checkbox"/> 18	<input type="checkbox"/> 18

<sup>1</sup> Self Propel Axle Plate selections limit tilt range  
<sup>2</sup> Not available with 12" Frame Depth  
<sup>3</sup> Not available with 12-16" Frame Depth  
<sup>4</sup> Not available with 17-20" Frame Depth  
<sup>5</sup> Not available in combination with Spoke wheels and Full Poly Rear Tire Type  
<sup>6</sup> Not available in combination with Spoke wheels, Full Poly Rear Tire Type and 17-20" Frame Depth  
<sup>7</sup> Not available in combination with Mag wheels and 17-20" Frame Depth  
<sup>8</sup> Not available in combination with Spoke wheels, Pneumatic tires and 17-20" Frame Depth  
<sup>9</sup> 17-20" Frame Depth requires Spoke wheels and Full Poly Rear Tire Type

Rear Wheel Size	Caster Size					
	5"	6"	6" x 2"	New! 7"	8"	8" x 2"
Seat Height						
16" Drum Brake	<input type="checkbox"/> 15	<input type="checkbox"/> 15	<input type="checkbox"/> 15			
	<input type="checkbox"/> 15.5	<input type="checkbox"/> 15.5	<input type="checkbox"/> 15.5	<input type="checkbox"/> 15.5		
	<input type="checkbox"/> 16	<input type="checkbox"/> 16	<input type="checkbox"/> 16	<input type="checkbox"/> 16	<input type="checkbox"/> 16	<input type="checkbox"/> 16
	<input type="checkbox"/> 16.5	<input type="checkbox"/> 16.5	<input type="checkbox"/> 16.5	<input type="checkbox"/> 16.5	<input type="checkbox"/> 16.5	<input type="checkbox"/> 16.5
	<input type="checkbox"/> 17	<input type="checkbox"/> 17	<input type="checkbox"/> 17	<input type="checkbox"/> 17	<input type="checkbox"/> 17	<input type="checkbox"/> 17
	<input type="checkbox"/> 17.5	<input type="checkbox"/> 17.5	<input type="checkbox"/> 17.5	<input type="checkbox"/> 17.5	<input type="checkbox"/> 17.5	<input type="checkbox"/> 17.5 <sup>3</sup>
	<input type="checkbox"/> 18	<input type="checkbox"/> 18	<input type="checkbox"/> 18	<input type="checkbox"/> 18	<input type="checkbox"/> 18	<input type="checkbox"/> 18
	<input type="checkbox"/> 18.5 <sup>4</sup>	<input type="checkbox"/> 18.5	<input type="checkbox"/> 18.5 <sup>4</sup>	<input type="checkbox"/> 18.5 <sup>4</sup>	<input type="checkbox"/> 18.5	<input type="checkbox"/> 18.5 <sup>3</sup>
	<input type="checkbox"/> 19	<input type="checkbox"/> 19	<input type="checkbox"/> 19	<input type="checkbox"/> 19	<input type="checkbox"/> 19	<input type="checkbox"/> 19
	<input type="checkbox"/> 19.5	<input type="checkbox"/> 19.5	<input type="checkbox"/> 19.5	<input type="checkbox"/> 19.5	<input type="checkbox"/> 19.5	<input type="checkbox"/> 19.5
20"	<input type="checkbox"/> 17	<input type="checkbox"/> 17	<input type="checkbox"/> 17	<input type="checkbox"/> 17	<input type="checkbox"/> 17	<input type="checkbox"/> 17 <sup>3</sup>
	<input type="checkbox"/> 17.5	<input type="checkbox"/> 17.5	<input type="checkbox"/> 17.5	<input type="checkbox"/> 17.5	<input type="checkbox"/> 17.5	<input type="checkbox"/> 17.5
	<input type="checkbox"/> 18	<input type="checkbox"/> 18	<input type="checkbox"/> 18	<input type="checkbox"/> 18	<input type="checkbox"/> 18	<input type="checkbox"/> 18
	<input type="checkbox"/> 18.5	<input type="checkbox"/> 18.5	<input type="checkbox"/> 18.5	<input type="checkbox"/> 18.5	<input type="checkbox"/> 18.5	<input type="checkbox"/> 18.5
	<input type="checkbox"/> 19	<input type="checkbox"/> 19	<input type="checkbox"/> 19	<input type="checkbox"/> 19	<input type="checkbox"/> 19	<input type="checkbox"/> 19
	<input type="checkbox"/> 19.5	<input type="checkbox"/> 19.5	<input type="checkbox"/> 19.5	<input type="checkbox"/> 19.5	<input type="checkbox"/> 19.5	<input type="checkbox"/> 19.5
	<input type="checkbox"/> 20	<input type="checkbox"/> 20 <sup>4</sup>	<input type="checkbox"/> 20 <sup>4</sup>	<input type="checkbox"/> 20	<input type="checkbox"/> 20	<input type="checkbox"/> 20
20" <sup>1</sup> Self Propel	<input type="checkbox"/> 15.5	<input type="checkbox"/> 15.5	<input type="checkbox"/> 15.5	<input type="checkbox"/> 15.5		
	<input type="checkbox"/> 16	<input type="checkbox"/> 16	<input type="checkbox"/> 16	<input type="checkbox"/> 16	<input type="checkbox"/> 16	<input type="checkbox"/> 16
	<input type="checkbox"/> 16.5	<input type="checkbox"/> 16.5	<input type="checkbox"/> 16.5	<input type="checkbox"/> 16.5	<input type="checkbox"/> 16.5	<input type="checkbox"/> 16.5
	<input type="checkbox"/> 17	<input type="checkbox"/> 17	<input type="checkbox"/> 17	<input type="checkbox"/> 17	<input type="checkbox"/> 17	<input type="checkbox"/> 17
	<input type="checkbox"/> 17.5	<input type="checkbox"/> 17.5	<input type="checkbox"/> 17.5	<input type="checkbox"/> 17.5	<input type="checkbox"/> 17.5	<input type="checkbox"/> 17.5
	<input type="checkbox"/> 18	<input type="checkbox"/> 18	<input type="checkbox"/> 18	<input type="checkbox"/> 18	<input type="checkbox"/> 18	<input type="checkbox"/> 18
	<input type="checkbox"/> 18.5 <sup>4</sup>	<input type="checkbox"/> 18.5	<input type="checkbox"/> 18.5	<input type="checkbox"/> 18.5 <sup>4</sup>	<input type="checkbox"/> 18.5	<input type="checkbox"/> 18.5
22" <sup>1</sup> Self Propel	<input type="checkbox"/> 16.5	<input type="checkbox"/> 16.5	<input type="checkbox"/> 16.5	<input type="checkbox"/> 16.5	<input type="checkbox"/> 16.5	<input type="checkbox"/> 16.5
	<input type="checkbox"/> 17	<input type="checkbox"/> 17	<input type="checkbox"/> 17	<input type="checkbox"/> 17	<input type="checkbox"/> 17	<input type="checkbox"/> 17
	<input type="checkbox"/> 17.5	<input type="checkbox"/> 17.5	<input type="checkbox"/> 17.5	<input type="checkbox"/> 17.5	<input type="checkbox"/> 17.5	<input type="checkbox"/> 17.5
	<input type="checkbox"/> 18	<input type="checkbox"/> 18	<input type="checkbox"/> 18	<input type="checkbox"/> 18	<input type="checkbox"/> 18	<input type="checkbox"/> 18
	<input type="checkbox"/> 19	<input type="checkbox"/> 19	<input type="checkbox"/> 19	<input type="checkbox"/> 19	<input type="checkbox"/> 19	<input type="checkbox"/> 19
24" <sup>1</sup> Self Propel	<input type="checkbox"/> 17	<input type="checkbox"/> 17	<input type="checkbox"/> 17	<input type="checkbox"/> 17	<input type="checkbox"/> 17	<input type="checkbox"/> 17
	<input type="checkbox"/> 17.5	<input type="checkbox"/> 17.5	<input type="checkbox"/> 17.5	<input type="checkbox"/> 17.5	<input type="checkbox"/> 17.5	<input type="checkbox"/> 17.5
	<input type="checkbox"/> 18	<input type="checkbox"/> 18	<input type="checkbox"/> 18	<input type="checkbox"/> 18	<input type="checkbox"/> 18	<input type="checkbox"/> 18
	<input type="checkbox"/> 18.5	<input type="checkbox"/> 18.5	<input type="checkbox"/> 18.5	<input type="checkbox"/> 18.5	<input type="checkbox"/> 18.5	<input type="checkbox"/> 18.5
	<input type="checkbox"/> 19	<input type="checkbox"/> 19	<input type="checkbox"/> 19	<input type="checkbox"/> 19	<input type="checkbox"/> 19	<input type="checkbox"/> 19
<input type="checkbox"/> 19.5 <sup>7,8</sup>	<input type="checkbox"/> 19.5	<input type="checkbox"/> 19.5	<input type="checkbox"/> 19.5 <sup>9</sup>	<input type="checkbox"/> 19.5 <sup>7,8</sup>	<input type="checkbox"/> 19.5	
<input type="checkbox"/> 20 <sup>5</sup>	<input type="checkbox"/> 20 <sup>6</sup>	<input type="checkbox"/> 20 <sup>6</sup>	<input type="checkbox"/> 20	<input type="checkbox"/> 20	<input type="checkbox"/> 20	

**Seat Height, Casters and Rear Wheels** continued

Caster Type	
<input type="checkbox"/> .75" Lighted Rollerblade	\$70
<i>Not available with 6", 7" and 8" Caster Sizes</i>	
<input type="checkbox"/> 1" Poly	N/C
<input type="checkbox"/> 1" Poly Aluminum	\$110
<i>Not available with 6", 7" and 8" Caster Sizes</i>	
<input type="checkbox"/> 1" Pneumatic	\$70
<i>Not available with 5" and 7" Caster Sizes</i>	
<input type="checkbox"/> 1.5" Poly	\$50
<i>Not available with 7" and 8" Caster Sizes</i>	
<input type="checkbox"/> 1.5" Soft Roll Aluminum	\$175
<i>Not available with 7" and 8" Caster Sizes</i>	
<input type="checkbox"/> 2" Poly	\$70
<i>Not available with 5", 7" and 8" Caster Sizes</i>	
<input type="checkbox"/> 2" Pneumatic w/ Foam Insert E2217	\$225
<i>Not available with 5", 6" and 7" Caster Sizes</i>	

Axle	
<b>No selection needed with Drum Brake</b>	
<input type="checkbox"/> Fixed	N/C
<i>Not available with Folding Base Frame and Standard Wheel Locks</i>	
<input type="checkbox"/> Quick Release	N/C

Rear Wheel Type	
<input type="checkbox"/> Mag (Black)	N/C
<input type="checkbox"/> Colored Mag	\$100
<i>Available only in 20" and 22" rear wheels</i>	
Color: <input type="checkbox"/> Florescent Red <input type="checkbox"/> White <input type="checkbox"/> Blue <input type="checkbox"/> Pink	
<input type="checkbox"/> Spoke	\$200
<i>Not available in 12" or 16" rear wheels</i>	

**Wheel Locks**

Wheel Lock	
<input type="checkbox"/> Push to Lock	N/C
<i>Rear wheels must be removed to fold chair with this wheel lock selection.</i>	
<input type="checkbox"/> Pull to Lock	N/C
<i>Not available with 12" or 16" Full Poly tires. Rear wheels must be removed to fold chair with this wheel lock selection.</i>	
<input type="checkbox"/> Attendant Foot Lock	N/C
<i>Available only with 12" and 16" Mag wheels.</i>	
<input type="checkbox"/> Drum Brake E2228	\$1200
<i>Not available with Hand Tilt Mechanism. Available only with 16" rear wheel. Axle is Quick Release.</i>	
<input type="checkbox"/> Omit Wheel Locks	N/C

**Armrests**

Armrests Please select Arm Pad option	
<input type="checkbox"/> Height Adjustable Pediatric Desk (6.5"- 9") E0973	\$185
<i>Best fit with Frame Depths less than 14"</i>	
<input type="checkbox"/> Height Adjustable Pediatric Desk Tall (8"-10.5") E0973	\$185
<i>Best fit with Frame Depths less than 14"</i>	
<input type="checkbox"/> Height Adjustable T-Arm (9.25"-14.25") E0973	\$265
<i>Not available with 12" Frame Depth</i>	
<input type="checkbox"/> Height Adjustable T-Arm Low (7.5"-11.5") E0973	\$265
<i>Not available with 12" Frame Depth</i>	
<input type="checkbox"/> Angle Adj. Locking Extendable Flip Up (9.5"-15") E0973	\$265
<i>Not available with Seating Dynamics Rocker Back</i>	

Rear Tire Type	
<input type="checkbox"/> Full Poly E2218	N/C
<input type="checkbox"/> Low Poly E2218 <b>New!</b>	N/C
<i>Available only with 16" rear wheel</i>	
<input type="checkbox"/> Pneumatic E2211, E2212	\$85
<input type="checkbox"/> Pneumatic with Airless Insert E2211, E2213	\$85

Handrim	
<b>Handrims not available on 12" or 16" rear wheels</b>	
Standard Handrims	
<input type="checkbox"/> Aluminum Anodized	Std
<input type="checkbox"/> Plastic Coated	\$95
<input type="checkbox"/> Projection (8 Vertical Projections)	\$295
<i>Available on 24" rear wheel</i>	
<input type="checkbox"/> Omit Handrims	N/C

Natural Fit Handrims	
<b>Not available on 20" rear wheel or 22" Mag wheels. Not recommended with Mag wheels.</b>	
<input type="checkbox"/> Natural Fit LT Regular Grip	\$445
<input type="checkbox"/> Natural Fit Standard Regular Grip	\$445
<input type="checkbox"/> Natural Fit LT Super Grip	\$495
<input type="checkbox"/> Natural Fit Standard Super Grip	\$495
<input type="checkbox"/> Natural Fit LT Omit Grip	\$345
<input type="checkbox"/> Natural Fit Standard Omit Grip	\$345

Handrim Option	
<b>Not available on Plastic Coated or Projection Handrims</b>	
<input type="checkbox"/> 3M Micro Grip	\$65

Wheel Lock continued	
<input type="checkbox"/> Dual - Attendant Foot Lock	\$70
<i>Available only with 12" and 16" Mag wheels. Rear wheels must be removed to fold chair with this wheel lock selection.</i>	
<input type="checkbox"/> Push to Lock <input type="checkbox"/> Pull to Lock	
<input type="checkbox"/> Dual - Drum Brake	\$1270
<i>Not available with Hand Tilt Mechanism. Available only with 16" rear wheel. Axle is Quick Release. Rear wheels must be removed to fold chair with this wheel lock selection.</i>	
<input type="checkbox"/> Push to Lock <input type="checkbox"/> Pull to Lock	

Wheel Lock Options	
<input type="checkbox"/> Extension Handles E0961	\$70
<i>Available only with Push to Lock or Pull to Lock Wheel Locks.</i>	

Arm Pad	
<b>Arm Pads not available with Tubular Flip Up Armrest. Full Arm Pad not available with Height Adjustable Pediatric Desk Armrests.</b>	
Standard Arm Pad	N/C <input type="checkbox"/> Desk <input type="checkbox"/> Full
Waterfall Arm Pad	\$25 <input type="checkbox"/> Desk <input type="checkbox"/> Full
<i>Not available with Angle Adjustable Locking Flip Up Extendable armrests</i>	
Foam Grip	N/C <input type="checkbox"/> Desk <input type="checkbox"/> Full
<i>Available only with Angle Adjustable Locking Extendable Flip Up Armrest</i>	

Armrests Continued	
<input type="checkbox"/> Tubular Flip Up	N/C
<i>Not available with Seating Dynamics Rocker Back</i>	
<input type="checkbox"/> Omit Armrests	N/C

## Accessories

### Removable Side Guard

*Not available with Height Adjustable Armrests*

- Adult Composite Side Guards (14" - 16") K0108 \_\_\_\_\_ \$140
- Pediatric Composite Side Guards (10" - 14") K0108 \_\_\_\_\_ \$140

### Positioning Belt

- 1" Auto Buckle E0978 \_\_\_\_\_ \$70
  - 1" Padded Auto Buckle E0978 \_\_\_\_\_ \$155
  - 1 1/2" Auto Buckle E0978 \_\_\_\_\_ \$70
  - 1 1/2" Padded Auto Buckle E0978 \_\_\_\_\_ \$155
  - 2" Airline Buckle E0978 \_\_\_\_\_ \$110
  - Bodypoint Evoflex - Pelvic Stabilizer E0978 \_\_\_\_\_ \$183
- Size:  Small (16-24")  Medium (19-29")  Large (25-37")

### Belt Mounting Kit

- Bodypoint 1" Band Clamp Belt Mounting Kit \_\_\_\_\_ \$105
- Seat Belt Mounting Hardware \_\_\_\_\_ \$30

### Anti-Tip

- Rear Anti-Tips E0971 \_\_\_\_\_ \$140

### Seat Pan

- Growing Solid Seat Pan E2231 \_\_\_\_\_ \$150

### Spoke Guards

- Clear Spoke Guards K0065 \_\_\_\_\_ \$115  
*Not available with Mag wheels*
- Black Spoke Guards K0065 \_\_\_\_\_ \$115  
*Available only with 22" and 24" Spoke wheels*
- Designer Spoke Guards K0065 \_\_\_\_\_ \$200  
*Available only with 20" and 22" Spoke wheels*

**Design:**

- Blue Monster
- Robo Fun
- Tutu Cute
- Fairy Garden
- Rockin' Reptiles
- Tuxedo Time
- Havin' a Hoot
- Super Puppies

### Other Accessories

- Tool Kit \_\_\_\_\_ \$55
- Removable Underseat Pouch \_\_\_\_\_ \$50
- Pediatric Ki Mobility Backpack **New!** \_\_\_\_\_ \$65
- Ki Mobility Backpack \_\_\_\_\_ \$65
- Black Neoprene Impact Guards \_\_\_\_\_ \$85
- Tilt Stop (Pair) K0108 \_\_\_\_\_ \$30

## Color

### Frame and Accent

*Select desired frame and accent colors. Accent color is used on side panel; if no accent color is selected, frame color will be used.*

Gloss Finish		Frame	Accent	Matte Finish		Frame	Accent
Black _____	N/C	<input type="checkbox"/>	<input type="checkbox"/>	Army Green _____	N/C	<input type="checkbox"/>	<input type="checkbox"/>
Blue Skies _____	N/C	<input type="checkbox"/>	<input type="checkbox"/>	Charcoal Gray _____	N/C	<input type="checkbox"/>	<input type="checkbox"/>
Bubblegum Pink _____	N/C	<input type="checkbox"/>	<input type="checkbox"/>	Emerald Green _____	N/C	<input type="checkbox"/>	<input type="checkbox"/>
Candy Blue _____	N/C	<input type="checkbox"/>	<input type="checkbox"/>	Iridescent Burgundy _____	N/C	<input type="checkbox"/>	<input type="checkbox"/>
Candy Purple _____	N/C	<input type="checkbox"/>	<input type="checkbox"/>	Iridescent Navy _____	N/C	<input type="checkbox"/>	<input type="checkbox"/>
Candy Red _____	N/C	<input type="checkbox"/>	<input type="checkbox"/>	Retro Red _____	N/C	<input type="checkbox"/>	<input type="checkbox"/>
Copper _____	N/C	<input type="checkbox"/>	<input type="checkbox"/>	Royal Blue _____	N/C	<input type="checkbox"/>	<input type="checkbox"/>
Fuchsia _____	N/C	<input type="checkbox"/>	<input type="checkbox"/>	Shadow Black _____	N/C	<input type="checkbox"/>	<input type="checkbox"/>
Granny Smith _____	N/C	<input type="checkbox"/>	<input type="checkbox"/>	Smooth White _____	N/C	<input type="checkbox"/>	<input type="checkbox"/>
Silver _____	N/C	<input type="checkbox"/>	<input type="checkbox"/>				
Teal _____	N/C	<input type="checkbox"/>	<input type="checkbox"/>				
Yellow _____	N/C	<input type="checkbox"/>	<input type="checkbox"/>				

## Cushions

### Matrix Seat Cushions

*Cushion size is determined by chair dimensions. Certain sizes will have higher price and longer lead time; visit [www.motionconcepts.com/matrix\\_international.html](http://www.motionconcepts.com/matrix_international.html) for size specifications.*

#### STEP 1: Choose Cushion

Kidabra Cushion \_\_\_\_\_ \$510  
*Weight capacity is 150 lbs*

#### STEP 2: Choose Cushion Size

Match to Chair Size  
 Other Size - Provide Part Number \_\_\_\_\_

## Backs

### Matrix Backs

*Matrix Backs come standard with quick release mounting hardware*

Mini Back \_\_\_\_\_ \$750  
*Not available with 16" Seat Width. Back Height of 14" not available with 10" and 11" Seat Widths. Comes standard with Kidabra cover.*

#### Back Height:

7"     10"     12"     14"

### Therafin Headrest

*Available with Matrix Back Systems*

Headrest, Multi-Axis Offset Hardware with Smooth/Soft Pad (E1028, E0955) 34187 \_\_\_\_\_ \$420  
 PAD SIZE:     5 x 8     5 x 10     5 x 14

### Other Matrix Mounting Hardware

*Must provide part number*

\_\_\_\_\_  
*For mounting hardware other than quick release*



**Vest**

- Full Shape, Comfort Strap (E0960) \_\_\_\_\_ \$167
  - Small 30589     Medium 30590     Large 30591
  - X-Large 30592     2X-Large 30593
- Full Shape, Extended Strap (E0960) \_\_\_\_\_ \$167
  - Small 30594     Medium 30595     Large 30596
  - X-Large 30597     2X-Large 30598
- Trim Shape, Comfort Strap (E0960) \_\_\_\_\_ \$169
  - Small 30618     Medium 30619     Large 30620
  - X-Large 30621     2X-Large 30622
- Trim Shape, Extended Strap (E0960) \_\_\_\_\_ \$169
  - Small 30604     Medium 30624     Large 30625
  - X-Large 30626     2X-Large 30627

**Shoulder Harness**

*With metal clips and cam locks*

- Small (E0960) 61403 \_\_\_\_\_ \$157
- Medium (E0960) 61503 \_\_\_\_\_ \$164
- Large (E0960) 61703 \_\_\_\_\_ \$170
- X-Large (E0960) 61803 \_\_\_\_\_ \$177

**Hip Belt**

- 4 Point Y-Style, Metal Buckle (E0978) \_\_\_\_\_ \$173
  - Small 60024     Medium 60025     Large 60026
- Dual Rear Pull, Metal Attaching Clips (E0978)
  - 1" Straps, 6.25"W x 1.75"H Pads 60161 \_\_\_\_\_ \$97
  - 1.5" Straps, 7.25"W x 2.25"H Pads 60171 \_\_\_\_\_ \$106
  - 1.5" Straps, 9.25"W x 2.5"H Pads 60181 \_\_\_\_\_ \$106
- Single Pull, Metal Buckle (E0978)
  - 1" Straps, 4.25"W x 1.75"H Pads 61001 \_\_\_\_\_ \$86
  - 1" Straps, 6.25"W x 1.75"H Pads 61201 \_\_\_\_\_ \$86
  - 1.5" Straps, 5.25"W x 2.25"H Pads 62001 \_\_\_\_\_ \$95
  - 1.5" Straps, 9.25"W x 2.5"H Pads 62201 \_\_\_\_\_ \$95
  - 2" Straps, 9.25"W x 3"H Pads 64001 \_\_\_\_\_ \$102
  - 2" Straps, 11.25"W x 3"H Pads 64101 \_\_\_\_\_ \$102

**Cane Mounted Swing Away Lateral Thoracic Support**

*Sold as pair*

- Swing Away Kits w/ 4"x6" Gel Pads (E0956, E1028) 50021 \_\_\_\_\_ \$478
- Swing Away Kits, No Pads (E0956, E1028) 50022 \_\_\_\_\_ \$309

**Foot / Shoe Holder**

*Compatible only with Aluminum Angle Adjustable and Composite Angle Adjustable Footplates*

- FootSure Ankle Support (K0108) \_\_\_\_\_ \$178  
*Sold in pairs*
  - X-Small w/ Buckle Closure, 6"-8" Circum. 60250
  - Small w/ Buckle Closure, 7"-8.5" Circum. 60251
  - Medium w/ Buckle Closure, 8"-11" Circum. 60252
  - Large w/ Buckle Closure, 9"-12" Circum. 60253
  - X-Large w/ Buckle Closure, 12.5"-15" Circum. 60254
  - X-Small w/ Hook-N-Loop Closure, 5"-6" Circum. 60255
  - Small w/ Hook-N-Loop Closure, 6"-8" Circum. 60256
  - Med w/ Hook-N-Loop Closure, 8"-10.5" Circum. 60257
  - Large w/ Hook-N-Loop Closure, 9"-12" Circum. 60258
  - X-Lrg w/ Hook-N-Loop Closure, 12"-15" Circum. 60259
- Padded Footbox (K0108)
  - Small, Full Design, Fits 14" Chair 30111 \_\_\_\_\_ \$287
  - Medium, Full Design, Fits 16" Chair 30112 \_\_\_\_\_ \$293
  - Large, Full Design, Fits 18" Chair 30113 \_\_\_\_\_ \$298
  - Small, Split Design, Fits 14" Chair 30152 \_\_\_\_\_ \$367
  - Medium, Split Design, Fits 16" Chair 30153 \_\_\_\_\_ \$373
  - Large, Split Design, Fits 18" Chair 30154 \_\_\_\_\_ \$378
- Molded Shoe Holders (K0108) \_\_\_\_\_ \$162  
*Sold in pairs*
  - Small with Plain Straps, 3.5" x 5.25" 30293
  - Medium with Plain Straps, 4" x 6.5" 30303
  - Large with Plain Straps, 4.75" x 8" 30309
  - X-Large with Plain Straps, 5" x 10" 30316





**Half Tray**

*Not compatible with Tubular Flip Up Armrest*

- Standard Shape, Clear, Flip-Away (E0956, E1028) \_\_\_\_\_ \$246
  - Right Side 41310       Left Side 41311
- Standard Shape, Padded, Flip-Away (E0956, E1028) \_\_\_\_\_ \$269
  - Right Side 41350       Left Side 41353
- Trimline Shape, Padded, Flip-Away (E0956, E1028) \_\_\_\_\_ \$269
  - Right Side 41352       Left Side 41351

**Clear Tray with Border and Pre-Drilled Holes**

**STEP 1: Choose Tray**

Chair Width	Tray Dimension (W x D)	Body Cutout (W x D)	HCPCS	Part Number	Price	
10"	17 x 15	7.5 x 4.5, H2.75	E0955	33203	\$209	<input type="checkbox"/>
11"	17 x 15	7.5 x 4.5	E0950	33204	\$209	<input type="checkbox"/>
	17 x 15	8 x 5	E0950	33205	\$220	<input type="checkbox"/>
12"	18 x 16	8.5 x 5	E0950	33206	\$220	<input type="checkbox"/>
	18 x 16	9.5 x 5.5	E0950	33207	\$220	<input type="checkbox"/>
13"	19 x 18	9.5 x 5.5	E0950	33208	\$220	<input type="checkbox"/>
	19 x 18	10.5 x 6	E0950	33209	\$220	<input type="checkbox"/>
14"	20 x 18	10.5 x 6	E0950	33210	\$220	<input type="checkbox"/>
	20 x 18	11.5 x 6.5	E0950	33211	\$220	<input type="checkbox"/>
15"	21 x 18	11.5 x 6.5	E0950	33212	\$220	<input type="checkbox"/>
	21 x 18	12.5 x 7	E0950	33213	\$220	<input type="checkbox"/>
16"	22 x 21	12.5 x 7.5	E0950	33214	\$255	<input type="checkbox"/>
	22 x 21	13.5 x 8.5	E0950	33215	\$255	<input type="checkbox"/>

**STEP 2: Choose Tray Hardware**

- E-Z Lock Clamps, Tube Arm (E1028) 31072 \_\_\_\_\_ \$150  
*Available only with Tubular Flip Up and Angle Adjustable Locking Extendable Flip Up with Foam Grip Armrests. Sold in Pairs.*
- E-Z Lock Clamps, Standard (E1028) 31098 \_\_\_\_\_ \$142  
*Not available with Tubular Flip Up and Angle Adjustable Locking Extendable Flip Up with Foam Grip Armrests. Sold in Pairs.*

**Communication Mount Hardware**

*Choose Device Holder AND Mount*

**STEP 1: What DEVICE will be used?**

- Smart Phone 32359 \_\_\_\_\_ \$150
- Large Smart Phone 32354 \_\_\_\_\_ \$150
- I-Pad 32355 \_\_\_\_\_ \$184
- I-Pad Mini 32356 \_\_\_\_\_ \$178
- 7" Tablet 32357 \_\_\_\_\_ \$173
- 10" Tablet 32358 \_\_\_\_\_ \$395

**STEP 2: Where does it need to MOUNT?**

- Armrest 32351 \_\_\_\_\_ \$395  
*Not compatible with Tubular Flip Up Armrest*
- Tray 32352 \_\_\_\_\_ \$212