



Phone (800) 981-1540 | Fax (715) 254-0996
 sales@kimobility.com | www.kimobility.com

Effective 01/01/2018

Company: _____

Ship To: _____

Account Number: _____ Date: _____

Address: _____

Quote PO Number: _____

City: _____ State: _____

RTS: _____

Zip: _____ Phone: _____

Marked For: _____

Email/Fax: _____

Mygo Stander

HCPC: E0641

Size 1 Weight Capacity: 110 lbs | Size 2 Weight Capacity: 132 lbs

Mygo Stander

Includes manual angle adjustment, polyurethane chest pad, hip pad and knee supports, two piece footplate and chassis

	Description	Item Code	Retail
<input type="checkbox"/>	Size 1 Mygo Stander	143-1600	\$4025
<input type="checkbox"/>	Size 2 Mygo Stander	143-2600	\$4025

Cover

Must match Stander size. Includes chest and hip pad covers and knee support covers.

	Description	Item Code	Retail
<input type="checkbox"/>	Size 1 Mygo Stander Covers - Blue	143-1770-08	N/C
<input type="checkbox"/>	Size 1 Mygo Stander Covers - Orange	143-1770-07	N/C
<input type="checkbox"/>	Size 1 Mygo Stander Covers - Pink	143-1770-09	N/C
<input type="checkbox"/>	Size 2 Mygo Stander Covers - Blue	143-2770-08	N/C
<input type="checkbox"/>	Size 2 Mygo Stander Covers - Orange	143-2770-07	N/C
<input type="checkbox"/>	Size 2 Mygo Stander Covers - Pink	143-2770-09	N/C

Shoulder Support

Cover selection required with Shoulder Pad

	Description	Item Code	Retail
<input type="checkbox"/>	Mygo Stander Shoulder Pad	143-1638	\$322
<input type="checkbox"/>	Mygo Stander Shoulder Pad Cover - Blue	143-1768-08	N/C
<input type="checkbox"/>	Mygo Stander Shoulder Pad Cover - Orange	143-1768-07	N/C
<input type="checkbox"/>	Mygo Stander Shoulder Pad Cover - Pink	143-1768-09	N/C
<input type="checkbox"/>	No Cover		N/C

HCPCS codes provided are based on code verification by the PDAC or our interpretation of the code definitions.

Proper use of HCPCS codes does not ensure coverage or payment. For coverage information, verify the appropriate payer's coverage policy.

DCN0286.1 Pricing and specifications subject to change without notice

Page 1 of 6

LECKEY^o

Head Support

Cushion selection required with hardware

	Description	Item Code	Retail
<input type="checkbox"/>	Contoured Headrest Support Hardware	120-859	\$317
<input type="checkbox"/>	Contoured Headrest Cushion Infection Control - Black	143-1787	\$149
<input type="checkbox"/>	Contoured Headrest Cushion - Black	137-800-05	\$145
<input type="checkbox"/>	Flat Headrest Support Hardware	117-697	\$194
<input type="checkbox"/>	Flat Headrest Cushion - Black	117-850-05	\$98
<input type="checkbox"/>	Flat Headrest Laterals with Cover - Black (Pair) <i>Requires Flat Headrest. Not available with Contoured Headrest.</i>	120-893-05	\$278

Chest Support

Chest Wraparound Harness requires Laterals

	Description	Item Code	Retail
<input type="checkbox"/>	Size 1 Fixed Chest Laterals (Pair)	143-1619	\$139
<input type="checkbox"/>	Size 2 Fixed Chest Laterals (Pair)	143-2619	\$139
<input type="checkbox"/>	Size 1 and 2 Flip-Away Chest Laterals (Pair)	143-1675	\$448
<input type="checkbox"/>	Size 1 Chest Wraparound Harness - Blue	143-1780-08	\$231
<input type="checkbox"/>	Size 1 Chest Wraparound Harness - Orange	143-1780-07	\$231
<input type="checkbox"/>	Size 1 Chest Wraparound Harness - Pink	143-1780-09	\$231
<input type="checkbox"/>	Size 1 Chest Wraparound Harness Infection Control - Black	143-1775-14	\$231
<input type="checkbox"/>	Size 2 Chest Wraparound Harness - Blue	143-2780-08	\$238
<input type="checkbox"/>	Size 2 Chest Wraparound Harness - Orange	143-2780-07	\$238
<input type="checkbox"/>	Size 2 Chest Wraparound Harness - Pink	143-2780-09	\$238
<input type="checkbox"/>	Size 2 Chest Wraparound Harness Infection Control - Black	143-2775-14	\$238
<input type="checkbox"/>	Size 3 Chest Wraparound Harness - Blue	143-3780-08	\$244
<input type="checkbox"/>	Size 3 Chest Wraparound Harness - Orange	143-3780-07	\$244
<input type="checkbox"/>	Size 3 Chest Wraparound Harness - Pink	143-3780-09	\$244
<input type="checkbox"/>	Size 3 Chest Wraparound Harness Infection Control - Black	143-3775-14	\$244
<input type="checkbox"/>	Size 1 and 2 Chest Complex Lateral Vertical with Cover Right <i>Not available with Wraparound Harness</i>	143-1674	\$262
<input type="checkbox"/>	Size 1 and 2 Chest Complex Lateral Vertical with Cover Left <i>Not available with Wraparound Harness</i>	143-1673	\$262

LECKEY^o

Hip Laterals

Hip Wraparound Harness requires Hip Laterals

	Description	Item Code	Retail
<input type="checkbox"/>	Hip Laterals (Pair) <i>Not available with Derotation Belt</i>	143-1663	\$139
<input type="checkbox"/>	Size 1 Hip Wraparound Harness - Blue	143-1780-08	\$231
<input type="checkbox"/>	Size 1 Hip Wraparound Harness - Orange	143-1780-07	\$231
<input type="checkbox"/>	Size 1 Hip Wraparound Harness - Pink	143-1780-09	\$231
<input type="checkbox"/>	Size 1 Hip Wraparound Harness Infection Control - Black	143-1775-14	\$231
<input type="checkbox"/>	Size 2 Hip Wraparound Harness - Blue	143-2780-08	\$238
<input type="checkbox"/>	Size 2 Hip Wraparound Harness - Orange	143-2780-07	\$238
<input type="checkbox"/>	Size 2 Hip Wraparound Harness - Pink	143-2780-09	\$238
<input type="checkbox"/>	Size 2 Hip Wraparound Harness Infection Control - Black	143-2775-14	\$238
<input type="checkbox"/>	Size 3 Hip Wraparound Harness - Blue	143-3780-08	\$244
<input type="checkbox"/>	Size 3 Hip Wraparound Harness - Orange	143-3780-07	\$244
<input type="checkbox"/>	Size 3 Hip Wraparound Harness - Pink	143-3780-09	\$244
<input type="checkbox"/>	Size 3 Hip Wraparound Harness Infection Control - Black	143-3775-14	\$244

Hip Belt

Not available with Hip Laterals or Hip Wraparound Harness

	Description	Item Code	Retail
<input type="checkbox"/>	Size 1 Prone Derotation Hip Belt with Hardware - Blue	143-1798-08	\$262
<input type="checkbox"/>	Size 1 Prone Derotation Hip Belt with Hardware - Orange	143-1798-07	\$262
<input type="checkbox"/>	Size 1 Prone Derotation Hip Belt with Hardware - Pink	143-1798-09	\$262
<input type="checkbox"/>	Size 2 Prone Derotation Hip Belt with Hardware - Blue	143-2798-08	\$296
<input type="checkbox"/>	Size 2 Prone Derotation Hip Belt with Hardware - Orange	143-2798-07	\$296
<input type="checkbox"/>	Size 2 Prone Derotation Hip Belt with Hardware - Pink	143-2798-09	\$296
<input type="checkbox"/>	Size 3 Prone Derotation Hip Belt with Hardware - Blue	143-3798-08	\$328
<input type="checkbox"/>	Size 3 Prone Derotation Hip Belt with Hardware - Orange	143-3798-07	\$328
<input type="checkbox"/>	Size 3 Prone Derotation Hip Belt with Hardware - Pink	143-3798-09	\$328
<input type="checkbox"/>	Size 1 Supine Derotation Hip Belt with Hardware - Blue	143-1783-08	\$262
<input type="checkbox"/>	Size 1 Supine Derotation Hip Belt with Hardware - Orange	143-1783-07	\$262
<input type="checkbox"/>	Size 1 Supine Derotation Hip Belt with Hardware - Pink	143-1783-09	\$262
<input type="checkbox"/>	Size 2 Supine Derotation Hip Belt with Hardware - Blue	143-2783-08	\$296
<input type="checkbox"/>	Size 2 Supine Derotation Hip Belt with Hardware - Orange	143-2783-07	\$296
<input type="checkbox"/>	Size 2 Supine Derotation Hip Belt with Hardware - Pink	143-2783-09	\$296
<input type="checkbox"/>	Size 3 Supine Derotation Hip Belt with Hardware - Blue	143-3783-08	\$328
<input type="checkbox"/>	Size 3 Supine Derotation Hip Belt with Hardware - Orange	143-3783-07	\$328
<input type="checkbox"/>	Size 3 Supine Derotation Hip Belt with Hardware - Pink	143-3783-09	\$328

LECKEY^o

Hip Infill

	Description	Item Code	Retail
<input type="checkbox"/>	Hip Infill Pad	143-1666	\$287
<input type="checkbox"/>	Size 1 Hip Infill Pad Cover - Blue	143-1799-08	\$120
<input type="checkbox"/>	Size 1 Hip Infill Pad Cover - Orange	143-1799-07	\$120
<input type="checkbox"/>	Size 1 Hip Infill Pad Cover - Pink	143-1799-09	\$120
<input type="checkbox"/>	Size 2 Hip Infill Pad Cover - Blue	143-2799-08	\$138
<input type="checkbox"/>	Size 2 Hip Infill Pad Cover - Orange	143-2799-09	\$138
<input type="checkbox"/>	Size 2 Hip Infill Pad Cover - Pink	143-2799-07	\$138

Posterior Support

Selection of Posterior Support and Pommel requires selection of both Posterior Support Cover AND Pommel Support Cover

	Description	Item Code	Retail
<input type="checkbox"/>	Size 1 Posterior Support and Pommel <i>Not available with Hip Laterals or Wraparound Harness</i>	143-1624	\$681
<input type="checkbox"/>	Posterior Support Cover - Blue	143-1790-08	\$127
<input type="checkbox"/>	Posterior Support Cover - Orange	143-1790-07	\$127
<input type="checkbox"/>	Posterior Support Cover - Pink	143-1790-09	\$127
<input type="checkbox"/>	Posterior Support Infection Control Cover - Black	143-1790-14	\$127
<input type="checkbox"/>	Pommel Support Cover - Blue	143-1792-08	\$106
<input type="checkbox"/>	Pommel Support Cover - Orange	143-1792-07	\$106
<input type="checkbox"/>	Pommel Support Cover - Pink	143-1792-09	\$106
<input type="checkbox"/>	Pommel Support Infection Control Cover - Black	143-1792-14	\$106
<input type="checkbox"/>	Size 2 Posterior Support and Pommel <i>Not available with Hip Laterals or Wraparound Harness</i>	143-2624	\$716
<input type="checkbox"/>	Posterior Support Cover - Blue	143-1790-08	\$127
<input type="checkbox"/>	Posterior Support Cover - Orange	143-1790-07	\$127
<input type="checkbox"/>	Posterior Support Cover - Pink	143-1790-09	\$127
<input type="checkbox"/>	Posterior Support Infection Control Cover - Black	143-1790-14	\$127
<input type="checkbox"/>	Pommel Support Cover - Blue	143-1792-08	\$106
<input type="checkbox"/>	Pommel Support Cover - Orange	143-1792-07	\$106
<input type="checkbox"/>	Pommel Support Cover - Pink	143-1792-09	\$106
<input type="checkbox"/>	Pommel Support Infection Control Cover - Black	143-1792-14	\$106

LECKEY^o

Knee

	Description	Item Code	Retail
<input type="checkbox"/>	Knee Cup Height Extension Brackets (Pair) <i>Increases distance from footplate to top of knee by 4"</i>	143-1690	\$202
<input type="checkbox"/>	Size 1 Basic Knee Straps - Blue (Pair)	143-1776-08	\$157
<input type="checkbox"/>	Size 1 Basic Knee Straps - Orange (Pair)	143-1776-07	\$157
<input type="checkbox"/>	Size 1 Basic Knee Straps - Pink (Pair)	143-1776-09	\$157
<input type="checkbox"/>	Size 1 Basic Knee Straps Infection Control - Black (Pair)	143-1776-14	\$157
<input type="checkbox"/>	Size 2 Basic Knee Straps - Blue (Pair)	143-2776-08	\$192
<input type="checkbox"/>	Size 2 Basic Knee Straps - Orange (Pair)	143-2776-07	\$192
<input type="checkbox"/>	Size 2 Basic Knee Straps - Pink (Pair)	143-2776-09	\$192
<input type="checkbox"/>	Size 2 Basic Knee Straps Infection Control - Black (Pair)	143-2776-14	\$192
<input type="checkbox"/>	Split Knee Straps - Blue (Pair)	143-1784-08	\$296
<input type="checkbox"/>	Split Knee Straps - Orange (Pair)	143-1784-07	\$296
<input type="checkbox"/>	Split Knee Straps - Pink (Pair)	143-1784-09	\$296

Accessories

	Description	Item Code	Retail
<input type="checkbox"/>	Size 1 Mygo Stander Tray	143-1609	\$380
<input type="checkbox"/>	Size 2 Mygo Stander Tray	143-2609	\$380
<input type="checkbox"/>	Size 1 Mygo Stander Tray with Bowl Insert	143-1633	\$392
<input type="checkbox"/>	Size 2 Mygo Stander Tray with Bowl Insert	143-2633	\$392
<input type="checkbox"/>	Size 1 Mygo Stander Tray with Clear Insert	143-1610	\$431
<input type="checkbox"/>	Size 2 Mygo Stander Tray with Clear Insert	143-2610	\$431
<input type="checkbox"/>	Grab Rail / Toy Bar for Tray	143-1649	\$244
<input type="checkbox"/>	Size 1 Sandals (Pair)	152-1600	\$225
<input type="checkbox"/>	Size 2 Sandals (Pair)	152-2600	\$245
<input type="checkbox"/>	Size 3 Sandals (Pair)	152-3600	\$265
<input type="checkbox"/>	Sandal Raiser (Each)	143-1615	\$110

Assembly

Product Assembly Fee not subject to discounts

<input type="checkbox"/>	Product Assembly Fee	\$100
--------------------------	----------------------	-------

LECKEY^o

Technical Specifications

Specification	Size 1	Size 2
Age (approximate)	4 - 10 years	8 - 14 years
User Height (Min / Max)	37" / 57.1"	53.1" / 66.9"
Maximum User Weight	110 lbs	132 lbs
Chest Support Height	29.9" - 42.5"	38.6" - 51.2"
Chest Width	6.3" - 11"	7.9" - 12.6"
Hip Support Height	18.1" - 28.7"	24.4" - 34.4"
Knee Cup Height	11" - 15"	12.25" - 16.25"
Knee Cup Height with Extension Bracket	15" - 19"	16.25" - 20.25"
Hip Width	8.7" - 13"	11" - 15.4"
Footplate Angle: Plantarflexion / Dorsiflexion	10° / 10°	
Maximum Knee Contracture	25°	
Tray Height Range	0" - 6"	
Tray Angle Adjustment: Prone / Supine	45° / 45°	
Product Weight	48 lbs	55 lbs
Product Angle Range: Prone / Supine	0° - 80°	
Stander Footprint (length x width)	38.5" x 24.8"	44.5" x 25.6"

Technical Specifications - Metric

Specification	Size 1	Size 2
Age (approximate)	4 - 10 years	8 - 14 years
User Height (Min / Max)	1080 mm / 1450 mm	1350 mm / 1700 mm
Maximum User Weight	50 kg	60 kg
Chest Support Height	760 mm - 1080 mm	980 mm - 1300 mm
Chest Width	160 mm - 280 mm	200 mm - 320 mm
Hip Support Height	460 mm - 730 mm	620 mm - 875 mm
Knee Cup Height	280 mm - 380 mm	310mm - 410 mm
Knee Cup Height with Extension Bracket	381 mm - 481 mm	411 mm - 511 mm
Hip Width	220 mm - 330 mm	280 mm - 390 mm
Footplate Angle: Plantarflexion / Dorsiflexion	10° / 10°	
Maximum Knee Contracture	25°	
Tray Height Range	0 mm - 150 mm	
Tray Angle Adjustment: Prone / Supine	45° / 45°	
Product Weight	22 kg	25 kg
Product Angle Range: Prone / Supine	0° - 80°	
Stander Footprint (length x width)	978 mm - 630 mm	1130 mm x 650 mm