



Phone (800) 981-1540 | Fax (715) 254-0996  
 sales@kimobility.com | www.kimobility.com

Company: \_\_\_\_\_ Ship To: \_\_\_\_\_  
 Account Number: \_\_\_\_\_ Date: \_\_\_\_\_ Address: \_\_\_\_\_  
 Quote  PO Number: \_\_\_\_\_ City: \_\_\_\_\_ State: \_\_\_\_\_  
 RTS: \_\_\_\_\_ Zip: \_\_\_\_\_ Phone: \_\_\_\_\_  
 Marked For: \_\_\_\_\_ Email/Fax: \_\_\_\_\_

Effective 03/19/2018

## Mygo Stander

**HCPC: E0641**

*Size 1 Weight Capacity: 110 lbs | Size 2 Weight Capacity: 132 lbs*

### New! Packages

*Package includes Mygo Stander, Stander Covers, Shoulder Pad and Cover, Contoured Headrest Hardware and Cushion - Black, Fixed Chest Laterals, Chest Wraparound Harness, Supine Derotation Hip Belt, Knee Cup Height Extension Brackets, Basic Knee Straps, Tray, Sandals with Straps and Product Assembly*

	Description	Item Code	HCPC	Retail
<input type="checkbox"/>	Size 1 Mygo Stander Package - Blue <i>Includes Medium Sandals</i>	PMSTB1	E0641	\$6545
<input type="checkbox"/>	Size 1 Mygo Stander Package - Pink <i>Includes Medium Sandals</i>	PMSTP1	E0641	\$6545
<input type="checkbox"/>	Size 2 Mygo Stander Package - Blue <i>Includes Large Sandals</i>	PMSTB2	E0641	\$6641
<input type="checkbox"/>	Size 2 Mygo Stander Package - Pink <i>Includes Large Sandals</i>	PMSTP2	E0641	\$6641

### Mygo Stander

*Selection required. Includes manual angle adjustment, polyurethane chest pad, hip pad and knee supports, two piece footplate and chassis*

	Description	Item Code	HCPC	Retail
<input type="checkbox"/>	Size 1 Mygo Stander	143-1600		\$4025
<input type="checkbox"/>	Size 2 Mygo Stander	143-2600		\$4025

### Cover

*Selection required. Must match Stander size. Includes chest and hip pad covers and knee support covers.*

	Description	Item Code	HCPC	Retail
<input type="checkbox"/>	Size 1 Mygo Stander Covers - Blue	143-1770-08		N/C
<input type="checkbox"/>	Size 1 Mygo Stander Covers - Orange	143-1770-07		N/C
<input type="checkbox"/>	Size 1 Mygo Stander Covers - Pink	143-1770-09		N/C
<input type="checkbox"/>	Size 2 Mygo Stander Covers - Blue	143-2770-08		N/C
<input type="checkbox"/>	Size 2 Mygo Stander Covers - Orange	143-2770-07		N/C
<input type="checkbox"/>	Size 2 Mygo Stander Covers - Pink	143-2770-09		N/C

*HCPCS codes provided are based on code verification by the PDAC or our interpretation of the code definitions.*

*Proper use of HCPCS codes does not ensure coverage or payment. For coverage information, verify the appropriate payer's coverage policy.*

## LECKEY<sup>o</sup>

### Shoulder Support

*Cover selection required with Shoulder Pad*

	Description	Item Code	HCPC	Retail
<input type="checkbox"/>	Mygo Stander Shoulder Pad	143-1638		\$322
<input type="checkbox"/>	Mygo Stander Shoulder Pad Cover - Blue	143-1768-08		N/C
<input type="checkbox"/>	Mygo Stander Shoulder Pad Cover - Orange	143-1768-07		N/C
<input type="checkbox"/>	Mygo Stander Shoulder Pad Cover - Pink	143-1768-09		N/C
<input type="checkbox"/>	No Cover			N/C

### Head Support

*Cushion selection required with hardware*

	Description	Item Code	HCPC	Retail
<input type="checkbox"/>	Contoured Headrest Support Hardware	120-859		\$317
<input type="checkbox"/>	Contoured Headrest Cushion Infection Control - Black	143-1787		\$149
<input type="checkbox"/>	Contoured Headrest Cushion - Black	137-800-05		\$145
<input type="checkbox"/>	Flat Headrest Support Hardware	117-697		\$194
<input type="checkbox"/>	Flat Headrest Cushion - Black	117-850-05		\$98
<input type="checkbox"/>	Flat Headrest Laterals with Covers - Black (Pair) <i>Requires Flat Headrest. Not available with Contoured Headrest.</i>	120-893-05		\$278

### Wraparound Harness

*Harnesses may be used as chest harness or hip harness. If used as chest harness, chest laterals are required. If used as hip harness, hip laterals are required.*

Description	Quantity 1	Quantity 2	Item Code	Retail (Each)
Size 1 Chest/Hip Wraparound Harness - Blue	<input type="checkbox"/>	<input type="checkbox"/>	143-1780-08	\$231
Size 1 Chest/Hip Wraparound Harness - Orange	<input type="checkbox"/>	<input type="checkbox"/>	143-1780-07	\$231
Size 1 Chest/Hip Wraparound Harness - Pink	<input type="checkbox"/>	<input type="checkbox"/>	143-1780-09	\$231
Size 1 Chest/Hip Wraparound Harness Infection Control - Black	<input type="checkbox"/>	<input type="checkbox"/>	143-1775-14	\$231
Size 2 Chest/Hip Wraparound Harness - Blue	<input type="checkbox"/>	<input type="checkbox"/>	143-2780-08	\$238
Size 2 Chest/Hip Wraparound Harness - Orange	<input type="checkbox"/>	<input type="checkbox"/>	143-2780-07	\$238
Size 2 Chest/Hip Wraparound Harness - Pink	<input type="checkbox"/>	<input type="checkbox"/>	143-2780-09	\$238
Size 2 Chest/Hip Wraparound Harness Infection Control - Black	<input type="checkbox"/>	<input type="checkbox"/>	143-2775-14	\$238
Size 3 Chest/Hip Wraparound Harness - Blue	<input type="checkbox"/>	<input type="checkbox"/>	143-3780-08	\$244
Size 3 Chest/Hip Wraparound Harness - Orange	<input type="checkbox"/>	<input type="checkbox"/>	143-3780-07	\$244
Size 3 Chest/Hip Wraparound Harness - Pink	<input type="checkbox"/>	<input type="checkbox"/>	143-3780-09	\$244
Size 3 Chest/Hip Wraparound Harness Infection Control - Black	<input type="checkbox"/>	<input type="checkbox"/>	143-3775-14	\$244

## LECKEY<sup>o</sup>

### Lateral Support

**Chest Laterals required if ordering Wraparound Harness to be used as chest harness. Hip Laterals required if ordering Wraparound Harness to be used as hip harness.**

	Description	Item Code	HCPC	Retail
<input type="checkbox"/>	Size 1 Fixed Chest Laterals (Pair)	143-1619		\$139
<input type="checkbox"/>	Size 2 Fixed Chest Laterals (Pair)	143-2619		\$139
<input type="checkbox"/>	Size 1 and 2 Flip-Away Chest Laterals (Pair)	143-1675		\$448
<input type="checkbox"/>	Size 1 and 2 Chest Complex Lateral Vertical with Cover Right <i>Not available if ordering Wraparound Harness to be used as chest harness</i>	143-1674		\$262
<input type="checkbox"/>	Size 1 and 2 Chest Complex Lateral Vertical with Cover Left <i>Not available if ordering Wraparound Harness to be used as chest harness</i>	143-1673		\$262
<input type="checkbox"/>	Hip Laterals (Pair) <i>Not available with Derotation Belt</i>	143-1663		\$139

### Hip Belt

**Not available with Hip Laterals or Hip Wraparound Harness**

	Description	Item Code	HCPC	Retail
<input type="checkbox"/>	Size 1 Prone Derotation Hip Belt with Hardware - Blue	143-1798-08		\$262
<input type="checkbox"/>	Size 1 Prone Derotation Hip Belt with Hardware - Orange	143-1798-07		\$262
<input type="checkbox"/>	Size 1 Prone Derotation Hip Belt with Hardware - Pink	143-1798-09		\$262
<input type="checkbox"/>	Size 2 Prone Derotation Hip Belt with Hardware - Blue	143-2798-08		\$296
<input type="checkbox"/>	Size 2 Prone Derotation Hip Belt with Hardware - Orange	143-2798-07		\$296
<input type="checkbox"/>	Size 2 Prone Derotation Hip Belt with Hardware - Pink	143-2798-09		\$296
<input type="checkbox"/>	Size 3 Prone Derotation Hip Belt with Hardware - Blue	143-3798-08		\$328
<input type="checkbox"/>	Size 3 Prone Derotation Hip Belt with Hardware - Orange	143-3798-07		\$328
<input type="checkbox"/>	Size 3 Prone Derotation Hip Belt with Hardware - Pink	143-3798-09		\$328
<input type="checkbox"/>	Size 1 Supine Derotation Hip Belt with Hardware - Blue	143-1783-08		\$262
<input type="checkbox"/>	Size 1 Supine Derotation Hip Belt with Hardware - Orange	143-1783-07		\$262
<input type="checkbox"/>	Size 1 Supine Derotation Hip Belt with Hardware - Pink	143-1783-09		\$262
<input type="checkbox"/>	Size 2 Supine Derotation Hip Belt with Hardware - Blue	143-2783-08		\$296
<input type="checkbox"/>	Size 2 Supine Derotation Hip Belt with Hardware - Orange	143-2783-07		\$296
<input type="checkbox"/>	Size 2 Supine Derotation Hip Belt with Hardware - Pink	143-2783-09		\$296
<input type="checkbox"/>	Size 3 Supine Derotation Hip Belt with Hardware - Blue	143-3783-08		\$328
<input type="checkbox"/>	Size 3 Supine Derotation Hip Belt with Hardware - Orange	143-3783-07		\$328
<input type="checkbox"/>	Size 3 Supine Derotation Hip Belt with Hardware - Pink	143-3783-09		\$328

# LECKEY<sup>o</sup>

## Hip Infill

	Description	Item Code	HCPC	Retail
<input type="checkbox"/>	Hip Infill Pad	143-1666		\$287
<input type="checkbox"/>	Size 1 Hip Infill Pad Cover - Blue	143-1799-08		\$120
<input type="checkbox"/>	Size 1 Hip Infill Pad Cover - Orange	143-1799-07		\$120
<input type="checkbox"/>	Size 1 Hip Infill Pad Cover - Pink	143-1799-09		\$120
<input type="checkbox"/>	Size 2 Hip Infill Pad Cover - Blue	143-2799-08		\$138
<input type="checkbox"/>	Size 2 Hip Infill Pad Cover - Orange	143-2799-09		\$138
<input type="checkbox"/>	Size 2 Hip Infill Pad Cover - Pink	143-2799-07		\$138

## Posterior Support

*Selection of Posterior Support and Pommel requires selection of both Posterior Support Cover AND Pommel Support Cover. Not available with Hip Laterals or Hip Wraparound Harness.*

	Description	Item Code	HCPC	Retail
<input type="checkbox"/>	Size 1 Posterior Support and Pommel	143-1624		\$681
<input type="checkbox"/>	Posterior Support Cover - Blue	143-1790-08		\$127
<input type="checkbox"/>	Posterior Support Cover - Orange	143-1790-07		\$127
<input type="checkbox"/>	Posterior Support Cover - Pink	143-1790-09		\$127
<input type="checkbox"/>	Posterior Support Infection Control Cover - Black	143-1790-14		\$127
<input type="checkbox"/>	Pommel Support Cover - Blue	143-1792-08		\$106
<input type="checkbox"/>	Pommel Support Cover - Orange	143-1792-07		\$106
<input type="checkbox"/>	Pommel Support Cover - Pink	143-1792-09		\$106
<input type="checkbox"/>	Pommel Support Infection Control Cover - Black	143-1792-14		\$106
<input type="checkbox"/>	Size 2 Posterior Support and Pommel	143-2624		\$716
<input type="checkbox"/>	Posterior Support Cover - Blue	143-1790-08		\$127
<input type="checkbox"/>	Posterior Support Cover - Orange	143-1790-07		\$127
<input type="checkbox"/>	Posterior Support Cover - Pink	143-1790-09		\$127
<input type="checkbox"/>	Posterior Support Infection Control Cover - Black	143-1790-14		\$127
<input type="checkbox"/>	Pommel Support Cover - Blue	143-1792-08		\$106
<input type="checkbox"/>	Pommel Support Cover - Orange	143-1792-07		\$106
<input type="checkbox"/>	Pommel Support Cover - Pink	143-1792-09		\$106
<input type="checkbox"/>	Pommel Support Infection Control Cover - Black	143-1792-14		\$106

## LECKEY<sup>o</sup>

### Knee

	Description	Item Code	HCPC	Retail
<input type="checkbox"/>	Knee Cup Height Extension Brackets (Pair) <i>Increases distance from footplate to top of knee by 4"</i>	143-1690		\$202
<input type="checkbox"/>	Size 1 Basic Knee Straps - Blue (Pair)	143-1776-08		\$157
<input type="checkbox"/>	Size 1 Basic Knee Straps - Orange (Pair)	143-1776-07		\$157
<input type="checkbox"/>	Size 1 Basic Knee Straps - Pink (Pair)	143-1776-09		\$157
<input type="checkbox"/>	Size 1 Basic Knee Straps Infection Control - Black (Pair)	143-1776-14		\$157
<input type="checkbox"/>	Size 2 Basic Knee Straps - Blue (Pair)	143-2776-08		\$192
<input type="checkbox"/>	Size 2 Basic Knee Straps - Orange (Pair)	143-2776-07		\$192
<input type="checkbox"/>	Size 2 Basic Knee Straps - Pink (Pair)	143-2776-09		\$192
<input type="checkbox"/>	Size 2 Basic Knee Straps Infection Control - Black (Pair)	143-2776-14		\$192
<input type="checkbox"/>	Split Knee Straps - Blue (Pair)	143-1784-08		\$296
<input type="checkbox"/>	Split Knee Straps - Orange (Pair)	143-1784-07		\$296
<input type="checkbox"/>	Split Knee Straps - Pink (Pair)	143-1784-09		\$296

### Accessories

	Description	Item Code	HCPC	Retail
<input type="checkbox"/>	Size 1 Mygo Stander Tray	143-1609		\$380
<input type="checkbox"/>	Size 2 Mygo Stander Tray	143-2609		\$380
<input type="checkbox"/>	Size 1 Mygo Stander Tray with Bowl Insert	143-1633		\$392
<input type="checkbox"/>	Size 2 Mygo Stander Tray with Bowl Insert	143-2633		\$392
<input type="checkbox"/>	Size 1 Mygo Stander Tray with Clear Insert	143-1610		\$431
<input type="checkbox"/>	Size 2 Mygo Stander Tray with Clear Insert	143-2610		\$431
<input type="checkbox"/>	Grab Rail/Toy Bar for Tray	143-1649		\$244
<input type="checkbox"/>	Small Sandals with Straps (Pair)	152-1600		\$225
<input type="checkbox"/>	Medium Sandals with Straps (Pair)	152-2600		\$245
<input type="checkbox"/>	Large Sandals with Straps (Pair)	152-3600		\$265
<input type="checkbox"/>	Sandal Raiser (Each)	<b>Quantity: _____ (1-6)</b>	143-1615	\$110

### Assembly

**Product Assembly Fee not subject to discounts**

<input type="checkbox"/>	Product Assembly Fee		\$120
--------------------------	----------------------	--	-------

# LECKEY<sup>o</sup>

## Technical Specifications

Specification	Size 1	Size 2
Age (approximate)	4 - 10 years	8 - 14 years
User Height (Min / Max)	37" / 57.1"	53.1" / 66.9"
Maximum User Weight	110 lbs	132 lbs
Chest Support Height	29.9" - 42.5"	38.6" - 51.2"
Chest Width	6.3" - 11"	7.9" - 12.6"
Hip Support Height	18.1" - 28.7"	24.4" - 34.4"
Knee Cup Height	11" - 15"	12.25" - 16.25"
Knee Cup Height with Extension Bracket	15" - 19"	16.25" - 20.25"
Hip Width	8.7" - 13"	11" - 15.4"
Footplate Angle: Plantarflexion / Dorsiflexion	10° / 10°	
Maximum Knee Contracture	25°	
Tray Height Range	0" - 6"	
Tray Angle Adjustment: Prone / Supine	45° / 45°	
Product Weight	48 lbs	55 lbs
Product Angle Range: Prone / Supine	0° - 80°	
Stander Footprint (length x width)	38.5" x 24.8"	44.5" x 25.6"

## Technical Specifications - Metric

Specification	Size 1	Size 2
Age (approximate)	4 - 10 years	8 - 14 years
User Height (Min / Max)	1080 mm / 1450 mm	1350 mm / 1700 mm
Maximum User Weight	50 kg	60 kg
Chest Support Height	760 mm - 1080 mm	980 mm - 1300 mm
Chest Width	160 mm - 280 mm	200 mm - 320 mm
Hip Support Height	460 mm - 730 mm	620 mm - 875 mm
Knee Cup Height	280 mm - 380 mm	310mm - 410 mm
Knee Cup Height with Extension Bracket	381 mm - 481 mm	411 mm - 511 mm
Hip Width	220 mm - 330 mm	280 mm - 390 mm
Footplate Angle: Plantarflexion / Dorsiflexion	10° / 10°	
Maximum Knee Contracture	25°	
Tray Height Range	0 mm - 150 mm	
Tray Angle Adjustment: Prone / Supine	45° / 45°	
Product Weight	22 kg	25 kg
Product Angle Range: Prone / Supine	0° - 80°	
Stander Footprint (length x width)	978 mm - 630 mm	1130 mm x 650 mm