



Phone (800) 981-1540 | Fax (715) 254-0996
sales@kimobility.com | www.kimobility.com



Company:

Ship To:

Account Number:

Date:

Address:

Quote: PO Number:

City:

State:

ATP:

Zip:

Phone:

Marked For:

Email/Fax:

Associated Quote/Order#:

Effective: 4/1/2024

Spark- Transit Included (E1238) \$2532

Spark - Transit Not Included (E1238) \$2204

165 lb weight capacity Not available with upholstery

Frame

Seat Width

10"	N/C
11"	N/C
12"	N/C
13"	N/C
14"	N/C
15"	N/C
16"	N/C

Seat Depth


12"	N/C
13"	N/C
14"	N/C
15"	N/C
16"	N/C
17"	N/C
18"	N/C

Transit Option

Transit System	\$328
No charge when selected with Spark - Transit Included. Not available with Removable Seat Pan.	
Omit Transit	N/C

Seat Pan

Removable Aluminum Seat Pan (E2231)	\$300
Pan depth 1/2" shorter than Frame Depth	
Omit	N/C



Frame Width

Spark is set up to allow 2" increase in frame width from initial frame width in 1" increments. Growth is achieved by expanding cross braces. 16" maximum frame width. 10" Frame Width expands only 1".



Frame Depth

Spark is set up to allow 4" increase in frame depth from initial frame depth in 1" increments. 19" maximum frame depth.

Wheelchair Drive Ready

Chair will ship from Ki Mobility with adaption kit installed for selected Alber system. Order system separately from Alber USA. Select Rear Frame Height (full profile selections only) in drive wheel size on Rear Frame Height chart (page 4). Requires Rear Anti-Tips. Not available with Spinerge wheel.

e-motion Ready	N/C
twion Ready	N/C



Phone (800) 981-1540 | Fax (715) 254-0996
sales@kimobility.com | www.kimobility.com



Backrest

Backrest

Back Posts are folding and angle adjustable in 4° increments.

Back Post	Price	Height Range	Back Height				
Stroller Handle Back Post (Depth Adjustable)	N/C	Short	20"				
		Tall	24"				
Height Adjustable Straight with Push Handle (Depth Adjustable)	\$77	Short ¹	13"	14"	15"	16"	17"
		Medium	16"	17"	18"	19"	20"
		Tall	20" ²	21"	22"	23"	24"

¹ Not available with Dynamic Back

² Not available with Angle Adjustable Locking Extendable Flip Up Armrest

Backrest Options

Removable Stroller Handle Extension \$80

Not available with Height Adjustable Straight with Push Handle Back Posts.

Dynamic Back (E2398) \$1395

Available only with T-Arms. Not available with 12" Frame Depths with Rear Frame Heights lower than 14.5".

Frame Type / Footrest

Frame Type

Swing Away N/C
Super Low - Swing Away \$51

Footrest Hanger Type

Hanger	60°	70°	80°	90°		Retail
Extension Tube Mount			1			N/C
Front Mount			1	1		N/C
Pediatric PRO ELR						\$318
Pro Elevating Leg Rests E0990 ²						\$318
Omit Hangers						(\$20)

¹ Not available with 7" or 8" Casters

² Not available with Seat Widths less than 14"

Minimum Front Seat Height with Front Mount Hangers						
Hanger Length	Short		Medium		Long	
Hanger Bend	70°	80/90°	70°	80/90°	70°	80°
Minimum Front Seat Height with Standard Selection	13"	13.5"	14"	15"	16"	17"
Minimum Front Seat Height with 3" Extension	14.5"	16"	16.5"	18"	18.5"	20"

PRO ELR Cover

Available only on PRO ELR.

Gel Pad Covers (Pair) K0108 \$436

Seat to Footrest Length

Select Extension and Seat to Footrest Length in row corresponding to desired Footrest. Footrest Length of at least 2.5" shorter than Front Seat Height recommended. Extension Tube Hanger measurements are based on 70° Hanger.

Footrest	Extension	Seat to Footrest Length											
Extension Tube Hangers (60°, 70°, 80°)	Short (11" - 15")	11"	11.5"	12"	12.5"	13"	13.5"	14"	14.5"	15"			
	Medium (13" - 17")	13"	13.5"	14"	14.5"	15"	15.5"	16"	16.5"	17"			
	Long (15.5" - 19.5")	15.5"	16"	16.5"	17"	17.5"	18"	18.5"	19"	19.5"			
Front Mount Hangers 1		Standard Selections									Requires 3" Extension		
	Short Front Mount Hanger	4" ⁵	4.5" ⁵	5"	5.5"	6"	6.5"				10"	10.5"	11"
		7"	7.5"	8"	8.5"	9"	9.5"				11.5"	12"	12.5"
	Medium Front Mount Hanger	4" ⁵	4.5" ⁵	5"	5.5"	6"	6.5"	7"	7.5"		12"	12.5"	13"
		8"	8.5"	9"	9.5"	10"	10.5"	11"	11.5"		13.5"	14"	14.5"
	Long Front Mount Hanger	5"	5.5"	6"	6.5"	7"	7.5"	8"	8.5"	9"	14"	14.5"	15"
		9.5"	10"	10.5"	11"	11.5"	12"	12.5"	13"	13.5"	15.5"	16"	16.5"
	Peds PRO ELR 2	Low Pivot Setting	7.5" ³	8" ³	8.5" ³	9" ⁴	9.5" ⁴	10"	10.5"	11"	11.5"	12"	12.5"
High Pivot Setting		7.5" ³	8" ³	8.5" ³	9" ⁴	9.5" ⁴	10"	10.5"	11"	11.5"			
PRO ELR Low Pivot Setting	PRO ELR Low Pivot Setting Short (12" - 17.5")	12" - 17.5"											
	PRO ELR Low Pivot Setting Long (15" - 21")	15" - 21"											
PRO ELR High Pivot Setting	PRO ELR High Pivot Setting Short (10.5" - 15.5")	10.5" - 15.5"											
	PRO ELR High Pivot Setting Long (13" - 19")	13" - 19"											

¹ See Minimum Front Seat Height chart above.

² Comes standard with 3" Extension

³ Not available on chairs 13" wide or more

⁴ Not available on chairs 15" wide or more

⁵ Available with 90° only



Phone (800) 981-1540 | Fax (715) 254-0996
sales@kimobility.com | www.kimobility.com



3" Extension

Available only with Front Mount Hangers.

3" Extension	\$56
3" Extension (Provided as Parts)	\$56

Footplates

Composite	N/C
Not available with Front Mount Hanger or Frame Widths less than 14".	
Composite Angle Adjustable (K0040)	\$200
Not available with Front Mount Hanger or Frame Widths less than 15".	
Aluminum Angle Adjustable (K0040)	\$205
Aluminum Locking Multi-Angle Adjustable (K0108)	\$308
Available only with Front Mount Hanger.	
One Piece Flip Up Angle Adjustable (K0040)	\$200
Not available with Front Mount Hanger, PRO ELR or Frame Widths less than 14".	

Hanger Release

If no selection is made Classic Hanger Release will be used.

Classic	N/C
4-Way - In-Line	N/C
Not available with PRO ELR or Peds PRO ELR.	
4-Way - Out	N/C
Not available with PRO ELR or Peds PRO ELR.	

Heel Loops

Not available with One Piece Flip Up Angle Adjustable Footplate.

Heel Loops (E0951)	\$41
Not available with Frame Widths less than 14".	
Adjustable Heel Loops (E0951)	\$51
Omit Heel Loops	N/C

Calf Straps

Velcro Adjustable (K0038)	\$41
Padded Velcro Adjustable (K0038)	\$72
Bodypoint Padded AeroMesh (K0038)	\$133

Residual Limb Support

	Limb Support	Price	Omit	Cover		Pad Size					
Left	Left Limb Support (E1020)	\$420	Omit Left Hanger	Smooth Cover Left	Air Mesh Cover Left	5X8 Pad Left Side	7X10 Pad Left Side	9x10 Pad Left Side	9x14 Pad Left Side	11x10 Pad Left Side	11x14 Pad Left Side
Right	Right Limb Support (E1020)	\$420	Omit Right Hanger	Smooth Cover Right	Air Mesh Cover Right	5X8 Pad Right Side	7X10 Pad Right Side	9x10 Pad Right Side	9x14 Pad Right Side	11x10 Pad Right Side	11x14 Pad Right Side

Front Frame Height and Casters

Front Seat Height

Select Front Frame Height in row corresponding to desired caster size. Selections in first three columns require Super Low Front Frame.

Caster Size	Front Seat Height												
	Uses Super Low Front Frame ¹												
3"	13.5"	14"		14.5"	15"	15.5"	16"	16.5"	17"	17.5"	18"	18.5"	
4"	13.5"	14"	14.5"	15"	15.5"	16"	16.5"	17"	17.5"	18"	18.5"	19"	19.5" 20"
5"	14.5"	15"	15.5"	16"	16.5"	17"	17.5"	18"	18.5"	19"	19.5"	20"	20.5"
6"	15.5"	16"	16.5"	17"	17.5"	18"	18.5"	19"	19.5"	20"	20.5"		
7"	17"	17.5"		18"	18.5"	19"	19.5"	20"	20.5"				
8"	17.5"	18"	18.5"	19"	19.5"	20"	20.5"						

¹ Price is \$51

Castors

		Caster Size						
		Price	3"	4"	5"	6"	7"	8"
Caster Type	.75" Rollerblade	N/C						
	.75" Lighted Rollerblade	\$77						
	1" Poly	N/C						
	1" Poly Aluminum	\$118						
	1" Pneumatic	\$92						
	1.5" Poly	\$56						
	1.5" Soft Roll Aluminum	\$200						
	2" Poly	\$97						
	2" Pneumatic w/Foam Insert	\$246						

Maximum Seat Dump Reference

(Front Seat Height minus Rear Seat Height)

Seat Depth	12"	13"	14"	15"	16"	17"	18"	19"
Maximum Dump Standard Caster Housing	1.5"	1.5"	1.5"	2"	2"	2"	2.5"	2.5"
Maximum Dump Performance Caster Housing	2.5"	2.5"	3"	3"	3"	3.5"	3.5"	4"



Phone (800) 981-1540 | Fax (715) 254-0996
sales@kimobility.com | www.kimobility.com



Fork Stem

Best Fit	N/C
<i>Ki Mobility will pick Fork and Stem combination to optimize performance based on Caster Size and Front Seat Height</i>	
Standard Fork Stem	N/C
3/4" Longer Fork Stem	N/C
1 1/2" Longer Fork Stem	N/C
Frog Legs Suspension Fork	\$461
<i>Available only with 3", 4", 5" or 6" Caster Sizes. Not available with Caster Pin Locks, 2" Poly Caster Type or 2" Pneumatic w/ Foam Insert Caster Type.</i>	

Fork Options

Frog Leg Options

User Weight Under 80 lbs	N/C
User Weight 80-120 lbs	N/C
User Weight 120-160 lbs	N/C

Caster Pin Locks (K0073)	\$210
<i>Not available with Performance Caster Housing, 1 1/2" Longer Fork Stem or Frog Legs Fork</i>	

Caster Housing

Standard Housing	N/C
Performance Housing	\$51
<i>Not available with 350 lb Weight Capacity</i>	

Wheel Locks

Wheel Lock

Push to Lock	N/C
Pull to Lock	N/C
Omit	N/C

Wheel Lock Options

Extension Handles (E0961)	\$77
<i>Available only with Push to Lock or Pull to Lock Wheel Locks.</i>	
Grade Aids (E0974)	\$154
<i>Available only with Push to Lock Wheel Locks and Pneumatic and Pneumatic with Airless Insert tires</i>	

Rear Frame Height and Wheels

Rear Seat Height

Select Rear Frame Height in row corresponding to Wheel Size and Tire Profile. 12"-14" Frame Depths not available with 24" wheel size.

Wheel Size	Tire Profile	Rear Seat Height										
		11" ²	11.5" ¹	12" ²	12.5" ¹	13" ²	13.5"	14" ²	14.5"	15" ²	15.5" ²	16" ²
18"	Low	11" ²	11.5" ¹	12" ²	12.5" ¹	13" ²	13.5"	14" ²	14.5"	15"	15.5" ²	16" ²
	Full	11.5" ²	12" ¹	12.5" ²	13" ¹	13.5" ²	14"	14.5" ²	15"	15.5" ²	16" ²	16.5" ²
20"	Low	11.5" ²	12" ¹	12.5" ²	13" ¹	13.5" ²	14"	14.5" ²	15"	15.5" ²	16" ²	16.5" ²
	Full	12" ²	12.5" ¹	13" ²	13.5" ¹	14" ²	14.5"	15" ²	15.5"	16" ²	16.5" ²	17" ²
22"	Low	13" ²	13.5" ¹	14" ²	14.5" ¹	15" ²	15.5"	16" ²	16.5"	17" ²	17.5" ²	18" ²
	Full	13.5" ²	14" ¹	14.5" ²	15" ¹	15.5" ²	16"	16.5" ²	17"	17.5" ²	18" ²	18.5" ²
24"	Low	13.5" ²	14" ¹	14.5" ²	15" ¹	15.5" ²	16"	16.5" ²	17"	17.5" ²	18" ²	18.5" ²
	Full	14" ²	14.5" ¹	15" ²	15.5" ¹	16" ²	16.5"	17" ²	17.5"	18" ²	18.5" ²	19" ²

1 Amputee Axle Plate not available with lowest two Rear Frame Height positions

One Arm Drive

One Arm Drive - Quick Release (E0958)	\$1122
<i>Available only with 20", 22" and 24" spoke wheels and Quick Release Axles. Spoke wheel style is exclusive to OAD. Not available with 4 degrees of camber. Standard handrim will be Aluminum Anodized.</i>	

Axle Plates

No selection needed with One Arm Drive.

Universal	N/C
Amputee (E0959)	\$187
<i>Not available with Wheelchair Drive Ready Option.</i>	

Drive Handrim Side

Left Side Drive	N/C
Right Side Drive	N/C

Drive Handrim Type

Aluminum Drive Handrim	N/C
Poly Coated Drive Handrim	\$205



Phone (800) 981-1540 | Fax (715) 254-0996
sales@kimobility.com | www.kimobility.com



Handrims

No selection needed with One Arm Drive.

Standard Handrims

Aluminum Anodized	N/C
Plastic Coated	\$108
Projection (8 Vertical Projections)	\$359
<i>Available only with full profile tires in combination with 22", 24", and 26" MAXX Spoke wheels, or 22" and 24" Mag wheels</i>	
Omit Handrim	N/C

Ergonomic Handrim

Not available with 20" or 26" wheels. Tab mount only available with 24" wheels.

Lite Regular Grip	\$477
Standard Regular Grip	\$456
Lite Extra Grip	\$528
Standard Extra Grip	\$507
Lite Omit Grip	\$405
Standard Omit Grip	\$354

Wheels

		Rear Wheel Size				
		Price	18"	20"	22"	24"
Rear Wheel Type	MAXX Spoke	\$128				
	Mag	N/C				
	Colored Mag ¹	\$108				
	Spinergy Spox	\$938				

Spinergy Option

Select spoke color and handrim tab for Spinergy wheels.

Black Spokes Red Spokes Blue Spokes White Spokes _____ N/C

Axles

Axles

No selection needed with One Arm Drive.

Quick Release	N/C
Quad Quick Release	\$190
<i>Not available with Mag or Spinergy wheels.</i>	

Handrim Options

3M Micro Grip	\$72
<i>Not available on Plastic Coated or Projection Handrims.</i>	

Handrim Mount

Close Mount	\$31
<i>If no selection is made, Standard Mount will be used</i>	

Rear Tire Type

Full Profile

Street Pneumatic 1 3/8" (E2211, E2212)	N/C
Full Poly 1 3/8" (E2218)	\$56
Pneumatic with Airless Insert 1 3/8" (E2211, E2213)	\$97

Low Profile

Low Poly 1" (E2218)	\$56
Sentinel Solid Tire 1"	\$97
<i>Not available on 18" wheel size or 20" and 22" Mag wheels.</i>	
High Pressure 1"	\$97
<i>Not available on Mag wheel or 18" wheel size.</i>	

Axle Position

Rear Wheel Camber

0° 2° 4°	N/C
--------------	-----

Center of Gravity

No selection needed with Amputee Axle Plate

Rear Center Forward	N/C
-------------------------	-----

Armrests

Armrest Type	Length		Price
Height Adjustable T-Arm Tall (11.5"-15") E0973	Desk		\$323
	Full		
Height Adjustable T-Arm (8.5"-13.5") E0973 ¹	Desk		\$323
	Full		
Height Adjustable T-Arm Low (6.5"-10.5") E0973 ²	Desk		\$323
	Full		

Armrest Type	Length		Price
Angle Adjustable Locking Extendable Flip Up E0973 ⁵	Desk		\$272
	Full		
Pediatric Height Adjustable T-Arm Desk (6.5"-10")			\$231
Tubular Flip Up ³			N/C
Omit Armrests			N/C

¹ Not available with Frame Depths less than 14".

² Not available with Frame Depths less than 14". If requested configuration is unavailable, Standard Height T-Arm will be used.

³ Not available with Arm Pads.

Arm Pads

Standard Arm Pad	N/C
Waterfall Arm Pad	\$26
Not available with Angle Adjustable Locking Extendable Flip Up Armrest	
Foam Grip	N/C
Available only with Angle Adjustable Locking Extendable Flip Up Armrest	

Accessories


Removable Side Guards

Not available with Height Adjustable T-Arm Armrests.

Adult Composite Side Guards (14" - 16") (K0108)	\$159
Pediatric Composite Side Guard (10" - 14") (K0108)	\$159

Pelvic Positioning Belts

1" Auto Buckle Pelvic Positioning Belt (E0978)	\$87
1" Padded Auto Buckle Pelvic Positioning Belt (E0978)	\$179
1 1/2" Auto Buckle Pelvic Positioning Belt (E0978)	\$87
1 1/2" Padded Auto Buckle Pelvic Positioning Belt (E0978)	\$179
2" Airline Buckle Pelvic Positioning Belt (E0978)	\$118
Small Bodypoint Evoflex (16-24") (E0978)	\$200
Medium Bodypoint Evoflex (19-29") (E0978)	\$200
Large Bodypoint Evoflex (25-37") (E0978)	\$200
Belt Mounting Band Clamp	\$36



Anti-tips

Rear Anti-Tips (E0971)	\$154
Omit	N/C

Spoke Guards

Clear Spoke Guards (K0065)	\$123
Not available with Mag wheel.	
Black Spoke Guards (K0065)	\$123
Not available with 18" or 20" wheel sizes or Mag wheel.	
Blue Monster Spoke Guards (K0065) (1)	\$261
Robo Fun Spoke Guards (K0065) (1)	\$261
Tutu Cute Spoke Guards (K0065) (1)	\$261
Tuxedo Time Spoke Guards (K0065) (1)	\$261

¹ Available only with 20" and 22" wheel sizes. Not available with Mag wheel.

Other Accessories

Tool Kit	\$46
Removable Underseat Pouch	\$51
Not available with 10" - 13" Seat Widths.	
Pediatric Ki Mobility Backpack	\$67
Ki Mobility Backpack	\$67
Black Neoprene Impact Guards	\$97
Cane and Crutch Holder	\$210
Canopy with Hardware (K0108)	\$287
See image to the left	



Phone (800) 981-1540 | Fax (715) 254-0996
sales@kimobility.com | www.kimobility.com



Frame Color

Gloss Finish		Matte Finish	
Black	N/C	Military Green	N/C
Blue Skies	N/C	Charcoal Gray	N/C
Bubblegum Pink	N/C	Emerald Green	N/C
Candy Blue	N/C	Iridescent Burgundy	N/C
Candy Purple	N/C	Iridescent Navy	N/C
Candy Red	N/C	Retro Red	N/C
Copper	N/C	Royal Blue	N/C
Granny Smith	N/C	Shadow Black	N/C
Metallic Orange	N/C	Smooth White	N/C
Yellow	N/C		

Hydrographics

AOR1 (Desert Digital Camo)	\$825
Pink Camo	\$825
Muddy Girl	\$825
NEXT Camo	\$825
Carbon Fiber	\$825
<i>Select Hydrographics Base Color.</i>	
Zebra	\$825
<i>Select Hydrographics Base Color. Not available in Black.</i>	

Hydrographics Base Color

Black	Blue	Granny Smith	Pink	Purple	Red	N/C
-------	------	--------------	------	--------	-----	-----

Mag Wheel Color

Gloss Finish

Black	Blue Skies	Bubblegum Pink	Candy Blue	Candy Purple	Candy Red	Copper	Granny Smith	N/C
Metallic Orange	Yellow							

Matte Finish

Military Green	Charcoal Gray	Emerald Green	Iridescent Burgundy	Iridescent Navy	Retro Red			
Royal Blue	Shadow Black	Smooth White						N/C

Spark Front Frame Height Matrix																														
3" Caster				4" Caster				5" Caster				6" Caster				7" Caster				8" Caster										
FSH	Super Low	Stem	Fork	FSH	Super Low	Stem	Fork	FSH	Super Low	Stem	Fork	FSH	Super Low	Stem	Fork	FSH	Super Low	Stem	Fork	FSH	Super Low	Stem	Fork							
14.5	13.5	Std	4"	15	13.5	Std	4"	16	14.5	Std	4"	17	15.5	Std	5"	18	16.5	Std	6"	19	17.5	Std	6"							
15	14			15.5	14			16.5	15.5	3/4		17.5	16.5	3/4		18.5	17			19.5	18.5	3/4								
15.5	14	3/4		15.5 ¹	14 ¹	3/4		17.5	16	1 1/2		18.5	17	1 1/2		18.5	17.5	3/4		20.5	n/a	1 1/2								
16	14.5			16	14.5			Std	17	15.5	3/4	17.5	16	Std	19	18	1 1/2			19	17.5	Std								
16	15	1 1/2		16.5 ¹	15 ¹	1 1/2			16.5	15		17.5	16		18	16.5		19.5		18	20		18.5							
16.5	15.5			17	15.5			16.5	15.5	3/4	18	16.5	5"	17.5	16.5	3/4	18.5	17		Std	20.5		19	7"						
15.5	14.5	Std	16	14.5	Std	17.5	16	1 1/2	18		17	3/4		19	17.5		Std	19.5	18.5		3/4									
16	15		3/4	16		15	3/4		18	16.5	Std		18.5	17.5	1 1/2	19.5		18	3/4	20		19								
17	15.5	1 1/2		16.5	15.5	1 1/2		16.5	15	Std		18.5	17	Std		19	17.5	Std		20.5		n/a	1 1/2							
17	16		17	15.5	17		15.5	18	16.5		19	17.5	19.5		18	20	18.5		20.5	n/a										
17.5	16.5	Std	5"	17.5	16	1 1/2	5"	17.5	16	Std	5"	19.5	18	1 1/2	6"	20	18.5	3/4	7"	21	n/a	1 1/2								
16.5	15.5			16.5	15			17.5	16			18	16.5			19	17.5			20	18.5		20.5	n/a						
17	16	17		15.5	18	16.5		19	17.5	20		18.5	20.5	n/a		21	n/a													
17.5	16.5	18		16.5	19	17.5		20	18.5	20.5		n/a	21	n/a		6" x 2" caster not available with 5" fork														
16.5	15.5	Std		6"	16.5	15		Std	6"	17.5		16	3/4	6"						18	16.5	Std	7"	20	19	1 1/2	7"	21	n/a	
17	16				17	15.5				18		16.5								19	17.5			20	18.5			20.5	n/a	
17.5	16	3/4	17		16	3/4	18	16.5		1 1/2	18.5	17	1 1/2		19					17.5	1 1/2	20		18.5	1 1/2					
18	16.5		17.5		16.5		18	16.5			19	17.5			20					18.5		20.5		19						
18	17	1 1/2	18		16.5	1 1/2	18.5	17		1 1/2	18.5	17.5	3/4		19					18	3/4	21		n/a						
18.5	17.5		18.5		17		19	17.5			19	18			20					18.5		20.5		19						
				18	16.5	Std	7"	18	16.5	Std	7"	19.5	18.5	1 1/2	7"															
				18.5	17			18.5	17			20	19																	
				18.5	17.5	3/4		19	17.5	3/4		19.5	18	1 1/2														20	18.5	
				19	18			18.5	17.5			20	18.5																	
				19.5	18.5	1 1/2		19	18	3/4		20.5	19	1 1/2														21	n/a	
				20	19			19.5	18.5			20	18.5			20.5	19													

¹ Not available with 4 x 1.5" Soft Roll Aluminum

Seat heights with "FSH" heading available with Standard Front Frame. Seat heights with "Super Low" heading available with Super Low Front Frame. Actual seat height may vary depending upon configuration.