



Phone (800) 981-1540 | Fax (715) 254-0996
 sales@kimobility.com | www.kimobility.com

Effective 05/23/2017

Company: _____

Ship To: _____

Account Number: _____ Date: _____

Address: _____

Quote PO Number: _____

City: _____ State: _____

RTS: _____

Zip: _____ Phone: _____

Marked For: _____

Email/Fax: _____

- Spark - Transit Included** E1238 **\$2320**
 165 lb weight capacity. Not available with upholstery.
- Spark - Transit Not Included** E1238 **\$2000**
 165 lb weight capacity. Not available with upholstery.

HCPCS codes provided are based on code verification by the PDAC or our interpretation of the code definitions.

Proper use of HCPCS codes does not ensure coverage or payment. For coverage information, verify the appropriate payer's coverage policy.

Frame

Frame Width / Frame Depth


Select box corresponding to desired Frame Width and Depth. See references at right.

		Frame Width						
		10"	11"	12"	13"	14"	15"	16"
Frame Depth	12"	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>
	13"	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>
	14"	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>
	15"	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>
	16"	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>
	17"	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>
	18"	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>



Frame Width

Spark is set up to allow 2" increase in frame width from initial frame width in 1" increments. Growth is achieved by expanding cross braces. 16" maximum frame width. 10" Frame Width expands only 1".



Frame Depth

Spark is set up to allow 4" increase in frame depth from initial frame depth in 1" increments. 19" maximum frame depth.

Transit Option

- Transit System** K0108 **\$320**
 No charge when selected with Spark - Transit Included. Not available with Removable Seat Pan.

Wheelchair Drive Ready

Chair will ship from Ki Mobility with adaption kit installed for selected Alber system. Order system separately from Alber USA. Select Rear Frame Height (full profile selections only) in drive wheel size on Rear Frame Height chart (page 4). Requires Rear Anti-Tips. Not available with Spinergy wheels.

- e-motion Ready _____ N/C
- twion Ready _____ N/C

Backrest

Back Post / Height Range / Back Height

Select Back Height in row corresponding to desired Back Post and Height Range. Both Back Post options are folding and angle adjustable in 4" increments.

Back Post	Price	Height Range	Back Height					
Standard Depth Adjustable	N/C	Short	<input type="checkbox"/> 20					
		Tall	<input type="checkbox"/> 24					
Height Adjustable Depth Adjustable	\$75	Short ¹	<input type="checkbox"/> 13	<input type="checkbox"/> 14	<input type="checkbox"/> 15	<input type="checkbox"/> 16	<input type="checkbox"/> 17	
		Medium	<input type="checkbox"/> 16	<input type="checkbox"/> 17	<input type="checkbox"/> 18	<input type="checkbox"/> 19	<input type="checkbox"/> 20	
		Tall	<input type="checkbox"/> 20	<input type="checkbox"/> 21	<input type="checkbox"/> 22	<input type="checkbox"/> 23	<input type="checkbox"/> 24	

¹ Not available with Seating Dynamics Rocker Back

Back Post Options

- Stroller Handle Extensions** _____ **\$75**
 Not available with Height Adjustable Depth Adjustable Back Posts
- Seating Dynamics Rocker Back** K0108 **\$1995**
 Available only with T-Arms. Not available with 12" Frame Depths with Rear Frame Heights lower than 14.5".

Footrest

Hanger / Front Frame Bend		Footplate	
<input type="checkbox"/> 60° _____ N/C	Available only with Extension Tube Hangers	<input type="checkbox"/> Composite _____ N/C	Not available with Front Mount Hangers, Frame Widths less than 14" or PRO ELR
<input type="checkbox"/> 70° _____ N/C		<input type="checkbox"/> Composite Angle Adjustable K0040 _____ \$185	Not available with Front Mount Hangers, ELR or Frame Widths less than 15"
<input type="checkbox"/> 80° _____ N/C	Not available with 7" or 8" Casters	<input type="checkbox"/> Aluminum Angle Adjustable K0040 _____ \$185	
<input type="checkbox"/> 90° _____ N/C	Available only with Front Mount Hangers. Not available with 7" or 8" Casters.	<input type="checkbox"/> One Piece Flip Up Angle Adjustable K0040 _____ \$190	Not available with Front Mount Hangers, ELR, PRO ELR or Frame Widths less than 14"
<input type="checkbox"/> Elevating Leg Rests (ELR) E0990 _____ \$295	Not available with Seat Widths less than 14"		
<input type="checkbox"/> Pro Elevating Leg Rests (PRO ELR) E0990 New! _____ \$350			
<input type="checkbox"/> Omit Hangers _____ - \$20			

Footrest / Extension / Seat to Footrest Length

Select Seat to Footrest Length in row corresponding to desired Footrest and Extension. Footrest Length of at least 2.5" shorter than Front Seat Height recommended. Extension Tube Hanger measurements are based on 70° Hanger / Front Frame Bend.

Footrest	Extension	Seat to Footrest Length												
Extension Tube Hangers (60°, 70°, 80°)	Short	<input type="checkbox"/> 11	<input type="checkbox"/> 11.5	<input type="checkbox"/> 12	<input type="checkbox"/> 12.5	<input type="checkbox"/> 13	<input type="checkbox"/> 13.5	<input type="checkbox"/> 14	<input type="checkbox"/> 14.5	<input type="checkbox"/> 15				
	Medium	<input type="checkbox"/> 13	<input type="checkbox"/> 13.5	<input type="checkbox"/> 14	<input type="checkbox"/> 14.5	<input type="checkbox"/> 15	<input type="checkbox"/> 15.5	<input type="checkbox"/> 16	<input type="checkbox"/> 16.5	<input type="checkbox"/> 17				
	Long	<input type="checkbox"/> 15.5	<input type="checkbox"/> 16	<input type="checkbox"/> 16.5	<input type="checkbox"/> 17	<input type="checkbox"/> 17.5	<input type="checkbox"/> 18	<input type="checkbox"/> 18.5	<input type="checkbox"/> 19	<input type="checkbox"/> 19.5				
Front Mount Hangers ¹	Standard Selections											Requires 3" Extension for \$55		
	Short (70°, 80°, 90°)	<input type="checkbox"/> 5	<input type="checkbox"/> 5.5	<input type="checkbox"/> 6	<input type="checkbox"/> 6.5	<input type="checkbox"/> 7					<input type="checkbox"/> 10	<input type="checkbox"/> 10.5	<input type="checkbox"/> 11	
		<input type="checkbox"/> 7.5	<input type="checkbox"/> 8	<input type="checkbox"/> 8.5	<input type="checkbox"/> 9	<input type="checkbox"/> 9.5					<input type="checkbox"/> 11.5	<input type="checkbox"/> 12	<input type="checkbox"/> 12.5	
	Medium (70°, 80°, 90°)	<input type="checkbox"/> 5	<input type="checkbox"/> 5.5	<input type="checkbox"/> 6	<input type="checkbox"/> 6.5	<input type="checkbox"/> 7	<input type="checkbox"/> 7.5	<input type="checkbox"/> 8			<input type="checkbox"/> 12	<input type="checkbox"/> 12.5	<input type="checkbox"/> 13	
		<input type="checkbox"/> 8.5	<input type="checkbox"/> 9	<input type="checkbox"/> 9.5	<input type="checkbox"/> 10	<input type="checkbox"/> 10.5	<input type="checkbox"/> 11	<input type="checkbox"/> 11.5			<input type="checkbox"/> 13.5	<input type="checkbox"/> 14	<input type="checkbox"/> 14.5	
	Long (70°, 80°)	<input type="checkbox"/> 5	<input type="checkbox"/> 5.5	<input type="checkbox"/> 6	<input type="checkbox"/> 6.5	<input type="checkbox"/> 7	<input type="checkbox"/> 7.5	<input type="checkbox"/> 8	<input type="checkbox"/> 8.5	<input type="checkbox"/> 9	<input type="checkbox"/> 14	<input type="checkbox"/> 14.5	<input type="checkbox"/> 15	
<input type="checkbox"/> 9.5		<input type="checkbox"/> 10	<input type="checkbox"/> 10.5	<input type="checkbox"/> 11	<input type="checkbox"/> 11.5	<input type="checkbox"/> 12	<input type="checkbox"/> 12.5	<input type="checkbox"/> 13	<input type="checkbox"/> 13.5	<input type="checkbox"/> 15.5	<input type="checkbox"/> 16	<input type="checkbox"/> 16.5		
Elevating Leg Rests (ELR)	Short	<input type="checkbox"/> 12	<input type="checkbox"/> 13	<input type="checkbox"/> 14	<input type="checkbox"/> 15	<input type="checkbox"/> 16								
	Long	<input type="checkbox"/> 13	<input type="checkbox"/> 14	<input type="checkbox"/> 15	<input type="checkbox"/> 16	<input type="checkbox"/> 17	<input type="checkbox"/> 18							
Extended Length ELR	Short	<input type="checkbox"/> 18	<input type="checkbox"/> 19	<input type="checkbox"/> 20	<input type="checkbox"/> 21	<input type="checkbox"/> 22								
	Long	<input type="checkbox"/> 18	<input type="checkbox"/> 19	<input type="checkbox"/> 20	<input type="checkbox"/> 21	<input type="checkbox"/> 22	<input type="checkbox"/> 23	<input type="checkbox"/> 24						
PRO ELR Low Pivot	Short	<input type="checkbox"/> 12 - 17.5												
	Long	<input type="checkbox"/> 15 - 21												
PRO ELR High Pivot	Short	<input type="checkbox"/> 10.5 - 15.5												
	Long	<input type="checkbox"/> 13 - 19												

¹ Available only with Aluminum Angle Adjustable Footplates. See Minimum Front Seat Height chart below.

Minimum Front Seat Height with Front Mount Hangers						
Hanger Length	Short		Medium		Long	
	70°	80/90°	70°	80/90°	70°	80°
Hanger Bend						
Minimum Front Seat Height with Standard Selection	13"	13.5"	14"	15"	16"	17"
Minimum Front Seat Height with 3" Extension	14.5"	16"	16.5"	18"	18.5"	20"

3" Extension



Available only with Front Mount Hangers. Provided as parts.

3" Extension for Front Mount Hanger _____ \$55

Hanger Release

If no selection is made Classic Hanger Release will be used

- Classic _____ Std
- 4-Way - In line **New!** _____ N/C
See reference A to the right. Not available with ELR or PRO ELR.
- 4-Way - Out **New!** _____ N/C
See reference B to the right



Footrest continued

- Heel Loop**
- Heel Loops E0951 \$30
Not available with One Piece Flip Up Angle Adjustable Footplate or Frame Widths less than 14"
 - Omit Heel Loops N/C

- Calf Strap**
- Velcro Adjustable K0038 \$35
 - Padded Velcro Adjustable K0038 \$68
 - Bodypoint Padded AeroMesh K0038 \$102

- Residual Limb Support**
- Residual Limb Support E1020
 - Left \$395
 - Right \$395
 - Omit Left Hanger - \$10
 - Omit Right Hanger - \$10
- Cover**
- Smooth
 - Air Mesh
- Pad Size**
- 5 x 8
 - 7 x 10
 - 9 x 10
 - 9 x 14
 - 11 x 10
 - 11 x 14

Front Frame Height and Casters

Caster Size / Front Frame Type / Front Frame Height

Select Front Frame Height in row corresponding to desired caster size. Selections in first three columns require Super Low Front Frame.

Caster Size	Front Frame Height													
	Requires Super Low Front Frame ¹													
3"	<input type="checkbox"/> 13.5	<input type="checkbox"/> 14		<input type="checkbox"/> 14.5	<input type="checkbox"/> 15	<input type="checkbox"/> 15.5	<input type="checkbox"/> 16	<input type="checkbox"/> 16.5	<input type="checkbox"/> 17	<input type="checkbox"/> 17.5	<input type="checkbox"/> 18	<input type="checkbox"/> 18.5		
4"	<input type="checkbox"/> 13.5 ²	<input type="checkbox"/> 14	<input type="checkbox"/> 14.5	<input type="checkbox"/> 15 ³	<input type="checkbox"/> 15.5	<input type="checkbox"/> 16	<input type="checkbox"/> 16.5	<input type="checkbox"/> 17	<input type="checkbox"/> 17.5	<input type="checkbox"/> 18	<input type="checkbox"/> 18.5	<input type="checkbox"/> 19	<input type="checkbox"/> 19.5	<input type="checkbox"/> 20
5"	<input type="checkbox"/> 14.5	<input type="checkbox"/> 15	<input type="checkbox"/> 15.5	<input type="checkbox"/> 16	<input type="checkbox"/> 16.5	<input type="checkbox"/> 17	<input type="checkbox"/> 17.5	<input type="checkbox"/> 18	<input type="checkbox"/> 18.5	<input type="checkbox"/> 19	<input type="checkbox"/> 19.5	<input type="checkbox"/> 20	<input type="checkbox"/> 20.5	
6"	<input type="checkbox"/> 15.5	<input type="checkbox"/> 16	<input type="checkbox"/> 16.5	<input type="checkbox"/> 17	<input type="checkbox"/> 17.5	<input type="checkbox"/> 18	<input type="checkbox"/> 18.5	<input type="checkbox"/> 19	<input type="checkbox"/> 19.5	<input type="checkbox"/> 20	<input type="checkbox"/> 20.5			
7"	<input type="checkbox"/> 17	<input type="checkbox"/> 17.5		<input type="checkbox"/> 18	<input type="checkbox"/> 18.5	<input type="checkbox"/> 19	<input type="checkbox"/> 19.5	<input type="checkbox"/> 20	<input type="checkbox"/> 20.5					
8"	<input type="checkbox"/> 17.5	<input type="checkbox"/> 18	<input type="checkbox"/> 18.5	<input type="checkbox"/> 19	<input type="checkbox"/> 19.5	<input type="checkbox"/> 20	<input type="checkbox"/> 20.5							

¹ Price is \$50
² Front Frame Height not available with 4 x 1" Poly Caster
³ Super Low Front Frame required when using 4 x 1" Poly Caster

- Caster Type**
- .75" Rollerblade N/C
Not available with 5", 6", 7" and 8" Caster Sizes
 - .75" Lighted Rollerblade \$50
Not available with 6", 7" and 8" Caster Sizes
 - 1" Poly N/C
Not available with 3" Caster Size
 - 1" Poly Aluminum \$110
Not available with 3", 6", 7" and 8" Caster Sizes
 - 1" Pneumatic \$80
Not available with 3", 4", 5" and 7" Caster Sizes
 - 1.5" Poly \$50
Not available with 3", 7" and 8" Caster Sizes
 - 1.5" Soft Roll Aluminum \$175
Not available with 7" and 8" Caster Sizes
 - 2" Poly \$70
Not available with 3", 4", 5", 7" and 8" Caster Sizes
 - 2" Pneumatic w/ Foam Insert \$225
Not available with 3", 4", 5", 6" and 7" Caster Sizes

- New! Caster Housing**
- Standard N/C
 - Performance \$35

Wheel Locks

- Wheel Lock**
- Push to Lock Std
 - Pull to Lock N/C
 - Omit Wheel Locks N/C

Maximum Seat Dump Reference
 (Front Seat Height minus Rear Seat Height)

Seat Depth	12"	13"	14"	15"	16"	17"	18"
Maximum Dump Standard Caster Housing	1.5"	1.5"	1.5"	2"	2"	2"	2.5"
Maximum Dump Performance Caster Housing	2"	2"	2.5"	2.5"	3"	3"	3"

- Fork Stem**
- Best Fit Std
Ki Mobility will pick optimal Fork Stem based on Caster Size and Front Frame Height
 - Standard Fork Stem N/C
 - 3/4" Longer Fork Stem N/C
 - 1 1/2" Longer Fork Stem N/C

- Fork Options**
- Frog Legs Ultra-Sport Fork - Black \$425
Available only with 3", 4", 5" or 6" Casters. Not available with Caster Pin Locks, 2" Poly Caster Type or 2" Pneumatic w/ Foam Insert Caster Type.
- End User Weight Required: _____ lbs**
- Caster Pin Locks K0073 \$195
Not available with Performance Caster Housing, 1 1/2" Longer Fork Stem or Frog Legs Fork

- Wheel Lock Options**
- Extension Handles E0961 \$67
Available only with Push to Lock or Pull to Lock Wheel Locks
 - Grade Aids E0974 \$150
Available only with Push to Lock Wheel Locks and Pneumatic and Pneumatic with Airless Insert Rear Tire Types

Rear Frame Height and Wheels

Wheel Size / Tire Profile / Rear Frame Height

Select Rear Frame Height in row corresponding to Wheel Size and Tire Profile. 14" Seat Depth not available with 24" Rear Wheels.

Requires curved axle plate Any axle plate can be used (Curved, Standard, Amputee, One Arm Drive)

Wheel Size	Tire Profile	Rear Frame Height									
18"	Low	<input type="radio"/> 11	<input type="checkbox"/> 11.5 ¹	<input type="radio"/> 12	<input type="checkbox"/> 12.5 ¹	<input type="radio"/> 13	<input type="checkbox"/> 13.5	<input type="radio"/> 14	<input type="checkbox"/> 14.5	<input type="radio"/> 15	<input type="radio"/> 15.5
	Full	<input type="radio"/> 11.5	<input type="checkbox"/> 12 ¹	<input type="radio"/> 12.5	<input type="checkbox"/> 13 ¹	<input type="radio"/> 13.5	<input type="checkbox"/> 14	<input type="radio"/> 14.5	<input type="checkbox"/> 15	<input type="radio"/> 15.5	<input type="radio"/> 16
20"	Low	<input type="radio"/> 11.5	<input type="checkbox"/> 12 ¹	<input type="radio"/> 12.5	<input type="checkbox"/> 13 ¹	<input type="radio"/> 13.5	<input type="checkbox"/> 14	<input type="radio"/> 14.5	<input type="checkbox"/> 15	<input type="radio"/> 15.5	<input type="radio"/> 16
	Full	<input type="radio"/> 12	<input type="checkbox"/> 12.5 ¹	<input type="radio"/> 13	<input type="checkbox"/> 13.5 ¹	<input type="radio"/> 14	<input type="checkbox"/> 14.5	<input type="radio"/> 15	<input type="checkbox"/> 15.5	<input type="radio"/> 16	<input type="radio"/> 16.5
22"	Low	<input type="radio"/> 12.5	<input type="checkbox"/> 13 ¹	<input type="radio"/> 13.5	<input type="checkbox"/> 14 ¹	<input type="radio"/> 14.5	<input type="checkbox"/> 15	<input type="radio"/> 15.5	<input type="checkbox"/> 16	<input type="radio"/> 16.5	<input type="radio"/> 17
	Full	<input type="radio"/> 13	<input type="checkbox"/> 13.5 ¹	<input type="radio"/> 14	<input type="checkbox"/> 14.5 ¹	<input type="radio"/> 15	<input type="checkbox"/> 15.5	<input type="radio"/> 16	<input type="checkbox"/> 16.5	<input type="radio"/> 17	<input type="radio"/> 17.5
24"	Low	<input type="radio"/> 13.5	<input type="checkbox"/> 14 ¹	<input type="radio"/> 14.5	<input type="checkbox"/> 15 ¹	<input type="radio"/> 15.5	<input type="checkbox"/> 16	<input type="radio"/> 16.5	<input type="checkbox"/> 17	<input type="radio"/> 17.5	<input type="radio"/> 18
	Full	<input type="radio"/> 14	<input type="checkbox"/> 14.5 ¹	<input type="radio"/> 15	<input type="checkbox"/> 15.5 ¹	<input type="radio"/> 16	<input type="checkbox"/> 16.5	<input type="radio"/> 17	<input type="checkbox"/> 17.5	<input type="radio"/> 18	<input type="radio"/> 18.5

¹ Amputee Axle Plate not available with lowest two Rear Frame Height positions

One Arm Drive

- One Arm Drive E0958 _____ \$1095
Available only with 24" Mag wheels. Not available with 14" Rear Seat Height in combination with 24" rear wheel, 4° camber or 10" Frame Width. Standard Handrim will be Aluminum Anodized. Axle not quick release.
- Drive Handrim Side: Left Right
- Drive Handrim Type: Aluminum _____ N/C
 Poly Coated _____ \$200

Axle Plate

No selection needed with One Arm Drive

- Standard _____ Std
 Curved _____ \$85
 Amputee _____ \$175
Not available with Wheelchair Drive Ready Option

Rear Wheel Type

- Spoke _____ \$100
 Mag (Black) _____ N/C
Not available in 18" rear wheel
- Colored Mag _____ \$100
Available only in 20" and 22" rear wheels
Color: Florescent Red White Blue Pink
- Spinergy Spox _____ \$795
Not available in 18" or 20" rear wheels
Color: Black Red White Blue Yellow
Handrim: Short Tab Long Tab
Available only with Aluminum Anodized and Plastic Coated Handrims

Rear Tire Type

Full Profile

- Street Pneumatic 1 3/8" E2211, E2212 _____ Std
 Full Poly 1 3/8" E2218 _____ \$55
 Pneumatic with Airless Insert 1 3/8" E2211, E2213 _____ \$95

Axles

Axle

No selection needed with Wheelchair Drive Ready Option or One Arm Drive

- Quick Release _____ Std
 Quad Quick Release _____ \$175
Not available with Mag wheels

Handrim

No selection needed with One Arm Drive

Standard Handrims

- Aluminum Anodized _____ Std
 Plastic Coated _____ \$95
 Projection (8 Vertical Projections) _____ \$295
Available only on 22" and 24" Spoke wheels & 24" Mag wheels w/full tires
 Omit Handrim _____ N/C

Natural Fit Handrims

Not available on 18" or 20" rear wheels or 22" Mag wheels. Not recommended with Mag wheels.

- Natural Fit LT Regular Grip _____ \$445
 Natural Fit Standard Regular Grip _____ \$445
 Natural Fit LT Super Grip _____ \$495
 Natural Fit Standard Super Grip _____ \$495
 Natural Fit LT Omit Grip _____ \$345
 Natural Fit Standard Omit Grip _____ \$345

Handrim Option

- 3M Micro Grip _____ \$65
Not available on Plastic Coated or Projection Handrims

Low Profile

- Low Poly 1" E2218 _____ \$55
 Sentinel Solid Tire 1" E2218 _____ \$95
Not available on 18" rear wheel
 High Pressure 1" E2211, E2212 _____ \$95
Not available on Mag or 18" rear wheel

Axle Position

- Rear Wheel Camber: 0° 2° 4°
Center of Gravity: Rear Center Forward

DCN0033.24 Pricing and specifications subject to change without notice

Armrests

- Armrest** Please select Arm Pad option
- Height Adjustable T-Arm (8.5" - 13.5") E0973 _____ \$265
Not available with Frame Depths less than 14"
 - Height Adjustable T-Arm Low (7.5" - 11.5") E0973 _____ \$265
 - Pediatric Height Adjustable T-Arm Desk (6.5" - 10") E0973 _____ \$185
 - Angle Adjustable Locking Extendable Flip Up E0973 _____ \$265

- Arm Pad**
- Not available with Tubular Flip Up Armrests. Full length Arm Pads not available with Pediatric Height Adjustable T-Arms.**
- Standard Arm Pad _____ N/C Desk Full
 - Waterfall Arm Pad _____ \$25 Desk Full
Not available with Angle Adjustable Locking Extendable Flip Up Armrest
 - Foam Grip _____ N/C Desk Full
Available only with Angle Adjustable Locking Extendable Flip Up Armrest

- Armrest continued**
- Tubular Flip Up _____ N/C
Not available with Arm Pads
 - Omit Armrests _____ N/C

Accessories

- Removable Side Guard**
- Not available with Height Adjustable T-Arm Armrests**
- Adult Composite Side Guards (14" - 16") K0108 _____ \$140
 - Pediatric Composite Side Guards (10" - 14") K0108 _____ \$140

- Anti-Tip**
- Rear Anti-Tips E0971 _____ \$140

- Positioning Belt**
- 1" Auto Buckle E0978 _____ \$70
 - 1" Padded Auto Buckle E0978 _____ \$155
 - 1 1/2" Auto Buckle E0978 _____ \$70
 - 1 1/2" Padded Auto Buckle E0978 _____ \$155
 - 2" Airline Buckle E0978 _____ \$110
 - Bodypoint Evoflex - Pelvic Stabilizer E0978 _____ \$183
Size: Small (16-24") Medium (19-29") Large (25-37")
 - Bodypoint Belt Mounting Kit _____ \$105
Style: Bracket Band Clamp 1"

- Seat Pan**
- Removable Aluminum Seat Pan E2231 _____ \$245
Not available with Transit System. Pan depth 1/2" shorter than Frame Depth.

- Spoke Guards**
- Clear Spoke Guards K0065 _____ \$115
Not available with Mag wheels
 - Black Spoke Guards K0065 _____ \$115
Not available with 18" or 20" rear wheels or Mag wheels
 - Designer Spoke Guards K0065 _____ \$200
Available only with 20" and 22" rear wheels. Not available with Mag wheels.
- Design:**
- Blue Monster Fairy Garden Havin' a Hoot
 - Robo Fun Rockin' Reptiles Super Puppies
 - Tutu Cute Tuxedo Time

- Other Accessories**
- Tool Kit _____ \$45
 - Removable Underseat Pouch _____ \$50
Not available with 10" - 13" Seat Widths
 - Pediatric Ki Mobility Backpack **New!** _____ \$65
 - Ki Mobility Backpack _____ \$65
 - Black Neoprene Impact Guards _____ \$85
 - Cane and Crutch Holder _____ \$195

Frame Color

- Gloss Finish**
- Black _____ N/C
 - Blue Skies _____ N/C
 - Bubblegum Pink _____ N/C
 - Candy Blue _____ N/C
 - Candy Purple _____ N/C
 - Candy Red _____ N/C
 - Copper _____ N/C
 - Fuchsia _____ N/C
 - Granny Smith _____ N/C
 - Yellow _____ N/C

- Matte Finish**
- Army Green _____ N/C
 - Charcoal Gray _____ N/C
 - Emerald Green _____ N/C
 - Iridescent Burgundy _____ N/C
 - Iridescent Navy _____ N/C
 - Retro Red _____ N/C
 - Royal Blue _____ N/C
 - Shadow Black _____ N/C
 - Smooth White _____ N/C

- Graphic Image Transfer**
- Visit www.kimobility.com for current graphic image transfer options**
Note selection below
- _____ \$595

Spark Front Frame Height Matrix																																																			
3" Caster				4" Caster				5" Caster				6" Caster				7" Caster				8" Caster																															
FSH	Super Low	Stem	Fork	FSH	Super Low	Stem	Fork	FSH	Super Low	Stem	Fork	FSH	Super Low	Stem	Fork	FSH	Super Low	Stem	Fork	FSH	Super Low	Stem	Fork																												
14.5	13.5	Std	4"	15 ¹	13.5 ¹	Std	4"	16	14.5	Std	4"	17	15.5	Std	5"	18	17	Std	6"	19	18	std	6"																												
15	14			15.5	14			16.5	15.5	3/4		18	16.5	3/4		18.5	17.5			20	18.5	3/4		6"																											
15.5	14	3/4		5"	16 ²	14.5 ²		3/4	5"	17.5	16	1 1/2	5"	18.5	17	1 1/2	6"	19		17.5	3/4	7"	20.5	n/a	1 1/2	7"																									
16	14.5				16.5	15				16	15	Std		17	16	Std		17.5		16.5			Std	19.5	18		19	17.5	19.5	18	std																				
16	15	1 1/2			6"	16.5 ²		15 ²		1 1/2	6"	16.5	15.5	3/4	6"	17.5	16.5	Std		7"	18		17	1 1/2	7"	20	19	Std	7"																						
16.5	15.5					17		15.5				17.5	16			18	17				19		18.5			17.5	20			19	20.5	19	19.5	18.5	20	19	std														
15.5	14.5	Std				7"		15.5		14.5		Std	7"	17.5		16	3/4	7"			18		16.5	3/4		7"	18.5	17		6"	7"	19	18	Std	7"																
16	15							16		15				17.5		16.5					1 1/2		18.5				17	19				17.5	19.5			18.5	20	19	19.5	18.5	20	19	3/4								
16.5	15	3/4						8"		16.5		15		3/4		8"	18				17		5"	8"			19	17.5		6"		8"	19.5	18.5		3/4	8"	20	19	Std	8"										
17	15.5									17		15.5					18.5				17.5						1 1/2	19.5					18.5	20				19	20.5			n/a	21	n/a	20.5	19	1 1/2				
17	16	1 1/2	9"				17			16		1 1/2		9"			17		16		Std		9"				19	18		1 1/2			9"	19		18		3/4	9"	21		n/a	Std	9"							
17.5	16.5						17.5			16.5							18.5		17.5		1 1/2						19.5	18.5						20		19				20.5		n/a			21	n/a	20.5	19	1 1/2		
16.5	15.5	Std		10"			17		15.5	Std		10"					17.5		16		3/4	10"					18	16.5		Std				10"		18		16.5		3/4		10"	20		19	3/4	10"				
17	16						17.5		16								18		16.5								3/4	18.5								17		19					17.5		20.5			n/a	21	n/a	20.5
17.5	16	3/4			11"		17.5		16.5	3/4	11"				18.5		17		6"	11"	19				17.5		Std	11"	19	17.5						1 1/2		11"		20			18.5		Std	11"					
18	16.5						18		17						19.5		18.5				1 1/2				20.5				19	21										n/a			20.5					19	21	n/a	20.5
18	17	1 1/2				12"	18.5		17	1 1/2			12"		18.5		17.5	1 1/2	12"		18.5				17.5	3/4	12"		18.5	17.5	7"				12"	21				n/a			Std		12"						
18.5	17.5						19		17.5						19		18				19				18				19.5	18.5						20				19								20.5	19	21	n/a
							18	17	Std	13"					18	16.5	Std	13"			19.5			18.5	3/4	13"			19.5	18	1 1/2	13"				8" x 2" caster not available with 6" fork															
							18.5	17.5							19	17.5					18.5			17					19.5	18.5										20	19		20.5					19	21	n/a	
							19	17.5	3/4					14"	19	17.5	7"				14"		19.5	18	1 1/2				14"	19.5	18		1 1/2							14"	6" x 2" caster not available with 5" fork										
							19.5	18							18.5	17.5							19	18						19.5	18.5																	20	18.5	20.5	19
				19.5			18.5	1 1/2	15"			19			18	3/4	15"					20.5	19	7"	15"					20.5	19		Std	15"													6" x 2" caster not available with 5" fork				
				20			19					19.5			18.5							20	18.5							20.5	19																				21
								19.5			18	1 1/2			16"	6" x 2" caster not available with 5" fork																																			
								20			18.5									20.5		19																													
								20.5			19	1 1/2	17"							6" x 2" caster not available with 5" fork																															
								20.5			19																																								

¹ Front Seat Height not available with 4 x 1" Poly Caster
² Not available with 4 x 1" Poly, 4 x 1.5" Poly or 4 x 1.5" Soft Roll Aluminum

Seat heights with "FSH" heading available with Standard Front Frame. Seat heights with "Super Low" heading available with Super Low Front Frame. Actual seat height may vary depending upon configuration.

Cushions

Matrix Seat Cushions

Cushion size is determined by chair dimensions. Certain sizes will have higher price and longer lead time; visit www.motionconcepts.com/matrix_international.html for size specifications.

STEP 1: Choose Cushion

- PS Cushion _____ \$510
Cushions 10" through 14" wide have a weight capacity of 150 lbs
- Kidabra Cushion _____ \$510
Weight capacity is 150 lbs

STEP 2: Choose Cushion Size

- Match to Chair Size
- Other Size - Provide Part Number _____

ROHO Seat Cushions

Cushion size is determined by chair dimensions. Certain sizes will have higher price and longer lead time; visit therohogroup.com for size specifications.

STEP 1: Choose Cushion

- HIGH PROFILE® Single Compartment (E2622) _____ \$603
- MID PROFILE™ Single Compartment (E2622) _____ \$603
- LOW PROFILE® Single Compartment (E2622) _____ \$603
- QUADTRO SELECT® High Profile® (E2624) _____ \$603
- QUADTRO SELECT® Mid Profile® (E2624) _____ \$603
- HYBRID ELITE® Single Comp Smart Check Ready E2622 _ \$623
- CONTOUR SELECT® (E2624) _____ \$603

STEP 2: Choose Cushion Size

- Match to Chair Size
- Other Size - Provide Part Number _____

Backs

<p>Matrix Backs</p> <p><input type="checkbox"/> Elite Back _____ \$750 <i>Available only with 14" through 16" Seat Widths</i></p> <p style="padding-left: 20px;">Back Height: <input type="checkbox"/> 10" <input type="checkbox"/> 12" <input type="checkbox"/> 14" <input type="checkbox"/> 16" <input type="checkbox"/> 18" <input type="checkbox"/> 20"</p> <p><input type="checkbox"/> Mini Back _____ \$750 <i>Not available with 16" Seat Width. Back Height of 14" not available with 10" and 11" Seat Widths. Comes standard with Kidabra cover.</i></p> <p style="padding-left: 20px;">Back Height: <input type="checkbox"/> 7" <input type="checkbox"/> 10" <input type="checkbox"/> 12" <input type="checkbox"/> 14"</p>	<p>Therafin Headrest for MATRIX Back Systems</p> <p><input type="checkbox"/> Headrest, Multi-Axis Offset Hardware with Smooth/Soft Pad 34187 _____ \$420</p> <p style="padding-left: 20px;">PAD SIZE: <input type="checkbox"/> 5 x 8 <input type="checkbox"/> 5 x 10 <input type="checkbox"/> 5 x 14</p>
	<p>Other Matrix Mounting Hardware</p> <p><i>Must provide part number</i></p> <p><input type="checkbox"/> _____ <i>For mounting hardware other than quick release</i></p>

ROHO Back Systems (Minimum chair width is 14")

STEP 1: Choose Back

<p><input type="checkbox"/> AGILITY™ Min. Contour w/ Foam Only (E2611) _____ \$549 BACKREST HEIGHT: <input type="checkbox"/> 10" <input type="checkbox"/> 13" <input type="checkbox"/> 16" <input type="checkbox"/> 18"</p> <p><input type="checkbox"/> AGILITY™ Min. Contour w/ Foam & Lumbar Pad (E2613) _____ \$659 BACKREST HEIGHT: <input type="checkbox"/> 10" <input type="checkbox"/> 13" <input type="checkbox"/> 16" <input type="checkbox"/> 18"</p> <p><input type="checkbox"/> AGILITY™ Min. Contour w/ Air Only (E2620) _____ \$769 BACKREST HEIGHT: <input type="checkbox"/> 10" <input type="checkbox"/> 13" <input type="checkbox"/> 16" <input type="checkbox"/> 18"</p> <p><input type="checkbox"/> AGILITY™ Min. Contour w/ Air & Lumbar Pad (E2615) _____ \$769 BACKREST HEIGHT: <input type="checkbox"/> 10" <input type="checkbox"/> 13" <input type="checkbox"/> 16" <input type="checkbox"/> 18"</p>	<p><input type="checkbox"/> AGILITY™ Mid Contour w/ Air Only (E2620) _____ \$781 BACKREST HEIGHT: <input type="checkbox"/> 10" <input type="checkbox"/> 13" <input type="checkbox"/> 16" <input type="checkbox"/> 18" <input type="checkbox"/> 20"</p> <p><input type="checkbox"/> AGILITY™ Mid Contour w/ Air & Lumbar Pad (E2615) _____ \$781 BACKREST HEIGHT: <input type="checkbox"/> 10" <input type="checkbox"/> 13" <input type="checkbox"/> 16" <input type="checkbox"/> 18" <input type="checkbox"/> 20"</p> <p><input type="checkbox"/> AGILITY™ Max Contour w/ Air Only (E2620) _____ \$920 BACKREST HEIGHT: <input type="checkbox"/> 16" <input type="checkbox"/> 18" <input type="checkbox"/> 20"</p> <p><input type="checkbox"/> AGILITY™ Max Contour w/ Air & Lumbar Pad (E2615) _____ \$920 BACKREST HEIGHT: <input type="checkbox"/> 16" <input type="checkbox"/> 18" <input type="checkbox"/> 20"</p>
---	---

STEP 2: Choose Back Hardware

Fixed _____ N/C
Not compatible with 18" and 20" tall backs

Quick Release _____ N/C

ROHO Laterals for AGILITY™ Back Systems

Sold in pairs

<p><input type="checkbox"/> Swing Away with Curved Pad (E1028/E0956) _____ \$1014</p> <p style="padding-left: 20px;">PAD SIZE: <input type="checkbox"/> Medium <input type="checkbox"/> Large <input type="checkbox"/> X-Large 4 T x 4 L 4 T x 5 L 5 T x 5 L</p> <p style="padding-left: 20px;">OFFSET TYPE: <input type="checkbox"/> Flush <input type="checkbox"/> Offset <input type="checkbox"/> 1" Offset <input type="checkbox"/> 2" Offset</p> <p><input type="checkbox"/> Fixed with Flat Pad (E0956) _____ \$452</p> <p style="padding-left: 20px;">PAD SIZE: <input type="checkbox"/> Medium <input type="checkbox"/> Large <input type="checkbox"/> X-Large 4 T x 4 L 4 T x 5 L 5 T x 5 L</p> <p style="padding-left: 20px;">OFFSET TYPE: <input type="checkbox"/> Flush <input type="checkbox"/> Offset</p>	<p>Therafin Headrest for AGILITY Back Systems</p> <p><input type="checkbox"/> Headrest, Multi-Axis Offset Hardware with Smooth/Soft Pad 34187 _____ \$420</p> <p style="padding-left: 20px;">PAD SIZE: <input type="checkbox"/> 5 x 8 <input type="checkbox"/> 5 x 10 <input type="checkbox"/> 5 x 14</p>
---	--



Vest

- Full Shape, Comfort Strap (E0960) _____ \$167
 - Small 30589 Medium 30590 Large 30591
 - X-Large 30592 2X-Large 30593
- Full Shape, Extended Strap (E0960) _____ \$167
 - Small 30594 Medium 30595 Large 30596
 - X-Large 30597 2X-Large 30598
- Trim Shape, Comfort Strap (E0960) _____ \$169
 - Small 30618 Medium 30619 Large 30620
 - X-Large 30621 2X-Large 30622
- Trim Shape, Extended Strap (E0960) _____ \$169
 - Small 30604 Medium 30624 Large 30625
 - X-Large 30626 2X-Large 30627

Shoulder Harness

With metal clips and cam locks

- Small (E0960) 61403 _____ \$157
- Medium (E0960) 61503 _____ \$164
- Large (E0960) 61703 _____ \$170
- X-Large (E0960) 61803 _____ \$177

Hip Belt

- 4 Point Y-Style, Metal Buckle (E0978) _____ \$173
 - Small 60024 Medium 60025 Large 60026
- Dual Rear Pull, Metal Attaching Clips (E0978)
 - 1" Straps, 6.25"W x 1.75"H Pads 60161 _____ \$97
 - 1.5" Straps, 7.25"W x 2.25"H Pads 60171 _____ \$106
 - 1.5" Straps, 9.25"W x 2.5"H Pads 60181 _____ \$106
- Single Pull, Metal Buckle (E0978)
 - 1" Straps, 4.25"W x 1.75"H Pads 61001 _____ \$86
 - 1" Straps, 6.25"W x 1.75"H Pads 61201 _____ \$86
 - 1.5" Straps, 5.25"W x 2.25"H Pads 62001 _____ \$95
 - 1.5" Straps, 9.25"W x 2.5"H Pads 62201 _____ \$95
 - 2" Straps, 9.25"W x 3"H Pads 64001 _____ \$102
 - 2" Straps, 11.25"W x 3"H Pads 64101 _____ \$102

Cane Mounted Swing Away Lateral Thoracic Support

Sold as pair

- Swing Away Kits w/ 4"x6" Gel Pads (E0956, E1028) 50021 _____ \$478
- Swing Away Kits, No Pads (E0956, E1028) 50022 _____ \$309

Knee Adductor

- Flip Down Hardware (E1028) _____ \$634
Sold in pairs, one right side and one left side
 - Short Arm, Fits 3/4", 7/8" and 1" Tubing 30248
 - Long Arm, Fits 3/4", 7/8" and 1" Tubing 30284
- Knee Adductor Pad with Neoprene Cover (E0956) _____ \$117
Sold single
 - Small Pad, 4.25" x 3.25" 50490 Quantity: _____
 - Medium Pad, 6.25" x 4.25" 50491 Quantity: _____

Foot / Shoe Holder

Compatible only with Aluminum Angle Adjustable and Composite Angle Adjustable Footplates

- FootSure Ankle Support (K0108) _____ \$178
Sold in pairs
 - X-Small w/ Buckle Closure, 6"-8" Circum. 60250
 - Small w/ Buckle Closure, 7"-8.5" Circum. 60251
 - Medium w/ Buckle Closure, 8"-11" Circum. 60252
 - Large w/ Buckle Closure, 9"-12" Circum. 60253
 - X-Large w/ Buckle Closure, 12.5"-15" Circum. 60254
 - X-Small w/ Hook-N-Loop Closure, 5"-6" Circum. 60255
 - Small w/ Hook-N-Loop Closure, 6"-8" Circum. 60256
 - Med w/ Hook-N-Loop Closure, 8"-10.5" Circum. 60257
 - Large w/ Hook-N-Loop Closure, 9"-12" Circum. 60258
 - X-Lrg w/ Hook-N-Loop Closure, 12"-15" Circum. 60259
- Padded Footbox (K0108)
 - Small, Full Design, Fits 14" Chair 30111 _____ \$287
 - Medium, Full Design, Fits 16" Chair 30112 _____ \$293
 - Large, Full Design, Fits 18" Chair 30113 _____ \$298
 - Small, Split Design, Fits 14" Chair 30152 _____ \$367
 - Medium, Split Design, Fits 16" Chair 30153 _____ \$373
 - Large, Split Design, Fits 18" Chair 30154 _____ \$378
- Molded Shoe Holders (K0108) _____ \$162
Sold in pairs
 - Small with Plain Straps, 3.5" x 5.25" 30293
 - Medium with Plain Straps, 4" x 6.5" 30303
 - Large with Plain Straps, 4.75" x 8" 30309
 - X-Large with Plain Straps, 5" x 10" 30316



Half Tray

Not compatible with Tubular Flip Up Armrests

- Standard Shape, Clear, Flip-Away (E0950, E1028) _____ \$246
 - Right Side 41310 Left Side 41311
- Standard Shape, Padded, Flip-Away (E0950, E1028) _____ \$269
 - Right Side 41350 Left Side 41353
- Trimline Shape, Padded, Flip-Away (E0950, E1028) _____ \$269
 - Right Side 41352 Left Side 41351

Clear Tray with Border and Pre-Drilled Holes

STEP 1: Choose Tray

Chair Width	Tray Dimension (W x D)	Body Cutout (W x D)	HCPCS	Part Number	Price	
10"	17 x 15	7.5 x 4.5, H2.75	E0950	33203	\$209	<input type="checkbox"/>
11"	17 x 15	7.5 x 4.5	E0950	33204	\$209	<input type="checkbox"/>
	17 x 15	8 x 5	E0950	33205	\$220	<input type="checkbox"/>
12"	18 x 16	8.5 x 5	E0950	33206	\$220	<input type="checkbox"/>
	18 x 16	9.5 x 5.5	E0950	33207	\$220	<input type="checkbox"/>
13"	19 x 18	9.5 x 5.5	E0950	33208	\$220	<input type="checkbox"/>
	19 x 18	10.5 x 6	E0950	33209	\$220	<input type="checkbox"/>
14"	20 x 18	10.5 x 6	E0950	33210	\$220	<input type="checkbox"/>
	20 x 18	11.5 x 6.5	E0950	33211	\$220	<input type="checkbox"/>
15"	21 x 18	11.5 x 6.5	E0950	33212	\$220	<input type="checkbox"/>
	21 x 18	12.5 x 7	E0950	33213	\$220	<input type="checkbox"/>
16"	22 x 21	12.5 x 7.5	E0950	33214	\$255	<input type="checkbox"/>
	22 x 21	13.5 x 8.5	E0950	33215	\$255	<input type="checkbox"/>

STEP 2: Choose Tray Hardware

- E-Z Lock Clamps, Tube Arm (E1028) 31072 _____ \$150
Available only with Tubular Flip Up and Angle Adjustable Locking Extendable Flip Up with Foam Grip Armrests. Sold in Pairs.
- E-Z Lock Clamps, Standard (E1028) 31098 _____ \$142
Not available with Tubular Flip Up and Angle Adjustable Locking Extendable Flip Up with Foam Grip Armrests. Sold in Pairs.

Communication Mount Hardware

Choose Device Holder AND Mount

STEP 1: What DEVICE will be used?

- Smart Phone 32359 _____ \$150
- Large Smart Phone 32354 _____ \$150
- I-Pad 32355 _____ \$184
- I-Pad Mini 32356 _____ \$178
- 7" Tablet 32357 _____ \$173
- 10" Tablet 32358 _____ \$395

STEP 2: Where does it need to MOUNT?

- Frame 32350 _____ \$512
- Armrest 32351 _____ \$395
Not compatible with Tubular Flip Up Armrests
- Tray 32352 _____ \$212