



Phone (800) 981-1540 | Fax (715) 254-0996
 sales@kimobility.com | www.kimobility.com

Company: _____

Ship To: _____

Account Number: _____ Date: _____

Address: _____

Quote PO Number: _____

City: _____ State: _____

RTS: _____

Zip: _____ Phone: _____

Marked For: _____

Email/Fax: _____

Effective 03/19/2018

Squiggles Seating System

Weight Capacity: 48.4 lbs

Seat Shell

Selection required

	Description	Item Code	HCPC	Retail
<input type="checkbox"/>	Squiggles Seat Shell	120-610		
	Seat Base with Depth Adjustment	120-610-1	E2294	\$924
	Seat Back - Fixed	120-610-2	E2293	\$924
	Seat Back Reclining Mechanism	120-610-3	E1014	\$404
	Seat Base Hip Guides / Support	120-610-4	E0956	\$112

Seat Cover

Selection required. Includes seat base and backrest cushion covers.

	Description	Item Code	HCPC	Retail
<input type="checkbox"/>	Squiggles Seat Covers - Blue	120-710-03		N/C
<input type="checkbox"/>	Squiggles Seat Covers - Green	120-710-01		N/C
<input type="checkbox"/>	Squiggles Seat Covers - Orange	120-710-02		N/C
<input type="checkbox"/>	Squiggles Seat Covers - Pink	120-710-04		N/C

Base

	Description	Item Code	HCPC	Retail
<input type="checkbox"/>	Size 1 Hi Low Chassis <i>Foot Pedal Operated. Includes Push Handle.</i>	120-1800	E1399	\$1895
<input type="checkbox"/>	Size 2 Hi Low Chassis <i>Foot Pedal Operated. Includes Push Handle.</i>	102-2753	E1399	\$1895
	Adapter to Mount Squiggles Seat on Size 2 Hi Low Chassis	120-818		\$188
<input type="checkbox"/>	Size 1 Kimba Neo - Black New!	470G7110000350	E1234	\$2205
	Transit System for Kimba	471X00ST007	E1399	\$51
<input type="checkbox"/>	Size 2 Kimba Neo - Black New!	470G7120000350	E1234	\$2205
	Transit System for Kimba	471X00ST007	E1399	\$51
<input type="checkbox"/>	Adapter to Mount Squiggles Seat on Size 2 Kimba Neo Base New!	471A71SK175	E1399	\$430

HCPCS codes provided are based on code verification by the PDAC or our interpretation of the code definitions.

Proper use of HCPCS codes does not ensure coverage or payment. For coverage information, verify the appropriate payer's coverage policy.

LECKEY^o

Foot Positioning

Sandals and Footplate Raisers require Squiggles Seat Adjustable Footrest

Description	Item Code	HCPC	Retail
<input type="checkbox"/> Squiggles Seat Adjustable Footrest	120-639	K0040	\$202
<input type="checkbox"/> Small Sandals with Straps (Pair)	152-1600	E0951, E0952	\$225
<input type="checkbox"/> Medium Sandals with Straps (Pair)	152-2600	E0951, E0952	\$245
<input type="checkbox"/> Footplate Raiser 1" (Each) Quantity: _____ (1-3)	120-816	E0951, E0952	\$76

Harness

Description	Item Code	HCPC	Retail
<input type="checkbox"/> Pelvic Harness Small - Blue	120-844-03	E0978	\$138
<input type="checkbox"/> Pelvic Harness Small - Green	120-844-01	E0978	\$138
<input type="checkbox"/> Pelvic Harness Small - Orange	120-844-02	E0978	\$138
<input type="checkbox"/> Pelvic Harness Small - Pink	120-844-04	E0978	\$138
<input type="checkbox"/> Pelvic Harness Medium - Blue	117-848-03	E0978	\$138
<input type="checkbox"/> Pelvic Harness Medium - Green	117-848-01	E0978	\$138
<input type="checkbox"/> Pelvic Harness Medium - Orange	117-848-02	E0978	\$138
<input type="checkbox"/> Pelvic Harness Medium - Pink	117-848-04	E0978	\$138
<input type="checkbox"/> Pelvic Harness Spacer Pads - Black (Pair)	120-789	K0108	\$78
<input type="checkbox"/> Small Trunk Harness	120-799	E0960	\$249
<input type="checkbox"/> Medium Trunk Harness	120-757	E0960	\$249

Lateral Support

Cover selection required with hardware

Description	Item Code	HCPC	Retail
<input type="checkbox"/> Flip-Away Lateral Support Hardware (Pair)	120-855	E1028	\$290
<input type="checkbox"/> Lateral Support Covers - Blue (Pair)	120-838-03	E0956	\$133
<input type="checkbox"/> Lateral Support Covers - Green (Pair)	120-838-01	E0956	\$133
<input type="checkbox"/> Lateral Support Covers - Orange (Pair)	120-838-02	E0956	\$133
<input type="checkbox"/> Lateral Support Covers - Pink (Pair)	120-838-04	E0956	\$133
<input type="checkbox"/> Rigid Lateral Support Hardware (Pair)	120-656	E1028	\$176
<input type="checkbox"/> Lateral Support Covers - Blue (Pair)	120-838-03	E0956	\$133
<input type="checkbox"/> Lateral Support Covers - Green (Pair)	120-838-01	E0956	\$133
<input type="checkbox"/> Lateral Support Covers - Orange (Pair)	120-838-02	E0956	\$133
<input type="checkbox"/> Lateral Support Covers - Pink (Pair)	120-838-04	E0956	\$133
<input type="checkbox"/> Laterals with Protraction Pads - Black (Pair)	120-954	E0956	\$372
<input type="checkbox"/> Chest Support - Blue <i>Requires laterals</i>	137-808-03	E0960	\$163
<input type="checkbox"/> Chest Support - Green <i>Requires laterals</i>	137-808-01	E0960	\$163
<input type="checkbox"/> Chest Support - Orange <i>Requires laterals</i>	137-808-02	E0960	\$163
<input type="checkbox"/> Chest Support - Pink <i>Requires laterals</i>	137-808-04	E0960	\$163

LECKEY^o

Tray

	Description	Item Code	HCPC	Retail
<input type="checkbox"/>	Tray - Gray	120-672	E0950	\$346
<input type="checkbox"/>	Extra High Tray + 2" Adjustment	120-688	E0950	\$437
<input type="checkbox"/>	Padded Tray Insert - Gray	120-270	K0108	\$280
<input type="checkbox"/>	Grab Rail/Toy Bar for Tray	117-769	K0108	\$251

Pommel

	Description	Item Code	HCPC	Retail
<input type="checkbox"/>	Pommel	120-641	E0957	\$98

Headrest

Cushion selection required with hardware

	Description	Item Code	HCPC	Retail
<input type="checkbox"/>	Contoured Headrest Support Hardware	120-859	E0128	\$317
<input type="checkbox"/>	<input type="checkbox"/> Contoured Headrest Cushion - Blue	137-800-03	E0955	\$145
	<input type="checkbox"/> Contoured Headrest Cushion - Green	137-800-01	E0955	\$145
	<input type="checkbox"/> Contoured Headrest Cushion - Orange	137-800-02	E0955	\$145
	<input type="checkbox"/> Contoured Headrest Cushion - Pink	137-800-04	E0955	\$145
<input type="checkbox"/>	Flat Headrest Support Hardware	117-697	E0128	\$194
<input type="checkbox"/>	<input type="checkbox"/> Flat Headrest Cushion - Blue	117-850-03	E0955	\$98
	<input type="checkbox"/> Flat Headrest Cushion - Green	117-850-01	E0955	\$98
	<input type="checkbox"/> Flat Headrest Cushion - Orange	117-850-02	E0955	\$98
	<input type="checkbox"/> Flat Headrest Cushion - Pink	117-850-04	E0955	\$98
<input type="checkbox"/>	Lateral Bracket Assembly Hardware (Pair) <i>Requires Flat Headrest Support Hardware and Cushion</i>	118-667	E1028	\$179
<input type="checkbox"/>	<input type="checkbox"/> Flat Headrest Lateral Cushions - Blue (Pair)	137-807-03	E0956	\$165
	<input type="checkbox"/> Flat Headrest Lateral Cushions - Green (Pair)	137-807-01	E0956	\$165
	<input type="checkbox"/> Flat Headrest Lateral Cushions - Orange (Pair)	137-807-02	E0956	\$165
	<input type="checkbox"/> Flat Headrest Lateral Cushions - Pink (Pair)	137-807-04	E0956	\$165
<input type="checkbox"/>	Removable Head Support Bracket for 1/2" Box Section	117-745	E1028	\$62

LECKEY^o

Accessories

Cover selection required with Shoulder Support Hardware

	Description	Item Code	HCPG	Retail
<input type="checkbox"/>	1" Back Spacer Pad and Flip-Away Lateral Supports - Blue	120-779-03	E0956, E1297	\$700
<input type="checkbox"/>	1" Back Spacer Pad and Flip-Away Lateral Supports - Green	120-779-01	E0956, E1297	\$700
<input type="checkbox"/>	1" Back Spacer Pad and Flip-Away Lateral Supports - Orange	120-779-02	E0956, E1297	\$700
<input type="checkbox"/>	1" Back Spacer Pad and Flip-Away Lateral Supports - Pink	120-779-04	E0956, E1297	\$700
<input type="checkbox"/>	Shoulder Support Hardware	120-702	K0108	\$237
<input type="checkbox"/>	Shoulder Support Cover - Blue	120-828-03	K0108	\$135
<input type="checkbox"/>	Shoulder Support Cover - Green	120-828-01	K0108	\$135
<input type="checkbox"/>	Shoulder Support Cover - Orange	120-828-02	K0108	\$135
<input type="checkbox"/>	Shoulder Support Cover - Pink	120-828-04	K0108	\$135
<input type="checkbox"/>	Long Sitting Extension	120-675	K0108	\$493

New! Kimba Neo Accessories

	Description	Item Code	HCPG	Retail
<input type="checkbox"/>	Kimba Neo Sun/Rain Canopy Complete Gray/Black <i>Requires Canopy Hardware</i>	471S71SK301	K0108	\$403
<input type="checkbox"/>	Kimba Neo Basic Canopy Black <i>Requires Canopy Hardware</i>	471S71SK505	K0108	\$175
<input type="checkbox"/>	Kimba Neo Size 1 Canopy Hardware <i>Must match base size. Requires Canopy.</i>	415S89SK023		\$84
<input type="checkbox"/>	Kimba Neo Size 2 Canopy Hardware <i>Must match base size. Requires Canopy.</i>	415S88SK024		\$84
<input type="checkbox"/>	Kimba Neo Vent Mount for Pulmonetics LTV	471C49NTKNLTVM	E1029/E1028	\$545
<input type="checkbox"/>	Kimba Neo Vent Mount for Respironics Trilogy	471C49NTKNWMTRI	E1029/E1028	\$545
<input type="checkbox"/>	Kimba Neo Size 1 Horizontal Vent Tray <i>Must match base size</i>	471S71SK060	E1029	\$464
<input type="checkbox"/>	Kimba Neo Size 2 Horizontal Vent Tray <i>Must match base size</i>	471S71SK160	E1029	\$464
<input type="checkbox"/>	Kimba Neo IV Pole	471C49NTKNIVP48	E0776	\$423
<input type="checkbox"/>	Kimba Neo O2 Tank Holder Side Mount	471C49NTKNO2H	E2208	\$261

Assembly

Product Assembly Fee not subject to discounts

<input type="checkbox"/>	Product Assembly Fee	\$120
--------------------------	----------------------	-------

LECKEY^o

Technical Specifications

Specification	Imperial	Metric
Age (Approximate)	1 - 5 years	
Stature	29.5" - 43.3"	750 mm - 1100 mm
Maximum User Weight	48.4 lbs	22 kg
Seat Depth	7.5" - 11.4"	190 mm - 290 mm
Seat Width	6.3" - 10.2"	160 mm - 260 mm
Hip Width	6.3" - 8.7"	160 mm - 220 mm
Backrest Height	12.6" - 15"	320 mm - 380 mm
Backrest Height with Shoulder Support	16.5" - 18.9"	420 mm - 480 mm
Minimum Chest Width	6.3"	160 mm
Maximum Knee Width	3.5"	90 mm
Minimum Distance Between Pommel and Side of Seat	7.5"	190 mm
Minimum Seat to Footplate	4.8"	122.5 mm
Seat to Floor	11.4" - 26.4"	290 mm - 670 mm
Backrest Angle (Prone / Recline)	10° / 25°	
Footplate Angle (Plantar / Dorsiflexion)	10° / 10°	
Tray Size	18.1" x 22"	460 mm x 560 mm
Seat Weight	13.2 lbs	6 kg
Tray Height	4.5" - 6.7"	115 mm x 170 mm

Squiggles Seat Hi-Lo Chassis Foot Operated (Excluding Push Bar)	Imperial	Metric
Weight	19.8 lbs	9 kg
Length	25.2"	640 mm
Width	21.3"	540 mm
Minimum Dimensions for Storage	12.4"	315 mm