



Phone (800) 981-1540 | Fax (715) 254-0996  
 sales@kimobility.com | www.kimobility.com

Effective 01/01/2018

Company: \_\_\_\_\_

Ship To: \_\_\_\_\_

Account Number: \_\_\_\_\_ Date: \_\_\_\_\_

Address: \_\_\_\_\_

Quote  PO Number: \_\_\_\_\_

City: \_\_\_\_\_ State: \_\_\_\_\_

RTS: \_\_\_\_\_

Zip: \_\_\_\_\_ Phone: \_\_\_\_\_

Marked For: \_\_\_\_\_

Email/Fax: \_\_\_\_\_

## Squiggles Stander

**HCPC: E0641**

*Weight Capacity: 48.4 lb*

### Standard Model

*Selection required. Includes hip and chest support with sternum prompt; chest and hip lateral supports; knee cups, footplate and tray attachment.*

	Description	Item Code	HCPC	Retail
<input type="checkbox"/>	Standard Model with Support Shell	118-600		\$1830
<input type="checkbox"/>	Pivot Base - Pneumatic Angle Adjustment for 3-in-1 Positioning	118-764		\$1337

### Cover

*Selection required. Includes chest and hip pad covers and knee support covers.*

	Description	Item Code	HCPC	Retail
<input type="checkbox"/>	Squiggles Stander Covers - Blue	118-762		\$517
<input type="checkbox"/>	Squiggles Stander Covers - Blue Vinyl	118-762-03VIN		\$517
<input type="checkbox"/>	Squiggles Stander Covers - Green	118-760		\$517
<input type="checkbox"/>	Squiggles Stander Covers - Orange	118-761		\$517
<input type="checkbox"/>	Squiggles Stander Covers - Pink	118-763		\$517

### Accessories

	Description	Item Code	HCPC	Retail
<input type="checkbox"/>	Size 1 Knee Cups	102-1640		\$107
<input type="checkbox"/>	Size 3 Knee Cups	102-3640		\$107
<input type="checkbox"/>	Tray - Gray	118-619		\$346
<input type="checkbox"/>	Padded Tray Insert - Gray	120-270		\$280
<input type="checkbox"/>	Grab Rail/Toy Bar for Tray	117-769		\$251
<input type="checkbox"/>	Sandals Small (Pair)	152-1600		\$225
<input type="checkbox"/>	Sandals Medium (Pair)	152-2600		\$245
<input type="checkbox"/>	Footplate Raiser 1"	120-816		\$76

*HCPCS codes provided are based on code verification by the PDAC or our interpretation of the code definitions.*

*Proper use of HCPCS codes does not ensure coverage or payment. For coverage information, verify the appropriate payer's coverage policy.*

# LECKEY<sup>o</sup>

## Headrest

*Cushion selection required with hardware*

	Description	Item Code	HCPC	Retail
<input type="checkbox"/>	Contoured Headrest Support Hardware	120-859	E0955	\$317
<input type="checkbox"/>	Contoured Headrest Cushion - Blue	137-800-03	E0955	\$145
<input type="checkbox"/>	Contoured Headrest Cushion - Blue Vinyl	137-800-03VIN	E0955	\$145
<input type="checkbox"/>	Contoured Headrest Cushion - Green	137-800-01	E0955	\$145
<input type="checkbox"/>	Contoured Headrest Cushion - Orange	137-800-02	E0955	\$145
<input type="checkbox"/>	Contoured Headrest Cushion - Pink	137-800-04	E0955	\$145
<input type="checkbox"/>	Flat Headrest Support Hardware	117-697	E0955	\$194
<input type="checkbox"/>	Flat Headrest Cushion - Blue	117-850-03	E0955	\$98
<input type="checkbox"/>	Flat Headrest Cushion - Blue Vinyl	117-850-03VIN	E0955	\$98
<input type="checkbox"/>	Flat Headrest Cushion - Green	117-850-01	E0955	\$98
<input type="checkbox"/>	Flat Headrest Cushion - Orange	117-850-02	E0955	\$98
<input type="checkbox"/>	Flat Headrest Cushion - Pink	117-850-04	E0955	\$98
<input type="checkbox"/>	Lateral Bracket Assembly Hardware	118-667	E1028	\$179
<input type="checkbox"/>	Lateral Headrest Cushion - Blue	137-807-03	E0956	\$165
<input type="checkbox"/>	Lateral Headrest Cushion - Blue Vinyl	137-807-03VIN	E0956	\$165
<input type="checkbox"/>	Lateral Headrest Cushion - Green	137-807-01	E0956	\$165
<input type="checkbox"/>	Lateral Headrest Cushion - Orange	137-807-02	E0956	\$165
<input type="checkbox"/>	Lateral Headrest Cushion - Pink	137-807-04	E0956	\$165

## Assembly

*Product Assembly Fee not subject to discounts*

<input type="checkbox"/>	Product Assembly Fee	\$100
--------------------------	----------------------	-------

# LECKEY<sup>o</sup>

## Technical Specifications

Specification	Imperial	Metric
Age (approximate)	1 - 5 years	
Maximum User Weight	48.4 lbs	22 kg
User Standing Height	29.5" - 43.5"	750 mm - 1110 mm
Chest Width Between Laterals	6.3" - 9"	160 mm - 230 mm
Chest Support Height	21.7" - 32"	550 mm - 820 mm
Chest Support Angle Adjustment	+/- 30°	
Chest Support Depth Adjustment	2"	50 mm
Hip Width Between Hip Guides	6.3" - 9"	160 mm - 230 mm
Hip Support Height from Footplate to Mid Buttocks	11.8" - 20.8"	300 mm - 530 mm
Footplate Angle	10° Plantarflexion / 10° Dorsiflexion	
Small Sandal Size	6.3" L x 3" W	175 x 75 mm
Medium Sandal Size	8.3" L x 3.2" W	210 x 80 mm
Tray Height: Prone	21.9" - 32.5"	555 mm - 825 mm
Tray Height: Supine	18.9" - 32.5"	480 mm - 825 mm
Tray Angle Adjustment	Prone 30° / Supine 40°	
Stander Support Frame Weight	15.4 lbs	7 kg
Stander Support Frame Dimensions	26.8" L x 15.4" W x 15.4" H	680 x 390 x 390 mm
Pivot Base Weight	12.1 lbs	5.5 kg
Pivot Base Dimensions (Unfolded)	32.3" L x 21.7" W x 20.3" H	820 x 550 x 515 mm
Pivot Base Dimensions (Folded)	27.6" L x 21.7" W x 11.8" H	700 x 550 x 300 mm
Knee Cup Opening Width	Size 1: 2.5"    Size 3: 3.5"	Size 1: 63 mm    Size 3: 89 mm