



Phone (800) 981-1540 | Fax (715) 254-0996  
 sales@kimobility.com | www.kimobility.com

Company: \_\_\_\_\_ Ship To: \_\_\_\_\_  
 Account Number: \_\_\_\_\_ Date: \_\_\_\_\_ Address: \_\_\_\_\_  
 Quote  PO Number: \_\_\_\_\_ City: \_\_\_\_\_ State: \_\_\_\_\_  
 RTS: \_\_\_\_\_ Zip: \_\_\_\_\_ Phone: \_\_\_\_\_  
 Marked For: \_\_\_\_\_ Email/Fax: \_\_\_\_\_

Effective 03/19/2018

## Squiggles Stander

HCPC: E0641

*Weight Capacity: 48.4 lb*

### New! Packages

*Package includes Squiggles Standard Model with Support Shell, Pivot Base, Positioning Pads, Flat Headrest Hardware, Cushion and Laterals, Tray, Small Sandals with Straps and Product Assembly*

	Description	Item Code	HCPC	Retail
<input type="checkbox"/>	Squiggles Stander Package - Blue	PSSTB	E0641	\$4984
<input type="checkbox"/>	Squiggles Stander Package - Blue Vinyl	PSSTBV	E0641	\$4984
<input type="checkbox"/>	Squiggles Stander Package - Green	PSSTG	E0641	\$4984
<input type="checkbox"/>	Squiggles Stander Package - Orange	PSSTO	E0641	\$4984
<input type="checkbox"/>	Squiggles Stander Package - Pink	PSSTP	E0641	\$4984

### Standard Model

*Selection required. Includes hip and chest support with sternum prompt; chest and hip lateral supports; size 2 knee cups, footplate and tray attachment.*

	Description	Item Code	HCPC	Retail
<input type="checkbox"/>	Standard Model with Support Shell	118-600		\$1830
<input type="checkbox"/>	Pivot Base - Pneumatic Angle Adjustment for 3-in-1 Positioning	118-764		\$1337

### Positioning Pads

*Selection required. Includes chest and hip pad covers and knee support covers.*

	Description	Item Code	HCPC	Retail
<input type="checkbox"/>	Squiggles Stander Positioning Pads - Blue	118-762		\$517
<input type="checkbox"/>	Squiggles Stander Positioning Pads - Blue Vinyl	118-762-03VIN		\$517
<input type="checkbox"/>	Squiggles Stander Positioning Pads - Green	118-760		\$517
<input type="checkbox"/>	Squiggles Stander Positioning Pads - Orange	118-761		\$517
<input type="checkbox"/>	Squiggles Stander Positioning Pads - Pink	118-763		\$517

HCPCS codes provided are based on code verification by the PDAC or our interpretation of the code definitions.

Proper use of HCPCS codes does not ensure coverage or payment. For coverage information, verify the appropriate payer's coverage policy.

# LECKEY<sup>o</sup>

## Accessories

Description				Item Code	HCPC	Retail
	Size 1 Knee Cup	<input type="checkbox"/> Quantity 1	<input type="checkbox"/> Quantity 2	102-1640		\$107
	Size 3 Knee Cup	<input type="checkbox"/> Quantity 1	<input type="checkbox"/> Quantity 2	102-3640		\$107
<input type="checkbox"/>	Tray - Gray			118-619		\$346
<input type="checkbox"/>	Padded Tray Insert - Gray			120-270		\$280
<input type="checkbox"/>	Grab Rail/Toy Bar for Tray			117-769		\$251
<input type="checkbox"/>	Small Sandals with Straps (Pair)			152-1600		\$225
<input type="checkbox"/>	Medium Sandals with Straps (Pair)			152-2600		\$245
<input type="checkbox"/>	Footplate Raiser 1" (Each)	Quantity: _____ (1-3)		120-816		\$76

## Headrest

**Cushion selection required with hardware**

Description				Item Code	HCPC	Retail
<input type="checkbox"/>	Contoured Headrest Support Hardware			120-859	E0128	\$317
<input type="checkbox"/>	Contoured Headrest Cushion - Blue			137-800-03	E0955	\$145
<input type="checkbox"/>	Contoured Headrest Cushion - Blue Vinyl			137-800-03VIN	E0955	\$145
<input type="checkbox"/>	Contoured Headrest Cushion - Green			137-800-01	E0955	\$145
<input type="checkbox"/>	Contoured Headrest Cushion - Orange			137-800-02	E0955	\$145
<input type="checkbox"/>	Contoured Headrest Cushion - Pink			137-800-04	E0955	\$145
<input type="checkbox"/>	Flat Headrest Support Hardware			117-697	E0128	\$194
<input type="checkbox"/>	Flat Headrest Cushion - Blue			117-850-03	E0955	\$98
<input type="checkbox"/>	Flat Headrest Cushion - Blue Vinyl			117-850-03VIN	E0955	\$98
<input type="checkbox"/>	Flat Headrest Cushion - Green			117-850-01	E0955	\$98
<input type="checkbox"/>	Flat Headrest Cushion - Orange			117-850-02	E0955	\$98
<input type="checkbox"/>	Flat Headrest Cushion - Pink			117-850-04	E0955	\$98
<input type="checkbox"/>	Lateral Bracket Assembly Hardware (Pair) <i>Requires Flat Headrest Support Hardware and Cushion</i>			118-667	E1028	\$179
<input type="checkbox"/>	Flat Headrest Lateral Cushions - Blue (Pair)			137-807-03	E0956	\$165
<input type="checkbox"/>	Flat Headrest Lateral Cushions - Blue Vinyl (Pair)			137-807-03VIN	E0956	\$165
<input type="checkbox"/>	Flat Headrest Lateral Cushions - Green (Pair)			137-807-01	E0956	\$165
<input type="checkbox"/>	Flat Headrest Lateral Cushions - Orange (Pair)			137-807-02	E0956	\$165
<input type="checkbox"/>	Flat Headrest Lateral Cushions - Pink (Pair)			137-807-04	E0956	\$165

## Assembly

**Product Assembly Fee not subject to discounts**

<input type="checkbox"/>	Product Assembly Fee	\$120
--------------------------	----------------------	-------

# LECKEY<sup>o</sup>

## Technical Specifications

Specification	Imperial	Metric
Age (approximate)	1 - 5 years	
Maximum User Weight	48.4 lbs	22 kg
User Standing Height	29.5" - 43.5"	750 mm - 1110 mm
Chest Width Between Laterals	6.3" - 9"	160 mm - 230 mm
Chest Support Height	21.7" - 32"	550 mm - 820 mm
Chest Support Angle Adjustment	+/- 30°	
Chest Support Depth Adjustment	2"	50 mm
Hip Width Between Hip Guides	6.3" - 9"	160 mm - 230 mm
Hip Support Height from Footplate to Mid Buttocks	11.8" - 20.8"	300 mm - 530 mm
Footplate Angle	10° Plantarflexion / 10° Dorsiflexion	
Small Sandal Size	6.3" L x 3" W	175 x 75 mm
Medium Sandal Size	8.3" L x 3.2" W	210 x 80 mm
Tray Height: Prone	21.9" - 32.5"	555 mm - 825 mm
Tray Height: Supine	18.9" - 32.5"	480 mm - 825 mm
Tray Angle Adjustment	Prone 30° / Supine 40°	
Stander Support Frame Weight	15.4 lbs	7 kg
Stander Support Frame Dimensions	26.8" L x 15.4" W x 15.4" H	680 x 390 x 390 mm
Pivot Base Weight	12.1 lbs	5.5 kg
Pivot Base Dimensions (Unfolded)	32.3" L x 21.7" W x 20.3" H	820 x 550 x 515 mm
Pivot Base Dimensions (Folded)	27.6" L x 21.7" W x 11.8" H	700 x 550 x 300 mm
Knee Cup Opening Width	Size 1: 2.5"    Size 3: 3.5"	Size 1: 63 mm    Size 3: 89 mm