

AXIOM

SEATING SOLUTIONS BY  Mobility




AIR | LTE

AXIOM

AIR | LTE

Axiom Air LTE, the lightweight and cost-effective air cushion option, is designed for general purpose applications. Delivering easy adjustment, comfort, and quality construction to suit a variety of care settings and environments.

Soft, durable and pliable materials create a comfortable and dependable solution for **general-use applications**.



Interconnected air cells adjust to a bodies shape and movement; improving **pressure distribution** and **reduced surface tension**.

Easy and **convenient** one-valve adjustment.

The large air cell design provides a **comfortable seating surface** for periods of prolonged sitting.

Specifications

Cover Type

Breathable Fabric

Height

3"

7.6 cm

Weight

1 lb

.5 kg

Sizes

16 x 16"
41 x 41 cm

16 x 18"
41 x 46 cm

18 x 16"
46 x 41 cm

18 x 18"
46 x 46 cm

