

CATALYST 5/4





INNOVATIVE AND PATENTED FRAME DESIGN

Patented frame concept uses intelligent geometry to create a significantly stronger, lighter and better performing chair.



VERTICAL AXLE PLATE

The Catalyst 5Vx, with all the performance you have come to expect from a Catalyst K5, offers three CG positions and five vertical seat height adjustments.



SWING AWAY FOOTRESTS

Enjoy the simplicity and strength of Catalyst's patented swing-away footrest hangers. The tapered engagement keeps the footrests snug and rattle free. Comes with a lifetime warranty.



The Catalyst 5Vx makes it a little easier to get reimbursement when CG adjustment is necessary.

SPECIFICATIONS

Transport Weight:	16 lbs.*	7.25 kg*
Seat Width:	14-22 in.	36-56 cm
Seat Depth:	14-20 in.	36-51 cm
Front Seat Height:	13.5-21 in.	34-53 cm
Rear Seat Height:	12.5-20 in.	32-51 cm
Top Frame Tubing Size:	1 in.	2.54 cm
Bottom Frame Tubing Size:	1.125 in.	2.86 cm
UPDATED! Weight Capacity:	300 lbs.	136 kg
Ansi Resna WC-19 Transit:	Yes	Yes

* 16"x16" Frame, 5"x1" Poly Casters, Height Adjustable T-Arm Receivers. Footrest and Rear Wheels removed.