

ROGUE²



FIND YOUR LIMIT... AND PUSH BEYOND.

The redefined style and superior design of Rogue² is the perfect balance of performance, functionality and self-expression. Developed to fit any lifestyle, every element of the chair has been enhanced to provide you with an optimal experience, wherever you go.



PERFECTLY FIT. AND PERFECTLY FITTING.

The lightest-in-class and engineered to be the perfect fit for your first chair or your next chair. Enhanced strength and durability paired with the right amount of adjustability lets you maximize your efficiency over the life of the wheelchair.



BACKREST DESIGN

Adjusting the backrest has never been easier, using one tool and no disassembly, change the backrest angle while the chair is occupied. Relax or adjust your clothing by simply engaging the relaxed position.



INDEXED DESIGN

Finally, a camber tube design that makes loading your chair easier. Effortlessly change center of gravity with the precisely aligned Dimple Track. Indexed rear seat height adjustment requires no disassembly.



WING DESIGN

Optimized wing design lets you get closer for transfers. Corrosion resistant components protect your chair from the elements. Easily get aligned with Ki's intuitive Caster Angle Adjustment System.



STYLE AND PERFORMANCE. REDEFINED.



“It is so important to have a chair that fits me and meets my needs. This is the first wheelchair that feels custom made for me. I can’t believe I’m saying this, but I feel so free in it! I’m so happy to be able to go about my day confidently.”

- Amanda Steijlen

*Founder,
Wheely Big Dreams*



“My passion for outdoor adventure takes me all over the world. Having a chair that can keep up and hold up in remote and rugged locations is essential.”

- Joe Stone

*Director of Mission,
Teton Adaptive Sports*



MAKE IT YOUR KIND OF... ONE OF A KIND.

Looking good and feeling good go hand in hand. The elegant, sculpted design will make you feel confident wherever you go. Minimalistic or edgy, make your chair one of a kind with the Color Visualizer.





SPECIFICATIONS

Transport Weight*	Seat Width
10.6 lbs. 4.8 kg	12-20" 30-51 cm
Camber	Tubing Size
0°, 2°, 4°, 6°, 8°	1 3/8" 3.5 cm
Seat Depth	Weight Capacity
14-20" 36-51 cm	275 lbs. 124.7 kg
Front Seat Height	Rear Seat Height
13.5-20.5" 34-52 cm	13.5-20" 34-51 cm
Ansi Resna WC-19 Transit Approved	

* 16"x16" Frame, Short Backrest Height, 4" Poly Casters, Tubular Open Footplate. Rear Wheels Removed.

Ki Mobility manufactures wheelchairs, seating and positioning products in the United States. Ki and Axiom products redefine their respective segments by raising the bar in innovation and quality. Ki Mobility continues to improve upon the experience for consumers and professionals by delivering on the "Better By Design" approach to products and services.

Ki Mobility has been a part of Etac since 2021.

Founded in 1973, Etac is a world-leading developer and provider of ergonomic assistive devices and patient handling equipment. Through specialized product brands, Etac offers state-of-the-art products for a wide range of daily needs and care settings, for people at all stages of life. Etac's product brands include Convoid, R82, Star, Molift, Immedia, Etac, Ki and Axiom.