



CATALYST 4 QUICK SHIP

Phone (800) 981-1540 | Fax (715) 254-0996
 sales@kimobility.com | www.kimobility.com

Company: _____ Ship To: _____
 Account Number: _____ Date: _____ Address: _____
 Quote PO Number: _____ City: _____ Province: _____
 RTS: _____ Postal Code: _____ Phone: _____
 Marked For: _____ Email/Fax: _____
 ADP Code: WMA2KI010 Effective 07/15/2020

Catalyst 4 \$1395
 300 lb weight capacity. Shadow Black frame color. Back upholstery not included.

<input type="checkbox"/> Configuration 1	<input type="checkbox"/> Configuration 2	<input type="checkbox"/> Configuration 3
Front Seat Height: 13.5" _____ N/C <i>Range: 13.5"-14"</i>	Front Seat Height: 15" _____ N/C <i>Range: 15"-16"</i>	Front Seat Height: 17" _____ N/C <i>Range: 17"-18.5"</i>
Rear Seat Height: 13.5" _____ N/C <i>Range: 13.5"-17.5"</i>	Rear Seat Height: 14" _____ N/C <i>Range: 14"-18"</i>	Rear Seat Height: 16" _____ N/C <i>Range: 16"-18.5"</i>
Frame Type: Super Low _____ N/C	Frame Type: Super Low _____ N/C	Frame Type: Standard _____ N/C
Width and Depth	Width and Depth	Width and Depth
<input type="checkbox"/> 16" x 16" _____ N/C	<input type="checkbox"/> 18" x 16" _____ N/C	<input type="checkbox"/> 16" x 18" _____ N/C
<input type="checkbox"/> 16" x 18" _____ N/C	<input type="checkbox"/> 18" x 18" _____ N/C	<input type="checkbox"/> 18" x 18" _____ N/C
<input type="checkbox"/> 18" x 18" _____ N/C	<input type="checkbox"/> 18" x 20" WAMB _____ \$225	<input type="checkbox"/> 18" x 20" WAMB _____ \$225
<input type="checkbox"/> 20" x 18" WAMA _____ \$295	<input type="checkbox"/> 20" x 20" WAMA, WAMB _____ \$520	<input type="checkbox"/> 20" x 18" WAMA _____ \$295
Wheels and Casters	Wheels and Casters	Wheels and Casters
Rear Wheel: 22" Mag _____ N/C	Rear Wheel: 22" Mag _____ N/C	Rear Wheel: 24" Mag _____ N/C
Tire: Low Poly 1" _____ N/C	Tire: Full Poly 1 3/8" _____ N/C	Tire: Full Poly 1 3/8" _____ N/C
Handrim: Aluminum Anodized _____ N/C	Handrim: Aluminum Anodized _____ N/C	Handrim: Aluminum Anodized _____ N/C
Caster: 4" x 1.5" Poly _____ N/C	Caster: 5" x 1" Poly _____ N/C	Caster: 6" x 1" Poly _____ N/C
Fork Size/Stem: 4"/Standard _____ N/C	Fork Size/Stem: 5"/Standard _____ N/C	Fork Size/Stem: 6"/Standard _____ N/C
Swing Away Footrest	Swing Away Footrest	Swing Away Footrest
Hanger: 60° Extension Mount _____ \$250	Hanger: 70° Extension Mount _____ \$250	Hanger: 70° Extension Mount _____ \$250
Footplate: Aluminum Angle Adj WAMJ _____ \$170	Footplate: Aluminum Angle Adj WAMJ _____ \$170	Footplate: Aluminum Angle Adj WAMJ _____ \$170
Heel Loop: Standard _____ N/C	Heel Loop: Standard _____ N/C	Heel Loop: Standard _____ N/C
Extension Tube: Short _____ N/C	Extension Tube: Medium _____ N/C	Extension Tube: Medium _____ N/C
Seat to Footrest Length: 11.5" _____ N/C <i>Range: 11.5"-15"</i>	Seat to Footrest Length: 13" _____ N/C <i>Range: 13"-17"</i>	Seat to Footrest Length: 14" _____ N/C <i>Range: 13"-17"</i>

All Configurations (Options on all configurations)

Backrest	Axle Plate	Other
Backrest: 8° Bend _____ N/C <i>Range: Tall, Set at 19.5"</i>	Axle Plate: Adjustable _____ N/C <i>Axle Position 3, 0° Camber</i>	Armrest: Height Adj Flip Back T-Arm _____ \$380 <i>Full Arm Pad</i>
Push Handle: Standard _____ N/C	Axle: Quick Release WAMU _____ N/C	Wheel Lock: Push to Lock w/ Extension _____ \$67
		Rear Anti-Tips _____ \$115

Quick Ship Seating, Parts and Accessories

All items selected below ship with the chair. Installation is required upon receipt.

		Axiom Posterior Back (PB)			
		Width			
		16"	18"	20"	
Height	Retail	\$911	\$911	\$911	
	ADP Code	SE0001624	SE0001624	SE0001624	
	Mid	16" <input type="checkbox"/> AXPB1616	18" <input type="checkbox"/> AXPB1816	20" <input type="checkbox"/> AXPB2016	
	18" <input type="checkbox"/> AXPB1618	18" <input type="checkbox"/> AXPB1818	20" <input type="checkbox"/> AXPB2018		

		Axiom Posterior Lateral Back (PL)			
		Width			
		16"	18"	20"	
Height	Retail	\$1045	\$1045	\$1045	
	ADP Code	SE0001626	SE0001626	SE0001626	
	Mid	16" <input type="checkbox"/> AXPL1616	18" <input type="checkbox"/> AXPL1816	20" <input type="checkbox"/> AXPL2016	
	18" <input type="checkbox"/> AXPL1618	18" <input type="checkbox"/> AXPL1818	20" <input type="checkbox"/> AXPL2018		

		Axiom S Cushion			
		Width			
		16"	18"	20"	
Depth	Retail	\$305	\$305	\$305	
	16"	<input type="checkbox"/> XS1616	<input type="checkbox"/> XS1816		
	17"	<input type="checkbox"/> XS1617	<input type="checkbox"/> XS1817		
	18"	<input type="checkbox"/> XS1618	<input type="checkbox"/> XS1818	<input type="checkbox"/> XS2018	
	19"	<input type="checkbox"/> XS1619	<input type="checkbox"/> XS1819	<input type="checkbox"/> XS2019	
	20"		<input type="checkbox"/> XS1820	<input type="checkbox"/> XS2020	
	21"			<input type="checkbox"/> XS2021	
	22"			<input type="checkbox"/> XS2022	

		Axiom P Cushion			
		Width			
		16"	18"	20"	
Depth	Retail	\$605	\$605	\$605	
	16"	<input type="checkbox"/> XP1616	<input type="checkbox"/> XP1816		
	17"	<input type="checkbox"/> XP1617	<input type="checkbox"/> XP1817		
	18"	<input type="checkbox"/> XP1618	<input type="checkbox"/> XP1818	<input type="checkbox"/> XP2018	
	19"	<input type="checkbox"/> XP1619	<input type="checkbox"/> XP1819	<input type="checkbox"/> XP2019	
	20"		<input type="checkbox"/> XP1820	<input type="checkbox"/> XP2020	
	21"			<input type="checkbox"/> XP2021	
	22"			<input type="checkbox"/> XP2022	

		Axiom SP Visco Cushion			
		Width			
		16"	18"	20"	
Depth	Retail	\$670	\$670	\$670	
	16"	<input type="checkbox"/> XSP1616V	<input type="checkbox"/> XSP1816V		
	17"	<input type="checkbox"/> XSP1617V	<input type="checkbox"/> XSP1817V		
	18"	<input type="checkbox"/> XSP1618V	<input type="checkbox"/> XSP1818V	<input type="checkbox"/> XSP2018V	
	19"	<input type="checkbox"/> XSP1619V	<input type="checkbox"/> XSP1819V	<input type="checkbox"/> XSP2019V	
	20"		<input type="checkbox"/> XSP1820V	<input type="checkbox"/> XSP2020V	
	21"			<input type="checkbox"/> XSP2021V	
	22"			<input type="checkbox"/> XSP2022V	

	Extra Outer Cushion Covers	Retail
<input type="checkbox"/>	Axiom S Extra Outer Cover	\$100
<input type="checkbox"/>	Axiom P Extra Outer Cover	\$100
<input type="checkbox"/>	Axiom SP Visco Extra Outer Cover	\$140

	Back Upholstery	Part Number	Retail
<input type="checkbox"/>	Standard 16" Wide Tall	CQS002	\$60
<input type="checkbox"/>	Standard 18" Wide Tall	CQS004	\$60
<input type="checkbox"/>	Standard 20" Wide Tall	CQS006	\$60
<input type="checkbox"/>	Tension Adj with Standard Cover 16" Wide (16.5"- 19.5")	CQS007	\$245
<input type="checkbox"/>	Tension Adj with Standard Cover 18" Wide (16.5"- 19.5")	CQS008	\$245
<input type="checkbox"/>	Tension Adj with Standard Cover 20" Wide (16.5"- 19.5")	CQS009	\$245

	Backrest Options	Part Number	Retail
<input type="checkbox"/>	Angle Adj with Std Push Handle (18"-20") (Pair)	CQS010	\$185

	Heel Loops	Part Number	Retail
<input type="checkbox"/>	Adjustable Heel Loops Short (14"-17") (Pair)	CQS023	\$45
<input type="checkbox"/>	Adjustable Heel Loops Medium (18"-20") (Pair)	CQS024	\$45

	Casters, Forks and Stems	Part Number	Retail
<input type="checkbox"/>	4" x 1.5" Poly Caster (Pair)	CQS029	\$63
<input type="checkbox"/>	5" x 1.5" Poly Caster (Pair)	CQS030	\$63
<input type="checkbox"/>	6" x 1.5" Poly Caster (Pair)	CQS031	\$63
<input type="checkbox"/>	7" x 1" Poly Caster (Pair)	CQS032	\$100
<input type="checkbox"/>	6" x 1" Pneumatic Caster (Pair)	CQS033	\$120
<input type="checkbox"/>	8" x 1" Pneumatic Caster (Pair)	CQS034	\$120
<input type="checkbox"/>	4" Standard Fork (Pair)	CQS035	\$60
<input type="checkbox"/>	5" Standard Fork (Pair)	CQS036	\$60
<input type="checkbox"/>	6" Standard Fork (Pair)	CQS037	\$60
<input type="checkbox"/>	7" Standard Fork (Pair)	CQS038	\$60
<input type="checkbox"/>	Standard Stem (Pair)	CQS039	\$20
<input type="checkbox"/>	3/4" Stem (Pair)	CQS040	\$20
<input type="checkbox"/>	1 1/2" Stem (Pair)	CQS041	\$20

	Handrims and Wheel Locks	Part Number	Retail
<input type="checkbox"/>	22" Plastic Coated Handrims (Pair)	CQS042	\$160
<input type="checkbox"/>	24" Plastic Coated Handrims (Pair)	CQS043	\$160
<input type="checkbox"/>	Pull to Lock (Pair)	CQS044	\$96
<input type="checkbox"/>	Pull to Lock with Extension Handles (Pair)	CQS045	\$180

	Accessories	Part Number	Retail
<input type="checkbox"/>	1 1/2" Auto Buckle 56"	CQS046	\$80
<input type="checkbox"/>	1 1/2" Auto Buckle 60"	CQS047	\$80
<input type="checkbox"/>	1 1/2" Padded Auto Buckle 56"	CQS048	\$175

Catalyst Front Seat Height Matrix																											
4" Caster				5" Caster				6" Caster				7" Caster				8" Caster											
FSH	Super Low	Stem	Fork	FSH	Super Low	Stem	Fork	FSH	Super Low	Stem	Fork	FSH	Super Low	Stem	Fork	FSH	Super Low	Stem	Fork								
15 ¹	13.5 ¹	Std	4"	16	14.5	Std	4"	17	15.5	Std	5"	18	17	Std	6"	19	18	Std	6"								
15.5	14			16.5	15.5	3/4		18	16.5	3/4		18.5	17.5			20	18.5	3/4		20.5	n/a	1 1/2					
16 ²	14.5 ²	3/4		17.5	16	1 1/2		18.5	17	1 1/2		19	17.5	3/4		19.5	18	19		17.5	Std						
16.5	15	1 1/2		16	15	Std	17	16	Std	19.5	18	1 1/2	19.5	18		20	18.5										
16.5 ²	15 ²			17	15.5		3/4	18		17	18.5		17.5	19.5		18.5	20	19	20.5	19							
15.5	14.5	Std		5"	17.5	16	3/4	5"	18	16.5	3/4	6"	18.5	17		Std	7"	19	18	3/4	7"						
16	15				17.5	16.5			1 1/2	18.5			17	19				17.5	19.5			18.5	20	19	20.5	19	
16.5	15	3/4			18	17			Std	18.5			17.5	1 1/2		19		18	3/4			19	18	20	19	20.5	n/a
17	15.5	1 1/2			16.5	15.5	17	16		19	18		19.5			18.5		20		19	20.5	19	21	n/a			
17	16	1 1/2			6"	17.5	16.5	Std	6"	18.5	17.5		1 1/2	7"		19		18	3/4	7"	19.5	18.5	1 1/2	7"			
17.5	16.5		17			16	19			18	20				19	20.5		n/a									
17	15.5	Std	17.5			16	3/4	18		16.5	Std		18.5		17	19.5		18.5	1 1/2		20	18.5	20.5		19	21	n/a
17.5	16		18			16.5		19		17.5			20		19	20.5		n/a									
17.5	16.5	3/4	18			17	1 1/2	18.5		17.5	3/4		19.5		18	1 1/2		20.5	19		1 1/2	21	n/a		8" x 2" caster not available with 6" fork		
18	17		18.5			17.5		19.5		18			20.5		19												
18.5	17	1 1/2	18.5	17.5		1 1/2	18.5	17.5		3/4	19	18	3/4		19	18	7"	19.5	18.5		1 1/2	7"					
19	17.5		18	16.5			19.5	18.5			20	19			20.5	19											
19.5	18	3/4	7"	19		17.5	Std	7"		19.5	18	1 1/2	7"		20	18.5	1 1/2	7"	20.5		19	1 1/2	7"				
19.5	18.5			18.5		17.5				19	18				20	18.5			20.5		19						
20	19	1 1/2		19.5	18.5	3/4	7"		20.5	19	21	n/a		6" x 2" caster not available with 5" fork													
				19.5	18	1 1/2																					
				20	18.5																						
				20.5	19																						

Wheel Size	Tire Profile	Rear Seat Height				
22"	Low	13.5	14.5	15.5	16.5	17.5
	Full	14	15	16	17	18
24"	Low	14.5	15.5	16.5	17.5	18.5
	Full	15	16	17	18	19

Maximum Seat Dump Reference (Front Seat Height minus Rear Seat Height)			
Seat Depth	16"	18"	20"
Maximum Dump	2"	2.5"	2.5"

¹ Front Seat Height not available with 4 x 1" Poly Caster
² Not available with 4 x 1" Poly, 4 x 1.5" Poly or 4 x 1.5" Soft Roll Aluminum

Seat heights with "FSH" heading available with Standard Front Frame. Seat heights with "Super Low" heading available with Super Low Front Frame. Actual seat height may vary depending upon configuration.