



Phone (800) 981-1540 | Fax (715) 254-0996
 sales@kimobility.com | www.kimobility.com

CATALYST 4C

AADL CAT # 151

Effective 10/01/2019

Company: _____

Ship To: _____

Account Number: _____ Date: _____

Address: _____

Quote PO Number: _____

City: _____ State: _____

RTS: _____

Zip: _____ Phone: _____

Marked For: _____

Email/Fax: _____

Catalyst 4C **\$1395**
 250 lb weight capacity

Frame

Seat Width / Seat Depth

Select box corresponding to desired Seat Width and Depth

		Seat Width									
		14"	15"	16"	17"	18"	19"	20"	21" ²	22" ²	
Price		N/C	N/C	N/C	N/C	N/C	\$188	\$395 ¹	N/C	N/C	
Seat Depth	14"	N/C	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	
	15"	\$80	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	
	16"	N/C	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	
	17"	\$113	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	
	18"	N/C	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	
	19"	\$244 ¹	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	
20"	\$395 ¹	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>		

¹ No charge with Heavy Duty Option
² Heavy Duty Option required

Weight Capacity Options

If no selection is made, 250 lb weight capacity will be used

- 300 lb Weight Capacity **New!** _____ N/C
Requires Adjustable Axle Plate
- 350 lb Weight Capacity _____ \$925
Requires minimum 6" caster height and Adjustable Axle Plate

Axle Options

- Offset Axle Bracket _____ Std
- Reverse Mount Offset Axle Bracket _____ N/C
Not available with Swing Away Armrests
- Adjustable Axle Plate **New!** _____ N/C
Provides center of gravity and finite Rear Seat Height adjustment. Comes standard with 300 lb Weight Capacity.

Cross Brace

- Standard Cross Braces _____ Std

Backrest

Back Post / Height Range / Back Height

Select Back Height in row corresponding to desired Back Post and Height Range

Back Post	Price	Height Range	Back Height			
Standard	N/C	Short	<input type="checkbox"/> 8.5	<input type="checkbox"/> 9.5	<input type="checkbox"/> 10.5	<input type="checkbox"/> 11.5
		Medium	<input type="checkbox"/> 12.5	<input type="checkbox"/> 13.5	<input type="checkbox"/> 14.5	<input type="checkbox"/> 15.5
		Tall	<input type="checkbox"/> 16.5	<input type="checkbox"/> 17.5	<input type="checkbox"/> 18.5	<input type="checkbox"/> 19.5
Angle Adjustable w/ Push Handle	\$231	Short	<input type="checkbox"/> 14	<input type="checkbox"/> 15	<input type="checkbox"/> 16	
		Medium	<input type="checkbox"/> 16	<input type="checkbox"/> 17	<input type="checkbox"/> 18	
		Tall	<input type="checkbox"/> 18	<input type="checkbox"/> 19	<input type="checkbox"/> 20	

Standard Back Post

Select option below only if a Standard Back Post option is selected on chart at left

- Straight with Push Handle _____ Std
- 8° Bend with Push Handle _____ N/C
Not available with Short Back Height Range

Back Upholstery

- Standard _____ Std
Not available with Short or Medium Back Height Ranges on 20" - 22" Seat Widths
- Tension Adjustable _____ \$175
- Omit Back Upholstery _____ -\$10

Footrest

Hanger Bend	
<input type="checkbox"/> 60°	\$195
<input type="checkbox"/> 70°	\$195
<input type="checkbox"/> 80°	\$195
<i>Not available with 7" or 8" Casters</i>	
<input type="checkbox"/> Elevating Leg Rests (ELR) (A2)	\$369
<input type="checkbox"/> Pro Elevating Leg Rests (PRO ELR) (A2)	\$500
<input type="checkbox"/> Omit Hangers	- \$25

Footplate	
<input type="checkbox"/> Composite	N/C
<i>Not available with PRO ELR</i>	
<input type="checkbox"/> Composite Angle Adjustable (A2)	\$195
<i>Not available with ELR or 14" Seat Width</i>	
<input type="checkbox"/> Aluminum Angle Adjustable (A2)	\$185

Footrest / Extension / Seat to Footrest Length
 Select Seat to Footrest Length in row corresponding to desired Footrest and Extension. Footrest Length of at least 2.5" shorter than Front Seat Height recommended. Extension Tube Hanger measurements are based on 70° Hanger / Front Frame Bend.

Footrest	Extension	Seat to Footrest Length									
		<input type="checkbox"/> 11	<input type="checkbox"/> 11.5	<input type="checkbox"/> 12	<input type="checkbox"/> 12.5	<input type="checkbox"/> 13	<input type="checkbox"/> 13.5	<input type="checkbox"/> 14	<input type="checkbox"/> 14.5	<input type="checkbox"/> 15	
Extension Tube Hangers (60°, 70°, 80°)	Short	<input type="checkbox"/> 11	<input type="checkbox"/> 11.5	<input type="checkbox"/> 12	<input type="checkbox"/> 12.5	<input type="checkbox"/> 13	<input type="checkbox"/> 13.5	<input type="checkbox"/> 14	<input type="checkbox"/> 14.5	<input type="checkbox"/> 15	
	Medium	<input type="checkbox"/> 13	<input type="checkbox"/> 13.5	<input type="checkbox"/> 14	<input type="checkbox"/> 14.5	<input type="checkbox"/> 15	<input type="checkbox"/> 15.5	<input type="checkbox"/> 16	<input type="checkbox"/> 16.5	<input type="checkbox"/> 17	
Elevating Leg Rests (ELR)	Short	<input type="checkbox"/> 12	<input type="checkbox"/> 13	<input type="checkbox"/> 14	<input type="checkbox"/> 15	<input type="checkbox"/> 16					
	Long	<input type="checkbox"/> 13	<input type="checkbox"/> 14	<input type="checkbox"/> 15	<input type="checkbox"/> 16	<input type="checkbox"/> 17	<input type="checkbox"/> 18				
Extended Length ELR	Short	<input type="checkbox"/> 18	<input type="checkbox"/> 19	<input type="checkbox"/> 20	<input type="checkbox"/> 21	<input type="checkbox"/> 22					
	Long	<input type="checkbox"/> 18	<input type="checkbox"/> 19	<input type="checkbox"/> 20	<input type="checkbox"/> 21	<input type="checkbox"/> 22	<input type="checkbox"/> 23	<input type="checkbox"/> 24			
PRO ELR Low Pivot	Short	<input type="checkbox"/> 12 - 17.5									
	Long	<input type="checkbox"/> 15 - 21									
PRO ELR High Pivot	Short	<input type="checkbox"/> 10.5 - 15.5									
	Long	<input type="checkbox"/> 13 - 19									

Heel Loop	
<input type="checkbox"/> Heel Loops	\$30
<input type="checkbox"/> Adjustable Heel Loops	\$45
<input type="checkbox"/> Omit Heel Loops	N/C

Calf Strap	
<input type="checkbox"/> Velcro Adjustable	\$44
<input type="checkbox"/> Padded Velcro Adjustable (A1)	\$85
<input type="checkbox"/> Bodypoint Padded AeroMesh (A1)	\$123

Residual Limb Support			
<input type="checkbox"/> Residual Limb Support (A2)			
<input type="checkbox"/> Left	\$450	<input type="checkbox"/> Right	\$450
<input type="checkbox"/> Omit Left Hanger	-\$14	<input type="checkbox"/> Omit Right Hanger	-\$14
Cover			
<input type="checkbox"/> Smooth	<input type="checkbox"/> Air Mesh	<input type="checkbox"/> Smooth	<input type="checkbox"/> Air Mesh
Pad Size			
<input type="checkbox"/> 5 x 8	<input type="checkbox"/> 9 x 14	<input type="checkbox"/> 5 x 8	<input type="checkbox"/> 9 x 14
<input type="checkbox"/> 7 x 10	<input type="checkbox"/> 11 x 10	<input type="checkbox"/> 7 x 10	<input type="checkbox"/> 11 x 10
<input type="checkbox"/> 9 x 10	<input type="checkbox"/> 11 x 14	<input type="checkbox"/> 9 x 10	<input type="checkbox"/> 11 x 14

Front Seat Height and Casters

Caster Size / Front Frame Type / Front Seat Height

Select Front Seat Height in row corresponding to desired caster size.

Note: Selections in first three Front Seat Height columns use Super Low Front Frame.

Caster Size	Front Seat Height													
	Uses Super Low Front Frame ³													
4"	<input type="checkbox"/> 13.5 ¹	<input type="checkbox"/> 14	<input type="checkbox"/> 14.5	<input type="checkbox"/> 15 ²	<input type="checkbox"/> 15.5	<input type="checkbox"/> 16	<input type="checkbox"/> 16.5	<input type="checkbox"/> 17	<input type="checkbox"/> 17.5	<input type="checkbox"/> 18	<input type="checkbox"/> 18.5	<input type="checkbox"/> 19	<input type="checkbox"/> 19.5	<input type="checkbox"/> 20
5"	<input type="checkbox"/> 14.5	<input type="checkbox"/> 15	<input type="checkbox"/> 15.5	<input type="checkbox"/> 16	<input type="checkbox"/> 16.5	<input type="checkbox"/> 17	<input type="checkbox"/> 17.5	<input type="checkbox"/> 18	<input type="checkbox"/> 18.5	<input type="checkbox"/> 19	<input type="checkbox"/> 19.5	<input type="checkbox"/> 20	<input type="checkbox"/> 20.5	
6"	<input type="checkbox"/> 15.5	<input type="checkbox"/> 16	<input type="checkbox"/> 16.5	<input type="checkbox"/> 17	<input type="checkbox"/> 17.5	<input type="checkbox"/> 18	<input type="checkbox"/> 18.5	<input type="checkbox"/> 19	<input type="checkbox"/> 19.5	<input type="checkbox"/> 20	<input type="checkbox"/> 20.5	<input type="checkbox"/> 21		
7"	<input type="checkbox"/> 17	<input type="checkbox"/> 17.5		<input type="checkbox"/> 18	<input type="checkbox"/> 18.5	<input type="checkbox"/> 19	<input type="checkbox"/> 19.5	<input type="checkbox"/> 20	<input type="checkbox"/> 20.5	<input type="checkbox"/> 21				
8"	<input type="checkbox"/> 17.5	<input type="checkbox"/> 18	<input type="checkbox"/> 18.5	<input type="checkbox"/> 19	<input type="checkbox"/> 19.5	<input type="checkbox"/> 20	<input type="checkbox"/> 20.5	<input type="checkbox"/> 21						

¹ Front Seat Height not available with 4 x 1" Poly Caster

² Super Low Front Frame required when using 4 x 1" Poly Caster

³ Super Low Front Frame \$100 retail

Caster Type

- 1" Poly _____ N/C
- 1" Poly Aluminum (A1) _____ \$100
Not available with 6", 7" and 8" Caster Sizes
- 1" Pneumatic _____ \$100
Not available with 4", 5" and 7" Caster Sizes
- 1.5" Poly _____ \$63
Not available with 7" and 8" Caster Sizes
- 1.5" Soft Roll Aluminum (A1) _____ \$150
Not available with 7" and 8" Caster Sizes
- 2" Poly _____ \$88
Not available with 4", 5", 7" and 8" Caster Sizes
- 2" Pneumatic w/ Foam Insert (A1) _____ \$285
Not available with 4", 5", 6" and 7" Caster Sizes

Maximum Seat Dump Reference

(Front Seat Height minus Rear Seat Height)

Seat Depth	14"	15"	16"	17"	18"	19"	20"
Maximum Dump	1.5"	2"	2"	2"	2.5"	2.5"	2.5"

Fork Stem

- Standard Fork Stem _____ Std
- 3/4" Longer Fork Stem _____ \$44
- 1 1/2" Longer Fork Stem _____ \$44

Wheel Locks

Wheel Lock

- Push to Lock _____ Std
- Push to Lock Flush Mount _____ N/C
- Pull to Lock _____ N/C
- Omit Wheel Locks _____ N/C

Wheel Lock Options

- Extension Handles _____ \$84
Available only with Push to Lock or Pull to Lock Wheel Locks
- Grade Aids (A2) _____ \$187
Available only with Push to Lock Wheel Locks and Pneumatic or Pneumatic with Airless Insert Rear Tire Types

Rear Seat Height and Wheels

Wheel Size / Tire Profile / Rear Seat Height

Select Rear Seat Height in row corresponding to desired Wheel Size and Tire Profile. 14" Seat Depth not available with 24" or 26" Rear Wheels.

Wheel Size	Tire Profile	Rear Seat Height				
20"	Low	<input type="checkbox"/> 12.5 ¹	<input type="checkbox"/> 13.5	<input type="checkbox"/> 14.5	<input type="checkbox"/> 15.5	<input type="checkbox"/> 16.5
	Full	<input type="checkbox"/> 13 ¹	<input type="checkbox"/> 14	<input type="checkbox"/> 15	<input type="checkbox"/> 16	<input type="checkbox"/> 17
22"	Low	<input type="checkbox"/> 13.5 ¹	<input type="checkbox"/> 14.5	<input type="checkbox"/> 15.5	<input type="checkbox"/> 16.5	<input type="checkbox"/> 17.5
	Full	<input type="checkbox"/> 14 ¹	<input type="checkbox"/> 15	<input type="checkbox"/> 16	<input type="checkbox"/> 17	<input type="checkbox"/> 18
24"	Low	<input type="checkbox"/> 14.5 ¹	<input type="checkbox"/> 15.5	<input type="checkbox"/> 16.5	<input type="checkbox"/> 17.5	<input type="checkbox"/> 18.5
	Full	<input type="checkbox"/> 15 ¹	<input type="checkbox"/> 16	<input type="checkbox"/> 17	<input type="checkbox"/> 18	<input type="checkbox"/> 19
26"	Low	<input type="checkbox"/> 15.5 ¹	<input type="checkbox"/> 16.5	<input type="checkbox"/> 17.5	<input type="checkbox"/> 18.5	<input type="checkbox"/> 19.5
	Full	<input type="checkbox"/> 16 ¹	<input type="checkbox"/> 17	<input type="checkbox"/> 18	<input type="checkbox"/> 19	<input type="checkbox"/> 20

¹ Not available with Flip Back Armrests

Rear Wheel Type

- Spoke _____ N/C
 Mag _____ N/C

Rear Tire Type

Full Profile

- Street Pneumatic 1 3/8" _____ Std
 Full Poly 1 3/8" _____ N/C
 Pneumatic with Airless Insert 1 3/8" (A1) _____ \$200

Low Profile

- Low Poly 1" _____ N/C
Not available on 26" rear wheels
 Sentinel Solid Tire 1" _____ N/C
Not available on 20" and 22" Mag wheels
 High Pressure 1" (A1) _____ \$119
Not available on Mag wheels
 Schwalbe Marathon Plus Evolution 1" (A1) _____ \$145
Not available on Mag wheel, 20" or 22" wheel sizes

Handrim

Standard Handrims

- Aluminum Anodized _____ Std
 Plastic Coated _____ \$119
 Projection (8 Vertical Projections) _____ \$299
Available on 22", 24", and 26" Spoke wheels & 24" Mag wheels w/full tires
 Omit Handrim _____ N/C

Handrim Option

- 3M Micro Grip (A1) _____ \$81
Not available on Plastic Coated or Projection Handrims

Axles

- Axle**
- Fixed Axle _____ Std
 - Quick Release _____ N/C
 - Quad Quick Release (A2) _____ \$219
Not available with Mag wheels

Armrests

- Armrest** Please select Arm Pad option
- Height Adjustable Flip Back T-Arm (9.5" - 13.5") _____ \$295

- Arm Pad**
- Not available with Swing Away Armrests*
- Standard Arm Pad _____ N/C Desk Full
 - Waterfall Arm Pad (A1) _____ \$35 Desk Full
Available only with Height Adjustable T-Arms

- Armrest** continued
- Swing Away (10.5" - 12.5") (A1) _____ \$175
Not available with Angle Adjustable Back Posts or Arm Pads
 - Omit Armrests _____ N/C

Accessories

- Removable Side Guard**
- Composite Side Guards (A1) _____ \$175
Not available with T-Arm or Flip Back Armrests

- Anti-Tip**
- Rear Anti-Tips _____ \$110

- Positioning Belt**
- 1 1/2" Auto Buckle _____ \$22

- Seat Pan**
- Removable Aluminum Seat Pan (A1) _____ \$245
Pan depth 1/2" shorter than Seat Depth

- Other Accessories**
- Tool Kit (A1) _____ \$60
 - Removable Underseat Pouch (A1) _____ \$75
 - Ki Mobility Backpack (A1) _____ \$90
 - Black Neoprene Impact Guards (A1) _____ \$95
 - Clear Spoke Guards (A1) _____ \$144
Not available with Mag wheels
 - Black Spoke Guards (A1) _____ \$144
Not available with 20" rear wheel or Mag wheels
 - Cane and Crutch Holder (A1) _____ \$110
 - IV Pole (A1) _____ \$175
Available only with Standard Straight Back Posts

- O₂ Holder**
- O₂ Holder _____ \$110
 Right Side Left Side

Frame Color

- Gloss Finish**
- Black _____ N/C
 - Blue _____ N/C

Catalyst Front Seat Height Matrix																														
4" Caster				5" Caster				6" Caster				7" Caster				8" Caster														
FSH	Super Low	Stem	Fork	FSH	Super Low	Stem	Fork	FSH	Super Low	Stem	Fork	FSH	Super Low	Stem	Fork	FSH	Super Low	Stem	Fork											
15 ¹	13.5 ¹	Std	4"	16	14.5	Std	4"	17	15.5	Std	5"	18	17	Std	6"	19	18	Std	6"	6"										
15.5	14			16.5	15.5	3/4		18	16.5	3/4		18.5	17.5			20	18.5	3/4												
16 ²	14.5 ²	17.5		16	1 1/2	18.5		17	1 1/2	19		17.5	20.5	n/a		1 1/2														
16.5	15	3/4		5"	16	15	Std	5"	17	16	Std	6"	19.5	18		3/4	7"	19	17.5		Std	7"	7"							
16.5 ²	15 ²				16.5	15.5			17	15.5			18	17				19.5	18.5					19.5	18	20	18.5			
17	15.5	1 1/2			17.5	16	3/4		18	16.5			18.5	17.5		20.5		19	20.5					19						
15.5	14.5	Std			5"	17.5	16	1 1/2	5"	18.5	17	3/4	6"	19		18		Std	7"		19.5	18.5		3/4	7"	7"				
16	15					18	17			19	17.5			20		19														
16.5	15	3/4				16.5	15.5			Std	18.5			17.5		19		18			19.5	18.5					20.5	n/a		
17	15.5	3/4				5"	17	16	Std	5"	19	18	1 1/2	6"		19		18			3/4	7"		20.5	19		1 1/2	7"	7"	
17.5	16.5		17.5				16.5	19			18	20.5			19															
17	16	1 1/2	17				16	Std	19		18	19.5			18.5	21		n/a		21	n/a									
17	15.5	Std	6"				17.5	16	3/4	6"	18	16.5	Std	7"	20	19		3/4		7"	8" x 2" caster not available with 6" fork				7"					
17.5	16			18			16.5	18.5			17	20.5			n/a															
17.5	16.5	3/4		18.5			17	19			17.5	20			18.5															
18	17	3/4		6"			18	17	1 1/2	6"	19.5	18.5	3/4	7"	20.5	19	1 1/2	7"									8" x 2" caster not available with 6" fork			
18.5	17				18.5		17.5	18.5			17.5	21			n/a															
19	17.5	1 1/2			19		18	19			18	20.5			19															
18	17	Std			7"		18	16.5	Std	7"	19.5	18.5	3/4	7"	8" x 2" caster not available with 6" fork				7"											
18.5	17.5					18.5	17	20			19																			
19	17.5	19				17.5	19.5	18																						
19.5	18	3/4				7"	18.5	17.5	3/4	7"	20	18.5	1 1/2	7"												8" x 2" caster not available with 6" fork				
19.5	18.5		19				18	20.5			19																			
20	19	1 1/2	19.5				18.5	21			n/a																			
							19.5	18	1 1/2	7"	6" x 2" caster not available with 5" fork									7"										
				20			18.5																							
				20.5			19																							

¹ Front Seat Height not available with 4 x 1" Poly Caster

² Not available with 4 x 1" Poly, 4 x 1.5" Poly or 4 x 1.5" Soft Roll Aluminum

Seat heights with "FSH" heading available with Standard Front Frame. Seat heights with "Super Low" heading available with Super Low Front Frame. Actual seat height may vary depending upon configuration.