



## Backrest

### Backrest

Select Back Height in row corresponding to desired Back Post and Height Range. Depth Adjustable not available with Titanium Frame.

Back Post	Price	Height Range	Back Height				
Standard [50BOP002]	N/C	Short <sup>3,4</sup>	8.5"	9.5"	10.5"	11.5"	
		Medium <sup>4</sup>	12.5" <sup>6</sup>	13.5"	14.5"	15.5"	
		Tall	16.5"	17.5"	18.5"	19.5"	
Angle Adjustable w/ Push Handle [50BOP004]	N/C	Short <sup>4</sup>	14"	15"	16"		
		Medium <sup>4</sup>	16"	17"	18"		
		Tall	18"	19"	20"		
Stroller Handle Back Post (Depth Adjustable) [50BOP032]	N/C	Short	20"				
		Tall	24"				
Height Adjustable Straight with Push Handle (Depth Adjustable) [50BOP033]	£73	Short <sup>1</sup>	13"	14"	15"	16"	17"
		Medium	16"	17"	18"	19"	20"
		Tall	20" <sup>8</sup>	21" <sup>5</sup>	22" <sup>5</sup>	23" <sup>5</sup>	24" <sup>5</sup>
Reclining w/ Adjustable Handle [50BOP028] <sup>1,7</sup>	N/C	Short	22"				
		Tall	26"				
Half-Folding [50BOP029] <sup>2,8</sup>	£83	Short	12.5" <sup>6,9</sup>	13.5"	14.5"	15.5"	
		Tall	16.5"	17.5"	18.5"	19.5"	

1 Not available with Dynamic Back

2 Not available with Standard Back Upholstery

3 Not available with Tension Adjustable Back Upholstery

4 Not available with Standard Back Upholstery in combination with 21" and 22" Seat Widths

5 Not available with Tension Adjustable Back Upholstery in combination with 14" and 15" Seat Widths

6 Not available with Tension Adjustable Back Upholstery in combination with 22" Seat Width

7 Not available with O2 Holder, Transit Option, Heavy Duty Option, Swing Away Armrests, Flip Back Armrest, Height Adjustable Flip Back T-Arm or IV Pole

8 Not available with Transit Option, Heavy Duty Option, Flip Back Armrest, Angle Adjustable Locking Extendable Flip Up Armrest or Fold Down Push Handle

9 Not available with 8° Bend with Push Handle Back Post

10 Not available with Angle Adjustable Locking Extendable Flip Up Armrest

11 Not available with Tension Adjustable Back Upholstery in combination with Angle Adjustable Locking Extendable Flip Up Armrest

### Reclining Back Options

Select options below if Reclining Back is selected on chart above. Recline angle range is 0° - 65°.

Standard Spring [50GSP001]	N/C
Heavy Duty Spring [50GSP002]	N/C

### Standard Back Post

Select option below only if a Standard Back Post option is selected on chart above

Straight with Push Handle [50BCK004]	N/C
Straight without Push Handle [50BCK008]	N/C
8 Degree Bend with Push Handle [50BOP003]	£37
<i>Not available with Short Back Height Range</i>	
Straight with Fold Down Push Handle [50BCK006]	£158
<i>Not available with 350 lb Weight Capacity</i>	

### Push Handle Options

Not available with Stroller Handle Back Post, Straight Back Post with Fold Down Push Handle, Reclining Back Post or 350 lb Weight Capacity

Low Profile Push Handle [50BCK007]	£67
<i>2" shorter than standard push handles</i>	

### Depth Adjustable Back Post Options

Options below may be selected only if a Depth Adjustable back post option is selected on the chart above

Stroller Handle Extensions [50BOP012]	£73
<i>Not available with Height Adjustable Straight with Push Handle Back Post</i>	
Dynamic Back (E2398) [50BOP014]	£1932
<i>Available only with T-Arms</i>	

### Back Upholstery Options

Not available with Reclining or Stroller Handle Back Post. If selecting both upholstery and Axiom Back, the upholstery will be installed.

Standard [50BUP001]	N/C
Tension Adjustable with Standard Cover [50BUP002]	N/C
Tension Adjustable with Air Mesh Breathable Cover [50BUP007]	£182
Omit [50BUP003]	N/C

## Frame Type / Footrest

Frame Type		PRO ELR Covers
Swing Away [50FRM001]	N/C	<i>Available only on PRO ELR</i> Gel Pad Covers (Pair) [50ELRPAD] £304
Super Low - Swing Away [50FRM002]	N/C	
70° Fixed Front [50FRM003]	N/C	
80° Fixed Front [50FRM004]	N/C	

### Footrest Hangers

Frame	Hanger	60°	70°	80°	90°	Retail
Swing Away	Extension Tube Mount	[50HAN001] <sup>2</sup>	[50HAN002]	[50HAN003] <sup>1</sup>		N/C
	Front Mount <sup>2</sup>		[50HAN016]	[50HAN017] <sup>1</sup>	[50HAN004] <sup>1</sup>	N/C
	Pro Elevating Leg Rests				[50HAN044]	£483
	Omit Hangers				[50HAN046]	N/C
Fixed	Extension Tube Mount		[50FRM003]	[50FRM004] <sup>1</sup>		N/C

<sup>1</sup> Not available with 7" or 8" casters  
<sup>2</sup> Not available with Titanium Frame

### Seat to Footrest Length

Select Extension and Seat to Footrest Length in row corresponding to desired Footrest. Footrest Length of at least 2.5" shorter than Front Seat Height recommended. Extension Tube Hangers and Fixed Front Frame measurements are based on 70° Hanger / Front Frame Bend.

Footrest	Extension	Seat to Footrest Length												
		11"	11.5"	12"	12.5"	13"	13.5"	14"	14.5"	15"	15.5"			
Extension Tube Hangers (60°,70°,80°)	Short (11" - 15")	11"	11.5"	12"	12.5"	13"	13.5"	14"	14.5"	15"				
	Medium (13" - 17")	13"	13.5"	14"	14.5"	15"	15.5"	16"	16.5"	17"				
	Long (15.5" - 19.5")	15.5"	16"	16.5"	17"	17.5"	18"	18.5"	19"	19.5"				
Front Mount Hangers <sup>1</sup>	Short Front Mount Hanger	Standard Selections										Requires 3" Extension		
		4" <sup>2</sup>	4.5" <sup>2</sup>	5"	5.5"	6"	6.5"					10"	10.5"	11"
	Medium Front Mount Hanger	7"	7.5"	8"	8.5"	9"	9.5"					11.5"	12"	12.5"
		4" <sup>2</sup>	4.5" <sup>2</sup>	5"	5.5"	6"	6.5"	7"	7.5"			12"	12.5"	13"
	Long Front Mount Hanger	8"	8.5"	9"	9.5"	10"	10.5"	11"	11.5"			13.5"	14"	14.5"
		5"	5.5"	6"	6.5"	7"	7.5"	8"	8.5"	9"		14"	14.5"	15"
		9.5"	10"	10.5"	11"	11.5"	12"	12.5"	13"	13.5"	15.5"	16"	16.5"	
PRO ELR Low Pivot Setting	PRO ELR Low Pivot Setting Short (12" - 17.5")	12 - 17.5												
	PRO ELR Low Pivot Setting Long (15" - 21")	15 - 21												
PRO ELR High Pivot Setting	PRO ELR High Pivot Setting Short (10.5" - 15.5")	10.5 - 15.5												
	PRO ELR High Pivot Setting Long (13" - 19")	13 - 19												
Fixed Front Frame (70°,80°)	Short (11" - 15")	11.5"	12"	12.5"	13"	13.5"	14"	14.5"	15"					
	Medium (13" - 17")	13"	13.5"	14"	14.5"	15"	15.5"	16"	16.5"	17"				
	Long (15.5" - 19.5")	15.5"	16"	16.5"	17"	17.5"	18"	18.5"	19"	19.5"				

<sup>1</sup> Not available with Titanium Frame. See Minimum Front Seat Height chart below  
<sup>2</sup> Available with 90° only

### 3" Extension

Available only with Front Mount Hangers

3" Extension [50EXT008]	£37
3" Extension (Provided as Parts) [50EXT018]	£37

### Footplates

Composite [50FPL001]	N/C
<i>Not available with Front Mount Hanger</i>	
Composite Angle Adjustable [50FPL014]	N/C
<i>Not available with Front Mount Hanger or 14" Seat Width</i>	
Aluminum Angle Adjustable [50FPL002]	£115
Aluminum Locking Multi-Angle Adjustable [50FPL020]	£290
<i>Available only with Front Mount Hanger</i>	
One Piece Flip Up Angle Adjustable [50FPL009]	£272
<i>Not available with Front Mount Hanger or PRO ELR</i>	

Hanger Length	Short		Medium		Long	
	70°	80/90°	70°	80/90°	70°	80°
Minimum Front Seat Height with Standard Selection	13.5"	13.5"	14"	15"	16"	17"
Minimum Front Seat Height with 3" Extension	14.5"	16"	16.5"	18"	18.5"	20"

### Hanger Release

If no selection is made Classic Hanger Release will be used

Classic [50REL001]	N/C
4-Way - In-Line [50REL002]	N/C
4-Way - Out [50REL003]	N/C

## Heel Loops

*Not available with One Piece Flip Up Angle Adjustable Footplate*

Heel Loops [50ACC008]	N/C
Adjustable Heel Loops [50ACC029]	£25
Omit Heel Loops [50ACC011]	N/C

## Calf Straps

Velcro Adjustable [50ACC002]	N/C
Padded Velcro Adjustable [50ACC003]	£67
Bodypoint Padded AeroMesh [50ACC010]	£103

## Residual Limb Support

	Limb Support	Price	Omit	Cover		Pad Size					
Left	Left Limb Support	£170	Omit Left Hanger [50RLS004]	Smooth [50RLS006]	Air Mesh [50RLS007]	5 x 8 [50RLS008]	7 x 10 [50RLS009]	9 x 10 [50RLS010]	9 x 14 [50RLS011]	11 x 10 [50RLS012]	11 x 14 [50RLS013]
Right	Right Limb Support	£170	Omit Right Hanger [50RLS005]	Smooth [50RLS020]	Air Mesh [50RLS021]	5 x 8 [50RLS014]	7 x 10 [50RLS015]	9 x 10 [50RLS016]	9 x 14 [50RLS017]	11 x 10 [50RLS018]	11 x 14 [50RLS019]

## Front Seat Height and Casters

### Front Seat Height

Select Front Seat Height in row corresponding to desired caster size.

*Note: Selections in first three Front Seat Height columns use Super Low Front Frame. Super Low Front Frame not available with Titanium Frame or Fixed Front Frame.*

Caster Size	Front Seat Height													
	Uses Super Low Front Frame													
4"	13.5"	14"	14.5"	15"	15.5"	16"	16.5"	17"	17.5"	18"	18.5"	19"	19.5"	20"
5"	14.5"	15"	15.5"	16"	16.5"	17"	17.5"	18"	18.5"	19"	19.5"	20"	20.5"	
6"	15.5"	16"	16.5"	17"	17.5"	18"	18.5"	19"	19.5"	20"	20.5"	21"		
7"	17"	17.5"		18"	18.5"	19"	19.5"	20"	20.5"	21"				
8"	17.5"	18"	18.5"	19"	19.5"	20"	20.5"	21"						

### Castors

		Caster Size					
		Price	4"	5"	6"	7"	8"
Caster Type	.75" Rollerblade	N/C	[50CAS004]				
	.75" Lighted Rollerblade	£81	[50CAS005]				
	1" Poly	N/C	[50CAS006]	[50CAS010]	[50CAS014]	[50CAS019]	[50CAS020]
	1" Poly Aluminum	£96	[50CAS007]	[50CAS011]			
	1" Pneumatic	£103			[50CAS017]		[50CAS021]
	1.5" Poly	N/C	[50CAS032]	[50CAS025]	[50CAS026]		
	1.5" Soft Roll Aluminum	£96	[50CAS009]	[50CAS013]	[50CAS016]		
	2" Poly	N/C			[50CAS024]		
2" Pneumatic w/Foam Insert	£242					[50CAS023]	

### Maximum Seat Dump Reference

(Front Seat Height minus Rear Seat Height)

Seat Depth	14"	15"	16"	17"	18"	19"	20"
Maximum Dump Standard Caster Housing	1.5"	2"	2"	2"	2.5"	2.5"	2.5"
Maximum Dump Performance Caster Housing	3"	3"	3"	3.5"	3.5"	4"	4"

*Depending on configuration, chairs with Fixed Front Frame Type may achieve up to an additional .75"*

## Fork Stems

Best Fit [50CST006]	N/C
<i>Ki Mobility will pick Fork and Stem combination to optimize performance based on Caster Size and Front Seat Height</i>	
Standard Fork Stem [50CST001]	N/C
3/4" Longer Fork Stem [50CST002]	N/C
Frog Legs Suspension Fork [50CST013]	£344
<i>Available only with 4", 5" or 6" Caster Sizes. Not available with Caster Pin Locks, 350 lb Weight Capacity, 2" Poly caster or 2" Pneumatic w/ Foam Insert caster.</i>	
1 1/2" Longer Fork Stem [50CST003]	N/C
Frog Legs II Suspension System [50CST014]	£424
<i>Lightweight shock absorbing fork. Available only with 4", 5" or 6" Caster Sizes. Not available with Caster Pin Locks, 350 lb Weight Capacity, 2" Poly caster or 2" Pneumatic w/ Foam Insert caster.</i>	

## Frog Legs Options

User Weight Under 80 lbs [50CST004]	N/C
User Weight 80-120 lbs [50CST008]	N/C
User Weight 120-160 lbs [50CST009]	N/C
User Weight 160-220 lbs [50CST010]	N/C
User Weight 220-260 lbs [50CST011]	N/C
User Weight 260-300 lbs [50CST012]	N/C

## Caster Pin Locks

Caster Pin Locks [50CST005]	£194
<i>Not available with Performance Caster Housing, 1 1/2" Longer Fork Stem or Frog Legs Fork</i>	

## Caster Housing

Standard Housing [50CSH001]	N/C
Performance Housing [50CSH002]	£37
<i>Not available with 350 lb Weight Capacity</i>	

## Wheel Locks

### Wheel Lock

Push to Lock [50WLK002]	N/C
Push to Lock Flush Mount [50WLK001]	£67
Pull to Lock [50WLK003]	N/C
Low Profile Scissor Lock [50WLK012]	£67
<i>Not available with Low Poly, Full Poly or Sentinel Solid tires</i>	
Drum Brake [50WLK015]	£532
<i>Available only with 20", 22" and 24" MAXX Spoke wheels. Not available with Reclining Backrest, Straight without Push Handle or Straight with Fold Down Push Handle Back Posts, or Low Profile Push Handle. Axle is Quick Release.</i>	
Dual - Drum Brake with Push to Lock [50WLK016]	£653
<i>Available only with 20", 22" and 24" MAXX Spoke wheels. Not available with Reclining Backrest, Straight without Push Handle or Straight with Fold Down Push Handle Back Posts, or Low Profile Push Handle. Axle is Quick Release.</i>	
Dual - Drum Brake with Pull to Lock [50WLK017]	£653
<i>Available only with 20", 22" and 24" MAXX Spoke wheels. Not available with Reclining Backrest, Straight without Push Handle or Straight with Fold Down Push Handle Back Posts, or Low Profile Push Handle. Axle is Quick Release.</i>	
Hemi Wheel Lock - Left [50WLK023]	£67
<i>Extension Handle included</i>	
Hemi Wheel Lock - Right [50WLK024]	£67
<i>Extension Handle included</i>	
Omit [50WLK007]	N/C

### Wheel Lock Options

Extension Handles [50WLK004]	£49
<i>Available only with Push to Lock or Pull to Lock Wheel Locks</i>	
Grade Aids [50WLK009]	£145
<i>Available only with Push to Lock Wheel Locks and Pneumatic or Pneumatic with Airless Insert tires</i>	

## Rear Seat Height and Wheels

### Rear Seat Height

Select Rear Seat Height in row corresponding to Wheel Size and Tire Profile. 14" Seat Depth not available with 24", 25" or 26" wheel sizes.

Wheel Size	Tire Profile	Rear Seat Height									
		11.5" <sup>3</sup>	12" <sup>1</sup>	12.5" <sup>3</sup>	13" <sup>1</sup>	13.5" <sup>3</sup>	14" <sup>3</sup>	14.5"	15" <sup>3</sup>	15.5"	16" <sup>3</sup>
20"	Low	11.5" <sup>3</sup>	12" <sup>1</sup>	12.5" <sup>3</sup>	13" <sup>1</sup>	13.5" <sup>3</sup>	14"	14.5" <sup>3</sup>	15"	15.5" <sup>3</sup>	16" <sup>3</sup>
	Full	12" <sup>3</sup>	12.5" <sup>1</sup>	13" <sup>3</sup>	13.5" <sup>1</sup>	14" <sup>3</sup>	14.5"	15" <sup>3</sup>	15.5"	16" <sup>3</sup>	16.5" <sup>3</sup>
22"	Low	12.5" <sup>3</sup>	13" <sup>1</sup>	13.5" <sup>3</sup>	14" <sup>1</sup>	14.5" <sup>3</sup>	15"	16"	16.5" <sup>3</sup>	17" <sup>3</sup>	
	Full	13" <sup>3</sup>	13.5" <sup>1</sup>	14" <sup>3</sup>	14.5" <sup>1</sup>	15" <sup>3</sup>	15.5"	16" <sup>3</sup>	16.5"	17" <sup>3</sup>	17.5" <sup>3</sup>
24"	Low	13.5" <sup>3</sup>	14" <sup>1</sup>	14.5" <sup>3</sup>	15" <sup>1</sup>	15.5" <sup>3</sup>	16"	16.5" <sup>3</sup>	17"	17.5" <sup>3</sup>	18" <sup>3</sup>
	Full	14" <sup>3</sup>	14.5" <sup>1</sup>	15" <sup>3</sup>	15.5" <sup>1</sup>	16" <sup>3</sup>	16.5"	17" <sup>3</sup>	17.5"	18" <sup>3</sup>	18" <sup>3</sup>
25"	Low	14" <sup>3</sup>	14.5" <sup>1</sup>	15" <sup>3</sup>	15.5" <sup>1</sup>	16" <sup>3</sup>	16.5"	17" <sup>3</sup>	17.5"	18" <sup>3</sup>	18.5" <sup>3</sup>
26"	Low	15" <sup>1,2</sup>	15.5" <sup>3</sup>	16" <sup>1</sup>	16.5" <sup>3</sup>	17"	17.5" <sup>3</sup>	18"	18.5" <sup>3</sup>	19" <sup>3</sup>	
	Full	15" <sup>2,3</sup>	15.5" <sup>1,2</sup>	16" <sup>3</sup>	16.5" <sup>1</sup>	17" <sup>3</sup>	17.5"	18" <sup>3</sup>	18.5"	19" <sup>3</sup>	19.5" <sup>3</sup>

## Axle Plates

No selection needed with One Arm Drive

Universal [50AXP009]	N/C
Amputee [50AXP003]	£170

*Not available with Swing Away Armrests or Wheelchair Drive Ready Option.  
Required with Reclining Backrest.*

## One Arm Drive

One Arm Drive - Quick Release [50DRV007]	£846
--	------

*Available only with 20", 22" and 24" spoke wheels and Quick Release Axles.  
Spoke wheel style is exclusive to OAD. Not available with 4 degrees of camber. Standard handrim will be Aluminum Anodized.*

## Drive Side

Left Side Drive [50DRV002]	N/C
Right Side Drive [50DRV003]	N/C

## Drive Handrim

Aluminum Drive Handrim [50DRV004]	N/C
Poly Coated Drive Handrim [50DRV005]	£194

## Wheels

No selection needed with One Arm Drive

		Rear Wheel Size					
		Price	20"	22"	24"	25"	26"
Rear Wheel Type	Maxx Spoke	N/C	[50WH1120]	[50WH1122]		[50WH1125]	[50WH1126]
	Maxx Spoke	N/C			[50WH1124]		
	Mag	N/C	[50WHE120]	[50WHE122]	[50WHE124]		
	Superlight Spoke <sup>1</sup>	N/C			[50WHE224]	[50WHE225] <sup>2</sup>	
	Spinergy Spox	£768		[50WHE422]	[50WHE424]	[50WHE425]	[50WHE426]
	Spinergy LX	£816		[50WHE522]	[50WHE524]	[50WHE525]	[50WHE526]

<sup>1</sup> Not available with 350 lb Weight Capacity

<sup>2</sup> Available only with Superlight and Plastic Coated Handrims

## Spinergy Options

Select spoke color and handrim tabs for Spinergy wheels

Black Spokes [50SSC001]	N/C
Red Spokes [50SSC002]	Blue Spokes [50SSC005]
White Spokes [50SSC003]	£91

## Rear Tire Type

### Full Profile

Street Pneumatic 1 3/8" [50TIR001]	N/C
<i>Not available on 25" wheel size</i>	
Full Poly 1 3/8" [50TIR002]	N/C
<i>Not available on Superlight wheel or 25" or 26" wheel sizes</i>	
Pneumatic with Airless Insert 1 3/8" [50TIR004]	£110
<i>Not available on Superlight wheel or 25" wheel size</i>	
Schwalbe HS-127 1 3/8" [50TIR017]	£110
<i>Available only on 24" wheel size. Not available on Mag wheel.</i>	
Schwalbe Downtown 1 3/8" [50TIR020]	£110

### Low Profile

Low Poly 1" [50TIR003]	£110
<i>Not available on Superlight wheel, 25" or 26" wheel sizes or 350 lb Weight Cap.</i>	
Sentinel Solid Tire 1" [50TIR005]	£110
<i>Not available on Superlight wheel or 20" and 22" Mag wheels</i>	
High Pressure 1" [50TIR006]	N/C
<i>Not available on Mag wheel</i>	
Schwalbe One 1" [50TIR018]	£110
<i>Not available on Mag wheel, 20", 22" or 26" wheel sizes</i>	
Schwalbe Marathon Plus Evolution 1" [50TIR010]	£110
<i>Not available on Mag wheel, 20" or 22" wheel sizes</i>	

## Handrims

No selection needed with One Arm Drive

### Standard Handrims

Aluminum Anodized [50HRM001]	N/C
Superlight Aluminum Anodized [50HRM010]	N/C
<i>Available only on 24" and 25" wheel sizes. Available only on MAXX Spoke, Superlight or Spinergy wheels.</i>	
Plastic Coated [50HRM002]	£110
Projection (8 Vertical Projections) [50HRM003]	£285
<i>Available only with full profile tires in combination with 22", 24", and 26" MAXX Spoke wheels, or 22" and 24" Mag wheels</i>	
Omit Handrim [50HRM012]	N/C

### Ergonomic Handrims

*Not available with 20" or 26" wheels. Tab mount only available with 24" wheels.*

Ergo Lite Regular Grip [50HRM005]	£272
Ergo Standard Regular Grip [50HRM008]	£272
Ergo Lite Extra Grip [50HRM006]	£272
Ergo Standard Extra Grip [50HRM009]	£272
Ergo Lite Omit Grip [50HRM004]	£272
Ergo Standard Omit Grip [50HRM007]	£272

## Handrim Options

### Handrim Mount

*If no selection is made, Standard Mount will be used.*

Close Mount [50TAB001]	N/C
3M Micro Grip [50ACC006]	£91

*Not available on Plastic Coated or Projection Handrims*

## Axles

*No selection needed with One Arm Drive or Drum Brake*

Quick Release [50AXO003]	N/C
Quad Quick Release [50AXO001]	£127

*Not available with Mag or Spinergy wheels*

## Axle Position

### Rear Wheel Camber

0° [50WC00]	2° [50WC02]	4° [50WC04]	N/C
-------------	-------------	-------------	-----

### Center of Gravity

Rear [50COG003]	Center [50COG002]	Forward [50COG001]	N/C
-----------------	-------------------	--------------------	-----

*No selection needed with Amputee Axle Plate*

## Armrests

Armrest Type	Length	Price
Height Adjustable T-Arm Tall (11.5"-15")	Desk [50ARM039]	N/C
	Full [50ARM040]	
Height Adjustable T-Arm (8.5"-13.5")	Desk [50ARM003]	£248
	Full [50ARM004]	
Height Adjustable T-Arm Low (6.5"-10.5") <sup>1</sup>	Desk [50ARM022]	£248
	Full [50ARM023]	
Height Adjustable Flip Back T-Arm Tall (11.5"-15") <sup>2</sup>	Desk [50ARM041]	£230
	Full [50ARM042]	
Height Adjustable Flip Back T-Arm (9.5"-13.5") <sup>2</sup>	Desk [50ARM029]	£230
	Full [50ARM030]	

<sup>1</sup> If requested configuration is unavailable, Standard Height T-Arm will be used

<sup>2</sup> Available only with Curved Rear Frame

<sup>3</sup> Not available with Seat Depths less than 16"

<sup>4</sup> Available only with Standard Back Posts

<sup>5</sup> Available only with Depth Adjustable and Reclining Back Posts

<sup>6</sup> Not available with Arm Pads

<sup>7</sup> Not Available with 15" Rear Seat Height and 24" wheels

Armrest Type	Length	Price
Flip back (10.5") <sup>3,4,7</sup>	Desk [50ARM009]	N/C
	Full [50ARM008]	
Angle Adjustable Locking Extendable Flip Up <sup>5</sup>	Desk [50ARM010]	£272
	Full [50ARM011]	
Swing Away (10.5"-12.5") <sup>4,6</sup>	[50ARM002]	£248
Swing Away with Side Guard (10.5"-12.5") <sup>4,6</sup>	[50ARM018]	£248
Omit Armrests	[50ARM007]	N/C

## Arm Pads

Standard Arm Pad [50ARM026]	N/C
Waterfall Arm Pad [50ARM020]	£25
<i>Available only with Height Adjustable T-Arm Tall, Height Adjustable T-Arm, and Height Adjustable T-Arm Low</i>	
Foam Grip [50ARM025]	N/C
<i>Available only with Angle Adjustable Locking Extendable Flip Up Armrest</i>	

## Accessories

### Removable Side Guard

*Not available with T-Arms, Flip Back or Swing Away with Side Guard Armrests*

Adult Composite Side Guards [50SDG002]	N/C
Pediatric Composite Side Guards [50SDG001]	N/C
Aluminum Straight - Black Matte [50SDG024]	£212
Aluminum Straight - Match Frame Color [50SDG025]	£272
<i>Not available with Hydrographics or Titanium Frame</i>	
Aluminum Adjustable Fender w/Foam Pad - Black Matte [50SDG014]	£303
Aluminum Adjustable Fender w/Foam Pad - Match Frame Color [50SDG015]	£303
<i>Not available with Hydrographics or Titanium Frame</i>	
Carbon Fiber Adjustable Fender [50SDG013]	£393

### Anti-tips

Rear Anti-Tips [50ATP001]	£91
Swing-In Anti-Tip Single - Left [50ATP007]	£105
<i>Not available with 350 lb Weight Capacity, Reclining Back Post or Depth Adjustable Rear Frame. Not available with left side O2 Holder.</i>	
Swing-In Anti-Tip Single - Right [50ATP008]	£105
<i>Not available with 350 lb Weight Capacity, Reclining Back Post or Depth Adjustable Rear Frame. Not available with right side O2 Holder.</i>	
Swing-In Anti-Tip Pair [50ATP009]	£213
<i>Not available with 350 lb Weight Capacity, Reclining Back Post or Depth Adjustable Rear Frame. Not available with Cane and Crutch Holder or O2 Holder.</i>	
Omit [50ATP000]	N/C

### Pelvic Positioning Belts

1 1/2" Auto Buckle Pelvic Positioning Belt [50POB001]	£68
1 1/2" Padded Auto Buckle Pelvic Positioning Belt [50POB002]	£151
2" Airline Buckle Pelvic Positioning Belt [50POB003]	£110
Small Bodypoint Evoflex (16-24") [50POB007]	£188
Medium Bodypoint Evoflex (19-29") [50POB008]	£188
Large Bodypoint Evoflex (25-37") [50POB009]	£188
Belt Mounting Band Clamp [50ACC054]	£103

### O2 Holder

*Not Available with Depth Adjustable Rear Frame*

O2 Holder - Left Side [50ACC032]	£303
O2 Holder - Right Side [50ACC033]	£303

### Other Accessories

Tool Kit [50ACC005]	£20
Removable Underseat Pouch [50ACC015]	£30
Pediatric Ki Mobility Backpack [50ACC026]	£100
Ki Mobility Backpack [50ACC009]	£100
Luggage Carrier [50ACC025]	£73
Black Neoprene Impact Guards [50IMG001]	£86
Clear Spoke Guards [50SPG001]	£151
Black Spoke Guards [50SPG002]	£151
Cane and Crutch Holder [50ACC007]	£91
IV Pole [50ACC022]	£260

## Color

### Gloss Finish

Black [50FCS001]	N/C
Blue Skies [50FCS018]	N/C
Bubblegum Pink [50FCS011]	N/C
Candy Blue [50FCS012]	N/C
Candy Purple [50FCS010]	N/C
Candy Red [50FCS003]	N/C
Copper [50FCS009]	N/C
Granny Smith [50FCS008]	N/C
Metallic Orange [50FCS032]	N/C
Yellow [50FCS029]	N/C

### Matte Finish

Military Green [50FCS025]	N/C
Charcoal Gray [50FCS031]	N/C
Emerald Green [50FCS019]	N/C
Iridescent Burgundy [50FCS021]	N/C
Iridescent Navy [50FCS016]	N/C
Retro Red [50FCS022]	N/C
Royal Blue [50FCS023]	N/C
Shadow Black [50FCS015]	N/C
Smooth White [50FCS024]	N/C

### Hydrographics

AOR1 (Desert Digital Camo) [50GIT002]	£725
NEXT Camo [50GIT003]	£725
Pink Camo [50GIT005]	£725
Carbon Fiber [50GIT008]	£725
<b>Select Hydrographics Base Color</b>	
Muddy Girl [50GIT010]	£725
Zebra [50GIT007]	£725
<i>Select Hydrographic Base Color. Zebra not available in Black.</i>	

### Hydrographics Base Color

*Required with Hydrographics*

Black [50GBC001]	Blue [50GBC002]	
Pink [50GBC005]	Granny Smith [50GBC004]	
Purple [50GBC006]	Red [50GBC007]	N/C

### Colored Anodizing Package

*Select desired color. If no selection is made, black will be used.*

	Option	Standard Caster Housing <sup>2</sup>	Caster Fork <sup>3</sup>	Axle Plate <sup>4</sup>
	Retail	£61	£61	£61
Color	Black <sup>1</sup>			
	Blue			
	Bright Green			
	Orange			
	Pewter			
	Pink			
	Purple			
	Red			

<sup>1</sup> Black is no charge

<sup>2</sup> Available only if corresponding option selection was made

<sup>3</sup> Not available with 2" wide casters, Frog Legs Suspension Fork or Frog Legs II Suspension System

<sup>4</sup> Color Anodized not available with amputee axle plate



Catalyst Front Seat Height Matrix																														
4" Caster				5" Caster				6" Caster				7" Caster				8" Caster														
FSH	Super Low	Stem	Fork	FSH	Super Low	Stem	Fork	FSH	Super Low	Stem	Fork	FSH	Super Low	Stem	Fork	FSH	Super Low	Stem	Fork											
15	13.5	Std	4"	16	14.5	Std	4"	17	15.5	Std	5"	18	16.5	Std	6"	19	17.5	Std	6"	6"										
15.5	14			16.5	15.5	3/4		17.5	16.5	3/4		18.5	17			19.5	18.5	3/4												
15.5 <sup>1</sup>	14 <sup>1</sup>	17.5		16	1 1/2	18.5		17	1 1/2	18.5		17.5	20.5	n/a		1 1/2														
16	14.5	3/4		5"	16	14.5	Std	5"	17	15.5	Std	6"	19	18		1 1/2	7"	19	17.5		Std	7"	7"							
16.5 <sup>1</sup>	15 <sup>1</sup>				16.5	15			17.5	16			18.5	17				19.5	18					20	18.5					
17	15.5	16.5			15.5	18	16.5		18.5	17.5			19.5	18		20.5		19												
15.5	14	Std			5"	17	16	3/4	5"	17.5	16.5	3/4	6"	18.5		17		Std	7"		20.5	19		Std	7"	7"				
16	14.5					17.5	16			18	17			19.5		17.5					19.5	18					20	19		
16	15	18				16.5	18.5	17.5		19.5	18			20		18.5		20.5			n/a									
16.5	15.5	3/4				5"	16.5	15	Std	5"	18.5	17	1 1/2	6"		20		18.5			3/4	7"		21	n/a		3/4	7"	7"	
17	15.5		17				15.5	19			17.5	19.5			18	20.5		19		21				n/a						
17.5	16	17.5	16				18.5	17	19.5		18	20			19	21		n/a												
16.5	15	Std	6"				17.5	16	3/4	6"	18	16.5	Std	7"	20	19		1 1/2		7"					7"					
17	15.5			18			16.5	18.5			17	19.5			18	20.5	n/a													
17	16	18.5		17			19	17.5	20		18.5	20			18.5															
17.5	16.5	3/4		6"			18	16.5	1 1/2	6"	18.5	17	3/4	7"	20.5	19	Std	7"												
18	16.5				18.5		17	19.5			18	20.5			19	21			n/a											
18.5	17	19			17.5		19.5	18	20		18.5	21			n/a															
18	16.5	Std			7"		18	16.5	Std	7"	19.5	18.5	1 1/2	7"					7"											
18.5	17					18.5	17	20			19	20.5														19				
18.5	17.5	19				17.5	19.5	18	20		18.5	21														n/a				
19	18	3/4				7"	18.5	17.5	3/4	7"	20	18.5	1 1/2	7"																
19.5	18.5		19				18	20.5			19	21								n/a										
20	19	19.5	18.5				20.5	19	21		n/a																			
							19.5	18	1 1/2	7"	6" x 2" caster not available with 5" fork									7"										
				20			18.5																							
				20.5			19																							

<sup>1</sup> Not available with 4 x 1.5" Soft Roll Aluminum

Seat heights with "FSH" heading available with Standard Front Frame. Seat heights with "Super Low" heading available with Super Low Front Frame. Actual seat height may vary depending upon configuration.

## Axiom G

		Width															
		14"		15"		16"		17"		18"		19"		20"		22"	
		Cushion	Outer Cover	Cushion	Outer Cover	Cushion	Outer Cover	Cushion	Outer Cover	Cushion	Outer Cover	Cushion	Outer Cover	Cushion	Outer Cover	Cushion	Outer Cover
Retail		£85	£70	£85	£70	£85	£70	£85	£70	£85	£70	£85	£70	£85	£70	£85	£70
Height	14"																
	15"																
	16"																
	17"																
	18"																
	19"																
	20"																
	21"																
	22"																

## Axiom S

		Width															
		14"		15"		16"		17"		18"		19"		20"		22"	
		Cushion	Outer Cover	Cushion	Outer Cover	Cushion	Outer Cover	Cushion	Outer Cover	Cushion	Outer Cover	Cushion	Outer Cover	Cushion	Outer Cover	Cushion	Outer Cover
Retail		£120	£70	£120	£70	£120	£70	£120	£70	£120	£70	£120	£70	£120	£70	£120	£70
Height	14"																
	15"																
	16"																
	17"																
	18"																
	19"																
	20"																
	21"																
	22"																

## Axiom P

		Width															
		14"		15"		16"		17"		18"		19"		20"		22"	
		Cushion	Outer Cover	Cushion	Outer Cover	Cushion	Outer Cover	Cushion	Outer Cover	Cushion	Outer Cover	Cushion	Outer Cover	Cushion	Outer Cover	Cushion	Outer Cover
Retail		£170	£70	£170	£70	£170	£70	£170	£70	£170	£70	£170	£70	£170	£70	£170	£70
Height	14"																
	15"																
	16"																
	17"																
	18"																
	19"																
	20"																
	21"																
	22"																

## Axiom SP - Visco

		Width																										
		14"			15"			16"			17"			18"			19"			20"			21"			22"		
		Cushion	Cover		Cushion	Cover		Cushion	Cover		Cushion	Cover		Cushion	Cover		Cushion	Cover		Cushion	Cover		Cushion	Cover				
Retail		£350	£70	£350	£70	£350	£70	£350	£70	£350	£70	£350	£70	£350	£70	£350	£70	£350	£70	£350	£70	£350	£70					
Height	14"																											
	15"																											
	16"																											
	17"																											
	18"																											
	19"																											
	20"																											
	21"																											
	22"																											

## Axiom SP - Fluid

		Width																										
		14"			15"			16"			17"			18"			19"			20"			21"			22"		
		Cushion	Cover		Cushion	Cover		Cushion	Cover		Cushion	Cover		Cushion	Cover		Cushion	Cover		Cushion	Cover		Cushion	Cover				
Retail		£320	£70	£320	£70	£320	£70	£320	£70	£320	£70	£320	£70	£320	£70	£320	£70	£320	£70	£320	£70	£320	£70					
Height	14"																											
	15"																											
	16"																											
	17"																											
	18"																											
	19"																											
	20"																											
	21"																											
	22"																											


**Solid Seat Insert**

*Weight capacity: 350 lbs*

		Width								
		14"	15"	16"	17"	18"	19"	20"	21"	22"
		Cushion	Cushion	Cushion	Cushion	Cushion	Cushion	Cushion	Cushion	Cushion
Retail		£55	£55	£55	£55	£55	£55	£55	£55	£55
Height	14"									
	15"									
	16"									
	17"									
	18"									
	19"									
	20"									
	21"									
	22"									

Axiom Back Width	Fits Back Frame Widths	
	Inner Slots	Outer Slots
12" (30.5cm)	11" (28.0cm)	12" -14" (30.5cm -35.5cm)
14" (35.5cm)	12" -13" (30.5cm -33.0cm)	14" -16" (35.5cm -40.5cm)
16" (40.5cm)	14" -15" (35.5cm -38.0cm)	16" -18" (40.5cm -45.5cm)
18" (45.5cm)	16" - 17" (40.5cm - 43.0cm)	18" - 20" (45.5cm - 50.5cm)
20" (51.0cm)	18" -19" (46.0cm -48.5cm)	20" -22" (50.5cm -55.5cm)
22" (56.0cm)	20" -21" (51.0cm -53.5cm)	22" -24" (56.0cm -61.0cm)

Hardware Description	Code
Improved DEX Hardware, removable adjustable hardware (300 lb max user weight)	<b>DEX</b>
New BTX Hardware, removable adjustable hardware not recommended for, or installed on Ultralight Rigid Frame Wheelchairs (400 lb max user weight)	<b>BTX</b>
New Fixed Hardware, non-removeable hardware up to (300 lb max user weight)	<b>FXT</b>



## Axiom Back

Hardware is included with back. Low Back selections are not available with headrests, modifications, or posture supports. Weight capacity: 300 lbs.

### Axiom Posterior Back (PB)

		Width						
		12"	14"	16"	18"	20"	22"	
Retail		£550	£550	£550	£550	£550	£550	
Height	Low	8"	[AXPB1208]	[AXPB1408]	[AXPB1608]	[AXPB1808]	[AXPB2008]	
		10"		[AXPB1410]	[AXPB1610]	[AXPB1810]	[AXPB2010]	[AXPB2210]
		12"	[AXPB1212]	[AXPB1412]	[AXPB1612]	[AXPB1812]	[AXPB2012]	
	Mid	14"		[AXPB1414]	[AXPB1614]	[AXPB1814]	[AXPB2014]	[AXPB2214]
		16"	[AXPB1216]	[AXPB1416]	[AXPB1616]	[AXPB1816]	[AXPB2016]	
		18"		[AXPB1418]	[AXPB1618]	[AXPB1818]	[AXPB2018]	[AXPB2218]

### Back Hardware

DEX Hardware [DEX]    BTX Hardware [BTX]    FXT Hardware [FXT] \_\_\_\_\_ N/C

### Lateral Support for AXPB Backs

Lateral Trunk Supports for AXPB Backs							
Pair	Left	Right	Description	Item Code	Pair	Each	
			3 x 4 Lateral Pads with Fixed Mounting Hardware for Axiom Posterior Back	PBLF3x4P	£295	£148	
			4 x 5 Lateral Pads with Fixed Mounting Hardware for Axiom Posterior Back	PBLF4x5P	£295	£148	
			3 x 4 Lateral Pads with Swing Away Mounting Hardware for Axiom Posterior Back	PBLS3x4P	£628	£314	
			4 x 5 Lateral Pads with Swing Away Mounting Hardware for Axiom Posterior Back	PBLS4x5P	£628	£314	
			3 x 4 Lateral Pads with Adjustable Swing Away Mounting Hardware	PBLAS3x4P	£628	£314	
			4 x 5 Lateral Pads with Adjustable Swing Away Mounting Hardware	PBLAS4x5P	£628	£314	

### Axiom Posterior Lateral Back (PL)

		Width						
		12"	14"	16"	18"	20"	22"	
Retail		£550	£550	£550	£550	£550	£550	
Height	Low	8"	[AXPL1208]	[AXPL1408]	[AXPL1608]	[AXPL1808]	[AXPL2008]	
		10"		[AXPL1410]	[AXPL1610]	[AXPL1810]	[AXPL2010]	[AXPL2210]
		12"	[AXPL1212]	[AXPL1412]	[AXPL1612]	[AXPL1812]	[AXPL2012]	
	Mid	14"		[AXPL1414]	[AXPL1614]	[AXPL1814]	[AXPL2014]	[AXPL2214]
		16"	[AXPL1216]	[AXPL1416]	[AXPL1616]	[AXPL1816]	[AXPL2016]	
		18"		[AXPL1418]	[AXPL1618]	[AXPL1818]	[AXPL2018]	[AXPL2218]
Tall	20"	[AXPL1220]	[AXPL1420]	[AXPL1620]	[AXPL1820]	[AXPL2020]		
	22"		[AXPL1422]	[AXPL1622]	[AXPL1822]	[AXPL2022]	[AXPL2222]	

### Back Hardware

DEX Hardware [DEX]    BTX Hardware [BTX]    FXT Hardware [FXT] \_\_\_\_\_ N/C

### Lateral Support for AXPL

Lateral Trunk Supports for AXPL & AXLB Backs						
Pair	Left	Right	Description	Item Code	Pair	Each
			3 x 4 Lateral Pads with Fixed Mounting Hardware for Axiom Posterior Lateral Back	PLLF3x4P	£295	£295
			4 x 5 Lateral Pads with Fixed Mounting Hardware for Axiom Posterior Lateral Back	PLLF4x5P	£295	£295
			3 x 4 Lateral Pads with Swing Away Mounting Hardware for Axiom Posterior Lateral Back	PLLS3x4P	£628	£314
			4 x 5 Lateral Pads with Swing Away Mounting Hardware for Axiom Posterior Lateral Back	PLLS4x5P	£628	£314
			3 x 4 Lateral Pads with Adjustable Swing Away Mounting Hardware	PLLAS3x4P	£628	£314
			4 x 5 Lateral Pads with Adjustable Swing Away Mounting Hardware	PLLAS4x5P	£628	£314

### Axiom Deep Lateral Back (DL)

		Width						
		12"	14"	16"	18"	20"	22"	
Retail		£550	£550	£550	£550	£550	£550	
HCPC		E2620	E2620	E2620	E2620	E2620	E2621	
Height	Mid	14"	[AXDL1214]	[AXDL1414]	[AXDL1614]	[AXDL1814]	[AXDL2014]	[AXDL2214]
		16"		[AXDL1416]	[AXDL1616]	[AXDL1816]	[AXDL2016]	
		18"	[AXDL1218]	[AXDL1418]	[AXDL1618]	[AXDL1818]	[AXDL2018]	[AXDL2218]
	Tall	20"	[AXDL1220]	[AXDL1420]	[AXDL1620]	[AXDL1820]	[AXDL2020]	
		22"		[AXDL1422]	[AXDL1622]	[AXDL1822]	[AXDL2022]	[AXDL2222]

### Back Hardware

DEX Hardware [DEX]    BTX Hardware [BTX]    FXT Hardware [FXT] \_\_\_\_\_ N/C

## Headrest

Description	Std. Straight Removable Adjustable Hardware		Multi-Axis Removable Adjustable Hardware		Posalinc Removable Adjustable Hardware		HCPCS
	Item Code	Retail	Item Code	Retail	Item Code	Retail	
Comfort Soft Headrest 8"	[CS08SS]	£268	[CS08MA]	£268	[CS08PL]	£268	E0955 / E1028
Comfort Soft Headrest 10"	[CS10SS]	£268	[CS10MA]	£268	[CS10PL]	£268	E0955 / E1028
Comfort Soft Headrest 14"	[CS14SS]	£268	[CS14MA]	£268	[CS14PL]	£268	E0955 / E1028
Contour Positioning Pad, 8" X 6.25"	[CP08SS]	£303	[CP08MA]	£303	[CP08PL]	£303	E0955 / E1028
Contour Positioning Pad, 9" X 6.5"	[CP09SS]	£303	[CP09MA]	£303	[CP09PL]	£303	E0955 / E1028
Contour Positioning Pad, 10" X 7.5"	[CP10SS]	£303	[CP10MA]	£303	[CP10PL]	£303	E0955 / E1028
Tri-soft Headrest 10"	[TS10SS]	£386	[TS10MA]	£386	[TS10PL]	£386	E0955 / E1028
Tri-soft Headrest 14"	[TS14SS]	£386	[TS14MA]	£386	[TS14PL]	£386	E0955 / E1028

Description	Item Code	Retail	HCPCS
Universal Headrest Mounting Plate	[UMP]	£43	K0108

*Not required for Axiom Headrests*

## Anterior Thoracic Supports

*Selection of a modification from the Axiom Back Modifications section required in order to mount posture supports. Not available with Low Backs.*

### Anterior Thoracic Supports

Description	Item Code	HCPCS	Retail
Dynaform Chest Support: 2-Piece with Hook&Loop, 3" High Pad, Small	[14941]	E0960	£100
Dynaform Chest Support: 2-Piece with Hook&Loop, 4" High Pad, Medium	[14942]	E0960	£100
Dynaform Chest Support: 2-Piece with Hook&Loop, 5" High Pad, Large	[14943]	E0960	£100
Dynaform Chest Support: 1-Piece with Buckle, 3" Height x 20" Width, Small	[14924]	E0960	£100
Dynaform Chest Support: 1-Piece with Buckle, 4" Height x 20" Width, Medium	[14927]	E0960	£100
Dynaform Chest Support: 1-Piece with Buckle, 4" Height x 22" Width, Large	[14928]	E0960	£100
Dynaform Anterior Trunk Support: Standard, Dynamic, X-Small	[14851]	E0960	£100
Dynaform Anterior Trunk Support: Standard, Dynamic, Small	[14852]	E0960	£100
Dynaform Anterior Trunk Support: Standard, Dynamic, Medium	[14853]	E0960	£100
Dynaform Anterior Trunk Support: Standard, Dynamic, Large	[14854]	E0960	£100
Dynaform Anterior Trunk Support: Slim Cut, Dynamic, X-Small	[14856]	E0960	£100
Dynaform Anterior Trunk Support: Slim Cut, Dynamic, Small	[14857]	E0960	£100
Dynaform Anterior Trunk Support: Slim Cut, Dynamic, Medium	[14858]	E0960	£100
Dynaform Anterior Trunk Support: Slim Cut, Dynamic, Large	[14859]	E0960	£100
AirLogic Anterior Trunk Support: Standard, Structure, Small	[14282]	E0960	£300
AirLogic Anterior Trunk Support: Standard, Structure, Medium	[14283]	E0960	£300
AirLogic Anterior Trunk Support: Standard, Structure, Large	[14284]	E0960	£300
AirLogic Anterior Trunk Support: Standard, Structure, X-Large	[14285]	E0960	£300
AirLogic Anterior Trunk Support: Slim Cut, Dynamic, Small	[14257]	E0960	£300
AirLogic Anterior Trunk Support: Slim Cut, Dynamic, Medium	[14258]	E0960	£300
AirLogic Anterior Trunk Support: Slim Cut, Dynamic, Large	[14259]	E0960	£300
AirLogic Anterior Trunk Support: Slim Cut, Dynamic, X-Large	[14260]	E0960	£300

## Accessories

Adjustable Shoulder Guides (K0108) [AXBSG] <i>Includes Anterior Thoracic Mounting Modification. Not available with Low Backs.</i>	£43
Extra Lumbar Pad (K0108) [LUMPAD] <i>Available only with Axiom Lateral Back and Axiom Deep Lateral back</i>	N/C
Privacy Flap (K0108) [PRFLAP]	£44

## Axiom Back Modifications

### Axiom Back Modifications

*Extended lead times will apply for all modifications. Makes back non-refundable.*

Anterior Thoracic Mounting Modification [AXSMM] <i>Provides holes for Anterior Thoracic Supports and Shoulder Guides</i>	N/C
---	-----

### Gloss Finish

*If no selection is made, standard Shadow Black will be used with no additional lead time*

Black Back Shell [AXMBLK]	N/C
Blue Skies Back Shell [AXMSKY]	N/C
Bubblegum Pink Back Shell [AXMPNK]	N/C
Candy Blue Back Shell [AXMCBL]	N/C
Candy Purple Back Shell [AXMPUR]	N/C
Candy Red Back Shell [AXMRED]	N/C
Copper Back Shell [AXMCOP]	N/C
Granny Smith Back Shell [AXMGSM]	N/C
Metallic Orange Back Shell [AXMMOR]	N/C
Yellow Back Shell [AXMYLW]	N/C
Silver Back Shell [AXMSIL]	N/C
Teal Back Shell [AXMTEA]	N/C

### Matte Finish

*If no selection is made, standard Shadow Black will be used with no additional lead time*

Military Green Back Shell [AXMAGM]	N/C
Charcoal Gray Back Shell [AXMCGM]	N/C
Emerald Green Back Shell [AXMEGM]	N/C
Iridescent Burgundy Back Shell [AXMIBM]	N/C
Iridescent Navy Back Shell [AXMINM]	N/C
Retro Red Back Shell [AXMRRM]	N/C
Royal Blue Back Shell [AXMRBM]	N/C
Smooth White Back Shell [AXMSWM]	N/C