



Phone (800) 981-1540 | Fax (715) 254-0996  
 sales@kimobility.com | www.kimobility.com



Company:

Ship To:

Account Number:

Date:

Address:

Quote: PO Number:

City:

State:

ATP:

Zip:

Phone:

Marked For:

Email/Fax:

Associated Quote/Order#:

Effective: 1/5/2026

Liberty Folding - Transit Not Included [LFOPTTRN] \_\_\_\_\_ £3300  
 250 lb weight capacity Tilt range: 0° - 20°

Liberty Rigid, Transit Not Included [LFRGDMT] \_\_\_\_\_ £3300  
 250 lb weight capacity Tilt range: 0° - 20°

## Frame

### Seat Width

Select Seat Width within desired Seat Width Range

Seat Width Range	Seat Width				
Small 14" - 18" [LFFRW001]	14 [LFSWT014] N/C	15 [LFSWT015] N/C	16 [LFSWT016] N/C	17 [LFSWT017] N/C	18 [LFSWT018] N/C
Medium 16" - 20" [LFFRW002]	16 [LFSWT016] N/C	17 [LFSWT017] N/C	18 [LFSWT018] N/C	19 [LFSWT019] N/C	20 [LFSWT020] N/C
Large 18" - 22" [LFFRW003]	18 [LFSWT018] N/C	19 [LFSWT019] N/C	20 [LFSWT020] N/C	21 [LFSWT021] £523	22 [LFSWT022] £523

### Seat Depth

Select Seat Depth within desired Seat Depth Range

Seat Depth Range	Seat Depth		
Short 14" - 16" [LFFDR001]	14 [LFSDE014] N/C	15 [LFSDE015] N/C	16 [LFSDE016] N/C
Medium 16" - 18" [LFFDR002]	16 [LFSDE016] N/C	17 [LFSDE017] N/C	18 [LFSDE018] N/C
Long 18" - 20" [LFFDR003]	18 [LFSDE018] N/C	19 [LFSDE019] N/C	20 [LFSDE020] N/C

### Assistive Gas Strut Set-Up

Client Weight > 200 lbs (Setting 1, 450N Spring) [LFGAS001] \_\_\_\_\_ £139  
 Client Weight 176-200 lbs (Setting 2, 450N Spring) [LFGAS002] \_\_\_\_\_ £139  
 Client Weight 151-175 lbs (Setting 3, 450N Spring) [LFGAS003] \_\_\_\_\_ £139  
*Standard on chairs less than 16" deep*  
 Client Weight 151-175 lbs (Setting 1, 300N Spring) [LFGAS004] \_\_\_\_\_ N/C  
*Standard on chairs 16" deep or larger*  
 Client Weight 125-150 lbs (Setting 2, 300N Spring) [LFGAS005] \_\_\_\_\_ N/C  
 Client Weight < 125 lbs (Setting 3, 300N Spring) [LFGAS006] \_\_\_\_\_ N/C

### Seat Option

Standard Seat Sling [LFGSS001] \_\_\_\_\_ N/C  
 Growing Adjustable Seat Pan [LFACT004] \_\_\_\_\_ £199  
 Omit Seat [LFCUS003] \_\_\_\_\_ N/C

### Transit Option

Transit System [LFRNPT] \_\_\_\_\_ N/C  
*No charge when selected with Liberty - Transit Included*  
 Omit [LFRNOP2] \_\_\_\_\_ N/C

## Backrest

### Back Post / Height Range / Back Height

Back Post	Price	Height Range	Back Height				
Stroller Handle Back Post <sup>1</sup> [LFBOP016]	N/C	Short [LFBHE002]	20 [LFBRH200]				
		Tall [LFBHE004]	24 [LFBRH240]				
Height Adjustable Straight with Push Handle [LFBOP001]	£29	Short [LFBHE002]	14 [LFBRH140]	15 [LFBRH150]	16 [LFBRH160]	17 [LFBRH170]	18 [LFBRH180]
		Medium [LFBHE003]	16 [LFBRH160]	17 [LFBRH170]	18 [LFBRH180]	19 [LFBRH190]	20 [LFBRH200]
		Tall [LFBHE004]	18 [LFBRH180]	19 [LFBRH190]	20 [LFBRH200]	21 [LFBRH210]	22 [LFBRH220]
Height Adjustable 8° Bend <sup>2</sup> with Push Handle [LFBOP003]	£29	Short [LFBHE002]	14 [LFBRH140]	15 [LFBRH150]	16 [LFBRH160]	17 [LFBRH170]	18 [LFBRH180]
		Medium [LFBHE003]	16 [LFBRH160]	17 [LFBRH170]	18 [LFBRH180]	19 [LFBRH190]	20 [LFBRH200]
		Tall [LFBHE004]	18 [LFBRH180]	19 [LFBRH190]	20 [LFBRH200]	21 [LFBRH210]	22 [LFBRH220]

<sup>1</sup> Not available with Back Upholstery

<sup>2</sup> Not available with IV Pole

### Stroller Handle Back Post Option

Available only with Stroller Handle Back Post

Removable Stroller Handle Extension [LFBOP012] £84

### Backrest Options

Fixed (Non-Folding) [LFFDB001] N/C

Angle Adjustable (Folding) [LFFDB002] £106

### Back Upholstery

If selecting both upholstery and Axiom Back, the upholstery will be installed

Standard [LFBUP001] N/C

Tension Adjustable [LFBUP002] £209

Not available on Seat Widths less than 16" combined with 21" and 22"

Backrest Heights

Omit [LFBUP003] N/C

## Footrest

### Hanger

Hanger	60°	70°	80°	90°		Retail
Extension Tube Mount	[LFHAN001]	[LFHAN002]	[LFHAN003] <sup>1</sup>			N/C
Front Mount		[LFHAN016]	[LFHAN017] <sup>1</sup>	[LFHAN004] <sup>2</sup>		N/C
Elevating Leg Rests					[LFHAN048]	N/C
Pro Elevating Leg Rests					[LFHAN044]	£390
Omit Hangers					[LFHAN046]	N/C

<sup>1</sup> Not available with 8" casters

<sup>2</sup> Not available with 8" casters in combination with 3" extension

### Footplates

Composite [LFFPLO01] N/C

Not available with Front Mount Hanger

Composite Angle Adjustable [LFFPLO14] £196

Not available with Front Mount Hanger or 14" Seat Width

Aluminum Angle Adjustable [LFFPLO02] £212

Aluminum Locking Multi-Angle Adjustable [LFFPLO20] £348

Available only with Front Mount Hanger

One Piece Flip Up Angle Adjustable [LFFPLO09] £216

Not available with Front Mount Hanger, ELR, or Pro ELR

### ELR Pad and Pro ELR Cover

Pro ELR Gel Pad Covers (Pair) [LFELRPAD] £300

Available only on Pro ELR

## Footrest / Extension / Seat to Footrest Length

Select Extension and Seat to Footrest Length in row corresponding to desired Footrest. Footrest Length of at least 2.5" shorter than Seat Height recommended. Extension Tube Hanger measurements based on 70° Hanger Bend.

Footrest	Extension	Seat to Footrest Length										
		10"	10.5"	11"	11.5"	12"	12.5"	13"	13.5"	15"	18"	
Extension Tube Hangers (60°, 70°, 80°)	Short [LFEXT002]	10"	10.5"	11"	11.5"	12"	12.5"	13"	13.5"	15"	18"	
	Medium [LFEXT003]	11.5"	12"	12.5"	13"	13.5"	14"	14.5"	15"	15.5"		
	Long [LFEXT004]	14"	14.5"	15"	15.5"	16"	16.5"	17"	17.5"	18"		
Front Mount Hangers <sup>1</sup>	Short Front Mount Hanger (70°, 80°, 90°)	Standard Selections								Requires 3" Extension		
		4" <sup>2</sup>	4.5" <sup>2</sup>	5"	5.5"	6"				8.5"	9"	9.5"
	Medium Front Mount Hanger (70°, 80°, 90°)	6.5"	7"	7.5"	8"					10"	10.5"	11"
		4" <sup>2</sup>	4.5" <sup>2</sup>	5"	5.5"	6"	6.5"	7"		10.5"	11"	11.5"
	Long Front Mount Hanger (70°, 80°)	7.5"	8"	8.5"	9"	9.5"	10"			12"	12.5"	13"
		5"	5.5"	6"	6.5"	7"	7.5"	8"	8.5"	12.5"	13"	13.5"
		9"	9.5"	10"	10.5"	11"	11.5"	12"	14"	14.5"	15"	
ELR	Standard ELR Extension Tube Short (10.25" - 15.25") [LFEXT025]											
	Standard ELR Extension Tube Medium (12.25" - 17.25") [LFEXT030]											
	Standard ELR Extension Tube Long (14.75" - 19.75") [LFEXT031]											
PRO ELR Low Pivot Setting	Pro ELR Low Pivot Setting Short (10.5" - 16") [LFEXT021]											
	Pro ELR Low Pivot Setting Long (13" - 19.5") [LFEXT022]											
PRO ELR High Pivot Setting	Pro ELR High Pivot Setting Short (8.5" - 14.5") [LFEXT023]											
	Pro ELR High Pivot Setting Long (11.5" - 18") [LFEXT024]											

<sup>1</sup> See Minimum Seat Height chart below

<sup>2</sup> Available with 90° front mount hangers

Minimum Seat Height with Front Mount Hangers						
Hanger Length	Short		Medium		Long	
	70°	80/90°	70°	80/90°	70°	80°
Hanger Bend	70°	80/90°	70°	80/90°	70°	80°
Minimum Seat Height with Standard Selection	13"	13"	13"	13.5"	14.5"	15.5"
Minimum Seat Height with 3" Extension	13.5"	14"	15.5"	16"	17"	18"

## Heel Loops

Not available with One Piece Flip Up Angle Adjustable Footplate

Heel Loops [LFACC008]	£33
Adjustable Heel Loops [LFACC029]	£56
Omit Heel Loops [LFACC011]	N/C

## Calf Straps

Velcro Adjustable [LFACC002]	£39
Padded Velcro Adjustable [LFACC003]	£76
Padded Airmesh Calf Strap [LFACC010]	£101

## Hanger Extension Tubes

Available only with Front Mount Hangers

3" Extension [LFEXT008]	£50
3" Extension (Provided as Parts) [LFEXT018]	£50

## Residual Limb Support

	Limb Support	Price	Omit	Cover		Pad Size				
Left	Left Limb Support	£209	Omit Left Hanger [LFRLS004]	Smooth [LFRLS006]	Air Mesh [LFRLS007]	5 x 8 [LFRLS008]	7 x 10 [LFRLS009]	9 x 10 [LFRLS010]	9 x 14 [LFRLS011]	11 x 10 [LFRLS012]
Right	Right Limb Support	£209	Omit Right Hanger [LFRLS005]	Smooth [LFRLS020]	Air Mesh [LFRLS021]	5 x 8 [LFRLS014]	7 x 10 [LFRLS015]	9 x 10 [LFRLS016]	9 x 14 [LFRLS017]	11 x 10 [LFRLS018]

## Hanger Release

4-Way - In-Line [LFREL002]	N/C
<i>Not available with ELR or Pro ELR</i>	
4-Way - Out [LFREL003]	N/C

## Seat Height, Casters and Rear Wheels

Select Seat Height corresponding to desired Rear Wheel and Caster Size

Rear Wheel Size	Caster Size			
	4"	5"	6"	8"
	Seat Height			
12"	13"	14"	15"	
	13.5"	14.5"	15.5"	17"
	14"	15"	16"	17.5"
	14.5"	15.5"	16.5"	18"
	15"	16"	17"	18.5"
	15.5"	16.5"	17.5"	19" <sup>2</sup>
	16"	17"	18"	
		17.5"	18.5"	
		18"	19" <sup>2</sup>	
		18.5"		
	19" <sup>2</sup>			
16"	13.5" <sup>1</sup>	14" <sup>1</sup>	15"	
	14" <sup>1</sup>	14.5"	15.5"	17"
	14.5"	15"	16"	17.5"
	15"	15.5"	16.5"	18"
	15.5"	16"	17"	18.5"
	16"	16.5"	17.5"	19"
		17"	18"	19.5"
		17.5"	18.5"	20"
		18"	19"	
		18.5"		
	19"			
20"	13"	14"	15"	
	13.5"	14.5"	15.5"	17"
	14"	15"	16"	17.5"
	14.5"	15.5"	16.5"	18"
	15"	16"	17"	18.5"
	15.5"	16.5"	17.5"	19"
	16"	17"	18"	19.5"
		17.5"	18.5"	20"
		18"	19"	
		18.5"		
	19"			

Rear Wheel Size	Caster Size			
	4"	5"	6"	8"
	Seat Height			
22"	13.5" <sup>1</sup>	14"	15"	
	14"	14.5"	15.5"	17"
	14.5"	15"	16"	17.5"
	15"	15.5"	16.5"	18"
	15.5"	16"	17"	18.5"
	16"	16.5"	17.5"	19"
		17"	18"	19.5"
		17.5"	18.5"	20"
		18"	19" <sup>2</sup>	
		18.5"		
	19" <sup>2</sup>			
24"	14.5" <sup>1</sup>	14.5" <sup>1</sup>	15"	
	15"	15"	15.5"	17"
	15.5"	15.5"	16"	17.5"
	16"	16"	16.5"	18"
		16.5"	17"	18.5"
		17"	17.5"	19"
		17.5"	18"	19.5"
		18"	18.5"	20"
		18.5"	19"	
		19"		

<sup>1</sup> Available only with Low Profile tires  
<sup>2</sup> Not available with Full Poly tire

### Casters

Caster Type		Caster Size				
		Price	4"	5"	6"	8"
Caster Type	1.25" Soft Roll	N/C	[LFCAS006]	[LFCAS010]	[LFCAS014]	[LFCAS020]
	1" Poly Aluminum	£122	[LFCAS007]	[LFCAS011]		
	1" Pneumatic	£105			[LFCAS017]	[LFCAS021]
	1.5" Soft Roll	N/C	[LFCAS032]	[LFCAS025]	[LFCAS026]	
	1.5" Soft Roll Aluminum	£209	[LFCAS009]	[LFCAS013]	[LFCAS016]	
	2" Poly	£105			[LFCAS024]	
	2" Pneumatic w/Foam Insert	£251				[LFCAS023]

## Wheels

Rear Wheel Type	Mag	Rear Wheel Size					
		Price	12"	16"	20"	22"	24"
		N/C	[LFWHE112]	[LFWHE116]	[LFWHE120]	[LFWHE122]	[LFWHE124]
	Spoke	N/C				[LFWH1122]	[LFWH1124]

## Rear Tire Type

### Full Profile

Street Pneumatic 1 3/8" [LFTIR001]	£127
Full Poly 1 3/8" [LFTIR002]	N/C
Pneumatic with Airless Insert 1 3/8" [LFTIR004]	£127

### Low Profile

Low Poly 1" [LFTIR003]	N/C
<i>Available only with 16" wheel size</i>	
Sentinel Solid Tire 1" [LFTIR005]	£127
Schwalbe Marathon Plus Evolution 1" [LFTIR010]	£167
<i>Available only on 24" Spoke</i>	

## Handrims

### Standard Handrims

*Handrims not available on 12" or 16" wheel sizes*

Aluminum Anodized [LFHRM001]	N/C
Plastic Coated [LFHRM002]	£139
Projection (8 Vertical Projections) [LFHRM003]	£308
<i>Available only with full profile tires in combination with 22" and 24" Spoke wheels, or 22" and 24" Mag wheels</i>	
Omit Handrim [LFHRM012]	N/C

## Standard Handrim Mount

*If no selection is made, Standard will be used*

Close Mount [LFHMT002]	N/C
<i>Reduces overall width by .75"</i>	

## Handrim Options

*Not available on Plastic Coated or Projection Handrims*

3M Micro Grip [LFACC006]	£106
--------------------------	------

## Axles

*No selection needed with Drum Brake*

Fixed [LFAXO002]	N/C
Quick Release [LFAXO003]	N/C

## Wheel Locks

### Wheel locks

Push to Lock [LFWLK002]	N/C
Pull to Lock [LFWLK003]	N/C
<i>Not available with 12" or 16" Full Poly tires</i>	
Hemi Wheel Lock - Left [LFWLK023]	£220
<i>Extension Handle included</i>	
Hemi Wheel Lock - Right [LFWLK024]	£220
<i>Extension Handle included</i>	

Drum Brake [LFWLK015]	£612
<i>Available only with Spoke wheels. Axle is Quick Release.</i>	
Dual - Drum Brake with Push to Lock [LFWLK016]	£751
<i>Available only with Spoke wheels. Axle is Quick Release.</i>	
Dual - Drum Brake with Pull to Lock [LFWLK017]	£751
<i>Available only with Spoke wheels. Axle is Quick Release.</i>	
Omit Wheel Locks [LFWLK008]	N/C

## Wheel Lock Options

Extension Handles [LFWLK004]	£75
<i>Available only with Push to Lock or Pull to Lock Wheel Locks</i>	
Grade Aids [LFWLK009]	£167
<i>Available only with Push to Lock Wheel Locks and Pneumatic or Pneumatic with Airless Insert tires</i>	

## Armrests

Armrest Type	Length	Price
Height Adjustable T-Arm Tall (11.5"-15")	Desk [LFARM039]	£244
	Full [LFARM040]	
Height Adjustable T-Arm (8.5"-13.5")	Desk [LFARM003]	£244
	Full [LFARM004]	
Angle Adjustable Locking Extendable Flip Up <sup>1</sup>	Desk [LFARM010]	£660
	Full [LFARM011]	
Omit Armrests	[LFARM007]	N/C

<sup>1</sup> Available only with 12" and 16" wheel sizes and Stroller Handle Back Post

## Arm Pads

Standard Arm Pad [LFARM026]	N/C
Waterfall Arm Pad [LFARM020]	£29
<i>Not available with Angle Adjustable Locking Extendable Flip Up Armrests</i>	
Foam Grip [LFARM025]	N/C
<i>Available only with Angle Adjustable Locking Extendable Flip Up Armrest</i>	

## Accessories

### Removable Side Guard

Composite Side Guards [LFSDG002]	£229
<i>Not available with T-Arms</i>	

### Pelvic Positioning Belt

1 1/2" Auto Buckle Pelvic Positioning Belt [LFPOB001]	£112
1 1/2" Padded Auto Buckle Pelvic Positioning Belt [LFPOB002]	£174
2" Airline Buckle Pelvic Positioning Belt [LFPOB003]	£122
Small Flex-Mount Positioning Belt (16-24") [LFPOB007]	£216
Medium Flex-Mount Positioning Belt (19-29") [LFPOB008]	£216
Large Flex-Mount Positioning Belt (25-37") [LFPOB009]	£216

### Belt Mounting Kit

Seat Belt Mounting Hardware [LFACC024]	£20
Belt Mounting Band Clamp [LFACC054]	£118
<i>Clamp includes hardware to mount Belt</i>	

### O2 Holder

O2 Holder [LFACC018]	£223
<i>Mounted to the left side. Not available with Fold Down Back.</i>	

### Upper Extremity Support Tray

*3/8" Polycarbonate with ABS ribs. Hardware ordered separately. BC = Body Cut Out. Armrest assembly is required for tray mounting.*

Tray	Fits Chair	Retail
Tray, Small 22" x 20" (BC 13" x 8") [LFTRY002]	10" - 14"	N/C
Tray, Medium 24" x 22" (BC 15" x 9") [LFTRY003]	12" - 16"	N/C
Tray, Large 26" x 24" (BC 17" x 10") [LFTRY004]	14" - 18"	N/C
Tray, X-Large 28" x 26" (BC 18" x 10") [LFTRY005]	16" - 20"	N/C
Tray Hardware Assy, Knob Release - Pair [LF123534]		N/C

## Frame Color

*Touch-up paint is available for purchase through the parts manual*

Black [LFFCS001]	N/C
Ruby Red [LFFCS063]	N/C
Radiant Orange [LFFCS056]	N/C
Rich Teal [LFFCS061]	N/C
Blue Steel [LFFCS051]	N/C
Deep Sapphire Blue [LFFCS065]	N/C
Purple Plum [LFFCS066]	N/C
Dark Pink [LFFCS053]	N/C
Shadow Black [LFFCS015]	N/C
Iridescent Navy [LFFCS016]	N/C
Bright Red [LFFCS062]	N/C
Tactical Green [LFFCS058]	N/C
Dry Earth Tan [LFFCS059]	N/C
Platinum [LFFCS050]	N/C
White [LFFCS060]	N/C

### Anti-Tip

Rear Anti-Tips [LFATP001]	£127
Omit [LFATP000]	N/C

### Transit Wheels

Transit wheels [LFACC035]	£348
<i>Not available with Long Depth Seat Frame, Seat Heights less than 14" or greater than 18", 12" or 16" rear wheels, Cane and Crutch Holder, and O2 Holder</i>	

### Other Accessories

Tool Kit [LFACC005]	£23
Removable Underseat Pouch [LFACC015]	£35
Ki Mobility Backpack [LFACC009]	£115
Black Neoprene Impact Guards [LFIMG001]	£94
Cane and Crutch Holder [LFACC007]	£216
IV Pole - Mounted Right Side [LFACC030]	£196
<i>Not available with Transit Option or 8° Bend with Push Handle Back Post</i>	
IV Pole - Mounted Left Side [LFACC031]	£196
<i>Not available with Transit Option or 8° Bend with Push Handle Back Post. Left side mount not available with O2 Holder.</i>	