



Etac UK 0(121)561-2222 | orders.uk@etac.com |
www.kimobility.com | www.etac.com



Company:

Ship To:

Account Number:

Date:

Address:

Quote: PO Number:

City:

State:

ATP:

Zip:

Phone:

Marked For:

Email/Fax:

Associated Quote/Order#:

Submit Form

Effective: 1/10/2024

Little Wave XP [TLXP] £2796

Little Wave Width Growable Frame 165 lb weight capacity Maximum Seat Width of 19"

Frame

Seat Width

10" [TLSWT010]	N/C
11" [TLSWT011]	N/C
12" [TLSWT012]	N/C
13" [TLSWT013]	N/C
14" [TLSWT014]	N/C
15" [TLSWT015]	N/C
16" [TLSWT016]	N/C

Seat Depth

10" [TLDSE010]	N/C
11" [TLDSE011]	N/C
12" [TLDSE012]	N/C
13" [TLDSE013]	N/C
14" [TLDSE014]	N/C
15" [TLDSE015]	N/C
16" [TLDSE016]	N/C

Sling Position

In Performance Position (0°), sling is at front of frame. The 1" or 2" sling position extends the frame in front of the sling by that amount. See reference C on page 7.

0" (Performance) [TLSP01]	N/C
1" [TLSP03]	N/C
2" (Classic) [TLSP02]	N/C

Transit Option

Transit System [TLTRNOPT]	N/C
125 lb maximum weight capacity	
Omit [TLTRNOP2]	N/C

Adjustable Depth

Adds the capability to grow seat in depth

No Rear Adjustment [TLAD05]	N/C
1/2" [TLAD06]	N/C
1" [TLAD01]	N/C
2" [TLAD02]	N/C
3" [TLAD03]	N/C

Long Sling

Available only with 1", 2" and 3" Adjustable Depth

Long Sling [TLAD04]	N/C
Sling will be maximum 2" longer than seat depth	

Seat Option

Standard Seat Sling [TLSS001]	N/C
Aluminum Solid Seat Pan [TLACC004]	£237
Seat Pan is 1" longer when Long Sling option is selected	

Front Frame Bend

Angle between front seat tube and ground. See reference D on page 7.

90° [TLFA01]	N/C
85° [TLFA03]	N/C
80° [TLFA02]	N/C
75° [TLFA04]	N/C
70° [TLFA05]	N/C

Backrest

Back Post / Height Range / Back Height							
Back	Price	Height	Back Height				
Height Adjustable	N/C	Short [TLBOP025]	11" ²	12"	13"	14"	
		Medium [TLBOP026]	13"	14"	15"	16"	17"
		Tall [TLBOP027]	17"	18"			
Fixed Height ¹	N/C	Fixed [TLBOP011]	9" ²	10" ²	11" ²	12" ²	13"
			14"	15"	16"	17"	18"
Offset ¹	N/C	Offset [TLBOP022]	11" ²	12" ²	13" ²	14" ²	15"
			16"	17"			

¹ Fixed Height. Adjustable Height Rigidizer Bar Included.

² Not available with Fold Down Push Handle

Offset Back Post Option

Selection required with Offset Back Posts. Not available with Height Adjustable or Fixed Height Back Posts. Each back post is offset 1" creating 2" overall offset.

2" Inward Offset [TLOBO001]	N/C
2" Outward Offset [TLOBO002]	N/C

Not available with Swing Away Armrests

Push Handle

Omit Push Handle [TLBCK008]	N/C
Standard Push Handle [TLBCK004]	N/C
<i>Not available on Offset Back Posts with 11" - 14" Back Heights</i>	
Bolt-On Push Handle [TLBOP005]	£182
<i>Handles will be 1" - 2" above back height</i>	
Fold Down Push Handle [TLBCK006]	£164

Stroller Handles

Stroller Handles are removable and height adjustable

Ergo Stroller Handle [TLBCK009]	£182
Single Stroller Handle [TLBOP007]	£164
Double Stroller Handle [TLBOP006]	£218

Not available with Push Handles or 10" Seat Width

Back Upholstery

Available only with Height Adjustable Back Posts. If selecting both upholstery and Axiom Back, the upholstery will be installed.

Padded Velcro Adjustable with Standard Cover [TLBUP001]	N/C
Padded Velcro Adjustable with Air Mesh Breathable Cover [TLBUP006]	£73
<i>Not available on chair widths less than 12"</i>	
Tension Adjustable with Standard Cover [TLBUP002]	N/C
Tension Adjustable with Air Mesh Breathable Cover [TLBUP007]	£182
<i>Not available on chair widths less than 12"</i>	
Omit [TLBUP003]	N/C

Rigidizer Bar

Adjustable Height Rigidizer Bar [TLBCK001]	N/C
<i>Included with Offset Back Posts</i>	
Non Adjustable Height Rigidizer Bar [TLBCK003]	N/C

Seat Back Angle

Standard is 90° (perpendicular to floor). Minus (-°) closes back angle. Plus (+°) opens back angle.

Standard [TLSBA001]	+3° [TLSBA002]	+6° [TLSBA003]	
+9° [TLSBA004]	-3° [TLSBA005]	-6° [TLSBA006]	
-9° [TLSBA007]	-12° [TLSBA008]		N/C

Front Seat Height and Casters

Caster Size / Footrest Taper / Front Seat Height

Select Front Seat Height and Footrest Taper in row corresponding to desired Caster Size. 6" Casters not available with Seat Depths less than 13".

Footrest Taper indicates bend of front frame creating footrest width. Taper is measured from outside of seat frame to outside of front tube. Taper is not available with 70° Front Frame Bend. See reference E on page 7.

Front Seat Height is measured from floor to top of seat tube at front of seat sling. Adjustable up or down by 1/2". See reference F on page 7.

Footrest Taper	Caster Size	Front Seat Height														
		11.5"	12"	12.5"	13"	13.5"	14"	14.5"	15"	15.5"	16"	16.5"	17"	17.5"	18"	18.5"
No Taper (Straight) [TLFT00]	3"	11.5"	12"	12.5"	13"	13.5"	14"	14.5"	15"	15.5"	16"	16.5"	17"	17.5"	18"	18.5"
	4"	12.5"	13"	13.5"	14"	14.5"	15"	15.5"	16"	16.5"	17"	17.5"	18"	18.5"	19"	19.5"
	5"	13.5"	14"	14.5"	15"	15.5"	16"	16.5"	17"	17.5"	18"	18.5"	19"	19.5"	20"	20.5"
	6"	14"	14.5"	15"	15.5"	16"	16.5"	17"	17.5"	18"	18.5"	19"	19.5"	20"	20.5"	
1" Taper [TLFT01]	3"	14.5"	15"	15.5"	16"	16.5"	17"	17.5"	18"	18.5"	19"	19.5"	20"	20.5"		
	4"	15.5"	16"	16.5"	17"	17.5"	18"	18.5"	19"	19.5"	20"	20.5"				
	5"	16.5"	17"	17.5"	18"	18.5"	19"	19.5"	20"	20.5"						
	6"	17"	17.5"	18"	18.5"	19"	19.5"	20"	20.5"							

Casters

		Caster Size				
		Price	3"	4"	5"	6"
Caster Type	.75" Rollerblade	N/C	[TLCAS002]	[TLCAS004]		
	.75" Lighted Rollerblade	£81	[TLCAS001]	[TLCAS005]	[TLCAS033]	
	1" Poly	N/C		[TLCAS006]	[TLCAS010]	[TLCAS014]
	1" Poly Aluminum	£97		[TLCAS007]	[TLCAS011]	
	1" Pneumatic	£103				[TLCAS017]
	1.5" Poly	N/C		[TLCAS032]	[TLCAS025]	[TLCAS026]
	1.5" Soft Roll Aluminum	£97	[TLCAS003]	[TLCAS009]	[TLCAS013]	[TLCAS016]

Fork Options

Standard Fork [TLFRK001]	N/C
Single Sided Fork [TLFRK002]	£121
<i>Not available with 6" Casters</i>	
Frog Legs II Suspension System [TLCST014]	£424
<i>Lightweight shock absorbing fork</i>	

Footrest

Footplates

Tubular "Open" Footplate [TLFPL004]	N/C
<i>Not available with footrest widths less than 9.5". See Reference G on page 7.</i>	
Tubular Footplate with ABS Cover [TLFPL005]	£86
<i>Not available with footrest widths less than 9.5". See Reference G on page 7.</i>	
Aluminum Angle Adjustable [TLFPL006]	£151
Aluminum Angle Adjustable with Riser [TLFPL007]	£151
High Mount Aluminum Angle Adjustable [TLFPL008]	£151
Angle Adjustable Flip Under [TLFPL012]	N/C
<i>Not available with footrest widths less than 9.5". See Reference G on page 7.</i>	
High Mount Angle Adjustable Flip Under [TLFPL011]	£272

Seat to Footrest Length

Measured from front edge of seat sling to top, rear of footrest. Footrest length at least 2.5" shorter than Front Seat Height recommended. See reference H on page 7.

5.5"	6"	6.5"	7"	7.5"	8"
8.5"	9"	9.5"	10"	10.5"	11"
11.5"	12"	12.5"	13"	13.5"	14"
14.5"	15"	15.5"	16"	16.5"	17"

Rear Seat Height and Wheels

Wheel Size / Rear Seat Height

Select Rear Seat Height in row corresponding to desired Wheel Size. For proper setup, selection is required when omitting rear wheels. Measured from floor to top of seat tube at front of back post. Rear Seat Height is custom to needs. See reference I on page 7.

Wheel Size	Rear Seat Height															
	11.5" ¹	12"	12.5"	13"	13.5"	14"	14.5"	15"	15.5"	16"	16.5"	17"	17.5"	18"	18.5"	19"
18"																
20"																
22"																
24"																

¹ Available only with Low Profile tires

Wheels

		Rear Wheel Size				
		Price	18"	20"	22"	24"
Rear Wheel	MAXX Spoke	N/C	[TLWH1118]	[TLWH1120]	[TLWH1122]	[TLWH1124]
	Mag	£121		[TLWH1220]	[TLWH1220]	[TLWH1224]
	Spinergy Spox	£665		[TLWHE420]	[TLWHE422]	[TLWHE424]

Spinergy Options

Select spoke color and handrim tab for Spinergy wheels

Black Spokes [TLSSC001]				N/C
Red Spokes [TLSSC002]	Blue Spokes [TLSSC005]	White Spokes [TLSSC003]		£91

Rear Tire Type

Full Poly 1 3/8" [TLTIR002]	N/C
Low Poly 1" [TLTIR003]	N/C
Pneumatic with Airless Insert [TLTIR004]	£110
Street Pneumatic 1 3/8" [TLTIR001]	N/C
High Pressure 1" [TLTIR006]	N/C
<i>Not available on Mag wheel or 18" wheel size</i>	
Schwalbe Marathon Plus Evolution 1" [TLTIR010]	£110
<i>Available only on 24" MAXX Spoke or Spinergy wheels</i>	

Handrims

Standard Handrims

Aluminum Anodized [TLHRM001]	N/C
Plastic Coated [TLHRM002]	£121
Projection (8 Vertical Projections) [TLHRM003]	£272
<i>Available only with full profile tires in combination with 22", 24", and 26" MAXX Spoke wheels, or 22" and 24" Mag wheels</i>	
Omit Handrim [TLHRM012]	N/C

Ergonomic Handrims

Not available with 20" or 26" wheels. Tab mount only available with 24" wheels.

Ergo Lite Regular Grip [TLHRM005]	£272
Ergo Standard Regular Grip [TLHRM008]	£272
Ergo Lite Extra Grip [TLHRM006]	£272
Ergo Standard Extra Grip [TLHRM009]	£272
Ergo Lite Omit Grip [TLHRM004]	£272
Ergo Standard Omit Grip [TLHRM007]	£272

Handrim Options

Handrim Mount

If no selection is made, Standard Mount will be used

Close Mount [TLTAB001]	N/C
------------------------	-----

Micro-Grip

3M Micro Grip [TLACC006]	£91
<i>Not available on Plastic Coated or Projection Handrims</i>	

Wheelchair Drive Ready

Chair will ship from Ki Mobility with adaption kit installed for selected Alber system. Order system separately from Alber USA. Select Rear Seat Height (full profile selections only) in drive wheel size on Rear Seat Height chart (page 3). Requires Rear Anti-Tips. Not available with 5th Wheel or Spinergy wheels.

e-motion Ready [TLWDR001]	£75
twion Ready [TLWDR002]	£75
0 Degree Camber [TLWDR003]	N/C
3 Degree Camber [TLWDR004]	N/C

Rear Wheel Spacing

Measured from outside of back to inside of rear tire. Adjustable out 1/2" from setting. 1.5" minimum Rear Wheel Spacing with Height Adjustable T-Arm Armrest, 1.25" minimum with Height Adjustable Pediatric Desk Armrest, 1" minimum with 0° Rear Wheel Camber. See reference K on page 7. No selection needed with Wheelchair Drive Ready option.

0.75" [TLWS075]	1" [TLWS100]	1.25" [TLWS125]	
1.5" [TLWS150]	1.75" [TLWS175]		N/C

Rear Wheel Camber

0° [TLWC00]	2° [TLWC02]	4° [TLWC04]	6° [TLWC06]	
8° [TLWC08]				N/C

Center of Gravity presets

Measured from front of back post to center of rear axle. 0° indicates axle will be directly under back post. See reference J on page 7.

Amputee Positions			
-1.5" [TLCGA150]	-1.25" [TLCGA125]	-1" [TLCGA100]	-.75" [TLCGA075]
-.5" [TLCGA050]	-.25" [TLCGA025]	0" [TLCGS000]	.25" [TLCGS025]
Standard Positions			
.5" [TLCGS050]	.75" [TLCGS075]	1" [TLCGS100]	1.25" [TLCGS125]
1.5" [TLCGS150]	1.75" [TLCGS175]	2" [TLCGS200]	2.25" [TLCGS225]
2.5" [TLCGS250]	2.75" [TLCGS275]	3" [TLCGS300]	3.25" [TLCGS325]
3.5" [TLCGS350]	3.75" [TLCGS375]	4" [TLCGS400]	4.25" [TLCGS425]

Axles

No selection needed with Drum Brake

Quick Release [TLAXO003]	N/C
Quad Quick Release [TLAXO001]	
<i>Not available with Mag or Spinergy wheels</i>	
	£182

Wheel Locks

Wheel Lock	Wheel Lock Options
Push to Lock [TLWLK002] N/C	Extension Handles [TLWLK004] £49 <i>Available only with Push to Lock or Pull to Lock Wheel Locks</i>
Pull to Lock [TLWLK003] N/C	Grade Aids [TLWLK009] £145 <i>Available only with Push to Lock Wheel Locks and Pneumatic or Pneumatic w/ Airless Insert tires</i>
Short Thro Scissor Lock [TLWLK006] £121	
Low Profile Scissor Lock [TLWLK012] £55	
<i>Not available with Low Poly or Full Poly tires</i>	
Drum Brake [TLWLK015] £532 <i>Available only with 20", 22" and 24" MAXX Spoke wheels. Not available with Fold Down Push Handles or Ergo and Single Stroller Handles without integrated push handles. Axle is Quick Release.</i>	
Dual - Drum Brake with Push to Lock [TLWLK016] £653 <i>Available only with 20", 22" and 24" MAXX Spoke wheels. Not available with Fold Down Push Handles or Ergo and Single Stroller Handles without integrated push handles. Axle is Quick Release.</i>	
Dual - Drum Brake with Pull to Lock [TLWLK017] £653 <i>Available only with 20", 22" and 24" MAXX Spoke wheels. Not available with Fold Down Push Handles or Ergo and Single Stroller Handles without integrated push handles. Axle is Quick Release.</i>	
Omit [TLWLK008] N/C	

Armrests

Armrest Type	Length	Price	Armrest Type	Price
Height Adjustable T-Arm Tall (11.5"-15") ^{1,3}	Desk [TLARM039]	£266	Height Adjustable Pediatric Desk (6.5" - 10") ^{7,8}	[TLARM014] £242
	Full [TLARM040]			
Height Adjustable T-Arm (8.5"-13.5") ^{1,3}	Desk [TLARM003]	£266	Tubular Flip Up ^{5,9,10,11}	[TLARM015] £61
	Full [TLARM004]			
Height Adjustable T-Arm Low (6.5"-10.5") ^{2,3,4}	Desk [TLARM022]	£266	Swing Away (8.5"-12.5") ^{9,11,12}	[TLARM002] £194
	Full [TLARM023]			
Angle Adjustable Locking Extendable Flip Up ^{5,6}	Desk [TLARM010]	N/C	Omit Armrests	[TLARM007] N/C
	Full [TLARM011]			

¹ Not available with Seat Depths less than 13"

² Not available with Seat Depths less than 12"

³ 1.5" minimum Rear Wheel Spacing

⁴ If requested configuration is unavailable, Standard Height T-Arm will be used

⁵ Requires Fixed Height Backrest and Adjustable Height Rigidizer Bar

⁶ Not available with back heights less than 12"

⁷ Best fit with Seat Depths less than 14"

⁸ 1.25" minimum Rear Wheel Spacing

⁹ Not available with Arm Pads

¹⁰ Not available with Back Heights less than 13"

¹¹ Not available with Offset

¹² Not available with Transit System AND more than 2" Adjustable Depth

Arm Pads

Not available with Tubular Flip Up or Swing Away Armrests. Full Arm Pads not available with Height Adjustable Pediatric Desk Armrests.

Standard Arm Pad [TLARM026]	N/C
Waterfall Arm Pad [TLARM020]	£30
<i>Not available with Angle Adjustable Locking Extendable Flip Up Armrests</i>	
Foam Grip [TLARM025]	N/C
<i>Available only with Angle Adjustable Locking Extendable Flip Up Armrest</i>	

Side Guards

Not available with T-Arms or Height Adjustable Pediatric Desk Armrests

Direct Mount Side Guards	Clamp-On Mount Side Guards
Composite Side Guards - Adult - Direct Mount Receiver [TLSDG002] £140 <i>12" minimum Seat Depth</i>	Aluminum Adjustable Fender w/Foam Pad - Clamp-On Receiver - Black [TLSDG017] £296
Composite Side Guards - Pediatric - Direct Mount Receiver [TLSDG001] £140	Aluminum Adjustable Fender w/Foam Pad - Clamp-On Receiver [TLSDG018] £357 <i>Not available with Hydrographics</i>

Accessories

5th Wheel

Rotating removable rear caster. See reference L on page 7.

Standard 5th Wheel [TLATP002]	£206
Dynamic 5th Wheel [TLATP004]	£206
Soft Elastomer [TLELA001]	N/C
<i>Both elastomers are included; the selected elastomer will be installed</i>	
Firm Elastomer [TLELA002]	N/C
<i>Both elastomers are included; the selected elastomer will be installed</i>	

Pelvic Positioning Belt

1" Auto Buckle Pelvic Positioning Belt [TLPOB004]	£55
1" Padded Auto Buckle Pelvic Positioning Belt [TLPOB005]	£151
1 1/2" Auto Buckle Pelvic Positioning Belt [TLPOB001]	£55
1 1/2" Padded Auto Buckle Pelvic Positioning Belt [TLPOB002]	£151
Small Bodypoint Evoflex (16-24") [TLPOB007]	£188
Medium Bodypoint Evoflex (19-29") [TLPOB008]	£188
Large Bodypoint Evoflex (25-37") [TLPOB009]	£188

Belt Mounting Kit

Ki Mobility Belt Mounting Clamp [TLACC023]	£37
Bodypoint Belt Mounting Bracket Kit [TLACC013]	£103

Anti-Tip

Rear Anti-Tips [TLATP001]	£212
Omit [TLATP000]	N/C

Calf Strap

Velcro Adjustable [TLACC002]	N/C
Padded Velcro Adjustable [TLACC003]	£67
Bodypoint Padded AeroMesh [TLACC010]	£103

Other Accessories

Tool Kit [TLACC005]	£20
Removable Underseat Pouch [TLACC015]	£30
Pediatric Ki Mobility Backpack [TLACC026]	£100
Ki Mobility Backpack [TLACC009]	£100
<i>Best fit on chairs more than 10" wide</i>	
Luggage Carrier [TLACC025]	£73
Black Neoprene Impact Guards [TLIMG001]	£86
Clear Spoke Guards [TLSPG001]	
<i>Not available with Mag wheel</i>	
Black Spoke Guards [TLSPG002]	
<i>Not available with 18" or 20" wheel sizes or Mag wheel</i>	
Cane and Crutch Holder [TLACC007]	£91

Frame Color

Gloss Finish

Black [TLFCS001]	N/C
Blue Skies [TLFCS018]	N/C
Bubblegum Pink [TLFCS011]	N/C
Candy Blue [TLFCS012]	N/C
Candy Purple [TLFCS010]	N/C
Candy Red [TLFCS003]	N/C
Copper [TLFCS009]	N/C
Granny Smith [TLFCS008]	N/C
Metallic Orange [TLFCS032]	N/C
Yellow [TLFCS029]	N/C

Matte Finish

Military Green [TLFCS025]	N/C
Charcoal Gray [TLFCS031]	N/C
Emerald Green [TLFCS019]	N/C
Iridescent Burgundy [TLFCS021]	N/C
Iridescent Navy [TLFCS016]	N/C
Retro Red [TLFCS022]	N/C
Royal Blue [TLFCS023]	N/C
Shadow Black [TLFCS015]	N/C
Smooth White [TLFCS024]	N/C

References

A. Seat Width *See Diagram 1*

Measured from outside of frame tube on one side to outside of frame tube on other side.

B. Seat Depth *See Diagram 3*

Measured from front of back posts to front edge of seat sling. Seat sling starts at beginning of bend at front of frame.

C. Sling Position *See Diagram 2*

In Performance Position (0°), sling is at front of frame to keep chair shorter and more maneuverable. The 1" or 2" seat sling position extends the frame in front of the sling by that amount. Adding more frame can improve stability and provide support to aid in transfers.

D. Front Frame Bend *See Diagram 3*

Angle between front seat tube and ground.

E. Footrest Taper *See Diagram 1*

Indicates bend of front frame creating footrest. Taper is measured from outside of seat frame to outside of front tube.

F. Front Seat Height *See Diagram 1*

Measured from floor to top of seat tube at front of seat sling.

G. Footrest Width *See Diagram 1*

Measured from inside of front tube to inside of front tube on other side. Footrest width is listed below for each seat width and taper.

Seat Width Outside Measurement	Straight	1" Taper
	2" Narrower than Seat Width	4" Narrower than Seat Width
10"	7.5"	
11"	8.5"	
12"	9.5"	7.5"
13"	10.5"	8.5"
14"	11.5"	9.5"
15"	12.5"	10.5"
16"	13.5"	11.5"

H. Seat to Footrest Length *See Diagram 3*

Measured from front edge of seat sling to top rear of footrest. Footrest length of at least 2.5" shorter than front seat height recommended.

I. Rear Seat Height *See Diagram 3*

Measured from floor to top of seat tube at front of back post. Rear seat height is custom to needs.

J. Center of Gravity Preset *See Diagram 3*

Measured from front of back post to center of rear axle. 0° indicates axle will be directly under back post.

K. Rear Wheel Spacing *See Diagram 1*

Measured from outside of seat back to inside of rear tire. Adjustable out .5" from setting.

L. 5th Wheel *See Diagram 4*

Rotating removable 5th wheel allows for aggressive center of gravity positioning of the rear wheels.

Diagram 1




Diagram 2




Diagram 3




Diagram 4



Axiom G

		Width															
		14"		15"		16"		17"		18"		19"		20"		22"	
		Cushion	Outer Cover	Cushion	Outer Cover	Cushion	Outer Cover	Cushion	Outer Cover	Cushion	Outer Cover	Cushion	Outer Cover	Cushion	Outer Cover	Cushion	Outer Cover
Retail		£85	£70	£85	£70	£85	£70	£85	£70	£85	£70	£85	£70	£85	£70	£85	£70
Height	14"																
	15"																
	16"																
	17"																
	18"																
	19"																
	20"																
	21"																
	22"																

Axiom S

		Width															
		14"		15"		16"		17"		18"		19"		20"		22"	
		Cushion	Outer Cover	Cushion	Outer Cover	Cushion	Outer Cover	Cushion	Outer Cover	Cushion	Outer Cover	Cushion	Outer Cover	Cushion	Outer Cover	Cushion	Outer Cover
Retail		£120	£70	£120	£70	£120	£70	£120	£70	£120	£70	£120	£70	£120	£70	£120	£70
Height	14"																
	15"																
	16"																
	17"																
	18"																
	19"																
	20"																
	21"																
	22"																

Axiom P

		Width															
		14"		15"		16"		17"		18"		19"		20"		22"	
		Cushion	Outer Cover	Cushion	Outer Cover	Cushion	Outer Cover	Cushion	Outer Cover	Cushion	Outer Cover	Cushion	Outer Cover	Cushion	Outer Cover	Cushion	Outer Cover
Retail		£170	£70	£170	£70	£170	£70	£170	£70	£170	£70	£170	£70	£170	£70	£170	£70
Height	14"																
	15"																
	16"																
	17"																
	18"																
	19"																
	20"																
	21"																
	22"																

Axiom SP - Visco

		Width																										
		14"			15"			16"			17"			18"			19"			20"			21"			22"		
		Cushion	Cover		Cushion	Cover		Cushion	Cover		Cushion	Cover		Cushion	Cover		Cushion	Cover		Cushion	Cover		Cushion	Cover				
Retail		£350	£70	£350	£70	£350	£70	£350	£70	£350	£70	£350	£70	£350	£70	£350	£70	£350	£70	£350	£70	£350	£70					
Height	14"																											
	15"																											
	16"																											
	17"																											
	18"																											
	19"																											
	20"																											
	21"																											
	22"																											

Axiom SP - Fluid

		Width																										
		14"			15"			16"			17"			18"			19"			20"			21"			22"		
		Cushion	Cover		Cushion	Cover		Cushion	Cover		Cushion	Cover		Cushion	Cover		Cushion	Cover		Cushion	Cover		Cushion	Cover				
Retail		£320	£70	£320	£70	£320	£70	£320	£70	£320	£70	£320	£70	£320	£70	£320	£70	£320	£70	£320	£70	£320	£70					
Height	14"																											
	15"																											
	16"																											
	17"																											
	18"																											
	19"																											
	20"																											
	21"																											
	22"																											


Solid Seat Insert

Weight capacity: 350 lbs

		Width								
		14"	15"	16"	17"	18"	19"	20"	21"	22"
		Cushion	Cushion	Cushion	Cushion	Cushion	Cushion	Cushion	Cushion	Cushion
Retail		£55	£55	£55	£55	£55	£55	£55	£55	£55
Height	14"									
	15"									
	16"									
	17"									
	18"									
	19"									
	20"									
	21"									
	22"									

Axiom Back Width	Fits Back Frame Widths	
	Inner Slots	Outer Slots
12" (30.5cm)	11" (28.0cm)	12" -14" (30.5cm -35.5cm)
14" (35.5cm)	12" -13" (30.5cm -33.0cm)	14" -16" (35.5cm -40.5cm)
16" (40.5cm)	14" -15" (35.5cm -38.0cm)	16" -18" (40.5cm -45.5cm)
18" (45.5cm)	16" - 17" (40.5cm - 43.0cm)	18" - 20" (45.5cm - 50.5cm)
20" (51.0cm)	18" -19" (46.0cm -48.5cm)	20" -22" (50.5cm -55.5cm)
22" (56.0cm)	20" -21" (51.0cm -53.5cm)	22" -24" (56.0cm -61.0cm)

Hardware Description	Code
Improved DEX Hardware, removable adjustable hardware (300 lb max user weight)	DEX
New BTX Hardware, removable adjustable hardware not recommended for, or installed on Ultralight Rigid Frame Wheelchairs (400 lb max user weight)	BTX
New Fixed Hardware, non-removeable hardware up to (300 lb max user weight)	FXT



Axiom Back

Hardware is included with back. Low Back selections are not available with headrests, modifications, or posture supports. Weight capacity: 300 lbs.

Axiom Posterior Back (PB)

		Width						
		12"	14"	16"	18"	20"	22"	
Retail		£550	£550	£550	£550	£550	£550	
Height	Low	8"	[AXPB1208]	[AXPB1408]	[AXPB1608]	[AXPB1808]	[AXPB2008]	
		10"		[AXPB1410]	[AXPB1610]	[AXPB1810]	[AXPB2010]	[AXPB2210]
		12"	[AXPB1212]	[AXPB1412]	[AXPB1612]	[AXPB1812]	[AXPB2012]	
	Mid	14"		[AXPB1414]	[AXPB1614]	[AXPB1814]	[AXPB2014]	[AXPB2214]
		16"	[AXPB1216]	[AXPB1416]	[AXPB1616]	[AXPB1816]	[AXPB2016]	
		18"		[AXPB1418]	[AXPB1618]	[AXPB1818]	[AXPB2018]	[AXPB2218]

Back Hardware

DEX Hardware [DEX] BTX Hardware [BTX] FXT Hardware [FXT] _____ N/C

Lateral Support for AXPB Backs

Lateral Trunk Supports for AXPB Backs							
Pair	Left	Right	Description	Item Code	Pair	Each	
			3 x 4 Lateral Pads with Fixed Mounting Hardware for Axiom Posterior Back	PBLF3x4P	£295	£148	
			4 x 5 Lateral Pads with Fixed Mounting Hardware for Axiom Posterior Back	PBLF4x5P	£295	£148	
			3 x 4 Lateral Pads with Swing Away Mounting Hardware for Axiom Posterior Back	PBLS3x4P	£628	£314	
			4 x 5 Lateral Pads with Swing Away Mounting Hardware for Axiom Posterior Back	PBLS4x5P	£628	£314	
			3 x 4 Lateral Pads with Adjustable Swing Away Mounting Hardware	PBLAS3x4P	£628	£314	
			4 x 5 Lateral Pads with Adjustable Swing Away Mounting Hardware	PBLAS4x5P	£628	£314	

Axiom Posterior Lateral Back (PL)

		Width						
		12"	14"	16"	18"	20"	22"	
Retail		£550	£550	£550	£550	£550	£550	
Height	Low	8"	[AXPL1208]	[AXPL1408]	[AXPL1608]	[AXPL1808]	[AXPL2008]	
		10"		[AXPL1410]	[AXPL1610]	[AXPL1810]	[AXPL2010]	[AXPL2210]
		12"	[AXPL1212]	[AXPL1412]	[AXPL1612]	[AXPL1812]	[AXPL2012]	
	Mid	14"		[AXPL1414]	[AXPL1614]	[AXPL1814]	[AXPL2014]	[AXPL2214]
		16"	[AXPL1216]	[AXPL1416]	[AXPL1616]	[AXPL1816]	[AXPL2016]	
		18"		[AXPL1418]	[AXPL1618]	[AXPL1818]	[AXPL2018]	[AXPL2218]
Tall	20"	[AXPL1220]	[AXPL1420]	[AXPL1620]	[AXPL1820]	[AXPL2020]		
	22"		[AXPL1422]	[AXPL1622]	[AXPL1822]	[AXPL2022]	[AXPL2222]	

Back Hardware

DEX Hardware [DEX] BTX Hardware [BTX] FXT Hardware [FXT] _____ N/C

Lateral Support for AXPL

Lateral Trunk Supports for AXPL & AXLB Backs						
Pair	Left	Right	Description	Item Code	Pair	Each
			3 x 4 Lateral Pads with Fixed Mounting Hardware for Axiom Posterior Lateral Back	PLLF3x4P	£295	£295
			4 x 5 Lateral Pads with Fixed Mounting Hardware for Axiom Posterior Lateral Back	PLLF4x5P	£295	£295
			3 x 4 Lateral Pads with Swing Away Mounting Hardware for Axiom Posterior Lateral Back	PLLS3x4P	£628	£314
			4 x 5 Lateral Pads with Swing Away Mounting Hardware for Axiom Posterior Lateral Back	PLLS4x5P	£628	£314
			3 x 4 Lateral Pads with Adjustable Swing Away Mounting Hardware	PLLAS3x4P	£628	£314
			4 x 5 Lateral Pads with Adjustable Swing Away Mounting Hardware	PLLAS4x5P	£628	£314

Axiom Deep Lateral Back (DL)

		Width						
		12"	14"	16"	18"	20"	22"	
Retail		£550	£550	£550	£550	£550	£550	
HCPC		E2620	E2620	E2620	E2620	E2620	E2621	
Height	Mid	14"	[AXDL1214]	[AXDL1414]	[AXDL1614]	[AXDL1814]	[AXDL2014]	[AXDL2214]
		16"		[AXDL1416]	[AXDL1616]	[AXDL1816]	[AXDL2016]	
		18"	[AXDL1218]	[AXDL1418]	[AXDL1618]	[AXDL1818]	[AXDL2018]	[AXDL2218]
	Tall	20"	[AXDL1220]	[AXDL1420]	[AXDL1620]	[AXDL1820]	[AXDL2020]	
		22"		[AXDL1422]	[AXDL1622]	[AXDL1822]	[AXDL2022]	[AXDL2222]

Back Hardware

DEX Hardware [DEX] BTX Hardware [BTX] FXT Hardware [FXT] _____ N/C

Headrest

Description	Std. Straight Removable Adjustable Hardware		Multi-Axis Removable Adjustable Hardware		Posalinc Removable Adjustable Hardware		HCPCS
	Item Code	Retail	Item Code	Retail	Item Code	Retail	
Comfort Soft Headrest 8"	[CS08SS]	£268	[CS08MA]	£268	[CS08PL]	£268	E0955 / E1028
Comfort Soft Headrest 10"	[CS10SS]	£268	[CS10MA]	£268	[CS10PL]	£268	E0955 / E1028
Comfort Soft Headrest 14"	[CS14SS]	£268	[CS14MA]	£268	[CS14PL]	£268	E0955 / E1028
Contour Positioning Pad, 8" X 6.25"	[CP08SS]	£303	[CP08MA]	£303	[CP08PL]	£303	E0955 / E1028
Contour Positioning Pad, 9" X 6.5"	[CP09SS]	£303	[CP09MA]	£303	[CP09PL]	£303	E0955 / E1028
Contour Positioning Pad, 10" X 7.5"	[CP10SS]	£303	[CP10MA]	£303	[CP10PL]	£303	E0955 / E1028
Tri-soft Headrest 10"	[TS10SS]	£386	[TS10MA]	£386	[TS10PL]	£386	E0955 / E1028
Tri-soft Headrest 14"	[TS14SS]	£386	[TS14MA]	£386	[TS14PL]	£386	E0955 / E1028

Description	Item Code	Retail	HCPCS
Universal Headrest Mounting Plate	[UMP]	£43	K0108

Not required for Axiom Headrests

Anterior Thoracic Supports

Selection of a modification from the Axiom Back Modifications section required in order to mount posture supports. Not available with Low Backs.

Anterior Thoracic Supports			
Description	Item Code	HCPCS	Retail
Dynaform Chest Support: 2-Piece with Hook&Loop, 3" High Pad, Small	[14941]	E0960	£100
Dynaform Chest Support: 2-Piece with Hook&Loop, 4" High Pad, Medium	[14942]	E0960	£100
Dynaform Chest Support: 2-Piece with Hook&Loop, 5" High Pad, Large	[14943]	E0960	£100
Dynaform Chest Support: 1-Piece with Buckle, 3" Height x 20" Width, Small	[14924]	E0960	£100
Dynaform Chest Support: 1-Piece with Buckle, 4" Height x 20" Width, Medium	[14927]	E0960	£100
Dynaform Chest Support: 1-Piece with Buckle, 4" Height x 22" Width, Large	[14928]	E0960	£100
Dynaform Anterior Trunk Support: Standard, Dynamic, X-Small	[14851]	E0960	£100
Dynaform Anterior Trunk Support: Standard, Dynamic, Small	[14852]	E0960	£100
Dynaform Anterior Trunk Support: Standard, Dynamic, Medium	[14853]	E0960	£100
Dynaform Anterior Trunk Support: Standard, Dynamic, Large	[14854]	E0960	£100
Dynaform Anterior Trunk Support: Slim Cut, Dynamic, X-Small	[14856]	E0960	£100
Dynaform Anterior Trunk Support: Slim Cut, Dynamic, Small	[14857]	E0960	£100
Dynaform Anterior Trunk Support: Slim Cut, Dynamic, Medium	[14858]	E0960	£100
Dynaform Anterior Trunk Support: Slim Cut, Dynamic, Large	[14859]	E0960	£100
AirLogic Anterior Trunk Support: Standard, Structure, Small	[14282]	E0960	£300
AirLogic Anterior Trunk Support: Standard, Structure, Medium	[14283]	E0960	£300
AirLogic Anterior Trunk Support: Standard, Structure, Large	[14284]	E0960	£300
AirLogic Anterior Trunk Support: Standard, Structure, X-Large	[14285]	E0960	£300
AirLogic Anterior Trunk Support: Slim Cut, Dynamic, Small	[14257]	E0960	£300
AirLogic Anterior Trunk Support: Slim Cut, Dynamic, Medium	[14258]	E0960	£300
AirLogic Anterior Trunk Support: Slim Cut, Dynamic, Large	[14259]	E0960	£300
AirLogic Anterior Trunk Support: Slim Cut, Dynamic, X-Large	[14260]	E0960	£300

Accessories

Adjustable Shoulder Guides (K0108) [AXBSG] <i>Includes Anterior Thoracic Mounting Modification. Not available with Low Backs.</i>	£43
Extra Lumbar Pad (K0108) [LUMPAD] <i>Available only with Axiom Lateral Back and Axiom Deep Lateral back</i>	N/C
Privacy Flap (K0108) [PRFLAP]	£44

Axiom Back Modifications

Axiom Back Modifications

Extended lead times will apply for all modifications. Makes back non-refundable.

Anterior Thoracic Mounting Modification [AXSMM] <i>Provides holes for Anterior Thoracic Supports and Shoulder Guides</i>	N/C
---	-----

Gloss Finish

If no selection is made, standard Shadow Black will be used with no additional lead time

Black Back Shell [AXMBLK]	N/C
Blue Skies Back Shell [AXMSKY]	N/C
Bubblegum Pink Back Shell [AXMPNK]	N/C
Candy Blue Back Shell [AXMCBL]	N/C
Candy Purple Back Shell [AXMPUR]	N/C
Candy Red Back Shell [AXMRED]	N/C
Copper Back Shell [AXMCOP]	N/C
Granny Smith Back Shell [AXMGSM]	N/C
Metallic Orange Back Shell [AXMMOR]	N/C
Yellow Back Shell [AXMYLW]	N/C
Silver Back Shell [AXMSIL]	N/C
Teal Back Shell [AXMTEA]	N/C

Matte Finish

If no selection is made, standard Shadow Black will be used with no additional lead time

Military Green Back Shell [AXMAGM]	N/C
Charcoal Gray Back Shell [AXMCGM]	N/C
Emerald Green Back Shell [AXMEGM]	N/C
Iridescent Burgundy Back Shell [AXMIBM]	N/C
Iridescent Navy Back Shell [AXMINM]	N/C
Retro Red Back Shell [AXMRRM]	N/C
Royal Blue Back Shell [AXMRBM]	N/C
Smooth White Back Shell [AXMSWM]	N/C