



Etac UK 0(121)561-2222 | orders.uk@etac.com |
www.kimobility.com | www.etac.com



Company:

Ship To:

Account Number:

Date:

Address:

Quote: PO Number:

City:

State:

ATP:

Zip:

Phone:

Marked For:

Email/Fax:

Associated Quote/Order#:

Submit Form

Effective: 1/10/2024

Rogue2 **£2796**
Rogue Standard Frame 275 lb Weight Capacity

Frame

Seat Width

See reference A on page 7

12" [RGSWT012]	N/C
13" [RGSWT013]	N/C
14" [RGSWT014]	N/C
15" [RGSWT015]	N/C
16" [RGSWT016]	N/C
17" [RGSWT017]	N/C
18" [RGSWT018]	N/C
19" [RGSWT019]	N/C
20" [RGSWT020]	£350

Seat Depth

See reference B on page 7

14" [RGSDE014]	N/C
15" [RGSDE015]	N/C
16" [RGSDE016]	N/C
17" [RGSDE017]	N/C
18" [RGSDE018]	N/C
19" [RGSDE019]	N/C
20" [RGSDE020]	£350

Transit Option

Transit System [RGTRNOPT]	£200
<i>Cannot be retrofitted. Requires Expanded Width Caster Position.</i>	
Omit [RGTRNOP2]	N/C

Adjustable Depth

Adds the capability to grow seat in depth

No Rear Adjustment [RGAD05]	N/C
1/2" [RGAD06]	N/C
1" [RGAD01]	£115
2" [RGAD02]	£115

Sling Position

In Performance Position (0"), sling is at front of frame. The 1" or 2" sling position extends the frame in front of the sling by that amount. See reference C on page 7.

0" (Performance) [RGSP01]	N/C
1" [RGSP03]	N/C
2" (Classic) [RGSP02]	N/C

Seat Option

Standard Seat Sling [RGSSS001]	N/C
Aluminum Solid Seat Pan [RGACC004]	£288

Front Frame Bend

Angle between front seat tube and ground. See reference D on page 7.

90° [RGFA01]	N/C
85° [RGFA03]	N/C
80° [RGFA02]	N/C
75° [RGFA04]	£173
70° [RGFA05]	£173

Backrest

Back Post / Height Range / Back Height

Back Post	Price	Height Range	Back Height				
			11" ²	12"	13"	14"	
Height Adjustable	N/C	Short [RGBOP025]					
		Medium [RGBOP026]	13"	14"	15"	16"	17"
		Tall [RGBOP027]	17"	18"	19"	20"	
Fixed Height ¹	N/C	Fixed [RGBOP011]	9" ²	10" ²	11" ²	12"	13"
			14"	15"	16"	17"	18"

¹ Not available with back upholstery

² Not available with Fold Down Push Handle

Push Handles

Omit Push Handle [RGBCK008]	N/C
Standard Push Handle [RGBCK004]	N/C
Low Profile Push Handle [RGBCK007]	£173
<i>2" shorter than standard push handles</i>	
Bolt-On Push Handle [RGBOP005]	£184
<i>Handles will be 1" - 2" above back height</i>	
Fold Down Push Handle [RGBCK006]	£173

Stroller Handle

Stroller Handles are removable and height adjustable

Single Stroller Handle [RGBOP007]	£184
Double Stroller Handle [RGBOP006]	£230

Not available with Push Handles

Rigidizer Bar

- Adjustable Height Rigidizer Bar [RGBCK001] _____ N/C
- Non-Adjustable Height Rigidizer Bar [RGBCK010] _____ N/C
- Not Available with Fixed Height Back Post*

Seat Back Angle

Standard is 90° (perpendicular to floor). Minus (-°) closes back angle. Plus (+°) opens back angle.

- Standard [RGSBA001] +3° [RGSBA002] +6° [RGSBA003]
- +9° [RGSBA004] -3° [RGSBA005] -6° [RGSBA006] -9°
- [RGSBA007]
- 12° [RGSBA008] _____ N/C

Back Upholstery Options

Available only with Height Adjustable Back Posts. If selecting both upholstery and Axiom Back, the upholstery will be installed.

- Padded Velcro Adjustable with Standard Cover [RGBUP001] _____ N/C
- Padded Velcro Adjustable with Air Mesh Breathable Cover [RGBUP006] _____ £92
- Tension Adjustable with Standard Cover [RGBUP002] _____ N/C
- Tension Adjustable with Air Mesh Breathable Cover [RGBUP007] _____ £196
- Omit [RGBUP003] _____ N/C

Front Seat Height and Casters

Caster Size / Footrest Taper / Front Seat Height

Select Front Seat Height and Footrest Taper in row corresponding to Caster Size. 6" casters not available with 85° and 90° Front Frame Bend. Footrest Taper indicates bend of front frame creating footrest width. Taper is measured from outside of seat frame to outside of front tube. 2" Taper not available with 75° Front Frame Bend. Taper is not available with 70° Front Frame Bend. See reference E on page 7. Front Seat Height is measured from floor to top of seat tube at front of seat sling. Adjustable up or down by 1/2". See reference F on page 7.

Footrest Taper	Price	Caster Size	Front Seat Height														
			13.5"	14"	14.5"	15"	15.5"	16"	16.5"	17"	17.5"	18"	18.5"	19"	19.5"	20"	20.5"
Straight [RGFT00]	N/C	3"	13.5"	14"	14.5"	15"	15.5"	16"	16.5"	17"	17.5"	18"	18.5"	19"	19.5"	20"	20.5"
		4"	14"	14.5"	15"	15.5"	16"	16.5"	17"	17.5"	18"	18.5"	19"	19.5"	20"	20.5"	
		5"	14.5"	15"	15.5"	16"	16.5"	17"	17.5"	18"	18.5"	19"	19.5"	20"	20.5"		
		6"	15"	15.5"	16"	16.5"	17"	17.5"	18"	18.5"	19"	19.5"	20"	20.5"			
1" Taper [RGFT01]	£60	3"	16.5"	17"	17.5"	18"	18.5"	19"	19.5"	20"	20.5"						
		4"	17"	17.5"	18"	18.5"	19"	19.5"	20"	20.5"							
		5"	17.5"	18"	18.5"	19"	19.5"	20"	20.5"								
		6"	18"	18.5"	19"	19.5"	20"	20.5"									
2" Taper [RGFT02]	£60	3"	16.5"	17"	17.5"	18"	18.5"	19"	19.5"	20"	20.5"						
		4"	17"	17.5"	18"	18.5"	19"	19.5"	20"	20.5"							
		5"	17.5"	18"	18.5"	19"	19.5"	20"	20.5"								
		6"	18"	18.5"	19"	19.5"	20"	20.5"									
3.5° V-Front [RGFT05]	£60	3"	13.5"	14"	14.5"	15"	15.5"	16"	16.5"	17"	17.5"	18"	18.5"	19"	19.5"	20"	20.5"
		4"	14"	14.5"	15"	15.5"	16"	16.5"	17"	17.5"	18"	18.5"	19"	19.5"	20"	20.5"	
		5"	14.5"	15"	15.5"	16"	16.5"	17"	17.5"	18"	18.5"	19"	19.5"	20"	20.5"		
		6"	15"	15.5"	16"	16.5"	17"	17.5"	18"	18.5"	19"	19.5"	20"	20.5"			
7° V-Front [RGFT06]	£60	3"	13.5"	14"	14.5"	15"	15.5"	16"	16.5"	17"	17.5"	18"	18.5"	19"	19.5"	20"	20.5"
		4"	14"	14.5"	15"	15.5"	16"	16.5"	17"	17.5"	18"	18.5"	19"	19.5"	20"	20.5"	
		5"	14.5"	15"	15.5"	16"	16.5"	17"	17.5"	18"	18.5"	19"	19.5"	20"	20.5"		
		6"	15"	15.5"	16"	16.5"	17"	17.5"	18"	18.5"	19"	19.5"	20"	20.5"			

Casters

Caster Type		Caster Size				
		Price	3"	4"	5"	6"
Caster Type	.75" Rollerblade	£104	[RGCAS002]	[RGCAS004]		
	.75" Lighted Rollerblade	£104	[RGCAS001]	[RGCAS005]		
	1" Poly	N/C		[RGCAS006]	[RGCAS010]	[RGCAS014]
	1" Poly Aluminum	£104		[RGCAS007]	[RGCAS011]	
	1" Pneumatic	£104				[RGCAS017]
	1.5" Poly	N/C		[RGCAS032]	[RGCAS025]	[RGCAS026]
	1.5" Soft Roll Aluminum	£104	[RGCAS003]	[RGCAS009]	[RGCAS013]	[RGCAS016]

Caster Position

See reference L on page 7

Optimized [RGCST015]	N/C
<i>Caster Placement minimizes the footprint of the chair</i>	
Expanded Width [RGCST016]	N/C
<i>Caster arm is positioned approximately .75" wider than Optimized position</i>	

Fork Options

Standard Fork [RGFRK001]	N/C
Single Sided Fork [RGFRK003]	£230
Frog Legs II Suspension System [RGCST014]	£518
<i>Lightweight shock absorbing fork. Not included in Colored Anodizing Package.</i>	

Footrest

Footplate

Tubular "Open" Footplate [RGFPL004]	N/C
<i>Not available with V-Taper</i>	
Tubular Footplate with ABS Cover [RGFPL005]	£104
<i>Not available w/footrest widths greater than 15.5". See Reference G on page 7. Not available with V-Taper.</i>	
Aluminum Angle Adjustable [RGFPL006]	£207
Aluminum Angle Adjustable with Riser [RGFPL007]	£207
High Mount Aluminum Angle Adjustable [RGFPL008]	£207
<i>Not available with V-Taper</i>	
Angle Adjustable Flip Under [RGFPL012]	£230
<i>Not available with footrest widths less than 9.5". Not available with V-Taper. See Reference G on page 7.</i>	
High Mount Angle Adjustable Flip Under [RGFPL011]	£288
<i>Not available w/footrest widths greater than 15.5". See Reference G on page 7. Not available with V-Taper.</i>	
Hybrid Angle Adjustable [RGFPL015]	£115
Hybrid Angle Adjustable with Wings [RGFPL016]	£115

Footrest Protector

Footrest Protector [RGACC020]	£23
<i>Not available with Angle Adjustable Flip Under, High Mount Angle Adjustable Flip Under or High Mount Aluminum Angle Adjustable Footplates</i>	

Seat to Footrest Length

Measured from front edge of seat sling to top, rear of footrest. Footrest length at least 2.5" shorter than Front Seat Height recommended. See reference H on page 7.

8"	8.5"	9"	9.5"	10"	10.5"
11"	11.5"	12"	12.5"	13"	13.5"
14"	14.5"	15"	15.5"	16"	16.5"
17"	17.5"	18"	18.5"	19"	19.5"

Rear Seat Height and Wheels

Wheel Size / Rear Seat Height

Select Rear Seat Height in row corresponding to desired Wheel Size. For proper setup, selection is required when omitting rear wheels. Measured from floor to top of seat tube at front of back post. Rear Seat Height is custom to needs. See reference I on page 7.

Wheel Size	Rear Seat Height															
	12.5" ¹	13"	13.5"	14"	14.5"	15"	15.5"	16"	16.5"	17"	17.5"	18"	18.5"	19"	19.5"	20"
20"	12.5" ¹	13"	13.5"	14"	14.5"	15"	15.5"	16"	16.5"	17"	17.5"	18"	18.5"	19"	19.5"	20"
22"	13.5" ¹	14"	14.5"	15"	15.5"	16"	16.5"	17"	17.5"	18"	18.5"	19"	19.5"	20"		
24"	14.5" ¹	15"	15.5"	16"	16.5"	17"	17.5"	18"	18.5"	19"	19.5"	20"				
25"	15" ¹	15.5"	16"	16.5"	17"	17.5"	18"	18.5"	19"	19.5"	20"					
26"	15.5" ¹	16"	16.5"	17"	17.5"	18"	18.5"	19"	19.5"	20"						

¹ Available only with Low Profile tires

One Arm Drive

One Arm Drive - Quick Release [RGDRV007]	£1150
<i>Available only with 20", 22" and 24" spoke wheels and quick release axles. Spoke wheel style provided is exclusive to OAD. Standard handrim will be Aluminum Anodized. Limitations exist based on rear wheel spacing, camber and armrest selections.</i>	

Left Side Drive [RGDRV002]	N/C
Right Side Drive [RGDRV003]	N/C
Aluminum Drive Handrim [RGDRV004]	N/C
Poly Coated Drive Handrim [RGDRV005]	£230

Wheels

No selection needed with One Arm Drive

		Rear Wheel Size					
		Price	20"	22"	24"	25"	26"
Rear Wheel Type	MAXX Spoke	N/C	[RGWH1120]	[RGWH1122]	[RGWH1124]	[RGWH1125]	[RGWH1126]
	Mag	£288	[RGWH1220]	[RGWH1222]	[RGWH1224]		
	Mag (Color) ²	£288	[RGWH1320]	[RGWH1322]	[RGWH1324]		
	Superlight Spoke	£230			[RGWHE224]	[RGWHE225] ¹	
	Spinergy Spox	£771		[RGWHE422]	[RGWHE424]	[RGWHE425]	[RGWHE426]
	Spinergy LX	£1001		[RGWHE522]	[RGWHE524]	[RGWHE525]	[RGWHE526]

¹ Available only with Superlight and Plastic Coated Handrims

² Color selection on page 6

Spinergy Options

Select spoke color and handrim tab for Spinergy wheels

Black Spokes [RGSSC001] Red Spokes [RGSSC002]
Blue Spokes [RGSSC005] White Spokes [RGSSC003] _____ N/C

Rear Tire Type

Full Profile

Street Pneumatic 1 3/8" [RGTIRO01] _____ N/C
Not available on 25" wheel size
Full Poly 1 3/8" [RGTIRO02] _____ £110
Not available on Superlight wheel or 25" or 26" wheel sizes
Pneumatic with Airless Insert 1 3/8" [RGTIRO04] _____ £110
Not available on Superlight wheel or 25" wheel size
Schwalbe HS-127 1 3/8" [RGTIRO17] _____ £110
Available only on 24" wheel size. Not available on Mag wheel
Mountain 2" Black [RGTIRO09] _____ £110
Available only on 22", 24" and 25" MAXX wheels and Spinergy wheels

Low Profile

Low Poly 1" [RGTIRO03] _____ N/C
Not available on Superlight wheel, 25" or 26" wheel sizes
Sentinel Solid Tire 1" [RGTIRO05] _____ £110
Not available on Superlight wheel
High Pressure 1" [RGTIRO06] _____ N/C
Not available on Mag wheel
Schwalbe One 1" [RGTIRO18] _____ £110
Available only on 24" and 25" wheel sizes. Not available on Mag wheel
Schwalbe Marathon Plus Evolution 1" [RGTIRO10] _____ £110
Not available on Mag wheel, 20" or 22" wheel sizes

Handrims

No selection needed with One Arm Drive

Standard Handrims

Aluminum Anodized [RGHRM001] _____ N/C
Superlight Aluminum Anodized [RGHRM010] _____ N/C
Available only on 24" and 25" wheel sizes. Available only on MAXX Spokes, Superlight or Spinergy wheels.
Plastic Coated [RGHRM002] _____ £127
Projection (8 Vertical Projections) [RGHRM003] _____ £345
Available only with full profile tires in combination with 22", 24", and 26" MAXX Spoke wheels, or 22" and 24" Mag wheels
Omit Handrim [RGHRM012] _____ N/C

Ergonomic Handrims

Not available with 20" or 26" wheels. Tab mount only available with 24" wheels.

Ergo Lite Regular Grip [RGHRM005] _____ £276
Ergo Standard Regular Grip [RGHRM008] _____ £276
Ergo Lite Extra Grip [RGHRM006] _____ £276
Ergo Standard Extra Grip [RGHRM009] _____ £276
Ergo Lite Omit Grip [RGHRM004] _____ £276
Ergo Standard Omit Grip [RGHRM007] _____ £276

Handrim Options

Handrim Mount

If no selection is made, Standard Mount will be used

Close Mount [RGTAB001] _____ N/C

Micro-Grip

3M Micro Grip [RGACC006] _____ £104
Not available on Plastic Coated or Projection Handrims

Axles

Quick Release [RGAXO003] _____ N/C
Quad Quick Release [RGAXO001] _____ £196
Not available with Mag or Spinergy wheels

Center of Gravity Presets

Measured from front of back post to center of rear axle. 0" indicates axle will be directly under back post. See reference J on page 7.

Amputee Positions				
-1.5" [RGCGA150]	-1" [RGCGA100]	-5" [RGCGA050]	0" [RGCGS000]	.25" [RGCGS025]
Standard Positions				
.5" [RGCGS050]	.75" [RGCGS075]	1" [RGCGS100]	1.25" [RGCGS125]	1.5" [RGCGS150]
1.75" [RGCGS175]	2" [RGCGS200]	2.25" [RGCGS225]	2.5" [RGCGS250]	2.75" [RGCGS275]
3" [RGCGS300]	3.25" [RGCGS325]	3.5" [RGCGS350]	3.75" [RGCGS375]	4" [RGCGS400]
4.25" [RGCGS425]	4.5" [RGCGS450]	4.75" [RGCGS475]	5" [RGCGS500]	

Rear Wheel Spacing

1" minimum with 0° camber. 1.5" minimum with T-Arms. The following apply if One Arm Drive is selected: 1.75" minimum with T-Arms. 1.25" minimum with select Rear Seat Heights. See Reference K on page 7.

0.75" [RGWS075]	1" [RGWS100]	1.25" [RGWS125]	
1.5" [RGWS150]	1.75" [RGWS175]		N/C

Rear Wheel Camber

The following apply if One Arm Drive is selected: 6° and 8° not available. 4° not available with T-Arms.

0° [RGWC00]	2° [RGWC02]	4° [RGWC04]	6° [RGWC06]	
8° [RGWC08]				N/C

Wheel Locks

Wheel Lock

Push to Lock [RGWLK002]	N/C
Push to Lock Flush Mount [RGWLK001]	£69
Pull to Lock [RGWLK003]	N/C
Under Seat Scissor Lock [RGWLK007]	£69
Low Profile Scissor Lock [RGWLK012]	£58
<i>Not available with Low Poly, Full Poly or Sentinel Solid tires</i>	
Short Thro Scissor Lock [RGWLK006]	£230
Omit [RGWLK008]	N/C

Wheel Lock Options

Extension Handles [RGWLK004]	£69
<i>Available only with Push to Lock or Pull to Lock Wheel Locks</i>	
Grade Aids [RGWLK009]	£173
<i>Available only with Push to Lock Wheel Locks and Pneumatic or Pneumatic with Airless Insert tires</i>	

Armrests

Armrest Type	Length	Price
Height Adjustable T-Arm Tall (11.5"-15")	Desk [RGARM039]	£230
	Full [RGARM040]	
Height Adjustable T-Arm (8.5"-13.5")	Desk [RGARM003]	£230
	Full [RGARM004]	
Height adjustable T-Arm Low (6.5"-10.5") ¹	Desk [RGARM022]	£230
	Full [RGARM023]	

¹ If requested configuration is unavailable, standard height T-Arm will be used. 14" minimum seat depth with Offset Frame.

² Includes direct mount receiver

Armrest Type	Price
Swing Away (8.5"-12.5") [RGARM002]	£230
Swing Away with Side Guard (8.5"-12.5") ² [RGARM018]	£288
Omit Armrests [RGARM007]	N/C

Side Guards

Direct Mount Side Guards

Not available with T-Arms or Swing Away with Side Guard Armrests

Composite Side Guards - Direct Mount Receiver [RGSDG002]	£58
Aluminum Straight - Direct Mount Receiver - Black Matte [RGSDG024]	£230
Aluminum Straight - Direct Mount Receiver - Match Frame Color [RGSDG025]	£288
<i>Not available with Hydrographics</i>	
Aluminum Adjustable Fender w/Foam Pad - Direct Mount Receiver [RGSDG014]	£230
<i>Not available with Adjustable Depth greater than 1/2"</i>	
Aluminum Adjustable Fender w/Foam Pad - Direct Mount Receiver [RGSDG015]	£288
<i>Not available with Adjustable Depth greater than 1/2". Not available with Hydrographics.</i>	
Carbon Fiber Adjustable Fender - Direct Mount Receiver [RGSDG013]	£288
<i>Not available with Adjustable Depth greater than 1/2"</i>	

Clamp-On Mount Side Guards

Not available with T-Arms or Swing Away with Side Guard Armrests

Composite Side Guards - Clamp-On Receiver [RGSDG004]	£58
Aluminum Straight - Clamp On Receiver - Black Matte [RGSDG020]	£230
Aluminum Straight - Clamp-On Receiver - Match Frame Color [RGSDG018]	£288
<i>Not available with Hydrographics</i>	
Aluminum Adjustable Fender w/Foam Pad - Clamp-On Receiver - Black [RGSDG017]	£230
Aluminum Adjustable Fender w/Foam Pad - Clamp-On Receiver [RGSDG018]	£288
<i>Not available with Hydrographics</i>	
Carbon Fiber Adjustable Fender - Clamp-On Receiver [RGSDG016]	£288

Accessories

Anti-Tips

Rear Anti-Tips [RGATP001]	£115
User Activated - Pair [RGATP012]	£345
User Activated - Left [RGATP010]	£173
User Activated - Right [RGATP011]	£173
Omit [RGATP000]	N/C

Calf Strap

Velcro Adjustable [RGACC002]	£41
Padded Velcro Adjustable [RGACC003]	£75
Bodypoint Padded AeroMesh [RGACC010]	£75

Pelvic Positioning Belt

1 1/2" Auto Buckle Pelvic Positioning Belt [RGPOB001]	£58
1 1/2" Padded Auto Buckle Pelvic Positioning Belt [RGPOB002]	£173
Small Bodypoint Evoflex (16-24") [RGPOB007]	£115
Medium Bodypoint Evoflex (19-29") [RGPOB008]	£115
Large Bodypoint Evoflex (25-37") [RGPOB009]	£115

Other Accessories

Tool Kit [RGACC005]	£20
Removable Underseat Pouch [RGACC015]	£30
Ki Mobility Backpack [RGACC009]	£100
Luggage Carrier [RGACC025]	£92
Black Neoprene Impact Guards [RGIMG001]	£58
Under Chair Carrier and Removable Bag [RGACC001]	£46
<i>Not available with Seat Widths less than 14" or scissor locks</i>	
Cane and Crutch Holder [RGACC007]	£115

Belt Mounting Kit

Ki Mobility Belt Mounting Clamp [RGACC023]	N/C
Bodypoint Belt Mounting Bracket Kit [RGACC013]	£115

Color

Options	Price	Colors																				
		Gloss												Matte								
		Black	Blue Skies	Bubblegum Pink	Candy Blue	Candy Purple	Candy Red	Copper	Granny Smith	MRGallic Orange	Silver	Teal	Yellow	Charcoal Gray	Emerald Green	Iridescent Burgundy	Iridescent Navy	Military Green	Retro Red	Royal Blue	Shadow Black	Smooth White
Frame	N/C																					
Backrest Release Bar	£58																					
	N/C	Black Anodized																				
MAXX Mag Color																						

Colored Anodizing Package

Select desired color. If no selection is made, black will be used.

	Color	Black ¹	Blue	Bright Green	Gold	Orange	Pewter	Pink	Purple	Red
Retail		N/C	N/C	N/C	N/C	N/C	N/C	N/C	N/C	N/C
Option	Back Plates ²	[RGANOB01]	[RGANOB02]	[RGANOB03]	[RGANOB09]	[RGANOB04]	[RGANOB05]	[RGANOB06]	[RGANOB07]	[RGANOB08]
	Camber Tube ³	[RGANOC01]	[RGANOC02]	[RGANOC03]	[RGANOC09]	[RGANOC04]	[RGANOC05]	[RGANOC06]	[RGANOC07]	[RGANOC08]
	Towers	[RGANOT01]	[RGANOT02]	[RGANOT03]	[RGANOT09]	[RGANOT04]	[RGANOT05]	[RGANOT06]	[RGANOT07]	[RGANOT08]
	Side Guard Receiver	[RGANOS01]	[RGANOS02]	[RGANOS03]	[RGANOS09]	[RGANOS04]	[RGANOS05]	[RGANOS07]	[RGANOS08]	[RGANOS09]
	Caster Fork ³	[RGANOF01]	[RGANOF02]	[RGANOF03]	[RGANOF09]	[RGANOF04]	[RGANOF05]	[RGANOF06]	[RGANOF07]	[RGANOF08]

¹ Black is no charge

² Not available with Transit Option

³ Available only with Standard or Single Sided

References

A. Seat Width See Diagram 1

Measured from outside of frame tube on one side to outside of frame tube on other side.

B. Seat Depth See Diagram 3

Measured from front of back posts to front edge of seat sling. Seat sling starts at beginning of bend at front of frame.

C. Sling Position See Diagram 2

In Performance Position (0°), sling is at front of frame to keep chair shorter and more maneuverable. The 1° or 2° seat sling position extends the frame in front of the sling by that amount. Adding more frame can improve stability and provide support to aid in transfers.

D. Front Frame Bend See Diagram 3

Angle between front seat tube and ground.

E. Footrest Taper See Diagram 1

Indicates bend of front frame creating footrest. Taper is measured from outside of seat frame to outside of front tube.

F. Front Seat Height See Diagram 1

Measured from floor to top of seat tube at front of seat sling.

G. Footrest Width See Diagram 1

Measured from inside of front tube to inside of front tube on other side. Footrest width is listed below for each seat width and taper.

Seat Width	Straight	1° Y Taper	2° Y Taper	3.5° V Taper	7° V Taper
Outside Measurement	2" Narrower than Seat Width	4" Narrower than Seat Width	6" Narrower than Seat Width	4" Narrower than Seat Width	6" Narrower than Seat Width
12"	9.5"	7.5"		7.5"	
13"	10.5"	8.5"		8.5"	
14"	11.5"	9.5"		9.5"	
15"	12.5"	10.5"		10.5"	
16"	13.5"	11.5"	9.5"	11.5"	9.5"
17"	14.5"	12.5"	10.5"	12.5"	10.5"
18"	15.5"	13.5"	11.5"	13.5"	11.5"
19"	16.5"	14.5"	12.5"	14.5"	12.5"
20"	17.5"	15.5"	13.5"	15.5"	13.5"

H. Seat to Footrest Length See Diagram 3

Measured from front edge of seat sling to top rear of footrest. Footrest length of at least 2.5" shorter than front seat height recommended.

I. Rear Seat Height See Diagram 3

Measured from floor to top of seat tube at front of back post. Rear seat height is custom to needs.

J. Center of Gravity Preset See Diagram 3

Measured from front of back post to center of rear axle. 0° indicates axle will be directly under back post.

K. Rear Wheel Spacing See Diagram 1

Measured from outside of seat back to inside of rear tire. Adjustable out .5" from setting.

L. Caster Position See Diagram 4

The Optimized Position keeps the casters tighter to the frame for enhanced maneuverability. Expanded Position pushes the casters out .75" past the Optimized Position for increased stability.

Diagram 1




Diagram 2




Diagram 3





Diagram 4



Axiom G

		Width															
		14"		15"		16"		17"		18"		19"		20"		22"	
		Cushion	Outer Cover	Cushion	Outer Cover	Cushion	Outer Cover	Cushion	Outer Cover	Cushion	Outer Cover	Cushion	Outer Cover	Cushion	Outer Cover	Cushion	Outer Cover
Retail		£85	£70	£85	£70	£85	£70	£85	£70	£85	£70	£85	£70	£85	£70	£85	£70
Height	14"																
	15"																
	16"																
	17"																
	18"																
	19"																
	20"																
	21"																
	22"																

Axiom S

		Width															
		14"		15"		16"		17"		18"		19"		20"		22"	
		Cushion	Outer Cover	Cushion	Outer Cover	Cushion	Outer Cover	Cushion	Outer Cover	Cushion	Outer Cover	Cushion	Outer Cover	Cushion	Outer Cover	Cushion	Outer Cover
Retail		£120	£70	£120	£70	£120	£70	£120	£70	£120	£70	£120	£70	£120	£70	£120	£70
Height	14"																
	15"																
	16"																
	17"																
	18"																
	19"																
	20"																
	21"																
	22"																

Axiom P

		Width															
		14"		15"		16"		17"		18"		19"		20"		22"	
		Cushion	Outer Cover	Cushion	Outer Cover	Cushion	Outer Cover	Cushion	Outer Cover	Cushion	Outer Cover	Cushion	Outer Cover	Cushion	Outer Cover	Cushion	Outer Cover
Retail		£170	£70	£170	£70	£170	£70	£170	£70	£170	£70	£170	£70	£170	£70	£170	£70
Height	14"																
	15"																
	16"																
	17"																
	18"																
	19"																
	20"																
	22"																

Axiom SP - Visco

		Width																										
		14"			15"			16"			17"			18"			19"			20"			21"			22"		
		Cushion	Cover		Cushion	Cover		Cushion	Cover		Cushion	Cover		Cushion	Cover		Cushion	Cover		Cushion	Cover		Cushion	Cover				
Retail		£350	£70	£350	£70	£350	£70	£350	£70	£350	£70	£350	£70	£350	£70	£350	£70	£350	£70	£350	£70	£350	£70					
Height	14"																											
	15"																											
	16"																											
	17"																											
	18"																											
	19"																											
	20"																											
	22"																											

Axiom SP - Fluid

		Width																										
		14"			15"			16"			17"			18"			19"			20"			21"			22"		
		Cushion	Cover		Cushion	Cover		Cushion	Cover		Cushion	Cover		Cushion	Cover		Cushion	Cover		Cushion	Cover		Cushion	Cover				
Retail		£320	£70	£320	£70	£320	£70	£320	£70	£320	£70	£320	£70	£320	£70	£320	£70	£320	£70	£320	£70	£320	£70					
Height	14"																											
	15"																											
	16"																											
	17"																											
	18"																											
	19"																											
	20"																											
	22"																											


Solid Seat Insert

Weight capacity: 350 lbs

		Width								
		14"	15"	16"	17"	18"	19"	20"	21"	22"
		Cushion	Cushion	Cushion	Cushion	Cushion	Cushion	Cushion	Cushion	Cushion
Retail		£55	£55	£55	£55	£55	£55	£55	£55	£55
Height	14"									
	15"									
	16"									
	17"									
	18"									
	19"									
	20"									
	21"									
	22"									

Axiom Back Width	Fits Back Frame Widths	
	Inner Slots	Outer Slots
12" (30.5cm)	11" (28.0cm)	12" -14" (30.5cm -35.5cm)
14" (35.5cm)	12" -13" (30.5cm -33.0cm)	14" -16" (35.5cm -40.5cm)
16" (40.5cm)	14" -15" (35.5cm -38.0cm)	16" -18" (40.5cm -45.5cm)
18" (45.5cm)	16" - 17" (40.5cm - 43.0cm)	18" - 20" (45.5cm - 50.5cm)
20" (51.0cm)	18" -19" (46.0cm -48.5cm)	20" -22" (50.5cm -55.5cm)
22" (56.0cm)	20" -21" (51.0cm -53.5cm)	22" -24" (56.0cm -61.0cm)

Hardware Description	Code
Improved DEX Hardware, removable adjustable hardware (300 lb max user weight)	DEX
New BTX Hardware, removable adjustable hardware not recommended for, or installed on Ultralight Rigid Frame Wheelchairs (400 lb max user weight)	BTX
New Fixed Hardware, non-removeable hardware up to (300 lb max user weight)	FXT



Axiom Back

Hardware is included with back. Low Back selections are not available with headrests, modifications, or posture supports. Weight capacity: 300 lbs.

Axiom Posterior Back (PB)

		Width						
		12"	14"	16"	18"	20"	22"	
Retail		£550	£550	£550	£550	£550	£550	
Height	Low	8"	[AXPB1208]	[AXPB1408]	[AXPB1608]	[AXPB1808]	[AXPB2008]	
		10"		[AXPB1410]	[AXPB1610]	[AXPB1810]	[AXPB2010]	[AXPB2210]
		12"	[AXPB1212]	[AXPB1412]	[AXPB1612]	[AXPB1812]	[AXPB2012]	
	Mid	14"		[AXPB1414]	[AXPB1614]	[AXPB1814]	[AXPB2014]	[AXPB2214]
		16"	[AXPB1216]	[AXPB1416]	[AXPB1616]	[AXPB1816]	[AXPB2016]	
		18"		[AXPB1418]	[AXPB1618]	[AXPB1818]	[AXPB2018]	[AXPB2218]

Back Hardware

DEX Hardware [DEX] BTX Hardware [BTX] FXT Hardware [FXT] _____ N/C

Lateral Support for AXPB Backs

Lateral Trunk Supports for AXPB Backs							
Pair	Left	Right	Description	Item Code	Pair	Each	
			3 x 4 Lateral Pads with Fixed Mounting Hardware for Axiom Posterior Back	PBLF3x4P	£295	£148	
			4 x 5 Lateral Pads with Fixed Mounting Hardware for Axiom Posterior Back	PBLF4x5P	£295	£148	
			3 x 4 Lateral Pads with Swing Away Mounting Hardware for Axiom Posterior Back	PBLS3x4P	£628	£314	
			4 x 5 Lateral Pads with Swing Away Mounting Hardware for Axiom Posterior Back	PBLS4x5P	£628	£314	
			3 x 4 Lateral Pads with Adjustable Swing Away Mounting Hardware	PBLAS3x4P	£628	£314	
			4 x 5 Lateral Pads with Adjustable Swing Away Mounting Hardware	PBLAS4x5P	£628	£314	

Axiom Posterior Lateral Back (PL)

		Width						
		12"	14"	16"	18"	20"	22"	
Retail		£550	£550	£550	£550	£550	£550	
Height	Low	8"	[AXPL1208]	[AXPL1408]	[AXPL1608]	[AXPL1808]	[AXPL2008]	
		10"		[AXPL1410]	[AXPL1610]	[AXPL1810]	[AXPL2010]	[AXPL2210]
		12"	[AXPL1212]	[AXPL1412]	[AXPL1612]	[AXPL1812]	[AXPL2012]	
	Mid	14"		[AXPL1414]	[AXPL1614]	[AXPL1814]	[AXPL2014]	[AXPL2214]
		16"	[AXPL1216]	[AXPL1416]	[AXPL1616]	[AXPL1816]	[AXPL2016]	
		18"		[AXPL1418]	[AXPL1618]	[AXPL1818]	[AXPL2018]	[AXPL2218]
Tall	20"	[AXPL1220]	[AXPL1420]	[AXPL1620]	[AXPL1820]	[AXPL2020]		
	22"		[AXPL1422]	[AXPL1622]	[AXPL1822]	[AXPL2022]	[AXPL2222]	

Back Hardware

DEX Hardware [DEX] BTX Hardware [BTX] FXT Hardware [FXT] _____ N/C

Lateral Support for AXPL

Lateral Trunk Supports for AXPL & AXLB Backs						
Pair	Left	Right	Description	Item Code	Pair	Each
			3 x 4 Lateral Pads with Fixed Mounting Hardware for Axiom Posterior Lateral Back	PLLF3x4P	£295	£295
			4 x 5 Lateral Pads with Fixed Mounting Hardware for Axiom Posterior Lateral Back	PLLF4x5P	£295	£295
			3 x 4 Lateral Pads with Swing Away Mounting Hardware for Axiom Posterior Lateral Back	PLLS3x4P	£628	£314
			4 x 5 Lateral Pads with Swing Away Mounting Hardware for Axiom Posterior Lateral Back	PLLS4x5P	£628	£314
			3 x 4 Lateral Pads with Adjustable Swing Away Mounting Hardware	PLLAS3x4P	£628	£314
			4 x 5 Lateral Pads with Adjustable Swing Away Mounting Hardware	PLLAS4x5P	£628	£314

Axiom Deep Lateral Back (DL)

		Width						
		12"	14"	16"	18"	20"	22"	
Retail		£550	£550	£550	£550	£550	£550	
HCPC		E2620	E2620	E2620	E2620	E2620	E2621	
Height	Mid	14"	[AXDL1214]	[AXDL1414]	[AXDL1614]	[AXDL1814]	[AXDL2014]	[AXDL2214]
		16"		[AXDL1416]	[AXDL1616]	[AXDL1816]	[AXDL2016]	
		18"	[AXDL1218]	[AXDL1418]	[AXDL1618]	[AXDL1818]	[AXDL2018]	[AXDL2218]
	Tall	20"	[AXDL1220]	[AXDL1420]	[AXDL1620]	[AXDL1820]	[AXDL2020]	
		22"		[AXDL1422]	[AXDL1622]	[AXDL1822]	[AXDL2022]	[AXDL2222]

Back Hardware

DEX Hardware [DEX] BTX Hardware [BTX] FXT Hardware [FXT] _____ N/C

Headrest

Description	Std. Straight Removable Adjustable Hardware		Multi-Axis Removable Adjustable Hardware		Posalinc Removable Adjustable Hardware		HCPCS
	Item Code	Retail	Item Code	Retail	Item Code	Retail	
Comfort Soft Headrest 8"	[CS08SS]	£268	[CS08MA]	£268	[CS08PL]	£268	E0955 / E1028
Comfort Soft Headrest 10"	[CS10SS]	£268	[CS10MA]	£268	[CS10PL]	£268	E0955 / E1028
Comfort Soft Headrest 14"	[CS14SS]	£268	[CS14MA]	£268	[CS14PL]	£268	E0955 / E1028
Contour Positioning Pad, 8" X 6.25"	[CP08SS]	£303	[CP08MA]	£303	[CP08PL]	£303	E0955 / E1028
Contour Positioning Pad, 9" X 6.5"	[CP09SS]	£303	[CP09MA]	£303	[CP09PL]	£303	E0955 / E1028
Contour Positioning Pad, 10" X 7.5"	[CP10SS]	£303	[CP10MA]	£303	[CP10PL]	£303	E0955 / E1028
Tri-soft Headrest 10"	[TS10SS]	£386	[TS10MA]	£386	[TS10PL]	£386	E0955 / E1028
Tri-soft Headrest 14"	[TS14SS]	£386	[TS14MA]	£386	[TS14PL]	£386	E0955 / E1028

Description	Item Code	Retail	HCPCS
Universal Headrest Mounting Plate	[UMP]	£43	K0108

Not required for Axiom Headrests

Anterior Thoracic Supports

Selection of a modification from the Axiom Back Modifications section required in order to mount posture supports. Not available with Low Backs.

Anterior Thoracic Supports			
Description	Item Code	HCPCS	Retail
Dynaform Chest Support: 2-Piece with Hook&Loop, 3" High Pad, Small	[14941]	E0960	£100
Dynaform Chest Support: 2-Piece with Hook&Loop, 4" High Pad, Medium	[14942]	E0960	£100
Dynaform Chest Support: 2-Piece with Hook&Loop, 5" High Pad, Large	[14943]	E0960	£100
Dynaform Chest Support: 1-Piece with Buckle, 3" Height x 20" Width, Small	[14924]	E0960	£100
Dynaform Chest Support: 1-Piece with Buckle, 4" Height x 20" Width, Medium	[14927]	E0960	£100
Dynaform Chest Support: 1-Piece with Buckle, 4" Height x 22" Width, Large	[14928]	E0960	£100
Dynaform Anterior Trunk Support: Standard, Dynamic, X-Small	[14851]	E0960	£100
Dynaform Anterior Trunk Support: Standard, Dynamic, Small	[14852]	E0960	£100
Dynaform Anterior Trunk Support: Standard, Dynamic, Medium	[14853]	E0960	£100
Dynaform Anterior Trunk Support: Standard, Dynamic, Large	[14854]	E0960	£100
Dynaform Anterior Trunk Support: Slim Cut, Dynamic, X-Small	[14856]	E0960	£100
Dynaform Anterior Trunk Support: Slim Cut, Dynamic, Small	[14857]	E0960	£100
Dynaform Anterior Trunk Support: Slim Cut, Dynamic, Medium	[14858]	E0960	£100
Dynaform Anterior Trunk Support: Slim Cut, Dynamic, Large	[14859]	E0960	£100
AirLogic Anterior Trunk Support: Standard, Structure, Small	[14282]	E0960	£300
AirLogic Anterior Trunk Support: Standard, Structure, Medium	[14283]	E0960	£300
AirLogic Anterior Trunk Support: Standard, Structure, Large	[14284]	E0960	£300
AirLogic Anterior Trunk Support: Standard, Structure, X-Large	[14285]	E0960	£300
AirLogic Anterior Trunk Support: Slim Cut, Dynamic, Small	[14257]	E0960	£300
AirLogic Anterior Trunk Support: Slim Cut, Dynamic, Medium	[14258]	E0960	£300
AirLogic Anterior Trunk Support: Slim Cut, Dynamic, Large	[14259]	E0960	£300
AirLogic Anterior Trunk Support: Slim Cut, Dynamic, X-Large	[14260]	E0960	£300

Accessories

Adjustable Shoulder Guides (K0108) [AXBSG] <i>Includes Anterior Thoracic Mounting Modification. Not available with Low Backs.</i>	£43
Extra Lumbar Pad (K0108) [LUMPAD] <i>Available only with Axiom Lateral Back and Axiom Deep Lateral back</i>	N/C
Privacy Flap (K0108) [PRFLAP]	£44

Axiom Back Modifications

Axiom Back Modifications

Extended lead times will apply for all modifications. Makes back non-refundable.

Anterior Thoracic Mounting Modification [AXSMM] <i>Provides holes for Anterior Thoracic Supports and Shoulder Guides</i>	N/C
---	-----

Gloss Finish

If no selection is made, standard Shadow Black will be used with no additional lead time

Black Back Shell [AXMBLK]	N/C
Blue Skies Back Shell [AXMSKY]	N/C
Bubblegum Pink Back Shell [AXMPNK]	N/C
Candy Blue Back Shell [AXMCBL]	N/C
Candy Purple Back Shell [AXMPUR]	N/C
Candy Red Back Shell [AXMRED]	N/C
Copper Back Shell [AXMCOP]	N/C
Granny Smith Back Shell [AXMGSM]	N/C
Metallic Orange Back Shell [AXMMOR]	N/C
Yellow Back Shell [AXMYLW]	N/C
Silver Back Shell [AXMSIL]	N/C
Teal Back Shell [AXMTEA]	N/C

Matte Finish

If no selection is made, standard Shadow Black will be used with no additional lead time

Military Green Back Shell [AXMAGM]	N/C
Charcoal Gray Back Shell [AXMCGM]	N/C
Emerald Green Back Shell [AXMEGM]	N/C
Iridescent Burgundy Back Shell [AXMIBM]	N/C
Iridescent Navy Back Shell [AXMINM]	N/C
Retro Red Back Shell [AXMRRM]	N/C
Royal Blue Back Shell [AXMRBM]	N/C
Smooth White Back Shell [AXMSWM]	N/C