



Etac UK 0(121)561-2222 | orders.uk@etac.com |
www.kimobility.com | www.etac.com



Company:

Ship To:

Account Number:

Date:

Address:

Quote: PO Number:

City:

State:

ATP:

Zip:

Phone:

Marked For:

Email/Fax:

Associated Quote/Order#:

Submit Form

Effective: 1/10/2024

Rogue ALX £1999
250 lb Weight Capacity

Frame

Seat Width

See reference A on page 6

14" [R2SWT014]	N/C
15" [R2SWT015]	N/C
16" [R2SWT016]	N/C
17" [R2SWT017]	N/C
18" [R2SWT018]	N/C
19" [R2SWT019]	N/C
20" [R2SWT020]	N/C
21" [R2SWT021]	N/C
22" [R2SWT022]	N/C

Seat Depth

See reference B on page 6

14" [R2SDE014]	N/C
15" [R2SDE015]	N/C
16" [R2SDE016]	N/C
17" [R2SDE017]	N/C
18" [R2SDE018]	N/C
19" [R2SDE019]	N/C
20" [R2SDE020]	N/C

Weight Capacity

250 lbs [R2STD001]	N/C
350 lb (Heavy Duty) [R2HVD001]	£303

Requires 6" casters

Transit Option

Transit System [R2TRNOPT]	N/C
125 lb maximum weight capacity	
Omit [R2TRNOP2]	N/C

Sling Position

In Performance Position (0"), sling is at front of frame. The 1" or 2" sling position extends the frame in front of the sling by that amount. See reference C on page 6.

0" (Performance) [R2SP01]	N/C
1" [R2SP03]	N/C
2" (Classic) [R2SP02]	N/C

Adjustable Depth

Adds the capability to grow seat in depth

No Rear Adjustment [R2AD05]	N/C
1/2" [R2AD06]	N/C
1" [R2AD01]	N/C
2" [R2AD02]	£164
3" [R2AD03]	£164

Long Sling

Available only with 1", 2" and 3" Adjustable Depth

Long Sling [R2AD04]	N/C
---------------------	-----

Sling will be maximum 2" longer than seat depth

Seat Option

Standard Seat Sling [R2SSS001]	N/C
Aluminum Solid Seat Pan [R2ACCO04]	£260

Seat Pan is 1" longer when Long Sling option is selected

Front Frame Bend

Angle between front seat tube and ground. See reference D on page 6.

90° [R2FA01]	N/C
85° [R2FA03]	N/C
80° [R2FA02]	N/C
75° [R2FA04]	£140
70° [R2FA05]	£140

Backrest

Back Post / Height Range / Back Height							
Back Post	Price	Height Range	Back Height				
Height Adjustable	N/C	Height Adjustable Short Back Post [RGBOP025]	11" ²	12"	13"	14"	
		Height Adjustable Medium Back Post [RGBOP026]	13"	14"	15"	16"	17"
		Height Adjustable Tall Back Post [RGBOP027]	17"	18"	19"	20"	
Fixed Height ¹	N/C	Fixed Height Straight Back Post [RGBOP011]	9" ²	10" ²	11" ²	12" ³	13"
			14"	15"	16"	17"	18"

¹ Not available with back upholstery

² Not available with Fold Down Push Handle

³ Fold Down Push Handle requires Height Adjustable Rigidizer Bar

Push Handles

Omit Push Handle [R2BCK008]	N/C
Standard Push Handle [R2BCK004]	N/C
Low Profile Push Handle [R2BCK007]	£134
<i>2" shorter than standard push handles</i>	
Bolt-On Push Handle [R2BOP005]	£182
<i>Handles will be 1" - 2" above back height</i>	
Fold Down Push Handle [R2BCK006]	£158
<i>Not available with 350 lb Weight Capacity</i>	

Rigidizer Bar

Adjustable Height Rigidizer Bar [R2BCK001]	N/C
Non-Adjustable Height Rigidizer Bar [R2BCK010]	N/C
<i>Not available with Fixed Height Back Post</i>	

Back Upholstery Options

Available only with Height Adjustable Back Posts. If selecting both upholstery and Axiom Back, the upholstery will be installed.

Padded Velcro Adjustable with Standard Cover [R2BUP001]	N/C
<i>Not available with 350 lb Weight Capacity</i>	
Padded Velcro Adjustable with Air Mesh Breathable Cover [R2BUP006]	£73
<i>Not available with 350 lb Weight Capacity</i>	
Tension Adjustable with Standard Cover [R2BUP002]	N/C
Tension Adjustable with Air Mesh Cover [R2BUP007]	£182
Omit [R2BUP003]	N/C

Locking Backrest

Locking Fold Down Back [R2BCK011]	£43
<i>Not available with 350 lb Weight Capacity</i>	

Seat Back Angle

Standard is 90° (perpendicular to floor). Minus (-°) closes back angle. Plus (+°) opens back angle.

Standard [R2SBA001]	+3° [R2SBA002]	+6° [R2SBA003]
+9° [R2SBA004]	-3° [R2SBA005]	-6° [R2SBA006]
-9° [R2SBA007]	-12° [R2SBA008]	N/C

Front Seat Height and Casters

Caster Size / Footrest Taper / Front Seat Height

Select Front Seat Height and Footrest Taper in row corresponding to Caster Size. 6" casters not available with 85° and 90° Front Frame Bend. Footrest Taper indicates bend of front frame creating footrest width. Taper is measured from outside of seat frame to outside of front tube. 2" Taper not available with 75° Front Frame Bend. Taper is not available with 70° Front Frame Bend. See reference E on page 6. Front Seat Height is measured from floor to top of seat tube at front of seat sling. Adjustable up or down by 1/2". See reference F on page 6.

Footrest Taper	Caster Size	Front Seat Height														
		13.5"	14"	14.5"	15"	15.5"	16"	16.5"	17"	17.5"	18"	18.5"	19"	19.5"	20"	20.5"
No Taper (Straight) [TXFT00]	3"	13.5"	14"	14.5"	15"	15.5"	16"	16.5"	17"	17.5"	18"	18.5"	19"	19.5"	20"	20.5"
	4"	14"	14.5"	15"	15.5"	16"	16.5"	17"	17.5"	18"	18.5"	19"	19.5"	20"	20.5"	
	5"	14.5"	15"	15.5"	16"	16.5"	17"	17.5"	18"	18.5"	19"	19.5"	20"	20.5"		
	6"	15"	15.5"	16"	16.5"	17"	17.5"	18"	18.5"	19"	19.5"	20"	20.5"			
1" [TXFT01]	3"	16.5"	17"	17.5"	18"	18.5"	19"	19.5"	20"	20.5"						
	4"	17"	17.5"	18"	18.5"	19"	19.5"	20"	20.5"							
	5"	17.5"	18"	18.5"	19"	19.5"	20"	20.5"								
	6"	18"	18.5"	19"	19.5"	20"	20.5"									
2" [TXFT02]	3"	16.5"	17"	17.5"	18"	18.5"	19"	19.5"	20"	20.5"						
	4"	17"	17.5"	18"	18.5"	19"	19.5"	20"	20.5"							
	5"	17.5"	18"	18.5"	19"	19.5"	20"	20.5"								
	6"	18"	18.5"	19"	19.5"	20"	20.5"									

Casters

		Caster Size				
		Price	3"	4"	5"	6"
Caster Type	.75" Rollerblade	N/C	[R2CAS002]	[R2CAS004]		
	.75" Lighted Rollerblade	£81	[R2CAS001]	[R2CAS005]		
	1" Poly	N/C		[R2CAS006]	[R2CAS010]	[R2CAS014]
	1" Poly Aluminum	£96		[R2CAS007]	[R2CAS011]	
	1" Pneumatic	£103				[R2CAS017]
	1.5" Poly	N/C		[R2CAS032]	[R2CAS025]	[R2CAS026]
	1.5" Soft Roll Aluminum	£96	[R2CAS003]	[R2CAS009]	[R2CAS013]	[R2CAS016]

Fork Options

Standard Fork [R2CST001]	N/C
Slotted Fork [R2FRK001]	N/C
Single Sided Fork [R2FRK002]	£206
<i>Not available with 6" Caster Size</i>	
Frog Legs II Suspension System [R2CST014]	£424
<i>Lightweight shock absorbing fork. Not available with 350 lb Weight Capacity.</i>	
Caster Pin Locks [R2CST005]	£194
<i>Not Available with Single Sided Fork</i>	

Footrest

Footplate

Tubular "Open" Footplate [R2FPL004]	N/C
<i>Not available w/footrest widths greater than 17.5". See Reference G on page 6.</i>	
Tubular Footplate with ABS Cover [R2FPL005]	£86
<i>Not available w/footrest widths greater than 15.5". See Reference G on page 6.</i>	
Aluminum Angle Adjustable [R2FPL006]	£182
Aluminum Angle Adjustable with Riser [R2FPL007]	£182
<i>Not available with 350 lb Weight Capacity</i>	
High Mount Aluminum Angle Adjustable [R2FPL008]	£182
<i>Not available with 350 lb Weight Capacity</i>	
Angle Adjustable Flip Under [R2FPL012]	£224
<i>Not available with footrest widths less than 9.5". Not available with 350 lb Weight Capacity</i>	
High Mount Angle Adjustable Flip Under [R2FPL011]	£272
<i>Not available with footrest widths greater than 15.5" or 350 lb Weight Capacity. See Reference G on page 6.</i>	

Footrest Protector

Footrest Protector [R2ACC020]	N/C
<i>Not available with Angle Adjustable Flip Under, High Mount Angle Adjustable Flip Under or High Mount Aluminum Angle Adjustable Footplates</i>	

Seat to Footrest Length

Measured from front edge of seat sling to top, rear of footrest. Footrest length at least 2.5" shorter than Front Seat Height recommended. See reference H on page 6.

8"	8.5"	9"	9.5"	10"	10.5"
11"	11.5"	12"	12.5"	13"	13.5"
14"	14.5"	15"	15.5"	16"	16.5"
17"	17.5"	18"	18.5"	19"	19.5"

Rear Seat Height and Wheels

Wheel Size / Rear Seat Height

Select Rear Seat Height in row corresponding to desired Wheel Size. For proper setup, selection is required when omitting rear wheels. Measured from floor to top of seat tube at front of back post. Rear Seat Height is custom to needs. See reference I on page 6

Wheel Size	Rear Seat Height													
	13.5" ¹	14"	14.5"	15"	15.5"	16"	16.5"	17"	17.5"	18"	18.5"	19"	19.5"	20"
22"														
24"														
25"														
26"														

¹ Available only with Low Profile tires

Camber Bar Tower

Standard Tower [R2CM01]	N/C
Mx Tower [R2CM02]	N/C

Standard Tower allows rear seat height adjustment up or down in 1/4" increments. 1" adjustment available; 1/2" in lowest end of range. Mx Tower required with select low seat heights or when large rear seat height adjustment range is required. Provides 3" adjustment range in 1/4" increments.

Rear Wheel

		Rear Wheel Size				
		Price	22"	24"	25"	26"
Rear Wheel Type	Spoke	N/C	[R2WHE022]	[R2WHE024]		
	Mag	N/C	[R2WHE122]	[R2WHE124]		
	Maxx Spoke	N/C	[R2WH1122]	[R2WH1124]	[R2WH1125]	[R2WH1126]
	Spinergy Spox	£768	[R2WHE422]	[R2WHE424]	[R2WHE425]	[R2WHE426]
	Spinergy LX	£816	[R2WHE522]	[R2WHE524]	[R2WHE525]	[R2WHE526]

Spinergy Options

Select spoke color and handrim tab for Spinergy wheels

Black Spokes [R2SSC001]	N/C
Red Spokes [R2SSC002]	Blue Spokes [R2SSC005]
White Spokes [R2SSC003]	£91

Handrims

Standard Handrims

Aluminum Anodized [R2HRM001]	N/C
Plastic Coated [R2HRM002]	£121
Projection (8 Vertical Projections) [R2HRM003]	£333

Available only with full profile tires in combination with 22", 24", and 26" MAXX Spoke wheels, or 22" and 24" Mag wheels

Ergonomic Handrims

Not available with 20" or 26" wheels. Tab mount only available with 24" wheels.

Ergo Lite Regular Grip [R2HRM005]	£272
Ergo Standard Regular Grip [R2HRM008]	£272
Ergo Lite Extra Grip [R2HRM006]	£272
Ergo Standard Extra Grip [R2HRM009]	£272
Ergo Lite Omit Grip [R2HRM004]	£272
Ergo Standard Omit Grip [R2HRM007]	£272

Handrim Options

Handrim Mount

If no selection is made, Standard Mount will be used

Close Mount [R2TAB001]	N/C
------------------------	-----

Micro-Grip

3M Micro Grip [R2ACC006]	£91
--------------------------	-----

Not available on Plastic Coated or Projection Handrims

Rear Tire Type

Full Poly 1 3/8" [R2TIR002]	N/C
<i>Not available on 25" or 26" wheel sizes</i>	
Low Poly 1" [R2TIR003]	N/C
<i>Not available on 25" or 26" wheel sizes or 350 lb Weight Capacity</i>	
Pneumatic with Airless Insert 1 3/8" [R2TIR004]	£110
<i>Not available on 25" wheel size</i>	
Sentinel Solid tire 1" [R2TIR005]	£110
<i>Not available on 22" Mag wheel</i>	
Street Pneumatic 1 3/8" [R2TIR001]	£110
<i>Not available on 25" wheel size</i>	
High Pressure 1" [R2TIR006]	N/C
<i>Not available on Mag wheel</i>	
Schwalbe Marathon Plus Evolution 1" [R2TIR010]	£110
<i>Not available on Mag wheel or 22" wheel size</i>	
Mountain 2" Black [R2TIR009]	£151
<i>Available only on 22", 24" and 25" MAXX Spoke and Spinergy wheels</i>	
Schwalbe HS-127 1 3/8" [R2TIR017]	£110
<i>Available only on 24" wheel size. Not available on Mag wheel.</i>	

Axle

Quick Release [R2AXO003]	N/C
Quad Quick Release [R2AXO001]	£182

Not available with Mag or Spinergy wheels

Center of Gravity Presets

Measured from front of back post to center of rear axle. 0" indicates axle will be directly under back post. See reference J on page 6.

Amputee Positions				
[TXCGA150]	[TXCGA100]	[TXCGA050]	[TXCGS000]	[TXCGS025]
Standard Positions				
[TXCGS050]	[TXCGS075]	[TXCGS100]	[TXCGS125]	[TXCGS150]
[TXCGS175]	[TXCGS200]	[TXCGS225]	[TXCGS250]	[TXCGS275]
[TXCGS300]	[TXCGS325]	[TXCGS350]	[TXCGS375]	[TXCGS400]
[TXCGS425]	[TXCGS450]	[TXCGS475]	[TXCGS500]	

Rear Wheel Spacing

Measured from outside of back to inside of rear tire. Adjustable out 1/2" from setting. 1.5" minimum Rear Wheel Spacing with Height Adjustable T-Arm Armrest. 1" minimum Rear Wheel Spacing with 0° Rear Wheel Camber. See reference K on page 6.

0.75" [R2WS075]	1" [R2WS100]	1.25" [R2WS125]	1.5" [R2WS150]
1.75" [R2WS175]			N/C

Rear Wheel Camber

0° [R2WC00]	2° [R2WC02]	4° [R2WC04]	6° [R2WC06]	
8° [R2WC08]				N/C

Wheel Locks

Wheel Lock

Push to Lock [R2WLK002]	N/C
Push to Lock Flush Mount [R2WLK001]	£67
Pull to Lock [R2WLK003]	N/C
Under Seat Scissor Lock [R2WLK007]	N/C
Low Profile Scissor Lock [R2WLK012]	N/C
<i>Not available with Low Poly, Full Poly or Sentinel Solid tires</i>	
Short Thro Scissor Lock [R2WLK006]	£121
Omit [R2WLK008]	N/C

Wheel Lock Options

Extension Handles [R2WLK004]	£49
<i>Available only with Push to Lock or Pull to Lock Wheel Locks</i>	
Grade Aids [R2WLK009]	£145
<i>Available only with Push to Lock Wheel Locks and Pneumatic or Pneumatic with Airless Insert tires</i>	

Armrests

Armrest Type	Length	Price
Height Adjustable T-Arm Tall (11.5"-15") ¹	Desk [R2ARM039]	£212
	Full [R2ARM040]	
Height Adjustable T-Arm (8.5"-13.5") ¹	Desk [R2ARM003]	£212
	Full [R2ARM004]	
Height Adjustable T-Arm Low (6.5"-10.5") ¹	Desk [R2ARM022]	£212
	Full [R2ARM023]	

¹ 1.5" minimum Rear Wheel Spacing
² Includes direct mount receiver

Armrest Type	Price
Swing Away (8.5"-12.5") [R2ARM002]	£140
Swing Away with Side Guard (8.5"-12.5") ² [R2ARM018]	£140
Omit Armrests [R2ARM007]	N/C

Side Guards

Direct Mount Side Guards

Not available with T-Arms or Swing Away with Side Guard Armrests

Composite Side Guards - Direct Mount Receiver [R2SDG002]	N/C
Aluminum Straight - Black Matte [R2SDG024]	£212
Aluminum Straight - Match Frame Color [R2SDG025]	£272
Aluminum Adjustable Fender with Foam Pad - Black Matte [R2SDG014]	£296
<i>Not available with Adjustable Depth greater than 1/2"</i>	
Aluminum Adjustable Fender with Foam Pad - Match Frame Color [R2SDG015]	£357
<i>Not available with Adjustable Depth greater than 1/2"</i>	
Carbon Fiber Adjustable Fender - Direct Mount Receiver [R2SDG013]	£393
<i>Not available with Adjustable Depth greater than 1/2"</i>	

Clamp-On Mount Side Guards

Not available with T-Arms or Swing Away with Side Guard Armrests

Composite Side Guards - Clamp-On Receiver [R2SDG004]	N/C
Aluminum Straight - Black Matte [R2SDG020]	£212
Aluminum Straight - Match Frame Color [R2SDG023]	£272
Aluminum Adjustable Fender with Foam Pad - Black Matte [R2SDG017]	£296
Aluminum Adjustable Fender w/Foam Pad - Clamp-On Receiver [R2SDG018]	£357
Carbon Fiber Adjustable Fender - Clamp-On Receiver [R2SDG016]	£393

Accessories

Anti-Tips

Rear Anti-Tips [R2ATP001]	£91
Omit [R2ATP000]	N/C

Calf Strap

Velcro Adjustable [R2ACC002]	N/C
Padded Velcro Adjustable [R2ACC003]	£67
Bodypoint Padded AeroMesh [R2ACC010]	£103

Pelvic Positioning Belt

1 1/2" Auto Buckle Pelvic Positioning Belt [R2POB001]	£55
1 1/2" Padded Auto Buckle Pelvic Positioning Belt [R2POB002]	£151
Small Bodypoint Evoflex (16-24") [R2POB007]	£188
Medium Bodypoint Evoflex (19-29") [R2POB008]	£188
Large Bodypoint Evoflex (25-37") [R2POB009]	£188

Other Accessories

Tool Kit [R2ACC005]	£20
Removable Underseat Pouch [R2ACC015]	£30
Ki Mobility Backpack [R2ACC009]	£100
Luggage Carrier [R2ACC025]	£73
Black Neoprene Impact Guards [R2IMG001]	£86
Under Chair Carrier and Removable Bag [R2ACC001]	£30
<i>Not available with scissor locks</i>	
Cane and Crutch Holder [R2ACC007]	£91

Belt Mounting Kit

Ki Mobility Belt Mounting Clamp [R2ACC023]	£37
Bodypoint Belt Mounting Bracket Kit [R2ACC013]	£103

Frame Color

Gloss Finish

Black [R2FCS001]	N/C
Blue Skies [R2FCS018]	N/C
Bubblegum Pink [R2FCS011]	N/C
Candy Blue [R2FCS012]	N/C
Candy Purple [R2FCS010]	N/C
Candy Red [R2FCS003]	N/C
Copper [R2FCS009]	N/C
Granny Smith [R2FCS008]	N/C
Metallic Orange [R2FCS032]	N/C
Silver [R2FCS026]	N/C
Yellow [R2FCS029]	N/C

Matte Finish

Charcoal Gray [R2FCS031]	£110
Emerald Green [R2FCS019]	£110
Iridescent Burgundy [R2FCS021]	£110
Iridescent Navy [R2FCS016]	£110
Military Green [R2FCS025]	£110
Retro Red [R2FCS022]	£110
Royal Blue [R2FCS023]	£110
Shadow Black [R2FCS015]	£110
Smooth White [R2FCS024]	£110

References

A. Seat Width *See Diagram 1*

Measured from outside of frame tube on one side to outside of frame tube on other side.

B. Seat Depth *See Diagram 3*

Measured from front of back posts to front edge of seat sling. Seat sling starts at beginning of bend at front of frame.

C. Sling Position *See Diagram 2*

In Performance Position (0°), sling is at front of frame to keep chair shorter and more maneuverable. The 1" or 2" seat sling position extends the frame in front of the sling by that amount. Adding more frame can improve stability and provide support to aid in transfers.

D. Front Frame Bend *See Diagram 3*

Angle between front seat tube and ground.

E. Footrest Taper *See Diagram 1*

Indicates bend of front frame creating footrest. Taper is measured from outside of seat frame to outside of front tube.

F. Front Seat Height *See Diagram 1*

Measured from floor to top of seat tube at front of seat sling.

G. Footrest Width *See Diagram 1*

Measured from inside of front tube to inside of front tube on other side. Footrest width is listed below for each seat width and taper.

Seat Width	Straight	1" Taper	2" Taper
Outside Measurement	2" Narrower than Seat Width	4" Narrower than Seat Width	6" Narrower than Seat Width
14"	11.5"	9.5"	
15"	12.5"	10.5"	
16"	13.5"	11.5"	9.5"
17"	14.5"	12.5"	10.5"
18"	15.5"	13.5"	11.5"
19"	16.5"	14.5"	12.5"
20"	17.5"	15.5"	13.5"
21"	18.5"	16.5"	14.5"
22"	19.5"	17.5"	15.5"

H. Seat to Footrest Length *See Diagram 3*

Measured from front edge of seat sling to top rear of footrest. Footrest length of at least 2.5" shorter than front seat height recommended.

I. Rear Seat Height *See Diagram 3*

Measured from floor to top of seat tube at front of back post. Rear seat height is custom to needs.

J. Center of Gravity Preset *See Diagram 3*

Measured from front of back post to center of rear axle. 0° indicates axle will be directly under back post.

K. Rear Wheel Spacing *See Diagram 1*

Measured from outside of seat back to inside of rear tire. Adjustable out .5" from setting.

Diagram 1




Diagram 2





Diagram 3



Axiom G

		Width															
		14"		15"		16"		17"		18"		19"		20"		22"	
		Cushion	Outer Cover	Cushion	Outer Cover	Cushion	Outer Cover	Cushion	Outer Cover	Cushion	Outer Cover	Cushion	Outer Cover	Cushion	Outer Cover	Cushion	Outer Cover
Retail		£85	£70	£85	£70	£85	£70	£85	£70	£85	£70	£85	£70	£85	£70	£85	£70
Height	14"																
	15"																
	16"																
	17"																
	18"																
	19"																
	20"																
	21"																
	22"																

Axiom S

		Width															
		14"		15"		16"		17"		18"		19"		20"		22"	
		Cushion	Outer Cover	Cushion	Outer Cover	Cushion	Outer Cover	Cushion	Outer Cover	Cushion	Outer Cover	Cushion	Outer Cover	Cushion	Outer Cover	Cushion	Outer Cover
Retail		£120	£70	£120	£70	£120	£70	£120	£70	£120	£70	£120	£70	£120	£70	£120	£70
Height	14"																
	15"																
	16"																
	17"																
	18"																
	19"																
	20"																
	21"																
	22"																

Axiom P

		Width															
		14"		15"		16"		17"		18"		19"		20"		22"	
		Cushion	Outer Cover	Cushion	Outer Cover	Cushion	Outer Cover	Cushion	Outer Cover	Cushion	Outer Cover	Cushion	Outer Cover	Cushion	Outer Cover	Cushion	Outer Cover
Retail		£170	£70	£170	£70	£170	£70	£170	£70	£170	£70	£170	£70	£170	£70	£170	£70
Height	14"																
	15"																
	16"																
	17"																
	18"																
	19"																
	20"																
	22"																

Axiom SP - Visco

		Width																										
		14"			15"			16"			17"			18"			19"			20"			21"			22"		
		Cushion	Cover		Cushion	Cover		Cushion	Cover		Cushion	Cover		Cushion	Cover		Cushion	Cover		Cushion	Cover		Cushion	Cover				
Retail		£350	£70	£350	£70	£350	£70	£350	£70	£350	£70	£350	£70	£350	£70	£350	£70	£350	£70	£350	£70	£350	£70					
Height	14"																											
	15"																											
	16"																											
	17"																											
	18"																											
	19"																											
	20"																											
	22"																											

Axiom SP - Fluid

		Width																										
		14"			15"			16"			17"			18"			19"			20"			21"			22"		
		Cushion	Cover		Cushion	Cover		Cushion	Cover		Cushion	Cover		Cushion	Cover		Cushion	Cover		Cushion	Cover		Cushion	Cover				
Retail		£320	£70	£320	£70	£320	£70	£320	£70	£320	£70	£320	£70	£320	£70	£320	£70	£320	£70	£320	£70	£320	£70					
Height	14"																											
	15"																											
	16"																											
	17"																											
	18"																											
	19"																											
	20"																											
	22"																											


Solid Seat Insert

Weight capacity: 350 lbs

		Width								
		14"	15"	16"	17"	18"	19"	20"	21"	22"
		Cushion	Cushion	Cushion	Cushion	Cushion	Cushion	Cushion	Cushion	Cushion
Retail		£55	£55	£55	£55	£55	£55	£55	£55	£55
Height	14"									
	15"									
	16"									
	17"									
	18"									
	19"									
	20"									
	21"									
	22"									

Axiom Back Width	Fits Back Frame Widths	
	Inner Slots	Outer Slots
12" (30.5cm)	11" (28.0cm)	12" -14" (30.5cm -35.5cm)
14" (35.5cm)	12" -13" (30.5cm -33.0cm)	14" -16" (35.5cm -40.5cm)
16" (40.5cm)	14" -15" (35.5cm -38.0cm)	16" -18" (40.5cm -45.5cm)
18" (45.5cm)	16" - 17" (40.5cm - 43.0cm)	18" - 20" (45.5cm - 50.5cm)
20" (51.0cm)	18" -19" (46.0cm -48.5cm)	20" -22" (50.5cm -55.5cm)
22" (56.0cm)	20" -21" (51.0cm -53.5cm)	22" -24" (56.0cm -61.0cm)

Hardware Description	Code
Improved DEX Hardware, removable adjustable hardware (300 lb max user weight)	DEX
New BTX Hardware, removable adjustable hardware not recommended for, or installed on Ultralight Rigid Frame Wheelchairs (400 lb max user weight)	BTX
New Fixed Hardware, non-removeable hardware up to (300 lb max user weight)	FXT



Axiom Back

Hardware is included with back. Low Back selections are not available with headrests, modifications, or posture supports. Weight capacity: 300 lbs.

Axiom Posterior Back (PB)

		Width						
		12"	14"	16"	18"	20"	22"	
Retail		£550	£550	£550	£550	£550	£550	
Height	Low	8"	[AXPB1208]	[AXPB1408]	[AXPB1608]	[AXPB1808]	[AXPB2008]	
		10"		[AXPB1410]	[AXPB1610]	[AXPB1810]	[AXPB2010]	[AXPB2210]
		12"	[AXPB1212]	[AXPB1412]	[AXPB1612]	[AXPB1812]	[AXPB2012]	
	Mid	14"		[AXPB1414]	[AXPB1614]	[AXPB1814]	[AXPB2014]	[AXPB2214]
		16"	[AXPB1216]	[AXPB1416]	[AXPB1616]	[AXPB1816]	[AXPB2016]	
		18"		[AXPB1418]	[AXPB1618]	[AXPB1818]	[AXPB2018]	[AXPB2218]

Back Hardware

DEX Hardware [DEX] BTX Hardware [BTX] FXT Hardware [FXT] _____ N/C

Lateral Support for AXPB Backs

Lateral Trunk Supports for AXPB Backs							
Pair	Left	Right	Description	Item Code	Pair	Each	
			3 x 4 Lateral Pads with Fixed Mounting Hardware for Axiom Posterior Back	PBLF3x4P	£295	£148	
			4 x 5 Lateral Pads with Fixed Mounting Hardware for Axiom Posterior Back	PBLF4x5P	£295	£148	
			3 x 4 Lateral Pads with Swing Away Mounting Hardware for Axiom Posterior Back	PBLS3x4P	£628	£314	
			4 x 5 Lateral Pads with Swing Away Mounting Hardware for Axiom Posterior Back	PBLS4x5P	£628	£314	
			3 x 4 Lateral Pads with Adjustable Swing Away Mounting Hardware	PBLAS3x4P	£628	£314	
			4 x 5 Lateral Pads with Adjustable Swing Away Mounting Hardware	PBLAS4x5P	£628	£314	

Axiom Posterior Lateral Back (PL)

		Width						
		12"	14"	16"	18"	20"	22"	
Retail		£550	£550	£550	£550	£550	£550	
Height	Low	8"	[AXPL1208]	[AXPL1408]	[AXPL1608]	[AXPL1808]	[AXPL2008]	
		10"		[AXPL1410]	[AXPL1610]	[AXPL1810]	[AXPL2010]	[AXPL2210]
		12"	[AXPL1212]	[AXPL1412]	[AXPL1612]	[AXPL1812]	[AXPL2012]	
	Mid	14"		[AXPL1414]	[AXPL1614]	[AXPL1814]	[AXPL2014]	[AXPL2214]
		16"	[AXPL1216]	[AXPL1416]	[AXPL1616]	[AXPL1816]	[AXPL2016]	
		18"		[AXPL1418]	[AXPL1618]	[AXPL1818]	[AXPL2018]	[AXPL2218]
Tall	20"	[AXPL1220]	[AXPL1420]	[AXPL1620]	[AXPL1820]	[AXPL2020]		
	22"		[AXPL1422]	[AXPL1622]	[AXPL1822]	[AXPL2022]	[AXPL2222]	

Back Hardware

DEX Hardware [DEX] BTX Hardware [BTX] FXT Hardware [FXT] _____ N/C

Lateral Support for AXPL

Lateral Trunk Supports for AXPL & AXLB Backs						
Pair	Left	Right	Description	Item Code	Pair	Each
			3 x 4 Lateral Pads with Fixed Mounting Hardware for Axiom Posterior Lateral Back	PLLF3x4P	£295	£295
			4 x 5 Lateral Pads with Fixed Mounting Hardware for Axiom Posterior Lateral Back	PLLF4x5P	£295	£295
			3 x 4 Lateral Pads with Swing Away Mounting Hardware for Axiom Posterior Lateral Back	PLLS3x4P	£628	£314
			4 x 5 Lateral Pads with Swing Away Mounting Hardware for Axiom Posterior Lateral Back	PLLS4x5P	£628	£314
			3 x 4 Lateral Pads with Adjustable Swing Away Mounting Hardware	PLLAS3x4P	£628	£314
			4 x 5 Lateral Pads with Adjustable Swing Away Mounting Hardware	PLLAS4x5P	£628	£314

Axiom Deep Lateral Back (DL)

		Width						
		12"	14"	16"	18"	20"	22"	
Retail		£550	£550	£550	£550	£550	£550	
HCPC		E2620	E2620	E2620	E2620	E2620	E2621	
Height	Mid	14"	[AXDL1214]	[AXDL1414]	[AXDL1614]	[AXDL1814]	[AXDL2014]	[AXDL2214]
		16"		[AXDL1416]	[AXDL1616]	[AXDL1816]	[AXDL2016]	
		18"	[AXDL1218]	[AXDL1418]	[AXDL1618]	[AXDL1818]	[AXDL2018]	[AXDL2218]
	Tall	20"	[AXDL1220]	[AXDL1420]	[AXDL1620]	[AXDL1820]	[AXDL2020]	
		22"		[AXDL1422]	[AXDL1622]	[AXDL1822]	[AXDL2022]	[AXDL2222]

Back Hardware

DEX Hardware [DEX] BTX Hardware [BTX] FXT Hardware [FXT] _____ N/C

Headrest

Description	Std. Straight Removable Adjustable Hardware		Multi-Axis Removable Adjustable Hardware		Posalinc Removable Adjustable Hardware		HCPCS
	Item Code	Retail	Item Code	Retail	Item Code	Retail	
Comfort Soft Headrest 8"	[CS08SS]	£268	[CS08MA]	£268	[CS08PL]	£268	E0955 / E1028
Comfort Soft Headrest 10"	[CS10SS]	£268	[CS10MA]	£268	[CS10PL]	£268	E0955 / E1028
Comfort Soft Headrest 14"	[CS14SS]	£268	[CS14MA]	£268	[CS14PL]	£268	E0955 / E1028
Contour Positioning Pad, 8" X 6.25"	[CP08SS]	£303	[CP08MA]	£303	[CP08PL]	£303	E0955 / E1028
Contour Positioning Pad, 9" X 6.5"	[CP09SS]	£303	[CP09MA]	£303	[CP09PL]	£303	E0955 / E1028
Contour Positioning Pad, 10" X 7.5"	[CP10SS]	£303	[CP10MA]	£303	[CP10PL]	£303	E0955 / E1028
Tri-soft Headrest 10"	[TS10SS]	£386	[TS10MA]	£386	[TS10PL]	£386	E0955 / E1028
Tri-soft Headrest 14"	[TS14SS]	£386	[TS14MA]	£386	[TS14PL]	£386	E0955 / E1028

Description	Item Code	Retail	HCPCS
Universal Headrest Mounting Plate	[UMP]	£43	K0108

Not required for Axiom Headrests

Anterior Thoracic Supports

Selection of a modification from the Axiom Back Modifications section required in order to mount posture supports. Not available with Low Backs.

Anterior Thoracic Supports

Description	Item Code	HCPCS	Retail
Dynaform Chest Support: 2-Piece with Hook&Loop, 3" High Pad, Small	[14941]	E0960	£100
Dynaform Chest Support: 2-Piece with Hook&Loop, 4" High Pad, Medium	[14942]	E0960	£100
Dynaform Chest Support: 2-Piece with Hook&Loop, 5" High Pad, Large	[14943]	E0960	£100
Dynaform Chest Support: 1-Piece with Buckle, 3" Height x 20" Width, Small	[14924]	E0960	£100
Dynaform Chest Support: 1-Piece with Buckle, 4" Height x 20" Width, Medium	[14927]	E0960	£100
Dynaform Chest Support: 1-Piece with Buckle, 4" Height x 22" Width, Large	[14928]	E0960	£100
Dynaform Anterior Trunk Support: Standard, Dynamic, X-Small	[14851]	E0960	£100
Dynaform Anterior Trunk Support: Standard, Dynamic, Small	[14852]	E0960	£100
Dynaform Anterior Trunk Support: Standard, Dynamic, Medium	[14853]	E0960	£100
Dynaform Anterior Trunk Support: Standard, Dynamic, Large	[14854]	E0960	£100
Dynaform Anterior Trunk Support: Slim Cut, Dynamic, X-Small	[14856]	E0960	£100
Dynaform Anterior Trunk Support: Slim Cut, Dynamic, Small	[14857]	E0960	£100
Dynaform Anterior Trunk Support: Slim Cut, Dynamic, Medium	[14858]	E0960	£100
Dynaform Anterior Trunk Support: Slim Cut, Dynamic, Large	[14859]	E0960	£100
AirLogic Anterior Trunk Support: Standard, Structure, Small	[14282]	E0960	£300
AirLogic Anterior Trunk Support: Standard, Structure, Medium	[14283]	E0960	£300
AirLogic Anterior Trunk Support: Standard, Structure, Large	[14284]	E0960	£300
AirLogic Anterior Trunk Support: Standard, Structure, X-Large	[14285]	E0960	£300
AirLogic Anterior Trunk Support: Slim Cut, Dynamic, Small	[14257]	E0960	£300
AirLogic Anterior Trunk Support: Slim Cut, Dynamic, Medium	[14258]	E0960	£300
AirLogic Anterior Trunk Support: Slim Cut, Dynamic, Large	[14259]	E0960	£300
AirLogic Anterior Trunk Support: Slim Cut, Dynamic, X-Large	[14260]	E0960	£300

Accessories

Adjustable Shoulder Guides (K0108) [AXBSG] <i>Includes Anterior Thoracic Mounting Modification. Not available with Low Backs.</i>	£43
Extra Lumbar Pad (K0108) [LUMPAD] <i>Available only with Axiom Lateral Back and Axiom Deep Lateral back</i>	N/C
Privacy Flap (K0108) [PRFLAP]	£44

Axiom Back Modifications

Axiom Back Modifications

Extended lead times will apply for all modifications. Makes back non-refundable.

Anterior Thoracic Mounting Modification [AXSMM] <i>Provides holes for Anterior Thoracic Supports and Shoulder Guides</i>	N/C
---	-----

Gloss Finish

If no selection is made, standard Shadow Black will be used with no additional lead time

Black Back Shell [AXMBLK]	N/C
Blue Skies Back Shell [AXMSKY]	N/C
Bubblegum Pink Back Shell [AXMPNK]	N/C
Candy Blue Back Shell [AXMCBL]	N/C
Candy Purple Back Shell [AXMPUR]	N/C
Candy Red Back Shell [AXMRED]	N/C
Copper Back Shell [AXMCOP]	N/C
Granny Smith Back Shell [AXMGSM]	N/C
Metallic Orange Back Shell [AXMMOR]	N/C
Yellow Back Shell [AXMYLW]	N/C
Silver Back Shell [AXMSIL]	N/C
Teal Back Shell [AXMTEA]	N/C

Matte Finish

If no selection is made, standard Shadow Black will be used with no additional lead time

Military Green Back Shell [AXMAGM]	N/C
Charcoal Gray Back Shell [AXMCGM]	N/C
Emerald Green Back Shell [AXMEGM]	N/C
Iridescent Burgundy Back Shell [AXMIBM]	N/C
Iridescent Navy Back Shell [AXMINM]	N/C
Retro Red Back Shell [AXMRRM]	N/C
Royal Blue Back Shell [AXMRBM]	N/C
Smooth White Back Shell [AXMSWM]	N/C