



Phone (800) 981-1540 | Fax (715) 254-0996
 sales@kimobility.com | www.kimobility.com

Company: _____

Ship To: _____

Account Number: _____ Date: _____

Address: _____

Quote PO Number: _____

City: _____ State: _____

RTS: _____

Zip: _____ Phone: _____

Marked For: _____

Email/Fax: _____

Select if chair quote has: Additional Ki Products with Different Quote/Order Number **Quote/Order Number:** _____ Effective 08/15/2022

- CR45 - Transit Included** E1161 **\$3295**
300 lb weight capacity. Tilt range: 0° anterior to 45° posterior.
- CR45 - Transit Not Included** E1161 **\$2975**
300 lb weight capacity. Tilt range: 0° anterior to 45° posterior.

HCPCS codes provided are based on code verification by the PDAC or our interpretation of the code definitions. Proper use of HCPCS codes does not ensure coverage or payment. For coverage information, verify the appropriate payer's coverage policy.

Frame Width

Frame Width

Select box corresponding to desired Frame Width

Frame Width								
<input type="checkbox"/> 14 N/C	<input type="checkbox"/> 15 N/C	<input type="checkbox"/> 16 N/C	<input type="checkbox"/> 17 N/C	<input type="checkbox"/> 18 N/C	<input type="checkbox"/> 19 \$500 E2201	<input type="checkbox"/> 20 \$500 E2201	<input type="checkbox"/> 21 \$500 E2201	<input type="checkbox"/> 22 \$500 E2201

Seating Setup

Make one selection from either Option 1 at right or Option 2 below based on desired seating setup

OPTION 2: Setup By Frame Depth

Select Frame Depth within desired Frame Depth Range

Frame Depth Range	Frame Depth				
<input type="checkbox"/> Short 14" - 18"	<input type="checkbox"/> 14 N/C	<input type="checkbox"/> 15 N/C	<input type="checkbox"/> 16 N/C	<input type="checkbox"/> 17 N/C	<input type="checkbox"/> 18 N/C
<input type="checkbox"/> Medium 16" - 20"	<input type="checkbox"/> 16 N/C	<input type="checkbox"/> 17 N/C	<input type="checkbox"/> 18 N/C	<input type="checkbox"/> 19 N/C	<input type="checkbox"/> 20 \$605 E2203
<input type="checkbox"/> Long 18" - 22"	<input type="checkbox"/> 18 N/C	<input type="checkbox"/> 19 N/C	<input type="checkbox"/> 20 \$605 E2203	<input type="checkbox"/> 21 \$605 E2203	<input type="checkbox"/> 22 \$1025 E2204


OPTION 1: Setup By Ki Mobility

Select Seat Depth corresponding to desired Frame Depth Range row and Backrest Offset column.

Note: All Seat Depths in a column and within the same frame depth range are achievable with future adjustment.

Frame Depth Range	Backrest Offset A						
	<input type="checkbox"/> -2	<input type="checkbox"/> -1	<input type="checkbox"/> Flush	<input type="checkbox"/> +1	<input type="checkbox"/> +2	<input type="checkbox"/> +3	<input type="checkbox"/> +4
<input type="checkbox"/> Short	Seat Depth						
	<input type="checkbox"/> 16	<input type="checkbox"/> 15	<input type="checkbox"/> 14				
	<input type="checkbox"/> 17	<input type="checkbox"/> 16	<input type="checkbox"/> 15	<input type="checkbox"/> 14			
	<input type="checkbox"/> 18	<input type="checkbox"/> 17	<input type="checkbox"/> 16	<input type="checkbox"/> 15	<input type="checkbox"/> 14		
	<input type="checkbox"/> 19	<input type="checkbox"/> 18	<input type="checkbox"/> 17	<input type="checkbox"/> 16	<input type="checkbox"/> 15	<input type="checkbox"/> 14	
	<input type="checkbox"/> 20 ¹	<input type="checkbox"/> 19	<input type="checkbox"/> 18	<input type="checkbox"/> 17	<input type="checkbox"/> 16	<input type="checkbox"/> 15	<input type="checkbox"/> 14
<input type="checkbox"/> Medium	<input type="checkbox"/> 18	<input type="checkbox"/> 17	<input type="checkbox"/> 16	<input type="checkbox"/> 15	<input type="checkbox"/> 14		
	<input type="checkbox"/> 19	<input type="checkbox"/> 18	<input type="checkbox"/> 17	<input type="checkbox"/> 16	<input type="checkbox"/> 15	<input type="checkbox"/> 14	
	<input type="checkbox"/> 20 ¹	<input type="checkbox"/> 19	<input type="checkbox"/> 18	<input type="checkbox"/> 17	<input type="checkbox"/> 16	<input type="checkbox"/> 15	<input type="checkbox"/> 14
	<input type="checkbox"/> 21 ¹	<input type="checkbox"/> 20 ¹	<input type="checkbox"/> 19	<input type="checkbox"/> 18	<input type="checkbox"/> 17	<input type="checkbox"/> 16	<input type="checkbox"/> 15
	<input type="checkbox"/> 22 ²	<input type="checkbox"/> 21 ¹	<input type="checkbox"/> 20 ¹	<input type="checkbox"/> 19	<input type="checkbox"/> 18	<input type="checkbox"/> 17	<input type="checkbox"/> 16
<input type="checkbox"/> Long	<input type="checkbox"/> 20 ¹	<input type="checkbox"/> 19	<input type="checkbox"/> 18	<input type="checkbox"/> 17	<input type="checkbox"/> 16	<input type="checkbox"/> 15	<input type="checkbox"/> 14
	<input type="checkbox"/> 21 ¹	<input type="checkbox"/> 20 ¹	<input type="checkbox"/> 19	<input type="checkbox"/> 18	<input type="checkbox"/> 17	<input type="checkbox"/> 16	<input type="checkbox"/> 15
	<input type="checkbox"/> 22 ²	<input type="checkbox"/> 21 ¹	<input type="checkbox"/> 20 ¹	<input type="checkbox"/> 19	<input type="checkbox"/> 18	<input type="checkbox"/> 17	<input type="checkbox"/> 16
		<input type="checkbox"/> 22 ²	<input type="checkbox"/> 21 ¹	<input type="checkbox"/> 20 ¹	<input type="checkbox"/> 19	<input type="checkbox"/> 18	<input type="checkbox"/> 17
			<input type="checkbox"/> 22 ²	<input type="checkbox"/> 21 ¹	<input type="checkbox"/> 20 ¹	<input type="checkbox"/> 19	<input type="checkbox"/> 18

¹ Price is \$605 E2203 ² Price is \$1025 E2204



Seat Pan

- Growing Adjustable Aluminum Seat Pan E2231 _____ \$205
- Omit Seat Pan _____ N/C

Frame Options

Transit

- Transit System _____ \$320
- Omit Transit _____ N/C

Tilt Mechanism

Hand

Hand Tilt - Dual _____ N/C

Foot

Foot Tilt - Right _____ N/C

Foot Tilt - Left _____ N/C

Backrest

Back Post

Stroller Handle Back Post (Folding)

Short 20" _____ Std

Tall 24" _____ N/C

Fixed Height with Adjustable Handle (Folding)

Short 22" _____ N/C

Tall 26" _____ N/C

Adjustable Height with Adjustable Handle (Folding)

Short 18"-25" _____ \$25

Not available with O₂ Holder

18" 19" 20" 21"

22" 23" 24" 25"

Tall 20"-27" _____ \$25

20" 21" 22" 23"

24" 25" 26" 27"

Reclining Back with Adjustable Handle (Non-Folding) E1225

Not available with O₂ Holder. If no Gas Spring selection is made, Standard Spring will be used.

Short 22" _____ \$1000

Tall 26" _____ \$1000

Gas Spring: Standard Spring

Heavy Duty Spring

Stroller Handle Back Post Option

Available only with Stroller Handle Back Post

Removable Stroller Handle Extension _____ \$85

Other Backrest Options

Non Folding Back _____ N/C
Standard with Reclining Back

Adjustable Height Rigidizer Bar _____ \$75

Footrest

Hanger / Front Frame Bend

Hanger	60°	70°	80°	90°		Retail
Extension Tube Mount	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/> ¹			Std
Front Mount		<input type="checkbox"/>	<input type="checkbox"/> ¹	<input type="checkbox"/>		N/C
Pro Elevating Leg Rests (PRO ELR) E0990					<input type="checkbox"/>	\$300
Pediatric Pro Elevating Leg Rests (Peds PRO ELR) E0990					<input type="checkbox"/>	\$300
Omit Hangers					<input type="checkbox"/>	-\$20

¹ 7" and 8" Casters may interfere with less than 5° tilt

PRO ELR Cover

Available only on PRO ELR

Gel Pad Covers (Pair) K0108 _____ \$390

3" Extension

Available only with Front Mount Hangers.

3" Extension _____ N/C

3" Extension (Provided as Parts) _____ N/C

Footrest continued

Footrest / Extension / Seat to Footrest Length

Select Extension and Seat to Footrest Length in row corresponding to desired Footrest. Footrest Length of at least 2.5" shorter than Front Seat Height recommended. Extension Tube Hanger measurements based on 70° Hanger Bend.

Footrest	Extension	Seat to Footrest Length													
Extension Tube Hangers (60°, 70°, 80°)	<input type="checkbox"/> Short	<input type="checkbox"/> 10	<input type="checkbox"/> 10.5	<input type="checkbox"/> 11	<input type="checkbox"/> 11.5	<input type="checkbox"/> 12	<input type="checkbox"/> 12.5	<input type="checkbox"/> 13	<input type="checkbox"/> 13.5						
	<input type="checkbox"/> Medium	<input type="checkbox"/> 11.5	<input type="checkbox"/> 12	<input type="checkbox"/> 12.5	<input type="checkbox"/> 13	<input type="checkbox"/> 13.5	<input type="checkbox"/> 14	<input type="checkbox"/> 14.5	<input type="checkbox"/> 15	<input type="checkbox"/> 15.5					
	<input type="checkbox"/> Long	<input type="checkbox"/> 14	<input type="checkbox"/> 14.5	<input type="checkbox"/> 15	<input type="checkbox"/> 15.5	<input type="checkbox"/> 16	<input type="checkbox"/> 16.5	<input type="checkbox"/> 17	<input type="checkbox"/> 17.5	<input type="checkbox"/> 18					
Front Mount Hangers ¹	Standard Selections										Requires 3" Extension				
	<input type="checkbox"/> Short (70°, 80°, 90°)	<input type="checkbox"/> 4 ⁴	<input type="checkbox"/> 4.5 ⁴	<input type="checkbox"/> 5	<input type="checkbox"/> 5.5	<input type="checkbox"/> 6						<input type="checkbox"/> 8.5	<input type="checkbox"/> 9	<input type="checkbox"/> 9.5	
		<input type="checkbox"/> 6.5	<input type="checkbox"/> 7	<input type="checkbox"/> 7.5	<input type="checkbox"/> 8						<input type="checkbox"/> 10	<input type="checkbox"/> 10.5	<input type="checkbox"/> 11		
	<input type="checkbox"/> Medium (70°, 80°, 90°)	<input type="checkbox"/> 4 ⁴	<input type="checkbox"/> 4.5 ⁴	<input type="checkbox"/> 5	<input type="checkbox"/> 5.5	<input type="checkbox"/> 6	<input type="checkbox"/> 6.5	<input type="checkbox"/> 7					<input type="checkbox"/> 10.5	<input type="checkbox"/> 11	<input type="checkbox"/> 11.5
		<input type="checkbox"/> 7.5	<input type="checkbox"/> 8	<input type="checkbox"/> 8.5	<input type="checkbox"/> 9	<input type="checkbox"/> 9.5	<input type="checkbox"/> 10						<input type="checkbox"/> 12	<input type="checkbox"/> 12.5	<input type="checkbox"/> 13
	<input type="checkbox"/> Long (70°, 80°)	<input type="checkbox"/> 5	<input type="checkbox"/> 5.5	<input type="checkbox"/> 6	<input type="checkbox"/> 6.5	<input type="checkbox"/> 7	<input type="checkbox"/> 7.5	<input type="checkbox"/> 8	<input type="checkbox"/> 8.5	<input type="checkbox"/> 12.5	<input type="checkbox"/> 13	<input type="checkbox"/> 13.5			
<input type="checkbox"/> 9		<input type="checkbox"/> 9.5	<input type="checkbox"/> 10	<input type="checkbox"/> 10.5	<input type="checkbox"/> 11	<input type="checkbox"/> 11.5	<input type="checkbox"/> 12					<input type="checkbox"/> 14	<input type="checkbox"/> 14.5	<input type="checkbox"/> 15	
Peds ² PRO ELR	<input type="checkbox"/> Low Pivot	<input type="checkbox"/> 7 ³	<input type="checkbox"/> 7.5 ³	<input type="checkbox"/> 8 ³	<input type="checkbox"/> 8.5	<input type="checkbox"/> 9	<input type="checkbox"/> 9.5	<input type="checkbox"/> 10	<input type="checkbox"/> 10.5	<input type="checkbox"/> 11	<input type="checkbox"/> 11.5				
	<input type="checkbox"/> High Pivot	<input type="checkbox"/> 7 ³	<input type="checkbox"/> 7.5 ³	<input type="checkbox"/> 8 ³	<input type="checkbox"/> 8.5	<input type="checkbox"/> 9	<input type="checkbox"/> 9.5	<input type="checkbox"/> 10							
PRO ELR Low Pivot Setting	<input type="checkbox"/> Short	10.5 - 15.5													
	<input type="checkbox"/> Long	13 - 19													
PRO ELR High Pivot Setting	<input type="checkbox"/> Short	8.5 - 17													
	<input type="checkbox"/> Long	11.5 - 17.5													

¹ See Minimum Front Seat Height chart below
² Comes standard with 3" Extension
³ Not available on chairs 15" wide or more
⁴ Available only with 90° Front Mount Hangers

Hanger Length	Short		Medium		Long	
	70°	80/90°	70°	80/90°	70°	80°
Hanger Bend	70°	80/90°	70°	80/90°	70°	80°
Minimum Front Seat Height with Standard Selection	13"	13"	13"	13.5"	14.5"	15.5"
Minimum Front Seat Height with 3" Extension	13.5"	14"	15.5"	16"	17"	18"


Footplate

- Composite _____ Std
Not available with Front Mount Hanger
- Composite Angle Adjustable k0040 _____ \$185
Not available with 14" Frame Width
- Aluminum Angle Adjustable k0040 _____ \$195
- One Piece Flip Up Angle Adjustable k0040 _____ \$190
Not available with Front Mount Hanger or PRO ELR

Hanger Release

If no selection is made 4-Way - In line Hanger Release will be used

- Classic _____ N/C
- 4-Way - In line _____ Std
See reference A below. Not available with PRO ELR or Peds PRO ELR.
- 4-Way - Out _____ N/C
See reference B below. Not available with PRO ELR or Peds PRO ELR.



Heel Loop

Not available with One Piece Flip Up Angle Adjustable Footplate

- Heel Loops E0951 _____ \$40
- Adjustable Heel Loops E0951 _____ \$50
- Omit Heel Loops _____ N/C

Calf Strap

- Velcro Adjustable k0038 _____ \$40
- Padded Velcro Adjustable k0038 _____ \$68
- Bodypoint Padded AeroMesh k0038 _____ \$118

Residual Limb Support

Residual Limb Support E1020

- Right and Left _____ \$810
- Left _____ \$405
- Right _____ \$405
- Omit Left Hanger _ -\$10
- Omit Right Hanger _ -\$10

Cover

- Smooth
- Air Mesh
- Smooth
- Air Mesh

Pad Size

- | | | | |
|---------------------------------|----------------------------------|---------------------------------|----------------------------------|
| <input type="checkbox"/> 5 x 8 | <input type="checkbox"/> 9 x 14 | <input type="checkbox"/> 5 x 8 | <input type="checkbox"/> 9 x 14 |
| <input type="checkbox"/> 7 x 10 | <input type="checkbox"/> 11 x 10 | <input type="checkbox"/> 7 x 10 | <input type="checkbox"/> 11 x 10 |
| <input type="checkbox"/> 9 x 10 | <input type="checkbox"/> 11 x 14 | <input type="checkbox"/> 9 x 10 | <input type="checkbox"/> 11 x 14 |

Seat Height, Casters and Rear Wheels

Seat Height / Caster Size / Rear Wheel Size

Select Seat Height corresponding to desired Rear Wheel and Caster Size

Rear Wheel Size	Caster Size					
	4"	5"	6"	6" x 2"	7"	8"
	Seat Height					
12" with Pneumatic or Pneumatic with Airless Insert Tire	<input type="checkbox"/> 15.5					
	<input type="checkbox"/> 16	<input type="checkbox"/> 16	<input type="checkbox"/> 16	<input type="checkbox"/> 16		
	<input type="checkbox"/> 16.5	<input type="checkbox"/> 16.5	<input type="checkbox"/> 16.5	<input type="checkbox"/> 16.5		
		<input type="checkbox"/> 17	<input type="checkbox"/> 17	<input type="checkbox"/> 17	<input type="checkbox"/> 17	
		<input type="checkbox"/> 17.5	<input type="checkbox"/> 17.5	<input type="checkbox"/> 17.5	<input type="checkbox"/> 17.5	
12" with Full Poly Tire	<input type="checkbox"/> 15	<input type="checkbox"/> 15				
	<input type="checkbox"/> 15.5	<input type="checkbox"/> 15.5				
	<input type="checkbox"/> 16	<input type="checkbox"/> 16	<input type="checkbox"/> 16	<input type="checkbox"/> 16		
	<input type="checkbox"/> 16.5	<input type="checkbox"/> 16.5	<input type="checkbox"/> 16.5	<input type="checkbox"/> 16.5		
16" with Pneumatic, Pneumatic with Airless Insert or Full Poly Tire	<input type="checkbox"/> 14					
	<input type="checkbox"/> 14.5					
	<input type="checkbox"/> 15	<input type="checkbox"/> 15				
	<input type="checkbox"/> 15.5	<input type="checkbox"/> 15.5 ¹				
		<input type="checkbox"/> 17		<input type="checkbox"/> 17 ²	<input type="checkbox"/> 17 ²	
		<input type="checkbox"/> 17.5	<input type="checkbox"/> 17.5	<input type="checkbox"/> 17.5	<input type="checkbox"/> 17.5	
		<input type="checkbox"/> 18	<input type="checkbox"/> 18	<input type="checkbox"/> 18	<input type="checkbox"/> 18	<input type="checkbox"/> 18
		<input type="checkbox"/> 18.5	<input type="checkbox"/> 18.5	<input type="checkbox"/> 18.5	<input type="checkbox"/> 18.5	<input type="checkbox"/> 18.5
		<input type="checkbox"/> 19	<input type="checkbox"/> 19	<input type="checkbox"/> 19	<input type="checkbox"/> 19	<input type="checkbox"/> 19
16" with Low Poly Tire	<input type="checkbox"/> 14					
	<input type="checkbox"/> 14.5					
		<input type="checkbox"/> 16				
	<input type="checkbox"/> 16.5	<input type="checkbox"/> 16.5	<input type="checkbox"/> 16.5	<input type="checkbox"/> 16.5		
		<input type="checkbox"/> 17	<input type="checkbox"/> 17	<input type="checkbox"/> 17	<input type="checkbox"/> 17	
		<input type="checkbox"/> 17.5	<input type="checkbox"/> 17.5	<input type="checkbox"/> 17.5	<input type="checkbox"/> 17.5	
	<input type="checkbox"/> 18	<input type="checkbox"/> 18	<input type="checkbox"/> 18	<input type="checkbox"/> 18	<input type="checkbox"/> 18	

Rear Wheel Size	Caster Size					
	4"	5"	6"	6" x 2"	7"	8"
	Seat Height					
20"	<input type="checkbox"/> 16	<input type="checkbox"/> 16	<input type="checkbox"/> 16	<input type="checkbox"/> 16		
	<input type="checkbox"/> 16.5	<input type="checkbox"/> 16.5	<input type="checkbox"/> 16.5	<input type="checkbox"/> 16.5		
		<input type="checkbox"/> 17	<input type="checkbox"/> 17	<input type="checkbox"/> 17	<input type="checkbox"/> 17	
		<input type="checkbox"/> 17.5	<input type="checkbox"/> 17.5	<input type="checkbox"/> 17.5	<input type="checkbox"/> 17.5	
			<input type="checkbox"/> 19 ¹			
		<input type="checkbox"/> 19.5	<input type="checkbox"/> 19.5	<input type="checkbox"/> 19.5	<input type="checkbox"/> 19.5	<input type="checkbox"/> 19.5
		<input type="checkbox"/> 20	<input type="checkbox"/> 20	<input type="checkbox"/> 20	<input type="checkbox"/> 20	<input type="checkbox"/> 20
22"	<input type="checkbox"/> 20.5	<input type="checkbox"/> 20.5	<input type="checkbox"/> 20.5	<input type="checkbox"/> 20.5	<input type="checkbox"/> 20.5	
		<input type="checkbox"/> 16.5		<input type="checkbox"/> 16.5		
		<input type="checkbox"/> 17	<input type="checkbox"/> 17	<input type="checkbox"/> 17	<input type="checkbox"/> 17	
		<input type="checkbox"/> 17.5	<input type="checkbox"/> 17.5	<input type="checkbox"/> 17.5	<input type="checkbox"/> 17.5	
		<input type="checkbox"/> 18	<input type="checkbox"/> 18	<input type="checkbox"/> 18	<input type="checkbox"/> 18	<input type="checkbox"/> 18
		<input type="checkbox"/> 18.5	<input type="checkbox"/> 18.5	<input type="checkbox"/> 18.5	<input type="checkbox"/> 18.5	<input type="checkbox"/> 18.5
		<input type="checkbox"/> 20.5	<input type="checkbox"/> 20.5	<input type="checkbox"/> 20.5	<input type="checkbox"/> 20.5	<input type="checkbox"/> 20.5
24"				<input type="checkbox"/> 17 ¹	<input type="checkbox"/> 17 ¹	
		<input type="checkbox"/> 17.5	<input type="checkbox"/> 17.5	<input type="checkbox"/> 17.5	<input type="checkbox"/> 17.5	
		<input type="checkbox"/> 18	<input type="checkbox"/> 18	<input type="checkbox"/> 18	<input type="checkbox"/> 18	<input type="checkbox"/> 18
		<input type="checkbox"/> 18.5	<input type="checkbox"/> 18.5	<input type="checkbox"/> 18.5	<input type="checkbox"/> 18.5	<input type="checkbox"/> 18.5
	<input type="checkbox"/> 19	<input type="checkbox"/> 19	<input type="checkbox"/> 19	<input type="checkbox"/> 19	<input type="checkbox"/> 19	

¹ Available only with Full Poly tire
² Not available with Full Poly tire

Casters and Rear Wheels

Casters

		Caster Size					
		Price	4"	5"	6"	7"	8"
Caster Type	1" Poly	Std	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>
	1.5" Poly	\$85	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>		
	2" Poly	\$90			<input type="checkbox"/>		

Axle

- Fixed _____ N/C
- Quick Release _____ N/C

Rear Wheel

		Rear Wheel Size				
		Price	12"	16"	20"	22"
Mag	Std	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>

Wheel Locks

Wheel Lock

- Push to Lock _____ Std
- Pull to Lock _____ N/C
Not available with 12" or 16" Full Poly tires
- Attendant Foot Lock _____ \$80
- Omit Wheel Locks _____ N/C

Armrests

Armrest

- Height Adjustable T-Arm Tall (11" - 15") E0973 **New!**
- Desk _____ \$310
 - Full _____ \$310
- Height Adjustable T-Arm (8.5"-13.5") E0973
- Desk _____ \$310
 - Full _____ \$310
- Height Adjustable T-Arm Low (6.5"-10.5") E0973
If requested configuration is unavailable, Standard Height T-Arm will be used
- Desk _____ \$310
 - Full _____ \$310
- Angle Adjustable Locking Extendable Flip Up E0973
Minimum armrest height with Reclining Backrest is 12".
- Desk _____ \$280
 - Full _____ \$280
- Tubular Flip Up _____ Std
Not available with Arm Pads
 - Omit Armrests _____ N/C

Rear Tire Type

- Full Poly E2218 _____ Std
- Low Poly E2218 _____ N/C
Available only with 16" wheel size
- Pneumatic E2211, E2212 _____ N/C
- Pneumatic with Airless Insert E2211, E2213 _____ \$90

Handrim

Handrims not available on 12" or 16" wheel sizes

Standard Handrims

- Aluminum Anodized _____ Std
- Plastic Coated _____ \$115
- Projection (8 Vertical Projections) _____ \$325
Available only with 24" wheels
- Omit Handrims _____ N/C

Wheel Lock Options

- Extension Handles E0961 _____ \$75
Available only with Push to Lock or Pull to Lock Wheel Locks.

Armrest Pad

- Standard Arm Pad _____ N/C
- Waterfall Arm Pad _____ \$25
Not available with Angle Adjustable Locking Extendable Flip Up Armrests
- Foam Grip _____ N/C
Available only with Angle Adjustable Locking Extendable Flip Up Armrest

Accessories

Removable Side Guard		O₂ Holder	
<input type="checkbox"/> Composite Side Guards K0108	\$140	<input type="checkbox"/> O ₂ Holder E2208	\$170
<i>Not available with T-Arm Armrests</i>		<i>Not available with Reclining Backrest</i>	
Anti-Tip		Other Accessories	
<input type="checkbox"/> Rear Anti-Tips E0971	\$155	<input type="checkbox"/> Tool Kit	\$55
<i>Required with Reclining Backrest</i>		<input type="checkbox"/> Removable Underseat Pouch	\$50
<input type="checkbox"/> Omit Rear Anti-Tips	N/C	<input type="checkbox"/> Ki Mobility Backpack	\$65
Pelvic Positioning Belt		<input type="checkbox"/> Black Neoprene Impact Guards	\$90
<input type="checkbox"/> 1 1/2" Auto Buckle E0978	\$80	<input type="checkbox"/> Tilt Stop (Pair) K0108	\$30
<input type="checkbox"/> 1 1/2" Padded Auto Buckle E0978	\$170	<input type="checkbox"/> IV Pole	\$275
<input type="checkbox"/> 2" Airline Buckle E0978	\$110	<input type="checkbox"/> Canopy with Hardware K0108	\$280
Bodypoint Evoflex - Pelvic Stabilizer E0978	\$188	<i>Not available on frame widths greater than 18". See image below.</i>	
Size: <input type="checkbox"/> Small (16-24") <input type="checkbox"/> Medium (19-29") <input type="checkbox"/> Large (25-37")			
Belt Mounting Kit			
<input type="checkbox"/> Bodypoint 1" Band Clamp Belt Mounting Kit	\$108		
Military Patch			
US Military Service Patch	\$90		
<input type="checkbox"/> Army <input type="checkbox"/> Navy <input type="checkbox"/> Marines <input type="checkbox"/> Air Force <input type="checkbox"/> Coast Guard			
<input type="checkbox"/> US Purple Heart Patch	\$90		
Purple Heart with Military Service Patch Combo	\$145		
<input type="checkbox"/> Army <input type="checkbox"/> Navy <input type="checkbox"/> Marines <input type="checkbox"/> Air Force <input type="checkbox"/> Coast Guard			

Frame Color

Standard Colors			
Gloss Finish		Matte Finish	
<input type="checkbox"/> Black	N/C	<input type="checkbox"/> Charcoal Gray	N/C
<input type="checkbox"/> Candy Blue	N/C	<input type="checkbox"/> Iridescent Burgundy	N/C
<input type="checkbox"/> Candy Purple	N/C	<input type="checkbox"/> Iridescent Navy	N/C
<input type="checkbox"/> Candy Red	N/C	<input type="checkbox"/> Military Green	N/C
<input type="checkbox"/> Copper	N/C	<input type="checkbox"/> Shadow Black	N/C
<input type="checkbox"/> Silver	N/C		
Optional Colors			
<i>Extended lead time</i>			
Gloss Finish		Matte Finish	
<input type="checkbox"/> Blue Skies	\$70	<input type="checkbox"/> Emerald Green	\$70
<input type="checkbox"/> Bubblegum Pink	\$70	<input type="checkbox"/> Retro Red	\$70
<input type="checkbox"/> Granny Smith	\$70	<input type="checkbox"/> Royal Blue	\$70
<input type="checkbox"/> Metallic Orange	\$70	<input type="checkbox"/> Smooth White	\$70
<input type="checkbox"/> Teal	\$70		
<input type="checkbox"/> Yellow	\$70		