



CATALYST 4

Phone (800) 981-1540 | Fax (715) 254-0996
 sales@kimobility.com | www.kimobility.com

Company: _____

Ship To: _____

Account Number: _____ Date: _____

Address: _____

Quote PO Number: _____

City: _____ State: _____

RTS: _____

Zip: _____ Phone: _____

Marked For: _____

Email/Fax: _____

Select if chair quote has: Additional Ki Products with Different Quote/Order Number **Quote/Order Number:** _____

Effective 07/18/2022

- Catalyst 4** K0004 **\$1195**
 15" Lowest Front Seat Height. 250 lb weight capacity.
- Catalyst 4 Super Low** K0004 **\$1285**
 13.5" Lowest Front Seat Height. 250 lb weight capacity.

Frame

- ### Seat Width
- 14" _____ \$100
 - 16" _____ N/C
 - 18" _____ N/C
 - 20" E2201 _____ \$455
 - 22" E2201 _____ \$500

- ### Weight Capacity
- 250 lb _____ Std
 - 350 lb (Heavy Duty) _____ \$500
Requires minimum 6" caster height

- ### Seat Depth
- 16" _____ N/C
 - 18" _____ N/C
 - 20" E2203 _____ \$495

Backrest

Back Post / Height Range / Back Height					
Back Post	Price	Back Height			
<input type="checkbox"/> Standard	N/C	<input type="checkbox"/> 16.5	<input type="checkbox"/> 17.5	<input type="checkbox"/> 18.5	<input type="checkbox"/> 19.5
<input type="checkbox"/> Angle Adjustable w/ Push Handle	\$185	<input type="checkbox"/> 16	<input type="checkbox"/> 17	<input type="checkbox"/> 18	

- ### Military Patch
- Attached only if ordered with Tension Adjustable Back Upholstery**
- US Military Service Patch _____ \$90
 Army Navy Marines Air Force Coast Guard
 - US Purple Heart Patch _____ \$90
 - Purple Heart with Military Service Patch Combo _____ \$145
 Army Navy Marines Air Force Coast Guard

- ### Back Upholstery
- If selecting both upholstery and Axiom Back, the upholstery will be installed.**
- Standard _____ Std
 - Tension Adjustable E2611
 - Standard Cover _____ \$275
 - Air Mesh Breathable Cover _____ \$350
 - Omit Back Upholstery _____ - \$10

HCPCS codes provided are based on code verification by the PDAC or our interpretation of the code definitions.
 Proper use of HCPCS codes does not ensure coverage or payment. For coverage information, verify the appropriate payer's coverage policy.

Footrest

Hanger				
Frame	Hanger	60°	70°	Retail
Swing Away	Extension Tube Mount	<input type="checkbox"/>	<input type="checkbox"/>	Std
	Pro Elevating Leg Rests (PRO ELR) E0990			<input type="checkbox"/> \$350
	Omit Hangers			<input type="checkbox"/> N/C

Footplate	
<input type="checkbox"/> Composite	N/C
<input type="checkbox"/> Composite Angle Adjustable K0040	\$190 <i>Not available with 14" Seat Width</i>
<input type="checkbox"/> Aluminum Angle Adjustable K0040	\$190

PRO ELR Covers

Available only on PRO ELR

Gel Pad Covers (Pair) K0108 \$390

Footrest / Extension / Seat to Footrest Length

Select Extension and Seat to Footrest Length in row corresponding to desired Footrest. Footrest Length of at least 2.5" shorter than Front Seat Height recommended. Extension Tube Hanger measurements are based on 70° Hanger.

Footrest	Extension	Seat to Footrest Length										
Extension Tube Hangers (60°, 70°)	<input type="checkbox"/> Short	<input type="checkbox"/> 11	<input type="checkbox"/> 11.5	<input type="checkbox"/> 12	<input type="checkbox"/> 12.5	<input type="checkbox"/> 13	<input type="checkbox"/> 13.5	<input type="checkbox"/> 14	<input type="checkbox"/> 14.5	<input type="checkbox"/> 15	<input type="checkbox"/> 16	<input type="checkbox"/> 17
	<input type="checkbox"/> Medium	<input type="checkbox"/> 13	<input type="checkbox"/> 13.5	<input type="checkbox"/> 14	<input type="checkbox"/> 14.5	<input type="checkbox"/> 15	<input type="checkbox"/> 15.5	<input type="checkbox"/> 16	<input type="checkbox"/> 16.5	<input type="checkbox"/> 17	<input type="checkbox"/> 18	<input type="checkbox"/> 19
	<input type="checkbox"/> Long	<input type="checkbox"/> 15.5	<input type="checkbox"/> 16	<input type="checkbox"/> 16.5	<input type="checkbox"/> 17	<input type="checkbox"/> 17.5	<input type="checkbox"/> 18	<input type="checkbox"/> 18.5	<input type="checkbox"/> 19	<input type="checkbox"/> 19.5		
PRO ELR Low Pivot Setting	<input type="checkbox"/> Short	12 - 17.5										
	<input type="checkbox"/> Long	15 - 21										
PRO ELR High Pivot Setting	<input type="checkbox"/> Short	10.5 - 15.5										
	<input type="checkbox"/> Long	13 - 19										

Heel Loop

- Heel Loops E0951 \$40
- Adjustable Heel Loops E0951 \$50
- Omit Heel Loops N/C

Calf Strap

- Velcro Adjustable K0038 \$40
- Padded Velcro Adjustable K0038 \$68
- Bodypoint Padded AeroMesh K0038 \$118

Residual Limb Support

Residual Limb Support E1020			
<input type="checkbox"/> Right and Left			\$790
<input type="checkbox"/> Left	\$395	<input type="checkbox"/> Right	\$395
<input type="checkbox"/> Omit Left Hanger	N/C	<input type="checkbox"/> Omit Right Hanger	N/C
Cover			
<input type="checkbox"/> Smooth	<input type="checkbox"/> Air Mesh	<input type="checkbox"/> Smooth	<input type="checkbox"/> Air Mesh
Pad Size			
<input type="checkbox"/> 5 x 8	<input type="checkbox"/> 9 x 14	<input type="checkbox"/> 5 x 8	<input type="checkbox"/> 9 x 14
<input type="checkbox"/> 7 x 10	<input type="checkbox"/> 11 x 10	<input type="checkbox"/> 7 x 10	<input type="checkbox"/> 11 x 10
<input type="checkbox"/> 9 x 10	<input type="checkbox"/> 11 x 14	<input type="checkbox"/> 9 x 10	<input type="checkbox"/> 11 x 14

Front Seat Height and Casters

Caster Size / Front Frame Type / Front Seat Height

Select Front Seat Height in row corresponding to desired caster size.
Note: Selections in first three Front Seat Height columns require Catalyst 4 Super Low.

Caster Size	Front Seat Height													
	Uses Super Low Front Frame													
4"	<input type="checkbox"/> 13.5 ¹	<input type="checkbox"/> 14	<input type="checkbox"/> 14.5	<input type="checkbox"/> 15 ²	<input type="checkbox"/> 15.5	<input type="checkbox"/> 16	<input type="checkbox"/> 16.5	<input type="checkbox"/> 17	<input type="checkbox"/> 17.5	<input type="checkbox"/> 18	<input type="checkbox"/> 18.5	<input type="checkbox"/> 19	<input type="checkbox"/> 19.5	<input type="checkbox"/> 20
5"	<input type="checkbox"/> 14.5	<input type="checkbox"/> 15	<input type="checkbox"/> 15.5	<input type="checkbox"/> 16	<input type="checkbox"/> 16.5	<input type="checkbox"/> 17	<input type="checkbox"/> 17.5	<input type="checkbox"/> 18	<input type="checkbox"/> 18.5	<input type="checkbox"/> 19	<input type="checkbox"/> 19.5	<input type="checkbox"/> 20	<input type="checkbox"/> 20.5	
6"	<input type="checkbox"/> 15.5	<input type="checkbox"/> 16	<input type="checkbox"/> 16.5	<input type="checkbox"/> 17	<input type="checkbox"/> 17.5	<input type="checkbox"/> 18	<input type="checkbox"/> 18.5	<input type="checkbox"/> 19	<input type="checkbox"/> 19.5	<input type="checkbox"/> 20	<input type="checkbox"/> 20.5	<input type="checkbox"/> 21		
7"	<input type="checkbox"/> 17	<input type="checkbox"/> 17.5		<input type="checkbox"/> 18	<input type="checkbox"/> 18.5	<input type="checkbox"/> 19	<input type="checkbox"/> 19.5	<input type="checkbox"/> 20	<input type="checkbox"/> 20.5	<input type="checkbox"/> 21				
8"	<input type="checkbox"/> 17.5	<input type="checkbox"/> 18	<input type="checkbox"/> 18.5	<input type="checkbox"/> 19	<input type="checkbox"/> 19.5	<input type="checkbox"/> 20	<input type="checkbox"/> 20.5	<input type="checkbox"/> 21						

Selection for Catalyst 4 Super Low is part of base frame selection on page 1

Front Seat Height and Casters Continued

Caster Type		Caster Size					
		Price	4"	5"	6"	7"	8"
Caster Type	1" Poly	N/C	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>
	1.5" Poly	\$55	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>		
	2" Poly	\$70			<input type="checkbox"/>		

Maximum Seat Dump Reference <small>(Front Seat Height minus Rear Seat Height)</small>			
Seat Depth	16"	18"	20"
Maximum Dump	2"	2.5"	2.5"

Fork Stem	
<input type="checkbox"/> Standard Fork Stem	Std
<input type="checkbox"/> 3/4" Longer Fork Stem	\$35
<input type="checkbox"/> 1 1/2" Longer Fork Stem	\$35

Rear Seat Height and Wheels

Wheel Size / Tire Profile / Rear Seat Height	
<i>Select Rear Seat Height in row corresponding to desired Wheel Size and Tire Profile</i>	

Wheel Size	Tire Profile	Rear Seat Height				
20"	Low	<input type="checkbox"/> 12.5 ¹	<input type="checkbox"/> 13.5	<input type="checkbox"/> 14.5	<input type="checkbox"/> 15.5	<input type="checkbox"/> 16.5
	Full	<input type="checkbox"/> 13 ¹	<input type="checkbox"/> 14	<input type="checkbox"/> 15	<input type="checkbox"/> 16	<input type="checkbox"/> 17
22"	Low	<input type="checkbox"/> 13.5 ¹	<input type="checkbox"/> 14.5	<input type="checkbox"/> 15.5	<input type="checkbox"/> 16.5	<input type="checkbox"/> 17.5
	Full	<input type="checkbox"/> 14 ¹	<input type="checkbox"/> 15	<input type="checkbox"/> 16	<input type="checkbox"/> 17	<input type="checkbox"/> 18
24"	Low	<input type="checkbox"/> 14.5 ¹	<input type="checkbox"/> 15.5	<input type="checkbox"/> 16.5	<input type="checkbox"/> 17.5	<input type="checkbox"/> 18.5
	Full	<input type="checkbox"/> 15 ¹	<input type="checkbox"/> 16	<input type="checkbox"/> 17	<input type="checkbox"/> 18	<input type="checkbox"/> 19

¹ Not available with Flip Back Armrests

Rear Wheel		Rear Wheel Size			
		Price	20"	22"	24"
Rear Wheel Type	MAXX Performance Spoke	\$75	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>
	Mag	N/C	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>

Rear Tire Type	
Low Profile	
<input type="checkbox"/> Low Poly 1" E2218	Std <small>Not available with 350 lb Weight Capacity</small>

Full Profile	
<input type="checkbox"/> Full Poly 1 3/8" E2218	\$40
<input type="checkbox"/> Street Pneumatic 1 3/8" E2211, E2212	\$80
<input type="checkbox"/> Pneumatic with Airless Insert 1 3/8" E2211, E2213	\$205

Handrim	
Standard Handrims	
<input type="checkbox"/> Aluminum Anodized	Std
<input type="checkbox"/> Plastic Coated	\$95
<input type="checkbox"/> Projection (8 Vertical Projections)	\$325 <small>Available only with full profile tires in combination with 22"-24" MAXX Performance or 24" Mag wheel</small>
<input type="checkbox"/> Omit Handrim	N/C

Natural Fit Handrims	
<i>Not available on 20" wheel size or 22" Mag wheel. Not recommended with Mag wheel.</i>	
<input type="checkbox"/> Natural Fit LT Regular Grip	\$465
<input type="checkbox"/> Natural Fit Standard Regular Grip	\$445
<input type="checkbox"/> Natural Fit LT Super Grip	\$515
<input type="checkbox"/> Natural Fit Standard Super Grip	\$495
<input type="checkbox"/> Natural Fit LT Omit Grip	\$395
<input type="checkbox"/> Natural Fit Standard Omit Grip	\$345

Handrim Option	
<input type="checkbox"/> 3M Micro Grip	\$65 <small>Not available on Plastic Coated or Projection Handrims</small>

Axles

Axle Bracket	
<input type="checkbox"/> Offset Axle Bracket	Std
<input type="checkbox"/> Reverse Mount Offset Axle Bracket	N/C <small>Not available with Swing Away Armrests</small>

Axle	
<input type="checkbox"/> Fixed Axle	Std
<input type="checkbox"/> Quick Release Axle	\$95

Wheel Locks

Wheel Lock	
<input type="checkbox"/> Push to Lock	Std
<input type="checkbox"/> Push to Lock Flush Mount	N/C
<input type="checkbox"/> Pull to Lock	N/C
<input type="checkbox"/> Omit Wheel Locks	N/C

Wheel Lock Options	
<input type="checkbox"/> Extension Handles E0961	\$75 <small>Available only with Push to Lock or Pull to Lock Wheel Locks</small>
<input type="checkbox"/> Grade Aids E0974	\$150 <small>Available only with Push to Lock Wheel Locks and Pneumatic or Pneumatic w/ Airless Insert tires</small>

Armrests

Armrest

- Height Adjustable T-Arm Tall (11.5" - 15") E0973 **New!**
- Desk _____ \$310
 - Full _____ \$310
- Height Adjustable T-Arm (8.5" - 13.5") E0973
- Desk _____ \$310
 - Full _____ \$310
- Height Adjustable T-Arm Low (6.5" - 10.5") E0973
If requested configuration is unavailable, Standard Height T-Arm will be used
- Desk _____ \$310
 - Full _____ \$310
- Height Adjustable Flip Back T-Arm Tall (11.5"-15") E0973 **New!**
- Desk _____ \$325
 - Full _____ \$325
- Height Adjustable Flip Back T-Arm (9.5"-13.5") E0973
- Desk _____ \$325
 - Full _____ \$325

Flip Back (10.5")

Not available with Angle Adjustable Back Posts

- Desk _____ \$35
- Full _____ \$35
- Swing Away (10.5" - 12.5") E0973 _____ \$225
Not available with Angle Adjustable Back Posts or Arm Pads
- Omit Armrests _____ N/C

Arm Pad

Not available with Swing Away Armrests

- Waterfall Arm Pad _____ \$25
Available only with Height Adjustable T-Arm and Height Adjustable T-Arm Low

Accessories

Removable Side Guard

- Composite Side Guards K0108 _____ \$145
Not available with T-Arm or Flip Back Armrests

Anti-Tip

- Rear Anti-Tips E0971 _____ \$150
- Omit Rear Anti-Tips _____ N/C

Pelvic Positioning Belt

- 1 1/2" Auto Buckle E0978 _____ \$80
- 1 1/2" Padded Auto Buckle E0978 _____ \$170
- 2" Airline Buckle E0978 _____ \$110
- Bodypoint Evoflex - Pelvic Stabilizer E0978 _____ \$188
 Size: Small (16-24") Medium (19-29") Large (25-37")
- Bodypoint Belt Mounting Kit _____ \$108
 Style: Bracket Band Clamp 1"

O₂ Holder

- O₂ Holder E2208 _____ \$295
- Right Side Left Side

Other Accessories

- Tool Kit _____ \$45
- Removable Underseat Pouch _____ \$50
- Ki Mobility Backpack _____ \$65
- Black Neoprene Impact Guards _____ \$90
- Clear Spoke Guards K0065 _____ \$115
Not available with Mag wheel
- Black Spoke Guards K0065 _____ \$115
Not available with 20" wheel size or Mag wheel
- Cane and Crutch Holder _____ \$195
- IV Pole _____ \$265
Straight Back Canes will be used when this option is selected

Frame Color

Matte Finish

- Charcoal Gray _____ N/C
- Emerald Green _____ N/C
- Iridescent Burgundy _____ N/C
- Iridescent Navy _____ N/C
- Military Green _____ N/C
- Retro Red _____ N/C
- Royal Blue _____ N/C
- Shadow Black _____ N/C
- Smooth White _____ N/C

Gloss Finish

- Black _____ N/C

Catalyst Front Seat Height Matrix																								
4" Caster				5" Caster				6" Caster				7" Caster				8" Caster								
FSH	Super Low	Stem	Fork	FSH	Super Low	Stem	Fork	FSH	Super Low	Stem	Fork	FSH	Super Low	Stem	Fork	FSH	Super Low	Stem	Fork					
15	13.5	Std	4"	16	14.5	Std	4"	17	15.5	Std	5"	18	16.5	Std	6"	19	17.5	Std	6"					
15.5	14			16.5	15.5	3/4		17.5	16.5	3/4		18.5	17			19.5	18.5	3/4						
15.5 ¹	14 ¹			17.5	16	1 1/2		18.5	17	1 1/2		18.5	17.5			20.5	n/a	1 1/2						
16	14.5	3/4		16	14.5	Std	5"	17	15.5	Std	6"	19	18	3/4		19	17.5	Std	7"					
16.5 ¹	15 ¹	1 1/2		16.5	15			17.5	16			19.5	18			19.5	18							
17	15.5	1 1/2		16.5	15.5			18	16.5			20	18.5			20	18.5							
15.5	14	Std		5"	17	16		3/4	5"	17.5		16.5	3/4	6"		18.5	17	1 1/2	20.5	19	Std	7"		
16	14.5				17.5	16				1 1/2		18				17	19		17.5	19.5			18.5	
16	15				18	16.5				1 1/2		18.5				17.5	19.5		18	20			19	
16.5	15.5	3/4			16.5	15		Std		6"		18.5	17			1 1/2	7"	20	18.5	3/4	20.5	n/a	3/4	7"
17	15.5	1 1/2	17		15.5	19						17.5	19		18			21	n/a					
17.5	16	17.5	16		19.5	18						19.5	18.5		21			n/a	1 1/2					
16.5	15	Std	6"		17.5	16		3/4				6"	18		16.5	Std		7"	20	19	3/4	7"	7"	7"
17	15.5				18	16.5	18.5				17		19		17.5				20.5	n/a				
17	16				18.5	17	20				18.5		20		18.5									
17.5	16.5	3/4			18	16.5	1 1/2	7"			19.5		18		1 1/2	7"			20.5	19	1 1/2			
18	16.5	1 1/2		18.5	17	18.5			17.5		19		18	21					n/a					
18.5	17	19		17.5	19	18			19		18													
18	16.5	Std		7"	18	16.5	Std		7"		19.5		18.5	3/4	7"				7"	7"				
18.5	17				18.5	17				20	19		20				19							
18.5	17.5				19	17.5				19.5	18		19.5				18							
19	18	3/4			18.5	17.5	3/4			7"	20		18.5	1 1/2			7"				7"			
19.5	18.5	1 1/2	19		18	20.5					19	20.5	19											
20	19	1 1/2	19.5		18.5	21					n/a	21	n/a											
					19.5	18	1 1/2				6" x 2" caster not available with 5" fork													
					20	18.5																		
					20.5	19																		

¹ Not available with 4 x 1.5" Soft Roll Aluminum

Seat heights with "FSH" heading available with Standard Front Frame. Seat heights with "Super Low" heading available with Super Low Front Frame. Actual seat height may vary depending upon configuration.