



Phone (800) 981-1540 | Fax (715) 254-0996  
 sales@kimobility.com | www.kimobility.com

Company: \_\_\_\_\_

Account Number: \_\_\_\_\_ Date: \_\_\_\_\_

Quote  PO Number: \_\_\_\_\_

RTS: \_\_\_\_\_

Marked For: \_\_\_\_\_

Select if chair quote has:  Additional Ki Products with Different Quote/Order Number **Quote/Order Number:** \_\_\_\_\_

**Catalyst 4C** K0004 **\$1595**  
 250 lb weight capacity

Effective 07/18/2022

# CATALYST 4C

Ship To: \_\_\_\_\_

Address: \_\_\_\_\_

City: \_\_\_\_\_ State: \_\_\_\_\_

Zip: \_\_\_\_\_ Phone: \_\_\_\_\_

Email/Fax: \_\_\_\_\_

## Frame

### Seat Width

- 14" \$45
- 15" \$45
- 16" N/C
- 17" \$50
- 18" N/C
- 19" \$155
- 20" E2201 \$455
- 21" E2201 \$500
- 22" E2201 \$500

### Seat Depth

- 14" N/C
- 15" \$70
- 16" N/C
- 17" \$90
- 18" N/C
- 19" \$195
- 20" E2203 \$505

### Weight Capacity

- 250 lb Std
- 350 lb (Heavy Duty) \$500  
*Requires minimum 6" casters*

### Cross Brace

- Standard Cross Braces Std
- Seating System Cross Braces N/C  
*Allows for easy attachment of seating system mounting hardware. Not available with Seat Sling.*

### Seat Pan

- Removable Aluminum Seat Pan E2231 \$275  
*Pan depth 1/2" shorter than Seat Depth*

## Backrest

### Back Post / Height Range / Back Height

Back Post	Price	Height Range	Back Height			
<input type="checkbox"/> Standard	N/C	<input type="checkbox"/> Short	<input type="checkbox"/> 8.5	<input type="checkbox"/> 9.5	<input type="checkbox"/> 10.5	<input type="checkbox"/> 11.5
		<input type="checkbox"/> Medium	<input type="checkbox"/> 12.5	<input type="checkbox"/> 13.5	<input type="checkbox"/> 14.5	<input type="checkbox"/> 15.5
		<input type="checkbox"/> Tall	<input type="checkbox"/> 16.5	<input type="checkbox"/> 17.5	<input type="checkbox"/> 18.5	<input type="checkbox"/> 19.5
<input type="checkbox"/> Angle Adjustable w/ Push Handle	\$195	<input type="checkbox"/> Short	<input type="checkbox"/> 14	<input type="checkbox"/> 15	<input type="checkbox"/> 16	
		<input type="checkbox"/> Medium	<input type="checkbox"/> 16	<input type="checkbox"/> 17	<input type="checkbox"/> 18	
		<input type="checkbox"/> Tall	<input type="checkbox"/> 18	<input type="checkbox"/> 19	<input type="checkbox"/> 20	

### Standard Back Post

Select option below only if a Standard Back Post option is selected on chart at left

- Straight with Push Handle Std
- Straight without Push Handle N/C
- 8° Bend with Push Handle \$25  
*Not available with Short Back Height Range*

HCPCS codes provided are based on code verification by the PDAC or our interpretation of the code definitions.

Proper use of HCPCS codes does not ensure coverage or payment. For coverage information, verify the appropriate payer's coverage policy.

**Backrest continued**

**Back Upholstery**

*If selecting both upholstery and Axiom Back, the upholstery will be installed.*

- Standard \_\_\_\_\_ Std  
Not available with Short or Medium Back Height Ranges on 21" and 22" Seat Widths
- Tension Adjustable E2611
  - Standard Cover \_\_\_\_\_ \$280
  - Air Mesh Breathable Cover \_\_\_\_\_ \$355
- Omit Back Upholstery \_\_\_\_\_ - \$10

**Footrest**

**Frame Type**

- Swing Away \_\_\_\_\_ N/C
- Super Low - Swing Away \_\_\_\_\_ N/C

Minimum Front Seat Height with Front Mount Hangers						
Hanger Length	Short		Medium		Long	
Hanger Bend	70°	80/90°	70°	80/90°	70°	80°
Minimum Front Seat Height with Standard Selection	13.5"	13.5"	14"	15"	16"	17"
Minimum Front Seat Height with 3" Extension	14.5"	16"	16.5"	18"	18.5"	20"

**Military Patch**

*Attached only if ordered with Tension Adjustable Back Upholstery*

- US Military Service Patch \_\_\_\_\_ \$90
  - Army  Navy  Marines  Air Force  Coast Guard
- US Purple Heart Patch \_\_\_\_\_ \$90
- Purple Heart with Military Service Patch Combo \_\_\_\_\_ \$145
  - Army  Navy  Marines  Air Force  Coast Guard

**Hanger**

Hanger	60°	70°	80°	90°		Retail
Extension Tube Mount	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/> <sup>1</sup>			Std
Front Mount		<input type="checkbox"/>	<input type="checkbox"/> <sup>1</sup>	<input type="checkbox"/> <sup>1</sup>		N/C
Pro Elevating Leg Rests (PRO ELR) E0990					<input type="checkbox"/>	\$350
Omit Hangers					<input type="checkbox"/>	-\$20

<sup>1</sup>Not available with 7" or 8" casters

**PRO ELR Covers**

- Available only on PRO ELR*
- Gel Pad Covers (Pair) K0108 \_\_\_\_\_ \$390

**Footrest / Extension / Seat to Footrest Length**

*Select Extension and Seat to Footrest Length in row corresponding to desired Footrest. Footrest Length of at least 2.5" shorter than Front Seat Height recommended. Extension Tube Hanger measurements are based on 70° Hanger.*

Footrest	Extension	Seat to Footrest Length											
Extension Tube Hangers (60°, 70°, 80°)	<input type="checkbox"/> Short	<input type="checkbox"/> 11	<input type="checkbox"/> 11.5	<input type="checkbox"/> 12	<input type="checkbox"/> 12.5	<input type="checkbox"/> 13	<input type="checkbox"/> 13.5	<input type="checkbox"/> 14	<input type="checkbox"/> 14.5	<input type="checkbox"/> 15			
	<input type="checkbox"/> Medium	<input type="checkbox"/> 13	<input type="checkbox"/> 13.5	<input type="checkbox"/> 14	<input type="checkbox"/> 14.5	<input type="checkbox"/> 15	<input type="checkbox"/> 15.5	<input type="checkbox"/> 16	<input type="checkbox"/> 16.5	<input type="checkbox"/> 17			
	<input type="checkbox"/> Long	<input type="checkbox"/> 15.5	<input type="checkbox"/> 16	<input type="checkbox"/> 16.5	<input type="checkbox"/> 17	<input type="checkbox"/> 17.5	<input type="checkbox"/> 18	<input type="checkbox"/> 18.5	<input type="checkbox"/> 19	<input type="checkbox"/> 19.5			
Front Mount Hangers <sup>1</sup>	<input type="checkbox"/> Short (70°, 80°, 90°)	Standard Selections									Requires 3" Extension for \$35		
		<input type="checkbox"/> 4 <sup>2</sup>	<input type="checkbox"/> 4.5 <sup>2</sup>	<input type="checkbox"/> 5	<input type="checkbox"/> 5.5	<input type="checkbox"/> 6	<input type="checkbox"/> 6.5		<input type="checkbox"/> 10	<input type="checkbox"/> 10.5	<input type="checkbox"/> 11		
	<input type="checkbox"/> Medium (70°, 80°, 90°)	<input type="checkbox"/> 7	<input type="checkbox"/> 7.5	<input type="checkbox"/> 8	<input type="checkbox"/> 8.5	<input type="checkbox"/> 9	<input type="checkbox"/> 9.5		<input type="checkbox"/> 11.5	<input type="checkbox"/> 12	<input type="checkbox"/> 12.5		
		<input type="checkbox"/> 4 <sup>2</sup>	<input type="checkbox"/> 4.5 <sup>2</sup>	<input type="checkbox"/> 5	<input type="checkbox"/> 5.5	<input type="checkbox"/> 6	<input type="checkbox"/> 6.5	<input type="checkbox"/> 7	<input type="checkbox"/> 7.5		<input type="checkbox"/> 12	<input type="checkbox"/> 12.5	<input type="checkbox"/> 13
	<input type="checkbox"/> Long (70°, 80°)	<input type="checkbox"/> 8	<input type="checkbox"/> 8.5	<input type="checkbox"/> 9	<input type="checkbox"/> 9.5	<input type="checkbox"/> 10	<input type="checkbox"/> 10.5	<input type="checkbox"/> 11	<input type="checkbox"/> 11.5		<input type="checkbox"/> 13.5	<input type="checkbox"/> 14	<input type="checkbox"/> 14.5
		<input type="checkbox"/> 5	<input type="checkbox"/> 5.5	<input type="checkbox"/> 6	<input type="checkbox"/> 6.5	<input type="checkbox"/> 7	<input type="checkbox"/> 7.5	<input type="checkbox"/> 8	<input type="checkbox"/> 8.5	<input type="checkbox"/> 9	<input type="checkbox"/> 14	<input type="checkbox"/> 14.5	<input type="checkbox"/> 15
PRO ELR Low Pivot Setting	<input type="checkbox"/> Short	12 - 17.5											
	<input type="checkbox"/> Long	15 - 21											
PRO ELR High Pivot Setting	<input type="checkbox"/> Short	10.5 - 15.5											
	<input type="checkbox"/> Long	13 - 19											

<sup>1</sup> See Minimum Front Seat Height chart above.

<sup>2</sup> Available with 90° only

**3" Extension**

*Available only with Front Mount Hangers*

- 3" Extension \_\_\_\_\_ \$35
- 3" Extension (Provided as Parts) \_\_\_\_\_ \$35

**Footplate**


- Composite \_\_\_\_\_ N/C  
Not available with Front Mount Hanger
- Composite Angle Adjustable K0040 \_\_\_\_\_ \$195  
Not available with Front Mount Hanger or 14" Seat Width
- Aluminum Angle Adjustable K0040 \_\_\_\_\_ \$195
- One Piece Flip Up Angle Adjustable K0040 \_\_\_\_\_ \$190  
Not available with Front Mount Hanger or PRO ELR

**Footrest continued**

**Hanger Release**

*If no selection is made Classic Hanger Release will be used*

- Classic \_\_\_\_\_ Std
- 4-Way - In line \_\_\_\_\_ N/C  
*See reference A below. Not available with PRO ELR.*
- 4-Way - Out \_\_\_\_\_ N/C  
*See reference B below. Not available with PRO ELR.*



**Calf Strap**

- Velcro Adjustable K0038 \_\_\_\_\_ \$40
- Padded Velcro Adjustable K0038 \_\_\_\_\_ \$68
- Bodypoint Padded AeroMesh K0038 \_\_\_\_\_ \$118

**Residual Limb Support**

- Residual Limb Support E1020
- Right and Left \_\_\_\_\_ \$790
  - Left \_\_\_\_\_ \$395
  - Omit Left Hanger \_ -\$10
  - Right \_\_\_\_\_ \$395
  - Omit Right Hanger \_ -\$10

**Cover**

- Smooth  Air Mesh
- Smooth  Air Mesh

**Pad Size**

- |                                 |                                  |                                 |                                  |
|---------------------------------|----------------------------------|---------------------------------|----------------------------------|
| <input type="checkbox"/> 5 x 8  | <input type="checkbox"/> 9 x 14  | <input type="checkbox"/> 5 x 8  | <input type="checkbox"/> 9 x 14  |
| <input type="checkbox"/> 7 x 10 | <input type="checkbox"/> 11 x 10 | <input type="checkbox"/> 7 x 10 | <input type="checkbox"/> 11 x 10 |
| <input type="checkbox"/> 9 x 10 | <input type="checkbox"/> 11 x 14 | <input type="checkbox"/> 9 x 10 | <input type="checkbox"/> 11 x 14 |

**Heel Loop**

*Not available with One Piece Flip Up Angle Adjustable Footplate*

- Heel Loops E0951 \_\_\_\_\_ \$40
- Adjustable Heel Loops E0951 \_\_\_\_\_ \$50
- Omit Heel Loops \_\_\_\_\_ N/C

**Front Seat Height and Casters**

**Caster Size / Front Frame Type / Front Seat Height**

*Select Front Seat Height in row corresponding to desired caster size.*

*Note: Selections in first three Front Seat Height columns use Super Low Front Frame.*

Caster Size	Front Seat Height																
	Uses Super Low Front Frame			13.5	14	14.5	15	15.5	16	16.5	17	17.5	18	18.5	19	19.5	20
4"	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>
5"	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>
6"	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>
7"	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>
8"	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>

**Caster Type**

Caster Type	Price	Caster Size				
		4"	5"	6"	7"	8"
1" Poly	N/C	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>
1" Pneumatic	\$90			<input type="checkbox"/>		<input type="checkbox"/>
1.5" Poly	\$55	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>		
1.5" Soft Roll Aluminum	\$185	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>		
2" Poly	\$70			<input type="checkbox"/>		
2" Pneumatic w/ Foam Insert	\$225					<input type="checkbox"/>

**Maximum Seat Dump Reference**

(Front Seat Height minus Rear Seat Height)

Seat Depth	14"	15"	16"	17"	18"	19"	20"
Maximum Dump	1.5"	2"	2"	2"	2.5"	2.5"	2.5"

**Fork Options**

- Frog Legs Suspension Fork \_\_\_\_\_ \$450  
*Available only with 4", 5" or 6" Caster Sizes. Not available with Caster Pin Locks, 350 lb Weight Capacity, 2" Poly caster or 2" Pneumatic w/ Foam Insert caster.*
  - User Weight Under 80 lbs
  - User Weight 80 - 120 lbs
  - User Weight 120 - 160 lbs
  - User Weight 160 - 220 lbs
  - User Weight 220 - 260 lbs
  - User Weight 260 - 300 lbs
- Caster Pin Locks K0073 \_\_\_\_\_ \$195  
*Not available with 1 1/2" Longer Fork Stem or Frog Legs Fork*

**Fork Stem**

- Standard Fork Stem \_\_\_\_\_ Std
- 3/4" Longer Fork Stem \_\_\_\_\_ \$35
- 1 1/2" Longer Fork Stem \_\_\_\_\_ \$35

**Wheel Locks**

- Wheel Lock**
- Push to Lock \_\_\_\_\_ Std
  - Push to Lock Flush Mount \_\_\_\_\_ N/C
  - Pull to Lock \_\_\_\_\_ N/C
  - Omit Wheel Locks \_\_\_\_\_ N/C

**Wheel Lock Options**

- Extension Handles E0961 \_\_\_\_\_ \$75  
*Available only with Push to Lock or Pull to Lock Wheel Locks*
- Grade Aids E0974 \_\_\_\_\_ \$150  
*Available only with Push to Lock Wheel Locks and Pneumatic or Pneumatic with Airless Insert tires*

**Rear Seat Height and Wheels**

**Wheel Size / Tire Profile / Rear Seat Height**

Select Rear Seat Height in row corresponding to desired Wheel Size and Tire Profile. 14" Seat Depth not available with 24", 25" or 26" wheel sizes.

Wheel Size	Tire Profile	Rear Seat Height				
20"	Low	<input type="checkbox"/> 12.5 <sup>1</sup>	<input type="checkbox"/> 13.5	<input type="checkbox"/> 14.5	<input type="checkbox"/> 15.5	<input type="checkbox"/> 16.5
	Full	<input type="checkbox"/> 13 <sup>1</sup>	<input type="checkbox"/> 14	<input type="checkbox"/> 15	<input type="checkbox"/> 16	<input type="checkbox"/> 17
22"	Low	<input type="checkbox"/> 13.5 <sup>1</sup>	<input type="checkbox"/> 14.5	<input type="checkbox"/> 15.5	<input type="checkbox"/> 16.5	<input type="checkbox"/> 17.5
	Full	<input type="checkbox"/> 14 <sup>1</sup>	<input type="checkbox"/> 15	<input type="checkbox"/> 16	<input type="checkbox"/> 17	<input type="checkbox"/> 18
24"	Low	<input type="checkbox"/> 14.5 <sup>1</sup>	<input type="checkbox"/> 15.5	<input type="checkbox"/> 16.5	<input type="checkbox"/> 17.5	<input type="checkbox"/> 18.5
	Full	<input type="checkbox"/> 15 <sup>1</sup>	<input type="checkbox"/> 16	<input type="checkbox"/> 17	<input type="checkbox"/> 18	<input type="checkbox"/> 19
25"	Low	<input type="checkbox"/> 15 <sup>1</sup>	<input type="checkbox"/> 16	<input type="checkbox"/> 17	<input type="checkbox"/> 18	<input type="checkbox"/> 19
26"	Low	<input type="checkbox"/> 15.5 <sup>1</sup>	<input type="checkbox"/> 16.5	<input type="checkbox"/> 17.5	<input type="checkbox"/> 18.5	<input type="checkbox"/> 19.5
	Full	<input type="checkbox"/> 16 <sup>1</sup>	<input type="checkbox"/> 17	<input type="checkbox"/> 18	<input type="checkbox"/> 19	<input type="checkbox"/> 20

<sup>1</sup> Not available with Flip Back Armrests

**Rear Wheel**

		Rear Wheel Size					
		Price	20"	22"	24"	25"	26"
Rear Wheel	MAXX Performance Spoke	\$125	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>
	MAXX Performance Spoke	N/C			<input type="checkbox"/>		
	Mag	N/C	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>		
	Superlight Spoke <sup>1</sup>	\$295			<input type="checkbox"/>	<input type="checkbox"/> <sup>2</sup>	
	Spinergy Spox	\$845		<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>

<sup>1</sup> Not available with 350 lb Weight Capacity

<sup>2</sup> Available only with Superlight and Plastic Coated Handrims

**Spinergy Options**

Select spoke color and handrim tab for Spinergy wheels

Spoke Color:

- Black
- Red
- White
- Blue
- Yellow

**Axles**

**Axle Bracket**

- Offset Axle Bracket \_\_\_\_\_ Std
- Reverse Mount Offset Axle Bracket \_\_\_\_\_ N/C  
*Not available with Swing Away Armrests*

**Axle**

- Quick Release \_\_\_\_\_ Std
- Quad Quick Release \_\_\_\_\_ \$195  
*Not available with Mag or Spinergy wheels*

**Rear Tire Type**

**Full Profile**

- Street Pneumatic 1 3/8" E2211, E2212 \_\_\_\_\_ Std  
*Not available on 25" wheel size*
- Full Poly 1 3/8" E2218 \_\_\_\_\_ N/C  
*Not available on Superlight wheel or 25" or 26" wheel sizes*
- Pneumatic with Airless Insert 1 3/8" E2211, E2213 \_\_\_\_\_ \$100  
*Not available on Superlight wheel or 25" wheel size*

**Low Profile**

- Low Poly 1" E2218 \_\_\_\_\_ N/C  
*Not available on Superlight wheel, 25" or 26" wheel sizes or 350 lb Weight Cap.*
- Sentinel Solid Tire 1" E2218 \_\_\_\_\_ \$95  
*Not available on Superlight wheel or 20" and 22" Mag wheels*
- High Pressure 1" E2211, E2212 \_\_\_\_\_ \$95  
*Not available on Mag wheel*

**Handrim**

**Standard Handrims**

- Aluminum Anodized \_\_\_\_\_ Std
- Superlight Aluminum Anodized \_\_\_\_\_ \$50  
*Available only on 24" and 25" wheel sizes. Available only on MAXX Performance, Superlight or Spinergy wheels. Not available with 350 lb Weight Capacity.*
- Plastic Coated \_\_\_\_\_ \$95
- Projection (8 Vertical Projections) \_\_\_\_\_ \$325  
*Available only with full profile tires in combination with 22", 24" or 26" MAXX Performance, or 24" Mag wheel*
- Omit Handrim \_\_\_\_\_ N/C

**Natural Fit Handrims**

*Not available on 20" wheel size or 22" Mag wheel. Not recommended with Mag wheel.*

- Natural Fit LT Regular Grip \_\_\_\_\_ \$465
- Natural Fit Standard Regular Grip \_\_\_\_\_ \$445
- Natural Fit LT Super Grip \_\_\_\_\_ \$515
- Natural Fit Standard Super Grip \_\_\_\_\_ \$495
- Natural Fit LT Omit Grip \_\_\_\_\_ \$395
- Natural Fit Standard Omit Grip \_\_\_\_\_ \$345

**Handrim Options**

- 3M Micro Grip \_\_\_\_\_ \$65  
*Not available on Plastic Coated or Projection Handrims*

Handrim Tab:

*Available only with Spinergy, Mag or Maxx Performance Spoke and Aluminum Anodized or Plastic Coated Handrims*

- Short Tab \_\_\_\_\_ N/C
- Long Tab \_\_\_\_\_ N/C

## Armrests

Armrest	
Height Adjustable T-Arm Tall (11.5" - 15") E0973 <b>New!</b>	
<input type="checkbox"/> Desk	\$310
<input type="checkbox"/> Full	\$310
Height Adjustable T-Arm (8.5" - 13.5") E0973	
<input type="checkbox"/> Desk	\$310
<input type="checkbox"/> Full	\$310
Height Adjustable T-Arm Low (6.5" - 10.5") E0973	
<i>If requested configuration is unavailable, Standard Height T-Arm will be used</i>	
<input type="checkbox"/> Desk	\$310
<input type="checkbox"/> Full	\$310
Height Adjustable Flip Back T-Arm Tall (11.5"-15") E0973 <b>New!</b>	
<input type="checkbox"/> Desk	\$325
<input type="checkbox"/> Full	\$325
Height Adjustable Flip Back T-Arm (9.5"-13.5") E0973	
<input type="checkbox"/> Desk	\$325
<input type="checkbox"/> Full	\$325
Flip Back (10.5")	
<i>Not available with Seat Depths less than 16" or Angle Adjustable Back Posts</i>	
<input type="checkbox"/> Desk	N/C
<input type="checkbox"/> Full	N/C
<input type="checkbox"/> Swing Away (10.5" - 12.5") E0973	\$225
<i>Not available with Angle Adjustable Back Posts or Arm Pads</i>	
<input type="checkbox"/> Omit Armrests	N/C

### Arm Pad

**Not available with Swing Away Armrests**

<input type="checkbox"/> Waterfall Arm Pad	\$25
<i>Available only with Height Adjustable T-Arm and Height Adjustable T-Arm Low</i>	

## Accessories

Removable Side Guard	
<input type="checkbox"/> Composite Side Guards K0108	\$145
<i>Not available with T-Arm or Flip Back Armrests</i>	
Anti-Tip	
<input type="checkbox"/> Rear Anti-Tips E0971	\$150
<input type="checkbox"/> Omit Rear Anti-Tips	N/C
Pelvic Positioning Belt	
<input type="checkbox"/> 1 1/2" Auto Buckle E0978	\$80
<input type="checkbox"/> 1 1/2" Padded Auto Buckle E0978	\$170
<input type="checkbox"/> 2" Airline Buckle E0978	\$110
Bodypoint Evoflex - Pelvic Stabilizer E0978	\$188
Size: <input type="checkbox"/> Small (16-24") <input type="checkbox"/> Medium (19-29") <input type="checkbox"/> Large (25-37")	
Bodypoint Belt Mounting Kit	\$108
Style: <input type="checkbox"/> Bracket <input type="checkbox"/> Band Clamp 1"	
O <sub>2</sub> Holder	
O <sub>2</sub> Holder E2208	\$295
<input type="checkbox"/> Right Side <input type="checkbox"/> Left Side	
Other Accessories	
<input type="checkbox"/> Tool Kit	\$45
<input type="checkbox"/> Removable Underseat Pouch	\$50
<input type="checkbox"/> Ki Mobility Backpack	\$65
<input type="checkbox"/> Black Neoprene Impact Guards	\$90
<input type="checkbox"/> Clear Spoke Guards K0065	\$115
<i>Not available with Mag wheel</i>	
<input type="checkbox"/> Black Spoke Guards K0065	\$115
<i>Not available with 20" wheel size or Mag wheel</i>	
<input type="checkbox"/> Cane and Crutch Holder	\$195
<input type="checkbox"/> IV Pole	\$275
<i>Available only with Standard Straight Back Posts</i>	

## Frame Color

Gloss Finish	
<input type="checkbox"/> Black	N/C
<input type="checkbox"/> Blue Skies	N/C
<input type="checkbox"/> Bubblegum Pink	N/C
<input type="checkbox"/> Candy Blue	N/C
<input type="checkbox"/> Candy Purple	N/C
<input type="checkbox"/> Candy Red	N/C
<input type="checkbox"/> Copper	N/C
<input type="checkbox"/> Granny Smith	N/C
<input type="checkbox"/> Metallic Orange	N/C
<input type="checkbox"/> Yellow	N/C
Matte Finish	
<input type="checkbox"/> Charcoal Gray	N/C
<input type="checkbox"/> Emerald Green	N/C
<input type="checkbox"/> Iridescent Burgundy	N/C
<input type="checkbox"/> Iridescent Navy	N/C
<input type="checkbox"/> Military Green	N/C
<input type="checkbox"/> Retro Red	N/C
<input type="checkbox"/> Royal Blue	N/C
<input type="checkbox"/> Shadow Black	N/C
<input type="checkbox"/> Smooth White	N/C
Hydrographics	
Hydrographics	\$775
<input type="checkbox"/> AOR1 (Desert Digital Camo)	
<input type="checkbox"/> Pink Camo	
<input type="checkbox"/> Muddy Girl	
<input type="checkbox"/> NEXT Camo	
<input type="checkbox"/> Carbon Fiber <sup>1</sup>	
<input type="checkbox"/> Zebra <sup>1</sup>	
<sup>1</sup> Select Hydrographics Base Color. Zebra not available in Black.	
Hydrographics Base Color	N/C
<input type="checkbox"/> Black <input type="checkbox"/> Blue <input type="checkbox"/> Granny Smith <input type="checkbox"/> Pink <input type="checkbox"/> Purple <input type="checkbox"/> Red	

Catalyst Front Seat Height Matrix																								
4" Caster				5" Caster				6" Caster				7" Caster				8" Caster								
FSH	Super Low	Stem	Fork	FSH	Super Low	Stem	Fork	FSH	Super Low	Stem	Fork	FSH	Super Low	Stem	Fork	FSH	Super Low	Stem	Fork					
15	13.5	Std	4"	16	14.5	Std	4"	17	15.5	Std	5"	18	16.5	Std	6"	19	17.5	Std	6"					
15.5	14			16.5	15.5	3/4		17.5	16.5	3/4		18.5	17			19.5	18.5	3/4						
15.5 <sup>1</sup>	14 <sup>1</sup>			17.5	16	1 1/2		18.5	17	1 1/2		18.5	17.5			20.5	n/a	1 1/2						
16	14.5	3/4		16	14.5	Std	5"	17	15.5	Std	6"	19	18	3/4		19	17.5	Std	7"					
16.5 <sup>1</sup>	15 <sup>1</sup>	1 1/2		16.5	15			17.5	16			19.5	18			19.5	18							
17	15.5	16.5		15.5	3/4			18	16.5			20	18.5			20	18.5							
15.5	14	Std		5"	17	16		3/4	5"	17.5		16.5	3/4	6"		18.5	17	1 1/2	20.5	19	Std	7"		
16	14.5				17.5	16				1 1/2		18				17	19		17.5	19.5			18.5	
16	15				18	16.5				1 1/2		18.5				17.5	19.5		18	20			19	
16.5	15.5	3/4			16.5	15		Std		6"		18.5	17			1 1/2	7"	20	18.5	3/4	20.5	n/a	3/4	7"
17	15.5	1 1/2	17		15.5	19						17.5	19		18			21	n/a					
17.5	16	17.5	16		19.5	18						19.5	18.5		21			n/a	1 1/2					
16.5	15	Std	6"		17.5	16		3/4				6"	18		16.5	Std		7"	20	19	3/4	7"	7"	7"
17	15.5				18	16.5	18.5				17		19		17.5				20.5	n/a				
17	16				18.5	17	20				18.5		20		18.5									
17.5	16.5	3/4			18	16.5	1 1/2	7"			19.5		18		1 1/2	7"			20.5	19	1 1/2			
18	16.5	1 1/2		18.5	17	18.5			17.5		19		18	21					n/a					
18.5	17	19		17.5	19	18			19		18													
18	16.5	Std		7"	18	16.5	Std		7"		19.5		18.5	3/4	7"				7"	7"				
18.5	17				18.5	17				20	19		20				19							
18.5	17.5				19	17.5				19.5	18		19.5				18							
19	18	3/4			18.5	17.5	3/4			7"	20		18.5	1 1/2			7"				7"			
19.5	18.5	1 1/2	19		18	20.5					19	20.5	19											
20	19	19.5	18.5		21	n/a					21	n/a												
					19.5	18	1 1/2				6" x 2" caster not available with 5" fork													
					20	18.5																		
					20.5	19																		

<sup>1</sup> Not available with 4 x 1.5" Soft Roll Aluminum

Seat heights with "FSH" heading available with Standard Front Frame. Seat heights with "Super Low" heading available with Super Low Front Frame. Actual seat height may vary depending upon configuration.