



CATALYST 4TTL

Phone (800) 981-1540 | Fax (715) 254-0996
 sales@kimobility.com | www.kimobility.com

Company: _____ Ship To: _____
 Account Number: _____ Date: _____ Address: _____
 Quote PO Number: _____ City: _____ State: _____
 RTS: _____ Zip: _____ Phone: _____
 Marked For: _____ Email/Fax: _____
 Select if chair quote has: Additional Ki Products with Different Quote/Order Number **Quote/Order Number:** _____

Catalyst 4TTL K0004 **\$1870** Effective 01/17/2022
 250 lb weight capacity

Frame

Seat Width

14" _____ N/C
 15" _____ N/C
 16" _____ N/C
 17" _____ N/C
 18" _____ N/C
 19" _____ N/C
 20" E2201 \$455
 21" E2201 \$500
 22" E2201 \$500

Weight Capacity

250 lb _____ Std
 350 lb (Heavy Duty) _____ \$500
Requires minimum 6" caster height

Cross Brace

Standard Cross Braces _____ Std
 Seating System Cross Braces _____ N/C
Allows for easy attachment of seating system mounting hardware. Not available with Seat Sling.

Seat Depth

14" _____ N/C
 15" _____ N/C
 16" _____ N/C
 17" _____ N/C
 18" _____ N/C
 19" _____ N/C
 20" E2203 \$495

Seat Pan

Removable Aluminum Seat Pan E2231 _____ \$175
Pan depth 1/2" shorter than Seat Depth

Backrest

Back Post / Height Range / Back Height						
Back Post	Price	Height Range	Back Height			
<input type="checkbox"/> Standard	N/C	<input type="checkbox"/> Short	<input type="checkbox"/> 8.5	<input type="checkbox"/> 9.5	<input type="checkbox"/> 10.5	<input type="checkbox"/> 11.5
		<input type="checkbox"/> Medium	<input type="checkbox"/> 12.5	<input type="checkbox"/> 13.5	<input type="checkbox"/> 14.5	<input type="checkbox"/> 15.5
		<input type="checkbox"/> Tall	<input type="checkbox"/> 16.5	<input type="checkbox"/> 17.5	<input type="checkbox"/> 18.5	<input type="checkbox"/> 19.5
<input type="checkbox"/> Angle Adjustable w/ Push Handle	\$185	<input type="checkbox"/> Short	<input type="checkbox"/> 14	<input type="checkbox"/> 15	<input type="checkbox"/> 16	
		<input type="checkbox"/> Medium	<input type="checkbox"/> 16	<input type="checkbox"/> 17	<input type="checkbox"/> 18	
		<input type="checkbox"/> Tall	<input type="checkbox"/> 18	<input type="checkbox"/> 19	<input type="checkbox"/> 20	

Standard Back Post

Select option below only if a Standard Back Post option is selected on chart above

Straight with Push Handle _____ Std
 Straight without Push Handle _____ N/C
 8° Bend with Push Handle _____ N/C
Not available with Short Back Height Range

HCPCS codes provided are based on code verification by the PDAC or our interpretation of the code definitions.

Proper use of HCPCS codes does not ensure coverage or payment. For coverage information, verify the appropriate payer's coverage policy.

Backrest continued

Back Upholstery

If selecting both upholstery and Axiom Back, the upholstery will be installed.

- Standard _____ Std
Not available with Short or Medium Back Height Ranges on 21" and 22" Seat Widths
- Tension Adjustable E2611
 - Standard Cover _____ \$275
 - Air Mesh Breathable Cover _____ \$350
- Omit Back Upholstery _____ - \$10

Military Patch

Attached only if ordered with Tension Adjustable Back Upholstery

- US Military Service Patch _____ \$90
 - Army Navy Marines Air Force Coast Guard
- US Purple Heart Patch _____ \$90
- Purple Heart with Military Service Patch Combo _____ \$145
 - Army Navy Marines Air Force Coast Guard

Footrest

Frame Type

- Swing Away _____ N/C
- Super Low - Swing Away _____ N/C

Minimum Front Seat Height with Front Mount Hangers						
Hanger Length	Short		Medium		Long	
	70°	80/90°	70°	80/90°	70°	80°
Hanger Bend	70°	80/90°	70°	80/90°	70°	80°
Minimum Front Seat Height with Standard Selection	13.5"	13.5"	14"	15"	16"	17"
Minimum Front Seat Height with 3" Extension	14.5"	16"	16.5"	18"	18.5"	20"

Hanger

Hanger	60°	70°	80°	90°		Retail
Extension Tube Mount	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/> ¹			Std
Front Mount		<input type="checkbox"/>	<input type="checkbox"/> ¹	<input type="checkbox"/> ¹		N/C
Pro Elevating Leg Rests (PRO ELR) E0990					<input type="checkbox"/>	\$295
Omit Hangers					<input type="checkbox"/>	N/C

¹Not available with 7" or 8" casters

PRO ELR Cover

Available only on PRO ELR

- Gel Pad Cover _____ \$390

Footrest / Extension / Seat to Footrest Length

Select Extension and Seat to Footrest Length in row corresponding to desired Footrest. Footrest Length of at least 2.5" shorter than Front Seat Height recommended. Extension Tube Hangers measurements are based on 70° Hanger.

Footrest	Extension	Seat to Footrest Length										
Extension Tube Hangers (60°, 70°, 80°)	<input type="checkbox"/> Short	<input type="checkbox"/> 11	<input type="checkbox"/> 11.5	<input type="checkbox"/> 12	<input type="checkbox"/> 12.5	<input type="checkbox"/> 13	<input type="checkbox"/> 13.5	<input type="checkbox"/> 14	<input type="checkbox"/> 14.5	<input type="checkbox"/> 15		
	<input type="checkbox"/> Medium	<input type="checkbox"/> 13	<input type="checkbox"/> 13.5	<input type="checkbox"/> 14	<input type="checkbox"/> 14.5	<input type="checkbox"/> 15	<input type="checkbox"/> 15.5	<input type="checkbox"/> 16	<input type="checkbox"/> 16.5	<input type="checkbox"/> 17		
	<input type="checkbox"/> Long	<input type="checkbox"/> 15.5	<input type="checkbox"/> 16	<input type="checkbox"/> 16.5	<input type="checkbox"/> 17	<input type="checkbox"/> 17.5	<input type="checkbox"/> 18	<input type="checkbox"/> 18.5	<input type="checkbox"/> 19	<input type="checkbox"/> 19.5		
Front Mount Hangers ¹	<input type="checkbox"/> Short (70°, 80°, 90°)	Standard Selections								Requires 3" Extension		
		<input type="checkbox"/> 5	<input type="checkbox"/> 5.5	<input type="checkbox"/> 6	<input type="checkbox"/> 6.5	<input type="checkbox"/> 7				<input type="checkbox"/> 10	<input type="checkbox"/> 10.5	<input type="checkbox"/> 11
	<input type="checkbox"/> Medium (70°, 80°, 90°)	<input type="checkbox"/> 7.5	<input type="checkbox"/> 8	<input type="checkbox"/> 8.5	<input type="checkbox"/> 9	<input type="checkbox"/> 9.5				<input type="checkbox"/> 11.5	<input type="checkbox"/> 12	<input type="checkbox"/> 12.5
		<input type="checkbox"/> 5	<input type="checkbox"/> 5.5	<input type="checkbox"/> 6	<input type="checkbox"/> 6.5	<input type="checkbox"/> 7	<input type="checkbox"/> 7.5	<input type="checkbox"/> 8		<input type="checkbox"/> 12	<input type="checkbox"/> 12.5	<input type="checkbox"/> 13
	<input type="checkbox"/> Long (70°, 80°)	<input type="checkbox"/> 8.5	<input type="checkbox"/> 9	<input type="checkbox"/> 9.5	<input type="checkbox"/> 10	<input type="checkbox"/> 10.5	<input type="checkbox"/> 11	<input type="checkbox"/> 11.5		<input type="checkbox"/> 13.5	<input type="checkbox"/> 14	<input type="checkbox"/> 14.5
		<input type="checkbox"/> 5	<input type="checkbox"/> 5.5	<input type="checkbox"/> 6	<input type="checkbox"/> 6.5	<input type="checkbox"/> 7	<input type="checkbox"/> 7.5	<input type="checkbox"/> 8	<input type="checkbox"/> 8.5	<input type="checkbox"/> 9	<input type="checkbox"/> 14	<input type="checkbox"/> 14.5
PRO ELR Low Pivot Setting	<input type="checkbox"/> Short	12 - 17.5										
	<input type="checkbox"/> Long	15 - 21										
PRO ELR High Pivot Setting	<input type="checkbox"/> Short	10.5 - 15.5										
	<input type="checkbox"/> Long	13 - 19										

¹ See Minimum Front Seat Height chart above

3" Extension

Available only with Front Mount Hangers

- 3" Extension _____ N/C
- 3" Extension (Provided as Parts) _____ N/C


Footplate

- Composite _____ N/C
Not available with Front Mount Hanger or PRO ELR
- Composite Angle Adjustable K0040 _____ \$190
Not available with Front Mount Hanger or 14" Seat Width
- Aluminum Angle Adjustable K0040 _____ \$190
- One Piece Flip Up Angle Adjustable K0040 _____ \$190
Not available with Front Mount Hanger or PRO ELR

Footrest continued

Hanger Release

- If no selection is made Classic Hanger Release will be used*
- Classic _____ Std
 - 4-Way - In line _____ N/C
See reference A below. Not available with PRO ELR.
 - 4-Way - Out _____ N/C
See reference B below. Not available with PRO ELR.



Heel Loop

- Not available with One Piece Flip Up Angle Adjustable Footplate*
- Heel Loops E0951 _____ \$40
 - Adjustable Heel Loops E0951 _____ \$50
 - Omit Heel Loops _____ N/C

Front Seat Height and Casters

Caster Size / Front Frame Type / Front Seat Height

Select Front Seat Height in row corresponding to desired caster size.
 Note: Selections in first three Front Seat Height columns use Super Low Front Frame.

Caster Size	Front Seat Height																
	Uses Super Low Front Frame			13.5	14	14.5	15	15.5	16	16.5	17	17.5	18	18.5	19	19.5	20
4"	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>
5"	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>
6"	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>
7"	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>
8"	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>

Caster

Caster Type		Price	Caster Size				
			4"	5"	6"	7"	8"
Caster Type	1" Poly	N/C	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>
	1" Pneumatic	\$90			<input type="checkbox"/>		<input type="checkbox"/>
	1.5" Poly	N/C	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>		
	1.5" Soft Roll Aluminum	\$65	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>		
	2" Poly	N/C			<input type="checkbox"/>		
	2" Pneumatic w/ Foam Insert	\$125				<input type="checkbox"/>	

Fork Stem

- Best Fit _____ Std
Ki Mobility will pick Fork and Stem combination to optimize performance based on Caster Size and Front Seat Height
- Standard Fork Stem _____ N/C
- 3/4" Longer Fork Stem _____ N/C
- 1 1/2" Longer Fork Stem _____ N/C

Calf Strap

- Velcro Adjustable K0038 _____ \$35
- Padded Velcro Adjustable K0038 _____ \$40
- Bodypoint Padded AeroMesh K0038 _____ \$118

Residual Limb Support

- Residual Limb Support E1020
- Right and Left _____ \$790
 - Left _____ \$395
 - Right _____ \$395
 - Omit Left Hanger _____ N/C
 - Omit Right Hanger _____ N/C
- Cover**
- Smooth Air Mesh
- Pad Size**
- | | | | |
|---------------------------------|----------------------------------|---------------------------------|----------------------------------|
| <input type="checkbox"/> 5 x 8 | <input type="checkbox"/> 9 x 14 | <input type="checkbox"/> 5 x 8 | <input type="checkbox"/> 9 x 14 |
| <input type="checkbox"/> 7 x 10 | <input type="checkbox"/> 11 x 10 | <input type="checkbox"/> 7 x 10 | <input type="checkbox"/> 11 x 10 |
| <input type="checkbox"/> 9 x 10 | <input type="checkbox"/> 11 x 14 | <input type="checkbox"/> 9 x 10 | <input type="checkbox"/> 11 x 14 |

Maximum Seat Dump Reference

(Front Seat Height minus Rear Seat Height)

Seat Depth	14"	15"	16"	17"	18"	19"	20"
Maximum Dump	1.5"	2"	2"	2"	2.5"	2.5"	2.5"

Fork Options

- Frog Legs Suspension Fork _____ \$450
Available only with 4", 5" or 6" Caster Sizes. Not available with Caster Pin Locks, 350 lb Weight Capacity, 2" Poly caster or 2" Pneumatic w/ Foam Insert caster.
 - User Weight Under 80 lbs
 - User Weight 80 - 120 lbs
 - User Weight 120 - 160 lbs
 - User Weight 160 - 220 lbs
 - User Weight 220 - 260 lbs
 - User Weight 260 - 300 lbs
- Caster Pin Locks K0073 _____ \$195
1 1/2" Longer Fork Stem or Frog Legs Fork

Wheel Locks

Wheel Lock	
<input type="checkbox"/> Push to Lock	Std
<input type="checkbox"/> Push to Lock Flush Mount	N/C
<input type="checkbox"/> Pull to Lock	N/C
<input type="checkbox"/> Omit Wheel Locks	N/C

Wheel Lock Options	
<input type="checkbox"/> Extension Handles E0961	\$75 <i>Available only with Push to Lock or Pull to Lock Wheel Locks</i>
<input type="checkbox"/> Grade Aids E0974	\$150 <i>Available only with Push to Lock Wheel Locks and Pneumatic or Pneumatic with Airless Insert tires</i>

Rear Seat Height and Wheels

Wheel Size / Tire Profile / Rear Seat Height
Select Rear Seat Height in row corresponding to desired Wheel Size and Tire Profile. 14" Seat Depth not available with 24", 25" or 26" wheel sizes.

Wheel Size	Tire Profile	Rear Seat Height				
20"	Low	<input type="checkbox"/> 12.5 ¹	<input type="checkbox"/> 13.5	<input type="checkbox"/> 14.5	<input type="checkbox"/> 15.5	<input type="checkbox"/> 16.5
	Full	<input type="checkbox"/> 13 ¹	<input type="checkbox"/> 14	<input type="checkbox"/> 15	<input type="checkbox"/> 16	<input type="checkbox"/> 17
22"	Low	<input type="checkbox"/> 13.5 ¹	<input type="checkbox"/> 14.5	<input type="checkbox"/> 15.5	<input type="checkbox"/> 16.5	<input type="checkbox"/> 17.5
	Full	<input type="checkbox"/> 14 ¹	<input type="checkbox"/> 15	<input type="checkbox"/> 16	<input type="checkbox"/> 17	<input type="checkbox"/> 18
24"	Low	<input type="checkbox"/> 14.5 ¹	<input type="checkbox"/> 15.5	<input type="checkbox"/> 16.5	<input type="checkbox"/> 17.5	<input type="checkbox"/> 18.5
	Full	<input type="checkbox"/> 15 ¹	<input type="checkbox"/> 16	<input type="checkbox"/> 17	<input type="checkbox"/> 18	<input type="checkbox"/> 19
25"	Low	<input type="checkbox"/> 15 ¹	<input type="checkbox"/> 16	<input type="checkbox"/> 17	<input type="checkbox"/> 18	<input type="checkbox"/> 19
26"	Low	<input type="checkbox"/> 15.5 ¹	<input type="checkbox"/> 16.5	<input type="checkbox"/> 17.5	<input type="checkbox"/> 18.5	<input type="checkbox"/> 19.5
	Full	<input type="checkbox"/> 16 ¹	<input type="checkbox"/> 17	<input type="checkbox"/> 18	<input type="checkbox"/> 19	<input type="checkbox"/> 20

¹ Not available with Flip Back Armrests

Rear Wheel		Rear Wheel Size					
		Price	20"	22"	24"	25"	26"
Rear Wheel Type	Assembled Six-Spoke Mag w/ Treaded Tire and Composite Handrim ^{1,2}	-\$25			<input type="checkbox"/>		
	MAXX Performance Spoke	\$125	<input type="checkbox"/>	<input type="checkbox"/>		<input type="checkbox"/>	<input type="checkbox"/>
	MAXX Performance Spoke	N/C			<input type="checkbox"/>		
	Mag	N/C	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>		
	Superlight Spoke ²	\$265			<input type="checkbox"/>	<input type="checkbox"/>	³
	Spinergy Spox	\$845		<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>

¹ No substitutions available. Do not select Rear Tire Type or Handrim. Measures as Full Profile tire.

² Not available with 350 lb Weight Capacity

³ Available only with Superlight and Plastic Coated Handrims

Spinergy Options	
Select spoke color and handrim tab for Spinergy wheels	
Spoke Color:	
<input type="checkbox"/> Black	<input type="checkbox"/> Red <input type="checkbox"/> White <input type="checkbox"/> Blue <input type="checkbox"/> Yellow
Handrim Tab:	
<i>Available only with Aluminum Anodized and Plastic Coated Handrims</i>	
<input type="checkbox"/> Short Tab	N/C
<input type="checkbox"/> Long Tab	N/C

Axles

Axle Bracket	
<input type="checkbox"/> Offset Axle Bracket	Std
<input type="checkbox"/> Reverse Mount Offset Axle Bracket	N/C <i>Not available with Swing Away Armrests</i>

Rear Tire Type	
Full Profile	
<input type="checkbox"/> Street Pneumatic 1 3/8" E2211, E2212	Std <i>Not available on 25" wheel size</i>
<input type="checkbox"/> Full Poly 1 3/8" E2218	N/C <i>Not available on Superlight wheel or 25" or 26" wheel sizes</i>
<input type="checkbox"/> Pneumatic with Airless Insert 1 3/8" E2211, E2213	\$100 <i>Not available on Superlight wheel or 25" wheel size</i>
Low Profile	
<input type="checkbox"/> Low Poly 1" E2218	N/C <i>Not available on Superlight wheel, 25" or 26" wheel sizes or v lb Weight Cap.</i>
<input type="checkbox"/> Sentinel Solid Tire 1" E2218	\$95 <i>Not available on Superlight wheel or 20" and 22" Mag wheels</i>
<input type="checkbox"/> High Pressure 1" E2211, E2212	\$95 <i>Not available on Mag wheel</i>

Handrim	
Standard Handrims	
<input type="checkbox"/> Aluminum Anodized	Std
<input type="checkbox"/> Superlight Aluminum Anodized	\$50 <i>Available only on 24" and 25" wheel sizes. Available only on MAXX Performance, Superlight or Spinergy wheels. Not available with 350 lb Weight Capacity.</i>
<input type="checkbox"/> Plastic Coated	\$50
<input type="checkbox"/> Projection (8 Vertical Projections)	\$325 <i>Available only with full profile tires in combination with 22"-26" MAXX Performance, or 24" Mag wheel</i>
<input type="checkbox"/> Omit Handrim	N/C

Natural Fit Handrims	
<i>Not available on 20" wheel size or 22" Mag wheel. Not recommended with Mag wheel.</i>	
<input type="checkbox"/> Natural Fit LT Regular Grip	\$465
<input type="checkbox"/> Natural Fit Standard Regular Grip	\$445
<input type="checkbox"/> Natural Fit LT Super Grip	\$515
<input type="checkbox"/> Natural Fit Standard Super Grip	\$495
<input type="checkbox"/> Natural Fit LT Omit Grip	\$395
<input type="checkbox"/> Natural Fit Standard Omit Grip	\$345

Handrim Option	
<input type="checkbox"/> 3M Micro Grip	\$65 <i>Not available on Plastic Coated or Projection Handrims</i>

Axle	
<input type="checkbox"/> Quick Release	Std
<input type="checkbox"/> Quad Quick Release	\$195 <i>Not available with Mag or Spinergy wheels</i>

Armrests

Armrest	
Height Adjustable T-Arm E0973	
<input type="checkbox"/> Desk (8.5" - 13.5")	\$295
<input type="checkbox"/> Full (8.5" - 13.5")	\$295
Height Adjustable T-Arm Low E0973	
<i>If selection is not available with rear wheel and rear seat height configuration, standard Height Adjustable T-Arm armrests will be used</i>	
<input type="checkbox"/> Desk (7.5" - 11.5")	\$295
<input type="checkbox"/> Full (7.5" - 11.5")	\$295
Height Adjustable Flip Back T-Arm E0973	
<input type="checkbox"/> Desk (9.5"-13.5")	\$305
<input type="checkbox"/> Full (9.5"-13.5")	\$305
Flip Back	
<i>Not available with Seat Depths less than 16" or Angle Adjustable Back Posts</i>	
<input type="checkbox"/> Desk (10.5")	N/C
<input type="checkbox"/> Full (10.5")	N/C
<input type="checkbox"/> Swing Away (10.5" - 12.5") E0973	\$225
<i>Not available with Angle Adjustable Back Posts or Arm Pads</i>	
<input type="checkbox"/> Omit Armrests	N/C
Arm Pad	
Not available with Swing Away Armrests	
<input type="checkbox"/> Waterfall Arm Pad	\$25
<i>Available only with Height Adjustable T-Arm and Height Adjustable T-Arm Low</i>	

Accessories

Removable Side Guard	
<input type="checkbox"/> Composite Side Guards K0108	\$145
<i>Not available with T-Arm or Flip Back Armrests</i>	
Anti-Tip	
<input type="checkbox"/> Rear Anti-Tips E0971	\$150
<input type="checkbox"/> Omit Rear Anti-Tips	N/C
Pelvic Positioning Belt	
<input type="checkbox"/> 1 1/2" Auto Buckle E0978	\$80
<input type="checkbox"/> 1 1/2" Padded Auto Buckle E0978	\$170
<input type="checkbox"/> 2" Airline Buckle E0978	\$110
Bodypoint Evoflex - Pelvic Stabilizer E0978	\$188
Size: <input type="checkbox"/> Small (16-24") <input type="checkbox"/> Medium (19-29") <input type="checkbox"/> Large (25-37")	
Bodypoint Belt Mounting Kit	\$108
Style: <input type="checkbox"/> Bracket <input type="checkbox"/> Band Clamp 1"	
O ₂ Holder	
O ₂ Holder E2208	\$295
<input type="checkbox"/> Right Side <input type="checkbox"/> Left Side	
Other Accessories	
<input type="checkbox"/> Tool Kit	\$45
<input type="checkbox"/> Removable Underseat Pouch	\$50
<input type="checkbox"/> Ki Mobility Backpack	\$65
<input type="checkbox"/> Black Neoprene Impact Guards	\$90
<input type="checkbox"/> Clear Spoke Guards K0065	\$115
<i>Not available with Mag wheel</i>	
<input type="checkbox"/> Black Spoke Guards K0065	\$115
<i>Not available with 20" wheel size or Mag wheel</i>	
<input type="checkbox"/> Cane and Crutch Holder	\$195
<input type="checkbox"/> IV Pole	\$275
<i>Available only with Standard Straight Back Posts</i>	

Frame Color

Gloss Finish	
<input type="checkbox"/> Black	N/C
<input type="checkbox"/> Blue Skies	N/C
<input type="checkbox"/> Bubblemum Pink	N/C
<input type="checkbox"/> Candy Blue	N/C
<input type="checkbox"/> Candy Purple	N/C
<input type="checkbox"/> Candy Red	N/C
<input type="checkbox"/> Copper	N/C
<input type="checkbox"/> Granny Smith	N/C
<input type="checkbox"/> Metallic Orange	N/C
<input type="checkbox"/> Yellow	N/C
Matte Finish	
<input type="checkbox"/> Charcoal Gray	N/C
<input type="checkbox"/> Emerald Green	N/C
<input type="checkbox"/> Iridescent Burgundy	N/C
<input type="checkbox"/> Iridescent Navy	N/C
<input type="checkbox"/> Military Green	N/C
<input type="checkbox"/> Retro Red	N/C
<input type="checkbox"/> Royal Blue	N/C
<input type="checkbox"/> Shadow Black	N/C
<input type="checkbox"/> Smooth White	N/C
Hydrographics	
Hydrographics	\$695
<input type="checkbox"/> AOR1 (Desert Digital Camo)	<input type="checkbox"/> Pink Camo
<input type="checkbox"/> Muddy Girl	<input type="checkbox"/> NEXT Camo
<input type="checkbox"/> Carbon Fiber ¹	<input type="checkbox"/> Zebra ¹
¹Select Hydrographics Base Color. Zebra not available in Black.	
Hydrographics Base Color	N/C
<input type="checkbox"/> Black <input type="checkbox"/> Blue <input type="checkbox"/> Granny Smith <input type="checkbox"/> Pink <input type="checkbox"/> Purple <input type="checkbox"/> Red	

Catalyst Front Seat Height Matrix																																					
4" Caster				5" Caster				6" Caster				7" Caster				8" Caster																					
FSH	Super Low	Stem	Fork	FSH	Super Low	Stem	Fork	FSH	Super Low	Stem	Fork	FSH	Super Low	Stem	Fork	FSH	Super Low	Stem	Fork																		
15	13.5	Std	4"	16	14.5	Std	4"	17	15.5	Std	5"	18	16.5	Std	6"	19	17.5	Std	6"																		
15.5	14			16.5	15.5	3/4		17.5	16.5	3/4		18.5	17			19.5	18.5	3/4																			
15.5 ¹	14 ¹			17.5	16	1 1/2		18.5	17	1 1/2		18.5	17.5			20.5	n/a	1 1/2																			
16	14.5	3/4		5"	16	14.5	Std	5"	17	15.5	Std	6"	19	18		3/4	7"	19	17.5	Std	7"																
16.5 ¹	15 ¹	1 1/2			16.5	15			17.5	16			18.5	17				19.5	18			19.5	18														
17	15.5	3/4			16.5	15.5			17.5	16.5			18	16.5				19.5	18.5			20	18.5														
15.5	14	Std			5"	17	16		3/4	5"	17.5		16.5	3/4		6"		18.5	17	1 1/2		7"	20.5	19	Std	7"											
16	14.5					1 1/2	17.5				16		18					17	18.5				17.5	19.5			18	19.5	18.5								
16	15					3/4	18				16.5		18.5					17.5	19.5				18	19.5			18	20	19								
16.5	15.5	1 1/2				5"	16.5		15		Std		5"	18.5				17	1 1/2	6"			20	18.5	Std		7"	20.5	n/a	3/4	7"						
17	15.5		17				15.5		19					17.5	19			18					20.5	n/a													
17.5	16		17.5				16		19.5					18	19.5			18.5					21	n/a													
16.5	15	Std	6"				17.5		16		3/4			6"	18			16.5	Std				7"	20	19			3/4	7"	20.5		n/a	1 1/2	7"			
17	15.5			18			16.5	18.5	17			19			17.5		20	18.5																			
17	16			18.5			17	19.5	18			20			18.5		20.5	19																			
17.5	16.5	3/4		6"			18	16.5	1 1/2		6"	18.5			17		3/4	7"	20.5		19			1 1/2	7"			21		n/a		1 1/2	7"				
18	16.5				18.5		17	19.5		18		20			18.5																						
18.5	17				19		17.5	19.5		18		20.5			19																						
18	16.5	Std			7"		18	16.5	Std	7"		19.5			18.5	3/4	7"		6" x 2" caster not available with 5" fork																		
18.5	17					18.5	17	20				19																									
18.5	17.5					19	17.5	19.5				18																									
19	18	3/4				7"	18.5	17.5	3/4			7"	20		18.5	1 1/2											7"	6" x 2" caster not available with 5" fork									
19.5	18.5		19				18	20.5					19																								
20	19		19.5				18.5	21					n/a																								
							19.5	18	1 1/2				7"	6" x 2" caster not available with 5" fork																							
				20			18.5																														
				20.5			19																														

¹ Not available with 4 x 1.5" Soft Roll Aluminum

Seat heights with "FSH" heading available with Standard Front Frame. Seat heights with "Super Low" heading available with Super Low Front Frame. Actual seat height may vary depending upon configuration.