



Phone (800) 981-1540 | Fax (715) 254-0996
 sales@kimobility.com | www.kimobility.com

CATALYST 5

Company: _____

Ship To: _____

Account Number: _____ Date: _____

Address: _____

Quote PO Number: _____

City: _____ State: _____

RTS: _____

Zip: _____ Phone: _____

Marked For: _____

Email/Fax: _____

Select if chair quote has: Additional Ki Products with Different Quote/Order Number **Quote/Order Number:** _____

Catalyst 5 K0005 **\$2145**
 300 lb weight capacity

Effective 01/17/2022

Frame

Seat Width

- 14" N/C
- 15" N/C
- 16" N/C
- 17" N/C
- 18" N/C
- 19" N/C
- 20" E2201 \$510
- 21" E2201 \$510
- 22" E2201 \$510

Seat Depth

- 14" N/C
- 15" N/C
- 16" N/C
- 17" N/C
- 18" N/C
- 19" \$75
- 20" E2203 \$560

Transit

- Transit System \$325
Not available with Removable Seat Pan or 350 lb Weight Capacity
- Omit Transit N/C

Wheelchair Drive Ready

Chair will ship from Ki Mobility with adaption kit installed for selected Alber system. Order system separately from Alber USA. Select Rear Seat Height (full profile selections only) in drive wheel size on Rear Seat Height chart (page 5). Requires Rear Anti-Tips. Not available with Spinergy wheels or Height Adjustable Flip Back T-Arm.

- e-motion Ready N/C
- twion Ready N/C

Titanium Option

- Titanium Frame K0108 \$1875
Includes hangers and base frame. Does not include Cross Braces.

Weight Capacity

- 300 lb Std
- 350 lb (Heavy Duty) \$500
Requires minimum 6" caster height

Cross Brace

- Standard Cross Braces Std
- Seating System Cross Braces N/C
Allows for easy attachment of seating system mounting hardware. Not available with Seat Sling.
- Titanium Cross Braces \$700
Not available with 21" and 22" Seat Widths

Seat Pan

- Removable Aluminum Seat Pan E2231 \$275
Not available with Transit System. Pan depth 1/2" shorter than Seat Depth.

Rear Frame Type

No selection needed with Titanium Frame

- Curved Std
Select from Standard or Angle Adjustable Back Posts on chart on the next page
- Depth Adjustable K0108 \$195
Select from depth adjustable back posts in chart on next page. Provides 2" forward and 2" rearward adjustment with minimum Seat Depth of 14" and maximum Seat Depth of 22". Available only with T-Arm or Flip Up Armrests.
- Reclining \$1095
Select from Reclining Back Post in chart on the next page. Not available with Seating Dynamics Rocker Back, O₂ Holder, Transit Option, 350 lb Weight Capacity, Swing Away Armrests, Flip Back Armrest, Height Adjustable Flip Back T-Arm, Back Upholstery or IV Pole.

HCPCS codes provided are based on code verification by the PDAC or our interpretation of the code definitions.

Proper use of HCPCS codes does not ensure coverage or payment. For coverage information, verify the appropriate payer's coverage policy.

Backrest

Back Post / Height Range / Back Height

Depth Adjustable not available with Titanium Frame.

Back Post	Price	Height Range	Back Height
<input type="checkbox"/> Standard	N/C	<input type="checkbox"/> Short ^{3,4}	<input type="checkbox"/> 8.5 <input type="checkbox"/> 9.5 <input type="checkbox"/> 10.5 <input type="checkbox"/> 11.5
		<input type="checkbox"/> Medium ⁴	<input type="checkbox"/> 12.5 ⁶ <input type="checkbox"/> 13.5 <input type="checkbox"/> 14.5 <input type="checkbox"/> 15.5
		<input type="checkbox"/> Tall	<input type="checkbox"/> 16.5 <input type="checkbox"/> 17.5 <input type="checkbox"/> 18.5 <input type="checkbox"/> 19.5
<input type="checkbox"/> Angle Adjustable with Push Handle	\$185	<input type="checkbox"/> Short ⁴	<input type="checkbox"/> 14 <input type="checkbox"/> 15 <input type="checkbox"/> 16
		<input type="checkbox"/> Medium ⁴	<input type="checkbox"/> 16 <input type="checkbox"/> 17 <input type="checkbox"/> 18
		<input type="checkbox"/> Tall	<input type="checkbox"/> 18 <input type="checkbox"/> 19 <input type="checkbox"/> 20
<input type="checkbox"/> Stroller Handle Back Post (Depth Adjustable)	N/C	<input type="checkbox"/> Short	<input type="checkbox"/> 20
		<input type="checkbox"/> Tall	<input type="checkbox"/> 24
<input type="checkbox"/> Height Adjustable ^{2,9} Straight with Push Handle (Depth Adjustable)	\$85	<input type="checkbox"/> Short ¹	<input type="checkbox"/> 13 <input type="checkbox"/> 14 <input type="checkbox"/> 15 <input type="checkbox"/> 16 <input type="checkbox"/> 17
		<input type="checkbox"/> Medium	<input type="checkbox"/> 16 <input type="checkbox"/> 17 <input type="checkbox"/> 18 <input type="checkbox"/> 19 <input type="checkbox"/> 20
		<input type="checkbox"/> Tall	<input type="checkbox"/> 20 ⁸ <input type="checkbox"/> 21 ⁵ <input type="checkbox"/> 22 ⁵ <input type="checkbox"/> 23 ⁵ <input type="checkbox"/> 24 ⁵
<input type="checkbox"/> Reclining with ^{1,7} Adjustable Handle E1226 (Removable)	N/C	<input type="checkbox"/> Short	<input type="checkbox"/> 22
		<input type="checkbox"/> Tall	<input type="checkbox"/> 26

¹ Not available with Seating Dynamics Rocker Back

² Not available with Standard Back Upholstery

³ Not available with Tension Adjustable Back Upholstery

⁴ Not available with Standard Back Upholstery in combination with 21" and 22" Seat Widths

⁵ Not available with Tension Adjustable Back Upholstery in combination with 14" and 15" Seat Widths

⁶ Not available with Tension Adjustable Back Upholstery in combination with 22" Seat Width

⁷ Not available with O₂ Holder, Transit Option, 350 lb Weight Capacity, Swing Away Armrests, Flip Back Armrest, Height Adjustable Flip Back T-Arm or IV Pole

⁸ Not available with Angle Adjustable Locking Extendable Flip Up Armrest

⁹ Not available with Tension Adjustable Back Upholstery in combination with Angle Adjustable Locking Extendable Flip Up Armrests

Reclining Back Post Option

Select options below if Reclining Back is selected on chart above. Recline angle range is 0° - 65°.

- Gas Spring: Standard Spring
 Heavy Duty Spring

Standard Back Post

Select option below only if a Standard Back Post option is selected on chart above

- Straight with Push Handle _____ Std
 Straight without Push Handle _____ N/C
 8° Bend with Push Handle _____ \$30
Not available with Short Back Height Range
 Straight with Fold Down Push Handle _____ \$125
Not available with 350 lb Weight Capacity

Push Handle Options

Not available with Stroller Handle Back Post, Straight Back Post with Fold Down Push Handle, Reclining Back Post or 350 lb Weight Capacity

- Low Profile Push Handle _____ \$85
2" shorter than standard push handles

Depth Adjustable Back Post Options

Options below may be selected only if a Depth Adjustable back post option is selected on the chart above

- Removable Stroller Handle Extension _____ \$85
Not available with Height Adjustable Straight with Push Handle Back Post
 Seating Dynamics Rocker Back K0108 _____ \$1995
Available only with T-Arms

Back Upholstery

Not available with Reclining or Stroller Handle Back Post. If selecting both upholstery and Axiom Back, the upholstery will be installed.

- Standard _____ Std
 Tension Adjustable E2611
 Standard Cover _____ \$315
 Air Mesh Breathable Cover _____ \$390
 Omit Back Upholstery _____ - \$10

Military Patch

Attached only if ordered with Tension Adjustable Back Upholstery

- US Military Service Patch _____ \$90
 Army Navy Marines Air Force Coast Guard
 US Purple Heart Patch _____ \$90
 Purple Heart with Military Service Patch Combo _____ \$145
 Army Navy Marines Air Force Coast Guard

Frame Type / Footrest

Frame Type

- Swing Away _____ N/C
- Super Low - Swing Away _____ N/C
- Fixed Front _____ N/C

Minimum Front Seat Height with Front Mount Hangers						
Hanger Length	Short		Medium		Long	
Hanger Bend	70°	80/90°	70°	80/90°	70°	80°
Minimum Front Seat Height with Standard Selection	13.5"	13.5"	14"	15"	16"	17"
Minimum Front Seat Height with 3" Extension	14.5"	16"	16.5"	18"	18.5"	20"

Hanger / Front Frame Bend

Frame	Hanger	60°	70°	80°	90°	Retail
Swing Away	Extension Tube Mount	<input type="checkbox"/> ²	<input type="checkbox"/>	<input type="checkbox"/> ¹	<input type="checkbox"/>	Std
	Front Mount		<input type="checkbox"/>	<input type="checkbox"/> ¹	<input type="checkbox"/> ^{1,2}	N/C
	Pro Elevating Leg Rests (PRO ELR) E0990				<input type="checkbox"/>	\$295
	Omit Hangers				<input type="checkbox"/>	-\$20
Fixed	Extension Tube Mount		<input type="checkbox"/>	<input type="checkbox"/> ¹		N/C

¹ Not available with 7" or 8" casters

² Not available with Titanium Frame

PRO ELR Cover

Available only on PRO ELR

- Gel Pad Cover _____ \$390

Footrest / Extension / Seat to Footrest Length

Select Extension and Seat to Footrest Length in row corresponding to desired Footrest. Footrest Length of at least 2.5" shorter than Front Seat Height recommended. Extension Tube Hangers and Fixed Front Frame measurements are based on 70° Hanger / Front Frame Bend.

Footrest	Extension	Seat to Footrest Length											
Extension Tube Hangers (60°, 70°, 80°)	<input type="checkbox"/> Short	<input type="checkbox"/> 11	<input type="checkbox"/> 11.5	<input type="checkbox"/> 12	<input type="checkbox"/> 12.5	<input type="checkbox"/> 13	<input type="checkbox"/> 13.5	<input type="checkbox"/> 14	<input type="checkbox"/> 14.5	<input type="checkbox"/> 15			
	<input type="checkbox"/> Medium	<input type="checkbox"/> 13	<input type="checkbox"/> 13.5	<input type="checkbox"/> 14	<input type="checkbox"/> 14.5	<input type="checkbox"/> 15	<input type="checkbox"/> 15.5	<input type="checkbox"/> 16	<input type="checkbox"/> 16.5	<input type="checkbox"/> 17			
	<input type="checkbox"/> Long	<input type="checkbox"/> 15.5	<input type="checkbox"/> 16	<input type="checkbox"/> 16.5	<input type="checkbox"/> 17	<input type="checkbox"/> 17.5	<input type="checkbox"/> 18	<input type="checkbox"/> 18.5	<input type="checkbox"/> 19	<input type="checkbox"/> 19.5			
Front Mount Hangers ¹	<input type="checkbox"/> Short (70°, 80°, 90°)	Standard Selections									Requires 3" Extension for \$35		
		<input type="checkbox"/> 5	<input type="checkbox"/> 5.5	<input type="checkbox"/> 6	<input type="checkbox"/> 6.5	<input type="checkbox"/> 7		<input type="checkbox"/> 10	<input type="checkbox"/> 10.5	<input type="checkbox"/> 11			
	<input type="checkbox"/> Medium (70°, 80°, 90°)	<input type="checkbox"/> 7.5	<input type="checkbox"/> 8	<input type="checkbox"/> 8.5	<input type="checkbox"/> 9	<input type="checkbox"/> 9.5		<input type="checkbox"/> 11.5	<input type="checkbox"/> 12	<input type="checkbox"/> 12.5			
		<input type="checkbox"/> 5	<input type="checkbox"/> 5.5	<input type="checkbox"/> 6	<input type="checkbox"/> 6.5	<input type="checkbox"/> 7	<input type="checkbox"/> 7.5	<input type="checkbox"/> 8	<input type="checkbox"/> 12	<input type="checkbox"/> 12.5	<input type="checkbox"/> 13		
	<input type="checkbox"/> Long (70°, 80°)	<input type="checkbox"/> 8.5	<input type="checkbox"/> 9	<input type="checkbox"/> 9.5	<input type="checkbox"/> 10	<input type="checkbox"/> 10.5	<input type="checkbox"/> 11	<input type="checkbox"/> 11.5	<input type="checkbox"/> 13.5	<input type="checkbox"/> 14	<input type="checkbox"/> 14.5		
		<input type="checkbox"/> 5	<input type="checkbox"/> 5.5	<input type="checkbox"/> 6	<input type="checkbox"/> 6.5	<input type="checkbox"/> 7	<input type="checkbox"/> 7.5	<input type="checkbox"/> 8	<input type="checkbox"/> 8.5	<input type="checkbox"/> 9	<input type="checkbox"/> 14	<input type="checkbox"/> 14.5	<input type="checkbox"/> 15
	<input type="checkbox"/> 9.5	<input type="checkbox"/> 10	<input type="checkbox"/> 10.5	<input type="checkbox"/> 11	<input type="checkbox"/> 11.5	<input type="checkbox"/> 12	<input type="checkbox"/> 12.5	<input type="checkbox"/> 13	<input type="checkbox"/> 13.5	<input type="checkbox"/> 15.5	<input type="checkbox"/> 16	<input type="checkbox"/> 16.5	
PRO ELR Low Pivot Setting	<input type="checkbox"/> Short	12 - 17.5											
	<input type="checkbox"/> Long	15 - 21											
PRO ELR High Pivot Setting	<input type="checkbox"/> Short	10.5 - 15.5											
	<input type="checkbox"/> Long	13 - 19											
Fixed Front Frame (70°, 80°)	<input type="checkbox"/> Short	<input type="checkbox"/> 11.5	<input type="checkbox"/> 12	<input type="checkbox"/> 12.5	<input type="checkbox"/> 13	<input type="checkbox"/> 13.5	<input type="checkbox"/> 14	<input type="checkbox"/> 14.5	<input type="checkbox"/> 15				
	<input type="checkbox"/> Medium	<input type="checkbox"/> 13	<input type="checkbox"/> 13.5	<input type="checkbox"/> 14	<input type="checkbox"/> 14.5	<input type="checkbox"/> 15	<input type="checkbox"/> 15.5	<input type="checkbox"/> 16	<input type="checkbox"/> 16.5	<input type="checkbox"/> 17			
	<input type="checkbox"/> Long	<input type="checkbox"/> 15.5	<input type="checkbox"/> 16	<input type="checkbox"/> 16.5	<input type="checkbox"/> 17	<input type="checkbox"/> 17.5	<input type="checkbox"/> 18	<input type="checkbox"/> 18.5	<input type="checkbox"/> 19	<input type="checkbox"/> 19.5			

¹ Not available with Titanium Frame. See Minimum Front Seat Height chart above.

3" Extension

Available only with Front Mount Hangers.

- 3" Extension _____ \$35
- 3" Extension (Provided as Parts) _____ \$35


Footplate

- Composite _____ N/C
Not available with Front Mount Hanger or PRO ELR
- Composite Angle Adjustable K0040 _____ \$190
Not available with Front Mount Hanger or 14" Seat Width
- Aluminum Angle Adjustable K0040 _____ \$190
- Aluminum Locking Multi-Angle Adjustable K0108 _____ \$300
Available only with Front Mount Hanger
- One Piece Flip Up Angle Adjustable K0040 _____ \$190
Not available with Front Mount Hanger or PRO ELR

Hanger Release

If no selection is made Classic Hanger Release will be used

- Classic _____ Std
- 4-Way - In line _____ N/C
See reference A below. Not available with PRO ELR.
- 4-Way - Out _____ N/C
See reference B below. Not available with PRO ELR.



Frame Type / Footrest continued

- Heel Loop**
- Not available with One Piece Flip Up Angle Adjustable Footplate*
- Heel Loops E0951 \$40
 - Adjustable Heel Loops E0951 \$50
 - Omit Heel Loops N/C
- Calf Strap**
- Velcro Adjustable K0038 \$40
 - Padded Velcro Adjustable K0038 \$68
 - Bodypoint Padded AeroMesh K0038 \$118

- Residual Limb Support**
- Residual Limb Support E1020
- Right and Left \$790
 - Left \$395 Right \$395
 - Omit Left Hanger -\$10 Omit Right Hanger -\$10
- Cover**
- Smooth Air Mesh Smooth Air Mesh
- Pad Size**
- | | | | |
|---------------------------------|----------------------------------|---------------------------------|----------------------------------|
| <input type="checkbox"/> 5 x 8 | <input type="checkbox"/> 9 x 14 | <input type="checkbox"/> 5 x 8 | <input type="checkbox"/> 9 x 14 |
| <input type="checkbox"/> 7 x 10 | <input type="checkbox"/> 11 x 10 | <input type="checkbox"/> 7 x 10 | <input type="checkbox"/> 11 x 10 |
| <input type="checkbox"/> 9 x 10 | <input type="checkbox"/> 11 x 14 | <input type="checkbox"/> 9 x 10 | <input type="checkbox"/> 11 x 14 |

Front Seat Height and Casters

Caster Size / Front Frame Type / Front Seat Height

Select Front Seat Height in row corresponding to desired caster size.
Note: Selections in first three Front Seat Height columns use Super Low Front Frame. Super Low Front Frame not available with Titanium Frame or Fixed Front Frame.

Caster Size	Front Seat Height													
	Uses Super Low Front Frame													
4"	<input type="checkbox"/> 13.5	<input type="checkbox"/> 14	<input type="checkbox"/> 14.5	<input type="checkbox"/> 15	<input type="checkbox"/> 15.5	<input type="checkbox"/> 16	<input type="checkbox"/> 16.5	<input type="checkbox"/> 17	<input type="checkbox"/> 17.5	<input type="checkbox"/> 18	<input type="checkbox"/> 18.5	<input type="checkbox"/> 19	<input type="checkbox"/> 19.5	<input type="checkbox"/> 20
5"	<input type="checkbox"/> 14.5	<input type="checkbox"/> 15	<input type="checkbox"/> 15.5	<input type="checkbox"/> 16	<input type="checkbox"/> 16.5	<input type="checkbox"/> 17	<input type="checkbox"/> 17.5	<input type="checkbox"/> 18	<input type="checkbox"/> 18.5	<input type="checkbox"/> 19	<input type="checkbox"/> 19.5	<input type="checkbox"/> 20	<input type="checkbox"/> 20.5	
6"	<input type="checkbox"/> 15.5	<input type="checkbox"/> 16	<input type="checkbox"/> 16.5	<input type="checkbox"/> 17	<input type="checkbox"/> 17.5	<input type="checkbox"/> 18	<input type="checkbox"/> 18.5	<input type="checkbox"/> 19	<input type="checkbox"/> 19.5	<input type="checkbox"/> 20	<input type="checkbox"/> 20.5	<input type="checkbox"/> 21		
7"	<input type="checkbox"/> 17	<input type="checkbox"/> 17.5		<input type="checkbox"/> 18	<input type="checkbox"/> 18.5	<input type="checkbox"/> 19	<input type="checkbox"/> 19.5	<input type="checkbox"/> 20	<input type="checkbox"/> 20.5	<input type="checkbox"/> 21				
8"	<input type="checkbox"/> 17.5	<input type="checkbox"/> 18	<input type="checkbox"/> 18.5	<input type="checkbox"/> 19	<input type="checkbox"/> 19.5	<input type="checkbox"/> 20	<input type="checkbox"/> 20.5	<input type="checkbox"/> 21						

Caster

Caster Type	Price	Caster Size				
		4"	5"	6"	7"	8"
.75" Rollerblade	N/C	<input type="checkbox"/>				
.75" Lighted Rollerblade	\$70	<input type="checkbox"/>				
1" Poly	N/C	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>
1" Poly Aluminum	\$110	<input type="checkbox"/>	<input type="checkbox"/>			
1" Pneumatic	\$85			<input type="checkbox"/>		<input type="checkbox"/>
1.5" Poly	\$50	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>		
1.5" Soft Roll Aluminum	\$185	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>		
2" Poly	\$70			<input type="checkbox"/>		
2" Pneumatic w/ Foam Insert	\$225					<input type="checkbox"/>

Maximum Seat Dump Reference
 (Front Seat Height minus Rear Seat Height)

Seat Depth	14"	15"	16"	17"	18"	19"	20"
Maximum Dump Standard Caster Housing	1.5"	2"	2"	2"	2.5"	2.5"	2.5"
Maximum Dump Performance Caster Housing	3"	3"	3"	3.5"	3.5"	4"	4"

Depending on configuration, chairs with Fixed Front Frame Type may achieve up to an additional .75"

- Fork Options**
- Frog Legs Suspension Fork \$450
Available only with 4", 5" or 6" Caster Sizes. Not available with Caster Pin Locks, 350 lb Weight Capacity, 2" Poly caster or 2" Pneumatic w/ Foam Insert caster.
 - User Weight Under 80 lbs
 - User Weight 80 - 120 lbs
 - User Weight 120 - 160 lbs
 - User Weight 160 - 220 lbs
 - User Weight 220 - 260 lbs
 - User Weight 260 - 300 lbs
 - Frog Legs II Suspension System \$450
Lightweight shock absorbing fork. Available only with 4", 5" or 6" Caster Sizes. Not available with Caster Pin Locks, 350 lb Weight Capacity, 2" Poly caster or 2" Pneumatic w/ Foam Insert caster.
 - Caster Pin Locks K0073 \$195
Not available with Performance Caster Housing, 1 1/2" Longer Fork Stem or Frog Legs Fork

- Fork Stem**
- Best Fit Std
Ki Mobility will pick Fork and Stem combination to optimize performance based on Caster Size and Front Seat Height
 - Standard Fork Stem N/C
 - 3/4" Longer Fork Stem N/C
 - 1 1/2" Longer Fork Stem N/C

- Caster Housing**
- Standard N/C
 - Performance \$40
Not available with 350 lb Weight Capacity

Wheel Locks

Wheel Lock		
<input type="checkbox"/> Push to Lock	Std	
<input type="checkbox"/> Push to Lock Flush Mount	N/C	
<input type="checkbox"/> Pull to Lock	N/C	
<input type="checkbox"/> Low Profile Scissor Lock	\$65	
<i>Not available with Low Poly, Full Poly or Sentinel Solid tires</i>		
<input type="checkbox"/> Short Thro Scissor Lock	\$110	
<input type="checkbox"/> Drum Brake E2228	\$560	
<i>Available only with 20", 22" and 24" MAXX Performance Spoke wheels. Not available with Reclining Backrest, Straight without Push Handle or Straight with Fold Down Push Handle Back Posts, or Low Profile Push Handle. Axle is Quick Release.</i>		
Dual - Drum Brake E2228		\$660
<i>Available only with 20", 22" and 24" MAXX Performance Spoke wheels. Not available with Reclining Backrest, Straight without Push Handle or Straight with Fold Down Push Handle Back Posts, or Low Profile Push Handle. Axle is Quick Release.</i>		
<input type="checkbox"/> Push to Lock	<input type="checkbox"/> Pull to Lock	
<input type="checkbox"/> Omit Wheel Locks		N/C

Wheel Lock Options

<input type="checkbox"/> Extension Handles E0961	\$75
<i>Available only with Push to Lock or Pull to Lock Wheel Locks</i>	
<input type="checkbox"/> Grade Aids E0974	\$150
<i>Available only with Push to Lock Wheel Locks and Pneumatic or Pneumatic with Airless Insert tires</i>	

Rear Seat Height and Wheels

Wheel Size / Tire Profile / Rear Seat Height

Select Rear Seat Height in row corresponding to Wheel Size and Tire Profile. 14" Seat Depth not available with 24", 25" or 26" wheel sizes.

Requires Universal Axle Plate Any axle plate can be used (Universal, Amputee, One Arm Drive)

Wheel Size	Tire Profile	Rear Seat Height									
		<input type="radio"/> 11.5	<input type="checkbox"/> 12 ¹	<input type="radio"/> 12.5	<input type="checkbox"/> 13 ¹	<input type="radio"/> 13.5	<input type="checkbox"/> 14	<input type="radio"/> 14.5	<input type="checkbox"/> 15	<input type="radio"/> 15.5	<input type="radio"/> 16
20"	Low	<input type="radio"/> 11.5	<input type="checkbox"/> 12 ¹	<input type="radio"/> 12.5	<input type="checkbox"/> 13 ¹	<input type="radio"/> 13.5	<input type="checkbox"/> 14	<input type="radio"/> 14.5	<input type="checkbox"/> 15	<input type="radio"/> 15.5	<input type="radio"/> 16
	Full	<input type="radio"/> 12	<input type="checkbox"/> 12.5 ¹	<input type="radio"/> 13	<input type="checkbox"/> 13.5 ¹	<input type="radio"/> 14	<input type="checkbox"/> 14.5	<input type="radio"/> 15	<input type="checkbox"/> 15.5	<input type="radio"/> 16	<input type="radio"/> 16.5
22"	Low	<input type="radio"/> 12.5	<input type="checkbox"/> 13 ¹	<input type="radio"/> 13.5	<input type="checkbox"/> 14 ¹	<input type="radio"/> 14.5	<input type="checkbox"/> 15	<input type="radio"/> 15.5	<input type="checkbox"/> 16	<input type="radio"/> 16.5	<input type="radio"/> 17
	Full	<input type="radio"/> 13	<input type="checkbox"/> 13.5 ¹	<input type="radio"/> 14	<input type="checkbox"/> 14.5 ¹	<input type="radio"/> 15	<input type="checkbox"/> 15.5	<input type="radio"/> 16	<input type="checkbox"/> 16.5	<input type="radio"/> 17	<input type="radio"/> 17.5
24"	Low	<input type="radio"/> 13.5	<input type="checkbox"/> 14 ¹	<input type="radio"/> 14.5	<input type="checkbox"/> 15 ¹	<input type="radio"/> 15.5	<input type="checkbox"/> 16	<input type="radio"/> 16.5	<input type="checkbox"/> 17	<input type="radio"/> 17.5	<input type="radio"/> 18
	Full	<input type="radio"/> 14	<input type="checkbox"/> 14.5 ¹	<input type="radio"/> 15	<input type="checkbox"/> 15.5 ¹	<input type="radio"/> 16	<input type="checkbox"/> 16.5	<input type="radio"/> 17	<input type="checkbox"/> 17.5	<input type="radio"/> 18	<input type="radio"/> 18.5
25"	Low	<input type="radio"/> 14	<input type="checkbox"/> 14.5 ¹	<input type="radio"/> 15	<input type="checkbox"/> 15.5 ¹	<input type="radio"/> 16	<input type="checkbox"/> 16.5	<input type="radio"/> 17	<input type="checkbox"/> 17.5	<input type="radio"/> 18	<input type="radio"/> 18.5
26"	Low		<input type="checkbox"/> 15 ^{1,2}	<input type="radio"/> 15.5	<input type="checkbox"/> 16 ¹	<input type="radio"/> 16.5	<input type="checkbox"/> 17	<input type="radio"/> 17.5	<input type="checkbox"/> 18	<input type="radio"/> 18.5	<input type="radio"/> 19
	Full	<input type="radio"/> 15 ²	<input type="checkbox"/> 15.5 ^{1,2}	<input type="radio"/> 16	<input type="checkbox"/> 16.5 ¹	<input type="radio"/> 17	<input type="checkbox"/> 17.5	<input type="radio"/> 18	<input type="checkbox"/> 18.5	<input type="radio"/> 19	<input type="radio"/> 19.5

¹ Not available with Amputee Axle Plate

² Not available with Height Adjustable Flip Back T-Arm

One Arm Drive

<input type="checkbox"/> One Arm Drive E0958	\$1095
<i>Available only with 20", 22" and 24" spoke wheels and Quick Release Axles. Spoke wheel style is exclusive to OAD. Not available with 4 degrees of camber. Standard handrim will be Aluminum Anodized.</i>	
Rear Wheel Size:	<input type="checkbox"/> 20" <input type="checkbox"/> 22" <input type="checkbox"/> 24"
Drive Handrim Side:	<input type="checkbox"/> Left <input type="checkbox"/> Right
Drive Handrim Type:	<input type="checkbox"/> Aluminum N/C <input type="checkbox"/> Poly Coated \$200

Axle Plate

No selection needed with One Arm Drive

<input type="radio"/> Universal	Std
<input type="checkbox"/> Amputee E0959	\$175
<i>Not available with Swing Away Armrests or Wheelchair Drive Ready Option. Required with Reclining Backrest.</i>	

Rear Wheel Type

No selection needed with One Arm Drive

Rear Wheel Type	Price	Rear Wheel Size				
		20"	22"	24"	25"	26"
MAXX Performance Spoke	\$200	<input type="checkbox"/>	<input type="checkbox"/>		<input type="checkbox"/>	<input type="checkbox"/>
MAXX Performance Spoke	N/C			<input type="checkbox"/>		
Mag	N/C	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>		
Superlight Spoke ¹	\$265			<input type="checkbox"/>	<input type="checkbox"/> ²	
Spinergy Spox	\$845		<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>
Spinergy LX	\$1050		<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>

¹ Not available with 350 lb Weight Capacity

² Available only with Superlight and Plastic Coated Handrims

Spinergy Options

Select spoke color and handrim tabs for Spinergy wheels

Spoke Color:

- Black Red White Blue Yellow
 Pink ¹ Green ¹ Purple ¹ Orange ¹

¹ Available only with Spinergy LX

Handrim Tabs:

Available only with Aluminum Anodized and Plastic Coated Handrims

- Short Tab \$25
 Long Tab N/C

Rear Seat Height and Wheels continued

Handrim

No selection needed with One Arm Drive

Standard Handrims

- Aluminum Anodized _____ Std
- Superlight Aluminum Anodized _____ \$50
Available only on 24" and 25" wheel sizes. Available only on MAXX Performance, Superlight or Spinerig wheels. Not available with 350 lb Weight Capacity.
- Plastic Coated _____ \$100
- Projection (8 Vertical Projections) _____ \$325
Available only with full profile tires in combination with 22"-26" MAXX Performance, or 24" Mag wheel
- Spinerig Flex Rim _____ \$700
Available only on 24" and 25" Spinerig LX wheels
- Omit Handrim _____ N/C

Natural Fit Handrims

Not available on 20" wheel size or 22" Mag wheel. Not recommended with Mag wheel.

- Natural Fit LT Regular Grip _____ \$465
- Natural Fit Standard Regular Grip _____ \$445
- Natural Fit LT Super Grip _____ \$515
- Natural Fit Standard Super Grip _____ \$495
- Natural Fit LT Omit Grip _____ \$395
- Natural Fit Standard Omit Grip _____ \$345

Handrim Option

- 3M Micro Grip _____ \$65
Not available on Plastic Coated or Projection Handrims

Rear Tire Type

Full Profile

- Street Pneumatic 1 3/8" E2211, E2212 _____ Std
Not available on 25" wheel size
- Full Poly 1 3/8" E2218 _____ N/C
Not available on Superlight wheel or 25" or 26" wheel sizes
- Pneumatic with Airless Insert 1 3/8" E2211, E2213 _____ \$85
Not available on Superlight wheel or 25" wheel size
- Schwalbe HS-127 1 3/8" E2211, E2212 _____ \$125
Available only on 24" wheel size. Not available on Mag wheel.
- Euro Trek Knobby 1 3/8" E2211, E2212 _____ \$95
Available only on 24" wheel size

Low Profile

- Low Poly 1" E2218 _____ N/C
Not available on Superlight wheel, 25" or 26" wheel sizes or 350 lb Weight Cap.
- Sentinel Solid Tire 1" E2218 _____ \$95
Not available on Superlight wheel or 20" and 22" Mag wheels
- High Pressure 1" E2211, E2212 _____ \$95
Not available on Mag wheel
- Schwalbe One 1" E2211, E2212 _____ \$135
Not available on Mag wheel, 20", 22" or 26" wheel sizes
- Schwalbe Marathon Plus Evolution 1" E2211, E2212 _____ \$145
Not available on Mag wheel, 20" or 22" wheel sizes

Axle

No selection needed with One Arm Drive or Drum Brake

- Quick Release _____ Std
- Quad Quick Release _____ \$135
Not available with Mag or Spinerig wheels

Axle Position

- Rear Wheel Camber: 0° 2° 4°
- Center of Gravity: Rear Center Forward
No selection needed with Amputee Axle Plate

Armrests

Armrest

- Height Adjustable T-Arm (8.5" - 13.5") E0973
- Desk _____ \$295
 - Full _____ \$295
- Height Adjustable T-Arm Low (7.5" - 11.5") E0973
If selection is not available with rear wheel and rear seat height configuration, standard Height Adjustable T-Arm armrests will be used
- Desk _____ \$295
 - Full _____ \$295
- Height Adjustable Flip Back T-Arm (9.5"-13.5") E0973
Available only with Standard Back Posts
- Desk _____ \$305
 - Full _____ \$305
- Flip Back (10.5")
Not available with Seat Depths less than 16". Available only with Standard Back Posts.
- Desk _____ N/C
 - Full _____ N/C

- Angle Adjustable Locking Extendable Flip Up E0973
Available only with Depth Adjustable and Reclining Back Posts
- Desk _____ \$275
 - Full _____ \$275
- Swing Away (10.5" - 12.5") E0973 _____ \$225
Available only with Standard Back Posts. Not available with Arm Pads.
 - Swing Away with Side Guard (10.5" - 12.5") E0973 _____ \$300
Available only with Standard Back Posts. Not available with Arm Pads.
 - Omit Armrests _____ N/C

Arm Pad

Not available with Swing Away Armrests

- Standard Arm Pad _____ N/C
- Waterfall Arm Pad _____ \$25
Available only with Height Adjustable T-Arm and Height Adjustable T-Arm Low
- Foam Grip _____ N/C
Available only with Angle Adjustable Locking Extendable Flip Up Armrest

Accessories

Removable Side Guard

- Not available with T-Arms, Flip Back or Swing Away with Side Guard Armrests**
- Composite Side Guards K0108 \$145
 - Aluminum Straight K0108
 - Black Matte \$195
 - Match Frame Color \$245
Not available with Hydrographics or Titanium Frame
 - Aluminum Adjustable Fender with Foam Pad K0108
 - Black Matte \$275
 - Match Frame Color \$325
Not available with Hydrographics or Titanium Frame
 - Carbon Fiber Adjustable Fender K0108 \$485

Pelvic Positioning Belt

- 1 1/2" Auto Buckle E0978 \$80
- 1 1/2" Padded Auto Buckle E0978 \$170
- 2" Airline Buckle E0978 \$110
- Bodypoint Evoflex - Pelvic Stabilizer E0978 \$188
 - Size: Small (16-24") Medium (19-29") Large (25-37")
- Bodypoint Belt Mounting Kit \$108
 - Style: Bracket Band Clamp 1"

O₂ Holder

- Not available with Depth Adjustable Rear Frame**
- O₂ Holder E2208 \$295
 - Right Side Left Side

Color

Frame Color not available with Titanium Frame Option

Gloss Frame Finish

- Black N/C
- Blue Skies N/C
- Bubblemum Pink N/C
- Candy Blue N/C
- Candy Purple N/C
- Candy Red N/C
- Copper N/C
- Granny Smith N/C
- Metallic Orange N/C
- Yellow N/C

Matte Frame Finish

- Charcoal Gray N/C
- Emerald Green N/C
- Iridescent Burgundy N/C
- Iridescent Navy N/C
- Military Green N/C
- Retro Red N/C
- Royal Blue N/C
- Shadow Black N/C
- Smooth White N/C

Anti-Tip

- Rear Anti-Tips E0971 \$150
- Swing-In Anti-Tip E0971
 - Not available with 350 lb Weight Capacity or Reclining Back Post*
 - Single - Left \$155
Not available with left side O₂ Holder
 - Single - Right \$155
Not available with right side O₂ Holder
 - Pair \$310
Not available with Cane and Crutch Holder or O₂ Holder
- Omit Rear Anti-Tips N/C

Other Accessories

- Tool Kit \$45
- Removable Underseat Pouch \$50
- Pediatric Ki Mobility Backpack \$50
- Ki Mobility Backpack \$65
- Luggage Carrier \$85
Not available with PRO ELR or Residual Limb Support
- Black Neoprene Impact Guards \$90
- Clear Spoke Guards K0065 \$115
Not available with Mag wheel
- Black Spoke Guards K0065 \$115
Not available with 20" wheel size or Mag wheel
- Cane and Crutch Holder \$195
- IV Pole \$275
Available only with Curved Rear Frame Type and Straight Back Posts

Hydrographics

- Hydrographics \$695
 - AOR1 (Desert Digital Camo)
 - Pink Camo
 - Muddy Girl
 - NEXT Camo
 - Carbon Fiber¹
 - Zebra¹
- ¹Select Hydrographics Base Color. Zebra not available in Black.
- Hydrographics Base Color N/C
- Black Blue Granny Smith Pink Purple Red

Colored Anodizing Package

Select desired color. If no selection is made, black will be used.

	Option	Standard Caster Housing ²	Caster Fork ³	Axle Plate ⁴
	Retail	\$75	\$75	\$75
Color	Black ¹	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>
	Blue	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>
	Bright Green	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>
	Gold	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>
	Orange	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>
	Pewter	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>
	Pink	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>
	Purple	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>
	Red	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>

¹Black is no charge
²Available only if corresponding option selection was made
³Not available with 2" wide casters, Frog Legs Suspension Fork or Frog Legs II Suspension System
⁴Color Anodized not available with amputee axle plate

Catalyst Front Seat Height Matrix																																												
4" Caster				5" Caster				6" Caster				7" Caster				8" Caster																												
FSH	Super Low	Stem	Fork	FSH	Super Low	Stem	Fork	FSH	Super Low	Stem	Fork	FSH	Super Low	Stem	Fork	FSH	Super Low	Stem	Fork																									
15	13.5	Std	4"	16	14.5	Std	4"	17	15.5	Std	5"	18	16.5	Std	6"	19	17.5	Std	6"																									
15.5	14			16.5	15.5	3/4		17.5	16.5	3/4		18.5	17			19.5	18.5	3/4																										
15.5 ¹	14 ¹			17.5	16	1 1/2		18.5	17	1 1/2		18.5	17.5			20.5	n/a	1 1/2																										
16	14.5	3/4		16	14.5	Std	5"	17	15.5	Std	6"	19	18	3/4		19	17.5	Std	7"																									
16.5 ¹	15 ¹	1 1/2		16.5	15			17.5	16			19.5	18			19.5	18																											
17	15.5	3/4		16.5	15.5			18	16.5			20	18.5			20	18.5																											
15.5	14	Std		5"	17	16		3/4	5"	17.5		16.5	3/4	6"		18.5	17	1 1/2	20.5	19	Std	7"																						
16	14.5				17.5	16				18		17				19	17.5		19.5	18																								
16	15				18	16.5				18.5		17.5				19.5	18		20	19																								
16.5	15.5	3/4			16.5	15		Std		5"		18.5	17			1 1/2	6"	20	18.5	3/4	20.5	n/a	3/4	7"																				
17	15.5	1 1/2	17		15.5	19						17.5	19		18			21	n/a																									
17.5	16	17.5	16		19.5	18						19.5	18.5		20.5			19	21		n/a																							
16.5	15	Std	6"		17.5	16		3/4				6"	18		16.5	Std		7"	20	19	3/4	7"	21	n/a	1 1/2	7"																		
17	15.5				18	16.5	18.5				17		20.5		n/a																													
17	16				18.5	17	19				17.5		20		18.5																													
17.5	16.5	3/4			18	16.5	1 1/2	6"			19.5		18		1 1/2	7"			20.5	19	1 1/2						7"	21	n/a	1 1/2	7"													
18	16.5	1 1/2		18.5	17	18.5			17.5		19		18	21					n/a																									
18.5	17	19		17.5	19	18			19.5		18.5		20	19																														
18	16.5	Std		7"	18	16.5	Std		7"		19.5		18.5	3/4	7"				21	n/a	1 1/2											7"	21	n/a	1 1/2	7"								
18.5	17				18.5	17				20	19		19.5				18																				20	18.5						
18.5	17.5				3/4	19				17.5	19.5		18				20																				18.5	20.5	19					
19	18	1 1/2			7"	18.5	17.5			3/4	7"		20.5	19			1 1/2																				7"	21	n/a	1 1/2	7"	21	n/a	1 1/2
19.5	18.5		19			18	20.5					19	19.5	18.5				20				18.5																						
20	19		19.5			18.5	21					n/a	6" x 2" caster not available with 5" fork																															
			19.5			18	1 1/2			7"		7"	21	n/a			1 1/2	7"				21	n/a	1 1/2	7"	21																		
			20			18.5		20.5								19											6" x 2" caster not available with 5" fork																	
			20.5			19		6" x 2" caster not available with 5" fork																																				

¹ Not available with 4 x 1.5" Soft Roll Aluminum

Seat heights with "FSH" heading available with Standard Front Frame. Seat heights with "Super Low" heading available with Super Low Front Frame. Actual seat height may vary depending upon configuration.