



# CATALYST 5 TTL

Phone (800) 981-1540 | Fax (715) 254-0996  
 sales@kimobility.com | www.kimobility.com

Company: \_\_\_\_\_

Ship To: \_\_\_\_\_

Account Number: \_\_\_\_\_ Date: \_\_\_\_\_

Address: \_\_\_\_\_

Quote  PO Number: \_\_\_\_\_

City: \_\_\_\_\_ State: \_\_\_\_\_

RTS: \_\_\_\_\_

Zip: \_\_\_\_\_ Phone: \_\_\_\_\_

Marked For: \_\_\_\_\_

Email/Fax: \_\_\_\_\_

Select if chair quote has:  Additional Ki Products with Different Quote/Order Number **Quote/Order Number:** \_\_\_\_\_

**Catalyst 5TTL K0005 \$2515**  
 300 lb weight capacity

Effective 01/17/2022

## Frame

### Seat Width

- 14" N/C
- 15" N/C
- 16" N/C
- 17" N/C
- 18" N/C
- 19" N/C
- 20" E2201 \$505
- 21" E2201 \$505
- 22" E2201 \$505

### Seat Depth

- 14" N/C
- 15" N/C
- 16" N/C
- 17" N/C
- 18" N/C
- 19" N/C
- 20" E2203 \$510

### Transit

- Transit System \$325  
*Not available with Removable Seat Pan or 350 lb Weight Capacity*
- Omit Transit N/C

### Wheelchair Drive Ready

**Chair will ship from Ki Mobility with adaption kit installed for selected Alber system. Order system separately from Alber USA. Select Rear Seat Height (full profile selections only) in drive wheel size on Rear Seat Height chart (page 5). Requires Rear Anti-Tips. Not available with Spinergy wheels or Height Adjustable Flip Back T-Arm.**

- e-motion Ready N/C
- twion Ready N/C

### Titanium Option

- Titanium Frame K0108 \$1775  
*Includes hangers and base frame. Does not include Cross Braces.*

### Weight Capacity

- 300 lb Std
- 350 lb (Heavy Duty) \$500  
*Requires minimum 6" caster height*

### Cross Brace

- Standard Cross Braces Std
- Seating System Cross Braces N/C  
*Allows for easy attachment of seating system mounting hardware. Not available with Seat Sling.*
- Titanium Cross Braces \$700  
*Not available with 21" and 22" Seat Widths. 250 lb maximum weight capacity.*

### Seat Pan

- Removable Aluminum Seat Pan E2231 \$275  
*Not available with Transit System. Pan depth 1/2" shorter than Seat Depth.*

### Rear Frame Type

**No selection needed with Titanium Frame**

- Curved Std  
*Select from Standard or Angle Adjustable Back Posts on chart on the next page*
- Depth Adjustable K0108 \$105  
*Select from depth adjustable back posts in chart on next page. Provides 2" forward and 2" rearward adjustment with minimum Seat Depth of 14" and maximum Seat Depth of 22". Available only with T-Arm or Flip Up Armrests.*
- Reclining \$1095  
*Select from Reclining Back Post in chart on the next page. Not available with Seating Dynamics Rocker Back, O<sub>2</sub> Holder, Transit Option, 350 lb Weight Capacity, Swing Away Armrests, Flip Back Armrests, Height Adjustable Flip Back Armrests, Back Upholstery or IV Pole.*

HCPCS codes provided are based on code verification by the PDAC or our interpretation of the code definitions.

Proper use of HCPCS codes does not ensure coverage or payment. For coverage information, verify the appropriate payer's coverage policy.

**Backrest**

**Back Post / Height Range / Back Height**

*Depth Adjustable not available with Titanium Frame.*

Back Post	Price	Height Range	Back Height				
<input type="checkbox"/> Standard	N/C	<input type="checkbox"/> Short <sup>3,4</sup>	<input type="checkbox"/> 8.5	<input type="checkbox"/> 9.5	<input type="checkbox"/> 10.5	<input type="checkbox"/> 11.5	
		<input type="checkbox"/> Medium <sup>4</sup>	<input type="checkbox"/> 12.5 <sup>6</sup>	<input type="checkbox"/> 13.5	<input type="checkbox"/> 14.5	<input type="checkbox"/> 15.5	
		<input type="checkbox"/> Tall	<input type="checkbox"/> 16.5	<input type="checkbox"/> 17.5	<input type="checkbox"/> 18.5	<input type="checkbox"/> 19.5	
<input type="checkbox"/> Angle Adjustable with Push Handle	\$185	<input type="checkbox"/> Short <sup>4</sup>	<input type="checkbox"/> 14	<input type="checkbox"/> 15	<input type="checkbox"/> 16		
		<input type="checkbox"/> Medium <sup>4</sup>	<input type="checkbox"/> 16	<input type="checkbox"/> 17	<input type="checkbox"/> 18		
		<input type="checkbox"/> Tall	<input type="checkbox"/> 18	<input type="checkbox"/> 19	<input type="checkbox"/> 20		
<input type="checkbox"/> Stroller Handle Back Post (Depth Adjustable)	N/C	<input type="checkbox"/> Short	<input type="checkbox"/> 20				
		<input type="checkbox"/> Tall	<input type="checkbox"/> 24				
<input type="checkbox"/> Height Adjustable <sup>2,9</sup> Straight with Push Handle (Depth Adjustable)	\$85	<input type="checkbox"/> Short <sup>1</sup>	<input type="checkbox"/> 13	<input type="checkbox"/> 14	<input type="checkbox"/> 15	<input type="checkbox"/> 16	<input type="checkbox"/> 17
		<input type="checkbox"/> Medium	<input type="checkbox"/> 16	<input type="checkbox"/> 17	<input type="checkbox"/> 18	<input type="checkbox"/> 19	<input type="checkbox"/> 20
		<input type="checkbox"/> Tall	<input type="checkbox"/> 20 <sup>8</sup>	<input type="checkbox"/> 21 <sup>5</sup>	<input type="checkbox"/> 22 <sup>5</sup>	<input type="checkbox"/> 23 <sup>5</sup>	<input type="checkbox"/> 24 <sup>5</sup>
<input type="checkbox"/> Reclining with <sup>1,7</sup> Adjustable Handle <small>E1226 (Removable)</small>	N/C	<input type="checkbox"/> Short	<input type="checkbox"/> 22				
		<input type="checkbox"/> Tall	<input type="checkbox"/> 26				

<sup>1</sup> Not available with Seating Dynamics Rocker Back

<sup>2</sup> Not available with Standard Back Upholstery

<sup>3</sup> Not available with Tension Adjustable Back Upholstery

<sup>4</sup> Not available with Standard Back Upholstery in combination with 21"-22" Seat Widths

<sup>5</sup> Not available with Tension Adjustable Back Upholstery in combination with 14" and 15" Seat Widths

<sup>6</sup> Not available with Tension Adjustable Back Upholstery in combination with 22" Seat Width

<sup>7</sup> Not available with O<sub>2</sub> Holder, Transit Option, 350 lb Weight Capacity, Swing Away Armrests, Flip Back Armrest, Height Adjustable Flip Back T-Arm or IV Pole

<sup>8</sup> Not available with Angle Adjustable Locking Extendable Flip Up Armrest

<sup>9</sup> Not available with Tension Adjustable Back Upholstery in combination with Angle Adjustable Locking Extendable Flip Up Armrests

**Reclining Back Post Option**

Select options below if Reclining Back is selected on chart above. Recline angle range is 0° - 65°.

- Gas Spring:  Standard Spring  
 Heavy Duty Spring

**Standard Back Post**

Select option below only if a Standard Back Post option is selected on chart above

- Straight with Push Handle \_\_\_\_\_ Std  
 Straight without Push Handle \_\_\_\_\_ N/C  
 8° Bend with Push Handle \_\_\_\_\_ N/C  
Not available with Short Back Height Range  
 Straight with Fold Down Push Handle \_\_\_\_\_ \$125  
Not available with 350 lb Weight Capacity

**Push Handle Options**

Not available with Stroller Handle Back Post, Straight Back Post with Fold Down Push Handle, Reclining Back Post or 350 Weight Capacity

- Low Profile Push Handle \_\_\_\_\_ \$85  
2" shorter than standard push handles

**Depth Adjustable Back Post Options**

Options below may be selected only if a Depth Adjustable back post option is selected on the chart above

- Removable Stroller Handle Extension \_\_\_\_\_ \$85  
Not available with Height Adjustable Straight with Push Handle Back Post  
 Seating Dynamics Rocker Back K0108 \_\_\_\_\_ \$1995  
Available only with T-Arms

**Back Upholstery**

Not available with Reclining or Stroller Handle back post. If selecting both upholstery and Axiom Back, the upholstery will be installed.

- Standard \_\_\_\_\_ Std  
 Tension Adjustable E2611  
 Standard Cover \_\_\_\_\_ \$315  
 Air Mesh Breathable Cover \_\_\_\_\_ \$390  
 Omit Back Upholstery \_\_\_\_\_ - \$10

**Military Patch**

Attached only if ordered with Tension Adjustable Back Upholstery

- US Military Service Patch \_\_\_\_\_ \$90  
 Army  Navy  Marines  Air Force  Coast Guard  
 US Purple Heart Patch \_\_\_\_\_ \$90  
 Purple Heart with Military Service Patch Combo \_\_\_\_\_ \$145  
 Army  Navy  Marines  Air Force  Coast Guard

**Frame Type / Footrest**

**Frame Type**

- Swing Away \_\_\_\_\_ N/C
- Super Low - Swing Away \_\_\_\_\_ N/C
- Fixed Front \_\_\_\_\_ N/C

Minimum Front Seat Height with Front Mount Hangers						
Hanger Length	Short		Medium		Long	
	70°	80/90°	70°	80/90°	70°	80°
Minimum Front Seat Height with Standard Selection	13.5"	13.5"	14"	15"	16"	17"
Minimum Front Seat Height with 3" Extension	14.5"	16"	16.5"	18"	18.5"	20"

**Hanger / Front Frame Bend**

Frame	Hanger	60°	70°	80°	90°	Retail
Swing Away	Extension Tube Mount	<input type="checkbox"/> <sup>2</sup>	<input type="checkbox"/>	<input type="checkbox"/> <sup>1</sup>		Std
	Front Mount		<input type="checkbox"/>	<input type="checkbox"/> <sup>1</sup>	<input type="checkbox"/> <sup>1</sup>	N/C
	Pro Elevating Leg Rests (PRO ELR) E0990					<input type="checkbox"/> \$295
	Omit Hangers					<input type="checkbox"/> N/C
Fixed	Extension Tube Mount		<input type="checkbox"/>	<input type="checkbox"/> <sup>1</sup>		N/C

<sup>1</sup> Not available with 7" or 8" casters  
<sup>2</sup> Not available with Titanium Frame

**PRO ELR Cover**

- Available only on PRO ELR
- Gel Pad Cover \_\_\_\_\_ \$390

**Footrest / Extension / Seat to Footrest Length**

Select Extension and Seat to Footrest Length in row corresponding to desired Footrest. Footrest Length of at least 2.5" shorter than Front Seat Height recommended. Extension Tube Hangers and Fixed Front Frame measurements are based on 70° Hanger / Front Frame Bend.

Footrest	Extension	Seat to Footrest Length													
Extension Tube Hangers (60° <sup>1</sup> , 70°, 80°)	<input type="checkbox"/> Short	<input type="checkbox"/> 11	<input type="checkbox"/> 11.5	<input type="checkbox"/> 12	<input type="checkbox"/> 12.5	<input type="checkbox"/> 13	<input type="checkbox"/> 13.5	<input type="checkbox"/> 14	<input type="checkbox"/> 14.5	<input type="checkbox"/> 15					
	<input type="checkbox"/> Medium	<input type="checkbox"/> 13	<input type="checkbox"/> 13.5	<input type="checkbox"/> 14	<input type="checkbox"/> 14.5	<input type="checkbox"/> 15	<input type="checkbox"/> 15.5	<input type="checkbox"/> 16	<input type="checkbox"/> 16.5	<input type="checkbox"/> 17					
	<input type="checkbox"/> Long	<input type="checkbox"/> 15.5	<input type="checkbox"/> 16	<input type="checkbox"/> 16.5	<input type="checkbox"/> 17	<input type="checkbox"/> 17.5	<input type="checkbox"/> 18	<input type="checkbox"/> 18.5	<input type="checkbox"/> 19	<input type="checkbox"/> 19.5					
Front Mount Hangers <sup>1</sup>	Standard Selections											Requires 3" Extension			
	<input type="checkbox"/> Short (70°, 80°, 90°)	<input type="checkbox"/> 5	<input type="checkbox"/> 5.5	<input type="checkbox"/> 6	<input type="checkbox"/> 6.5	<input type="checkbox"/> 7							<input type="checkbox"/> 10	<input type="checkbox"/> 10.5	<input type="checkbox"/> 11
		<input type="checkbox"/> 7.5	<input type="checkbox"/> 8	<input type="checkbox"/> 8.5	<input type="checkbox"/> 9	<input type="checkbox"/> 9.5							<input type="checkbox"/> 11.5	<input type="checkbox"/> 12	<input type="checkbox"/> 12.5
	<input type="checkbox"/> Medium (70°, 80°, 90°)	<input type="checkbox"/> 5	<input type="checkbox"/> 5.5	<input type="checkbox"/> 6	<input type="checkbox"/> 6.5	<input type="checkbox"/> 7	<input type="checkbox"/> 7.5	<input type="checkbox"/> 8					<input type="checkbox"/> 12	<input type="checkbox"/> 12.5	<input type="checkbox"/> 13
		<input type="checkbox"/> 8.5	<input type="checkbox"/> 9	<input type="checkbox"/> 9.5	<input type="checkbox"/> 10	<input type="checkbox"/> 10.5	<input type="checkbox"/> 11	<input type="checkbox"/> 11.5					<input type="checkbox"/> 13.5	<input type="checkbox"/> 14	<input type="checkbox"/> 14.5
	<input type="checkbox"/> Long (70°, 80°)	<input type="checkbox"/> 5	<input type="checkbox"/> 5.5	<input type="checkbox"/> 6	<input type="checkbox"/> 6.5	<input type="checkbox"/> 7	<input type="checkbox"/> 7.5	<input type="checkbox"/> 8	<input type="checkbox"/> 8.5	<input type="checkbox"/> 9			<input type="checkbox"/> 14	<input type="checkbox"/> 14.5	<input type="checkbox"/> 15
<input type="checkbox"/> 9.5		<input type="checkbox"/> 10	<input type="checkbox"/> 10.5	<input type="checkbox"/> 11	<input type="checkbox"/> 11.5	<input type="checkbox"/> 12	<input type="checkbox"/> 12.5	<input type="checkbox"/> 13	<input type="checkbox"/> 13.5	<input type="checkbox"/> 14	<input type="checkbox"/> 14.5	<input type="checkbox"/> 15.5	<input type="checkbox"/> 16	<input type="checkbox"/> 16.5	
PRO ELR Low Pivot Setting	<input type="checkbox"/> Short	12 - 17.5													
	<input type="checkbox"/> Long	15 - 21													
PRO ELR High Pivot Setting	<input type="checkbox"/> Short	10.5 - 15.5													
	<input type="checkbox"/> Long	13 - 19													
Fixed Front Frame (70°, 80°)	<input type="checkbox"/> Short	<input type="checkbox"/> 11.5	<input type="checkbox"/> 12	<input type="checkbox"/> 12.5	<input type="checkbox"/> 13	<input type="checkbox"/> 13.5	<input type="checkbox"/> 14	<input type="checkbox"/> 14.5	<input type="checkbox"/> 15						
	<input type="checkbox"/> Medium	<input type="checkbox"/> 13	<input type="checkbox"/> 13.5	<input type="checkbox"/> 14	<input type="checkbox"/> 14.5	<input type="checkbox"/> 15	<input type="checkbox"/> 15.5	<input type="checkbox"/> 16	<input type="checkbox"/> 16.5	<input type="checkbox"/> 17					
	<input type="checkbox"/> Long	<input type="checkbox"/> 15.5	<input type="checkbox"/> 16	<input type="checkbox"/> 16.5	<input type="checkbox"/> 17	<input type="checkbox"/> 17.5	<input type="checkbox"/> 18	<input type="checkbox"/> 18.5	<input type="checkbox"/> 19	<input type="checkbox"/> 19.5					

<sup>1</sup> Not available with Titanium Frame. See Minimum Front Seat Height chart above.

**3" Extension**

Available only with Front Mount Hangers

- 3" Extension \_\_\_\_\_ N/C
- 3" Extension (Provided as Parts) \_\_\_\_\_ N/C


**Footplate**

- Composite \_\_\_\_\_ N/C  
*Not available with Front Mount Hanger or PRO ELR*
- Composite Angle Adjustable K0040 \_\_\_\_\_ \$190  
*Not available with Front Mount Hanger or 14" Seat Width*
- Aluminum Angle Adjustable K0040 \_\_\_\_\_ \$190
- Aluminum Locking Multi-Angle Adjustable K0108 \_\_\_\_\_ \$300  
*Available only with Front Mount Hanger*
- One Piece Flip Up Angle Adjustable K0040 \_\_\_\_\_ \$190  
*Not available with Front Mount Hanger or PRO ELR*

**Hanger Release**

If no selection is made Classic Hanger Release will be used

- Classic \_\_\_\_\_ Std
- 4-Way - In line \_\_\_\_\_ N/C  
*See reference A below. Not available with PRO ELR.*
- 4-Way - Out \_\_\_\_\_ N/C  
*See reference B below. Not available with PRO ELR.*



**Frame Type / Footrest continued**

- Heel Loop**
- Not available with One Piece Flip Up Angle Adjustable Footplate*
- Heel Loops E0951 \$40
  - Adjustable Heel Loops E0951 \$50
  - Omit Heel Loops N/C

- Calf Strap**
- Velcro Adjustable K0038 \$40
  - Padded Velcro Adjustable K0038 \$68
  - Bodypoint Padded AeroMesh K0038 \$118

**Front Seat Height and Casters**

**Caster Size / Front Frame Type / Front Seat Height**

Select Front Seat Height in row corresponding to desired caster size.

Note: Selections in first three Front Seat Height columns use Super Low Front Frame. Super Low Front Frame not available with Titanium Frame or Fixed Front Frame.

Caster Size	Front Seat Height													
	Uses Super Low Front Frame													
4"	<input type="checkbox"/> 13.5	<input type="checkbox"/> 14	<input type="checkbox"/> 14.5	<input type="checkbox"/> 15	<input type="checkbox"/> 15.5	<input type="checkbox"/> 16	<input type="checkbox"/> 16.5	<input type="checkbox"/> 17	<input type="checkbox"/> 17.5	<input type="checkbox"/> 18	<input type="checkbox"/> 18.5	<input type="checkbox"/> 19	<input type="checkbox"/> 19.5	<input type="checkbox"/> 20
5"	<input type="checkbox"/> 14.5	<input type="checkbox"/> 15	<input type="checkbox"/> 15.5	<input type="checkbox"/> 16	<input type="checkbox"/> 16.5	<input type="checkbox"/> 17	<input type="checkbox"/> 17.5	<input type="checkbox"/> 18	<input type="checkbox"/> 18.5	<input type="checkbox"/> 19	<input type="checkbox"/> 19.5	<input type="checkbox"/> 20	<input type="checkbox"/> 20.5	
6"	<input type="checkbox"/> 15.5	<input type="checkbox"/> 16	<input type="checkbox"/> 16.5	<input type="checkbox"/> 17	<input type="checkbox"/> 17.5	<input type="checkbox"/> 18	<input type="checkbox"/> 18.5	<input type="checkbox"/> 19	<input type="checkbox"/> 19.5	<input type="checkbox"/> 20	<input type="checkbox"/> 20.5	<input type="checkbox"/> 21		
7"	<input type="checkbox"/> 17	<input type="checkbox"/> 17.5		<input type="checkbox"/> 18	<input type="checkbox"/> 18.5	<input type="checkbox"/> 19	<input type="checkbox"/> 19.5	<input type="checkbox"/> 20	<input type="checkbox"/> 20.5	<input type="checkbox"/> 21				
8"	<input type="checkbox"/> 17.5	<input type="checkbox"/> 18	<input type="checkbox"/> 18.5	<input type="checkbox"/> 19	<input type="checkbox"/> 19.5	<input type="checkbox"/> 20	<input type="checkbox"/> 20.5	<input type="checkbox"/> 21						

**Caster**

Caster Type	Price	Caster Size				
		4"	5"	6"	7"	8"
.75" Rollerblade	N/C	<input type="checkbox"/>				
.75" Lighted Rollerblade	\$25	<input type="checkbox"/>				
1" Poly	N/C	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>
1" Poly Aluminum	\$110	<input type="checkbox"/>	<input type="checkbox"/>			
1" Pneumatic	\$85			<input type="checkbox"/>		<input type="checkbox"/>
1.5" Poly	N/C	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>		
1.5" Soft Roll Aluminum	\$65	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>		
2" Poly	N/C			<input type="checkbox"/>		
2" Pneumatic w/ Foam Insert	\$125					<input type="checkbox"/>

**Fork Stem**

- Best Fit Std  
*Ki Mobility will pick Fork and Stem combination to optimize performance based on Caster Size and Front Seat Height*
- Standard Fork Stem N/C
- 3/4" Longer Fork Stem N/C
- 1 1/2" Longer Fork Stem N/C

**Caster Housing**

- Standard N/C
- Performance \$35  
*Not available with 350 lb Weight Capacity*

**Residual Limb Support**

- Residual Limb Support E1020
- Right and Left \$810
  - Left \$405  Right \$405
  - Omit Left Hanger N/C  Omit Right Hanger N/C
- Cover**
- Smooth  Air Mesh  Smooth  Air Mesh
- Pad Size**
- |                                 |                                  |                                 |                                  |
|---------------------------------|----------------------------------|---------------------------------|----------------------------------|
| <input type="checkbox"/> 5 x 8  | <input type="checkbox"/> 9 x 14  | <input type="checkbox"/> 5 x 8  | <input type="checkbox"/> 9 x 14  |
| <input type="checkbox"/> 7 x 10 | <input type="checkbox"/> 11 x 10 | <input type="checkbox"/> 7 x 10 | <input type="checkbox"/> 11 x 10 |
| <input type="checkbox"/> 9 x 10 | <input type="checkbox"/> 11 x 14 | <input type="checkbox"/> 9 x 10 | <input type="checkbox"/> 11 x 14 |

**Maximum Seat Dump Reference**

(Front Seat Height minus Rear Seat Height)

Seat Depth	14"	15"	16"	17"	18"	19"	20"
Maximum Dump Standard Caster Housing	1.5"	2"	2"	2"	2.5"	2.5"	2.5"
Maximum Dump Performance Caster Housing	3"	3"	3"	3.5"	3.5"	4"	4"

Depending on configuration, chairs with Fixed Front Frame Type may achieve up to an additional .75"

**Fork Options**

- Frog Legs Suspension Fork \$450  
*Available only with 4", 5" or 6" Caster Sizes. Not available with Caster Pin Locks, 350 lb Weight Capacity, 2" Poly caster or 2" Pneumatic w/ Foam Insert caster.*
  - User Weight Under 80 lbs
  - User Weight 80 - 120 lbs
  - User Weight 120 - 160 lbs
  - User Weight 160 - 220 lbs
  - User Weight 220 - 260 lbs
  - User Weight 260 - 300 lbs
- Frog Legs II Suspension System \$450  
*Lightweight shock absorbing fork. Available only with 4", 5" or 6" Caster Sizes. Not available with Caster Pin Locks, 350 lb Weight Capacity, 2" Poly caster or 2" Pneumatic w/ Foam Insert caster.*
- Caster Pin Locks K0073 \$195  
*Not available with Performance Caster Housing, 1 1/2" Longer Fork Stem or Frog Legs Fork*

## Wheel Locks

<b>Wheel Lock</b>		Dual - Drum Brake _____ \$660
<input type="checkbox"/> Push to Lock _____ Std		<i>Available only with 20", 22" and 24" MAXX Performance Spoke wheels. Not available with Reclining Backrest, Straight without Push Handle or Straight with Fold Down Push Handle Back Posts, or Low Profile Push Handle. Axle is Quick Release.</i>
<input type="checkbox"/> Push to Lock Flush Mount _____ N/C		<input type="checkbox"/> Push to Lock <input type="checkbox"/> Pull to Lock
<input type="checkbox"/> Pull to Lock _____ N/C		
<input type="checkbox"/> Low Profile Scissor Lock _____ N/C		
<i>Not available with Low Poly, Full Poly or Sentinel Solid tires</i>		
<input type="checkbox"/> Short Thro Scissor Lock _____ \$110		<input type="checkbox"/> Omit Wheel Locks _____ N/C
<input type="checkbox"/> Drum Brake _____ \$560		
<i>Available only with 20", 22" and 24" MAXX Performance Spoke wheels. Not available with Reclining Backrest, Straight without Push Handle or Straight with Fold Down Push Handle Back Posts, or Low Profile Push Handle. Axle is Quick Release.</i>		

<b>Wheel Lock Options</b>	
<input type="checkbox"/> Extension Handles E0961 _____ \$75	<i>Available only with Push to Lock or Pull to Lock Wheel Locks</i>
<input type="checkbox"/> Grade Aids E0974 _____ \$150	<i>Available only with Push to Lock Wheel Locks and Pneumatic or Pneumatic with Airless Insert tires</i>

## Rear Seat Height and Wheels

**Wheel Size / Tire Profile / Rear Seat Height**

Select **Rear Seat Height** in row corresponding to **Wheel Size and Tire Profile**. 14" Seat Depth not available with 24", 25" or 26" wheel sizes.

Requires Universal Axle Plate     Any axle plate can be used (Universal, Amputee, One Arm Drive)

Wheel Size	Tire Profile	Rear Seat Height									
20"	Low	<input type="radio"/> 11.5	<input type="checkbox"/> 12 <sup>1</sup>	<input type="radio"/> 12.5	<input type="checkbox"/> 13 <sup>1</sup>	<input type="radio"/> 13.5	<input type="checkbox"/> 14	<input type="radio"/> 14.5	<input type="checkbox"/> 15	<input type="radio"/> 15.5	<input type="radio"/> 16
	Full	<input type="radio"/> 12	<input type="checkbox"/> 12.5 <sup>1</sup>	<input type="radio"/> 13	<input type="checkbox"/> 13.5 <sup>1</sup>	<input type="radio"/> 14	<input type="checkbox"/> 14.5	<input type="radio"/> 15	<input type="checkbox"/> 15.5	<input type="radio"/> 16	<input type="radio"/> 16.5
22"	Low	<input type="radio"/> 12.5	<input type="checkbox"/> 13 <sup>1</sup>	<input type="radio"/> 13.5	<input type="checkbox"/> 14 <sup>1</sup>	<input type="radio"/> 14.5	<input type="checkbox"/> 15	<input type="radio"/> 15.5	<input type="checkbox"/> 16	<input type="radio"/> 16.5	<input type="radio"/> 17
	Full	<input type="radio"/> 13	<input type="checkbox"/> 13.5 <sup>1</sup>	<input type="radio"/> 14	<input type="checkbox"/> 14.5 <sup>1</sup>	<input type="radio"/> 15	<input type="checkbox"/> 15.5	<input type="radio"/> 16	<input type="checkbox"/> 16.5	<input type="radio"/> 17	<input type="radio"/> 17.5
24"	Low	<input type="radio"/> 13.5	<input type="checkbox"/> 14 <sup>1</sup>	<input type="radio"/> 14.5	<input type="checkbox"/> 15 <sup>1</sup>	<input type="radio"/> 15.5	<input type="checkbox"/> 16	<input type="radio"/> 16.5	<input type="checkbox"/> 17	<input type="radio"/> 17.5	<input type="radio"/> 18
	Full	<input type="radio"/> 14	<input type="checkbox"/> 14.5 <sup>1</sup>	<input type="radio"/> 15	<input type="checkbox"/> 15.5 <sup>1</sup>	<input type="radio"/> 16	<input type="checkbox"/> 16.5	<input type="radio"/> 17	<input type="checkbox"/> 17.5	<input type="radio"/> 18	<input type="radio"/> 18.5
25"	Low	<input type="radio"/> 14	<input type="checkbox"/> 14.5 <sup>1</sup>	<input type="radio"/> 15	<input type="checkbox"/> 15.5 <sup>1</sup>	<input type="radio"/> 16	<input type="checkbox"/> 16.5	<input type="radio"/> 17	<input type="checkbox"/> 17.5	<input type="radio"/> 18	<input type="radio"/> 18.5
26"	Low		<input type="checkbox"/> 15 <sup>1,2</sup>	<input type="radio"/> 15.5	<input type="checkbox"/> 16 <sup>1</sup>	<input type="radio"/> 16.5	<input type="checkbox"/> 17	<input type="radio"/> 17.5	<input type="checkbox"/> 18	<input type="radio"/> 18.5	<input type="radio"/> 19
	Full	<input type="radio"/> 15 <sup>2</sup>	<input type="checkbox"/> 15.5 <sup>1,2</sup>	<input type="radio"/> 16	<input type="checkbox"/> 16.5 <sup>1</sup>	<input type="radio"/> 17	<input type="checkbox"/> 17.5	<input type="radio"/> 18	<input type="checkbox"/> 18.5	<input type="radio"/> 19	<input type="radio"/> 19.5

<sup>1</sup> Not available with Amputee Axle Plate  
<sup>2</sup> Not available with Height Adjustable Flip Back T-Arm  
<sup>3</sup> Not available with Depth Adjustable Rear Frame Type, Reclining Rear Frame or Height Adjustable Flip Back T-Arms

**One Arm Drive**

One Arm Drive E0958 \_\_\_\_\_ \$1095

*Available only with 20", 22" and 24" spoke wheels and Quick Release Axles. Spoke wheel style is exclusive to OAD. Not available with 4 degrees of camber. Standard handrim will be Aluminum Anodized.*

Rear Wheel Size:     20"     22"     24"

Drive Handrim Side:     Left     Right

Drive Handrim Type:     Aluminum \_\_\_\_\_ N/C  
                                    Poly Coated \_\_\_\_\_ \$200

**Axle Plate**

*No selection needed with One Arm Drive*

Universal \_\_\_\_\_ Std

Amputee \_\_\_\_\_ \$175

*Not available with Swing Away Armrests or Wheelchair Drive Ready Option. Required with Reclining Backrest.*

**Rear Wheel**

*No selection needed with One Arm Drive*

		Rear Wheel Size					
		Price	20"	22"	24"	25"	26"
Rear Wheel Type	MAXX Performance Spoke	N/C	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>
	Mag	N/C	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>		
	Superlight Spoke <sup>1</sup>	\$295			<input type="checkbox"/>	<input type="checkbox"/> <sup>2</sup>	
	Spinergy Spox	\$845		<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>
	Spinergy LX	\$1050		<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>

<sup>1</sup> Not available with 350 lb Weight Capacity  
<sup>2</sup> Available only with Superlight and Plastic Coated Handrims

**Spinergy Options**

Select Spoke Color and Handrim Tab for Spinergy wheels

Spoke Color:

Black     Red     White     Blue     Yellow  
 Pink <sup>1</sup>     Green <sup>1</sup>     Purple <sup>1</sup>     Orange <sup>1</sup>

<sup>1</sup> Available only with Spinergy LX

Handrim Tab:

*Available only with Aluminum Anodized and Plastic Coated Handrims*

Short Tab \_\_\_\_\_ \$25

Long Tab \_\_\_\_\_ N/C

## Rear Seat Height and Wheels continued

### Handrim

**No selection needed with One Arm Drive**

#### Standard Handrims

- Aluminum Anodized \_\_\_\_\_ Std
- Superlight Aluminum Anodized \_\_\_\_\_ \$50  
*Available only on 24" and 25" wheel sizes. Available only on MAXX Performance, Superlight or Spinergy wheels. Not available with 350 lb Weight Capacity.*
- Plastic Coated \_\_\_\_\_ \$50
- Projection (8 Vertical Projections) \_\_\_\_\_ \$325  
*Available only with full profile tires in combination with 22"-26" MAXX Performance or 24" Mag wheel*
- Spinergy Flex Rim \_\_\_\_\_ \$700  
*Available only on 24" and 25" Spinergy LX wheels*
- Omit Handrim \_\_\_\_\_ N/C

#### Natural Fit Handrims

**Not available on 20" wheel size or 22" Mag wheel. Not recommended with Mag wheel.**

- Natural Fit LT Regular Grip \_\_\_\_\_ \$465
- Natural Fit Standard Regular Grip \_\_\_\_\_ \$445
- Natural Fit LT Super Grip \_\_\_\_\_ \$515
- Natural Fit Standard Super Grip \_\_\_\_\_ \$495
- Natural Fit LT Omit Grip \_\_\_\_\_ \$395
- Natural Fit Standard Omit Grip \_\_\_\_\_ \$345

#### Handrim Option

- 3M Micro Grip \_\_\_\_\_ \$25  
*Not available on Plastic Coated or Projection Handrims*

### Rear Tire Type

#### Full Profile

- Street Pneumatic 1 3/8" E2211, E2212 \_\_\_\_\_ Std  
*Not available on 25" wheel size*
- Full Poly 1 3/8" E2218 \_\_\_\_\_ N/C  
*Not available on Superlight wheel or 25" or 26" wheel sizes*
- Pneumatic with Airless Insert 1 3/8" E2211, E2213 \_\_\_\_\_ \$85  
*Not available on Superlight wheel or 25" wheel size*
- Schwalbe HS-127 1 3/8" E2211, E2212 \_\_\_\_\_ \$125  
*Available only on 24" wheel size. Not available on Mag wheel.*
- Euro Trek Knobby 1 3/8" E2211, E2212 \_\_\_\_\_ \$95  
*Available only on 24" wheel size*

#### Low Profile

- Low Poly 1" E2218 \_\_\_\_\_ N/C  
*Not available on Superlight wheel, 25" or 26" wheel sizes or 350 lb Weight Cap.*
- Sentinel Solid Tire 1" E2218 \_\_\_\_\_ \$25  
*Not available on Superlight wheels or 20" and 22" Mag wheels*
- High Pressure 1" E2211, E2212 \_\_\_\_\_ N/C  
*Not available on Mag wheel*
- Schwalbe One 1" E2211, E2212 \_\_\_\_\_ \$115  
*Not available on Mag wheel, 20", 22" or 26" wheel sizes*
- Schwalbe Marathon Plus Evolution 1" E2211, E2212 \_\_\_\_\_ \$125  
*Not available on Mag wheel, 20" or 22" wheel sizes*

### Axle

**No selection needed with One Arm Drive**

- Quick Release \_\_\_\_\_ Std
- Quad Quick Release \_\_\_\_\_ \$125  
*Not available with Mag or Spinergy wheels*

### Axle Position

- Rear Wheel Camber:  0°       2°       4°
- Center of Gravity:  Rear       Center       Forward
- No selection needed with Amputee Axle Plate*

### Armrests

#### Armrest

- Height Adjustable T-Arm (8.5" - 13.5") E0973
- Desk \_\_\_\_\_ \$295
  - Full \_\_\_\_\_ \$295
- Height Adjustable T-Arm Low (7.5" - 11.5") E0973  
*If selection is not available with rear wheel and rear seat height configuration, standard Height Adjustable T-Arm armrests will be used*
- Desk \_\_\_\_\_ \$295
  - Full \_\_\_\_\_ \$295
- Height Adjustable Flip Back T-Arm (9.5"-13.5") E0973  
*Available only with Standard Back Post*
- Desk \_\_\_\_\_ \$305
  - Full \_\_\_\_\_ \$305
- Flip Back (10.5")  
*Not available with Seat Depths less than 16". Available only with Standard Back Post.*
- Desk \_\_\_\_\_ N/C
  - Full \_\_\_\_\_ N/C

- Angle Adjustable Locking Extendable Flip Up E0973  
*Available only with Depth Adjustable and Reclining Back Posts*
- Desk \_\_\_\_\_ \$265
  - Full \_\_\_\_\_ \$265
- Swing Away (10.5" - 12.5") E0973 \_\_\_\_\_ \$225  
*Not available with Angle Adjustable Back Posts or Arm Pads*
  - Swing Away with Side Guard (10.5" - 12.5") E0973 \_\_\_\_\_ \$295  
*Not available with Angle Adjustable Back Posts or Arm Pads*
  - Omit Armrests \_\_\_\_\_ N/C

#### Arm Pad

**Not available with Swing Away Armrests**

- Standard Arm Pad \_\_\_\_\_ N/C
- Waterfall Arm Pad \_\_\_\_\_ \$25  
*Available only with Height Adjustable T-Arm and Height Adjustable T-Arm Low*
- Foam Grip \_\_\_\_\_ N/C  
*Available only with Angle Adjustable Locking Extendable Flip Up Armrest*

**Accessories**

- Removable Side Guard**
- Not available with T-Arms, Flip Back or Swing Away with Side Guard Armrests*
- Composite Side Guards K0108 \$140
  - Aluminum Straight K0108
    - Black Matte \$195
    - Match Frame Color \$245  
*Not available with Hydrographics or Titanium Frame*
  - Aluminum Adjustable Fender with Foam Pad K0108
    - Black Matte \$275
    - Match Frame Color \$325  
*Not available with Hydrographics or Titanium Frame*
  - Carbon Fiber Adjustable Fender K0108 \$485

- Pelvic Positioning Belt**
- 1 1/2" Auto Buckle E0978 \$80
  - 1 1/2" Padded Auto Buckle E0978 \$170
  - 2" Airline Buckle E0978 \$110
  - Bodypoint Evoflex - Pelvic Stabilizer E0978 \$188  
Size:  Small (16-24")  Medium (19-29")  Large (25-37")
  - Bodypoint Belt Mounting Kit \$108  
Style:  Bracket  Band Clamp 1"

- O<sub>2</sub> Holder**
- Not available with Depth Adjustable Rear Frame*
- O<sub>2</sub> Holder E2208 \$295  
 Right Side  Left Side

**Color**

- Frame Color not available with Titanium Frame Option*
- Gloss Frame Finish**
- Black N/C
  - Blue Skies N/C
  - Bubblegum Pink N/C
  - Candy Blue N/C
  - Candy Purple N/C
  - Candy Red N/C
  - Copper N/C
  - Granny Smith N/C
  - Metallic Orange N/C
  - Yellow N/C

- Matte Frame Finish**
- Charcoal Gray N/C
  - Emerald Green N/C
  - Iridescent Burgundy N/C
  - Iridescent Navy N/C
  - Military Green N/C
  - Retro Red N/C
  - Royal Blue N/C
  - Shadow Black N/C
  - Smooth White N/C

- Anti-Tip**
- Rear Anti-Tips E0971 \$145
  - Swing-In Anti-Tip E0971  
*Not available with 350 lb Weight Capacity or Reclining Back Post*
    - Single - Left \$155  
*Not available with left side O<sub>2</sub> Holder*
    - Single - Right \$155  
*Not available with right side O<sub>2</sub> Holder*
    - Pair \$310  
*Not available with Cane and Crutch Holder or O<sub>2</sub> Holder*
  - Omit Rear Anti-Tips N/C

- Other Accessories**
- Tool Kit \$45
  - Removable Underseat Pouch \$50
  - Pediatric Ki Mobility Backpack \$50
  - Ki Mobility Backpack \$65
  - Luggage Carrier \$85  
*Not available with PRO ELR or Residual Limb Support*
  - Black Neoprene Impact Guards \$90
  - Clear Spoke Guards K0065 \$115  
*Not available with Mag wheel*
  - Black Spoke Guards K0065 \$115  
*Not available with 20" wheel size or Mag wheel*
  - Cane and Crutch Holder \$195
  - IV Pole \$275  
*Available only with Curved Rear Frame Type and Straight Back Posts*

- Hydrographics**
- Hydrographics \$695
    - AOR1 (Desert Digital Camo)
    - Pink Camo
    - Muddy Girl
    - NEXT Camo
    - Carbon Fiber<sup>1</sup>
    - Zebra<sup>1</sup>
- <sup>1</sup>Select Hydrographics Base Color. Zebra not available in Black.

- Hydrographics Base Color N/C
- Black  Blue  Granny Smith  Pink  Purple  Red

**Colored Anodizing Package**

Select desired color. If no selection is made, black will be used.

	Option	Standard Caster Housing <sup>2</sup>	Caster Fork <sup>3</sup>	Axle Plate <sup>4</sup>
	Retail	\$75	\$75	\$75
Color	Black <sup>1</sup>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>
	Blue	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>
	Bright Green	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>
	Gold	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>
	Orange	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>
	Pewter	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>
	Pink	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>
	Purple	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>
Red	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	

<sup>1</sup>Black is no charge  
<sup>2</sup>Available only if corresponding option selection was made  
<sup>3</sup>Not available with 2" wide casters, Frog Legs Suspension Fork or Frog Legs II Suspension System  
<sup>4</sup>Color Anodized not available with amputee axle plate

Catalyst Front Seat Height Matrix																																	
4" Caster				5" Caster				6" Caster				7" Caster				8" Caster																	
FSH	Super Low	Stem	Fork	FSH	Super Low	Stem	Fork	FSH	Super Low	Stem	Fork	FSH	Super Low	Stem	Fork	FSH	Super Low	Stem	Fork														
15	13.5	Std	4"	16	14.5	Std	4"	17	15.5	Std	5"	18	16.5	Std	6"	19	17.5	Std	6"														
15.5	14			16.5	15.5	3/4		17.5	16.5	3/4		18.5	17			19.5	18.5	3/4															
15.5 <sup>1</sup>	14 <sup>1</sup>			17.5	16	1 1/2		18.5	17	1 1/2		18.5	17.5			20.5	n/a	1 1/2															
16	14.5	3/4		5"	16	14.5	Std	5"	17	15.5	Std	6"	19	18		3/4	7"	19	17.5	Std	7"												
16.5 <sup>1</sup>	15 <sup>1</sup>	1 1/2			16.5	15			17.5	16			18	16.5				19.5	18			19.5	18										
17	15.5	3/4			16.5	15.5			17.5	16.5			18	16.5				19.5	18.5			20	18.5										
15.5	14	Std			5"	17	16		3/4	5"	17.5		16.5	3/4		6"		18.5	17	Std		7"	20.5	19	Std	7"							
16	14.5					1 1/2	17.5				16		18					17	18.5				17.5	19			17.5	19.5	18				
16	15					3/4	18				16.5		18.5					17.5	19.5				18	19.5			18	20	19				
16.5	15.5	1 1/2				5"	16.5		15		Std		5"	18.5				17	1 1/2	6"			20	18.5	3/4		7"	20.5	n/a	3/4	7"		
17	15.5		17				15.5		19					17.5	19			18					21	n/a									
17.5	16		17.5				16		19.5					18	19.5			18.5					20.5	19				21	n/a			1 1/2	
16.5	15	Std	6"				17.5		16		3/4			6"	18			16.5	Std				7"	20	19			1 1/2	7"				
17	15.5			18			16.5	18.5	17			19			17.5		20.5	n/a															
17	16			18.5			17	19.5	18			20			18.5		20.5	19															
17.5	16.5	3/4		6"			18	16.5	1 1/2		6"	18.5			17		3/4	7"	21		n/a			1 1/2	7"								
18	16.5				18.5		17	19.5		18		20			18.5																		
18.5	17				19		17.5	19		18		20.5			19																		
18	16.5	Std			7"		18	16.5	Std	7"		19.5			18.5	1 1/2	7"																
18.5	17					18.5	17	19.5				18	20		19																		
18.5	17.5					19	17.5	20				18.5	20.5		19																		
19	18	3/4				7"	18.5	17.5	3/4			7"	20.5		19	1 1/2								7"									
19.5	18.5		19				18	20.5					19																				
20	19		19.5				18.5	21					n/a																				
							19.5	18	1 1/2				6" x 2" caster not available with 5" fork																				
				20			18.5																										
				20.5			19																										

<sup>1</sup> Not available with 4 x 1.5" Soft Roll Aluminum

Seat heights with "FSH" heading available with Standard Front Frame. Seat heights with "Super Low" heading available with Super Low Front Frame. Actual seat height may vary depending upon configuration.