



Phone (800) 981-1540 | Fax (715) 254-0996
 sales@kimobility.com | www.kimobility.com

Company: _____

Ship To: _____

Account Number: _____ Date: _____

Address: _____

Quote PO Number: _____

City: _____ State: _____

RTS: _____

Zip: _____ Phone: _____

Marked For: _____

Email/Fax: _____

Select if chair quote has: Additional Ki Products with Different Quote/Order Number **Quote/Order Number:** _____

Catalyst 5Vx K0005 _____ **\$1915**
 300 lb weight capacity

Effective 01/17/2022

HCPCS codes provided are based on code verification by the PDAC or our interpretation of the code definitions. Proper use of HCPCS codes does not ensure coverage or payment. For coverage information, verify the appropriate payer's coverage policy.

Frame

Seat Width

- 14" _____ N/C
- 15" _____ N/C
- 16" _____ N/C
- 17" _____ N/C
- 18" _____ N/C
- 19" _____ N/C
- 20" E2201 _____ \$500
- 21" E2201 _____ \$500
- 22" E2201 _____ \$500

Seat Depth

- 14" _____ N/C
- 15" _____ N/C
- 16" _____ N/C
- 17" _____ N/C
- 18" _____ N/C
- 19" _____ N/C
- 20" E2203 _____ \$550

Weight Capacity

- 300 lb _____ Std
- 350 lb (Heavy Duty) _____ \$500
Requires minimum 6" caster height

Backrest

Back Post / Height Range / Back Height

| Back Post | Price | Height Range | Back Height | | | |
|--|-------|---------------------------------|-------------------------------|-------------------------------|-------------------------------|-------------------------------|
| | | | <input type="checkbox"/> 8.5 | <input type="checkbox"/> 9.5 | <input type="checkbox"/> 10.5 | <input type="checkbox"/> 11.5 |
| <input type="checkbox"/> Standard | N/C | <input type="checkbox"/> Short | <input type="checkbox"/> 12.5 | <input type="checkbox"/> 13.5 | <input type="checkbox"/> 14.5 | <input type="checkbox"/> 15.5 |
| | | <input type="checkbox"/> Medium | <input type="checkbox"/> 16.5 | <input type="checkbox"/> 17.5 | <input type="checkbox"/> 18.5 | <input type="checkbox"/> 19.5 |
| | | <input type="checkbox"/> Tall | | | | |
| <input type="checkbox"/> Angle Adjustable w/ Push Handle | \$185 | <input type="checkbox"/> Short | <input type="checkbox"/> 14 | <input type="checkbox"/> 15 | <input type="checkbox"/> 16 | |
| | | <input type="checkbox"/> Medium | <input type="checkbox"/> 16 | <input type="checkbox"/> 17 | <input type="checkbox"/> 18 | |
| | | <input type="checkbox"/> Tall | <input type="checkbox"/> 18 | <input type="checkbox"/> 19 | <input type="checkbox"/> 20 | |

Transit

- Transit System _____ \$320
Not available with Removable Seat Pan or 350 lb Weight Capacity
- Omit Transit _____ N/C

Cross Brace

- Standard Cross Braces _____ Std
- Seating System Cross Braces _____ N/C
Allows for easy attachment of seating system mounting hardware. Not available with Seat Sling.

Seat Option

- Removable Aluminum Seat Pan E2231 _____ \$275
Not available with Transit System. Pan depth 1/2" shorter than Seat Depth.

Wheelchair Drive Ready

Chair will ship from Ki Mobility with adaption kit installed for selected Alber system. Order system separately from Alber USA. Select Rear Seat Height (full profile selections only) in drive wheel size on Rear Seat Height chart (page 5). Requires Rear Anti-Tips. Not available with Height Adjustable Flip Back T-Arm.

- e-motion Ready _____ N/C
- twion Ready _____ N/C

Standard Back Post

Select option below only if a Standard Back Post option is selected on chart to the left

- Straight with Push Handle _____ Std
- Straight without Push Handle _____ N/C
- 8° Bend with Push Handle _____ \$25
Not available with Short Back Height Range
- Straight with Fold Down Push Handle _____ \$125
Not available with 350 lb Weight Capacity

Push Handle Options

Not available with Straight Back Post with Fold Down Push Handle or 350 lb Weight Capacity

- Low Profile Push Handle _____ \$85
2" shorter than standard push handles

Backrest continued

Back Upholstery

If selecting both upholstery and Axiom Back, the upholstery will be installed.

- Standard _____ Std
Not available with Short or Medium Back Height Ranges on 21" - 22" Seat Widths
- Tension Adjustable E2611
 - Standard Cover _____ \$315
 - Air Mesh Breathable Cover _____ \$390
- Omit Back Upholstery _____ - \$15

Military Patch

Attached only if ordered with Tension Adjustable Back Upholstery

- US Military Service Patch _____ \$90
 - Army Navy Marines Air Force Coast Guard
- US Purple Heart Patch _____ \$90
- Purple Heart with Military Service Patch Combo _____ \$145
 - Army Navy Marines Air Force Coast Guard

Footrest / Frame Type

Frame Type

- Swing Away _____ N/C
- Super Low - Swing Away _____ N/C
- Fixed Front _____ N/C

| Minimum Front Seat Height with Front Mount Hangers | | | | | | |
|--|-------|--------|--------|--------|-------|-----|
| Hanger Length | Short | | Medium | | Long | |
| | 70° | 80/90° | 70° | 80/90° | 70° | 80° |
| Minimum Front Seat Height with Standard Selection | 13.5" | 13.5" | 14" | 15" | 16" | 17" |
| Minimum Front Seat Height with 3" Extension | 14.5" | 16" | 16.5" | 18" | 18.5" | 20" |

Hanger / Front Frame Bend

| Frame | Hanger | 60° | 70° | 80° | 90° | Retail |
|------------|----------------------|--------------------------|--------------------------|---------------------------------------|---------------------------------------|--------|
| Swing Away | Extension Tube Mount | <input type="checkbox"/> | <input type="checkbox"/> | <input type="checkbox"/> ¹ | | Std |
| | Front Mount | | <input type="checkbox"/> | <input type="checkbox"/> ¹ | <input type="checkbox"/> ¹ | N/C |
| | PRO ELR | | | | <input type="checkbox"/> | \$295 |
| | Omit Hangers | | | | <input type="checkbox"/> | -\$20 |
| Fixed | Extension Tube Mount | | <input type="checkbox"/> | <input type="checkbox"/> ¹ | | N/C |

¹ Not available with 7" or 8" casters

PRO ELR Cover

- Available only on PRO ELR*
- Gel Pad Cover _____ \$390

Footrest / Extension / Seat to Footrest Length

Select Extension and Seat to Footrest Length in row corresponding to desired Footrest. Footrest Length of at least 2.5" shorter than Front Seat Height recommended. Extension Tube Hangers and Fixed Front Frame measurements are based on 70° Hanger / Front Frame Bend.

| Footrest | Extension | Seat to Footrest Length | | | | | | | | | | |
|---|--|-------------------------------|-------------------------------|-------------------------------|-------------------------------|-------------------------------|-------------------------------|-------------------------------|-------------------------------|-------------------------------|-------------------------------|-------------------------------|
| Extension Tube Hangers (60°, 70°, 80°) | <input type="checkbox"/> Short | <input type="checkbox"/> 11 | <input type="checkbox"/> 11.5 | <input type="checkbox"/> 12 | <input type="checkbox"/> 12.5 | <input type="checkbox"/> 13 | <input type="checkbox"/> 13.5 | <input type="checkbox"/> 14 | <input type="checkbox"/> 14.5 | <input type="checkbox"/> 15 | | |
| | <input type="checkbox"/> Medium | <input type="checkbox"/> 13 | <input type="checkbox"/> 13.5 | <input type="checkbox"/> 14 | <input type="checkbox"/> 14.5 | <input type="checkbox"/> 15 | <input type="checkbox"/> 15.5 | <input type="checkbox"/> 16 | <input type="checkbox"/> 16.5 | <input type="checkbox"/> 17 | | |
| | <input type="checkbox"/> Long | <input type="checkbox"/> 15.5 | <input type="checkbox"/> 16 | <input type="checkbox"/> 16.5 | <input type="checkbox"/> 17 | <input type="checkbox"/> 17.5 | <input type="checkbox"/> 18 | <input type="checkbox"/> 18.5 | <input type="checkbox"/> 19 | <input type="checkbox"/> 19.5 | | |
| Front Mount Hangers ¹ | Standard Selections | | | | | | | | | | | |
| | <input type="checkbox"/> Short (70°, 80°, 90°) | <input type="checkbox"/> 5 | <input type="checkbox"/> 5.5 | <input type="checkbox"/> 6 | <input type="checkbox"/> 6.5 | <input type="checkbox"/> 7 | | | | <input type="checkbox"/> 10 | <input type="checkbox"/> 10.5 | <input type="checkbox"/> 11 |
| | | <input type="checkbox"/> 7.5 | <input type="checkbox"/> 8 | <input type="checkbox"/> 8.5 | <input type="checkbox"/> 9 | <input type="checkbox"/> 9.5 | | | | <input type="checkbox"/> 11.5 | <input type="checkbox"/> 12 | <input type="checkbox"/> 12.5 |
| | <input type="checkbox"/> Medium (70°, 80°, 90°) | <input type="checkbox"/> 5 | <input type="checkbox"/> 5.5 | <input type="checkbox"/> 6 | <input type="checkbox"/> 6.5 | <input type="checkbox"/> 7 | <input type="checkbox"/> 7.5 | <input type="checkbox"/> 8 | | <input type="checkbox"/> 12 | <input type="checkbox"/> 12.5 | <input type="checkbox"/> 13 |
| | | <input type="checkbox"/> 8.5 | <input type="checkbox"/> 9 | <input type="checkbox"/> 9.5 | <input type="checkbox"/> 10 | <input type="checkbox"/> 10.5 | <input type="checkbox"/> 11 | <input type="checkbox"/> 11.5 | | <input type="checkbox"/> 13.5 | <input type="checkbox"/> 14 | <input type="checkbox"/> 14.5 |
| | <input type="checkbox"/> Long (70°, 80°) | <input type="checkbox"/> 5 | <input type="checkbox"/> 5.5 | <input type="checkbox"/> 6 | <input type="checkbox"/> 6.5 | <input type="checkbox"/> 7 | <input type="checkbox"/> 7.5 | <input type="checkbox"/> 8 | <input type="checkbox"/> 8.5 | <input type="checkbox"/> 9 | <input type="checkbox"/> 14 | <input type="checkbox"/> 14.5 |
| <input type="checkbox"/> 9.5 | | <input type="checkbox"/> 10 | <input type="checkbox"/> 10.5 | <input type="checkbox"/> 11 | <input type="checkbox"/> 11.5 | <input type="checkbox"/> 12 | <input type="checkbox"/> 12.5 | <input type="checkbox"/> 13 | <input type="checkbox"/> 13.5 | <input type="checkbox"/> 15.5 | <input type="checkbox"/> 16 | <input type="checkbox"/> 16.5 |
| PRO ELR Low Pivot Setting | <input type="checkbox"/> Short | 12 - 17.5 | | | | | | | | | | |
| | <input type="checkbox"/> Long | 15 - 21 | | | | | | | | | | |
| PRO ELR High Pivot Setting | <input type="checkbox"/> Short | 10.5 - 15.5 | | | | | | | | | | |
| | <input type="checkbox"/> Long | 13 - 19 | | | | | | | | | | |
| Fixed Front Frame (70°, 80°) | <input type="checkbox"/> Short | <input type="checkbox"/> 11.5 | <input type="checkbox"/> 12 | <input type="checkbox"/> 12.5 | <input type="checkbox"/> 13 | <input type="checkbox"/> 13.5 | <input type="checkbox"/> 14 | <input type="checkbox"/> 14.5 | <input type="checkbox"/> 15 | | | |
| | <input type="checkbox"/> Medium | <input type="checkbox"/> 13 | <input type="checkbox"/> 13.5 | <input type="checkbox"/> 14 | <input type="checkbox"/> 14.5 | <input type="checkbox"/> 15 | <input type="checkbox"/> 15.5 | <input type="checkbox"/> 16 | <input type="checkbox"/> 16.5 | <input type="checkbox"/> 17 | | |
| | <input type="checkbox"/> Long | <input type="checkbox"/> 15.5 | <input type="checkbox"/> 16 | <input type="checkbox"/> 16.5 | <input type="checkbox"/> 17 | <input type="checkbox"/> 17.5 | <input type="checkbox"/> 18 | <input type="checkbox"/> 18.5 | <input type="checkbox"/> 19 | <input type="checkbox"/> 19.5 | | |

¹ See Minimum Front Seat Height chart above.

Footrest / Frame Type continued

3" Extension

Available only with Front Mount Hangers

- 3" Extension \$35
- 3" Extension (Provided as Parts) \$35


Footplate

- Composite N/C
Not available with Front Mount Hanger or PRO ELR
- Composite Angle Adjustable K0040 \$190
Not available with Front Mount Hanger or 14" Seat Width
- Aluminum Angle Adjustable K0040 \$190
- One Piece Flip Up Angle Adjustable K0040 \$190
Not available with Front Mount Hanger or PRO ELR

Hanger Release

If no selection is made Classic Hanger Release will be used

- Classic Std
- 4-Way - In line N/C
See reference A to the right. Not available with PRO ELR.
- 4-Way - Out N/C
See reference B to the right. Not available with PRO ELR.



Heel Loop

Not available with One Piece Flip Up Angle Adjustable Footplate

- Heel Loops E0951 \$40
- Adjustable Heel Loops E0951 \$50
- Omit Heel Loops N/C

Calf Strap

- Velcro Adjustable K0038 \$40
- Padded Velcro Adjustable K0038 \$68
- Bodypoint Padded AeroMesh K0038 \$118

Residual Limb Support

- Residual Limb Support E1020
- Right and Left \$810
 - Left \$405 Right \$405
 - Omit Left Hanger -\$10 Omit Right Hanger -\$10

Cover

- Smooth Air Mesh Smooth Air Mesh

Pad Size

- | | | | |
|---------------------------------|----------------------------------|---------------------------------|----------------------------------|
| <input type="checkbox"/> 5 x 8 | <input type="checkbox"/> 9 x 14 | <input type="checkbox"/> 5 x 8 | <input type="checkbox"/> 9 x 14 |
| <input type="checkbox"/> 7 x 10 | <input type="checkbox"/> 11 x 10 | <input type="checkbox"/> 7 x 10 | <input type="checkbox"/> 11 x 10 |
| <input type="checkbox"/> 9 x 10 | <input type="checkbox"/> 11 x 14 | <input type="checkbox"/> 9 x 10 | <input type="checkbox"/> 11 x 14 |

Front Seat Height and Casters

Caster Size / Front Frame Type / Front Seat Height

Select Front Seat Height in row corresponding to desired caster size.

Note: Selections in first three Front Seat Height columns use Super Low Front Frame. Super Low Front Frame not available with Fixed Front Frame.

| Caster Size | Front Seat Height | | | | | | | | | | | | | | | | |
|-------------|----------------------------|--------------------------|--------------------------|--------------------------|--------------------------|--------------------------|--------------------------|--------------------------|--------------------------|--------------------------|--------------------------|--------------------------|--------------------------|--------------------------|--------------------------|--------------------------|--------------------------|
| | Uses Super Low Front Frame | | | 13.5 | 14 | 14.5 | 15 | 15.5 | 16 | 16.5 | 17 | 17.5 | 18 | 18.5 | 19 | 19.5 | 20 |
| 4" | <input type="checkbox"/> | <input type="checkbox"/> | <input type="checkbox"/> | <input type="checkbox"/> | <input type="checkbox"/> | <input type="checkbox"/> | <input type="checkbox"/> | <input type="checkbox"/> | <input type="checkbox"/> | <input type="checkbox"/> | <input type="checkbox"/> | <input type="checkbox"/> | <input type="checkbox"/> | <input type="checkbox"/> | <input type="checkbox"/> | <input type="checkbox"/> | <input type="checkbox"/> |
| 5" | <input type="checkbox"/> | <input type="checkbox"/> | <input type="checkbox"/> | <input type="checkbox"/> | <input type="checkbox"/> | <input type="checkbox"/> | <input type="checkbox"/> | <input type="checkbox"/> | <input type="checkbox"/> | <input type="checkbox"/> | <input type="checkbox"/> | <input type="checkbox"/> | <input type="checkbox"/> | <input type="checkbox"/> | <input type="checkbox"/> | <input type="checkbox"/> | <input type="checkbox"/> |
| 6" | <input type="checkbox"/> | <input type="checkbox"/> | <input type="checkbox"/> | <input type="checkbox"/> | <input type="checkbox"/> | <input type="checkbox"/> | <input type="checkbox"/> | <input type="checkbox"/> | <input type="checkbox"/> | <input type="checkbox"/> | <input type="checkbox"/> | <input type="checkbox"/> | <input type="checkbox"/> | <input type="checkbox"/> | <input type="checkbox"/> | <input type="checkbox"/> | <input type="checkbox"/> |
| 7" | <input type="checkbox"/> | <input type="checkbox"/> | <input type="checkbox"/> | <input type="checkbox"/> | <input type="checkbox"/> | <input type="checkbox"/> | <input type="checkbox"/> | <input type="checkbox"/> | <input type="checkbox"/> | <input type="checkbox"/> | <input type="checkbox"/> | <input type="checkbox"/> | <input type="checkbox"/> | <input type="checkbox"/> | <input type="checkbox"/> | <input type="checkbox"/> | <input type="checkbox"/> |
| 8" | <input type="checkbox"/> | <input type="checkbox"/> | <input type="checkbox"/> | <input type="checkbox"/> | <input type="checkbox"/> | <input type="checkbox"/> | <input type="checkbox"/> | <input type="checkbox"/> | <input type="checkbox"/> | <input type="checkbox"/> | <input type="checkbox"/> | <input type="checkbox"/> | <input type="checkbox"/> | <input type="checkbox"/> | <input type="checkbox"/> | <input type="checkbox"/> | <input type="checkbox"/> |

Caster

| Caster Type | | Price | Caster Size | | | | |
|-----------------------------|---------------------------|-------|--------------------------|--------------------------|--------------------------|--------------------------|--------------------------|
| | | | 4" | 5" | 6" | 7" | 8" |
| Caster Type | .75" Roller Blade | N/C | <input type="checkbox"/> | | | | |
| | .75" Lighted Roller Blade | \$70 | <input type="checkbox"/> | | | | |
| | 1" Poly | N/C | <input type="checkbox"/> | <input type="checkbox"/> | <input type="checkbox"/> | <input type="checkbox"/> | <input type="checkbox"/> |
| | 1" Poly Aluminum | \$160 | <input type="checkbox"/> | <input type="checkbox"/> | | | |
| | 1" Pneumatic | \$85 | | | <input type="checkbox"/> | | <input type="checkbox"/> |
| | 1.5" Poly | \$50 | <input type="checkbox"/> | <input type="checkbox"/> | <input type="checkbox"/> | | |
| | 1.5" Soft Roll Aluminum | \$185 | <input type="checkbox"/> | <input type="checkbox"/> | <input type="checkbox"/> | | |
| | 2" Poly | \$70 | | | <input type="checkbox"/> | | |
| 2" Pneumatic w/ Foam Insert | \$225 | | | | | <input type="checkbox"/> | |

Maximum Seat Dump Reference

| Seat Depth | 14" | 15" | 16" | 17" | 18" | 19" | 20" |
|---|------|-----|-----|------|------|------|------|
| Maximum Dump Standard Caster Housing | 1.5" | 2" | 2" | 2" | 2.5" | 2.5" | 2.5" |
| Maximum Dump Performance Caster Housing | 3" | 3" | 3" | 3.5" | 3.5" | 4" | 4" |

Depending on configuration, chairs with Fixed Front Frame Type may achieve up to an additional .75"

Front Seat Height and Casters continued

Fork Stem

- Best Fit _____ Std
Ki Mobility will pick Fork and Stem combination to optimize performance based on Caster Size and Front Seat Height
- Standard Fork Stem _____ N/C
- 3/4" Longer Fork Stem _____ N/C
- 1 1/2" Longer Fork Stem _____ N/C

Caster Housing

- Standard _____ N/C
- Performance _____ \$40
Not available with 350 lb Weight Capacity

Fork Options

- Frog Legs Suspension Fork _____ \$450
Available only with 4", 5" or 6" Caster Sizes. Not available with Caster Pin Locks, 350 lb Weight Capacity, 2" Poly caster or 2" Pneumatic w/ Foam Insert caster.
 - User Weight Under 80 lbs
 - User Weight 80 - 120 lbs
 - User Weight 120 - 160 lbs
 - User Weight 160 - 220 lbs
 - User Weight 220 - 260 lbs
 - User Weight 260 - 300 lbs
- Frog Legs II Suspension System _____ \$450
Lightweight shock absorbing fork. Available only with 4", 5" or 6" Caster Sizes. Not available with Caster Pin Locks, 350 lb Weight Capacity, 2" Poly caster or 2" Pneumatic w/ Foam Insert caster.
- Caster Pin Locks K0073 _____ \$195
Not available with Performance Caster Housing, 1 1/2" Longer Fork Stem or Frog Legs Fork

Wheel Locks

Wheel Lock

- Push to Lock _____ Std
- Push to Lock Flush Mount _____ N/C
- Pull to Lock _____ N/C
- Low Profile Scissor Lock _____ \$65
Not available with Low Poly, Full Poly or Sentinel Solid tires
- Short Thro Scissor Lock _____ \$110
- Drum Brake _____ \$560
Available only with 20", 22" and 24" MAXX Performance Spoke wheels. Not available with Straight without Push Handle or Straight with Fold Down Push Handle Back Posts, or Low Profile Push Handle. Axle is Quick Release.

- Dual - Drum Brake _____ \$660
Available only with 20", 22" and 24" MAXX Performance Spoke wheels. Not available with Straight without Push Handle or Straight with Fold Down Push Handle Back Posts, or Low Profile Push Handle. Axle is Quick Release.
 - Push to Lock
 - Pull to Lock

- Omit Wheel Locks _____ N/C

Wheel Lock Options

- Extension Handles E0961 _____ \$75
Available only with Push to Lock or Pull to Lock Wheel Locks
- Grade Aids E0974 _____ \$150
Available only with Push to Lock Wheel Locks and Pneumatic or Pneumatic with Airless Insert tires

Rear Seat Height and Wheels

Wheel Size / Tire Profile / Rear Seat Height

Select Rear Seat Height in row corresponding to desired Wheel Size and Tire Profile. 14" Seat Depth not available with 24", 25" or 26" wheel sizes.

| Wheel Size | Tire Profile | Rear Seat Height | | | | | |
|------------|--------------|--------------------------------|--------------------------------|--------------------------------|--------------------------------|--------------------------------|--------------------------------|
| 20" | Low | <input type="checkbox"/> 11.5" | <input type="checkbox"/> 12.5" | <input type="checkbox"/> 13.5" | <input type="checkbox"/> 14.5" | <input type="checkbox"/> 15.5" | <input type="checkbox"/> 16.5" |
| | Full | <input type="checkbox"/> 12" | <input type="checkbox"/> 13" | <input type="checkbox"/> 14" | <input type="checkbox"/> 15" | <input type="checkbox"/> 16" | <input type="checkbox"/> 17" |
| 22" | Low | <input type="checkbox"/> 12.5" | <input type="checkbox"/> 13.5" | <input type="checkbox"/> 14.5" | <input type="checkbox"/> 15.5" | <input type="checkbox"/> 16.5" | <input type="checkbox"/> 17.5" |
| | Full | <input type="checkbox"/> 13" | <input type="checkbox"/> 14" | <input type="checkbox"/> 15" | <input type="checkbox"/> 16" | <input type="checkbox"/> 17" | <input type="checkbox"/> 18" |
| 24" | Low | <input type="checkbox"/> 13.5" | <input type="checkbox"/> 14.5" | <input type="checkbox"/> 15.5" | <input type="checkbox"/> 16.5" | <input type="checkbox"/> 17.5" | <input type="checkbox"/> 18.5" |
| | Full | <input type="checkbox"/> 14" | <input type="checkbox"/> 15" | <input type="checkbox"/> 16" | <input type="checkbox"/> 17" | <input type="checkbox"/> 18" | <input type="checkbox"/> 19" |
| 25" | Low | <input type="checkbox"/> 14" | <input type="checkbox"/> 15" | <input type="checkbox"/> 16" | <input type="checkbox"/> 17" | <input type="checkbox"/> 18" | <input type="checkbox"/> 19" |
| 26" | Low | <input type="checkbox"/> 14.5" | <input type="checkbox"/> 15.5" | <input type="checkbox"/> 16.5" | <input type="checkbox"/> 17.5" | <input type="checkbox"/> 18.5" | <input type="checkbox"/> 19.5" |
| | Full | <input type="checkbox"/> 15" | <input type="checkbox"/> 16" | <input type="checkbox"/> 17" | <input type="checkbox"/> 18" | <input type="checkbox"/> 19" | <input type="checkbox"/> 20" |

Rear Seat Height and Wheels continued

Rear Wheel

| Rear Wheel Type | Price | Rear Wheel Size | | | | |
|--|-------|--------------------------|--------------------------|--------------------------|---------------------------------------|--------------------------|
| | | 20" | 22" | 24" | 25" | 26" |
| Assembled Six-Spoke Mag w/ Treaded Tire and Composite Handrim ^{1,2} | -\$25 | | | <input type="checkbox"/> | | |
| MAXX Performance Spoke | \$200 | <input type="checkbox"/> | <input type="checkbox"/> | | <input type="checkbox"/> | <input type="checkbox"/> |
| MAXX Performance Spoke | N/C | | | <input type="checkbox"/> | | |
| Mag | N/C | <input type="checkbox"/> | <input type="checkbox"/> | <input type="checkbox"/> | | |
| Superlight Spoke ² | \$295 | | | <input type="checkbox"/> | <input type="checkbox"/> ³ | |
| Spinergy Spox | \$845 | | <input type="checkbox"/> | <input type="checkbox"/> | <input type="checkbox"/> | <input type="checkbox"/> |

¹ No substitutions available. Do not select Rear Tire Type or Handrim. Measures as Full Profile tire.

² Not available with 350 lb Weight Capacity

³ Available only with Superlight and Plastic Coated Handrims

Rear Tire Type

Full Profile

- Street Pneumatic 1 3/8" E2211, E2212 _____ Std
Not available on 25" wheel size
- Full Poly 1 3/8" E2218 _____ N/C
Not available on Superlight wheel or 25" or 26" wheel sizes
- Pneumatic with Airless Insert 1 3/8" E2211, E2213 _____ \$100
Not available on Superlight wheel or 25" wheel size
- Schwalbe HS-127 1 3/8" E2211, E2212 _____ \$125
Available only on 24" wheel size. Not available on Mag wheel.
- Euro Trek Knobby 1 3/8" E2211, E2212 _____ \$95
Available only on 24" wheel size

Low Profile

- Low Poly 1" E2218 _____ \$50
Not available on Superlight wheel, 25" or 26" wheel sizes or 350 lb Weight Cap.
- Sentinel Solid Tire 1" E2218 _____ \$95
Not available on Superlight wheel or 20" and 22" Mag wheels
- High Pressure 1" E2211, E2212 _____ \$95
Not available on Mag wheel
- Schwalbe Marathon Plus Evolution 1" E2211, E2212 _____ \$135
Not available on Mag wheel, 20" or 22" wheel sizes

Axles

Axle

No selection needed with Drum Brake

- Quick Release _____ Std
- Quad Quick Release _____ \$175
Not available with Mag or Spinergy wheels

Camber

Combination of Standard and Adjustable Camber cannot exceed 4°. Selection will default to a max of 4° combined if exceeded.

- Standard Camber _____ N/C
Chair will be set up with this camber selection
- 0° 1° 2° 3° 4°

- Adjustable Camber _____ N/C
Extra camber spacers mounted to the inside of the rear frame for later use
- 1° 2° 3° 4°

Spinergy Options

Select spoke color and handrim tab for Spinergy wheels

- Spoke Color:
 Black Red White Blue Yellow

- Handrim Tab:
 Available only with Aluminum Anodized and Plastic Coated Handrims
 Short Tab _____ \$25
 Long Tab _____ N/C

Handrim

Standard Handrims

- Aluminum Anodized _____ Std
- Superlight Aluminum Anodized _____ \$50
Available only on 24" and 25" wheel sizes. Available only on MAXX Performance or Superlight wheels. Not available with 350 lb Weight Capacity.
- Plastic Coated _____ \$100
- Projection (8 Vertical Projections) _____ \$325
Available only with full profile tires in combination with 22"-26" MAXX Performance Spoke, or 24" Mag wheel
- Omit Handrim _____ N/C

Natural Fit Handrims

Not available on 20" wheel size or 22" Mag wheel. Not recommended with Mag wheel.

- Natural Fit LT Regular Grip _____ \$465
- Natural Fit Standard Regular Grip _____ \$445
- Natural Fit LT Super Grip _____ \$515
- Natural Fit Standard Super Grip _____ \$495
- Natural Fit LT Omit Grip _____ \$395
- Natural Fit Standard Omit Grip _____ \$345

Handrim Option


- 3M Micro Grip _____ \$65
Not available on Plastic Coated or Projection Handrims

Axle Position

If no selection is made position 3 will be used. Position 4 is 7/8" forward of the back cane.

Center of Gravity

- 1
- 2
- 3
- 4



Armrests

| Armrest | |
|---|-------|
| Height Adjustable T-Arm E0973 | |
| <input type="checkbox"/> Desk (8.5" - 13.5") | \$295 |
| <input type="checkbox"/> Full (8.5" - 13.5") | \$295 |
| Height Adjustable T-Arm Low E0973 | |
| <i>If selection is not available with rear wheel and rear seat height configuration, standard Height Adjustable T-Arm armrests will be used</i> | |
| <input type="checkbox"/> Desk (7.5" - 11.5") | \$295 |
| <input type="checkbox"/> Full (7.5" - 11.5") | \$295 |
| Height Adjustable Flip Back T-Arm E0973 | |
| <input type="checkbox"/> Desk (9.5"-13.5") | \$305 |
| <input type="checkbox"/> Full (9.5"-13.5") | \$305 |
| Flip Back Desk (10.5") | |
| <i>Not available with Seat Depths less than 16" or Angle Adjustable Back Posts</i> | |
| <input type="checkbox"/> Desk (10.5") | \$50 |
| <input type="checkbox"/> Full (10.5") | \$50 |
| <input type="checkbox"/> Swing Away (10.5" - 12.5") E0973 | \$225 |
| <i>Not available with Angle Adjustable Back Posts or Arm Pads</i> | |
| <input type="checkbox"/> Swing Away with Side Guard (10.5" - 12.5") E0973 | \$300 |
| <i>Not available with Angle Adjustable Back Posts or Arm Pads</i> | |
| <input type="checkbox"/> Omit Armrests | N/C |

| Arm Pad | |
|--|------|
| Not available with Swing Away Armrests | |
| <input type="checkbox"/> Waterfall Arm Pad | \$25 |
| <i>Available only with Height Adjustable T-Arm and Height Adjustable T-Arm Low</i> | |

Accessories

| Removable Side Guard | |
|--|-------|
| Not available with T-Arms, Flip Back or Swing Away with Side Guard Armrests | |
| <input type="checkbox"/> Composite Side Guards K0108 | \$145 |
| Aluminum Straight K0108 | |
| <input type="checkbox"/> Black Matte | \$195 |
| <input type="checkbox"/> Match Frame Color | \$245 |
| <i>Not available with Hydrographics</i> | |
| Fixed Mount Fender Side Guards K0108 | |
| <i>Not available with 24" wheel size in combination with 14" Rear Seat Height</i> | |
| <input type="checkbox"/> Aluminum | \$150 |
| <input type="checkbox"/> Carbon Fiber | \$475 |
| Aluminum Adjustable Fender with Foam Pad K0108 | |
| <input type="checkbox"/> Black Matte | \$275 |
| <input type="checkbox"/> Match Frame Color | \$325 |
| <i>Not available with Hydrographics</i> | |
| <input type="checkbox"/> Carbon Fiber Adjustable Fender K0108 | \$485 |

| Anti-Tip | |
|--|-------|
| <input type="checkbox"/> Rear Anti-Tips E0971 | \$150 |
| Swing-In Anti-Tip E0971 | |
| <i>Not available with 350 lb Weight Capacity</i> | |
| <input type="checkbox"/> Single - Left | \$155 |
| <i>Not available with left side O₂ Holder</i> | |
| <input type="checkbox"/> Single - Right | \$155 |
| <i>Not available with right side O₂ Holder</i> | |
| <input type="checkbox"/> Pair | \$310 |
| <i>Not available with Cane and Crutch Holder or O₂ Holder</i> | |
| <input type="checkbox"/> Omit Rear Anti-Tips | N/C |

| Pelvic Positioning Belt | |
|--|-------|
| <input type="checkbox"/> 1 1/2" Auto Buckle E0978 | \$80 |
| <input type="checkbox"/> 1 1/2" Padded Auto Buckle E0978 | \$170 |
| <input type="checkbox"/> 2" Airline Buckle E0978 | \$110 |
| Bodypoint Evoflex - Pelvic Stabilizer E0978 | \$188 |
| Size: <input type="checkbox"/> Small (16-24") <input type="checkbox"/> Medium (19-29") <input type="checkbox"/> Large (25-37") | |
| Bodypoint Belt Mounting Kit | \$108 |
| Style: <input type="checkbox"/> Bracket <input type="checkbox"/> Band Clamp 1" | |

| Other Accessories | |
|--|-------|
| <input type="checkbox"/> Tool Kit | \$45 |
| <input type="checkbox"/> Removable Underseat Pouch | \$50 |
| <input type="checkbox"/> Pediatric Ki Mobility Backpack | \$50 |
| <input type="checkbox"/> Ki Mobility Backpack | \$65 |
| <input type="checkbox"/> Luggage Carrier | \$85 |
| <i>Not available with PRO ELR or Residual Limb Support</i> | |
| <input type="checkbox"/> Black Neoprene Impact Guards | \$90 |
| <input type="checkbox"/> Clear Spoke Guards K0065 | \$115 |
| <i>Not available with Mag wheel</i> | |
| <input type="checkbox"/> Black Spoke Guards K0065 | \$115 |
| <i>Not available with 20" wheel size or Mag wheel</i> | |
| <input type="checkbox"/> Cane and Crutch Holder | \$195 |
| <input type="checkbox"/> IV Pole | \$275 |
| <i>Available only with Standard Straight Back Posts</i> | |
| Step Tube | \$50 |
| <i>Not available with Transit Wheels or 350 lb Weight Capacity</i> | |
| <input type="checkbox"/> Right Side <input type="checkbox"/> Left Side | |

| O ₂ Holder | |
|--|-------|
| O ₂ Holder E2208 | \$295 |
| <input type="checkbox"/> Right Side <input type="checkbox"/> Left Side | |

| Transit Wheels | |
|--|-------|
| <input type="checkbox"/> Transit Wheels K0108 | \$145 |
| <i>Not available with Step Tube, Rear Seat Heights less than 14" or 350 lb Weight Capacity</i> | |

Color

Gloss Finish

- Black N/C
- Blue Skies N/C
- Bubblegum Pink N/C
- Candy Blue N/C
- Candy Purple N/C
- Candy Red N/C
- Copper N/C
- Granny Smith N/C
- Metallic Orange N/C
- Yellow N/C

Matte Finish

- Charcoal Gray N/C
- Emerald Green N/C
- Iridescent Burgundy N/C
- Iridescent Navy N/C
- Military Green N/C
- Retro Red N/C
- Royal Blue N/C
- Shadow Black N/C
- Smooth White N/C

Hydrographics

Hydrographics _____ \$695

- AOR1 (Desert Digital Camo)
- Pink Camo
- Muddy Girl
- NEXT Camo
- Carbon Fiber¹
- Zebra¹

¹Select Hydrographics Base Color. Zebra not available in Black.

Hydrographics Base Color _____ N/C

- Black Blue Granny Smith Pink Purple Red

Colored Anodizing Package

Select desired color. If no selection is made, black will be used.

| | Option | Standard Caster Housing ² | Caster Fork ³ | Axle Plate |
|--------------|--------------------|--------------------------------------|--------------------------|--------------------------|
| | Retail | \$75 | \$75 | \$75 |
| Color | Black ¹ | <input type="checkbox"/> | <input type="checkbox"/> | <input type="checkbox"/> |
| | Blue | <input type="checkbox"/> | <input type="checkbox"/> | <input type="checkbox"/> |
| | Bright Green | <input type="checkbox"/> | <input type="checkbox"/> | <input type="checkbox"/> |
| | Gold | <input type="checkbox"/> | <input type="checkbox"/> | <input type="checkbox"/> |
| | Orange | <input type="checkbox"/> | <input type="checkbox"/> | <input type="checkbox"/> |
| | Pewter | <input type="checkbox"/> | <input type="checkbox"/> | <input type="checkbox"/> |
| | Pink | <input type="checkbox"/> | <input type="checkbox"/> | <input type="checkbox"/> |
| | Purple | <input type="checkbox"/> | <input type="checkbox"/> | <input type="checkbox"/> |
| | Red | <input type="checkbox"/> | <input type="checkbox"/> | <input type="checkbox"/> |

¹Black is no charge

²Available only if corresponding option selection was made

³Not available with 2" wide casters, Frog Legs Suspension Fork or Frog Legs II Suspension System

| Catalyst Front Seat Height Matrix | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | |
|-----------------------------------|-----------------|-------|------|-----------|-----------|-------|------|-----------|-----------|-------|------|-----------|-----------|---|------|-----------|-----------|-------|-------|-------|-----|------|------|-------|-----|----|------|-------|------|-----|----|
| 4" Caster | | | | 5" Caster | | | | 6" Caster | | | | 7" Caster | | | | 8" Caster | | | | | | | | | | | | | | | |
| FSH | Super Low | Stem | Fork | FSH | Super Low | Stem | Fork | FSH | Super Low | Stem | Fork | FSH | Super Low | Stem | Fork | FSH | Super Low | Stem | Fork | | | | | | | | | | | | |
| 15 | 13.5 | Std | 4" | 16 | 14.5 | Std | 4" | 17 | 15.5 | Std | 5" | 18 | 16.5 | Std | 6" | 19 | 17.5 | Std | 6" | | | | | | | | | | | | |
| 15.5 | 14 | | | 16.5 | 15.5 | 3/4 | | 17.5 | 16.5 | 3/4 | | 18.5 | 17 | | | 19.5 | 18.5 | 3/4 | | | | | | | | | | | | | |
| 15.5 ¹ | 14 ¹ | | | 17.5 | 16 | 1 1/2 | | 18.5 | 17 | 1 1/2 | | 18.5 | 17.5 | | | 20.5 | n/a | 1 1/2 | | | | | | | | | | | | | |
| 16 | 14.5 | 3/4 | | 5" | 16 | 14.5 | Std | 5" | 17 | 15.5 | Std | 6" | 19 | 18 | | 3/4 | 7" | 19 | 17.5 | Std | 7" | | | | | | | | | | |
| 16.5 ¹ | 15 ¹ | 1 1/2 | | | 16.5 | 15 | | | 17.5 | 16 | | | 18.5 | 17 | | | | 19.5 | 18 | | | 19.5 | 18 | | | | | | | | |
| 17 | 15.5 | 3/4 | | | 16.5 | 15.5 | | | 17.5 | 16.5 | | | 18 | 16.5 | | | | 18.5 | 17 | | | 19.5 | 18.5 | | | | | | | | |
| 15.5 | 14 | Std | | | 5" | 17 | 16 | | 3/4 | 5" | 17.5 | | 16.5 | 3/4 | | 6" | | 18.5 | 17 | 1 1/2 | | 7" | 20.5 | 19 | Std | 7" | | | | | |
| 16 | 14.5 | | | | | 1 1/2 | 17.5 | | | | 16 | | 18 | | | | | 17 | 18.5 | | | | 17.5 | 19 | | | 17.5 | 19.5 | 18.5 | | |
| 16 | 15 | | | | | 3/4 | 18 | | | | 16.5 | | 18.5 | | | | | 17.5 | 19 | | | | 17.5 | 19.5 | | | 18 | 20 | 19 | | |
| 16.5 | 15.5 | 1 1/2 | | | | 5" | 16.5 | | 15 | | Std | | 5" | 18.5 | | | | 17 | 1 1/2 | 6" | | | 20 | 18.5 | 3/4 | | 7" | 20.5 | n/a | 3/4 | 7" |
| 17 | 15.5 | | 17 | | | | 15.5 | | 19 | | | | | 17.5 | 19 | | | 18 | | | | | 21 | n/a | | | | | | | |
| 17.5 | 16 | | 17.5 | | | | 16 | | 19.5 | | | | | 18 | 19.5 | | | 18.5 | | | | | 20.5 | 19 | | | | 21 | n/a | | |
| 16.5 | 15 | Std | 6" | | | | 17.5 | | 16 | | 3/4 | | | 6" | 18 | | | 16.5 | Std | | | | 7" | 20 | 19 | | | 1 1/2 | 7" | 21 | |
| 17 | 15.5 | | | 18 | | | 16.5 | 18.5 | 17 | | | 19 | | | 17.5 | | 20.5 | n/a | | | 20 | | | 18.5 | | | | | | | |
| 17 | 16 | | | 18.5 | | | 17 | 19.5 | 18 | | | 20 | | | 18.5 | | 20.5 | 19 | | | 21 | | | n/a | | | | | | | |
| 17.5 | 16.5 | 3/4 | | 6" | | | 18 | 16.5 | 1 1/2 | | 6" | 18.5 | | | 17 | | 3/4 | 7" | 21 | | n/a | | | 1 1/2 | 7" | | | | | | |
| 18 | 16.5 | | | | 18.5 | | 17 | 19.5 | | 18 | | 20 | | | 18.5 | 20.5 | | | 19 | | 21 | n/a | | | | | | | | | |
| 18.5 | 17 | | | | 19 | | 17.5 | 19 | | 18 | | 20 | | | 18.5 | 20.5 | | | 19 | | 21 | n/a | | | | | | | | | |
| 18 | 16.5 | Std | | | 7" | | 18 | 16.5 | Std | 7" | | 19.5 | | | 18.5 | 1 1/2 | 7" | | | | | | | 7" | | | | | | | |
| 18.5 | 17 | | | | | 18.5 | 17 | 19.5 | | | | 18 | 20 | | 19 | | | | 20.5 | 19 | 21 | n/a | | | | | | | | | |
| 18.5 | 17.5 | | | | | 19 | 17.5 | 19.5 | | | | 18 | 20 | | 18.5 | | | | 20.5 | 19 | 21 | n/a | | | | | | | | | |
| 19 | 18 | 3/4 | | | | 7" | 18.5 | 17.5 | 3/4 | | | 7" | 20 | | 18.5 | 1 1/2 | | | 7" | | | | | | | 7" | | | | | |
| 19.5 | 18.5 | | 19 | | | | 18 | 20.5 | | | | | 19 | 21 | n/a | | | | | | | | | | | | | | | | |
| 20 | 19 | | 19.5 | | | | 18.5 | 20.5 | | | | | 19 | 21 | n/a | | | | | | | | | | | | | | | | |
| | | | | | | | 19.5 | 18 | 1 1/2 | | | | 7" | 6" x 2" caster not available with 5" fork | | | | | | 7" | | | | | | | | | | | |
| | | | | 20 | | | 18.5 | 20.5 | | | 19 | | | | | | | | | | | | | | | | | | | | |
| | | | | 20.5 | | | 19 | | | | | | | | | | | | | | | | | | | | | | | | |

¹ Not available with 4 x 1.5" Soft Roll Aluminum

Seat heights with "FSH" heading available with Standard Front Frame. Seat heights with "Super Low" heading available with Super Low Front Frame. Actual seat height may vary depending upon configuration.