



# Corner Seat

Phone (800) 981-1540 | Fax (715) 254-0996  
 sales@kimobility.com | www.kimobility.com

Company: \_\_\_\_\_ Ship To: \_\_\_\_\_  
 Account Number: \_\_\_\_\_ Date: \_\_\_\_\_ Address: \_\_\_\_\_  
 Quote  PO Number: \_\_\_\_\_ City: \_\_\_\_\_ State: \_\_\_\_\_  
 RTS: \_\_\_\_\_ Zip: \_\_\_\_\_ Phone: \_\_\_\_\_  
 Marked For: \_\_\_\_\_ Email/Fax: \_\_\_\_\_

Effective 08/27/2018

## Corner Seat

### Standard Model

**Selection Required. Includes Angle Adjustable Backrest, Adjustable Lateral Supports, Abductor and Lap Belt.**

	Description	Item Code	Retail
<input type="checkbox"/>	Size 1 Corner Sitter - Meadow	LCS/1-MW	\$1031
<input type="checkbox"/>	Size 2 Corner Sitter - Mondeo Blue	LCS/2-MB	\$1086
<input type="checkbox"/>	Size 3 Corner Sitter - Mondeo Blue	LCS/3-MB	\$1138

### Accessories

	Description	Item Code	Retail
<input type="checkbox"/>	Size 1 Multi-Adjustable Contoured Head Support - Meadow	LCS/1/104-MW	\$281
<input type="checkbox"/>	Size 2 Multi-Adjustable Contoured Head Support - Mondeo Blue	LCS/2/104-MB	\$296
<input type="checkbox"/>	Size 3 Multi-Adjustable Contoured Head Support - Mondeo Blue	LCS/3/104-MB	\$310
<input type="checkbox"/>	Corner Sitter Tray Table <i>One size fits all</i>	LCT/1/105	\$323

## LECKEY<sup>®</sup>

### Technical Specifications

Specification	Size 1	Size 2	Size 3
Age (approximate)	1 1/2 - 6 years	5 - 9 years	8 - 14 years
Overall User Height	45.9"	52.4"	64.2"
Maximum User Weight	66 lbs	88 lbs	110 lbs
Seat Base Length	20.5"	25.6"	29"
Seat Base Width	18.5"	20.9"	23.6"
Height in Folding Position	7"		
Weight of Seat	33 lbs	35.5 lbs	37 lbs
Backrest Height	11.8"	15.7"	19.7"
Head Support	4.7"	5.5"	6.3"
Head Support Adjustment	4"	5.3"	6.7"
Pommel Adjustment	2.7"		
Tray Height (From Ground)	8.8" - 13.8"		

### Technical Specifications - Metric

Specification	Size 1	Size 2	Size 3
Age (approximate)	1 1/2 - 6 years	5 - 9 years	8 - 14 years
Overall User Height	1165 mm	1330 mm	1630 mm
Maximum User Weight	30 kg	40 kg	50 kg
Seat Base Length	520 mm	650 mm	740 mm
Seat Base Width	470 mm	530 mm	600 mm
Height in Folding Position	180 mm		
Weight of Seat	15 kg	16 kg	17 kg
Backrest Height	300 mm	400 mm	500 mm
Head Support	120 mm	140 mm	160 mm
Head Support Adjustment	100 mm	135 mm	170 mm
Pommel Adjustment	70 mm		
Tray Height (From Ground)	225 mm - 350 mm		