



Easy Seat

Phone (800) 981-1540 | Fax (715) 254-0996
 sales@kimobility.com | www.kimobility.com

Company: _____ Ship To: _____
 Account Number: _____ Date: _____ Address: _____
 Quote PO Number: _____ City: _____ State: _____
 RTS: _____ Zip: _____ Phone: _____
 Marked For: _____ Email/Fax: _____

Effective 08/27/2018

Easy Seat

HCPC: E0638

Easy Seat Complete

Selection required if Easy Seat Basic is not selected. Includes Hip Belt, 1" Side Pads, Tray, Footplate, Sandal Straps, 2 Swivel and 2 Fixed Casters.

	Description	Item Code	Retail
<input type="checkbox"/>	Size 1 Easy Seat Complete - Meadow	LES/1/COMP-MW	\$1315
<input type="checkbox"/>	Size 2 Easy Seat Complete - Mondeo Blue	LES/2/COMP-MB	\$1360
<input type="checkbox"/>	Size 3 Easy Seat Complete - Mondeo Blue	LES/3/COMP-MB	\$1475
<input type="checkbox"/>	Size 4 Easy Seat Complete - Mondeo Blue	LES/4/COMP-MB	\$1575

Easy Seat Basic

Selection required if Easy Seat Complete is not selected. Includes Non-Wheeled Legs, Hip Belt and 1" Side Pads.

	Description	Item Code	Retail
<input type="checkbox"/>	Size 1 Easy Seat Basic - Meadow	LES/1/FL-MW	\$855
<input type="checkbox"/>	Size 2 Easy Seat Basic - Mondeo Blue	LES/2/FL-MB	\$895
<input type="checkbox"/>	Size 3 Easy Seat Basic - Mondeo Blue	LES/3/FL-MB	\$995
<input type="checkbox"/>	Size 4 Easy Seat Basic - Mondeo Blue	LES/4/FL-MB	\$1078

Head Support

Selection of Lateral Supports are required with Flat High Back Headrest. Lateral selection must match Headrest size.

	Description	Item Code	Retail
<input type="checkbox"/>	Size 1 Flat High Back Headrest - Meadow	LES/1/L-MW	\$256
<input type="checkbox"/>	Size 2 Flat High Back Headrest - Mondeo Blue	LES/2/L-MB	\$263
<input type="checkbox"/>	Size 3 Flat High Back Headrest - Mondeo Blue	LES/3/L-MB	\$275
<input type="checkbox"/>	Size 4 Flat High Back Headrest - Mondeo Blue	LES/4/L-MB	\$278
<input type="checkbox"/>	Size 1 Lateral Supports for Flat High Back Headrest - Meadow	LES/1/M-MW	\$102
<input type="checkbox"/>	Size 2 Lateral Supports for Flat High Back Headrest - Mondeo Blue	LES/2/M-MB	\$109
<input type="checkbox"/>	Size 3 Lateral Supports for Flat High Back Headrest - Mondeo Blue	LES/3/M-MB	\$120
<input type="checkbox"/>	Size 4 Lateral Supports for Flat High Back Headrest - Mondeo Blue	LES/4/M-MB	\$125

HCPCS codes provided are based on code verification by the PDAC or our interpretation of the code definitions.

Proper use of HCPCS codes does not ensure coverage or payment. For coverage information, verify the appropriate payer's coverage policy.

LECKEY[®]

Trunk Support

Select only Flexible Lateral Supports or Rigid Lateral Supports

	Description	Item Code	Retail
<input type="checkbox"/>	Size 1 Flexible Lateral Supports - Meadow	LES/1/D-MW	\$169
<input type="checkbox"/>	Size 1 Flexible Lateral Supports - Mondeo Blue	LES/1/D-MB	\$169
<input type="checkbox"/>	Size 2 Flexible Lateral Supports - Mondeo Blue	LES/2/D-MB	\$175
<input type="checkbox"/>	Size 3 Flexible Lateral Supports - Mondeo Blue	LES/3/D-MB	\$180
<input type="checkbox"/>	Size 4 Flexible Lateral Supports - Mondeo Blue	LES/4/D-MB	\$184
<input type="checkbox"/>	Size 1 and 2 Rigid Lateral Supports - Meadow	LES/2/E-MW	\$210
<input type="checkbox"/>	Size 1 and 2 Rigid Lateral Supports - Mondeo Blue	LES/2/E-MB	\$210
<input type="checkbox"/>	Size 3 and 4 Rigid Lateral Supports - Mondeo Blue	LES/3/E-MB	\$233

Additional Side Pads

Must match Easy Seat size

	Description	Item Code	Retail
<input type="checkbox"/>	Size 1 Easy Seat Side Pads 1/2" Thick - Meadow (Pair)	LES/1/B1-MW	\$89
<input type="checkbox"/>	Size 2 Easy Seat Side Pads 1/2" Thick - Mondeo Blue (Pair)	LES/2/B1-MB	\$98
<input type="checkbox"/>	Size 3 Easy Seat Side Pads 1/2" Thick - Mondeo Blue (Pair)	LES/3/B1-MB	\$120
<input type="checkbox"/>	Size 4 Easy Seat Side Pads 1/2" Thick - Mondeo Blue (Pair)	LES/4/B1-MB	\$132
<input type="checkbox"/>	Size 1 Easy Seat Side Pads 1" Thick - Meadow (Pair)	LES/1/B2-MW	\$108
<input type="checkbox"/>	Size 2 Easy Seat Side Pads 1" Thick - Mondeo Blue (Pair)	LES/2/B2-MB	\$120
<input type="checkbox"/>	Size 3 Easy Seat Side Pads 1" Thick - Mondeo Blue (Pair)	LES/3/B2-MB	\$130
<input type="checkbox"/>	Size 4 Easy Seat Side Pads 1" Thick - Mondeo Blue (Pair)	LES/4/B2-MB	\$132
<input type="checkbox"/>	Size 1 Easy Seat Side Pads 2" Thick - Meadow (Pair)	LES/1/B3-MW	\$108
<input type="checkbox"/>	Size 2 Easy Seat Side Pads 2" Thick - Mondeo Blue (Pair)	LES/2/B3-MB	\$117
<input type="checkbox"/>	Size 3 Easy Seat Side Pads 2" Thick - Mondeo Blue (Pair)	LES/3/B3-MB	\$139
<input type="checkbox"/>	Size 4 Easy Seat Side Pads 2" Thick - Mondeo Blue (Pair)	LES/4/B3-MB	\$147

Pommel

	Description	Item Code	Retail
<input type="checkbox"/>	Size 1 Easy Seat Pommel - Meadow	LES/1/C-MW	\$108
<input type="checkbox"/>	Size 1 Easy Seat Pommel - Mondeo Blue	LES/1/C-MB	\$108
<input type="checkbox"/>	Size 2 Easy Seat Pommel - Mondeo Blue	LES/2/C-MB	\$120
<input type="checkbox"/>	Size 3 and 4 Easy Seat Pommel - Mondeo Blue	LES/3/C-MB	\$129

Potty Trainer Seat Insert

Must match Easy Seat size

	Description	Item Code	Retail
<input type="checkbox"/>	Size 1 Easy Potty Trainer Seat Insert	LES/1/S	\$72
<input type="checkbox"/>	Size 2 Easy Potty Trainer Seat Insert	LES/2/S	\$72
<input type="checkbox"/>	Size 3 Easy Potty Trainer Seat Insert	LES/3/S	\$103
<input type="checkbox"/>	Size 4 Easy Potty Trainer Seat Insert	LES/4/S	\$103

LECKEY[®]

Padded Seat Insert

Used with Poty Trainer Seat Insert. Must match Easy Seat size.

	Description	Item Code	Retail
<input type="checkbox"/>	Size 1 Padded Seat Insert - Meadow	LES/1/R-MW	\$79
<input type="checkbox"/>	Size 2 Padded Seat Insert - Mondeo Blue	LES/2/R-MB	\$90
<input type="checkbox"/>	Size 3 Padded Seat Insert - Mondeo Blue	LES/3/R-MB	\$98
<input type="checkbox"/>	Size 4 Padded Seat Insert - Mondeo Blue	LES/4/R-MB	\$109

Foot Positioning

	Description	Item Code	Retail
<input type="checkbox"/>	Size 1 Footplate with Sandal Straps	LES/1/N	\$124
<input type="checkbox"/>	Size 2 Footplate with Sandal Straps	LES/2/N	\$130
<input type="checkbox"/>	Size 3 Footplate with Sandal Straps	LES/3/N	\$141
<input type="checkbox"/>	Size 4 Footplate with Sandal Straps	LES/4/N	\$151
<input type="checkbox"/>	Size 1 Easy Seat Footplate Raiser 1"	Quantity: _____ (1-3) LES/1/K	\$113
<input type="checkbox"/>	Size 2 Easy Seat Footplate Raiser 1"	Quantity: _____ (1-3) LES/2/K	\$124
<input type="checkbox"/>	Size 3 Easy Seat Footplate Raiser 1"	Quantity: _____ (1-3) LES/3/K	\$135
<input type="checkbox"/>	Size 4 Easy Seat Footplate Raiser 1"	Quantity: _____ (1-3) LES/4/K	\$147

Accessories

Must match Easy Seat size

	Description	Item Code	Retail
<input type="checkbox"/>	Size 1 Easy Seat Activity Tray	LES/1/A	\$264
<input type="checkbox"/>	Size 2 Easy Seat Activity Tray	LES/2/A	\$264
<input type="checkbox"/>	Size 3 Easy Seat Activity Tray	LES/3/A	\$264
<input type="checkbox"/>	Size 4 Easy Seat Activity Tray	LES/4/A	\$267
<input type="checkbox"/>	Size 1 Easy Seat Grab/Suspension Rail	LES/1/I	\$244
<input type="checkbox"/>	Size 2 Easy Seat Grab/Suspension Rail	LES/2/I	\$244
<input type="checkbox"/>	Size 3 Easy Seat Grab/Suspension Rail	LES/3/I	\$244
<input type="checkbox"/>	Size 4 Easy Seat Grab/Suspension Rail	LES/4/I	\$244
<input type="checkbox"/>	Size 1 Easy Seat Grab Posts	LES/1/P	\$244
<input type="checkbox"/>	Size 2 Easy Seat Grab Posts	LES/2/P	\$244
<input type="checkbox"/>	Size 3 Easy Seat Grab Posts	LES/3/P	\$244
<input type="checkbox"/>	Size 4 Easy Seat Grab Posts	LES/4/P	\$244
<input type="checkbox"/>	Size 1 Easy Seat Push Handle	LES/1/Q	\$143
<input type="checkbox"/>	Size 2 Easy Seat Push Handle	LES/2/Q	\$143
<input type="checkbox"/>	Size 3 Easy Seat Push Handle	LES/3/Q	\$143
<input type="checkbox"/>	Size 4 Easy Seat Push Handle	LES/4/Q	\$143
<input type="checkbox"/>	Size 1 Caster Legs (2 Fixed and 2 Swivel)	LES/1/H	\$286
<input type="checkbox"/>	Size 2 Caster Legs (2 Fixed and 2 Swivel)	LES/2/H	\$286
<input type="checkbox"/>	Size 3 Caster Legs (2 Fixed and 2 Swivel)	LES/3/H	\$286
<input type="checkbox"/>	Size 4 Caster Legs (2 Fixed and 2 Swivel)	LES/4/H	\$286

LECKEY[®]

Technical Specifications

Specification	Size 1	Size 2	Size 3	Size 4
User Height (Min / Max)	30" / 40"	37" / 46"	40" / 52"	46" / 58"
Maximum User Weight	66 lbs	66 lbs	88 lbs	88 lbs
Seat Width (Min / Max)	5" / 9"	6" / 10"	9" / 13"	11" / 15"
Seat Depth (Min / Max)	7" / 10"	9" / 12"	10" / 14"	12" / 16"
Backrest Height (Min / Max)	9" / 11"	10" / 14"	13" / 17"	16" / 20"
Seat to Footrest (Min / Max)	6" / 9"	9" / 12"	10" / 14"	12" / 17"
Seat to Armrest (Min / Max)	4" / 6"	5" / 7"	6" / 8"	7" / 10"
Tray Dimensions (Length x Width)	18" L x 16.5" W		23" L x 20.5" W	

Technical Specifications - Metric

Specification	Size 1	Size 2	Size 3	Size 4
User Height (Min / Max)	760 mm / 1000 mm	930 mm / 1150 mm	1000 mm / 1300 mm	1150 mm / 1450 mm
Maximum User Weight	30 Kg	30 Kg	40 Kg	40 Kg
Seat Width (Min / Max)	125 mm / 225 mm	150 mm / 250 mm	225 mm / 325 mm	275 mm / 375 mm
Seat Depth (Min / Max)	175 mm / 250 mm	225 mm / 330 mm	250 mm / 350 mm	300 mm / 400 mm
Backrest Height (Min / Max)	225 mm / 285 mm	265 mm / 340 mm	325 mm / 425 mm	400 mm / 500 mm
Seat to Footrest (Min / Max)	150 mm / 230 mm	220 mm / 300 mm	250 mm / 350 mm	300 mm / 420 mm
Seat to Armrest (Min / Max)	110 mm / 150 mm	100 mm / 180 mm	160 mm / 210 mm	180 mm / 240 mm
Tray Dimensions (Length x Width)	450 mm L x 420 mm W		590 mm L x 520 mm W	