



Everyday Activity Seat

Phone (800) 981-1540 | Fax (715) 254-0996
 sales@kimobility.com | www.kimobility.com

Company: _____ Ship To: _____
 Account Number: _____ Date: _____ Address: _____
 Quote PO Number: _____ City: _____ State: _____
 RTS: _____ Zip: _____ Phone: _____
 Marked For: _____ Email/Fax: _____

Effective 03/19/2018

Everyday Activity Seat

HCPC: E1399

Size 1 Weight Capacity: 66 lbs | Size 2 Weight Capacity: 110 lbs | Size 3 weight Capacity: 220 lbs

Base Frame

Selection required

	Description	Item Code	HCPC	Retail
<input type="checkbox"/>	Size 1 Seat Base, Backrest and Hi-Low Chassis with Gas Spring	151-6102	E1399	\$2625
<input type="checkbox"/>	Size 2 Seat Base, Backrest and Hi-Low Chassis with Gas Spring	151-6202	E1399	\$2725
<input type="checkbox"/>	Size 3 Seat Base, Backrest and Hi-Low Hydraulic Assist Chassis	151-6302	E1399	\$5050

Chair Cover

Selection required. Must match base frame size.

	Description	Item Code	HCPC	Retail
<input type="checkbox"/>	Size 1 Microfiber Seat and Back Covers - Blue	151-6104-03		\$315
<input type="checkbox"/>	Size 1 Microfiber Seat and Back Covers - Green	151-6104-01		\$315
<input type="checkbox"/>	Size 1 Microfiber Seat and Back Covers - Raspberry	151-6104-04		\$315
<input type="checkbox"/>	Size 2 Microfiber Seat and Back Covers - Blue	151-6206-03		\$355
<input type="checkbox"/>	Size 2 Microfiber Seat and Back Covers - Green	151-6206-01		\$355
<input type="checkbox"/>	Size 2 Microfiber Seat and Back Covers - Raspberry	151-6206-04		\$355
<input type="checkbox"/>	Size 3 Microfiber Seat and Back Covers - Blue	151-6304-03		\$590
<input type="checkbox"/>	Size 3 Microfiber Seat and Back Covers - Green	151-6304-01		\$590
<input type="checkbox"/>	Size 3 Microfiber Seat and Back Covers - Raspberry	151-6304-04		\$590
<input type="checkbox"/>	Size 1 Vinyl Seat and Back Covers - Blue	151-6104-03VIN		\$315
<input type="checkbox"/>	Size 1 Vinyl Seat and Back Covers - Green	151-6104-01VIN		\$315
<input type="checkbox"/>	Size 1 Vinyl Seat and Back Covers - Raspberry	151-6104-04VIN		\$315
<input type="checkbox"/>	Size 2 Vinyl Seat and Back Covers - Blue	151-6206-03VIN		\$355
<input type="checkbox"/>	Size 2 Vinyl Seat and Back Covers - Green	151-6206-01VIN		\$355
<input type="checkbox"/>	Size 2 Vinyl Seat and Back Covers - Raspberry	151-6206-04VIN		\$355
<input type="checkbox"/>	Size 3 Vinyl Seat and Back Covers - Blue	151-6304-03VIN		\$590
<input type="checkbox"/>	Size 3 Vinyl Seat and Back Covers - Green	151-6304-01VIN		\$590
<input type="checkbox"/>	Size 3 Vinyl Seat and Back Covers - Raspberry	151-6304-04VIN		\$590

HCPCS codes provided are based on code verification by the PDAC or our interpretation of the code definitions.

Proper use of HCPCS codes does not ensure coverage or payment. For coverage information, verify the appropriate payer's coverage policy.

LECKEY[®]

Lateral Support

Laterals must match base frame size.

	Description	Item Code	HCPC	Retail
<input type="checkbox"/>	Size 1 Rigid Chest Laterals - Black Vinyl (Pair)	151-6106-07VIN		\$330
<input type="checkbox"/>	Size 2 Rigid Chest Laterals - Black Vinyl (Pair)	151-6210-07VIN		\$330
<input type="checkbox"/>	Size 3 Rigid Chest Laterals - Black Vinyl (Pair)	151-6306-07VIN		\$450
<input type="checkbox"/>	Size 1 Swing-Away Chest Laterals - Black Vinyl (Pair)	151-6108-07VIN		\$410
<input type="checkbox"/>	Size 2 Swing-Away Chest Laterals - Black Vinyl (Pair)	151-6212-07VIN		\$415
<input type="checkbox"/>	Size 3 Swing-Away Chest Laterals - Black Vinyl (Pair)	151-6308-07VIN		\$545
<input type="checkbox"/>	Size 1 Basic Laterals Fixed - Black Vinyl (Pair)	151-6111-07VIN		\$305
<input type="checkbox"/>	Size 2 Basic Laterals Fixed - Black Vinyl (Pair)	151-6216-07VIN		\$310
<input type="checkbox"/>	Size 3 Basic Laterals Fixed - Black Vinyl (Pair)	151-6312-07VIN		\$370
<input type="checkbox"/>	Size 1 and 2 Seats Protraction Pads - Black Vinyl (Pair)	151-6218-07VIN		\$400
<input type="checkbox"/>	Size 3 Seat Protraction Pads - Black Vinyl (Pair)	151-6314-07VIN		\$430
<input type="checkbox"/>	Chest Support - Black <i>Requires laterals</i>	137-808-05		\$157

Upper Extremity Positioning

Trays require Armrests. Selections must match base frame size.

	Description	Item Code	HCPC	Retail
<input type="checkbox"/>	Size 1 Armrests (Pair)	151-6117		\$343
<input type="checkbox"/>	Size 2 Armrests (Pair)	151-6226		\$343
<input type="checkbox"/>	Size 3 Armrests (Pair)	151-3747		\$482
<input type="checkbox"/>	Size 1 Tray - Plastic Gray	151-1644		\$343
<input type="checkbox"/>	Size 2 Tray - Plastic Gray	151-2644		\$343
<input type="checkbox"/>	Size 3 Tray - Plastic Black	151-3767		\$429

Head Support

	Description	Item Code	HCPC	Retail
<input type="checkbox"/>	Size 1 and 2 Flat Headrest Support and Cushion - Black Vinyl <i>Must match base frame size</i>	151-6208-07VIN		\$320
<input type="checkbox"/>	Size 1 and 2 Laterals for Flat Headrest - Black Vinyl (Pair) <i>Must match base frame size. Includes width adjustable lateral pads and covers. Requires Flat Headrest Support and Cushion.</i>	151-6209-07VIN		\$270
<input type="checkbox"/>	Size 3 Contoured Headrest - Black Neoprene <i>Must match base frame size</i>	151-3751-05		\$417
	Contoured Headrest Adapter Bracket	151-1756		\$71
<input type="checkbox"/>	Contoured Headrest Support Hardware	120-859		\$317
	Contoured Headrest Cushion - Black	137-800-05		\$145

LECKEY^o

Foot Positioning

	Description	Item Code	HCPC	Retail
<input type="checkbox"/>	Size 1 One Piece Footplate - Metal	151-1745		\$343
<input type="checkbox"/>	Size 2 One Piece Footplate - Metal	151-2745		\$343
<input type="checkbox"/>	Size 3 One Piece Footplate - Metal	151-3745		\$583
<input type="checkbox"/>	Small Sandals with Straps (Pair)	152-1600		\$225
<input type="checkbox"/>	Medium Sandals with Straps (Pair)	152-2600		\$245
<input type="checkbox"/>	Large Sandals with Straps (Pair)	152-3600		\$265
<input type="checkbox"/>	Size 1 and 2 Sandal Raiser 1" (Each) <i>Must match Sandal size</i>	Quantity: _____ (1-6)	151-1760	\$50
<input type="checkbox"/>	Size 3 Sandal Raiser 1" (Each) <i>Must match Sandal size</i>	Quantity: _____ (1-6)	151-3760	\$58

Pelvic Fixation

	Description	Item Code	HCPC	Retail
<input type="checkbox"/>	Size 1 Pelvic Belt 4-Point	151-1660		\$126
<input type="checkbox"/>	Size 2 Pelvic Belt 4-Point	151-2660		\$135
<input type="checkbox"/>	Size 3 Pelvic Belt 4-Point	151-3660		\$145
<input type="checkbox"/>	Size 1 Side Pads - Black Vinyl (Pair)	151-6115-07VIN		\$445
<input type="checkbox"/>	Size 2 Side Pads - Black Vinyl (Pair)	151-6222-07VIN		\$530
<input type="checkbox"/>	Size 3 Side Pads - Black Vinyl (Pair)	151-6317-07VIN		\$680
<input type="checkbox"/>	Size 2 Adduction Side Pad - Black Vinyl <i>Not available with size 1 base frame</i>	151-6224-07VIN		\$255
<input type="checkbox"/>	Size 3 Adduction Side Pad - Black Vinyl <i>Not available with size 1 base frame</i>	151-6320-07VIN		\$255
<input type="checkbox"/>	Size 1 Wedge Shape Pommel - Black Vinyl	151-6120-07VIN		\$180
<input type="checkbox"/>	Size 2 Wedge Shape Pommel - Black Vinyl	151-6229-07VIN		\$215
<input type="checkbox"/>	Size 3 Wedge Shape Pommel - Black Vinyl	151-6334-07VIN		\$225

Harnesses

	Description	Item Code	HCPC	Retail
<input type="checkbox"/>	Small Trunk Harness	120-799		\$249
<input type="checkbox"/>	Medium Trunk Harness	120-757		\$249
<input type="checkbox"/>	Large Trunk Harness	117-761		\$249
<input type="checkbox"/>	Harness Extension Strap - 4" Additional Length	117-750		\$62

Assembly

Product Assembly Fee not subject to discounts

<input type="checkbox"/>	Product Assembly Fee			\$120
--------------------------	----------------------	--	--	-------

LECKEY[®]

Technical Specifications

Specification	Size 1	Size 2	Size 3
Max Load	66 lbs	110 lbs	220 lbs
Weight (Seat and Hi-Low Chassis)	65 lbs	75 lbs	108 lbs
Seat Height (Hi-Low Chassis)	9.8" - 23.6"	12.6" - 23.6"	17.7" - 23.6"
Seat Width	6.3" - 11"	10.2" - 15"	13.4" - 19.7"
Seat Depth	5.9" - 11"	9" - 15"	13" - 18.9"
Lateral Distance (Chest)	5.9" - 8.7"	6.7" - 9.4"	9.4" - 14.6"
Back Height	11.8" - 15.7"	15.7" - 20.5"	17.7" - 24.8"
Back Angle	30° Supine / 10° Prone		
Armrest Height	5.9" - 9.8"	7.9" - 11.8"	7" - 12.7"
Lower Leg Length	5.9" - 9.8"	9.8" - 17.7"	9.8" - 17.7"
Lower Leg Angle	-10° to +50°		
Footplate Angle	-15° to +10°		
Tilt-In-Space	30° Supine / 10° Prone		

Technical Specifications - Metric

Specification	Size 1	Size 2	Size 3
Max Load	30 kg	50 kg	100 kg
Weight (Seat and Hi-Low Chassis)	29.5 kg	34 kg	49 kg
Seat Height (Hi-Low Chassis)	250 mm - 600 mm	320 mm - 600 mm	450 mm - 600 mm
Seat Width	160 mm - 280 mm	260 mm - 380 mm	340 mm - 500 mm
Seat Depth	150 mm - 280 mm	230 mm - 380 mm	330 mm - 480 mm
Lateral Distance (Chest)	150 mm - 220 mm	170 mm - 240 mm	240 mm - 370 mm
Back Height	300 mm - 400 mm	400 mm - 520 mm	450 mm - 630 mm
Back Angle	30° Supine / 10° Prone		
Armrest Height	150 mm - 250 mm	200 mm - 300 mm	180 mm - 325 mm
Lower Leg Length	150 mm - 250 mm	250 mm - 450 mm	250 mm - 450 mm
Lower Leg Angle	-10° to +50°		
Footplate Angle	-15° to +10°		
Tilt-In-Space	30° Supine / 10° Prone		