

Company:

Ship To:

Account Number:

Date:

Address:

Quote: PO Number:

City:

State:

ATP:

Zip:

Phone:

Marked For:

Email/Fax:

Associated Quote/Order#:

Effective: 9/6/2023

Little Wave Flip for x:panda (Folding Base)

E1234, E1161 _____ **\$3140**

129 lb weight capacity. Tilt range: 0° - 45°. Hand Tilt only.

Base does not come with seating, armrests, or footrests.

Frame / Backrest

Frame Width

Will accept x:panda sizes 1, 2, and 3

16" Frame for x:panda _____ N/C

Transit Option

Transit System _____ \$240
Omit Transit _____ N/C

Seating System Position

Select Seat Depth of x:panda System if known. If no selection is made, Position 1 will be used.
Compatible with x:panda Size 1, Size 2, & Size 3 and x:panda shape Size 1 w/infant seating Size 1 & Size 2.

Flip

Position 1 (6.25" to 10.25") _____ N/C
Position 2 (10.5" to 13") _____ N/C
Position 3 (13.25" to 19.25") _____ N/C

Backrest

Fixed Height Back with Adjustable Handle (Short 22") _____ N/C
Fixed Height Back with Adjustable Handle (Tall 26") _____ N/C

Frame Height, Casters and Rear Wheels

Frame Height

x:panda seat height will be 3.5" higher above selected Frame Height

Rear Wheel	Frame Height		
	6" Caster	7" Caster	8" Caster
12" Mag	15.5"	15.5"	
	16.5"	16.5"	16.5"
	17.5" ¹	17.5" ¹	17.5" ¹
16" Mag	15.5"	15.5"	
	16.5"	16.5"	16.5"
	17.5"	17.5"	17.5"

¹ Not available with Full Poly Tire

Axles

Quick Release _____ N/C

Casters

Caster Type		Caster Size			
		Price	6"	7"	8"
Caster Type	1" Poly	N/C			
	1" Pneumatic	\$85			
	1.5" Poly	\$85			
	2" Poly	\$90			

Rear Tire Type

Low Poly (E2218) _____ N/C
Available only with 16" wheel size
Full Poly (E2218) _____ N/C
Pneumatic (E2211, E2212) _____ \$90
Pneumatic with Airless Insert (E2211, E2213) _____ \$95

Wheel Locks

Push to Lock _____ N/C
Rear wheels must be removed to fold chair with this wheel lock selection
Attendant Foot Lock _____ N/C

HCPCS codes provided are based on code verification by the PDAC or our interpretation of the code definitions.

Proper use of HCPCS codes does not ensure coverage or payment. For coverage information, verify the appropriate payer's coverage policy.

DCN1068.0 Pricing and specifications subject to change without notice.

Accessories

Anti-tips		Other Accessories	
Rear Anti-Tips (E0971)	\$195	Tool Kit	\$55
Omit	N/C	Pediatric Ki Mobility Backpack	\$70
		Ki Mobility Backpack	\$70
		Tilt Stop Posterior (Pair) (K0108)	\$30
		Tilt Stop Anterior (Pair) (K0108)	\$30
		<i>Use prevents folding</i>	
		IV Pole	\$285

Color

Frame and Accent

Select desired frame and accent colors. Accent color is used on side panel; if no accent color is selected, frame color will be used.

Options	Price	Colors																				
		Gloss													Matte							
		Black	Blue Skies	Bubblegum Pink	Candy Blue	Candy Purple	Candy Red	Copper	Granny Smith	Metallic Orange	Silver	Teal	Yellow	Charcoal Gray	Emerald Green	Iridescent Burgundy	Iridescent Navy	Military Green	Retro Red	Royal Blue	Shadow Black	Smooth White
Frame	N/C																					
Accent	N/C																					