



Freestander

Phone (800) 981-1540 | Fax (715) 254-0996
 sales@kimobility.com | www.kimobility.com

Company: _____ Ship To: _____
 Account Number: _____ Date: _____ Address: _____
 Quote PO Number: _____ City: _____ State: _____
 RTS: _____ Zip: _____ Phone: _____
 Marked For: _____ Email/Fax: _____

Effective 08/27/2018

Freestander

HCPC: E0638

Freestander

Selection required. Includes Chest and Hip Pad, Flexible Lateral, Lap Belt, PU Knee Supports, Footplate, Sandals, Tray Cover and Activity Bowl.

	Description	Item Code	Retail
<input type="checkbox"/>	Size 1 Freestander - Meadow	LFS/1-MW-EXP	\$2326
<input type="checkbox"/>	Size 2 Freestander - Mondeo Blue	LFS/2-MB-EXP	\$2520
<input type="checkbox"/>	Size 3 Freestander - Mondeo Blue	LFS/3-MB-EXP	\$2850

Additional Positioning

	Description	Item Code	Retail
<input type="checkbox"/>	Size 1 Lap Belt - Meadow	LFS/1/B-MW	\$169
<input type="checkbox"/>	Size 1 Lap Belt - Mondeo Blue	LFS/1/B-MB	\$155
<input type="checkbox"/>	Size 2 Lap Belt - Meadow	LFS/2/B-MW	\$214
<input type="checkbox"/>	Size 2 Lap Belt - Mondeo Blue	LFS/2/B-MB	\$175
<input type="checkbox"/>	Size 3 Lap Belt - Mondeo Blue	LFS/3/B-MB	\$190
<input type="checkbox"/>	Size 1 High Back Strap - Meadow	LFS/1/D-MW	\$105
<input type="checkbox"/>	Size 2 High Back Strap - Mondeo Blue	LFS/2/D-MB	\$112
<input type="checkbox"/>	Size 3 High Back Strap - Mondeo Blue	LFS/3/D-MB	\$125
<input type="checkbox"/>	Size 1 Solid Pelvic Plate - Meadow <i>Lap Belt selection required</i>	LFS/1/L-MW	\$273
<input type="checkbox"/>	Size 2 Solid Pelvic Plate - Mondeo Blue <i>Lap Belt selection required</i>	LFS/2/L-MB	\$309
<input type="checkbox"/>	Size 3 Solid Pelvic Plate - Mondeo Blue <i>Lap Belt selection required</i>	LFS/3/L-MB	\$345

HCPCS codes provided are based on code verification by the PDAC or our interpretation of the code definitions.

Proper use of HCPCS codes does not ensure coverage or payment. For coverage information, verify the appropriate payer's coverage policy.

LECKEY[®]

Additional Knee Support

Description			Item Code	Retail (Each)
<input type="checkbox"/>	Size 1 Upholstered Knee Cups - Meadow		LPS/1/K-MW	\$109
<input type="checkbox"/>	Size 2 Upholstered Knee Cups - Mondeo Blue		LPS/2/K-MB	\$109
<input type="checkbox"/>	Size 3 Upholstered Knee Cups - Mondeo Blue		LPS/3/K-MB	\$109
<input type="checkbox"/>	Size 4 Upholstered Knee Cups - Mondeo Blue		LPS/4/K-MB	\$109
<input type="checkbox"/>	Size 1 Knee Spacer (Each)	<input type="checkbox"/> Quantity 1 <input type="checkbox"/> Quantity 2	LPS/1/F	\$55
<input type="checkbox"/>	Size 2 Knee Spacer (Each)	<input type="checkbox"/> Quantity 1 <input type="checkbox"/> Quantity 2	LPS/2/F	\$55
<input type="checkbox"/>	Size 3 Knee Spacer (Each)	<input type="checkbox"/> Quantity 1 <input type="checkbox"/> Quantity 2	LPS/3/F	\$55

Additional Foot Positioning

Description			Item Code	Retail (Each)
<input type="checkbox"/>	Small Sandals with Straps (Pair)		152-1600	\$225
<input type="checkbox"/>	Medium Sandals with Straps (Pair)		152-2600	\$245
<input type="checkbox"/>	Large Sandals with Straps (Pair)		152-3600	\$265
<input type="checkbox"/>	Size 1 Sandal Raiser (Each)	Quantity: _____ (1-6)	LFS/1/J	\$64
<input type="checkbox"/>	Size 2 Sandal Raiser (Each)	Quantity: _____ (1-6)	LFS/2/J	\$64
<input type="checkbox"/>	Size 3 Sandal Raiser (Each)	Quantity: _____ (1-6)	LFS/3/J	\$64

Headrest

Description			Item Code	Retail
<input type="checkbox"/>	Headrest Interface Frame <i>Selection required with Contoured Headrest</i>		LIF/1	\$470
<input type="checkbox"/>	Size 1 Contoured Headrest - Meadow		LIF/2/A-MW	\$319
<input type="checkbox"/>	Size 1 Contoured Headrest - Mondeo Blue		LIF/2/A-MB	\$319
<input type="checkbox"/>	Size 1 Contoured Headrest - Vinyl		LIF/2/A	\$319
<input type="checkbox"/>	Size 2 Contoured Headrest - Mondeo Blue		LIF/2/B-MB	\$319
<input type="checkbox"/>	Size 2 Contoured Headrest - Vinyl		LIF/2/B	\$319
<input type="checkbox"/>	Size 3 Contoured Headrest - Mondeo Blue		LIF/2/C-MB	\$319
<input type="checkbox"/>	Size 3 Contoured Headrest - Vinyl		LIF/2/C	\$319

Accessories

Description			Item Code	Retail
<input type="checkbox"/>	Size 1 Grab/Suspension Rail		LPS/1/C	\$179
<input type="checkbox"/>	Size 2 Grab/Suspension Rail		LPS/2/C	\$179
<input type="checkbox"/>	Size 3 Grab/Suspension Rail		LPS/3/C	\$179

LECKEY[®]

Technical Specifications

Specification	Size 1	Size 2	Size 3
Age (approximate)	1 1/2 - 5 years	4 - 10 years	9 - 18 years
User Body Height (Min / Max)	31.5" / 45.3"	43.3" / 57"	55.1" / 72.8"
Maximum User Weight	132 lbs	176 lbs	198 lbs
Chest Pad Height (Min / Max)	20.7" / 33.5"	26.8" / 43.3"	38.4" / 56.1"
Width Between Chest Supports (Min / Max)	6" / 10"	8" / 12"	10" / 14"
Activity Tray	16" x 21"	18" x 23"	20" x 25"

Technical Specifications - Metric

Specification	Size 1	Size 2	Size 3
Age (approximate)	1 1/2 - 5 years	4 - 10 years	9 - 18 years
User Body Height (Min / Max)	800 mm / 1151 mm	1100 mm / 1448 mm	1400 mm / 1849 mm
Maximum User Weight	60 kg	80 kg	90 kg
Chest Pad Height (Min / Max)	470 mm / 680 mm	680 mm / 975 mm	975 mm / 1350 mm
Width Between Chest Supports (Min / Max)	152 mm / 254 mm	203 mm / 305 mm	254 mm / 356 mm
Activity Tray	400 mm x 505 mm	450 mm x 575 mm	500 mm x 625 mm