



Horizon Stander

Phone (800) 981-1540 | Fax (715) 254-0996
 sales@kimobility.com | www.kimobility.com

Company: _____ Ship To: _____
 Account Number: _____ Date: _____ Address: _____
 Quote PO Number: _____ City: _____ State: _____
 RTS: _____ Zip: _____ Phone: _____
 Marked For: _____ Email/Fax: _____

Effective 08/27/2018

Horizon Stander

HCPC: E0641

Frame

Selection Required. Includes Spine Assembly, Hip Pad, Chest Pad, Flexible Chest Laterals, Standard Knees and Basic Footplate.

	Description	Item Code	HCPC	Retail
<input type="checkbox"/>	Size 1 Horizon Support Frame - Blue	LHZ/1/178-BB-EXP	E1399	\$2236
<input type="checkbox"/>	Size 1 Horizon Support Frame - Pink	LHZ/1/178-PB-EXP	E1399	\$2236
<input type="checkbox"/>	Size 2 Horizon Support Frame - Blue	LHZ/2/178-BB-EXP	E1399	\$2369
<input type="checkbox"/>	Size 2 Horizon Support Frame - Pink	LHZ/2/178-PB-EXP	E1399	\$2369
<input type="checkbox"/>	Size 3 Horizon Support Frame - Blue	LHZ/3/178-BB-EXP	E1399	\$2929
<input type="checkbox"/>	Size 3 Horizon Support Frame - Pink	LHZ/3/178-PB-EXP	E1399	\$2929

Base

Selection Required. Must match Frame size. Includes Worldwide Adapters.

	Description	Item Code	HCPC	Retail
<input type="checkbox"/>	Size 1 Horizon Base/Chassis	LHZ/1/177-EXP	E0641	\$2744
<input type="checkbox"/>	Size 2 Horizon Base/Chassis	LHZ/2/177-EXP	E0641	\$2911
<input type="checkbox"/>	Size 3 Horizon Base/Chassis	LHZ/3/177-EXP	E0641	\$3634

HCPCS codes provided are based on code verification by the PDAC or our interpretation of the code definitions.

Proper use of HCPCS codes does not ensure coverage or payment. For coverage information, verify the appropriate payer's coverage policy.

LECKEY^o

Headrest

Head Support Assembly or Extended Head Support Assembly selection required if Head Lateral Supports is selected. Head Support Assembly size must match Head Lateral Supports size.

	Description	Item Code	HCPC	Retail
<input type="checkbox"/>	Size 1 Head Support Assembly - Blue	LHZ/1/169-BB	E1399	\$215
<input type="checkbox"/>	Size 1 Head Support Assembly - Pink	LHZ/1/169-PB	E1399	\$215
<input type="checkbox"/>	Size 2 and 3 Head Support Assembly - Blue	LHZ/2/169-BB	E1399	\$237
<input type="checkbox"/>	Size 2 and 3 Head Support Assembly - Pink	LHZ/2/169-PB	E1399	\$237
<input type="checkbox"/>	Size 1 Head Lateral Supports - Blue	LHZ/1/117-BB	E1399	\$202
	Size 1 Head Lateral Supports - Pink	LHZ/1/117-PB	E1399	\$202
	Size 2 and 3 Head Lateral Supports - Blue	LHZ/2/117-BB	E1399	\$204
	Size 2 and 3 Head Lateral Supports - Pink	LHZ/2/117-PB	E1399	\$204
<input type="checkbox"/>	Size 1 Extended 4" Head Support Assembly - Blue	LHZ/1/187-BB	E1399	\$317
<input type="checkbox"/>	Size 1 Extended 4" Head Support Assembly - Pink	LHZ/1/187-PB	E1399	\$317
<input type="checkbox"/>	Size 2 and 3 Extended 4" Head Support Assembly - Blue	LHZ/2/187-BB	E1399	\$317
<input type="checkbox"/>	Size 2 and 3 Extended 4" Head Support Assembly - Pink	LHZ/2/187-PB	E1399	\$317
<input type="checkbox"/>	Size 1 Head Lateral Supports - Blue	LHZ/1/117-BB	E1399	\$202
	Size 1 Head Lateral Supports - Pink	LHZ/1/117-PB	E1399	\$202
	Size 2 and 3 Head Lateral Supports - Blue	LHZ/2/117-BB	E1399	\$204
	Size 2 and 3 Head Lateral Supports - Pink	LHZ/2/117-PB	E1399	\$204
<input type="checkbox"/>	Size 1-3 Contoured Head Support	LHZ/1/180	E1399	\$519

Additional Chest Support

	Description	Item Code	HCPC	Retail
<input type="checkbox"/>	Size 1 Flexible Chest Laterals - Blue	LHZ/1/168-BB		\$363
<input type="checkbox"/>	Size 1 Flexible Chest Laterals - Pink	LHZ/1/168-PB		\$363
<input type="checkbox"/>	Size 2 Flexible Chest Laterals - Blue	LHZ/2/168-BB		\$426
<input type="checkbox"/>	Size 2 Flexible Chest Laterals - Pink	LHZ/2/168-PB		\$426
<input type="checkbox"/>	Size 3 Flexible Chest Laterals - Blue	LHZ/3/168-BB		\$475
<input type="checkbox"/>	Size 3 Flexible Chest Laterals - Pink	LHZ/3/168-PB		\$475
<input type="checkbox"/>	Size 3 Flip Away Chest Laterals - Blue	LHZ/3/154-BB		\$523
<input type="checkbox"/>	Size 3 Flip Away Chest Laterals - Pink	LHZ/3/154-PB		\$523

LECKEY^o

Hip Laterals

Selection required if Hip Belt is not selected

	Description	Item Code	HCPC	Retail
<input type="checkbox"/>	Size 1 Horizon Hip Laterals - Blue	LHZ/1/164-BB		\$271
<input type="checkbox"/>	Size 1 Horizon Hip Laterals - Pink	LHZ/1/164-PB		\$271
<input type="checkbox"/>	Size 2 Horizon Hip Laterals - Blue	LHZ/2/164-BB		\$304
<input type="checkbox"/>	Size 2 Horizon Hip Laterals - Pink	LHZ/2/164-PB		\$304
<input type="checkbox"/>	Size 3 Horizon Hip Laterals - Blue	LHZ/3/164-BB		\$315
<input type="checkbox"/>	Size 3 Horizon Hip Laterals - Pink	LHZ/3/164-PB		\$315
<input type="checkbox"/>	Size 3 Horizon Hip Laterals for Size 2 Frame - Blue	LHZ/2/162-BB		\$315
<input type="checkbox"/>	Size 3 Horizon Hip Laterals for Size 2 Frame - Pink	LHZ/2/162-PB		\$315
<input type="checkbox"/>	Size 3 Flip Away Hip Laterals - Blue	LHZ/3/153-BB		\$423
<input type="checkbox"/>	Size 3 Flip Away Hip Laterals - Pink	LHZ/3/153-PB		\$423

Hip Belts

Selection required if Hip Laterals are not selected . Multiple selections can be made.

	Description	Item Code	HCPC	Retail
<input type="checkbox"/>	Size 1 Supine Belt - Blue	LHZ/1/181-BB		\$167
<input type="checkbox"/>	Size 1 Supine Belt - Pink	LHZ/1/181-PB		\$167
<input type="checkbox"/>	Size 2 Supine Belt - Blue	LHZ/2/181-BB		\$182
<input type="checkbox"/>	Size 2 Supine Belt - Pink	LHZ/2/181-PB		\$182
<input type="checkbox"/>	Size 3 Supine Belt - Blue	LHZ/3/181-BB		\$200
<input type="checkbox"/>	Size 3 Supine Belt - Pink	LHZ/3/181-PB		\$200
<input type="checkbox"/>	Size 1 Prone Belt - Blue	LHZ/1/182-BB		\$167
<input type="checkbox"/>	Size 1 Prone Belt - Pink	LHZ/1/182-PB		\$167
<input type="checkbox"/>	Size 2 Prone Belt - Blue	LHZ/2/182-BB		\$182
<input type="checkbox"/>	Size 2 Prone Belt - Pink	LHZ/2/182-PB		\$182
<input type="checkbox"/>	Size 3 Prone Belt - Blue	LHZ/3/182-BB		\$200
<input type="checkbox"/>	Size 3 Prone Belt - Pink	LHZ/3/182-PB		\$200

LECKEY[®]

Additional Knee Support

Standard Knee Assembly has 1"-3" adjustment or 1"-5" adjustment if flipped. Medium Extended Knee Assembly has 5" - 7" adjustment and Long Extended Knee Assembly has 7" - 9" adjustment.

	Description	Item Code	HCPC	Retail
<input type="checkbox"/>	Size 1 Standard Knee Assembly - Blue	LHZ/1/167-BB		\$426
<input type="checkbox"/>	Size 1 Standard Knee Assembly - Pink	LHZ/1/167-PB		\$426
<input type="checkbox"/>	Size 2 and 3 Standard Knee Assembly - Blue	LHZ/2/167-BB		\$426
<input type="checkbox"/>	Size 2 and 3 Standard Knee Assembly - Pink	LHZ/2/167-PB		\$426
<input type="checkbox"/>	Size 1 Medium Extended Knee Assembly - Blue	LHZ/1/183-BB		\$560
<input type="checkbox"/>	Size 1 Medium Extended Knee Assembly - Pink	LHZ/1/183-PB		\$560
<input type="checkbox"/>	Size 2 and 3 Medium Extended Knee Assembly - Blue	LHZ/2/183-BB		\$560
<input type="checkbox"/>	Size 2 and 3 Medium Extended Knee Assembly - Pink	LHZ/2/183-PB		\$560
<input type="checkbox"/>	Size 1 Long Extended Knee Assembly - Blue	LHZ/1/184-BB		\$593
<input type="checkbox"/>	Size 1 Long Extended Knee Assembly - Pink	LHZ/1/184-PB		\$593
<input type="checkbox"/>	Size 2 and 3 Long Extended Knee Assembly - Blue	LHZ/2/184-BB		\$593
<input type="checkbox"/>	Size 2 and 3 Long Extended Knee Assembly - Pink	LHZ/2/184-PB		\$593
<input type="checkbox"/>	Size 1-3 Medium Extended Knee Brackets Only (Pair)	LHZ/2/185		\$471
<input type="checkbox"/>	Size 1 Long Extended Knee Brackets Only (Pair)	LHZ/2/186		\$475

Lower Leg Positioning

	Description	Item Code	HCPC	Retail
<input type="checkbox"/>	Small Sandals with Straps (Pair)	152-1600		\$225
<input type="checkbox"/>	Medium Sandals with Straps (Pair)	152-2600		\$245
<input type="checkbox"/>	Large Sandals with Straps (Pair)	152-3600		\$265
<input type="checkbox"/>	Size 1 Flip Up Footplate Assembly	LHZ/1/161		\$672
<input type="checkbox"/>	Size 2 and 3 Flip Up Footplate Assembly	LHZ/2/161		\$671
<input type="checkbox"/>	Size 1-3 3D Footplate Assembly <i>One size fits all</i>	LHZ/3/160		\$1098

Accessories

Must match Horizon Stander size

	Description	Item Code	HCPC	Retail
<input type="checkbox"/>	Size 1 Supine Tray	LHZ/1/143		\$330
<input type="checkbox"/>	Size 2 and 3 Supine Tray	LHZ/2/143		\$356
<input type="checkbox"/>	Size 1 Prone Tray	LHZ/1/158		\$404
<input type="checkbox"/>	Size 2 and 3 Prone Tray	LHZ/2/158		\$471
<input type="checkbox"/>	Size 1 Horizon Grab Rail	LHZ/1/192		\$234
<input type="checkbox"/>	Size 2 and 3 Horizon Grab Rail	LHZ/2/192		\$234
<input type="checkbox"/>	Size 1 Horizon Grab Posts	LHZ/1/193		\$234
<input type="checkbox"/>	Size 2 and 3 Horizon Grab Posts	LHZ/2/193		\$234

LECKEY[®]

Technical Specifications

Specification	Size 1	Size 2	Size 3
Age Range	4 - 10	9 - 18	14 - Adult
Horizontal Position L x W x H	44.7" x 22.8" x 30.7"	50.4" x 25.2" x 31.7"	
Weight of Product	97 lbs	110.6 lbs	130 lbs
Maximum User Weight	110 lbs	176 lbs	220 lbs
User Standing Height (Min / Max)	39.4" / 55.1"	55.1" x 70.9"	61" / 70.9"
Shoulder Height - Footplate to Top of Pad (Min / Max)	32.3" / 45.3"	39.4" / 57"	
Hip Pad Height - Footplate to Top of Pad (Min / Max)	17.7" / 29.5"	29.5" / 37.4"	29.5" / 45.3"
Hip Pad Width (Min / Max)	6.9" / 10.6"	8.3" / 15.3"	12.2" / 18.5"
Knee Height - Footplate to Middle of Pad (Min / Max)	9.8" / 17.7"	15.7" / 21.6"	15.7" / 25.6"
Width Between Chest Laterals (Min / Max)	7.9" / 11.8"	9.4" / 14.9"	13.4" / 18.1"
Footplate Angle Adjustable + / -	10°		
Tray Sizes	16.5" x 17.3"	20.5" x 23.2"	
Tray Angle Adjustment - Prone Version	0° to 40°		
Tray Angle Adjustment - Supine Version	-20° to +30°		
Frame Clearance - Bottom of Ski to Ground	6.1"	7"	

Technical Specifications - Metric

Specification	Size 1	Size 2	Size 3
Age Range	4 - 10	9 - 18	14 - Adult
Horizontal Position L x W x H	1135 mmx580 mm x 780 mm	1280 mm x 640 mm x 805 mm	
Weight of Product	44 kgs	50 kgs	95 kg
Maximum User Weight	50 kgs	80 kgs	100 kgs
User Standing Height (Min / Max)	1000 mm / 1400 mm	1400 mm / 1800 mm	1550 mm / 1800 mm
Shoulder Height - Footplate to Top of Pad (Min / Max)	820 mm / 1150 mm	1000 mm / 1450 mm	
Hip Pad Height - Footplate to Top of Pad (Min / Max)	450 mm / 750 mm	750 mm / 950 mm	750 mm / 1150 mm
Hip Pad Width (Min / Max)	175 mm / 270 mm	210 mm / 390 mm	310 mm / 470 mm
Knee Height - Footplate to Middle of Pad (Min / Max)	250 mm / 450 mm	400 mm / 550 mm	400 mm / 650 mm
Width Between Chest Laterals (Min / Max)	200 mm / 300 mm	240 mm / 380 mm	340 mm / 460 mm
Footplate Angle Adjustable + / -	10°		
Tray Sizes	420 mm x 440 mm	520 mm x 590 mm	
Tray Angle Adjustment - Prone Version	0° to 40°		
Tray Angle Adjustment - Supine Version	-20° to +30°		
Frame Clearance - Bottom of Ski to Ground	155 mm	180 mm	