



# LIBERTY FT Basic

Phone (800) 981-1540 | Fax (715) 254-0996  
 sales@kimobility.com | www.kimobility.com

Company: \_\_\_\_\_

Ship To: \_\_\_\_\_

Account Number: \_\_\_\_\_ Date: \_\_\_\_\_

Address: \_\_\_\_\_

Quote  PO Number: \_\_\_\_\_

City: \_\_\_\_\_ State: \_\_\_\_\_

RTS: \_\_\_\_\_

Zip: \_\_\_\_\_ Phone: \_\_\_\_\_

Marked For: \_\_\_\_\_

Email/Fax: \_\_\_\_\_

Select if chair quote has:  Additional Ki Products with Different Quote/Order Number **Quote/Order Number:** \_\_\_\_\_

Effective 02/14/2022

**Liberty FT Basic** E1161 **\$2795**  
 250 lb weight capacity. Tilt range: 0° - 20°. Standard seat upholstery.

**Liberty FT Basic - Rigid** E1161 **\$2795**  
 250 lb weight capacity. Tilt range: 0° - 20°. Standard seat upholstery.

## Frame

### Seat Width

**14"-16" Seat Widths ship with 300 N Assistive Gas Strut, 18"-22" Seat Widths ship with 450 N Assistive Gas Strut**

- 14" \_\_\_\_\_ N/C
- 16" \_\_\_\_\_ N/C
- 18" \_\_\_\_\_ N/C
- 20" E2201 \$470
- 22" E2201 \$470

### Seat Depth

- 14" \_\_\_\_\_ N/C
- 16" \_\_\_\_\_ N/C
- 18" \_\_\_\_\_ N/C
- 20" E2203 \$605

### Seat Option

- Standard Seat Sling \_\_\_\_\_ Std
- Growing Adjustable Seat Pan E2231 **New!** \$200
- Omit Seat \_\_\_\_\_ N/C

### Transit

- Transit System \_\_\_\_\_ \$320
- Omit Transit \_\_\_\_\_ N/C

## Backrest

### Back Post / Height Range / Back Height

Back Post	Price	Height Range	Back Height				
<input type="checkbox"/> Height Adjustable 8° Bend with Push Handle	\$25	<input type="checkbox"/> Short	<input type="checkbox"/> 14	<input type="checkbox"/> 15	<input type="checkbox"/> 16	<input type="checkbox"/> 17	<input type="checkbox"/> 18
		<input type="checkbox"/> Medium	<input type="checkbox"/> 16	<input type="checkbox"/> 17	<input type="checkbox"/> 18	<input type="checkbox"/> 19	<input type="checkbox"/> 20
		<input type="checkbox"/> Tall	<input type="checkbox"/> 18	<input type="checkbox"/> 19	<input type="checkbox"/> 20	<input type="checkbox"/> 21	<input type="checkbox"/> 22

### Backrest Options

- Fixed (Non-Folding) \_\_\_\_\_ Std
- Angle Adjustable (Folding) \_\_\_\_\_ \$95

### Back Upholstery

**If selecting both upholstery and Axiom Back, the upholstery will be installed**

- Standard \_\_\_\_\_ Std
- Tension Adjustable E2611 \$245  
Not available on Seat Widths less than 16" combined with 21" and 22" Backrest Heights
- Omit Back Upholstery \_\_\_\_\_ - \$15

## Footrest

### Hanger / Front Frame Bend

Hanger	60°	70°		Retail
Extension Tube Mount	<input type="checkbox"/>	<input type="checkbox"/>		Std
Pro Elevating Leg Rests (PRO ELR) E0990, Rental K0195			<input type="checkbox"/>	\$295
Omit Hangers			<input type="checkbox"/>	N/C

### Footplate

- Composite \_\_\_\_\_ Std
- Aluminum Angle Adjustable K0040 \$195

HCPCS codes provided are based on code verification by the PDAC or our interpretation of the code definitions. Proper use of HCPCS codes does not ensure coverage or payment. For coverage information, verify the appropriate payer's coverage policy.

**Footrest continued**

**Footrest / Extension / Seat to Footrest Length**

Select Extension and Seat to Footrest Length in row corresponding to desired Footrest. Footrest Length of at least 2.5" shorter than Seat Height recommended. Extension Tube Hanger measurements based on 70° Hanger Bend.

Footrest	Extension	Seat to Footrest Length									
Extension Tube Hangers (60°, 70°)	<input type="checkbox"/> Short	<input type="checkbox"/> 10	<input type="checkbox"/> 10.5	<input type="checkbox"/> 11	<input type="checkbox"/> 11.5	<input type="checkbox"/> 12	<input type="checkbox"/> 12.5	<input type="checkbox"/> 13	<input type="checkbox"/> 13.5		
	<input type="checkbox"/> Medium	<input type="checkbox"/> 11.5	<input type="checkbox"/> 12	<input type="checkbox"/> 12.5	<input type="checkbox"/> 13	<input type="checkbox"/> 13.5	<input type="checkbox"/> 14	<input type="checkbox"/> 14.5	<input type="checkbox"/> 15	<input type="checkbox"/> 15.5	
	<input type="checkbox"/> Long	<input type="checkbox"/> 14	<input type="checkbox"/> 14.5	<input type="checkbox"/> 15	<input type="checkbox"/> 15.5	<input type="checkbox"/> 16	<input type="checkbox"/> 16.5	<input type="checkbox"/> 17	<input type="checkbox"/> 17.5	<input type="checkbox"/> 18	
PRO ELR Low Pivot Setting	<input type="checkbox"/> Short	10.5 - 15.5									
	<input type="checkbox"/> Long	13 - 19									
PRO ELR High Pivot Setting	<input type="checkbox"/> Short	8.5 - 17									
	<input type="checkbox"/> Long	11.5 - 17.5									

**PRO ELR Cover**

Available only on PRO ELR

Gel Pad Cover K0108 \_\_\_\_\_ \$390

**Heel Loop**

Heel Loops E0951 \_\_\_\_\_ \$40  
 Omit Heel Loops \_\_\_\_\_ N/C

**Calf Strap**

Velcro Adjustable K0038 \_\_\_\_\_ \$35

**Seat Height, Casters and Rear Wheels**

**Seat Height / Caster Size / Rear Wheel Size**


Select Seat Height corresponding to desired Rear Wheel and Caster Size

Rear Wheel Size	Caster Size		
	4"	5"	6"
16"	Seat Height		
	<input type="checkbox"/> 13.5 <sup>1</sup>	<input type="checkbox"/> 14	<input type="checkbox"/> 15
	<input type="checkbox"/> 14 <sup>1</sup>	<input type="checkbox"/> 14.5	<input type="checkbox"/> 15.5
	<input type="checkbox"/> 14.5	<input type="checkbox"/> 15	<input type="checkbox"/> 16
	<input type="checkbox"/> 15	<input type="checkbox"/> 15.5	<input type="checkbox"/> 16.5
	<input type="checkbox"/> 15.5	<input type="checkbox"/> 16	<input type="checkbox"/> 17
	<input type="checkbox"/> 16	<input type="checkbox"/> 16.5	<input type="checkbox"/> 17.5
		<input type="checkbox"/> 17	<input type="checkbox"/> 18
		<input type="checkbox"/> 17.5	<input type="checkbox"/> 18.5
		<input type="checkbox"/> 18	<input type="checkbox"/> 19
		<input type="checkbox"/> 18.5	
		<input type="checkbox"/> 19	
20"	<input type="checkbox"/> 13	<input type="checkbox"/> 14	<input type="checkbox"/> 15
	<input type="checkbox"/> 13.5	<input type="checkbox"/> 14.5	<input type="checkbox"/> 15.5
	<input type="checkbox"/> 14	<input type="checkbox"/> 15	<input type="checkbox"/> 16
	<input type="checkbox"/> 14.5	<input type="checkbox"/> 15.5	<input type="checkbox"/> 16.5
	<input type="checkbox"/> 15	<input type="checkbox"/> 16	<input type="checkbox"/> 17
	<input type="checkbox"/> 15.5	<input type="checkbox"/> 16.5	<input type="checkbox"/> 17.5
	<input type="checkbox"/> 16	<input type="checkbox"/> 17	<input type="checkbox"/> 18
		<input type="checkbox"/> 17.5	<input type="checkbox"/> 18.5
		<input type="checkbox"/> 18	<input type="checkbox"/> 19
		<input type="checkbox"/> 18.5	
		<input type="checkbox"/> 19	

<sup>1</sup> Available only with Low Profile tires

**Hanger Release**

- 4-Way - In line \_\_\_\_\_ Std  
See reference A below. Not available with PRO ELR.
- 4-Way - Out \_\_\_\_\_ N/C  
See reference B below.



Rear Wheel Size	Caster Size		
	4"	5"	6"
22"	Seat Height		
	<input type="checkbox"/> 14	<input type="checkbox"/> 14	<input type="checkbox"/> 15
	<input type="checkbox"/> 14.5	<input type="checkbox"/> 14.5	<input type="checkbox"/> 15.5
	<input type="checkbox"/> 15	<input type="checkbox"/> 15	<input type="checkbox"/> 16
	<input type="checkbox"/> 15.5	<input type="checkbox"/> 15.5	<input type="checkbox"/> 16.5
	<input type="checkbox"/> 16	<input type="checkbox"/> 16	<input type="checkbox"/> 17
		<input type="checkbox"/> 16.5	<input type="checkbox"/> 17.5
		<input type="checkbox"/> 17	<input type="checkbox"/> 18
		<input type="checkbox"/> 17.5	<input type="checkbox"/> 18.5
		<input type="checkbox"/> 18	<input type="checkbox"/> 19
		<input type="checkbox"/> 18.5	
		<input type="checkbox"/> 19	
24"	<input type="checkbox"/> 15	<input type="checkbox"/> 14.5	<input type="checkbox"/> 15
	<input type="checkbox"/> 15.5	<input type="checkbox"/> 15	<input type="checkbox"/> 15.5
	<input type="checkbox"/> 16	<input type="checkbox"/> 15.5	<input type="checkbox"/> 16
		<input type="checkbox"/> 16	<input type="checkbox"/> 16.5
		<input type="checkbox"/> 16.5	<input type="checkbox"/> 17
		<input type="checkbox"/> 17	<input type="checkbox"/> 17.5
		<input type="checkbox"/> 17.5	<input type="checkbox"/> 18
		<input type="checkbox"/> 18	<input type="checkbox"/> 18.5
		<input type="checkbox"/> 18.5	<input type="checkbox"/> 19
		<input type="checkbox"/> 19	

## Seat Height, Casters and Rear Wheels continued

### Caster

		Caster Size			
		Price	4"	5"	6"
Caster Type	1" Poly	Std	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>
	1.5" Poly	\$55	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>

### Rear Wheel

		Rear Wheel Size				
		Price	16"	20"	22"	24"
Rear Wheel Type	Mag	Std	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>
	MAXX Performance Spoke	N/C				<input type="checkbox"/>
	MAXX Performance Spoke	\$180		<input type="checkbox"/>	<input type="checkbox"/>	

## Wheel Locks

### Wheel Lock

- Push to Lock \_\_\_\_\_ N/C  
 Hemi Wheel Lock K0108 \_\_\_\_\_ \$300  
*Extension Handle included*  
 Left  Right

- Omit Wheel Locks \_\_\_\_\_ N/C

### Wheel Lock Options

- Extension Handles E0961 \_\_\_\_\_ \$70

## Accessories

### Anti-Tip

- Rear Anti-Tips E0971 \_\_\_\_\_ \$150  
 Omit Rear Anti-Tips \_\_\_\_\_ N/C

### Transit Wheels

- Transit Wheels K0108 \_\_\_\_\_ \$300  
*Not available with Seat Heights less than 14" or greater than 18", Seat Depths greater than 18", or 12" and 16" rear wheels*

### Other Accessories

- 1.5" Auto Buckle E0978 \_\_\_\_\_ \$80  
 Seat Belt Mounting Hardware \_\_\_\_\_ \$16  
 O<sub>2</sub> Holder E2208 \_\_\_\_\_ \$295  
*Mounted to the left side. Not available with Fold Down Back.*  
 Cane and Crutch Holder \_\_\_\_\_ \$195

### Rear Tire Type

#### Full Profile

- Street Pneumatic E2211, E2212 \_\_\_\_\_ Std  
 Full Poly E2218 \_\_\_\_\_ N/C  
 Pneumatic with Airless Insert E2211, E2213 \_\_\_\_\_ \$90

#### Low Profile

- Low Poly E2218 \_\_\_\_\_ N/C  
*Available only with 16" Mag Wheels*

### Handrim

**Not Available with 16" Rear Wheels**

#### Standard Handrims

- Aluminum Anodized \_\_\_\_\_ Std  
 Plastic Coated \_\_\_\_\_ \$115

#### Standard Handrim Mount

**If no selection is made, Standard will be used**

- Handrim Mount \_\_\_\_\_ N/C  
 Standard Mount  Close Mount  
*Reduces overall width by .75"*

### Axle

- Quick Release \_\_\_\_\_ Std

## Armrests

### Armrest

**Selection Required**

Height Adjustable T-Arm (8.5" - 13.5") E0973

- Desk \_\_\_\_\_ \$270  
 Full \_\_\_\_\_ \$270

### Arm Pad

- Standard Arm Pad \_\_\_\_\_ N/C

## Frame Color

### Gloss Colors

- Black \_\_\_\_\_ N/C  
 Candy Blue \_\_\_\_\_ N/C  
 Candy Red \_\_\_\_\_ N/C

### Matte Colors

- Charcoal Gray \_\_\_\_\_ N/C  
 Shadow Black \_\_\_\_\_ N/C