



Phone (800) 981-1540 | Fax (715) 254-0996
 sales@kimobility.com | www.kimobility.com

Company: _____ Ship To: _____
 Account Number: _____ Date: _____ Address: _____
 Quote PO Number: _____ City: _____ State: _____
 RTS: _____ Zip: _____ Phone: _____
 Marked For: _____ Email/Fax: _____

Select if chair quote has: Additional Ki Products with Different Quote/Order Number **Quote/Order Number:** _____
 Effective 01/17/2022

- Little Wave Clik XPe - Transit Included**
 K0005, E1235, E1237 **\$2695**
Width Growable Frame. Price does not include one time free growth kit. 165 lb weight capacity. Maximum Seat Width of 18".
- Little Wave Clik XPe - Transit Not Included**
 K0005, E1235, E1237 **\$2475**
Width Growable Frame. Price does not include one time free growth kit. 165 lb weight capacity. Maximum Seat Width of 18".

HCPCS codes provided are based on code verification by the PDAC or our interpretation of the code definitions.
 Proper use of HCPCS codes does not ensure coverage or payment. For coverage information, verify the appropriate payer's coverage policy.

Frame

Seat Width

- See reference A on page 7*
- 8" _____ N/C
Not available with 14"-16" Seat Depth
 - 9" _____ N/C
Not available with 15"-16" Seat Depth
 - 10" _____ N/C
Not available with 16" Seat Depth
 - 11" _____ N/C
 - 12" _____ N/C
 - 13" _____ N/C
 - 14" _____ N/C
Not available with 8" Seat Depth
 - 15" _____ N/C
Not available with 8"-9" Seat Depth
 - 16" _____ N/C
Not available with 8"-10" Seat Depth

Seat Depth

- See reference B on page 7*
- 8" _____ N/C
 - 9" _____ N/C
 - 10" _____ N/C
 - 11" _____ N/C
 - 12" _____ N/C
 - 13" _____ N/C
 - 14" _____ N/C
 - 15" _____ N/C
 - 16" _____ N/C

Adjustable Depth

- Adds the capability to grow seat in depth. For restrictions see Chart 2 on page 7.*
- No Rear Adjustment _____ N/C
 - 1" _____ N/C
 - 2" _____ N/C
 - 3" _____ N/C

Transit

- Transit System _____ \$220
No charge when selected with Little Wave Clik XPe - Transit Included
- Omit Transit _____ N/C

Offset Option

- See Chart 1 on page 7*
- Straight Frame _____ Std
 - Offset Frame _____ N/C
Not available with 0 and -1 Sling Positions

Seat Option

- Growing Seat Sling _____ N/C
- Growing Solid Seat Pan E2231 _____ \$265
- Omit Seat _____ N/C

Sling Position

- With 0 Sling Position, sling is approximately 4" from front down tube. The 1 and 2 Sling Positions extend the frame in front of the sling by that amount. -1 pushes the seat 1" farther forward. See reference D on page 7 for further detail and restrictions.*
- 2 _____ N/C
 - 1 _____ N/C
 - 0 _____ N/C
Not available with Offset Frame
 - 1 _____ N/C
Not available with Offset Frame

Front Caster Wing

- Certain combinations of Seat Depth, Center of Gravity Position and Wheel Size may require a forward facing caster wing*
- Forward Facing Wing _____ N/C

Front Frame Bend

- Angle between front seat tube and ground. See reference C on page 7.*
- 90° _____ N/C
 - 85° _____ N/C
 - 80° _____ N/C

Backrest

Back Post / Height Range / Back Height						
Back Post	Height Range	Back Height				
Height Adjustable (Standard or 1" Inward Offset)	<input type="checkbox"/> Ex. Short ¹	<input type="checkbox"/> 9	<input type="checkbox"/> 10	<input type="checkbox"/> 11	<input type="checkbox"/> 12	
	<input type="checkbox"/> Short	<input type="checkbox"/> 11	<input type="checkbox"/> 12	<input type="checkbox"/> 13	<input type="checkbox"/> 14 ¹	
	<input type="checkbox"/> Medium	<input type="checkbox"/> 13	<input type="checkbox"/> 14	<input type="checkbox"/> 15	<input type="checkbox"/> 16	
	<input type="checkbox"/> Tall	<input type="checkbox"/> 15	<input type="checkbox"/> 16	<input type="checkbox"/> 17	<input type="checkbox"/> 18	
Fixed Height (Straight or 1" Inward Offset)	<input type="checkbox"/> Fixed	<input type="checkbox"/> 9 ¹	<input type="checkbox"/> 10 ¹	<input type="checkbox"/> 11 ¹	<input type="checkbox"/> 12	<input type="checkbox"/> 13
		<input type="checkbox"/> 14	<input type="checkbox"/> 15	<input type="checkbox"/> 16	<input type="checkbox"/> 17	<input type="checkbox"/> 18
Fixed Height 2" Offset	<input type="checkbox"/> Fixed Offset 2"	<input type="checkbox"/> 11 ¹	<input type="checkbox"/> 12 ¹	<input type="checkbox"/> 13 ¹	<input type="checkbox"/> 14 ¹	
		<input type="checkbox"/> 15	<input type="checkbox"/> 16	<input type="checkbox"/> 17		

¹ Not available with Fold Down Push Handle

Height Adjustable Back Post Option

Selection required with Height Adjustable Back Post. Each Offset Back Post is offset 1/2" creating 1" overall offset.

Standard _____ N/C

1" Inward Offset _____ N/C
Not available with 8" Seat Width

Fixed Height Back Post Option

Selection required with Fixed Height Back Post. Uses 1" tubes compatible with after-market backs. Not available with back upholstery. 2" Inward or Outward Offset Back Posts are offset 1" creating 2" overall offset. 1" Offset Back Post is offset 1/2" creating 1" overall offset.

Straight _____ N/C

1" Inward Offset _____ N/C
Not available with 8" Seat Widths

2" Inward Offset _____ \$95
Not available with 8" or 9" Seat Widths

2" Outward Offset _____ \$95
Not available with Swing Away Armrests or 8" and 9" Seat Widths

Push Handle

Omit Push Handle _____ N/C

Standard Push Handle _____ N/C

Fold Down Push Handle _____ \$110
Not available on Height Adjustable Back Posts with 9"-11" Back Heights or Offset Back Posts with 11" - 14" Back Heights

Bolt-On Push Handle _____ \$185
Handles will be 1" - 2" above back height

Stroller Handle

Stroller Handles are removable and height adjustable

Ergo Stroller Handle _____ \$185

Single Straight Stroller Handle _____ \$170

Double Straight Stroller Handle _____ \$330
Not available with Push Handles or 8"-10" Seat Widths

Back Upholstery

Available only with Height Adjustable Back Posts. If selecting both upholstery and Axiom Back, the upholstery will be installed.

Padded Velcro Adjustable
Not available with 17" or 18" Back Heights in combination with 8" or 9" Seat Widths

Standard Cover _____ Std

Air Mesh Breathable Cover _____ \$75
Not available on chair widths less than 12"

Tension Adjustable E2611
Not available with 8" and 9" Seat Widths

Standard Cover _____ \$320
Not available on 10" Seat Width and 1" Inward Offset

Air Mesh Breathable Cover _____ \$390
Not available on chair widths less than 12"

Omit Back Upholstery _____ N/C

Seat Back Angle

Standard is 90° (perpendicular to floor). Minus (-°) closes back angle. Plus (+°) opens back angle.

Std +2° +4° +6° +8° +10°
 -2° -4° -6° -8° -10°

Spinergy Backrest Release Cable Color

Select desired color. If no selection is made, black will be used.

Black Red White Blue
 Yellow Orange Green Pink

Front Seat Height and Casters

Caster Size / Front Seat Height

Select Front Seat Height in row corresponding to desired Caster Size.

Front Seat Height is measured from floor to top of seat tube at front of seat sling. Adjustable up or down by 1/2". See reference E on page 7.

Caster Size	Front Seat Height														
3"	<input type="checkbox"/> 11.5	<input type="checkbox"/> 12	<input type="checkbox"/> 12.5	<input type="checkbox"/> 13	<input type="checkbox"/> 13.5	<input type="checkbox"/> 14	<input type="checkbox"/> 14.5	<input type="checkbox"/> 15	<input type="checkbox"/> 15.5	<input type="checkbox"/> 16	<input type="checkbox"/> 16.5	<input type="checkbox"/> 17	<input type="checkbox"/> 17.5	<input type="checkbox"/> 18	<input type="checkbox"/> 18.5
4"	<input type="checkbox"/> 12.5	<input type="checkbox"/> 13	<input type="checkbox"/> 13.5	<input type="checkbox"/> 14	<input type="checkbox"/> 14.5	<input type="checkbox"/> 15	<input type="checkbox"/> 15.5	<input type="checkbox"/> 16	<input type="checkbox"/> 16.5	<input type="checkbox"/> 17	<input type="checkbox"/> 17.5	<input type="checkbox"/> 18	<input type="checkbox"/> 18.5	<input type="checkbox"/> 19	<input type="checkbox"/> 19.5
5"	<input type="checkbox"/> 13.5	<input type="checkbox"/> 14	<input type="checkbox"/> 14.5	<input type="checkbox"/> 15	<input type="checkbox"/> 15.5	<input type="checkbox"/> 16	<input type="checkbox"/> 16.5	<input type="checkbox"/> 17	<input type="checkbox"/> 17.5	<input type="checkbox"/> 18	<input type="checkbox"/> 18.5	<input type="checkbox"/> 19	<input type="checkbox"/> 19.5	<input type="checkbox"/> 20	<input type="checkbox"/> 20.5
6"	<input type="checkbox"/> 14.5	<input type="checkbox"/> 15	<input type="checkbox"/> 15.5	<input type="checkbox"/> 16	<input type="checkbox"/> 16.5	<input type="checkbox"/> 17	<input type="checkbox"/> 17.5	<input type="checkbox"/> 18	<input type="checkbox"/> 18.5	<input type="checkbox"/> 19	<input type="checkbox"/> 19.5	<input type="checkbox"/> 20	<input type="checkbox"/> 20.5		

Caster

Caster Type	Caster Size	Caster Size			
		Price	3"	4"	5"
.75" Rollerblade	N/C	<input type="checkbox"/>	<input type="checkbox"/>		
.75" Lighted Rollerblade	\$80	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	
1" Poly	Std		<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>
1" Poly Aluminum	\$125		<input type="checkbox"/>	<input type="checkbox"/>	
1" Pneumatic	\$85				<input type="checkbox"/>
1.5" Poly	\$50		<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>
1.5" Soft Roll Aluminum	\$190	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>

Fork Options

- Standard Fork _____ Std
- Single Sided Fork _____ \$170
- Frog Legs Suspension Fork E1015 _____ \$450
Not included in Colored Anodizing Package
 - User Weight Under 80 lbs
 - User Weight 80 - 120 lbs
 - User Weight 120 - 160 lbs

Footrest

Footplate

- Aluminum Angle Adjustable K0040 _____ \$95
Not available with Front Seat Heights less than 14.5"
- Aluminum Angle Adjustable with Riser K0040, K0037 _____ \$185
Not available with Front Seat Heights less than 14.5"
- High Mount Aluminum Angle Adjustable K0040, K0037 _____ \$160
- Angle Adjustable Flip Under K0040, K0037 _____ \$185
Not available with footrest widths less than 9" or with Front Seat Heights less than 15". See reference F on page 7.
- High Mount Angle Adjustable Flip Under K0040, K0037 _____ \$270
Not available with 16" Seat Width in combination with Offset Frame

Seat to Footrest Length

Measured from front edge of seat sling to top, rear of footrest. Footrest length at least 2.5" shorter than Front Seat Height recommended. See reference G on page 7.

<input type="checkbox"/> 5.5"	<input type="checkbox"/> 6"	<input type="checkbox"/> 6.5"	<input type="checkbox"/> 7"
<input type="checkbox"/> 7.5"	<input type="checkbox"/> 8"	<input type="checkbox"/> 8.5"	<input type="checkbox"/> 9"
<input type="checkbox"/> 9.5"	<input type="checkbox"/> 10"	<input type="checkbox"/> 10.5"	<input type="checkbox"/> 11"
<input type="checkbox"/> 11.5"	<input type="checkbox"/> 12"	<input type="checkbox"/> 12.5"	<input type="checkbox"/> 13"
<input type="checkbox"/> 13.5"	<input type="checkbox"/> 14"	<input type="checkbox"/> 14.5"	<input type="checkbox"/> 15"
<input type="checkbox"/> 15.5"	<input type="checkbox"/> 16"	<input type="checkbox"/> 16.5"	<input type="checkbox"/> 17"
<input type="checkbox"/> 17.5"	<input type="checkbox"/> 18"		

Rear Seat Height and Wheels

Wheel Size / Rear Seat Height

Select Rear Seat Height in row corresponding to desired Wheel Size. For proper setup, selection is required when omitting rear wheels. Measured from floor to top of seat tube at front of back post. Rear Seat Height is custom to needs. See reference H on page 7.

Wheel Size	Rear Seat Height															
18"	<input type="checkbox"/> 11.5 ¹	<input type="checkbox"/> 12	<input type="checkbox"/> 12.5	<input type="checkbox"/> 13	<input type="checkbox"/> 13.5	<input type="checkbox"/> 14	<input type="checkbox"/> 14.5	<input type="checkbox"/> 15	<input type="checkbox"/> 15.5	<input type="checkbox"/> 16	<input type="checkbox"/> 16.5	<input type="checkbox"/> 17	<input type="checkbox"/> 17.5	<input type="checkbox"/> 18	<input type="checkbox"/> 18.5	<input type="checkbox"/> 19
20"	<input type="checkbox"/> 12.5 ¹	<input type="checkbox"/> 13	<input type="checkbox"/> 13.5	<input type="checkbox"/> 14	<input type="checkbox"/> 14.5	<input type="checkbox"/> 15	<input type="checkbox"/> 15.5	<input type="checkbox"/> 16	<input type="checkbox"/> 16.5	<input type="checkbox"/> 17	<input type="checkbox"/> 17.5	<input type="checkbox"/> 18	<input type="checkbox"/> 18.5	<input type="checkbox"/> 19	<input type="checkbox"/> 19.5	<input type="checkbox"/> 20
22"	<input type="checkbox"/> 13.5 ¹	<input type="checkbox"/> 14	<input type="checkbox"/> 14.5	<input type="checkbox"/> 15	<input type="checkbox"/> 15.5	<input type="checkbox"/> 16	<input type="checkbox"/> 16.5	<input type="checkbox"/> 17	<input type="checkbox"/> 17.5	<input type="checkbox"/> 18	<input type="checkbox"/> 18.5	<input type="checkbox"/> 19	<input type="checkbox"/> 19.5	<input type="checkbox"/> 20		
24"	<input type="checkbox"/> 14.5 ¹	<input type="checkbox"/> 15	<input type="checkbox"/> 15.5	<input type="checkbox"/> 16	<input type="checkbox"/> 16.5	<input type="checkbox"/> 17	<input type="checkbox"/> 17.5	<input type="checkbox"/> 18	<input type="checkbox"/> 18.5	<input type="checkbox"/> 19	<input type="checkbox"/> 19.5	<input type="checkbox"/> 20				

¹ Available only with Low Profile tires

One Arm Drive

- One Arm Drive E0958 _____ \$1495
Available only with 20", 22" and 24" spoke wheels and Quick Release axles. Spoke wheel style is exclusive to OAD. Standard handrim will be Aluminum anodized. Rear Wheel Spacing limitations exist based on camber and armrest selections.
 - Rear Wheel Size: 20" 22" 24"
 - Drive Handrim Side: Left Right
 - Drive Handrim Type: Aluminum _____ N/C
 Poly Coated _____ \$200

Rear Seat Height and Wheels

Rear Wheel

No selection needed with One Arm Drive

		Rear Wheel Size				
		Price	18"	20"	22"	24"
Rear Wheel Type	MAXX Performance Spoke	Std	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>
	Mag (Black)	\$55		<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>
	Colored Mag	\$100		<input type="checkbox"/>	<input type="checkbox"/>	
	Spinergy Spox	\$845		<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>
	Spinergy LX	\$1050			<input type="checkbox"/>	<input type="checkbox"/>

Spinergy Options

Select spoke color and handrim tab for Spinergy wheels

Spoke Color:

- Black Red White Blue Yellow
 Pink ¹ Green ¹ Purple ¹ Orange ¹

¹Available only with Spinergy LX

Handrim Tab:

Available only with Aluminum Anodized and Plastic Coated Handrims

- Short Tab _____ N/C
 Long Tab _____ N/C

Colored Mag

Mag Color: Florescent Red White Blue Pink

Rear Tire Type

Full Profile

- Street Pneumatic 1 3/8" E2211, E2212 _____ Std
 Full Poly 1 3/8" E2218 _____ \$75
 Pneumatic with Airless Insert 1 3/8" E2211, E2213 _____ \$95

Low Profile

- Low Poly 1" E2218 _____ \$75
 Sentinel Solid Tire 1" E2218 _____ \$95
Not available on 18" wheel size or 20" and 22" Mag wheels
 High Pressure 1" E2211, E2212 _____ \$97
Not available on Mag wheel or 18" wheel size
 Schwalbe Marathon Plus Evolution 1" E2211, E2212 _____ \$135
Available only on 24" MAXX Performance or Spinergy wheels

Handrim

No selection needed with One Arm Drive

Standard Handrims

- Aluminum Anodized _____ Std
 Plastic Coated _____ \$112
 Projection (8 Vertical Projections) _____ \$325
Available only with full profile tires in combination with 22"-24" MAXX Performance or 24" Mag wheel

Omit Handrim

Natural Fit Handrims

Not available on 18" or 20" wheel sizes or 22" Mag wheel. Not recommended with Mag wheel.

- Natural Fit LT Regular Grip _____ \$465
 Natural Fit Standard Regular Grip _____ \$445
 Natural Fit LT Super Grip _____ \$515
 Natural Fit Standard Super Grip _____ \$495
 Natural Fit LT Omit Grip _____ \$395
 Natural Fit Standard Omit Grip _____ \$345

Handrim Option

- 3M Micro Grip _____ \$65
Not available on Plastic Coated or Projection Handrims

Axle

No selection needed with Drum Brake

- Quick Release _____ Std
 Quad Quick Release _____ \$185
Not available with Mag or Spinergy wheels

Rear Wheel Spacing

Adjustable up to .5" outward. Minimum spacing with offset frame is 1.25" with 0° camber, 1" with 2° or 4° camber and .75" with 6° or 8° camber. .75" minimum with 0° camber on standard frame. 0° Camber and Mag wheel combination requires 1" minimum Rear Wheel Spacing. No selection needed with Wheelchair Drive Ready option. .5" and .75" not available with One Arm Drive. See Reference J on page 7.

- .5" .75" 1" 1.25" 1.5" 1.75"

Wheelchair Drive Ready

Chair will ship from Ki Mobility with adaption kit installed for selected Alber system. Order system separately from Alber USA. Select Rear Seat Height (full profile selections only) in drive wheel size on Rear Seat Height chart above. Requires Rear Anti-Tips. Not available with 5th Wheel or Spinergy wheels.

- e-motion Ready _____ N/C
 twion Ready _____ N/C

Camber Selection: 0° 3°

Center of Gravity Preset

Measured from front of back post to center of rear axle. 0" indicates axle will be directly under back post. Center of Gravity setup may be limited with certain configurations. 2" minimum Center of Gravity with 5th Wheel and One Arm Drive and 0" minimum with Transit Option. See reference I on page 7.

Amputee Positions

<input type="checkbox"/> -1.5"	<input type="checkbox"/> -1.25"	<input type="checkbox"/> -1"	<input type="checkbox"/> -.75"
<input type="checkbox"/> -.5"	<input type="checkbox"/> -.25"	<input type="checkbox"/> 0"	<input type="checkbox"/> .25"

Standard Positions

<input type="checkbox"/> .5"	<input type="checkbox"/> .75"	<input type="checkbox"/> 1"	<input type="checkbox"/> 1.25"
<input type="checkbox"/> 1.5"	<input type="checkbox"/> 1.75"	<input type="checkbox"/> 2"	<input type="checkbox"/> 2.25"
<input type="checkbox"/> 2.5"	<input type="checkbox"/> 2.75"	<input type="checkbox"/> 3"	<input type="checkbox"/> 3.25"
<input type="checkbox"/> 3.5"	<input type="checkbox"/> 3.75"	<input type="checkbox"/> 4"	<input type="checkbox"/> 4.25"

Rear Wheel Camber

The following apply if One Arm Drive is selected: 6° and 8° not available. 4° not available with T-Arms.

- 0° 2° 4° 6° 8°

Wheel Locks

Wheel Lock	
<input type="checkbox"/> Push to Lock _____ Std	Dual - Drum Brake E2228 _____ \$660
<input type="checkbox"/> Push to Lock Flush Mount _____ N/C	<i>Available only with 20", 22" and 24" MAXX Performance Spoke wheels. Not available with Fold Down Push Handles or Ergo and Single Stroller Handles without integrated push handles. Axle is Quick Release.</i>
<i>Not available with Seat Depths less than 11" or Offset Frame</i>	<input type="checkbox"/> Push to Lock <input type="checkbox"/> Pull to Lock
<input type="checkbox"/> Pull to Lock _____ N/C	<input type="checkbox"/> Omit Wheel Locks _____ N/C
<input type="checkbox"/> Low Profile Scissor Lock _____ \$65	
<i>Not available with Low Poly, Full Poly, Sentinel tires or Offset Hanger</i>	
<input type="checkbox"/> Short Thro Scissor Lock _____ \$125	
<input type="checkbox"/> Drum Brake E2228 _____ \$560	
<i>Available only with 20", 22" and 24" MAXX Performance Spoke wheels. Not available with Fold Down Push Handles or Ergo and Single Stroller Handles without integrated push handles. Axle is Quick Release.</i>	

Armrests

Armrest	
Height Adjustable T-Arm (8.5" - 13.5") E0973	<input type="checkbox"/> Tubular Flip Up _____ Std
<i>Not available with Seat Depths less than 13". 16" minimum Seat Depth with Offset Frame. 1.25" minimum Rear Wheel Spacing.</i>	<i>Requires Fixed Height Backrest. Not available with Arm Pads, Back Heights less than 13" or Back Heights less than 14" with Offset.</i>
<input type="checkbox"/> Desk _____ \$280	<input type="checkbox"/> Short Swing Away (7" - 9") E0973 _____ \$230
<input type="checkbox"/> Full _____ \$280	<i>Not available with Arm Pads or Outward Offset Back Posts. If selection is not available with rear wheel and rear seat height configuration, then Tall Swing Away armrest will be used.</i>
Height Adjustable T-Arm Low (7.5" - 11.5") E0973	<input type="checkbox"/> Tall Swing Away (9.25" - 11.25") E0973 _____ \$230
<i>Not available with Seat Depths less than 12". 13" minimum Seat Depth with Offset Frame. 1.25" minimum Rear Wheel Spacing.</i>	<i>Not available with Arm Pads or Outward Offset Back Posts</i>
<input type="checkbox"/> Desk _____ \$280	<input type="checkbox"/> Omit Armrests _____ N/C
<input type="checkbox"/> Full _____ \$280	
Angle Adjustable Locking Extendable Flip Up E0973	
<i>Requires Fixed Height Backrest. Not available with backrest heights less than 12".</i>	
<input type="checkbox"/> Desk _____ \$265	
<input type="checkbox"/> Full _____ \$265	
<input type="checkbox"/> Height Adjustable Pediatric Desk (6.5" - 10") E0973 _____ \$230	
<i>Seat Depth plus Sling Position must be at least 10". Seat Depth plus Sling Position must be at least 12" with Offset Frame 1.25" minimum Rear Wheel Spacing.</i>	

Side Guards

Not available with T-Arms or Height Adjustable Pediatric Desk Armrests

Direct Mount		Clamp-On Mount	
<input type="checkbox"/> Composite Side Guards - Adult (14" - 16") K0108 _____ \$143	<i>Seat Depth plus Sling Position must be at least 12"</i>	Aluminum Adjustable Fender with Foam Pad K0108	
<input type="checkbox"/> Composite Side Guards - Pediatric (10" - 14") K0108 _____ \$143	<i>Seat Depth plus Sling Position must be at least 10"</i>	<i>10" minimum Seat Depth</i>	
		<input type="checkbox"/> Black Matte _____ \$305	
		<input type="checkbox"/> Match Frame Color _____ \$355	
		<i>Not available with Hydrographics</i>	
		<input type="checkbox"/> Carbon Fiber Adjustable Fender K0108 _____ \$515	
		<i>10" minimum Seat Depth</i>	

Accessories

5th Wheel		Spoke Guards	
Rotating removable rear caster. Not available with Anti-Tips.		<input type="checkbox"/> Clear Spoke Guards K0065 _____ \$115	
<input type="checkbox"/> Standard 5th Wheel E0971 _____ \$250		<i>Not available with Mag wheel</i>	
Dynamic 5th Wheel K0108 _____ \$250		<input type="checkbox"/> Black Spoke Guards K0065 _____ \$115	
Elastomer: <input type="checkbox"/> Soft <input type="checkbox"/> Firm		<i>Not available with 18" or 20" wheel sizes or Mag wheel</i>	
<i>Both elastomers are included; the selected elastomer will be installed.</i>		Designer Spoke Guards K0065 _____ \$250	
		<i>Available only with 20" and 22" wheel sizes. Not available with Mag wheel.</i>	
		Design:	
		<input type="checkbox"/> Blue Monster <input type="checkbox"/> Fairy Garden <input type="checkbox"/> Havin' a Hoot	
		<input type="checkbox"/> Robo Fun <input type="checkbox"/> Rockin' Reptiles <input type="checkbox"/> Super Puppies	
		<input type="checkbox"/> Tutu Cute <input type="checkbox"/> Tuxedo Time	
<input type="checkbox"/> Rear Anti-Tips E0971 _____ \$150			
<input type="checkbox"/> Omit Rear Anti-Tips _____ N/C			

Accessories

Calf Strap

Not available with 8" or 9" wide chairs

- Velcro Adjustable K0038 \$40
- Padded Velcro Adjustable K0038 \$68
- Bodypoint Padded AeroMesh K0038 \$118

Pelvic Positioning Belt

- 1" Auto Buckle E0978 \$80
- 1" Padded Auto Buckle E0978 \$170
- 1 1/2" Auto Buckle E0978 \$80
- 1 1/2" Padded Auto Buckle E0978 \$170
- 2" Airline Buckle E0978 \$110


Bodypoint Evoflex - Pelvic Stabilizer E0978 \$188
 Size: Small (16-24") Medium (19-29") Large (25-37")

Belt Mounting Kit

- Ki Mobility Belt Mounting Clamp \$50

Other Accessories

- Tool Kit \$45
- Removable Underseat Pouch \$50
- Pediatric Ki Mobility Backpack \$65
Best fit on chairs more than 8" wide
- Ki Mobility Backpack \$65
Best fit on chairs more than 10" wide
- Luggage Carrier \$85
- Cane and Crutch Holder \$195
- Canopy with Hardware K0108 \$280
See image below



Color

Frame and Accent

Select desired frame and accent colors. Accent color is used on backrest rigidizer bar and front caster wings; if no accent color is selected, frame color will be used.

Gloss Finish

	Frame	Accent
Black	<input type="checkbox"/>	<input type="checkbox"/>
Blue Skies	<input type="checkbox"/>	<input type="checkbox"/>
Bubblegum Pink	<input type="checkbox"/>	<input type="checkbox"/>
Candy Blue	<input type="checkbox"/>	<input type="checkbox"/>
Candy Purple	<input type="checkbox"/>	<input type="checkbox"/>
Candy Red	<input type="checkbox"/>	<input type="checkbox"/>
Copper	<input type="checkbox"/>	<input type="checkbox"/>
Granny Smith	<input type="checkbox"/>	<input type="checkbox"/>
Metallic Orange	<input type="checkbox"/>	<input type="checkbox"/>
Silver	<input type="checkbox"/>	<input type="checkbox"/>
Teal	<input type="checkbox"/>	<input type="checkbox"/>
Yellow	<input type="checkbox"/>	<input type="checkbox"/>

Matte Finish

	Frame	Accent
Charcoal Gray	<input type="checkbox"/>	<input type="checkbox"/>
Emerald Green	<input type="checkbox"/>	<input type="checkbox"/>
Iridescent Burgundy	<input type="checkbox"/>	<input type="checkbox"/>
Iridescent Navy	<input type="checkbox"/>	<input type="checkbox"/>
Military Green	<input type="checkbox"/>	<input type="checkbox"/>
Retro Red	<input type="checkbox"/>	<input type="checkbox"/>
Royal Blue	<input type="checkbox"/>	<input type="checkbox"/>
Shadow Black	<input type="checkbox"/>	<input type="checkbox"/>
Smooth White	<input type="checkbox"/>	<input type="checkbox"/>

Hydrographics

- Hydrographics \$695
- AOR1 (Desert Digital Camo)
- Muddy Girl
- Carbon Fiber¹
- Pink Camo
- NEXT Camo
- Zebra¹

¹Select Hydrographics Base Color. Zebra not available in Black.

- Hydrographics Base Color N/C
- Black Blue Granny Smith Pink Purple Red

Colored Anodizing Package

Select desired color. If no selection is made, black will be used. Included in Colored Anodizing Package: Backrest plates, camber tube, towers, CG clamps, side guard receivers, front caster forks (standard or unitine) and Dynamic 5th Wheel split clamp. Camber tube not included if Wheelchair Drive Ready is selected.

- Black N/C
- Blue \$165
- Bright Green \$165
- Orange \$165
- Pewter \$165
- Pink \$165
- Purple \$165
- Red \$165
- Gold \$165

References

A. Seat Width See Diagram 1

Measured from outside of frame tube on one side to outside of frame tube on other side.

B. Seat Depth See Diagram 3

Measured from front of back posts to front edge of seat sling. Seat sling starts at beginning of bend at front of frame.

C. Front Frame Bend See Diagram 3

Angle between front seat tube and ground.

D. Sling Position See Diagram 2

With 0 Sling Position, sling is approximately 4" from the front tube. The 1" and 2" Sling Positions extend the frame in front of the sling by that amount. -1 pushes the seat 1" farther forward. See Charts 1 and 2 for more detail.

E. Front Seat Height See Diagram 1

Measured from floor to top of seat tube at front of seat sling.

F. Footrest Width See Diagram 1

Measured from inside of front tube to inside of front tube on other side. Footrest width is listed below for each seat width. See Chart 1 for more detail.

G. Seat to Footrest Length See Diagram 3

Measured from front edge of seat sling to top of rear of footrest. Footrest length of at least 2.5" shorter than front seat height recommended.

H. Rear Seat Height See Diagram 3

Measured from floor to top of seat tube at front of back post. Rear seat height is custom to needs.

I. Center of Gravity Preset See Diagram 3

Measured from front of back post to center of rear axle. 0" indicates axle will be directly under back post.

J. Rear Wheel Spacing See Diagram 1

Measured from outside of seat back to inside of rear tire. Adjustable up to .5" outward.

Chart 1

F. Footrest Width		
Seat Width	Straight Frame (Diagram 1)	Offset Frame
8"	5.75"	9"
9"	6.75"	10"
10"	7.75"	11"
11"	8.75"	12"
12"	9.75"	13"
13"	10.75"	14"
14"	11.75"	15"
15"	12.75"	16"
16"	13.75"	17"

Diagram 1




Diagram 2




Diagram 3




Chart 2

Sling Position / Adjustable Depth / Seat Depth Compatibility									
Sling Position	-1", 0"			1"			2"		
	1"	2"	3"	1"	2"	3"	1"	2"	3"
Seat Depth	14"								X
	15"					X		X	X
	16"			X		X	X	X	X

"X" indicates Sling Position, Adjustable Depth and Seat Depth combination is not available.
 -1" Sling Position, 0" Adjustable Depth and 8" - 13" Seat Depths have no restrictions.